

---

# HVL Common Code Base Documentation

*Release 0.10.1*

**Mikolaj Rybiński, David Graber, Henrik Menne, Alise Chachereau,**

**Jan 24, 2022**



## CONTENTS:

<b>1</b>	<b>HVL Common Code Base</b>	<b>1</b>
1.1	Features . . . . .	1
1.2	Documentation . . . . .	5
1.3	Credits . . . . .	5
<b>2</b>	<b>Installation</b>	<b>7</b>
2.1	Stable release . . . . .	7
2.2	From sources . . . . .	7
2.3	Additional system libraries . . . . .	8
<b>3</b>	<b>Usage</b>	<b>9</b>
<b>4</b>	<b>API Documentation</b>	<b>11</b>
4.1	hvl_ccb . . . . .	11
<b>5</b>	<b>Contributing</b>	<b>183</b>
5.1	Types of Contributions . . . . .	183
5.2	Get Started! . . . . .	184
5.3	Merge Request Guidelines . . . . .	185
5.4	Tips . . . . .	185
5.5	Deploying . . . . .	186
<b>6</b>	<b>Credits</b>	<b>187</b>
6.1	Maintainers . . . . .	187
6.2	Authors . . . . .	187
6.3	Contributors . . . . .	187
<b>7</b>	<b>History</b>	<b>189</b>
7.1	0.10.1 (2022-01-24) . . . . .	189
7.2	0.10.0 (2022-01-17) . . . . .	189
7.3	0.9.0 (2022-01-07) . . . . .	190
7.4	0.8.5 (2021-11-05) . . . . .	190
7.5	0.8.4 (2021-10-22) . . . . .	191
7.6	0.8.3 (2021-09-27) . . . . .	191
7.7	0.8.2 (2021-08-27) . . . . .	191
7.8	0.8.1 (2021-08-13) . . . . .	191
7.9	0.8.0 (2021-07-02) . . . . .	191
7.10	0.7.1 (2021-06-04) . . . . .	192
7.11	0.7.0 (2021-05-25) . . . . .	192
7.12	0.6.1 (2021-05-08) . . . . .	192
7.13	0.6.0 (2021-04-23) . . . . .	193

7.14	0.5.0 (2020-11-11) . . . . .	193
7.15	0.4.0 (2020-07-16) . . . . .	194
7.16	0.3.5 (2020-02-18) . . . . .	194
7.17	0.3.4 (2019-12-20) . . . . .	195
7.18	0.3.3 (2019-05-08) . . . . .	195
7.19	0.3.2 (2019-05-08) . . . . .	195
7.20	0.3.1 (2019-05-02) . . . . .	195
7.21	0.3 (2019-05-02) . . . . .	195
7.22	0.2.1 (2019-04-01) . . . . .	196
7.23	0.2.0 (2019-03-31) . . . . .	196
7.24	0.1.0 (2019-02-06) . . . . .	196
<b>8</b>	<b>Indices and tables</b>	<b>197</b>
	<b>Python Module Index</b>	<b>199</b>
	<b>Index</b>	<b>201</b>

---

**CHAPTER  
ONE**

---

## **HVL COMMON CODE BASE**

Python common code base (CCB) to control devices, which are used in high-voltage research. All implemented devices are used and tested in the High Voltage Laboratory ([HVL](#)) of the Federal Institute of Technology Zurich (ETH Zurich).

- Free software: GNU General Public License v3
- Copyright (c) 2019-2022 ETH Zurich, SIS ID and HVL D-ITET

### **1.1 Features**

For managing multi-device experiments instantiate the `ExperimentManager` utility class.

#### **1.1.1 Devices**

The device wrappers in `hvl_ccb` provide a standardised API with configuration dataclasses, various settings and options, as well as start/stop methods. Currently wrappers are available to control the following devices:

Function/Type	Devices
Data acquisition	LabJack (T4, T7, T7-PRO; requires <a href="#">LJM Library</a> ) Pico Technology PT-104 Platinum Resistance Data Logger (requires <a href="#">PicoSDK/libusbpt104</a> )
Digital Delay Generator	Highland T560
Digital IO	LabJack (T4, T7, T7-PRO; requires <a href="#">LJM Library</a> )
Experiment control	HVL Cube with and without Power Inverter
Gas Analyser	MBW 973-SF6 gas dew point mirror analyzer Pfeiffer Vacuum TPG (25x, 26x and 36x) controller for compact pressure gauges SST Luminox oxygen sensor
I2C host	TiePie (HS5, WS5; requires <a href="#">LibTiePie SDK</a> )
Laser	CryLaS pulsed laser CryLaS laser attenuator
Oscilloscope	Rhode & Schwarz RTO 1024 TiePie (HS5, HS6, WS5; requires <a href="#">LibTiePie SDK</a> )
Power supply	Elektro-Automatik PSI9000 FuG Elektronik Heinzinger PNC Technix capacitor charger
Stepper motor drive	Newport SMC100PP Schneider Electric ILS2T
Temperature control	Lauda PRO RP 245 E circulation thermostat
Waveform generator	TiePie (HS5, WS5; requires <a href="#">LibTiePie SDK</a> )
2	<b>Chapter 1. HVL Common Code Base</b>

Each device uses at least one standardised communication protocol wrapper.

### 1.1.2 Communication protocols

In hvl\_ccb by “communication protocol” we mean different levels of communication standards, from the low level actual communication protocols like serial communication to application level interfaces like VISA TCP standard. There are also devices in hvl\_ccb that use a dummy communication protocol; this is because these devices are build on proprietary manufacturer libraries that communicate with the corresponding devices, as in the case of TiePie or LabJack devices.

The communication protocol wrappers in hvl\_ccb provide a standardised API with configuration dataclasses, as well as open/close and read/write/query methods. Currently, wrappers for the following communication protocols are available:

Communication protocol	Devices using
Modbus TCP	Schneider Electric ILS2T stepper motor drive
OPC UA	HVL Cube with and without Power Inverter
Serial	CryLaS pulsed laser and laser attenuator FuG Elektronik power supply (e.g. capacitor charger HCK) using the Probus V protocol Heinzinger PNC power supply using Heinzinger Digital Interface I/II SST Luminox oxygen sensor MBW 973-SF6 gas dew point mirror analyzer Newport SMC100PP single axis driver for 2-phase stepper motors Pfeiffer Vacuum TPG (25x, 26x and 36x) controller for compact pressure gauges Technix capacitor charger
TCP	Lauda PRO RP 245 E circulation thermostat
Telnet	Technix capacitor charger
VISA TCP	Elektro-Automatik PSI9000 DC power supply Rhode & Schwarz RTO 1024 oscilloscope
<i>proprietary</i>	LabJack (T4, T7, T7-PRO) devices, which communicate via <a href="#">LJM Library</a> Pico Technology PT-104 Platinum Resistance Data Logger, which communicate via <a href="#">PicoSDK/libusbpt104</a> TiePie (HS5, HS6, WS5) oscilloscopes, generators and I2C hosts, which communicate via <a href="#">LibTiePie SDK</a>

### 1.1.3 Sensor and Unit Conversion Utility

The Sensor and Unit Conversion Utility is a submodule that allows on the one hand a unified implementation of hardware-sensors and on the other hand provides a unified way to convert units. The utility can be used with single numbers (`int`, `float`) as well as array-like structures containing single numbers (`np.array()`, `list`, `dict`, `tuple`).

Currently the following sensors are implemented:

- LEM LT 4000S
- LMT 70A

The following unit conversion classes are implemented:

- Temperature (Kelvin, Celsius, Fahrenheit)
- Pressure (Pascal, Bar, Atmosphere, Psi, Torr, Millimeter Mercury)

## 1.2 Documentation

Note: if you're planning to contribute to the `hvl_ccb` project read the **Contributing** section in the HVL CCB documentation.

Do either:

- read [HVL CCB documentation at RTD](#),

or

- build and read HVL CCB documentation locally; install first [Graphviz](#) (make sure to have the `dot` command in the executable search path) and the Python build requirements for documentation:

```
$ pip install docs/requirements.txt
```

and then either on Windows in Git BASH run:

```
$ ./make.sh docs
```

or from any other shell with GNU Make installed run:

```
$ make docs
```

The target index HTML ("`docs/_build/html/index.html`") should open automatically in your Web browser.

## 1.3 Credits

This package was created with [Cookiecutter](#) and the [audreyr/cookiecutter-pypackage](#) project template.



## INSTALLATION

### 2.1 Stable release

To install HVL Common Code Base, run this command in your terminal:

```
$ pip install hvl_ccb
```

To install HVL Common Code Base with optional Python libraries that require manual installations of additional system libraries, you need to specify on installation extra requirements corresponding to these controllers. For instance, to install Python requirements for LabJack and TiePie devices, run:

```
$ pip install "hvl_ccb[tiepie,labjack]"
```

See below for the info about additional system libraries and the corresponding extra requirements.

To install all extra requirements run:

```
$ pip install "hvl_ccb[all]"
```

This is the preferred method to install HVL Common Code Base, as it will always install the most recent stable release.

If you don't have `pip` installed, this [Python installation guide](#) can guide you through the process.

### 2.2 From sources

The sources for HVL Common Code Base can be downloaded from the [GitLab](#) repo.

You can either clone the repository:

```
$ git clone git@gitlab.com:ethz_hvl/hvl_ccb.git
```

Or download the [tarball](#):

```
$ curl -OL https://gitlab.com/ethz_hvl/hvl_ccb/-/archive/master/hvl_ccb.tar.gz
```

Once you have a copy of the source, you can install it with:

```
$ pip install .
```

## 2.3 Additional system libraries

If you have installed *hvl\_ccb* with any of the extra features corresponding to device controllers, you must additionally install respective system library; these are:

Extra feature	Additional system library
labjack	<a href="#">LJM Library</a>
picotech	<a href="#">PicoSDK</a> (Windows) / <a href="#">libusbpt104</a> (Ubuntu/Debian)
tiepie	<a href="#">LibTiePie SDK</a>

For more details on installation of the libraries see docstrings of the corresponding `hvl_ccb` modules.

---

**CHAPTER  
THREE**

---

**USAGE**

To use HVL Common Code Base in a project:

```
import hvl_ccb
```



## API DOCUMENTATION


### 4.1 hvl\_ccb

#### 4.1.1 Subpackages

[hvl\\_ccb.comm](#)

**Submodules**

[hvl\\_ccb.comm.base](#)



Module with base classes for communication protocols.

**class** `AsyncCommunicationProtocol(config)`

Bases: `hvl_ccb.comm.base.CommunicationProtocol`

Abstract base class for asynchronous communication protocols

**static config\_cls()** → Type[`hvl_ccb.comm.base.AsyncCommunicationProtocolConfig`]

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**read()** → str

Read a single line of text as `str` from the communication.

**Returns** text as `str` including the terminator, which can also be empty “”

**read\_all**(*n\_attempts\_max*: *Optional[int]* = *None*, *attempt\_interval\_sec*: *Optional[Union[int, float]]* = *None*) → *Optional[str]*

Read all lines of text from the connection till nothing is left to read.

#### Parameters

- **n\_attempts\_max** – Amount of attempts how often a non-empty text is tried to be read
- **attempt\_interval\_sec** – time between the reading attempts

**Returns** A multi-line *str* including the terminator internally

**abstract** **read\_bytes**() → *bytes*

Read a single line as *bytes* from the communication.

This method uses *self.access\_lock* to ensure thread-safety.

**Returns** a single line as *bytes* containing the terminator, which can also be empty b””

**read\_nonempty**(*n\_attempts\_max*: *Optional[int]* = *None*, *attempt\_interval\_sec*: *Optional[Union[int, float]]* = *None*) → *Optional[str]*

Try to read a non-empty single line of text as *str* from the communication. If the host does not reply or reply with white space only, it will return *None*.

**Returns** a non-empty text as a *str* or *None* in case of an empty string

#### Parameters

- **n\_attempts\_max** – Amount of attempts how often a non-empty text is tried to be read
- **attempt\_interval\_sec** – time between the reading attempts

**read\_text**() → *str*

Read one line of text from the serial port. The input buffer may hold additional data afterwards, since only one line is read.

NOTE: backward-compatibility proxy for *read* method; to be removed in v1.0

**Returns** String read from the serial port; “ if there was nothing to read.

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**read\_text\_nonempty**(*n\_attempts\_max*: *Optional[int]* = *None*, *attempt\_interval\_sec*: *Optional[Union[int, float]]* = *None*) → *Optional[str]*

Reads from the serial port, until a non-empty line is found, or the number of attempts is exceeded.

NOTE: backward-compatibility proxy for *read* method; to be removed in v1.0

Attention: in contrast to *read\_text*, the returned answer will be stripped of a whitespace newline terminator at the end, if such terminator is set in the initial configuration (default).

#### Parameters

- **n\_attempts\_max** – maximum number of read attempts
- **attempt\_interval\_sec** – time between the reading attempts

**Returns** String read from the serial port; “ if number of attempts is exceeded or serial port is not opened.

**write**(*text*: *str*)

Write text as *str* to the communication.

**Parameters** **text** – test as a *str* to be written

```
abstract write_bytes(data: bytes) → int
    Write data as bytes to the communication.

    This method uses self.access_lock to ensure thread-safety.

    Parameters data – data as bytes-string to be written

    Returns number of bytes written

write_text(text: str)
    Write text to the serial port. The text is encoded and terminated by the configured terminator.

    NOTE: backward-compatibility proxy for read method; to be removed in v1.0

    Parameters text – Text to send to the port.

    Raises SerialCommunicationIOError – when communication port is not opened

class AsyncCommunicationProtocolConfig(terminator: bytes = b'\r\n', encoding: str = 'utf-8',
                                             encoding_error_handling: str = 'strict',
                                             wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                             default_n_attempts_read_text_nonempty: int = 10)

Bases: object
Base configuration data class for asynchronous communication protocols

clean_values()

default_n_attempts_read_text_nonempty: int = 10
    default number of attempts to read a non-empty text

encoding: str = 'utf-8'
    Standard encoding of the connection. Typically this is utf-8, but can also be latin-1 or something from here: https://docs.python.org/3/library/codecs.html#standard-encodings

encoding_error_handling: str = 'strict'
    Encoding error handling scheme as defined here: https://docs.python.org/3/library/codecs.html#error-handlers By default strict error handling that raises UnicodeError.

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.

    Parameters
        • fieldname – name of the field
        • value – value to assign

is_configdataclass = True

classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.
```

**classmethod required\_keys()** → Sequence[str]  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**terminator: bytes = b'\r\n'**  
The terminator character. Typically this is b'\r\n' or b'\n', but can also be b'\r' or other combinations.  
This defines the end of a single line.

**wait\_sec\_read\_text\_nonempty: Union[int, float] = 0.5**  
time to wait between attempts of reading a non-empty text

**exception CommunicationException**

Bases: *hvl\_ccb.exception.CCBException*

**class CommunicationProtocol(config)**

Bases: *hvl\_ccb.configuration.ConfigurationMixin, abc.ABC*

Communication protocol abstract base class.

Specifies the methods to implement for communication protocol, as well as implements some default settings and checks.

**access\_lock**

Access lock to use with context manager when accessing the communication protocol (thread safety)

**abstract close()**

Close the communication protocol

**abstract open()**

Open communication protocol

**class NullCommunicationProtocol(config)**

Bases: *hvl\_ccb.comm.base.CommunicationProtocol*

Communication protocol that does nothing.

**close() → None**

Void close function.

**static config\_cls()** → Type[*hvl\_ccb.configuration.EmptyConfig*]

Empty configuration

**Returns** EmptyConfig

**open() → None**

Void open function.

**class SyncCommunicationProtocol(config)**

Bases: *hvl\_ccb.comm.base.AsyncCommunicationProtocol, abc.ABC*

Abstract base class for synchronous communication protocols with *query()*

**static config\_cls()** → Type[*hvl\_ccb.comm.base.SyncCommunicationProtocolConfig*]

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query(command: str) → Optional[str]**


Send a command to the interface and handle the status message. Eventually raises an exception.

**Parameters** **command** – Command to send

**Returns** Answer from the interface, which can be None instead of an empty reply

```
class SyncCommunicationProtocolConfig(terminator: bytes = b'\r\n', encoding: str = 'utf-8',
                                      encoding_error_handling: str = 'strict',
                                      wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                      default_n_attempts_read_text_nonempty: int = 10)
Bases: hvl_ccb.comm.base.AsyncCommunicationProtocolConfig
```

## hvl\_ccb.comm.labjack\_ljm



Communication protocol for LabJack using the LJM Library. Originally developed and tested for LabJack T7-PRO.

Makes use of the LabJack LJM Library Python wrapper. This wrapper needs an installation of the LJM Library for Windows, Mac OS X or Linux. Go to: <https://labjack.com/support/software/installers/ljm> and <https://labjack.com/support/software/examples/ljm/python>

```
class LJMCommunication(configuration)
Bases: hvl_ccb.comm.base.CommunicationProtocol
```

Communication protocol implementing the LabJack LJM Library Python wrapper.

**close()** → *None*

Close the communication port.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**property is\_open: bool**

Flag indicating if the communication port is open.

**Returns** *True* if the port is open, otherwise *False*

**open()** → *None*

Open the communication port.

```
read_name(*names: str, return_num_type: typing.Type[numbers.Real] = <class 'float'>) →
    Union[numbers.Real, Sequence[numbers.Real]]
```

Read one or more input numeric values by name.

### Parameters

- **names** – one or more names to read out from the LabJack
- **return\_num\_type** – optional numeric type specification for return values; by default *float*.

**Returns** answer of the LabJack, either single number or multiple numbers in a sequence, respectively, when one or multiple names to read were given

**Raises** **TypeError** – if read value of type not compatible with *return\_num\_type*

**write\_name**(*name*: str, *value*: numbers.Real) → *None*

Write one value to a named output.

### Parameters

- **name** – String or with name of LabJack IO
- **value** – is the value to write to the named IO port

**write\_names**(*name\_value\_dict*: Dict[str, numbers.Real]) → *None*

Write more than one value at once to named outputs.

**Parameters** **name\_value\_dict** – is a dictionary with string names of LabJack IO as keys and corresponding numeric values

```
class LJMCommunicationConfig(device_type: Union[str, hvl_ccb._dev.labjack.DeviceType] = 'ANY',
                               connection_type: Union[str,
                               hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig.ConnectionType] =
                               'ANY', identifier: str = 'ANY')
```

Bases: object

Configuration dataclass for [LJMCommunication](#).

```
class ConnectionType(value=<no_arg>, names=None, module=None, type=None, start=1,
                      boundary=None)
```

Bases: [hvl\\_ccb.utils.enum.AutoNumberNameEnum](#)

LabJack connection type.

**ANY** = 1

**ETHERNET** = 4

**TCP** = 3

**USB** = 2

**WIFI** = 5

```
class DeviceType(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
```

Bases: [hvl\\_ccb.utils.enum.AutoNumberNameEnum](#)

LabJack device types.

Can be also looked up by ambiguous Product ID (*p\_id*) or by instance name: `python LabJackDeviceType(4)` is `LabJackDeviceType('T4')`

**ANY** = 1

**T4** = 2

**T7** = 3

**T7\_PRO** = 4

```
classmethod get_by_p_id(p_id: int) → Union[hvl_ccb._dev.labjack.DeviceType,  
List[hvl_ccb._dev.labjack.DeviceType]]
```

Get LabJack device type instance via LabJack product ID.

Note: Product ID is not unambiguous for LabJack devices.

**Parameters** **p\_id** – Product ID of a LabJack device

**Returns** Instance or list of instances of *LabJackDeviceType*

**Raises** **ValueError** – when Product ID is unknown

```
clean_values() → None
```

Performs value checks on device\_type and connection\_type.

```
connection_type: Union[str,
```

```
hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig.ConnectionType] = 'ANY'
```

Can be either string or of enum *ConnectionType*.

```
device_type: Union[str, hvl_ccb._dev.labjack.DeviceType] = 'ANY'
```

Can be either string ‘ANY’, ‘T7\_PRO’, ‘T7’, ‘T4’, or of enum *DeviceType*.

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
identifier: str = 'ANY'
```

The identifier specifies information for the connection to be used. This can be an IP address, serial number, or device name. See the LabJack docs (<https://labjack.com/support/software/api/ljm/function-reference/ljmopens/identifier-parameter>) for more information.

```
is_configdataclass = True
```

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.


**Returns** a list of strings containing all required keys.

```
exception LJMCommunicationError
```

Bases: *hvl\_ccb.comm.base.CommunicationException*

Errors coming from LJMCommunication.

## `hvl_ccb.comm.modbus_tcp`



Communication protocol for modbus TCP ports. Makes use of the [pymodbus](#) library.

**class** `ModbusTcpCommunication(configuration)`  
*Bases:* `hvl_ccb.comm.base.CommunicationProtocol`

Implements the Communication Protocol for modbus TCP.

**close()**

Close the Modbus TCP connection.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**open() → None**

Open the Modbus TCP connection.

**Raises** `ModbusTcpConnectionFailedException` – if the connection fails.

**read\_holding\_registers(address: int, count: int) → List[int]**

Read specified number of register starting with given address and return the values from each register.

**Parameters**

- **address** – address of the first register
- **count** – count of registers to read

**Returns** list of `int` values

**read\_input\_registers(address: int, count: int) → List[int]**

Read specified number of register starting with given address and return the values from each register in a list.

**Parameters**

- **address** – address of the first register
- **count** – count of registers to read

**Returns** list of `int` values

**write\_registers(address: int, values: Union[List[int], int])**

Write values from the specified address forward.

**Parameters**

- **address** – address of the first register
- **values** – list with all values

```
class ModbusTcpCommunicationConfig(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address], unit:  
int, port: int = 502)
```

Bases: object

Configuration dataclass for [ModbusTcpCommunication](#).

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address]**

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port: int = 502**

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**unit: int**


```
exception ModbusTcpConnectionFailedException(string="")
```

Bases: pymodbus.exceptions.ConnectionException,

[hvl\\_ccb.comm.base.CommunicationException](#)

Exception raised when the connection failed.

## hvl\_ccb.comm.opc



Communication protocol implementing an OPC UA connection. This protocol is used to interface with the “Super-cube” PLC from Siemens.

**class Client(url: str, timeout: int = 4)**

Bases: `asyncua.sync.Client`

**disconnect()**

**get\_objects\_node()**

Get Objects node of client. Returns a Node object.

**property is\_open**

**send\_hello(\*args, \*\*kwargs)**

**class OpcUaCommunication(config)**

Bases: `hvl_ccb.comm.base.CommunicationProtocol`

Communication protocol implementing an OPC UA connection. Makes use of the package `python-opcua`.

**close() → None**

Close the connection to the OPC UA server.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**init\_monitored\_nodes(node\_id: Union[object, Iterable], ns\_index: int) → None**

Initialize monitored nodes.

### Parameters

- **node\_id** – one or more strings of node IDs; node IDs are always casted via `str()` method here, hence do not have to be strictly string objects.
- **ns\_index** – the namespace index the nodes belong to.

**Raises** `OpcUaCommunicationIOError` – when protocol was not opened or can't communicate with a OPC UA server

**property is\_open: bool**

Flag indicating if the communication port is open.

**Returns** `True` if the port is open, otherwise `False`

**open() → None**

Open the communication to the OPC UA server.

**Raises** `OpcUaCommunicationIOError` – when communication port cannot be opened.

### `read(node_id, ns_index)`

Read a value from a node with id and namespace index.

#### Parameters

- `node_id` – the ID of the node to read the value from
- `ns_index` – the namespace index of the node

**Returns** the value of the node object.

**Raises** `OpcUaCommunicationIOError` – when protocol was not opened or can't communicate with a OPC UA server

### `write(node_id, ns_index, value) → None`

Write a value to a node with name `name`.

#### Parameters

- `node_id` – the id of the node to write the value to.
- `ns_index` – the namespace index of the node.
- `value` – the value to write.

**Raises** `OpcUaCommunicationIOError` – when protocol was not opened or can't communicate with a OPC UA server

```
class OpcUaCommunicationConfig(host: typing.Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address],
                               endpoint_name: str, port: int = 4840, sub_handler:
                               hvl_ccb.comm.opc.OpcUaSubHandler =
                               <hvl_ccb.comm.opc.OpcUaSubHandler object>, update_parameter:
                               asyncua.ua.uaproto_col_auto.CreateSubscriptionParameters =
                               CreateSubscriptionParameters(RequestedPublishingInterval=1000,
                               RequestedLifetimeCount=300, RequestedMaxKeepAliveCount=22,
                               MaxNotificationsPerPublish=10000, PublishingEnabled=True, Priority=0),
                               wait_timeout_retry_sec: typing.Union[int, float] = 1,
                               max_timeout_retry_nr: int = 5)
```

Bases: `object`

Configuration dataclass for OPC UA Communciation.

### `clean_values()`

### `endpoint_name: str`

Endpoint of the OPC server, this is a path like ‘OPCUA/SimulationServer’

### `force_value(fieldname, value)`

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

#### Parameters

- `fieldname` – name of the field
- `value` – value to assign

### `host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address]`

Hostname or IP-Address of the OPC UA server.

### `is_configdataclass = True`

```
classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.

max_timeout_retry_nr: int = 5
    Maximal number of call re-tries on underlying OPC UA client timeout error

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.

port: int = 4840
    Port of the OPC UA server to connect to.

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

sub_handler: hvl_ccb.comm.opc.OpcUaSubHandler = <hvl_ccb.comm.opc.OpcUaSubHandler
object>
    object to use for handling subscriptions.

update_parameter: asyncua.ua.uaprotoauto.CreateSubscriptionParameters =
CreateSubscriptionParameters(RequestedPublishingInterval=1000,
RequestedLifetimeCount=300, RequestedMaxKeepAliveCount=22,
MaxNotificationsPerPublish=10000, PublishingEnabled=True, Priority=0)
    Values are given as a ua.CreateSubscriptionParameters as these parameters are requested by the OPC server.
    Other values will lead to an automatic revision of the parameters and a warning in the opc-logger, cf. MR !173

wait_timeout_retry_sec: Union[int, float] = 1
    Wait time between re-trying calls on underlying OPC UA client timeout error

exception OpcUaCommunicationIOError
    Bases: OSError, hvl_ccb.comm.base.CommunicationException
    OPC-UA communication I/O error.

exception OpcUaCommunicationTimeoutError
    Bases: hvl_ccb.comm.opc.OpcUaCommunicationIOError
    OPC-UA communication timeout error.

class OpcUaSubHandler
    Bases: object
    Base class for subscription handling of OPC events and data change events. Override methods from this class to add own handling capabilities.


    To receive events from server for a subscription data_change and event methods are called directly from receiving thread. Do not do expensive, slow or network operation there. Create another thread if you need to do such a thing.

        datachange_notification(node, val, data)
        event_notification(event)

class Server(shelf_file=None, tloop=None)
    Bases: asyncua.sync.Server
```

**get\_objects\_node()**  
Get Objects node of server. Returns a Node object.

## hvl\_ccb.comm.serial



Communication protocol for serial ports. Makes use of the [pySerial](#) library.

**class SerialCommunication(configuration)**  
Bases: [hvl\\_ccb.comm.base.AsyncCommunicationProtocol](#)

Implements the Communication Protocol for serial ports.

**close()**

Close the serial connection.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**property is\_open: bool**

Flag indicating if the serial port is open.

**Returns** *True* if the serial port is open, otherwise *False*

**open()**

Open the serial connection.

**Raises** [SerialCommunicationIOError](#) – when communication port cannot be opened.

**read\_bytes() → bytes**

Read the bytes from the serial port till the terminator is found. The input buffer may hold additional lines afterwards.

This method uses *self.access\_lock* to ensure thread-safety.

**Returns** Bytes read from the serial port; *b”* if there was nothing to read.

**Raises** [SerialCommunicationIOError](#) – when communication port is not opened

**read\_single\_bytes(size: int = 1) → bytes**

Read the specified number of bytes from the serial port. The input buffer may hold additional data afterwards.

**Returns** Bytes read from the serial port; *b”* if there was nothing to read.

```
write_bytes(data: bytes) → int
    Write bytes to the serial port.

    This method uses self.access_lock to ensure thread-safety.

    Parameters data – data to write to the serial port

    Returns number of bytes written

    Raises SerialCommunicationIOError – when communication port is not opened

class SerialCommunicationBytesize(value=<no_arg>, names=None, module=None, type=None, start=1,
                                         boundary=None)
Bases: hvl_ccb.utils.enum.ValueEnum

Serial communication bytesize.

EIGHTBITS = 8
FIVEBITS = 5
SEVENBITS = 7
SIXBITS = 6

class SerialCommunicationConfig(terminator: bytes = b'\r\n', encoding: str = 'utf-8',
                                         encoding_error_handling: str = 'strict', wait_sec_read_text_nonempty:
                                         Union[int, float] = 0.5, default_n_attempts_read_text_nonempty: int = 10,
                                         port: Optional[str] = None, baudrate: int = 9600, parity: Union[str,
                                         hvl_ccb.comm.serial.SerialCommunicationParity] =
                                         SerialCommunicationParity.NONE, stopbits: Union[int, float,
                                         hvl_ccb.comm.serial.SerialCommunicationStopbits] =
                                         SerialCommunicationStopbits.ONE, bytesize: Union[int,
                                         hvl_ccb.comm.serial.SerialCommunicationBytesize] =
                                         SerialCommunicationBytesize.EIGHTBITS, timeout: Union[int, float] = 2)
Bases: hvl_ccb.comm.base.AsyncCommunicationProtocolConfig

Configuration dataclass for SerialCommunication.

Bytesize
    alias of hvl_ccb.comm.serial.SerialCommunicationBytesize

Parity
    alias of hvl_ccb.comm.serial.SerialCommunicationParity

Stopbits
    alias of hvl_ccb.comm.serial.SerialCommunicationStopbits

baudrate: int = 9600
    Baudrate of the serial port

bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize] = 8
    Size of a byte, 5 to 8

clean_values()

create_serial_port() → serial.serialposix.Serial
    Create a serial port instance according to specification in this configuration

    Returns Closed serial port instance

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.
```

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity:** Union[str, [hvl\\_ccb.comm.serial.SerialCommunicationParity](#)] = 'N'

Parity to be used for the connection.

**port:** Optional[str] = None

Port is a string referring to a COM-port (e.g. 'COM3') or a URL. The full list of capabilities is found [on the pyserial documentation](#).

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits:** Union[int, float, [hvl\\_ccb.comm.serial.SerialCommunicationStopbits](#)] = 1

Stopbits setting, can be 1, 1.5 or 2.

**terminator\_str()** → str

**timeout:** Union[int, float] = 2

Timeout in seconds for the serial port

**exception SerialCommunicationIOError**

Bases: OSError, [hvl\\_ccb.comm.base.CommunicationException](#)

Serial communication related I/O errors.

**class SerialCommunicationParity**(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Serial communication parity.

**EVEN** = 'E'

**MARK** = 'M'

**NAMES** = {'E': 'Even', 'M': 'Mark', 'N': 'None', 'O': 'Odd', 'S': 'Space'}

**NONE** = 'N'

**ODD** = 'O'

**SPACE** = 'S'

```
class SerialCommunicationStopbits(value=<no_arg>, names=None, module=None, type=None, start=1,
                                     boundary=None)
```

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)


Serial communication stopbits.

**ONE** = 1

**ONE\_POINT\_FIVE** = 1.5

**TWO** = 2

## [hvl\\_ccb.comm.tcp](#)



Makes use of the socket library.

TCP communication protocol.

**class Tcp**(configuration)

Bases: [hvl\\_ccb.comm.base.CommunicationProtocol](#)

Tcp Communication Protocol.

**close()** → [None](#)

Close TCP connection.

**static config\_cls()** → Type[[hvl\\_ccb.comm.tcp.TcpCommunicationConfig](#)]

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**open()** → [None](#)

Open TCP connection.

**read()** → str

TCP read function :return: information read from TCP buffer formatted as string

```
write(command: str = "") → None
    TCP write function :param command: command string to be sent :return: none

class TcpCommunicationConfig(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address], port: int =
    54321, bufsize: int = 1024)
Bases: object
Configuration dataclass for TcpCommunication.

bufsize: int = 1024
clean_values()
force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing
    after value has been forced on fieldname.

Parameters
    • fieldname – name of the field
    • value – value to assign

host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address]
is_configdataclass = True
classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.


    Returns a list of strings containing all keys.

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified
    on instantiation.

    Returns a list of strings containing all optional keys.

port: int = 54321
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on
    instantiation.

    Returns a list of strings containing all required keys.
```

**hvl\_ccb.comm.telnet**

Communication protocol for telnet. Makes use of the [telnetlib](#) library.

**class TelnetCommunication(configuration)**  
Bases: [hvl\\_ccb.comm.base.AsyncCommunicationProtocol](#)

Implements the Communication Protocol for telnet.

**close()**

Close the telnet connection unless it is not closed.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**property is\_open: bool**

Is the connection open?

**Returns** True for an open connection

**open()**

Open the telnet connection unless it is not yet opened.

**read\_bytes() → bytes**

Read data as *bytes* from the telnet connection.

**Returns** data from telnet connection

**Raises** [TelnetError](#) – when connection is not open, raises an Error during the communication

**write\_bytes(data: bytes)**

Write the data as *bytes* to the telnet connection.

**Parameters** **data** – Data to be sent.

**Raises** [TelnetError](#) – when connection is not open, raises an Error during the communication

```
class TelnetCommunicationConfig(terminator: bytes = b'\r\n', encoding: str = 'utf-8',
                                encoding_error_handling: str = 'strict', wait_sec_read_text_nonempty:
                                Union[int, float] = 0.5, default_n_attempts_read_text_nonempty: int = 10,
                                host: Optional[Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address]] =
                                None, port: int = 0, timeout: Union[int, float] = 0.2)
```

Bases: `hvl_ccb.comm.base.AsyncCommunicationProtocolConfig`

Configuration dataclass for `TelnetCommunication`.

**clean\_values()**

**create\_telnet()** → `Optional[telnetlib.Telnet]`

Create a telnet client :return: Opened Telnet object or None if connection is not possible

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**host:** `Optional[Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address]] = None`

Host to connect to can be localhost or

**classmethod keys()** → `Sequence[str]`

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → `Dict[str, object]`

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port:** `int = 0`

Port at which the host is listening

**classmethod required\_keys()** → `Sequence[str]`

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**timeout:** `Union[int, float] = 0.2`


Timeout for reading a line

**exception TelnetError**

Bases: `OSSError, hvl_ccb.comm.base.CommunicationException`

Telnet communication related errors.

## `hvl_ccb.comm.visa`



Communication protocol for VISA. Makes use of the pyvisa library. The backend can be NI-Visa or pyvisa-py.

Information on how to install a VISA backend can be found here: [https://pyvisa.readthedocs.io/en/master/getting\\_nivisa.html](https://pyvisa.readthedocs.io/en/master/getting_nivisa.html)

So far only TCPIP SOCKET and TCPIP INSTR interfaces are supported.

**class VisaCommunication(configuration)**  
Bases: `hvl_ccb.comm.base.CommunicationProtocol`

Implements the Communication Protocol for VISA / SCPI.

**MULTI\_COMMANDS\_MAX = 5**

The maximum of commands that can be sent in one round is 5 according to the VISA standard.

**MULTI\_COMMANDS\_SEPARATOR = ';'**

The character to separate two commands is ; according to the VISA standard.

**WAIT\_AFTER\_WRITE = 0.08**

Small pause in seconds to wait after write operations, allowing devices to really do what we tell them before continuing with further tasks.

**close() → None**

Close the VISA connection and invalidates the handle.

**static config\_cls() → Type[hvl\_ccb.comm.visa.VisaCommunicationConfig]**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**open() → None**

Open the VISA connection and create the resource.

**query(\*commands: str) → Union[str, Tuple[str, ...]]**

A combination of write(message) and read.

**Parameters** `commands` – list of commands

**Returns** list of values

**Raises** `VisaCommunicationError` – when connection was not started, or when trying to issue too many commands at once.

**spoll()** → int

Execute serial poll on the device. Reads the status byte register STB. This is a fast function that can be executed periodically in a polling fashion.

**Returns** integer representation of the status byte

**Raises** `VisaCommunicationError` – when connection was not started

**write(\*commands: str)** → `None`

Write commands. No answer is read or expected.

**Parameters** `commands` – one or more commands to send

**Raises** `VisaCommunicationError` – when connection was not started

```
class VisaCommunicationConfig(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address],
                               interface_type: Union[str,
                               hvl_ccb.comm.visa.VisaCommunicationConfig.InterfaceType], board: int =
                               0, port: int = 5025, timeout: int = 5000, chunk_size: int = 204800,
                               open_timeout: int = 1000, write_termination: str = '\n', read_termination: str
                               = '\n', visa_backend: str = '')
```

Bases: `object`

*VisaCommunication* configuration dataclass.

**class InterfaceType**(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)

Bases: `hvl_ccb.utils.enum.AutoNumberNameEnum`

Supported VISA Interface types.

**TCPIP\_INSTR** = 2

VXI-11 protocol

**TCPIP\_SOCKET** = 1

VISA-RAW protocol

**address**(host: str, port: Optional[int] = None, board: Optional[int] = None) → str

Address string specific to the VISA interface type.

**Parameters**

- **host** – host IP address
- **port** – optional TCP port
- **board** – optional board number

**Returns** address string

**property address: str**

Address string depending on the VISA protocol's configuration.

**Returns** address string corresponding to current configuration

**board: int = 0**

Board number is typically 0 and comes from old bus systems.

**chunk\_size: int = 204800**

Chunk size is the allocated memory for read operations. The standard is 20kB, and is increased per default here to 200kB. It is specified in bytes.

**clean\_values()**

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**host:** Union[str, [ipaddress.IPv4Address](#), [ipaddress.IPv6Address](#)]

**interface\_type:** Union[str, [hvl\\_ccb.comm.visa.VisaCommunicationConfig.InterfaceType](#)]

Interface type of the VISA connection, being one of [InterfaceType](#).

**is\_configdataclass = True**

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**open\_timeout: int = 1000**

Timeout for opening the connection, in milli seconds.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port: int = 5025**

TCP port, standard is 5025.

**read\_termination: str = '\n'**

Read termination character.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**timeout: int = 5000**

Timeout for commands in milli seconds.

**visa\_backend: str = ''**

Specifies the path to the library to be used with PyVISA as a backend. Defaults to None, which is NI-VISA (if installed), or pyvisa-py (if NI-VISA is not found). To force the use of pyvisa-py, specify '@py' here.

**write\_termination: str = '\n'**

Write termination character.

**exception VisaCommunicationError**

Bases: OSError, [hvl\\_ccb.comm.base.CommunicationException](#)

Base class for VisaCommunication errors.

## Module contents

Communication protocols subpackage.


### `hvl_ccb.dev`

#### Subpackages

##### `hvl_ccb.dev.cube`

#### Submodules

##### `hvl_ccb.dev.cube.base`



Classes for the BaseCube device.

```

class BaseCube(com, dev_config=None)
    Bases: hvl_ccb.dev.base.SingleCommDevice

    Base class for Cube variants.

    OPC_MAX_YEAR = 2089
    OPC_MIN_YEAR = 1990

    active_alarms(human_readable: bool = True) → List[Union[int, str]]
        Displays all active alarms / messages.

        Parameters human_readable – True for human readable message, False for corresponding integer

        Returns list with active alarms
  
```

**property breakdown\_detection\_active: bool**

Get the state of the breakdown detection functionality. Returns True if it is enabled, False otherwise.

**Returns** state of the breakdown detection functionality

**breakdown\_detection\_reset() → None**

Reset the breakdown detection circuitry so that it is ready to detect breakdowns again.

**property breakdown\_detection\_triggered: bool**

See if breakdown detection unit has been triggered. Returns True if it is triggered, False otherwise.

**Returns** trigger status of the breakdown detection unit

**property cee16\_socket**

Read the on-state of the IEC CEE16 three-phase power socket.

**Returns** the on-state of the CEE16 power socket

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**classmethod datetime\_to\_opc(time\_dt: datetime.datetime) → List[int]**

Converts python datetime format into opc format (list of 8 integers) as defined in the following link: <https://support.industry.siemens.com/cs/mdm/109798671?c=133950752267&lc=de-WW> Each byte corresponds to one list entry. [yy, MM, dd, hh, mm, ss, milliseconds, weekday] Milliseconds and Weekday are not used, as this precision / information is not needed. The conversion of the numbers is special. Each decimal number is treated as it would be a hex-number and then converted back to decimal. This is tested with the used PLC in the BaseCube. yy: 0 to 99 (0 -> 2000, 89 -> 2089, 90 -> 1990, 99 -> 1999) MM: 1 to 12 dd: 1 to 31 hh: 0 to 23 mm: 0 to 59 ss: 0 to 59

**Parameters** **time\_dt** – time to be converted

**Returns** time in opc list format

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**display\_message\_board() → None**

Display 15 newest messages

**display\_status\_board() → None**

Display status board.

**door\_1\_status**

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**door\_2\_status**

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**door\_3\_status**

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**earthing\_rod\_1\_status**

Get the status of a earthing rod. See `constants.EarthingRodStatus` for possible returned earthing rod statuses.

**earthing\_rod\_2\_status**

Get the status of a earthing rod. See `constants.EarthingRodStatus` for possible returned earthing rod statuses.

**earthing\_rod\_3\_status**

Get the status of a earthing rod. See `constants.EarthingRodStatus` for possible returned earthing rod statuses.

**property operate: Optional[bool]**

Indicates if ‘operate’ is activated. ‘operate’ means locket safety circuit, red lamps, high voltage on and locked safety switches.

**Returns** `True` if `operate` is activated (RED\_OPERATE), `False` if `ready` is deactivated (RED\_READY), `None` otherwise

**quit\_error() → None**

Quits errors that are active on the Cube.

**read(node\_id: str)**

Local wrapper for the OPC UA communication protocol read method.

**Parameters** `node_id` – the id of the node to read.

**Returns** the value of the variable

**property ready: Optional[bool]**

Indicates if ‘ready’ is activated. ‘ready’ means locket safety circuit, red lamps, but high voltage still off.

**Returns** `True` if `ready` is activated (RED\_READY), `False` if `ready` is deactivated (GREEN\_READY), `None` otherwise

**set\_message\_board(msgs: List[str], display\_board: bool = True) → None**

Fills messages into message board that display that 15 newest messages with a timestamp.

**Parameters**

- `msgs` – list of strings
- `display_board` – display 15 newest messages if `True` (default)

**Raises** `ValueError` – if there are too many messages or the positions indices are invalid.

**set\_status\_board(msgs: List[str], pos: Optional[List[int]] = None, clear\_board: bool = True, display\_board: bool = True) → None**

Sets and displays a status board. The messages and the position of the message can be defined.

**Parameters**

- `msgs` – list of strings
- `pos` – list of integers [0...14]
- `clear_board` – clear unspecified lines if `True` (default), keep otherwise
- `display_board` – display new status board if `True` (default)

**Raises** `ValueError` – if there are too many messages or the positions indices are invalid.

**start() → None**

Starts the device. Sets the root node for all OPC read and write commands to the Siemens PLC object node which holds all our relevant objects and variables.

**property status: hvl\_ccb.dev.cube.constants.SafetyStatus**

Get the safety circuit status of the Cube. This methods is for the user.

**Returns** the safety status of the Cube’s state machine.

**stop() → None**

Stop the Cube device. Deactivates the remote control and closes the communication protocol.

**Raises** `CubeStopError` – when the cube is not in the correct status to stop the operation

**t13\_socket\_1**

Set and get the state of a SEV T13 power socket.

**t13\_socket\_2**

Set and get the state of a SEV T13 power socket.

**t13\_socket\_3**

Set and get the state of a SEV T13 power socket.

**write(node\_id, value) → None**

Local wrapper for the OPC UA communication protocol write method.

**Parameters**

- **node\_id** – the id of the node to write
- **value** – the value to write to the variable

```
class BaseCubeConfiguration(namespace_index: int = 3, polling_delay_sec: Union[int, float] = 5.0,
                            polling_interval_sec: Union[int, float] = 1.0, timeout_status_change:
                            Union[int, float] = 6, timeout_interval: Union[int, float] = 0.1,
                            noise_level_measurement_channel_1: Union[int, float] = 100,
                            noise_level_measurement_channel_2: Union[int, float] = 100,
                            noise_level_measurement_channel_3: Union[int, float] = 100,
                            noise_level_measurement_channel_4: Union[int, float] = 100)
```

Bases: `object`

Configuration dataclass for the BaseCube devices.

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**namespace\_index: int = 3**

Namespace of the OPC variables, typically this is 3 (coming from Siemens)

**noise\_level\_measurement\_channel\_1: Union[int, float] = 100**

**noise\_level\_measurement\_channel\_2: Union[int, float] = 100**

**noise\_level\_measurement\_channel\_3: Union[int, float] = 100**

**noise\_level\_measurement\_channel\_4: Union[int, float] = 100**

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**polling\_delay\_sec:** Union[int, float] = 5.0

**polling\_interval\_sec:** Union[int, float] = 1.0

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**timeout\_interval:** Union[int, float] = 0.1

**timeout\_status\_change:** Union[int, float] = 6

**class BaseCubeOpcUaCommunication(config)**

Bases: *hvl\_ccb.comm.opc.OpcUaCommunication*

Communication protocol specification for BaseCube devices.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class BaseCubeOpcUaCommunicationConfig(host: typing.Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address], endpoint\_name: hvl\_ccb.dev.cube.constants.\_CubeOpcEndpoint = \_CubeOpcEndpoint.BASE\_CUBE, port: int = 4840, sub\_handler: hvl\_ccb.comm.opc.OpcUaSubHandler = <hvl\_ccb.dev.cube.base.\_BaseCubeSubscriptionHandler object>, update\_parameter: asyncua.ua.uaproto\_col\_auto.CreateSubscriptionParameters = CreateSubscriptionParameters(RequestedPublishingInterval=1000, RequestedLifetimeCount=300, RequestedMaxKeepAliveCount=22, MaxNotificationsPerPublish=10000, PublishingEnabled=True, Priority=0), wait\_timeout\_retry\_sec: typing.Union[int, float] = 1, max\_timeout\_retry\_nr: int = 5)**

Bases: *hvl\_ccb.comm.opc.OpcUaCommunicationConfig*

Communication protocol configuration for OPC UA, specifications for the BaseCube devices.

**endpoint\_name: hvl\_ccb.dev.cube.constants.\_CubeOpcEndpoint = 'BaseCube'**

Endpoint of the OPC server, this is a path like 'OPCUA/SimulationServer'

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]


Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**sub\_handler:** *hvl\_ccb.comm.opc.OpcUaSubHandler* =  
<hvl\_ccb.dev.cube.base.\_BaseCubeSubscriptionHandler object>

Subscription handler for data change events

## [hvl\\_ccb.dev.cube.constants](#)



Constants, variable names for the BaseCube OPC-connected devices.

**exception CubeEarthlingStickOperationError**

Bases: *hvl\_ccb.dev.cube.constants.CubeError*

**exception CubeError**

Bases: *hvl\_ccb.dev.base.DeviceException*

**exception CubeRemoteControlError**

Bases: *hvl\_ccb.dev.cube.constants.CubeError*

**exception CubeStatusChangeError**

Bases: *hvl\_ccb.dev.cube.constants.CubeError*

**exception CubeStopError**

Bases: *hvl\_ccb.dev.cube.constants.CubeError*

**class DoorStatus**(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)

Bases: *aenum.IntEnum*

Possible status values for doors.

**CLOSED = 2**

Door is closed, but not locked.

**ERROR = 4**

Door has an error or was opened in locked state (either with emergency stop or from the inside).

**INACTIVE = 0**

not enabled in BaseCube HMI setup, this door is not supervised.

**LOCKED = 3**

Door is closed and locked (safe state).

---

```

OPEN = 1
Door is open.

class EarthingRodStatus(value=<no_arg>, names=None, module=None, type=None, start=1,
boundary=None)
Bases: aenum.IntEnum

Possible status values for earthing rods.

EXPERIMENT_BLOCKED = 0
earthing rod is somewhere in the experiment and blocks the start of the experiment

EXPERIMENT_READY = 1
earthing rod is hanging next to the door, experiment is ready to operate

class EarthingStickOperatingStatus(value=<no_arg>, names=None, module=None, type=None, start=1,
boundary=None)
Bases: aenum.IntEnum

Operating Status for an earthing stick. Stick can be used in auto or manual mode.

AUTO = 0

MANUAL = 1

class EarthingStickOperation(value=<no_arg>, names=None, module=None, type=None, start=1,
boundary=None)
Bases: hvl_ccb.utils.enum.BoolEnum

Operation of the earthing stick in manual operating mode. Can be closed or opened.

CLOSE = True

OPEN = False

class EarthingStickStatus(value=<no_arg>, names=None, module=None, type=None, start=1,
boundary=None)
Bases: aenum.IntEnum

Status of an earthing stick. These are the possible values in the status integer e.g. in _EarthingStick.status.

CLOSED = 1

ERROR = 3

INACTIVE = 0

OPEN = 2

class MessageBoard(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: hvl_ccb.dev.cube.constants._LineEnumBase

Variable NodeID strings for message board lines.

LINE_1 = '"DB_OPC_Connection"."Is_status_Line_1"
LINE_10 = '"DB_OPC_Connection"."Is_status_Line_10"'
LINE_11 = '"DB_OPC_Connection"."Is_status_Line_11"'
LINE_12 = '"DB_OPC_Connection"."Is_status_Line_12"'
LINE_13 = '"DB_OPC_Connection"."Is_status_Line_13"'
LINE_14 = '"DB_OPC_Connection"."Is_status_Line_14"'

```

```
LINE_15 = '"DB_OPC_Connection"."Is_status_Line_15"'
LINE_2 = '"DB_OPC_Connection"."Is_status_Line_2"'
LINE_3 = '"DB_OPC_Connection"."Is_status_Line_3"'
LINE_4 = '"DB_OPC_Connection"."Is_status_Line_4"'
LINE_5 = '"DB_OPC_Connection"."Is_status_Line_5"'
LINE_6 = '"DB_OPC_Connection"."Is_status_Line_6"'
LINE_7 = '"DB_OPC_Connection"."Is_status_Line_7"'
LINE_8 = '"DB_OPC_Connection"."Is_status_Line_8"'
LINE_9 = '"DB_OPC_Connection"."Is_status_Line_9"'

exception PICubeTestParameterError
    Bases: hvl\_ccb.dev\(cube.constants.CubeError

class Polarity(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
    Bases: aenum.IntEnum

    An enumeration.

    NEGATIVE = 0
    POSITIVE = 1

class PowerSetup(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
    Bases: aenum.IntEnum

    Possible power setups corresponding to the value of variable Power.setup. The values for slope_min are experimentally defined, below these values the slope is more like a staircase

    The name of the first argument needs to be ‘value’, otherwise the IntEnum is not working correctly.

    AC_100KV = 3
    AC_150KV = 4
    AC_200KV = 5
    AC_50KV = 2
    DC_140KV = 7
    DC_280KV = 8
    EXTERNAL_SOURCE = 1
    IMPULSE_140KV = 9
    NO_SOURCE = 0
    POWER_INVERTER_220V = 6

STOP_SAFETY_STATUSES: Tuple[hvl\_ccb.dev\(cube.constants.SafetyStatus, ...] =
(<SafetyStatus.GREEN_NOT_READY: 1>, <SafetyStatus.GREEN_READY: 2>)

    BaseCube’s safety statuses required to close the connection to the device.

class SafetyStatus(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
    Bases: aenum.IntEnum

    Safety status values that are possible states returned from hvl\_ccb.dev\(cube.base.BaseCube.status\(\). These values correspond to the states of the BaseCube’s safety circuit statemachine.


    ERROR = 6
```

```

GREEN_NOT_READY = 1
GREEN_READY = 2
INITIALIZING = 0
QUICK_STOP = 5
RED_OPERATE = 4
RED_READY = 3

```

## `hvl_ccb.dev(cube.picube)`



A PICube is a BaseCube with build in Power Inverter

```

class PICube(com, dev_config=None)
    Bases: hvl_ccb.dev(cube.base.BaseCube)

```

Variant of the BaseCube with build in Power Inverter

```

static config_cls()
    Return the default configdataclass class.

```

**Returns** a reference to the default configdataclass class

```

property current_primary: float

```

Read the current primary current at the output of the frequency converter (before transformer).

**Returns** primary current in A

```

static default_com_cls()

```

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

```

property frequency: float

```

Read the electrical frequency of the current PICube setup.

**Returns** the frequency in Hz

**property operate: Optional[bool]**

Indicates if ‘operate’ is activated. ‘operate’ means locket safety circuit, red lamps, high voltage on and locked safety switches.

**Returns** *True* if operate is activated (RED\_OPERATE), *False* if ready is deactivated (RED\_READY), *None* otherwise

**property polarity: Optional[hvl\_ccb.dev.cube.constants.Polarity]**

Polarity of a DC setup. :return: if a DC setup is programmed the polarity is returned, else None.

**property power\_setup: hvl\_ccb.dev.cube.constants.PowerSetup**

Return the power setup selected in the PICube’s settings.

**Returns** the power setup

**property voltage\_actual: float**

Reads the actual measured voltage and returns the value in V.

**Returns** the actual voltage of the setup in V.

**property voltage\_max: float**

Reads the maximum voltage of the setup and returns in V.

**Returns** the maximum voltage of the setup in V.

**property voltage\_primary: float**

Read the current primary voltage at the output of the frequency converter (before transformer).

**Returns** primary voltage in V

```
class PICubeConfiguration(namespace_index: int = 3, polling_delay_sec: Union[int, float] = 5.0,
                           polling_interval_sec: Union[int, float] = 1.0, timeout_status_change: Union[int,
                           float] = 6, timeout_interval: Union[int, float] = 0.1,
                           noise_level_measurement_channel_1: Union[int, float] = 100,
                           noise_level_measurement_channel_2: Union[int, float] = 100,
                           noise_level_measurement_channel_3: Union[int, float] = 100,
                           noise_level_measurement_channel_4: Union[int, float] = 100,
                           timeout_test_parameters: 'Number' = 2.0)
```

Bases: *hvl\_ccb.dev.cube.base.BaseCubeConfiguration*

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**timeout\_test\_parameters: Union[int, float] = 2.0**

**class PICubeOpcUaCommunication(config)**

Bases: *hvl\_ccb.dev.cube.base.BaseCubeOpcUaCommunication*

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class PICubeOpcUaCommunicationConfig(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address],**

```
endpoint_name: '_CubeOpcEndpoint' =
    <_CubeOpcEndpoint.PI_CUBE: 'PICube'>, port: int = 4840,
    sub_handler: hvl_ccb.comm.opc.OpcUaSubHandler =
        <hvl_ccb.dev.cube.base._BaseCubeSubscriptionHandler object at
        0x7f295edb9310>, update_parameter:
            asynccua.ua.uaprotoauto.CreateSubscriptionParameters =
                CreateSubscriptionParameters(RequestedPublishingInterval=1000,
                RequestedLifetimeCount=300, RequestedMaxKeepAliveCount=22,
                MaxNotificationsPerPublish=10000, PublishingEnabled=True,
                Priority=0), wait_timeout_retry_sec: Union[int, float] = 1,
                max_timeout_retry_nr: int = 5)
```

Bases: *hvl\_ccb.dev.cube.base.BaseCubeOpcUaCommunicationConfig*

**endpoint\_name: \_CubeOpcEndpoint = 'PICube'**

Endpoint of the OPC server, this is a path like ‘OPCUA/SimulationServer’

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.


## Module contents

Cube package with implementation for system versions from 2019 on (new concept with hard-PLC Siemens S7-1500 as CPU).

### [hvl\\_ccb.dev.highland\\_t560](#)

#### Submodules

##### [hvl\\_ccb.dev.highland\\_t560.base](#)



Module containing base device and communication classes and enums.

Communication with device is performed via its ethernet port and a Telnet connection.

**class AutoInstallMode**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)  
Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Modes for installing configuration settings to the device.

**INSTALL** = 1

**OFF** = 0

**QUEUE** = 2

**class GateMode**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)  
Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Available T560 gate modes

**INPUT** = 'INP'

**OFF** = 'OFF'

**OUTPUT** = 'OUT'

**class Polarity**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)  
Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Possible channel polarity states

```

ACTIVE_HIGH = 'POS'
ACTIVE_LOW = 'NEG'

class T560Communication(configuration)
    Bases:           hvl_ccb.comm.base.SyncCommunicationProtocol,      hvl_ccb.comm.telnet.
                  TelnetCommunication

    Communication class for T560. It uses a TelnetCommunication with the SyncCommunicationProtocol

    static config_cls()
        Return the default configdataclass class.

        Returns a reference to the default configdataclass class

    query(command: str) → str
        Send a command to the device and handle the response.

        For device setting queries, response will be ‘OK’ if successful, or ‘??’ if setting cannot be carried out,
        raising an error.

        Parameters command – Command string to be sent
        Raises T560Error – if no response is received, or if the device responds with an error message.
        Returns Response from the device.

class T560CommunicationConfig(terminator: bytes = b'\r', encoding: str = 'utf-8', encoding_error_handling:
                                str = 'strict', wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                default_n_attempts_read_text_nonempty: int = 10, host: Union[str,
                                ipaddress.IPv4Address, ipaddress.IPv6Address, NoneType] = None, port: int
                                = 2000, timeout: Union[int, float] = 0.2)
    Bases:           hvl_ccb.comm.base.SyncCommunicationProtocolConfig,      hvl_ccb.comm.telnet.
                  TelnetCommunicationConfig

    force_value(fieldname, value)
        Forces a value to a dataclass field despite the class being frozen.

        NOTE: you can define post_force_value method with same signature as this method to do extra processing
        after value has been forced on fieldname.

        Parameters
            • fieldname – name of the field
            • value – value to assign

    classmethod keys() → Sequence[str]
        Returns a list of all configdataclass fields key-names.

        Returns a list of strings containing all keys.

    classmethod optional_defaults() → Dict[str, object]
        Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified
        on instantiation.

        Returns a list of strings containing all optional keys.

    port: int = 2000
        Port at which the host is listening

    classmethod required_keys() → Sequence[str]
        Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on
        instantiation.

```

**Returns** a list of strings containing all required keys.

**terminator: bytes = b'\r'**

The terminator character. Typically this is b'\r\n' or b'\n', but can also be b'\r' or other combinations.  
This defines the end of a single line.

**exception T560Error**

Bases: *hvl\_ccb.dev.base.DeviceException*

T560 related errors.

**class TriggerMode(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**

Bases: *hvl\_ccb.utils.enum.ValueEnum*

Available T560 trigger modes

**COMMAND = 'REM'**

**EXT\_FALLING\_EDGE = 'NEG'**


**EXT\_RISING\_EDGE = 'POS'**

**INT\_SYNTHESIZER = 'SYN'**

**OFF = 'OFF'**

***hvl\_ccb.dev.highland\_t560.channel***

Module for controlling pulse output channels A, B, C and D.

**hvl\_ccb.dev.highland\_t560.device**

Module for controlling device, including TRIG, CLOCK and GATE I/Os.

```

class T560(com, dev_config=None)
  Bases: hvl_ccb.dev.base.SingleCommDevice

  activate_clock_output()
    Outputs 10 MHz clock signal

  property auto_install_mode: hvl_ccb.dev.highland_t560.base.AutoInstallMode
    Check the autoinstall settings of the T560. The autoinstall mode sets how changes to device settings are applied. See manual section 4.7.2 for more information about these modes.

  property ch_a: hvl_ccb.dev.highland_t560.channel._Channel
    Channel A of T560

  property ch_b: hvl_ccb.dev.highland_t560.channel._Channel
    Channel B of T560

  property ch_c: hvl_ccb.dev.highland_t560.channel._Channel
    Channel C of T560

  property ch_d: hvl_ccb.dev.highland_t560.channel._Channel
    Channel D of T560

  static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class
  
```

```
static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

disarm_trigger()
    Disarm DDG by disabling all trigger sources.

fire_trigger()
    Fire a software trigger.

property frequency: float
    The frequency of the timing cycle in Hz.

property gate_mode: hvl_ccb.dev.highland_t560.base.GateMode
    Check the mode setting of the GATE I/O port.

property gate_polarity: hvl_ccb.dev.highland_t560.base.Polarity
    Check the polarity setting of the GATE I/O port.

load_device_configuration()
    Load the settings saved in nonvolatile memory.

property period: float
    The period of the timing cycle (time between triggers) in seconds.

save_device_configuration()
    Save the current settings to nonvolatile memory.

property trigger_level
    Get external trigger level.

property trigger_mode
    Get device trigger source.

use_external_clock()
    Finds and accepts an external clock signal to the CLOCK input

class T560Config
    Bases: object

    auto_install_mode = 1

    clean_values()
        Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

    force_value(fieldname, value)
        Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.

    Parameters
        • fieldname – name of the field
        • value – value to assign

    is_configdataclass = True

    classmethod keys() → Sequence[str]
        Returns a list of all configdataclass fields key-names.

        Returns a list of strings containing all keys.
```

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## Module contents

This module establishes methods for interfacing with the Highland Technology T560-2 via its ethernet adapter with a telnet communication protocol.

The T560 is a small digital delay & pulse generator. It outputs up to four individually timed pulses with 10-ps precision, given an internal or external trigger.


This module introduces methods for configuring channels, gating, and triggering. Further documentation and a more extensive command list may be obtained from:

<https://www.highlandtechnology.com/DSS/T560DS.shtml>

## hvl\_ccb.dev.supercube2015

### Submodules

#### hvl\_ccb.dev.supercube2015.base



Base classes for the Supercube device.

**exception InvalidSupercubeStatusError**

Bases: `Exception`

Exception raised when supercube has invalid status.

**class Supercube2015Base**(*com, dev\_config=None*)

Bases: `hvl_ccb.dev.base.SingleCommDevice`

Base class for Supercube variants.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**get\_cee16\_socket() → bool**

Read the on-state of the IEC CEE16 three-phase power socket.

**Returns** the on-state of the CEE16 power socket

**get\_door\_status(door: int) → hvl\_ccb.dev.supercube2015.constants.DoorStatus**

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**Parameters** `door` – the door number (1..3)

**Returns** the door status

**get\_earthing\_manual(number: int) → bool**

Get the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters** `number` – number of the earthing stick (1..6)

**Returns** earthing stick manual status

**get\_earthing\_status(number: int) → int**

Get the status of an earthing stick, whether it is closed, open or undefined (moving).

**Parameters** `number` – number of the earthing stick (1..6)

**Returns** earthing stick status; see `constants.EarthlingStickStatus`

**get\_measurement\_ratio(channel: int) → float**

Get the set measurement ratio of an AC/DC analog input channel. Every input channel has a divider ratio assigned during setup of the Supercube system. This ratio can be read out.

**Attention:** Supercube 2015 does not have a separate ratio for every analog input. Therefore there is only one ratio for `channel = 1`.

**Parameters** `channel` – number of the input channel (1..4)

**Returns** the ratio

**get\_measurement\_voltage(channel: int) → float**

Get the measured voltage of an analog input channel. The voltage read out here is already scaled by the configured divider ratio.

**Attention:** In contrast to the *new* Supercube, the old one returns here the input voltage read at the ADC. It is not scaled by a factor.

**Parameters** `channel` – number of the input channel (1..4)

**Returns** measured voltage

**get\_status() → int**

Get the safety circuit status of the Supercube.

**Returns** the safety status of the supercube's state machine; see `constants.SafetyStatus`.

**get\_support\_input(port: int, contact: int) → bool**

Get the state of a support socket input.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital input read state**get\_support\_output**(*port: int, contact: int*) → bool

Get the state of a support socket output.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital output read state**get\_t13\_socket**(*port: int*) → bool

Read the state of a SEV T13 power socket.

**Parameters** **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS***Returns** on-state of the power socket**horn**(*state: bool*) → *None*

Turns acoustic horn on or off.

**Parameters** **state** – Turns horn on (True) or off (False)**operate**(*state: bool*) → *None*

Set operate state. If the state is RedReady, this will turn on the high voltage and close the safety switches.

**Parameters** **state** – set operate state**quit\_error**() → *None*

Quits errors that are active on the Supercube.

**read**(*node\_id: str*)

Local wrapper for the OPC UA communication protocol read method.

**Parameters** **node\_id** – the id of the node to read.**Returns** the value of the variable**ready**(*state: bool*) → *None*

Set ready state. Ready means locket safety circuit, red lamps, but high voltage still off.

**Parameters** **state** – set ready state**set\_cee16\_socket**(*state: bool*) → *None*

Switch the IEC CEE16 three-phase power socket on or off.

**Parameters** **state** – desired on-state of the power socket**Raises** **TypeError** – if state is not of type bool**set\_earthing\_manual**(*number: int, manual: bool*) → *None*

Set the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters**

- **number** – number of the earthing stick (1..6)
- **manual** – earthing stick manual status (True or False)

**set\_remote\_control**(state: bool) → *None*

Enable or disable remote control for the Supercube. This will effectively display a message on the touch-screen HMI.

**Parameters** **state** – desired remote control state

**set\_support\_output**(port: int, contact: int, state: bool) → *None*

Set the state of a support output socket.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **state** – is the desired state of the support output

**Raises** **TypeError** – when state is not of type bool

**set\_support\_output\_impulse**(port: int, contact: int, duration: float = 0.2, pos\_pulse: bool = True) → *None*

Issue an impulse of a certain duration on a support output contact. The polarity of the pulse (On-wait-Off or Off-wait-On) is specified by the pos\_pulse argument.

This function is blocking.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **duration** – is the length of the impulse in seconds
- **pos\_pulse** – is True, if the pulse shall be HIGH, False if it shall be LOW

**Raises** **TypeError** – when state is not of type bool

**set\_t13\_socket**(port: int, state: bool) → *None*

Set the state of a SEV T13 power socket.

**Parameters**

- **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*
- **state** – is the desired on-state of the socket

**Raises**

- **ValueError** – when port is not valid
- **TypeError** – when state is not of type bool

**start()** → *None*

Starts the device. Sets the root node for all OPC read and write commands to the Siemens PLC object node which holds all our relevant objects and variables.

**stop()** → *None*

Stop the Supercube device. Deactivates the remote control and closes the communication protocol.

**write(node\_id, value)** → *None*

Local wrapper for the OPC UA communication protocol write method.

**Parameters**

- **node\_id** – the id of the node to read
- **value** – the value to write to the variable

```
class SupercubeConfiguration(namespace_index: int = 7)
Bases: object

    Configuration dataclass for the Supercube devices.

    clean_values()
        Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

    force_value(fieldname, value)
        Forces a value to a dataclass field despite the class being frozen.

        NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.

    Parameters
        • fieldname – name of the field
        • value – value to assign

    is_configdataclass = True

    classmethod keys() → Sequence[str]
        Returns a list of all configdataclass fields key-names.

        Returns a list of strings containing all keys.

    namespace_index: int = 7
        Namespace of the OPC variables, typically this is 3 (coming from Siemens)

    classmethod optional_defaults() → Dict[str, object]
        Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

        Returns a list of strings containing all optional keys.

    classmethod required_keys() → Sequence[str]
        Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

        Returns a list of strings containing all required keys.

class SupercubeOpcUaCommunication(config)
Bases: hvl_ccb.comm.opc.OpcUaCommunication

    Communication protocol specification for Supercube devices.

    static config_cls()
        Return the default configdataclass class.

        Returns a reference to the default configdataclass class
```

```
class SupercubeOpcUaCommunicationConfig(host: typing.Union[str, ipaddress.IPv4Address,
                                                ipaddress.IPv6Address], endpoint_name: str, port: int = 4845,
                                                sub_handler: hvl_ccb.comm.opc.OpcUaSubHandler =
                                                <hvl_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler
                                                object>, update_parameter:
                                                asynccua.ua.uaprotoauto.CreateSubscriptionParameters =
                                                CreateSubscriptionParameters(
                                                RequestedPublishingInterval=1000,
                                                RequestedLifetimeCount=300,
                                                RequestedMaxKeepAliveCount=22,
                                                MaxNotificationsPerPublish=10000, PublishingEnabled=True,
                                                Priority=0), wait_timeout_retry_sec: typing.Union[int, float] =
                                                1, max_timeout_retry_nr: int = 5)
```

Bases: [hvl\\_ccb.comm.opc.OpcUaCommunicationConfig](#)

Communication protocol configuration for OPC UA, specifications for the Supercube devices.

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port: int = 4845**

Port of the OPC UA server to connect to.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**sub\_handler: hvl\_ccb.comm.opc.OpcUaSubHandler =**  
**<hvl\_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler object>**

Subscription handler for data change events

**class SupercubeSubscriptionHandler**

Bases: [hvl\\_ccb.comm.opc.OpcUaSubHandler](#)

OPC Subscription handler for datachange events and normal events specifically implemented for the Supercube devices.


**datachange\_notification**(*node*: asynccua.sync.SyncNode, *val*, *data*)

In addition to the standard operation (debug logging entry of the datachange), alarms are logged at INFO level using the alarm text.

## Parameters

- **node** – the node object that triggered the datachange event
- **val** – the new value
- **data** –

## `hvl_ccb.dev.supercube2015.constants`



Constants, variable names for the Supercube OPC-connected devices.

```
class AlarmText(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: hvl_ccb.utils.enum.ValueEnum
```

This enumeration contains textual representations for all error classes (stop, warning and message) of the Supercube system. Use the `AlarmText.get()` method to retrieve the enum of an alarm number.

```

Alarm0 = 'No Alarm.'
Alarm1 = 'STOP Safety switch 1 error'
Alarm10 = 'STOP Earthing stick 2 error'
Alarm11 = 'STOP Earthing stick 3 error'
Alarm12 = 'STOP Earthing stick 4 error'
Alarm13 = 'STOP Earthing stick 5 error'
Alarm14 = 'STOP Earthing stick 6 error'
Alarm17 = 'STOP Source switch error'
Alarm19 = 'STOP Fence 1 error'
Alarm2 = 'STOP Safety switch 2 error'
Alarm20 = 'STOP Fence 2 error'
Alarm21 = 'STOP Control error'
Alarm22 = 'STOP Power outage'
Alarm3 = 'STOP Emergency Stop 1'
Alarm4 = 'STOP Emergency Stop 2'
Alarm5 = 'STOP Emergency Stop 3'
Alarm6 = 'STOP Door 1 lock supervision'
Alarm7 = 'STOP Door 2 lock supervision'
Alarm8 = 'STOP Door 3 lock supervision'
Alarm9 = 'STOP Earthing stick 1 error'
```

```
classmethod get(alarm: int)
    Get the attribute of this enum for an alarm number.

    Parameters alarm – the alarm number

    Returns the enum for the desired alarm number

not_defined = 'NO ALARM TEXT DEFINED'

class BreakdownDetection(value=<no_arg>, names=None, module=None, type=None, start=1,
                           boundary=None)
Bases: hvl\_ccb.utils.enum.ValueEnum

Node ID strings for the breakdown detection.

activated = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.Breakdowndetection.connect'
    Boolean read-only variable indicating whether breakdown detection and fast switchoff is enabled in the system or not.

reset = 'hvl-ipc.WINAC.Support6OutA'
    Boolean writable variable to reset the fast switch-off. Toggle to re-enable.

triggered = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.Breakdowndetection.triggered'
    Boolean read-only variable telling whether the fast switch-off has triggered. This can also be seen using the safety circuit state, therefore no method is implemented to read this out directly.

class DoorStatus(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: aenum.IntEnum

Possible status values for doors.

closed = 2
    Door is closed, but not locked.

error = 4
    Door has an error or was opened in locked state (either with emergency stop or from the inside).

inactive = 0
    not enabled in Supercube HMI setup, this door is not supervised.

locked = 3
    Door is closed and locked (safe state).

open = 1
    Door is open.

class EarthingStick(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: hvl\_ccb.utils.enum.ValueEnum

Variable NodeID strings for all earthing stick statuses (read-only integer) and writable booleans for setting the earthing in manual mode.

classmethod manual(number: int)
    Get the manual enum attribute for an earthing stick number.

    Parameters number – the earthing stick (1..6)

    Returns the manual enum

manual_1 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.MANUAL'
manual_2 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.MANUAL'
manual_3 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.MANUAL'
manual_4 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.MANUAL'
```

```

manual_5 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.MANUAL'
manual_6 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.MANUAL'
status_1_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.CLOSE'
status_1_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.CONNECT'
status_1_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.OPEN'
status_2_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.CLOSE'
status_2_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.CONNECT'
status_2_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.OPEN'
status_3_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.CLOSE'
status_3_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.CONNECT'
status_3_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.OPEN'
status_4_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.CLOSE'
status_4_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.CONNECT'
status_4_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.OPEN'
status_5_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.CLOSE'
status_5_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.CONNECT'
status_5_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.OPEN'
status_6_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.CLOSE'
status_6_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.CONNECT'
status_6_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.OPEN'

classmethod status_closed(number: int)
    Get the status enum attribute for an earthing stick number.

    Parameters number – the earthing stick (1..6)

    Returns the status enum

classmethod status_connected(number: int)
    Get the status enum attribute for an earthing stick number.

    Parameters number – the earthing stick (1..6)

    Returns the status enum

classmethod status_open(number: int)
    Get the status enum attribute for an earthing stick number.

    Parameters number – the earthing stick (1..6)

    Returns the status enum

class EarthingStickStatus(value=<no_arg>, names=None, module=None, type=None, start=1,
                           boundary=None)
    Bases: aenum.IntEnum

    Status of an earthing stick. These are the possible values in the status integer e.g. in EarthingStick.status_1.

```

**closed = 1**

Earthing is closed (safe).

**error = 3**

Earthing is in error, e.g. when the stick did not close correctly or could not open.

**inactive = 0**

Earthing stick is deselected and not enabled in safety circuit. To get out of this state, the earthing has to be enabled in the Supercube HMI setup.

**open = 2**

Earthing is open (not safe).

**class Errors**(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Variable NodeID strings for information regarding error, warning and message handling.

**quit = 'hvl-ipc.WINAC.SYSTEMSTATE.Faultconfirmation'**

Writable boolean for the error quit button.

**stop = 'hvl-ipc.WINAC.SYSTEMSTATE.ERROR'**

Boolean read-only variable telling if a stop is active.

**stop\_number = 'hvl-ipc.WINAC.SYSTEMSTATE.Errornumber'**

**class GeneralSockets**(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

NodeID strings for the power sockets (3x T13 and 1xCEE16).

**cee16 = 'hvl-ipc.WINAC.SYSTEM\_COMPONENTS.CEE16'**

CEE16 socket (writeable boolean).

**t13\_1 = 'hvl-ipc.WINAC.SYSTEM\_COMPONENTS.T13\_1'**

SEV T13 socket No. 1 (writable boolean).

**t13\_2 = 'hvl-ipc.WINAC.SYSTEM\_COMPONENTS.T13\_2'**

SEV T13 socket No. 2 (writable boolean).

**t13\_3 = 'hvl-ipc.WINAC.SYSTEM\_COMPONENTS.T13\_3'**

SEV T13 socket No. 3 (writable boolean).

**class GeneralSupport**(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

NodeID strings for the support inputs and outputs.

**in\_1\_1 = 'hvl-ipc.WINAC.Support1InA'**

**in\_1\_2 = 'hvl-ipc.WINAC.Support1InB'**

**in\_2\_1 = 'hvl-ipc.WINAC.Support2InA'**

**in\_2\_2 = 'hvl-ipc.WINAC.Support2InB'**

**in\_3\_1 = 'hvl-ipc.WINAC.Support3InA'**

**in\_3\_2 = 'hvl-ipc.WINAC.Support3InB'**

**in\_4\_1 = 'hvl-ipc.WINAC.Support4InA'**

**in\_4\_2 = 'hvl-ipc.WINAC.Support4InB'**

**in\_5\_1 = 'hvl-ipc.WINAC.Support5InA'**

**in\_5\_2 = 'hvl-ipc.WINAC.Support5InB'**

```
in_6_1 = 'hvl-ipc.WINAC.Support6InA'  
in_6_2 = 'hvl-ipc.WINAC.Support6InB'  
classmethod input(port, contact)  
    Get the NodeID string for a support input.
```

**Parameters**

- **port** – the desired port (1..6)
- **contact** – the desired contact at the port (1..2)

**Returns** the node id string

```
out_1_1 = 'hvl-ipc.WINAC.Support1OutA'  
out_1_2 = 'hvl-ipc.WINAC.Support1OutB'  
out_2_1 = 'hvl-ipc.WINAC.Support2OutA'  
out_2_2 = 'hvl-ipc.WINAC.Support2OutB'  
out_3_1 = 'hvl-ipc.WINAC.Support3OutA'  
out_3_2 = 'hvl-ipc.WINAC.Support3OutB'  
out_4_1 = 'hvl-ipc.WINAC.Support4OutA'  
out_4_2 = 'hvl-ipc.WINAC.Support4OutB'  
out_5_1 = 'hvl-ipc.WINAC.Support5OutA'  
out_5_2 = 'hvl-ipc.WINAC.Support5OutB'  
out_6_1 = 'hvl-ipc.WINAC.Support6OutA'  
out_6_2 = 'hvl-ipc.WINAC.Support6OutB'
```

**classmethod** **output**(*port, contact*)

Get the NodeID string for a support output.

**Parameters**

- **port** – the desired port (1..6)
- **contact** – the desired contact at the port (1..2)

**Returns** the node id string

```
class MeasurementsDividerRatio(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
```

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Variable NodeID strings for the measurement input scaling ratios. These ratios are defined in the Supercube HMI setup and are provided in the python module here to be able to read them out, allowing further calculations.

**classmethod** **get**(*channel: int*)

Get the attribute for an input number.

**Parameters** **channel** – the channel number (1..4)**Returns** the enum for the desired channel.

```
input_1 = 'hvl-ipc.WINAC.SYSTEM_INTERN.DivididerRatio'
```

```
class MeasurementsScaledInput(value=<no_arg>, names=None, module=None, type=None, start=1,
                                boundary=None)
```

Bases: `hvl_ccb.utils.enum.ValueEnum`

Variable NodeID strings for the four analog BNC inputs for measuring voltage. The voltage returned in these variables is already scaled with the set ratio, which can be read using the variables in `MeasurementsDividerRatio`.

```
classmethod get(channel: int)
```

Get the attribute for an input number.

**Parameters** `channel` – the channel number (1..4)

**Returns** the enum for the desired channel.

```
input_1 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI1Volt'
```

```
input_2 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI2Volt'
```

```
input_3 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI3Volt'
```

```
input_4 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI4Volt'
```

```
class Power(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
```

Bases: `hvl_ccb.utils.enum.ValueEnum`

Variable NodeID strings concerning power data.

```
current_primary = 'hvl-ipc.WINAC.SYSTEM_INTERN.FUCurrentprim'
```

Primary current in ampere, measured by the frequency converter. (read-only)

```
frequency = 'hvl-ipc.WINAC.FU.Frequency'
```

Frequency converter output frequency. (read-only)

```
setup = 'hvl-ipc.WINAC.FU.TrafoSetup'
```

Power setup that is configured using the Supercube HMI. The value corresponds to the ones in `PowerSetup`.  
(read-only)

```
voltage_max = 'hvl-ipc.WINAC.FU.maxVoltagekV'
```

Maximum voltage allowed by the current experimental setup. (read-only)

```
voltage_primary = 'hvl-ipc.WINAC.SYSTEM_INTERN.FUVoltageprim'
```

Primary voltage in volts, measured by the frequency converter at its output. (read-only)

```
voltage_slope = 'hvl-ipc.WINAC.FU.dUdt_-1'
```

Voltage slope in V/s.

```
voltage_target = 'hvl-ipc.WINAC.FU.SOLL'
```

Target voltage setpoint in V.

```
class PowerSetup(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
```

Bases: `aenum.IntEnum`

Possible power setups corresponding to the value of variable `Power.setup`.

```
AC_DoubleStage_150kV = 3
```

AC voltage with two MWB transformers, one at 100kV and the other at 50kV, resulting in a total maximum voltage of 150kV.

```
AC_DoubleStage_200kV = 4
```

AC voltage with two MWB transformers both at 100kV, resulting in a total maximum voltage of 200kV

```
AC_SingleStage_100kV = 2
```

AC voltage with MWB transformer set to 100kV maximum voltage.

**AC\_SingleStage\_50kV = 1**

AC voltage with MWB transformer set to 50kV maximum voltage.

**DC\_DoubleStage\_280kV = 7**

DC voltage with two AC transformers set to 100kV AC each, resulting in 280kV DC in total (or a single stage transformer with Greinacher voltage doubling rectifier)

**DC\_SingleStage\_140kV = 6**

DC voltage with one AC transformer set to 100kV AC, resulting in 140kV DC

**External = 0**

External power supply fed through blue CEE32 input using isolation transformer and safety switches of the Supercube, or using an external safety switch attached to the Supercube Type B.

**Internal = 5**

Internal usage of the frequency converter, controlling to the primary voltage output of the supercube itself (no measurement transformer used)

**class Safety(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**

Bases: `hvl_ccb.utils.enum.ValueEnum`

NodeID strings for the basic safety circuit status and green/red switches “ready” and “operate”.

**horn = 'hvl-ipc.WINAC.SYSTEM\_INTERN.hornen'**

Writeable boolean to manually turn on or off the horn

**status\_error = 'hvl-ipc.WINAC.SYSTEMSTATE.ERROR'**

**status\_green = 'hvl-ipc.WINAC.SYSTEMSTATE.GREEN'**

**status\_ready\_for\_red = 'hvl-ipc.WINAC.SYSTEMSTATE.ReadyForRed'**

Status is a read-only integer containing the state number of the supercube-internal state machine. The values correspond to numbers in `SafetyStatus`.

**status\_red = 'hvl-ipc.WINAC.SYSTEMSTATE.RED'**

**switchto\_green = 'hvl-ipc.WINAC.SYSTEMSTATE.GREEN\_REQUEST'**

**switchto\_operate = 'hvl-ipc.WINAC.SYSTEMSTATE.switchon'**

Writable boolean for switching to Red Operate (locket, HV on) state.

**switchto\_ready = 'hvl-ipc.WINAC.SYSTEMSTATE.RED\_REQUEST'**

Writable boolean for switching to Red Ready (locked, HV off) state.

**class SafetyStatus(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**

Bases: `aenum.IntEnum`

Safety status values that are possible states returned from `hvl_ccb.dev.supercube.base.Supercube.get_status()`. These values correspond to the states of the Supercube’s safety circuit statemachine.

**Error = 6**

System is in error mode.

**GreenNotReady = 1**

System is safe, lamps are green and some safety elements are not in place such that it cannot be switched to red currently.

**GreenReady = 2**

System is safe and all safety elements are in place to be able to switch to *ready*.

**Initializing = 0**

System is initializing or booting.

**QuickStop = 5**

Fast turn off triggered and switched off the system. Reset FSO to go back to a normal state.

**RedOperate = 4**

System is locked in red state and in *operate* mode, i.e. high voltage on.

**RedReady = 3**

System is locked in red state and *ready* to go to *operate* mode.

```
class SupercubeOpcEndpoint(value=<no_arg>, names=None, module=None, type=None, start=1,
                             boundary=None)
```

Bases: *hvl\_ccb.utils.enum.ValueEnum*

OPC Server Endpoint strings for the supercube variants.


A = 'OPC.SimaticNET.S7'

B = 'OPC.SimaticNET.S7'

T13\_SOCKET\_PORTS = (1, 2, 3)

Port numbers of SEV T13 power socket

## ***hvl\_ccb.dev.supercube2015.typ\_a***



Supercube Typ A module.

```
class Supercube2015WithFU(com, dev_config=None)
```

Bases: *hvl\_ccb.dev.supercube2015.base.Supercube2015Base*

Variant A of the Supercube with frequency converter.

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**fso\_reset()** → *None*

Reset the fast switch off circuitry to go back into normal state and allow to re-enable operate mode.

**get\_frequency()** → float

Read the electrical frequency of the current Supercube setup.

**Returns** the frequency in Hz

**get\_fso\_active()** → bool

Get the state of the fast switch off functionality. Returns True if it is enabled, False otherwise.

**Returns** state of the FSO functionality

**get\_max\_voltage()** → float

Reads the maximum voltage of the setup and returns in V.

**Returns** the maximum voltage of the setup in V.

**get\_power\_setup()** → *hvl\_ccb.dev.supercube2015.constants.PowerSetup*

Return the power setup selected in the Supercube's settings.

**Returns** the power setup

**get\_primary\_current()** → float

Read the current primary current at the output of the frequency converter ( before transformer).

**Returns** primary current in A

**get\_primary\_voltage()** → float

Read the current primary voltage at the output of the frequency converter ( before transformer).

**Returns** primary voltage in V

**get\_target\_voltage()** → float

Gets the current setpoint of the output voltage value in V. This is not a measured value but is the corresponding function to *set\_target\_voltage()*.

**Returns** the setpoint voltage in V.

**set\_slope(*slope: float*)** → *None*

Sets the dV/dt slope of the Supercube frequency converter to a new value in V/s.

**Parameters** **slope** – voltage slope in V/s (0..15'000)

**set\_target\_voltage(*volt\_v: float*)** → *None*

Set the output voltage to a defined value in V.

**Parameters** **volt\_v** – the desired voltage in V

**class SupercubeA0pcUaCommunication(*config*)**

Bases: *hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaCommunication*

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class SupercubeA0pcUaConfiguration(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address],  
                                    endpoint_name: str = 'OPC.SimaticNET.S7', port: int = 4845,  
                                    sub_handler: hvl_ccb.comm.opc.OpcUaSubHandler =  
                                    <hvl_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler  
                                    object at 0x7f29600e3710>, update_parameter:  
                                    asynccua.ua.uaproto_col_auto.CreateSubscriptionParameters =  
                                    CreateSubscriptionParameters(RequestedPublishingInterval=1000,  
                                    RequestedLifetimeCount=300, RequestedMaxKeepAliveCount=22,  
                                    MaxNotificationsPerPublish=10000, PublishingEnabled=True,  
                                    Priority=0), wait_timeout_retry_sec: Union[int, float] = 1,  
                                    max_timeout_retry_nr: int = 5)
```

Bases: `hvl_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig`

**endpoint\_name:** `str = 'OPC.SimaticNET.S7'`

Endpoint of the OPC server, this is a path like ‘OPCUA/SimulationServer’

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.


## Module contents

Supercube package with implementation for the old system version from 2015 based on Siemens WinAC soft-PLC on an industrial 32bit Windows computer.

## `hvl_ccb.dev.technix`

### Submodules

#### `hvl_ccb.dev.technix.base`



Communication and auxiliary classes for Technix

#### `exception TechnixError`

Bases: `hvl_ccb.dev.base.DeviceException`

Technix related errors.

#### `exception TechnixFaultError`

Bases: `hvl_ccb.dev.technix.base.TechinxError`

Raised when the fault flag was detected while the interlock is closed

#### `class TechnixSerialCommunication(configuration)`

Bases: `hvl_ccb.dev.technix.base._TechnixCommunication, hvl_ccb.comm.serial.SerialCommunication`

Serial communication for Technix

##### `static config_cls()`

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class TechnixSerialCommunicationConfig(terminator: bytes = b'\r', encoding: str = 'utf-8',
                                         encoding_error_handling: str = 'strict',
                                         wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                         default_n_attempts_read_text_nonempty: int = 10, port:
                                         Optional[str] = None, baudrate: int = 9600, parity: Union[str,
                                         hvl_ccb.comm.serial.SerialCommunicationParity] =
                                         SerialCommunicationParity.NONE, stopbits: Union[int, float,
                                         hvl_ccb.comm.serial.SerialCommunicationStopbits] =
                                         SerialCommunicationStopbits.ONE, bytesize: Union[int,
                                         hvl_ccb.comm.serial.SerialCommunicationBytesize] =
                                         SerialCommunicationBytesize.EIGHTBITS, timeout: Union[int,
                                         float] = 2)
```

Bases: `hvl_ccb.dev.technix.base._TechnixCommunicationConfig`, `hvl_ccb.comm.serial.SerialCommunicationConfig`

Configuration for the serial communication for Technix

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**class TechnixTelnetCommunication**(*configuration*)

Bases: `hvl_ccb.comm.telnet.TelnetCommunication`, `hvl_ccb.dev.technix.base._TechnixCommunication`

Telnet communication for Technix

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class TechnixTelnetCommunicationConfig(terminator: bytes = b'\r', encoding: str = 'utf-8',
                                         encoding_error_handling: str = 'strict',
                                         wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                         default_n_attempts_read_text_nonempty: int = 10, host:
                                         Optional[Union[str, ipaddress.IPv4Address,
                                         ipaddress.IPv6Address]] = None, port: int = 4660, timeout:
                                         Union[int, float] = 0.2)
Bases:      hvl_ccb.dev.technix.base._TechnixCommunicationConfig,  hvl_ccb.comm.telnet.
                                         TelnetCommunicationConfig
```

Configuration for the telnet communication for Technix

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.


**port: int = 4660**

Port at which Technix is listening

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**hvl\_ccb.dev.technix.device**

its corresponding configuration class

```

class Technix(com, dev_config)
  Bases: hvl_ccb.dev.base.SingleCommDevice
  Device class to control capacitor chargers from Technix

  static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

  property current: Optional[Union[int, float]]
    Actual current of the output in A

  default_com_cls() → Union[Type[hvl_ccb.dev.technix.base.TechnixSerialCommunication],
    Type[hvl_ccb.dev.technix.base.TechnixTelnetCommunication]]
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

  property inhibit: Optional[bool]
    Is the output of the voltage inhibited? The output stage can still be active.

  property is_started: bool
    Is the device started?

  property max_current: Union[int, float]
    Maximal output current of the hardware in A
  
```

The device class *Technix* and

```

property max_voltage: Union[int, float]
    Maximal output voltage of the hardware in V

property open_interlock: Optional[bool]
    Is the interlock open? (in safe mode)

property output: Optional[bool]
    State of the high voltage output

query_status(*, _retry: bool = False)
    Query the status of the device.

Returns This function returns nothing

property remote: Optional[bool]
    Is the device in remote control mode?

start()
    Start the device and set it into the remote controllable mode. The high voltage is turn off, and the status poller is started.

property status: Optional[hvl_ccb.dev.technix.base._Status]
    The status of the device with the different states as sub-fields

stop()
    Stop the device. The status poller is stopped and the high voltage output is turn off.

property voltage: Optional[Union[int, float]]
    Actual voltage at the output in V

property voltage_regulation: Optional[bool]
    Status if the output is in voltage regulation mode (or current regulation)

class TechnixConfig(communication_channel:
    Union[Type[hvl_ccb.dev.technix.base.TechnixSerialCommunication],
    Type[hvl_ccb.dev.technix.base.TechnixTelnetCommunication]], max_voltage: Union[int,
    float], max_current: Union[int, float], polling_interval_sec: Union[int, float] = 4,
    post_stop_pause_sec: Union[int, float] = 1, register_pulse_time: Union[int, float] = 0.1,
    read_output_while_polling: bool = False)

Bases: object

clean_values()
    Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

communication_channel:
Union[Type[hvl_ccb.dev.technix.base.TechnixSerialCommunication],
Type[hvl_ccb.dev.technix.base.TechnixTelnetCommunication]]
    communication channel between computer and Technix

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.

Parameters

- fieldname – name of the field
- value – value to assign

is_configdataclass = True

```

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**max\_current: Union[int, float]**

Maximal Output current

**max\_voltage: Union[int, float]**

Maximal Output voltage

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**polling\_interval\_sec: Union[int, float] = 4**

Polling interval in s to maintain to watchdog of the device

**post\_stop\_pause\_sec: Union[int, float] = 1**

Time to wait after stopping the device

**read\_output\_while\_polling: bool = False**

Read output voltage and current within the polling event

**register\_pulse\_time: Union[int, float] = 0.1**

Time for pulsing a register

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.


## Module contents

Device classes for “RS 232” and “Ethernet” interfaces, which are used to control power supplies from Technix. Manufacturer homepage: <https://www.technix-hv.com>

The regulated power supplies series and capacitor chargers series from Technix are series of low and high voltage direct current power supplies as well as capacitor chargers. The class *Technix* is tested with a CCR10KV-7,5KJ via an ethernet connection as well as a CCR15-P-2500-OP via a serial connection. Check the code carefully before using it with other devices or device series

**This Python package may support the following interfaces from Technix:**

- Remote Interface RS232
- Ethernet Remote Interface
- Optic Fiber Remote Interface

**hvl\_ccb.dev.tiepie****Submodules****hvl\_ccb.dev.tiepie.base**

```
class TiePieDeviceConfig(serial_number: int, require_block_measurement_support: bool = True,
                        n_max_try_get_device: int = 10, wait_sec_retry_get_device: Union[int, float] =
                        1.0, is_data_ready_polling_interval_sec: Union[int, float] = 0.01)
```

Bases: `object`

Configuration dataclass for TiePie

`clean_values()`

`force_value(fieldname, value)`

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

#### Parameters

- `fieldname` – name of the field
- `value` – value to assign

`is_configdataclass = True`

`is_data_ready_polling_interval_sec: Union[int, float] = 0.01`

`classmethod keys() → Sequence[str]`

Returns a list of all configdataclass fields key-names.

`Returns` a list of strings containing all keys.

```
n_max_try_get_device: int = 10
classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.

require_block_measurement_support: bool = True
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

serial_number: int
wait_sec_retry_get_device: Union[int, float] = 1.0
class TiePieDeviceType(value=<no_arg>, names=None, module=None, type=None, start=1,
                       boundary=None)
Bases: hvl_ccb.utils.enum.NameEnum
TiePie device type.

GENERATOR = 2
I2C = 4
OSCILLOSCOPE = 1
exception TiePieError
Bases: hvl_ccb.dev.base.DeviceException
Error of the class TiePie

get_device_by_serial_number(serial_number: int, device_type: Union[str, Tuple[int,
                           hvl_ccb.dev.tiepie.base._LtpDeviceReturnType]], n_max_try_get_device: int =
                           10, wait_sec_retry_get_device: float = 1.0) →
                           hvl_ccb.dev.tiepie.base._LtpDeviceReturnType
Open and return handle of TiePie device with a given serial number
```

#### Parameters

- **serial\_number** – int serial number of the device
- **device\_type** – a *TiePieDeviceType* instance containing device identifier (int number) and its corresponding class, both from *libtiepie*, or a string name of such instance
- **n\_max\_try\_get\_device** – maximal number of device list updates (int number)
- **wait\_sec\_retry\_get\_device** – waiting time in seconds between retries (int number)

**Returns** Instance of a *libtiepie* device class according to the specified *device\_type*

#### Raises

- **TiePieError** – when there is no device with given serial number
- **ValueError** – when *device\_type* is not an instance of *TiePieDeviceType*

```
wrap_libtiepie_exception(func: Callable) → Callable
```


Decorator wrapper for *libtiepie* methods that use *libtiepie.library.check\_last\_status\_raise\_on\_error()* calls.

**Parameters** **func** – Function or method to be wrapped

**Raises** `TiePieError` – instead of `LibTiePieException` or one of its subtypes.

**Returns** whatever `func` returns

## hvl\_ccb.dev.tiepie.channel



### class SafeGround

Bases: `object`

Class that dynamically adds the `safe_ground_enabled` attribute getter/setter if the bound oscilloscope has the safe ground option.

```

class TiePieOscilloscopeChannelConfig(ch_number: int, channel:
    libtiepie.oscilloscopechannel.OscilloscopeChannel)
Bases: hvl_ccb.dev.tiepie.utils.PublicPropertiesReprMixin

Oscilloscope's channel configuration, with cleaning of values in properties setters as well as setting and reading them on and from the device's channel.

static clean_coupling(coupling: Union[str,
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeChannelCoupling]) →
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeChannelCoupling

static clean_enabled(enabled: bool) → bool

clean_input_range(input_range: Union[float, hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeRange]) →
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeRange

clean_probe_offset(probe_offset: float) → float

static clean_trigger_enabled(trigger_enabled)

clean_trigger_hysteresis(trigger_hysteresis: float) → float

static clean_trigger_kind(trigger_kind: Union[str,
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeTriggerKind]) →
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeTriggerKind

clean_trigger_level(trigger_level: Union[int, float]) → float

static clean_trigger_level_mode(level_mode: Union[str,
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeTriggerLevelMode]) →
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeTriggerLevelMode

property coupling: hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeChannelCoupling
property enabled: bool
  
```

```
property has_safe_ground: bool
    Check whether bound oscilloscope device has “safe ground” option

        Returns bool: 1=safe ground available

property input_range: hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeRange
property probe_offset: float
property trigger_enabled: bool
property trigger_hysteresis: float
property trigger_kind: hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeTriggerKind
property trigger_level: float
property trigger_level_mode:
    hvl_ccb.dev.tiepie.channel.TiePieOscilloscopeTriggerLevelMode

class TiePieOscilloscopeChannelConfigLimits(osc_channel:
                                             libtiepie.oscilloscopechannel.OscilloscopeChannel)
Bases: object
Default limits for oscilloscope channel parameters.

class TiePieOscilloscopeChannelCoupling(value=<no_arg>, names=None, module=None, type=None,
                                         start=1, boundary=None)
Bases: hvl_ccb.utils.enum.NameEnum
An enumeration.

ACA = 8
ACV = 2
DCA = 4
DCV = 1

class TiePieOscilloscopeRange(value=<no_arg>, names=None, module=None, type=None, start=1,
                               boundary=None)
Bases: hvl_ccb.utils.enum.NameEnum
An enumeration.

EIGHTY_VOLT = 80
EIGHT_HUNDRED_MILLI_VOLT = 0.8
EIGHT_VOLT = 8
FORTY_VOLT = 40
FOUR_HUNDRED_MILLI_VOLT = 0.4
FOUR_VOLT = 4
TWENTY_VOLT = 20
TWO_HUNDRED_MILLI_VOLT = 0.2
TWO_VOLT = 2

static suitable_range(value)
```

```
class TiePieOscilloscopeTriggerKind(value=<no_arg>, names=None, module=None, type=None, start=1,
                                     boundary=None)
```

Bases: `hvl_ccb.utils.enum.NameEnum`

An enumeration.

**ANY** = 16

**FALLING** = 2

**RISING** = 1

**RISING\_OR\_FALLING** = 16

```
class TiePieOscilloscopeTriggerLevelMode(value=<no_arg>, names=None, module=None, type=None,
                                         start=1, boundary=None)
```

Bases: `hvl_ccb.utils.enum.NameEnum`


An enumeration.

**ABSOLUTE** = 2

**RELATIVE** = 1

**UNKNOWN** = 0

## [hvl\\_ccb.dev.tiepie.device](#)



TiePie devices.

```
class TiePieHS5(com, dev_config)
Bases:      hvl_ccb.dev.tiepie.i2c.TiePieI2CHostMixin,    hvl_ccb.dev.tiepie.generator.
TiePieGeneratorMixin, hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscope

TiePie HS5 device.


class TiePieHS6(com, dev_config)
Bases: hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscope

TiePie HS6 DIFF device.

class TiePieWS5(com, dev_config)
Bases:      hvl_ccb.dev.tiepie.i2c.TiePieI2CHostMixin,    hvl_ccb.dev.tiepie.generator.
TiePieGeneratorMixin, hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscope

TiePie WS5 device.
```

## hvl\_ccb.dev.tiepie.generator



```
class TiePieGeneratorConfig(dev_gen: libtiepie.generator.Generator)
Bases: hvl_ccb.dev.tiepie.utils.PublicPropertiesReprMixin

Generator's configuration with cleaning of values in properties setters.

property amplitude: float
clean_amplitude(amplitude: float) → float
static clean_enabled(enabled: bool) → bool
clean_frequency(frequency: float) → float
clean_offset(offset: float) → float
static clean_signal_type(signal_type: Union[int,
                                              hvl_ccb.dev.tiepie.generator.TiePieGeneratorSignalType]) →
                                              hvl_ccb.dev.tiepie.generator.TiePieGeneratorSignalType
clean_waveform(waveform: numpy.ndarray[Any, numpy.dtype[numpy.typing._generic_alias.ScalarType]]) →
                                              numpy.ndarray[Any, numpy.dtype[numpy.typing._generic_alias.ScalarType]]
property enabled: bool
```

```
property frequency: float
property offset: float
property signal_type: hvl_ccb.dev.tiepie.generator.TiePieGeneratorSignalType
property waveform: Optional[numumpy.ndarray[Any,
numpy.dtype[numpy.typing._generic_alias.ScalarType]]]

class TiePieGeneratorConfigLimits(dev_gen: libtiepie.generator.Generator)
Bases: object

Default limits for generator parameters.

class TiePieGeneratorMixin(com, dev_config)
Bases: object

TiePie Generator sub-device.

A wrapper for the libtiepie.generator.Generator class. To be mixed in with TiePieOscilloscope base class.

config_gen: Optional[hvl_ccb.dev.tiepie.generator.TiePieGeneratorConfig]
    Generator's dynamical configuration.

generator_start()
    Start signal generation.

generator_stop()
    Stop signal generation.

start() → None
    Start the Generator.


stop() → None
    Stop the generator.

class TiePieGeneratorSignalType(value=<no_arg>, names=None, module=None, type=None, start=1,
                                boundary=None)
Bases: hvl_ccb.utils.enum.NameEnum

An enumeration.

ARBITRARY = 32
DC = 8
NOISE = 16
PULSE = 64
SINE = 1
SQUARE = 4
TRIANGLE = 2
UNKNOWN = 0
```

## hvl\_ccb.dev.tiepie.i2c



```
class TiePieI2CHostConfig(dev_i2c: libtiepie.i2chost.I2CHost)
    Bases: hvl_ccb.dev.tiepie.utils.PublicPropertiesReprMixin
    I2C Host's configuration with cleaning of values in properties setters.

class TiePieI2CHostConfigLimits(dev_i2c: libtiepie.i2chost.I2CHost)
    Bases: object
    Default limits for I2C host parameters.


class TiePieI2CHostMixin(com, dev_config)
    Bases: object
    TiePie I2CHost sub-device.
    A wrapper for the libtiepie.i2chost.I2CHost class. To be mixed in with TiePieOscilloscope base class.

config_i2c: Optional[hvl_ccb.dev.tiepie.i2c.TiePieI2CHostConfig]
    I2C host's dynamical configuration.

start() → None
    Start the I2C Host.

stop() → None
    Stop the I2C host.
```

## hvl\_ccb.dev.tiepie.oscilloscope



**class TiePieOscilloscope(*com, dev\_config*)**  
Bases: *hvl\_ccb.dev.base.SingleCommDevice*

TiePie oscilloscope.

A wrapper for TiePie oscilloscopes, based on the class *libtiepie.oscilloscope.Oscilloscope* with simplifications for starting of the device (using serial number) and managing mutable configuration of both the device and its channels, including extra validation and typing hints support for configurations.

Note that, in contrast to *libtiepie* library, since all physical TiePie devices include an oscilloscope, this is the base class for all physical TiePie devices. The additional TiePie sub-devices: “Generator” and “I2CHost”, are mixed-in to this base class in subclasses.

The channels use *1..N* numbering (not *0..N-1*), as in, e.g., the Multi Channel software.

**property channels\_enabled: Generator[int, None, None]**

Yield numbers of enabled channels.

**Returns** Numbers of enabled channels

**collect\_measurement\_data(*timeout: Optional[Union[int, float]] = 0*) → Optional[numpy.ndarray[Any, numpy.dtype[numpy.typing.\_generic\_alias.ScalarType]]]**

Try to collect the data from TiePie; return *None* if data is not ready.

**Parameters** **timeout** – The timeout to wait until data is available. This option makes this function blocking the code. *timeout = None* blocks the code infinitely till data will be available. Per default, the *timeout* is set to 0: The function will not block.

**Returns** Measurement data of only enabled channels and time vector in a 2D-*numpy.ndarray* with float sample data; or *None* if there is no data available.

**static config\_cls() → Type[hvl\_ccb.dev.tiepie.base.TiePieDeviceConfig]**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**config\_osc: Optional[hvl\_ccb.dev.tiepie.oscilloscope.TiePieOscilloscopeConfig]**

Oscilloscope’s dynamical configuration.

**config\_osc\_channel\_dict: Dict[int, hvl\_ccb.dev.tiepie.channel.TiePieOscilloscopeChannelConfig]**

Channel configuration. A *dict* mapping actual channel number, numbered *1..N*, to channel configuration. The channel info is dynamically read from the device only on the first *start()*; beforehand the *dict* is empty.

**static default\_com\_cls()** → Type[*hvl\_ccb.comm.base.NullCommunicationProtocol*]

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**force\_trigger()** → *None*

Forces the TiePie to trigger with a software sided trigger event.

**Return None**

**Raises TiePieError** – when device is not started or status of underlying device gives an error

**is\_measurement\_data\_ready()** → bool

Reports if TiePie has data which is ready to collect

**Returns** if the data is ready to collect.

**Raises TiePieError** – when device is not started or status of underlying device gives an error

**is\_triggered()** → bool

Reports if TiePie has triggered. Maybe data is not yet available. One can check with the function *is\_measurement\_data\_ready()*.

**Returns** if a trigger event occurred

**static list\_devices()** → libtiepie.devicelist.DeviceList

List available TiePie devices.

**Returns** libtiepie up to date list of devices

**property n\_channels**

Number of channels in the oscilloscope.

**Returns** Number of channels.

**start()** → *None*

Start the oscilloscope.

**start\_measurement()** → *None*

Start a measurement using set configuration.

**Raises TiePieError** – when device is not started or status of underlying device gives an error

**stop()** → *None*

Stop the oscilloscope.

**class TiePieOscilloscopeAutoResolutionModes**(*value=<no\_arg>*, *names=None*, *module=None*,  
*type=None*, *start=1*, *boundary=None*)

Bases: *hvl\_ccb.utils.enum.NameEnum*

An enumeration.

**ALL = 4**

**DISABLED = 1**

**NATIVEONLY = 2**

**UNKNOWN = 0**

**class TiePieOscilloscopeConfig**(*dev\_osc: libtiepie.oscilloscope.Oscilloscope*)

Bases: *hvl\_ccb.dev.tiepie.utils.PublicPropertiesReprMixin*

Oscilloscope's configuration with cleaning of values in properties setters.

**property auto\_resolution\_mode:**

*hvl\_ccb.dev.tiepie.oscilloscope.TiePieOscilloscopeAutoResolutionModes*

```

static clean_auto_resolution_mode(auto_resolution_mode: Union[int,
    hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscopeAutoResolutionModes]) →
    hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscopeAutoResolutionModes

clean_pre_sample_ratio(pre_sample_ratio: float) → float
clean_record_length(record_length: Union[int, float]) → int
static clean_resolution(resolution: Union[int,
    hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscopeResolution]) →
    hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscopeResolution

clean_sample_frequency(sample_frequency: float) → float
clean_trigger_timeout(trigger_timeout: Optional[Union[int, float]]) → float
property pre_sample_ratio: float
property record_length: int
property resolution: hvl_ccb.dev.tiepie.oscilloscope.TiePieOscilloscopeResolution
property sample_frequency: float
property trigger_timeout: Optional[float]

class TiePieOscilloscopeConfigLimits(dev_osc: libtiepie.oscilloscope.Oscilloscope)
    Bases: object
    Default limits for oscilloscope parameters.

class TiePieOscilloscopeResolution(value=<no_arg>, names=None, module=None, type=None, start=1,
    boundary=None)
    Bases: aenum.IntEnum
    An enumeration.

EIGHT_BIT = 8
FOURTEEN_BIT = 14
SIXTEEN_BIT = 16
TWELVE_BIT = 12

```

## **hvl\_ccb.dev.tiepie.utils**

### PublicPropertiesReprMixin

```

class PublicPropertiesReprMixin
    Bases: object

```

General purpose utility mixin that overwrites object representation to a one analogous to *dataclass* instances, but using public properties and their values instead of *fields*.

## Module contents

This module is a wrapper around LibTiePie SDK devices; see <https://www.tiepie.com/en/libtiepie-sdk>.

The device classes adds simplifications for starting of the device (using serial number) and managing mutable configuration of both the device and oscilloscope's channels. This includes extra validation and typing hints support.

## Extra installation

LibTiePie SDK library is available only on Windows and on Linux.

To use this LibTiePie SDK devices wrapper:

1. install the hvl\_ccb package with a tiepie extra feature:

```
$ pip install "hvl_ccb[tiepie]"
```

this will install the Python bindings for the library.

2. install the library

- on Linux: follow instructions in <https://www.tiepie.com/en/libtiepie-sdk/linux>;
- on Windows: the additional DLL is included in Python bindings package.

## Troubleshooting

On a Windows system, if you encounter an OSError like this:


```
...
self._handle = _dlopen(self._name, mode)
OSSError: [WinError 126] The specified module could not be found
```

most likely the python-libtiepie package was installed in your site-packages/ directory as a python-libtiepie-\*.egg file via python setup.py install or python setup.py develop command. In such case uninstall the library and re-install it using pip:

```
$ pip uninstall python-libtiepie
$ pip install python-libtiepie
```

This should create libtiepie/ folder. Alternatively, manually move the folder libtiepie/ from inside of the .egg archive file to the containing it site-packages/ directory (PyCharm's Project tool window supports reading and extracting from .egg archives).

## Submodules

**hvl\_ccb.dev.base**

Module with base classes for devices.

```

class Device(dev_config=None)
    Bases: hvl_ccb.configuration.ConfigurationMixin, abc.ABC

    Base class for devices. Implement this class for a concrete device, such as measurement equipment or voltage sources.

    Specifies the methods to implement for a device.

static config_cls()
    Return the default configdataclass class

    Returns a reference to the default configdataclass class

abstract start() → None
    Start or restart this Device. To be implemented in the subclass.

abstract stop() → None
    Stop this Device. To be implemented in the subclass.

exception DeviceException
    Bases: hvl_ccb.exception.CCBException

exception DeviceExistingException
    Bases: hvl_ccb.dev.base.DeviceException

    Exception to indicate that a device with that name already exists.

exception DeviceFailuresException(failures: Dict[str, Exception], *args)
    Bases: hvl_ccb.dev.base.DeviceException

    Exception to indicate that one or several devices failed.

    failures: Dict[str, Exception]
        A dictionary of named devices failures (exceptions).

class DeviceSequenceMixin(devices: Dict[str, hvl_ccb.dev.base.Device])
    Bases: abc.ABC

    Mixin that can be used on a device or other classes to provide facilities for handling multiple devices in a sequence.

add_device(name: str, device: hvl_ccb.dev.base.Device) → None
    Add a new device to the device sequence.

```

### Parameters

- **name** – is the name of the device.
- **device** – is the instantiated Device object.

**Raises** *DeviceExistingException* –

**devices\_failed\_start:** Dict[str, *hvl\_ccb.dev.base.Device*]

Dictionary of named device instances from the sequence for which the most recent *start()* attempt failed.

Empty if *stop()* was called last; cf. *devices\_failed\_stop*.

**devices\_failed\_stop:** Dict[str, *hvl\_ccb.dev.base.Device*]

Dictionary of named device instances from the sequence for which the most recent *stop()* attempt failed.

Empty if *start()* was called last; cf. *devices\_failed\_start*.

**get\_device(name: str) → hvl\_ccb.dev.base.Device**

Get a device by name.

**Parameters** **name** – is the name of the device.

**Returns** the device object from this sequence.

**get\_devices() → List[Tuple[str, hvl\_ccb.dev.base.Device]]**

Get list of name, device pairs according to current sequence.

**Returns** A list of tuples with name and device each.

**remove\_device(name: str) → hvl\_ccb.dev.base.Device**

Remove a device from this sequence and return the device object.

**Parameters** **name** – is the name of the device.

**Returns** device object or *None* if such device was not in the sequence.

**Raises** **ValueError** – when device with given name was not found

**start() → None**

Start all devices in this sequence in their added order.

**Raises** *DeviceFailuresException* – if one or several devices failed to start

**stop() → None**

Stop all devices in this sequence in their reverse order.

**Raises** *DeviceFailuresException* – if one or several devices failed to stop

**class EmptyConfig**

Bases: object

Empty configuration dataclass that is the default configuration for a Device.

**clean\_values()**

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

### Parameters

- **fieldname** – name of the field

- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**class SingleCommDevice(com, dev\_config=None)**

Bases: *hvl\_ccb.dev.base.Device*, *abc.ABC*

Base class for devices with a single communication protocol.

**property com**

Get the communication protocol of this device.

**Returns** an instance of *CommunicationProtocol* subtype

**abstract static default\_com\_cls() → Type[hvl\_ccb.comm.base.CommunicationProtocol]**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device


**start() → None**

Open the associated communication protocol.

**stop() → None**

Close the associated communication protocol.

## **hvl\_ccb.dev.crylas**



Device classes for a CryLas pulsed laser controller and a CryLas laser attenuator, using serial communication.

There are three modes of operation for the laser 1. Laser-internal hardware trigger (default): fixed to 20 Hz and max energy per pulse. 2. Laser-internal software trigger (for diagnosis only). 3. External trigger: required for arbitrary pulse energy or repetition rate. Switch to “external” on the front panel of laser controller for using option 3.

After switching on the laser with `laser_on()`, the system must stabilize for some minutes. Do not apply abrupt changes of pulse energy or repetition rate.

Manufacturer homepage: [https://www.crylas.de/products/pulsed\\_laser.html](https://www.crylas.de/products/pulsed_laser.html)

```
class CryLasAttenuator(com, dev_config=None)
    Bases: hvl_ccb.dev.base.SingleCommDevice

    Device class for the CryLas laser attenuator.

    property attenuation: Union[int, float]

    static config_cls()
        Return the default configdataclass class.

        Returns a reference to the default configdataclass class

    static default_com_cls()
        Get the class for the default communication protocol used with this device.

        Returns the type of the standard communication protocol for this device

    set_attenuation(percent: Union[int, float]) → None
        Set the percentage of attenuated light (inverse of set_transmission). :param percent: percentage of attenuation, number between 0 and 100 :raises ValueError: if param percent not between 0 and 100 :raises SerialCommunicationIOError: when communication port is not opened :raises CryLasAttenuatorError: if the device does not confirm success

    set_init_attenuation()
        Sets the attenuation to its configured initial/default value

        Raises SerialCommunicationIOError – when communication port is not opened

    set_transmission(percent: Union[int, float]) → None
        Set the percentage of transmitted light (inverse of set_attenuation). :param percent: percentage of transmitted light :raises ValueError: if param percent not between 0 and 100 :raises SerialCommunicationIOError: when communication port is not opened :raises CryLasAttenuatorError: if the device does not confirm success

    start() → None
        Open the com, apply the config value ‘init_attenuation’

        Raises SerialCommunicationIOError – when communication port cannot be opened

    property transmission: Union[int, float]

class CryLasAttenuatorConfig(init_attenuation: Union[int, float] = 0, response_sleep_time: Union[int, float] = 1)
    Bases: object

    Device configuration dataclass for CryLas attenuator.

    clean_values()

    force_value(fieldname, value)
        Forces a value to a dataclass field despite the class being frozen.

        NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.

    Parameters
        • fieldname – name of the field
        • value – value to assign

    init_attenuation: Union[int, float] = 0
    is_configdataclass = True
```

```
classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

response_sleep_time: Union[int, float] = 1

exception CryLasAttenuatorError
    Bases: hvl\_ccb.dev.base.DeviceException
    General error with the CryLas Attenuator.

class CryLasAttenuatorSerialCommunication(configuration)
    Bases: hvl\_ccb.comm.serial.SerialCommunication
    Specific communication protocol implementation for the CryLas attenuator. Already predefines device-specific protocol parameters in config.

    static config_cls()
        Return the default configdataclass class.

        Returns a reference to the default configdataclass class

class CryLasAttenuatorSerialCommunicationConfig(terminator: bytes = b'', encoding: str = 'utf-8', encoding_error_handling: str = 'strict', wait_sec_read_text_nonempty: Union[int, float] = 0.5, default_n_attempts_read_text_nonempty: int = 10, port: Union[str, NoneType] = None, baudrate: int = 9600, parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity] = <SerialCommunicationParity.NONE: 'N'>, stopbits: Union[int, hvl_ccb.comm.serial.SerialCommunicationStopbits] = <SerialCommunicationStopbits.ONE: 1>, bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize] = <SerialCommunicationBytesize.EIGHTBITS: 8>, timeout: Union[int, float] = 3)
        Bases: hvl\_ccb.comm.serial.SerialCommunicationConfig

baudrate: int = 9600
    Baudrate for CryLas attenuator is 9600 baud

bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize] = 8
    One byte is eight bits long

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.
```

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity:** Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = 'N'

CryLas attenuator does not use parity

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits:** Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = 1

CryLas attenuator uses one stop bit

**terminator:** bytes = b''

No terminator

**timeout:** Union[int, float] = 3

use 3 seconds timeout as default

**class CryLasLaser(com, dev\_config=None)**

Bases: hvl\_ccb.dev.base.SingleCommDevice

CryLas laser controller device class.

**class AnswersShutter(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**

Bases: aenum.Enum

Standard answers of the CryLas laser controller to ‘Shutter’ command passed via *com*.

**CLOSED** = ‘Shutter inaktiv’

**OPENED** = ‘Shutter aktiv’

**class AnswersStatus(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**

Bases: aenum.Enum

Standard answers of the CryLas laser controller to ‘STATUS’ command passed via *com*.

**ACTIVE** = ‘STATUS: Laser active’

**HEAD** = ‘STATUS: Head ok’

**INACTIVE** = ‘STATUS: Laser inactive’

**READY** = ‘STATUS: System ready’

```

TEC1 = 'STATUS: TEC1 Regulation ok'
TEC2 = 'STATUS: TEC2 Regulation ok'

class LaserStatus(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
    Bases: aenum.Enum

    Status of the CryLas laser

    READY_ACTIVE = 2
    READY_INACTIVE = 1
    UNREADY_INACTIVE = 0
    property is_inactive
    property is_ready

class RepetitionRates(value=<no_arg>, names=None, module=None, type=None, start=1,
                      boundary=None)
    Bases: aenum.IntEnum

    Repetition rates for the internal software trigger in Hz

    HARDWARE = 0
    SOFTWARE_INTERNAL_SIXTY = 60
    SOFTWARE_INTERNAL_TEN = 10
    SOFTWARE_INTERNAL_TWENTY = 20

ShutterStatus
    alias of hvl\_ccb.dev.crylas.CryLasLaserShutterStatus

close_shutter() → None
    Close the laser shutter.

    Raises
        • SerialCommunicationIOError – when communication port is not opened
        • CryLasLaserError – if success is not confirmed by the device

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

get_pulse_energy_and_rate() → Tuple[int, int]
    Use the debug mode, return the measured pulse energy and rate.

    Returns (energy in micro joule, rate in Hz)

    Raises
        • SerialCommunicationIOError – when communication port is not opened
        • CryLasLaserError – if the device does not answer the query

laser_off() → None
    Turn the laser off.

```

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**laser\_on()** → *None*

Turn the laser on.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserNotReadyError* – if the laser is not ready to be turned on
- *CryLasLaserError* – if success is not confirmed by the device

**open\_shutter()** → *None*

Open the laser shutter.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**set\_init\_shutter\_status()** → *None*

Open or close the shutter, to match the configured shutter\_status.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**set\_pulse\_energy(energy: int)** → *None*

Sets the energy of pulses (works only with external hardware trigger). Proceed with small energy steps, or the regulation may fail.

**Parameters** **energy** – energy in micro joule

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if the device does not confirm success

**set\_repetition\_rate(rate: Union[int, hvl\_ccb.dev.crylas.CryLasLaser.RepetitionRates])** → *None*

Sets the repetition rate of the internal software trigger.

**Parameters** **rate** – frequency (Hz) as an integer

**Raises**

- *ValueError* – if rate is not an accepted value in RepetitionRates Enum
- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**start()** → *None*

Opens the communication protocol and configures the device.

**Raises** *SerialCommunicationIOError* – when communication port cannot be opened

**stop()** → *None*

Stops the device and closes the communication protocol.

**Raises**

- **SerialCommunicationIOError** – if com port is closed unexpectedly
- **CryLasLaserError** – if laser\_off() or close\_shutter() fail

**property target\_pulse\_energy**

**update\_laser\_status()** → *None*

Update the laser status to *LaserStatus.NOT\_READY* or *LaserStatus.INACTIVE* or *LaserStatus.ACTIVE*.

Note: laser never explicitly says that it is not ready (*LaserStatus.NOT\_READY*) in response to ‘STATUS’ command. It only says that it is ready (heated-up and implicitly inactive/off) or active (on). If it’s not either of these then the answer is *Answers.HEAD*. Moreover, the only time the laser explicitly says that its status is inactive (*Answers.INACTIVE*) is after issuing a ‘LASER OFF’ command.

**Raises** **SerialCommunicationIOError** – when communication port is not opened

**update\_repetition\_rate()** → *None*

Query the laser repetition rate.

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **CryLasLaserError** – if success is not confirmed by the device

**update\_shutter\_status()** → *None*

Update the shutter status (OPENED or CLOSED)

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **CryLasLaserError** – if success is not confirmed by the device

**update\_target\_pulse\_energy()** → *None*

Query the laser pulse energy.

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **CryLasLaserError** – if success is not confirmed by the device

**wait\_until\_ready()** → *None*

Block execution until the laser is ready

**Raises** **CryLasLaserError** – if the polling thread stops before the laser is ready

```
class CryLasLaserConfig(calibration_factor: Union[int, float] = 4.35, polling_period: Union[int, float] = 12,
                        polling_timeout: Union[int, float] = 300, auto_laser_on: bool = True,
                        init_shutter_status: Union[int, hvl_ccb.dev.crylas.CryLasLaserShutterStatus] =
                        CryLasLaserShutterStatus.CLOSED)
```

Bases: object

Device configuration dataclass for the CryLas laser controller.

**ShutterStatus**

alias of *hvl\_ccb.dev.crylas.CryLasLaserShutterStatus*

**auto\_laser\_on:** bool = True

**calibration\_factor:** Union[int, float] = 4.35

**clean\_values()**

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**init\_shutter\_status:** Union[int, *hvl\_ccb.dev.crylas.CryLasLaserShutterStatus*] = 0

**is\_configdataclass** = True

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**polling\_period:** Union[int, float] = 12

**polling\_timeout:** Union[int, float] = 300

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception CryLasLaserError**

Bases: *hvl\_ccb.dev.base.DeviceException*

General error with the CryLas Laser.

**exception CryLasLaserNotReadyError**

Bases: *hvl\_ccb.dev.crylas.CryLasLaserError*

Error when trying to turn on the CryLas Laser before it is ready.

**class CryLasLaserPoller**(*spoll\_handler*: Callable, *check\_handler*: Callable, *check\_laser\_status\_handler*: Callable, *polling\_delay\_sec*: Union[int, float] = 0, *polling\_interval\_sec*: Union[int, float] = 1, *polling\_timeout\_sec*: Optional[Union[int, float]] = None)

Bases: *hvl\_ccb.dev.utils.Poller*

Poller class for polling the laser status until the laser is ready.

**Raises**

- **CryLasLaserError** – if the timeout is reached before the laser is ready
- **SerialCommunicationIOError** – when communication port is closed.

**class CryLasLaserSerialCommunication**(*configuration*)

Bases: *hvl\_ccb.comm.serial.SerialCommunication*

Specific communication protocol implementation for the CryLas laser controller. Already predefines device-specific protocol parameters in config.

**READ\_TEXT\_SKIP\_PREFIXES** = ('>', 'MODE:')

Prefixes of lines that are skipped when read from the serial port.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query**(*cmd*: str, *prefix*: str, *post\_cmd*: Optional[str] = None) → str

Send a command, then read the com until a line starting with prefix, or an empty line, is found. Returns the line in question.

#### Parameters

- **cmd** – query message to send to the device
- **prefix** – start of the line to look for in the device answer
- **post\_cmd** – optional additional command to send after the query

**Returns** line in question as a string

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**query\_all**(*cmd*: str, *prefix*: str)

Send a command, then read the com until a line starting with prefix, or an empty line, is found. Returns a list of successive lines starting with prefix.

#### Parameters

- **cmd** – query message to send to the device
- **prefix** – start of the line to look for in the device answer

**Returns** line in question as a string

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**read()** → str

Read first line of text from the serial port that does not start with any of *self.READ\_TEXT\_SKIP\_PREFIXES*.

**Returns** String read from the serial port; '' if there was nothing to read.

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**class CryLasLaserSerialCommunicationConfig**(*terminator*: bytes = b'\n', *encoding*: str = 'utf-8', *encoding\_error\_handling*: str = 'strict', *wait\_sec\_read\_text\_nonempty*: Union[int, float] = 0.5, *default\_n\_attempts\_read\_text\_nonempty*: int = 10, *port*: Union[str, NoneType] = None, *baudrate*: int = 19200, *parity*: Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = <SerialCommunicationParity.NONE: 'N'>, *stopbits*: Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = <SerialCommunicationStopbits.ONE: 1>, *bytesize*: Union[int, hvl\_ccb.comm.serial.SerialCommunicationBytesize] = <SerialCommunicationBytesize.EIGHTBITS: 8>, *timeout*: Union[int, float] = 10)

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

```
baudrate: int = 19200
    Baudrate for CryLas laser is 19200 baud

bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize] = 8
    One byte is eight bits long

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing
    after value has been forced on fieldname.

Parameters

- fieldname – name of the field
- value – value to assign

classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified
    on instantiation.

    Returns a list of strings containing all optional keys.

parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity] = 'N'
    CryLas laser does not use parity

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on
    instantiation.

    Returns a list of strings containing all required keys.

stopbits: Union[int, hvl_ccb.comm.serial.SerialCommunicationStopbits] = 1
    CryLas laser uses one stop bit

terminator: bytes = b'\n'
    The terminator is LF


timeout: Union[int, float] = 10
    use 10 seconds timeout as default (a long timeout is needed!)

class CryLasLaserShutterStatus(value=<no_arg>, names=None, module=None, type=None, start=1,
                                boundary=None)
    Bases: aenum.Enum

    Status of the CryLas laser shutter

    CLOSED = 0
    OPENED = 1
```

## `hvl_ccb.dev.ea_psi9000`



Device class for controlling a Elektro Automatik PSI 9000 power supply over VISA.

It is necessary that a backend for pyvisa is installed. This can be NI-Visa oder pyvisa-py (up to now, all the testing was done with NI-Visa)

```

class PSI9000(com: Union[hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunication,
hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig, dict], dev_config:
Optional[Union[hvl_ccb.dev.ea_psi9000.PSI9000Config, dict]] = None)
  
```

Bases: *hvl\_ccb.dev.visa.VisaDevice*

Elektro Automatik PSI 9000 power supply.

**MS\_NOMINAL\_CURRENT = 2040**

**MS\_NOMINAL\_VOLTAGE = 80**

**SHUTDOWN\_CURRENT\_LIMIT = 0.1**

**SHUTDOWN\_VOLTAGE\_LIMIT = 0.1**

**check\_master\_slave\_config()** → *None*

Checks if the master / slave configuration and initializes if successful

**Returns** *PSI9000Error* – if master-slave configuration failed

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Return the default communication protocol for this device type, which is VisaCommunication.

**Returns** the VisaCommunication class

**get\_output()** → *bool*

Reads the current state of the DC output of the source. Returns True, if it is enabled, false otherwise.

**Returns** the state of the DC output

**get\_system\_lock()** → bool

Get the current lock state of the system. The lock state is true, if the remote control is active and false, if not.

**Returns** the current lock state of the device

**get\_ui\_lower\_limits()** → Tuple[float, float]

Get the lower voltage and current limits. A lower power limit does not exist.

**Returns** Umin in V, Imin in A

**get\_uip\_upper\_limits()** → Tuple[float, float, float]

Get the upper voltage, current and power limits.

**Returns** Umax in V, Imax in A, Pmax in W

**get\_voltage\_current\_setpoint()** → Tuple[float, float]

Get the voltage and current setpoint of the current source.

**Returns** Uset in V, Iset in A

**measure\_voltage\_current()** → Tuple[float, float]

Measure the DC output voltage and current

**Returns** Umeas in V, Imeas in A

**set\_lower\_limits(voltage\_limit: Optional[float] = None, current\_limit: Optional[float] = None)** → None

Set the lower limits for voltage and current. After writing the values a check is performed if the values are set correctly.

#### Parameters

- **voltage\_limit** – is the lower voltage limit in V
- **current\_limit** – is the lower current limit in A

**Raises** *PSI9000Error* – if the limits are out of range

**set\_output(target\_onstate: bool)** → None

Enables / disables the DC output.

**Parameters** **target\_onstate** – enable or disable the output power

**Raises** *PSI9000Error* – if operation was not successful

**set\_system\_lock(lock: bool)** → None

Lock / unlock the device, after locking the control is limited to this class unlocking only possible when voltage and current are below the defined limits

**Parameters** **lock** – True: locking, False: unlocking

**set\_upper\_limits(voltage\_limit: Optional[float] = None, current\_limit: Optional[float] = None, power\_limit: Optional[float] = None)** → None

Set the upper limits for voltage, current and power. After writing the values a check is performed if the values are set. If a parameter is left blank, the maximum configurable limit is set.

#### Parameters

- **voltage\_limit** – is the voltage limit in V
- **current\_limit** – is the current limit in A
- **power\_limit** – is the power limit in W

**Raises** *PSI9000Error* – if limits are out of range

**set\_voltage\_current**(*volt*: float, *current*: float) → *None*

Set voltage and current setpoints.

After setting voltage and current, a check is performed if writing was successful.

#### Parameters

- **volt** – is the setpoint voltage: 0.81.6 V (1.02 \* 0-80 V) (absolute max, can be smaller if limits are set)
- **current** – is the setpoint current: 0..2080.8 A (1.02 \* 0 - 2040 A) (absolute max, can be smaller if limits are set)

**Raises** *PSI9000Error* – if the desired setpoint is out of limits

**start()** → *None*

Start this device.

**stop()** → *None*

Stop this device. Turns off output and lock, if enabled.

```
class PSI9000Config(spoll_interval: Union[int, float] = 0.5, spoll_start_delay: Union[int, float] = 2,
                     power_limit: Union[int, float] = 43500, voltage_lower_limit: Union[int, float] = 0.0,
                     voltage_upper_limit: Union[int, float] = 10.0, current_lower_limit: Union[int, float] = 0.0,
                     current_upper_limit: Union[int, float] = 2040.0, wait_sec_system_lock: Union[int, float]
                     = 0.5, wait_sec_settings_effect: Union[int, float] = 1, wait_sec_initialisation: Union[int,
                     float] = 2)
```

Bases: *hvl\_ccb.dev.visa.VisaDeviceConfig*

Elektro Automatik PSI 9000 power supply device class. The device is communicating over a VISA TCP socket.

Using this power supply, DC voltage and current can be supplied to a load with up to 2040 A and 80 V (using all four available units in parallel). The maximum power is limited by the grid, being at 43.5 kW available through the CEE63 power socket.

**clean\_values()** → *None*

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**current\_lower\_limit: Union[int, float] = 0.0**

Lower current limit in A, depending on the experimental setup.

**current\_upper\_limit: Union[int, float] = 2040.0**

Upper current limit in A, depending on the experimental setup.

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**power\_limit: Union[int, float] = 43500**

Power limit in W depending on the experimental setup. With 3x63A, this is 43.5kW. Do not change this value, if you do not know what you are doing. There is no lower power limit.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**voltage\_lower\_limit: Union[int, float] = 0.0**

Lower voltage limit in V, depending on the experimental setup.

**voltage\_upper\_limit: Union[int, float] = 10.0**

Upper voltage limit in V, depending on the experimental setup.

**wait\_sec\_initialisation: Union[int, float] = 2**

**wait\_sec\_settings\_effect: Union[int, float] = 1**

**wait\_sec\_system\_lock: Union[int, float] = 0.5**

**exception PSI9000Error**

Bases: *hvl\_ccb.dev.base.DeviceException*

Base error class regarding problems with the PSI 9000 supply.

**class PSI9000VisaCommunication(configuration)**

Bases: *hvl\_ccb.comm.visa.VisaCommunication*

Communication protocol used with the PSI 9000 power supply.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class PSI9000VisaCommunicationConfig(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address], interface\_type: Union[str, hvl\_ccb.comm.visa.VisaCommunicationConfig.InterfaceType] = InterfaceType.TCPIP\_SOCKET, board: int = 0, port: int = 5025, timeout: int = 5000, chunk\_size: int = 204800, open\_timeout: int = 1000, write\_termination: str = '\n', read\_termination: str = '\n', visa\_backend: str = '')**

Bases: *hvl\_ccb.comm.visa.VisaCommunicationConfig*

Visa communication protocol config dataclass with specification for the PSI 9000 power supply.

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**interface\_type:** Union[str, *hvl\_ccb.comm.visa.VisaCommunicationConfig.InterfaceType*]  
**= 1**

Interface type of the VISA connection, being one of InterfaceType.

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.


**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## hvl\_ccb.dev.fug



Device classes for “Probus V - ADDAT30” Interfaces which are used to control power supplies from FuG Elektronik GmbH

This interface is used for many FuG power units. Manufacturer homepage: <https://www.fug-elektronik.de>

The Professional Series of Power Supplies from FuG is a series of low, medium and high voltage direct current power supplies as well as capacitor chargers. The class FuG is tested with a HCK 800-20 000 in Standard Mode. The addressable mode is not implemented. Check the code carefully before using it with other devices. Manufacturer homepage: <https://www.fug-elektronik.de/netzgeraete/professional-series/>

The documentation of the interface from the manufacturer can be found here: [https://www.fug-elektronik.de/wp-content/uploads/download/de/SOFTWARE/Probus\\_V.zip](https://www.fug-elektronik.de/wp-content/uploads/download/de/SOFTWARE/Probus_V.zip)

The provided classes support the basic and some advanced commands. The commands for calibrating the power supplies are not implemented, as they are only for very special porpoises and should not be used by “normal” customers.

**class FuG(com, dev\_config=None)**  
Bases: *hvl\_ccb.dev.fug.FuGProbusV*

FuG power supply device class.

The power supply is controlled over a FuG ADDA Interface with the PROBUS V protocol

**property config\_status:** *hvl\_ccb.dev.fug.FuGProbusVConfigRegisters*  
Returns the registers for the registers with the configuration and status values

**Returns** FuGProbusVConfigRegisters

**property current:** *hvl\_ccb.dev.fug.FuGProbusVSetRegisters*

Returns the registers for the current output

**Returns**

**property current\_monitor:** *hvl\_ccb.dev.fug.FuGProbusVMonitorRegisters*

Returns the registers for the current monitor.

A typically usage will be “self.current\_monitor.value” to measure the output current

**Returns**

**property di:** *hvl\_ccb.dev.fug.FuGProbusVDIRegisters*

Returns the registers for the digital inputs

**Returns** FuGProbusVDIRegisters

**identify\_device()** → *None*

Identify the device nominal voltage and current based on its model number.

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**property max\_current:** *Union[int, float]*

Returns the maximal current which could provided within the test setup

**Returns**

**property max\_current\_hardware:** *Union[int, float]*

Returns the maximal current which could provided with the power supply

**Returns**

**property max\_voltage:** *Union[int, float]*

Returns the maximal voltage which could provided within the test setup

**Returns**

**property max\_voltage\_hardware:** *Union[int, float]*

Returns the maximal voltage which could provided with the power supply

**Returns**

**property on:** *hvl\_ccb.dev.fug.FuGProbusVDORegisters*

Returns the registers for the output switch to turn the output on or off

**Returns** FuGProbusVDORegisters

**property outX0:** *hvl\_ccb.dev.fug.FuGProbusVDORegisters*

Returns the registers for the digital output X0

**Returns** FuGProbusVDORegisters

**property outX1:** *hvl\_ccb.dev.fug.FuGProbusVDORegisters*

Returns the registers for the digital output X1

**Returns** FuGProbusVDORegisters

**property outX2:** *hvl\_ccb.dev.fug.FuGProbusVDORegisters*

Returns the registers for the digital output X2

**Returns** FuGProbusVDORegisters

**property outXCMD:** *hvl\_ccb.dev.fug.FuGProbusVDORegisters*

Returns the registers for the digital outputX-CMD

**Returns** FuGProbusVDORegisters

**start**(*max\_voltage=0, max\_current=0*) → *None*  
Opens the communication protocol and configures the device.

#### Parameters

- **max\_voltage** – Configure here the maximal permissible voltage which is allowed in the given experimental setup
- **max\_current** – Configure here the maximal permissible current which is allowed in the given experimental setup

**property voltage:** *hvl\_ccb.dev.fug.FuGProbusVSetRegisters*  
Returns the registers for the voltage output

#### Returns

**property voltage\_monitor:** *hvl\_ccb.dev.fug.FuGProbusVMonitorRegisters*  
Returns the registers for the voltage monitor.  
A typically usage will be “self.voltage\_monitor.value” to measure the output voltage

#### Returns

**class FuGConfig**(*wait\_sec\_stop\_commands: Union[int, float] = 0.5*)  
Bases: *object*

Device configuration dataclass for FuG power supplies.

**clean\_values()**

**force\_value**(*fieldname, value*)  
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys()** → Sequence[str]  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**wait\_sec\_stop\_commands: Union[int, float] = 0.5**  
Time to wait after subsequent commands during stop (in seconds)

```
class FuGDigitalVal(value)
    Bases: enum.IntEnum
    An enumeration.

    NO = 0
    OFF = 0
    ON = 1
    YES = 1

exception FuGError(*args, **kwargs)
    Bases: hvl_ccb.dev.base.DeviceException
    Error with the FuG voltage source.

    errorcode: str
        Errorcode from the Probus, see documentation of Probus V chapter 5. Errors with three-digit errorcodes
        are thrown by this python module.

class FuGErrorcodes(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
    Bases: hvl_ccb.utils.enum.NameEnum

    The power supply can return an errorcode. These errorcodes are handled by this class. The original errorcodes
    from the source are with one or two digits, see documentation of Probus V chapter 5. All three-digit errorcodes
    are from this python module.

    E0 = ('no error', 'standard response on each command')
    E1 = ('no data available', 'Customer tried to read from GPIB but there were no data
          prepared. (IBIG50 sent command ~T2 to ADDA)')
    E10 = ('unknown SCPI command', 'This SCPI command is not implemented')
    E100 = ('Command is not implemented', 'You tried to execute a command, which is not
           implemented or does not exist')
    E106 = ('The rampstate is a read-only register', 'You tried to write data to the
           register, which can only give you the status of the ramping.')
    E11 = ('not allowed Trigger-on-Talk', 'Not allowed attempt to Trigger-on-Talk (~T1)
           while ADDA was in addressable mode.')
    E115 = ('The given index to select a digital value is out of range', 'Only integer
           values between 0 and 1 are allowed.')
    E12 = ('invalid argument in ~Tn command', 'Only ~T1 and ~T2 is implemented.')
    E125 = ('The given index to select a ramp mode is out of range', 'Only integer
           values between 0 and 4 are allowed.')
    E13 = ('invalid N-value', 'Register > K8 contained an invalid value. Error code is
           output on an attempt to query data with ? or ~T1')
    E135 = ('The given index to select the readback channel is out of range', 'Only
           integer values between 0 and 6 are allowed.')
    E14 = ('register is write only', 'Some registers can only be written to (i.e. > H0)')
    E145 = ('The given value for the AD-conversion is unknown', 'Valid values for the
           ad-conversion are integer values from "0" to "7".')
    E15 = ('string too long', 'i.e. serial number string too long during calibration')
```

```

E155 = ('The given value to select a polarity is out range.', 'The value should be 0 or 1.')
E16 = ('wrong checksum', 'checksum over command string was not correct, refer also to 4.4 of the Probus V documentation')
E165 = ('The given index to select the terminator string is out of range', '')
E2 = ('unknown register type', "No valid register type after '>'")
E206 = ('This status register is read-only', 'You tried to write data to this register, which can only give you the actual status of the corresponding digital output.')
E306 = ('The monitor register is read-only', 'You tried to write data to a monitor, which can only give you measured data.')
E4 = ('invalid argument', 'The argument of the command was rejected .i.e. malformed number')
E5 = ('argument out of range', 'i.e. setvalue higher than type value')
E504 = ('Empty string as response', 'The connection is broken.')
E505 = ('The returned register is not the requested.', 'Maybe the connection is overburden.')
E6 = ('register is read only', 'Some registers can only be read but not written to. (i.e. monitor registers)')
E666 = ('You cannot overwrite the most recent error in the interface of the power supply. But, well: You created an error anyway...', '')
E7 = ('Receive Overflow', 'Command string was longer than 50 characters.')
E8 = ('EEPROM is write protected', 'Write attempt to calibration data while the write protection switch was set to write protected.')
E9 = ('address error', 'A non addressed command was sent to ADDA while it was in addressable mode (and vice versa.)')
raise_()

```

**class FuGMonitorModes(*value*)**

Bases: enum.IntEnum

An enumeration.

**T1MS = 1**  
15 bit + sign, 1 ms integration time

**T200MS = 6**  
typ. 19 bit + sign, 200 ms integration time

**T20MS = 3**  
17 bit + sign, 20 ms integration time

**T256US = 0**  
14 bit + sign, 256 us integration time

**T40MS = 4**  
17 bit + sign, 40 ms integration time

**T4MS = 2**  
15 bit + sign, 4 ms integration time

```
T800MS = 7
typ. 20 bit + sign, 800 ms integration time

T80MS = 5
typ. 18 bit + sign, 80 ms integration time

class FuGPolarities(value)
Bases: enum.IntEnum

An enumeration.

NEGATIVE = 1
POSITIVE = 0

class FuGProbusIV(com, dev_config=None)
Bases: hvl_ccb.dev.base.SingleCommDevice, abc.ABC

FuG Probus IV device class

Sends basic SCPI commands and reads the answer. Only the special commands and PROBUS IV instruction set is implemented.

command(command: hvl_ccb.dev.fug.FuGProbusIVCommands, value=None) → str

Parameters

- command – one of the commands given within FuGProbusIVCommands
- value – an optional value, depending on the command

Returns a String if a query was performed

static config_cls()
Return the default configdataclass class.

Returns a reference to the default configdataclass class

static default_com_cls()
Get the class for the default communication protocol used with this device.

Returns the type of the standard communication protocol for this device

output_off() → None
Switch DC voltage output off.

reset() → None
Reset of the interface: All setvalues are set to zero

abstract start()
Open the associated communication protocol.

stop() → None
Close the associated communication protocol.

class FuGProbusIVCommands(value=<no_arg>, names=None, module=None, type=None, start=1,
                           boundary=None)
Bases: hvl_ccb.utils.enum.NameEnum

An enumeration.

ADMODE = ('S', (<enum 'FuGMonitorModes'>, <class 'int'>))
CURRENT = ('I', (<class 'int'>, <class 'float'>))
EXECUTE = ('X', None)
```

```

EXECUTEONX = ('G', (<enum 'FuGDigitalVal'>, <class 'int'>))
    Wait for "X" to execute pending commands

ID = ('*IDN?', None)

OUTPUT = ('F', (<enum 'FuGDigitalVal'>, <class 'int'>))

POLARITY = ('P', (<enum 'FuGPolarities'>, <class 'int'>))

QUERY = ('?', None)

READBACKCHANNEL = ('N', (<enum 'FuGReadbackChannels'>, <class 'int'>))

RESET = ('=', None)

TERMINATOR = ('Y', (<enum 'FuGTerminators'>, <class 'int'>))

VOLTAGE = ('U', (<class 'int'>, <class 'float'>))

XOUTPUTS = ('R', <class 'int'>)
    TODO: the possible values are limited to 0..13

```

**class FuGProbusV**(*com, dev\_config=None*)  
Bases: *hvl\_ccb.dev.fug.FuGProbusIV*

FuG Probus V class which uses register based commands to control the power supplies

**get\_register**(*register: str*) → *str*  
get the value from a register

**Parameters** **register** – the register from which the value is requested

**Returns** the value of the register as a String

**set\_register**(*register: str, value: Union[int, float, str]*) → *None*  
generic method to set value to register

**Parameters**

- **register** – the name of the register to set the value
- **value** – which should be written to the register

**class FuGProbusVConfigRegisters**(*fug, super\_register: hvl\_ccb.dev.fug.FuGProbusVRegisterGroups*)  
Bases: *object*

Configuration and Status values, acc. 4.2.5

**property execute\_on\_x:** *hvl\_ccb.dev.fug.FuGDigitalVal*  
status of Execute-on-X

**Returns** FuGDigitalVal of the status

**property most\_recent\_error:** *hvl\_ccb.dev.fug.FuGErrorcodes*  
Reads the Error-Code of the most recent command

**Return FuGError**

**Raises** *FuGError* – if code is not “E0”

**property readback\_data:** *hvl\_ccb.dev.fug.FuGReadbackChannels*  
Preselection of readout data for Trigger-on-Talk

**Returns** index for the readback channel

**property srq\_mask:** *int*  
SRQ-Mask, Service-Request Enable status bits for SRQ 0: no SRQ Bit 2: SRQ on change of status to CC  
Bit 1: SRQ on change to CV

**Returns** representative integer value

**property srq\_status: str**  
SRQ-Statusbyte output as a decimal number: Bit 2: PS is in CC mode Bit 1: PS is in CV mode

**Returns** representative string

**property status: str**  
Statusbyte as a string of 0/1. Combined status (compatibel to Probus IV), MSB first: Bit 7: I-REG Bit 6: V-REG Bit 5: ON-Status Bit 4: 3-Reg Bit 3: X-Stat (polarity) Bit 2: Cal-Mode Bit 1: unused Bit 0: SEL-D

**Returns** string of 0/1

**property terminator: hvl\_ccb.dev.fug.FuGTerminators**  
Terminator character for answer strings from ADDA

**Returns** FuGTerminators

**class FuGProbusVDIRegisters(fug, super\_register: hvl\_ccb.dev.fug.FuGProbusVRegisterGroups)**  
Bases: object

Digital Inputs acc. 4.2.4

**property analog\_control: hvl\_ccb.dev.fug.FuGDigitalVal**

**Returns** shows 1 if power supply is controlled by the analog interface

**property calibration\_mode: hvl\_ccb.dev.fug.FuGDigitalVal**

**Returns** shows 1 if power supply is in calibration mode

**property cc\_mode: hvl\_ccb.dev.fug.FuGDigitalVal**

**Returns** shows 1 if power supply is in CC mode

**property cv\_mode: hvl\_ccb.dev.fug.FuGDigitalVal**

**Returns** shows 1 if power supply is in CV mode

**property digital\_control: hvl\_ccb.dev.fug.FuGDigitalVal**

**Returns** shows 1 if power supply is digitally controlled

**property on: hvl\_ccb.dev.fug.FuGDigitalVal**

**Returns** shows 1 if power supply ON

**property reg\_3: hvl\_ccb.dev.fug.FuGDigitalVal**  
For special applications.

**Returns** input from bit 3-REG

**property x\_stat: hvl\_ccb.dev.fug.FuGPolarities**

**Returns** polarity of HVPS with polarity reversal

**class FuGProbusVDORegisters(fug, super\_register: hvl\_ccb.dev.fug.FuGProbusVRegisterGroups)**  
Bases: object

Digital outputs acc. 4.2.2

---

**property out: Union[int, hvl\_ccb.dev.fug.FuGDigitalVal]**  
 Status of the output according to the last setting. This can differ from the actual state if output should only pulse.

**Returns** FuGDigitalVal

**property status: hvl\_ccb.dev.fug.FuGDigitalVal**  
 Returns the actual value of output. This can differ from the set value if pulse function is used.

**Returns** FuGDigitalVal

**class FuGProbusVMonitorRegisters(fug, super\_register: hvl\_ccb.dev.fug.FuGProbusVRegisterGroups)**  
 Bases: object

Analog monitors acc. 4.2.3

**property adc\_mode: hvl\_ccb.dev.fug.FuGMonitorModes**  
 The programmed resolution and integration time of the AD converter

**Returns** FuGMonitorModes

**property value: float**  
 Value from the monitor.

**Returns** a float value in V or A

**property value\_raw: float**  
 uncalibrated raw value from AD converter

**Returns** float value from ADC

**class FuGProbusVRegisterGroups(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**

Bases: *hvl\_ccb.utils.enum.NameEnum*

An enumeration.

**CONFIG = 'K'**

**INPUT = 'D'**

**MONITOR\_I = 'M1'**

**MONITOR\_V = 'M0'**

**OUTPUTONCMD = 'BON'**

**OUTPUTX0 = 'B0'**

**OUTPUTX1 = 'B1'**

**OUTPUTX2 = 'B2'**

**OUTPUTXCMD = 'BX'**

**SETCURRENT = 'S1'**

**SETVOLTAGE = 'S0'**

**class FuGProbusVSetRegisters(fug, super\_register: hvl\_ccb.dev.fug.FuGProbusVRegisterGroups)**  
 Bases: object

Setvalue control acc. 4.2.1 for the voltage and the current output

**property actualsetvalue: float**  
 The actual valid set value, which depends on the ramp function.

**Returns** actual valid set value

**property high\_resolution:** *hvl\_ccb.dev.fug.FuGDigitalVal*  
Status of the high resolution mode of the output.

**Return 0** normal operation

**Return 1** High Res. Mode

**property rampmode:** *hvl\_ccb.dev.fug.FuGRampModes*  
The set ramp mode to control the setvalue.

**Returns** the mode of the ramp as instance of FuGRampModes

**property ramprate:** *float*  
The set ramp rate in V/s.

**Returns** ramp rate in V/s

**property rampstate:** *hvl\_ccb.dev.fug.FuGDigitalVal*  
Status of ramp function.

**Return 0** if final setvalue is reached

**Return 1** if still ramping up

**property setvalue:** *float*  
For the voltage or current output this setvalue was programmed.

**Returns** the programmed setvalue

**class FuGRampModes(*value*)**

Bases: enum.IntEnum

An enumeration.

**FOLLOWRAMP = 1**

Follow the ramp up- and downwards

**IMMEDIATELY = 0**

Standard mode: no ramp

**ONLYUPWARDSOFFTOZERO = 4**

Follow the ramp up- and downwards, if output is OFF set value is zero

**RAMPUPWARDS = 2**

Follow the ramp only upwards, downwards immediately

**SPECIALRAMPUPWARDS = 3**

Follow a special ramp function only upwards

**class FuGReadbackChannels(*value*)**

Bases: enum.IntEnum

An enumeration.

**CURRENT = 1**

**FIRMWARE = 5**

**RATEDCURRENT = 4**

**RATEDVOLTAGE = 3**

**SN = 6**

**STATUSBYTE = 2**

**VOLTAGE = 0**

```
class FuGSerialCommunication(configuration)
Bases: hvl_ccb.comm.serial.SerialCommunication

Specific communication protocol implementation for FuG power supplies. Already predefines device-specific
protocol parameters in config.

static config_cls()
Return the default configdataclass class.

Returns a reference to the default configdataclass class

query(command: str) → str
Send a command to the interface and handle the status message. Eventually raises an exception.

Parameters command – Command to send

Raises FuGError – if the connection is broken or the error from the power source itself

Returns Answer from the interface or empty string

class FuGSerialCommunicationConfig(terminator: bytes = b'\n', encoding: str = 'utf-8',
encoding_error_handling: str = 'strict', wait_sec_read_text_nonempty:
Union[int, float] = 0.5, default_n_attempts_read_text_nonempty: int =
10, port: Union[str, NoneType] = None, baudrate: int = 9600, parity:
Union[str, hvl_ccb.comm.serial.SerialCommunicationParity] =
<SerialCommunicationParity.NONE: 'N'>, stopbits: Union[int,
hvl_ccb.comm.serial.SerialCommunicationStopbits] =
<SerialCommunicationStopbits.ONE: 1>, bytesize: Union[int,
hvl_ccb.comm.serial.SerialCommunicationBytesize] =
<SerialCommunicationBytesize.EIGHTBITS: 8>, timeout: Union[int,
float] = 3)
Bases: hvl_ccb.comm.serial.SerialCommunicationConfig

baudrate: int = 9600
Baudrate for FuG power supplies is 9600 baud

bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize] = 8
One byte is eight bits long

default_n_attempts_read_text_nonempty: int = 10
default number of attempts to read a non-empty text

force_value(fieldname, value)
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define post_force_value method with same signature as this method to do extra processing
after value has been forced on fieldname.

Parameters

- fieldname – name of the field
- value – value to assign

classmethod keys() → Sequence[str]
Returns a list of all configdataclass fields key-names.

Returns a list of strings containing all keys.

classmethod optional_defaults() → Dict[str, object]
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified
on instantiation.

Returns a list of strings containing all optional keys.
```

```

parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity] = 'N'
    FuG does not use parity

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

        Returns a list of strings containing all required keys.

stopbits: Union[int, hvl_ccb.comm.serial.SerialCommunicationStopbits] = 1
    FuG uses one stop bit

terminator: bytes = b'\n'
    The terminator is LF

timeout: Union[int, float] = 3
    use 3 seconds timeout as default

wait_sec_read_text_nonempty: Union[int, float] = 0.5
    default time to wait between attempts of reading a non-empty text


class FuGTerminators(value)
    Bases: enum.IntEnum

        An enumeration.

            CR = 3
            CRLF = 0
            LF = 2
            LFCR = 1

```

## [hvl\\_ccb.dev.heinzinger](#)



Device classes for Heinzinger Digital Interface I/II and Heinzinger PNC power supply.

The Heinzinger Digital Interface I/II is used for many Heinzinger power units. Manufacturer homepage: <https://www.heinzinger.com/products/accessories-and-more/digital-interfaces/>

The Heinzinger PNC series is a series of high voltage direct current power supplies. The class HeinzingerPNC is tested with two PNChp 60000-1neg and a PNChp 1500-1neg. Check the code carefully before using it with other PNC devices, especially PNC3p or PNCCap. Manufacturer homepage: <https://www.heinzinger.com/products/high-voltage/universal-high-voltage-power-supplies/>

```
class HeinzingerConfig(default_number_of_recordings: Union[int,  
                                hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum] = 1,  
                                number_of_decimals: int = 6, wait_sec_stop_commands: Union[int, float] = 0.5)
```

Bases: object

Device configuration dataclass for Heinzinger power supplies.

```
class RecordingsEnum(value)
```

Bases: enum.IntEnum

An enumeration.

```
EIGHT = 8
```

```
FOUR = 4
```

```
ONE = 1
```

```
SIXTEEN = 16
```

```
TWO = 2
```

```
clean_values()
```

```
default_number_of_recordings: Union[int,  
                                hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum] = 1
```

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
is_configdataclass = True
```

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
number_of_decimals: int = 6
```

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
wait_sec_stop_commands: Union[int, float] = 0.5
```

Time to wait after subsequent commands during stop (in seconds)

```
class HeinzingerDI(com, dev_config=None)
```

Bases: hvl\_ccb.dev.base.SingleCommDevice, abc.ABC

Heinzinger Digital Interface I/II device class

Sends basic SCPI commands and reads the answer. Only the standard instruction set from the manual is implemented.

**class OutputStatus(*value*)**

Bases: enum.IntEnum

Status of the voltage output

**OFF = 0**

**ON = 1**

**UNKNOWN = -1**

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**get\_current() → float**

Queries the set current of the Heinzinger PNC (not the measured current!).

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_interface\_version() → str**

Queries the version number of the digital interface.

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_number\_of\_recordings() → int**

Queries the number of recordings the device is using for average value calculation.

**Returns** int number of recordings

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_serial\_number() → str**

Ask the device for its serial number and returns the answer as a string.

**Returns** string containing the device serial number

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_voltage() → float**

Queries the set voltage of the Heinzinger PNC (not the measured voltage!).

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**measure\_current() → float**

Ask the Device to measure its output current and return the measurement result.

**Returns** measured current as float

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**measure\_voltage() → float**

Ask the Device to measure its output voltage and return the measurement result.

**Returns** measured voltage as float

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**output\_off()** → *None*

Switch DC voltage output off and updates the output status.

**Raises *SerialCommunicationIOError*** – when communication port is not opened

**output\_on()** → *None*

Switch DC voltage output on and updates the output status.

**Raises *SerialCommunicationIOError*** – when communication port is not opened

**property output\_status: *hvl\_ccb.dev.heinzinger.HeinzingerDI.OutputStatus*****reset\_interface()** → *None*

Reset of the digital interface; only Digital Interface I: Power supply is switched to the Local-Mode (Manual operation)

**Raises *SerialCommunicationIOError*** – when communication port is not opened

**set\_current(*value: Union[int, float]*)** → *None*

Sets the output current of the Heinzinger PNC to the given value.

**Parameters *value*** – current expressed in *self.unit\_current*

**Raises *SerialCommunicationIOError*** – when communication port is not opened

**set\_number\_of\_recordings(*value: Union[int,***

***hvl\_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum*)** → *None*

Sets the number of recordings the device is using for average value calculation. The possible values are 1, 2, 4, 8 and 16.

**Raises *SerialCommunicationIOError*** – when communication port is not opened

**set\_voltage(*value: Union[int, float]*)** → *None*

Sets the output voltage of the Heinzinger PNC to the given value.

**Parameters *value*** – voltage expressed in *self.unit\_voltage*

**Raises *SerialCommunicationIOError*** – when communication port is not opened

**abstract start()**

Opens the communication protocol.

**Raises *SerialCommunicationIOError*** – when communication port cannot be opened.

**stop()** → *None*

Stop the device. Closes also the communication protocol.

**class HeinzingerPNC(*com, dev\_config=None*)**

Bases: *hvl\_ccb.dev.heinzinger.HeinzingerDI*

Heinzinger PNC power supply device class.

The power supply is controlled over a Heinzinger Digital Interface I/II

**class UnitCurrent(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)**

Bases: *hvl\_ccb.utils.enum.AutoNumberNameEnum*

An enumeration.

**A = 3**

**UNKNOWN = 1**

**mA = 2**

```
class UnitVoltage(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: hvl_ccb.utils.enum.AutoNumberNameEnum

An enumeration.

UNKNOWN = 1
V = 2
kV = 3

identify_device() → None
Identify the device nominal voltage and current based on its serial number.

    Raises SerialCommunicationIOError – when communication port is not opened

property max_current: Union[int, float]
property max_current_hardware: Union[int, float]
property max_voltage: Union[int, float]
property max_voltage_hardware: Union[int, float]

set_current(value: Union[int, float]) → None
Sets the output current of the Heinzinger PNC to the given value.

    Parameters value – current expressed in self.unit_current
    Raises SerialCommunicationIOError – when communication port is not opened

set_voltage(value: Union[int, float]) → None
Sets the output voltage of the Heinzinger PNC to the given value.

    Parameters value – voltage expressed in self.unit_voltage
    Raises SerialCommunicationIOError – when communication port is not opened

start() → None
Opens the communication protocol and configures the device.

property unit_current: hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitCurrent
property unit_voltage: hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitVoltage

exception HeinzingerPNCDeviceNotRecognizedException
Bases: hvl\_ccb.dev.heinzinger.HeinzingerPNCError

Error indicating that the serial number of the device is not recognized.

exception HeinzingerPNCError
Bases: hvl\_ccb.dev.base.DeviceException

General error with the Heinzinger PNC voltage source.

exception HeinzingerPNCMaxCurrentExceededException
Bases: hvl\_ccb.dev.heinzinger.HeinzingerPNCError

Error indicating that program attempted to set the current to a value exceeding ‘max_current’.

exception HeinzingerPNCMaxVoltageExceededException
Bases: hvl\_ccb.dev.heinzinger.HeinzingerPNCError

Error indicating that program attempted to set the voltage to a value exceeding ‘max_voltage’.

class HeinzingerSerialCommunication(configuration)
Bases: hvl\_ccb.comm.serial.SerialCommunication
```

Specific communication protocol implementation for Heinzinger power supplies. Already predefines device-specific protocol parameters in config.

### **static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class HeinzingerSerialCommunicationConfig(terminator: bytes = b'\n', encoding: str = 'utf-8',
                                         encoding_error_handling: str = 'strict',
                                         wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                         default_n_attempts_read_text_nonempty: int = 40, port:
                                         Union[str, NoneType] = None, baudrate: int = 9600, parity:
                                         Union[str, hvl_ccb.comm.serial.SerialCommunicationParity]
                                         = <SerialCommunicationParity.NONE: 'N'>, stopbits:
                                         Union[int,
                                         hvl_ccb.comm.serial.SerialCommunicationStopbits] =
                                         <SerialCommunicationStopbits.ONE: 1>, bytesize:
                                         Union[int,
                                         hvl_ccb.comm.serial.SerialCommunicationBytesize] =
                                         <SerialCommunicationBytesize.EIGHTBITS: 8>, timeout:
                                         Union[int, float] = 3)
```

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

#### **baudrate: int = 9600**

Baudrate for Heinzinger power supplies is 9600 baud

#### **bytesize: Union[int, hvl\_ccb.comm.serial.SerialCommunicationBytesize] = 8**

One byte is eight bits long

#### **default\_n\_attempts\_read\_text\_nonempty: int = 40**

increased to 40 default number of attempts to read a non-empty text

#### **force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### **Parameters**

- **fieldname** – name of the field
- **value** – value to assign

#### **classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

#### **classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

#### **parity: Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = 'N'**

Heinzinger does not use parity

#### **classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.


**stopbits:** Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = 1  
Heinzinger uses one stop bit

**terminator:** bytes = b'\n'  
The terminator is LF

**timeout:** Union[int, float] = 3  
use 3 seconds timeout as default

**wait\_sec\_read\_text\_nonempty:** Union[int, float] = 0.5  
default time to wait between attempts of reading a non-empty text

## hvl\_ccb.dev.labjack



A  
LabJack T-series devices wrapper around the LabJack's LJM Library; see <https://labjack.com/ljm> . The wrapper was originally developed and tested for a LabJack T7-PRO device.

## Extra installation

To use this LabJack T-series devices wrapper:

1. install the hvl\_ccb package with a labjack extra feature:

```
$ pip install "hvl_ccb[labjack]"
```

this will install the Python bindings for the library.

2. install the library - follow instruction in <https://labjack.com/support/software/installers/ljm> .

**class LabJack**(*com*, *dev\_config=None*)  
Bases: *hvl\_ccb.dev.base.SingleCommDevice*

LabJack Device.

This class is tested with a LabJack T7-Pro and should also work with T4 and T7 devices communicating through the LJM Library. Other or older hardware versions and variants of LabJack devices are not supported.

**class AInRange**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)  
Bases: *hvl\_ccb.utils.enum.StrEnumBase*

An enumeration.

**ONE = 1.0**

**ONE\_HUNDREDTH = 0.01**

**ONE\_TENTH = 0.1**

**TEN = 10.0**

**property value: float**

**class BitLimit**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)  
Bases: *aenum.IntEnum*

Maximum integer values for clock settings

**THIRTY\_TWO\_BIT = 4294967295**

**class CalMicroAmpere**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*,  
*boundary=None*)

Bases: *aenum.Enum*

Pre-defined microampere (uA) values for calibration current source query.

**TEN = '10uA'**

**TWO\_HUNDRED = '200uA'**

**class CjcType**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)  
Bases: *hvl\_ccb.utils.enum.NameEnum*

CJC slope and offset

**internal = (1, 0)**

**lm34 = (55.56, 255.37)**

**class ClockFrequency**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*,  
*boundary=None*)

Bases: *aenum.IntEnum*

Available clock frequencies, in Hz

```
FIVE_MHZ = 5000000
FORTY_MHZ = 40000000
MAXIMUM = 80000000
MINIMUM = 312500
TEN_MHZ = 10000000
TWELVE_HUNDRED_FIFTY_KHZ = 1250000
TWENTY_FIVE_HUNDRED_KHZ = 2500000
TWENTY_MHZ = 20000000

DIOChannel
alias of hvl_ccb._dev.labjack.TSeriesDIOChannel

class DIOStatus(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: aenum.IntEnum
State of a digital I/O channel.

HIGH = 1
LOW = 0

class DeviceType(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: hvl_ccb.utils.enum.AutoNumberNameEnum
LabJack device types.

Can be also looked up by ambiguous Product ID (p_id) or by instance name: `python
LabJackDeviceType(4) is LabJackDeviceType('T4')` 

ANY = 1
T4 = 2
T7 = 3
T7_PRO = 4

classmethod get_by_p_id(p_id: int) → Union[hvl_ccb._dev.labjack.DeviceType,
                                              List[hvl_ccb._dev.labjack.DeviceType]]
Get LabJack device type instance via LabJack product ID.

Note: Product ID is not unambiguous for LabJack devices.
Parameters p_id – Product ID of a LabJack device
Returns Instance or list of instances of LabJackDeviceType
Raises ValueError – when Product ID is unknown

class TemperatureUnit(value=<no_arg>, names=None, module=None, type=None, start=1,
                      boundary=None)
Bases: hvl_ccb.utils.enum.NameEnum
Temperature unit (to be returned)

C = 1
F = 2
K = 0
```

---

```
class ThermocoupleType(value=<no_arg>, names=None, module=None, type=None, start=1,
                           boundary=None)
```

Bases: `hvl_ccb.utils.enum.NameEnum`

Thermocouple type; NONE means disable thermocouple mode.

**C = 30**

**E = 20**

**J = 21**

**K = 22**

**NONE = 0**

**PT100 = 40**

**PT1000 = 42**

**PT500 = 41**

**R = 23**

**S = 25**

**T = 24**

**config\_high\_pulse**(address: *Union[str, hvl\_ccb.\_dev.labjack.TSeriesDIOChannel]*, t\_start: *Union[int, float]*, t\_width: *Union[int, float]*, n\_pulses: int = 1) → `None`

Configures one FIO channel to send a timed HIGH pulse. Configure multiple channels to send pulses with relative timing accuracy. Times have a maximum resolution of 1e-7 seconds @ 10 MHz. :param address: FIO channel: [T7] FIO0;2;3;4;5. [T4] FIO6;7. :raises LabJackError if address is not supported. :param t\_start: pulse start time in seconds. :raises ValueError: if t\_start is negative or would exceed the clock period. :param t\_width: duration of high pulse, in seconds. :raises ValueError: if t\_width is negative or would exceed the clock period. :param n\_pulses: number of pulses to be sent; single pulse default. :raises TypeError if n\_pulses is not of type int. :raises ValueError if n\_pulses is negative or exceeds the 32bit limit.

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**disable\_pulses**(\*addresses: *Optional[Union[str, hvl\_ccb.\_dev.labjack.TSeriesDIOChannel]]*) → `None`

Disable previously configured pulse channels. :param addresses: tuple of FIO addresses. All channels disabled if no argument is passed.

**enable\_clock**(clock\_enabled: bool) → `None`

Enable/disable LabJack clock to configure or send pulses. :param clock\_enabled: True -> enable, False -> disable. :raises TypeError: if clock\_enabled is not of type bool

**get\_ain**(\*channels: int) → *Union[float, Sequence[float]]*

Read currently measured value (voltage, resistance, ...) from one or more of analog inputs.

**Parameters** **channels** – AIN number or numbers (0..254)

**Returns** the read value (voltage, resistance, ...) as `float` or `tuple` of them in case multiple channels given

**get\_cal\_current\_source**(name: *Union[str, CalMicroAmpere]*) → float

This function will return the calibration of the chosen current source, this ist not a measurement!

The value was stored during fabrication.

**Parameters** `name` – ‘200uA’ or ‘10uA’ current source

**Returns** calibration of the chosen current source in ampere

**get\_clock()** → Dict[str, Union[int, float]]

Return clock settings read from LabJack.

**get\_digital\_input**(*address*: Union[str, hvl\_ccb.\_dev.labjack.TSeriesDIOChannel]) →

*hvl\_ccb.dev.labjack.LabJack.DIOStatus*

Get the value of a digital input.

allowed names for T7 (Pro): FIO0 - FIO7, EIO0 - EIO 7, CIO0- CIO3, MIO0 - MIO2 :param *address*: name of the output -> ‘FIO0’ :return: HIGH when *address* DIO is high, and LOW when *address* DIO is low

**get\_product\_id()** → int

This function returns the product ID reported by the connected device.

Attention: returns 7 for both T7 and T7-Pro devices!

**Returns** integer product ID of the device

**get\_product\_name**(*force\_query\_id=False*) → str

This function will return the product name based on product ID reported by the device.

Attention: returns “T7” for both T7 and T7-Pro devices!

**Parameters** `force_query_id` – boolean flag to force `get_product_id` query to device instead of using cached device type from previous queries.

**Returns** device name string, compatible with *LabJack.DeviceType*

**get\_product\_type**(*force\_query\_id: bool = False*) → *hvl\_ccb.\_dev.labjack.DeviceType*

This function will return the device type based on reported device type and in case of unambiguity based on configuration of device’s communication protocol (e.g. for “T7” and “T7\_PRO” devices), or, if not available first matching.

**Parameters** `force_query_id` – boolean flag to force `get_product_id` query to device instead of using cached device type from previous queries.

**Returns** *DeviceType* instance

**Raises** *LabJackIdentifierDIOError* – when read Product ID is unknown

**get\_sbus\_rh**(*number: int*) → float

Read the relative humidity value from a serial SBUS sensor.

**Parameters** `number` – port number (0..22)

**Returns** relative humidity in %RH

**get\_sbus\_temp**(*number: int*) → float

Read the temperature value from a serial SBUS sensor.

**Parameters** `number` – port number (0..22)

**Returns** temperature in Kelvin

**get\_serial\_number()** → int

Returns the serial number of the connected LabJack.

**Returns** Serial number.

**read\_resistance**(*channel: int*) → float

Read resistance from specified channel.

**Parameters** `channel` – channel with resistor

**Returns** resistance value with 2 decimal places

**read\_thermocouple**(`pos_channel: int`) → `float`

Read the temperature of a connected thermocouple.

**Parameters** `pos_channel` – is the AIN number of the positive pin

**Returns** temperature in specified unit

**send\_pulses**(\*`addresses: Union[str, hvl_ccb._dev.labjack.TSeriesDIOChannel]`) → `None`

Sends pre-configured pulses for specified addresses. :param `addresses`: tuple of FIO addresses :raises LabJackError if an address has not been configured.

**set\_ain\_differential**(`pos_channel: int, differential: bool`) → `None`

Sets an analog input to differential mode or not. T7-specific: For base differential channels, positive must be even channel from 0-12 and negative must be positive+1. For extended channels 16-127, see Mux80 datasheet.

**Parameters**

- `pos_channel` – is the AIN number (0..12)

- `differential` – True or False

**Raises** `LabJackError` – if parameters are unsupported

**set\_ain\_range**(`channel: int, vrange: Union[Real, AInRange]`) → `None`

Set the range of an analog input port.

**Parameters**

- `channel` – is the AIN number (0..254)

- `vrange` – is the voltage range to be set

**set\_ain\_resistance**(`channel: int, vrange: Union[Real, AInRange], resolution: int`) → `None`

Set the specified channel to resistance mode. It utilized the 200uA current source of the LabJack.

**Parameters**

- `channel` – channel that should measure the resistance

- `vrange` – voltage range of the channel

- `resolution` – resolution index of the channel T4: 0-5, T7: 0-8, T7-Pro 0-12

**set\_ain\_resolution**(`channel: int, resolution: int`) → `None`

Set the resolution index of an analog input port.

**Parameters**

- `channel` – is the AIN number (0..254)

- `resolution` – is the resolution index within 0...`get\_product\_type().ain\_max\_resolution` range; 0 will set the resolution index to default value.

**set\_ain\_thermocouple**(`pos_channel: int, thermocouple: Union[None, str, ThermocoupleType], cjc_address: int = 60050, cjc_type: Union[str, CjcType] = CjcType.internal, vrange: Union[Real, AInRange] = AInRange.ONE_HUNDREDTH, resolution: int = 10, unit: Union[str, TemperatureUnit] = TemperatureUnit.K`) → `None`

Set the analog input channel to thermocouple mode.

**Parameters**

- `pos_channel` – is the analog input channel of the positive part of the differential pair

- **thermocouple** – None to disable thermocouple mode, or string specifying the thermocouple type
- **cjc\_address** – modbus register address to read the CJC temperature
- **cjc\_type** – determines cjc slope and offset, ‘internal’ or ‘lm34’
- **vrange** – measurement voltage range
- **resolution** – resolution index (T7-Pro: 0-12)
- **unit** – is the temperature unit to be returned (‘K’, ‘C’ or ‘F’)

**Raises** `LabJackError` – if parameters are unsupported

**set\_analog\_output**(*channel: int, value: Union[int, float]*) → *None*

Set the voltage of a analog output port

#### Parameters

- **channel** – DAC channel number 1/0
- **value** – The output voltage value 0-5 Volts int/float

**set\_clock**(*clock\_frequency: Union[Number, ClockFrequency] = 10000000, clock\_period: Number = 1*) → *None*

Configure LabJack clock for pulse out feature. :param `clock_frequency`: clock frequency in Hz; default 10 MHz for base 10. :raises `ValueError`: if `clock_frequency` is not allowed (see `ClockFrequency`). :param `clock_period`: clock roll time in seconds; default 1s, 0 for max. :raises `ValueError`: if `clock_period` exceeds the 32bit tick limit. Clock period determines pulse spacing when using multi-pulse settings. Ensure period exceeds maximum intended pulse end time.

**set\_digital\_output**(*address: str, state: Union[int, DIOStatus]*) → *None*

Set the value of a digital output.

#### Parameters

- **address** – name of the output -> ‘FIO0’
- **state** – state of the output -> `DIOStatus` instance or corresponding `int` value

**start()** → *None*

Start the Device.

**stop()** → *None*

Stop the Device.

### exception `LabJackError`


Bases: `hvl_ccb.dev.base.DeviceException`

General Exception for the LabJack device.

### exception `LabJackIdentifierDIOError`

Bases: `hvl_ccb.dev.labjack.LabJackError`

Error indicating a wrong DIO identifier

**[hvl\\_ccb.dev.lauda](#)**

Device class for controlling a Lauda PRO RP245E, circulation chiller over TCP.

```

class LaudaProRp245e(com, dev_config=None)
  Bases: hvl\_ccb.dev.base.SingleCommDevice

  Lauda RP245E circulation chiller class.

  static config_cls() → Type[hvl\_ccb.dev.lauda.LaudaProRp245eConfig]
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

  continue_ramp() → str
    Continue current ramp program.

    Returns reply of the device to the last call of “query”

  static default_com_cls() → Type[hvl\_ccb.dev.lauda.LaudaProRp245eTcpCommunication]
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

  get_bath_temp() → float
    :return : float value of measured lauda bath temp in °C

  get_device_type() → str
    :return : Connected Lauda device type (for connection/com test)

  pause() → str
    Stop temperature control and pump.

    Returns reply of the device to the last call of “query”

  pause_ramp() → str
    Pause current ramp program.

    Returns reply of the device to the last call of “query”

  reset_ramp() → str
    Delete all segments from current ramp program.

    Returns reply of the device to the last call of “query”

  run() → str
    Start temperature control & pump.

    Returns reply of the device to the last call of “query”
  
```

**set\_control\_mode**(*mod: Union[int, hvl\_ccb.dev.lauda.LaudaProRp245eConfig.ExtControlModeEnum] = ExtControlModeEnum.INTERNAL*) → str

Define control mode. 0 = INTERNAL (control bath temp), 1 = EXPT100 (pt100 attached to chiller), 2 = ANALOG, 3 = SERIAL, 4 = USB, 5 = ETH (to be used when passing the ext. temp. via ethernet) (temperature then needs to be passed every second, when not using options 3, 4, or 5)

**Parameters** **mod** – temp control mode (control internal temp or external temp).

**Returns** reply of the device to the last call of “query” (“OK”, if command was recognized”)

**set\_external\_temp**(*external\_temp: float = 20.0*) → str

Pass value of external controlled temperature. Should be done every second, when control of external temperature is active. Has to be done right before control of external temperature is activated.

**Parameters** **external\_temp** – current value of external temperature to be controlled.

**Returns** reply of the device to the last call of “query”

**set\_pump\_level**(*pump\_level: int = 6*) → str

Set pump level Raises ValueError, if pump level is invalid.

**Parameters** **pump\_level** – pump level.

**Returns** reply of the device to the last call of “query”

**set\_ramp\_iterations**(*num: int = 1*) → str

Define number of ramp program cycles.

**Parameters** **num** – number of program cycles to be performed.

**Returns** reply of the device to the last call of “query”

**set\_ramp\_program**(*program: int = 1*) → str

Define ramp program for following ramp commands. Raises ValueError if maximum number of ramp programs (5) is exceeded.

**Parameters** **program** – Number of ramp program to be activated for following commands.

**Returns** reply of the device to the last call of “query”

**set\_ramp\_segment**(*temp: float = 20.0, dur: int = 0, tol: float = 0.0, pump: int = 6*) → str

Define segment of current ramp program - will be attached to current program. Raises ValueError, if pump level is invalid.

#### Parameters

- **temp** – target temperature of current ramp segment
- **dur** – duration in minutes, in which target temperature should be reached
- **tol** – tolerance at which target temperature should be reached (for 0.00, next segment is started after dur has passed).
- **pump** – pump level to be used for this program segment.

**Returns** reply of the device to the last call of “query”

**set\_temp\_set\_point**(*temp\_set\_point: float = 20.0*) → str

Define temperature set point

**Parameters** **temp\_set\_point** – temperature set point.

**Returns** reply of the device to the last call of “query”

**start()** → *None*

Start this device.

---

**start\_ramp()** → str  
Start current ramp program.

**Returns** reply of the device to the last call of “query”

**stop()** → *None*  
Stop this device. Disables access and closes the communication protocol.

**stop\_ramp()** → str  
Stop current ramp program.

**Returns** reply of the device to the last call of “query”

**validate\_pump\_level(*level: int*)**  
Validates pump level. Raises ValueError, if pump level is incorrect. :param level: pump level, integer

```
class LaudaProRp245eCommand(value=<no_arg>, names=None, module=None, type=None, start=1,
                                boundary=None)
```

Bases: *hvl\_ccb.utils.enum.ValueEnum*

Commands for Lauda PRO RP245E Chiller Command strings most often need to be complimented with a parameter (attached as a string) before being sent to the device. Commands implemented as defined in “Lauda Betriebsanleitung fuer PRO Badthermostate und Umwaelzthermostate” pages 42 - 49

**BATH\_TEMP = 'IN\_PV\_00'**  
Request internal bath temperature

**COM\_TIME\_OUT = 'OUT\_SP\_08\_'**  
Define communication time out

**CONT\_MODE = 'OUT\_MODE\_01\_'**  
Set control mode 1=internal, 2=ext. analog, 3=ext. serial, 4=USB, 5=ethernet

**DEVICE\_TYPE = 'TYPE'**  
Request device type

**EXTERNAL\_TEMP = 'OUT\_PV\_05\_'**  
Pass on external controlled temperature

**LOWER\_TEMP = 'OUT\_SP\_05\_'**  
Define lower temp limit

**OPERATION\_MODE = 'OUT\_SP\_02\_'**  
Define operation mode

**PUMP\_LEVEL = 'OUT\_SP\_01\_'**  
Define pump level 1-8

**RAMP\_CONTINUE = 'RMP\_CONT'**  
Continue a paused ramp program

**RAMP\_DELETE = 'RMP\_RESET'**  
Reset a selected ramp program

**RAMP\_ITERATIONS = 'RMP\_OUT\_02\_'**  
Define how often a ramp program should be iterated

**RAMP\_PAUSE = 'RMP\_PAUSE'**  
Pause a selected ramp program

**RAMP\_SELECT = 'RMP\_SELECT\_'**  
Select a ramp program (target for all further ramp commands)

```
RAMP_SET = 'RMP_OUT_00_'
Define parameters of a selected ramp program

RAMP_START = 'RMP_START'
Start a selected ramp program

RAMP_STOP = 'RMP_STOP'
Stop a running ramp program

START = 'START'
Start temp control (pump and heating/cooling)

STOP = 'STOP'
Stop temp control (pump and heating/cooling)

TEMP_SET_POINT = 'OUT_SP_00_'
Define temperature set point

UPPER_TEMP = 'OUT_SP_04_'
Define upper temp limit

build_str(param: str = "", terminator: str = '\r\n')
Build a command string for sending to the device
```

#### Parameters

- **param** – Command's parameter given as string
- **terminator** – Command's terminator

**Returns** Command's string with a parameter and terminator

### exception LaudaProRp245eCommandError

Bases: `hvl_ccb.dev.base.DeviceException`

Exception raised when an error is returned upon a command.

```
class LaudaProRp245eConfig(temp_set_point_init: Union[int, float] = 20.0, pump_init: int = 6, upper_temp:
    Union[int, float] = 80.0, lower_temp: Union[int, float] = -55.0, com_time_out:
    Union[int, float] = 0, max_pump_level: int = 8, max_pr_number: int = 5,
    operation_mode: Union[int,
        hvl_ccb.dev.lauda.LaudaProRp245eConfig.OperationModeEnum] =
        OperationModeEnum.AUTO, control_mode: Union[int,
            hvl_ccb.dev.lauda.LaudaProRp245eConfig.ExtControlModeEnum] =
            ExtControlModeEnum.INTERNAL)
```

Bases: `object`

Configuration for the Lauda RP245E circulation chiller.

### class ExtControlModeEnum(value)

Bases: `enum.IntEnum`

Source for definition of external, controlled temperature (option 2, 3 and 4 are not available with current configuration of the Lauda RP245E, add-on hardware would required)

```
ANALOG = 2
ETH = 5
EXPT100 = 1
INTERNAL = 0
SERIAL = 3
```

```

USB = 4
class OperationModeEnum(value)
    Bases: enum.IntEnum
        Operation Mode (Cooling OFF/Cooling On/AUTO - set to AUTO)
AUTO = 2
    Automatically select heating/cooling
COOLOFF = 0
COOLON = 1
clean_values() → None
com_time_out: Union[int, float] = 0
    Communication time out (0 = OFF)
control_mode: Union[int, hvl_ccb.dev.lauda.LaudaProRp245eConfig.ExtControlModeEnum] = 0
force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.
    NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.
Parameters

- fieldname – name of the field
- value – value to assign

is_configdataclass = True
classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.
    Returns a list of strings containing all keys.
lower_temp: Union[int, float] = -55.0
    Lower temperature limit (safe for Galden HT135 cooling liquid)
max_pr_number: int = 5
    Maximum number of ramp programs that can be stored in the memory of the chiller
max_pump_level: int = 8
    Highest pump level of the chiller
operation_mode: Union[int, hvl_ccb.dev.lauda.LaudaProRp245eConfig.OperationModeEnum] = 2
classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.
    Returns a list of strings containing all optional keys.
pump_init: int = 6
    Default pump Level
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

```

**Returns** a list of strings containing all required keys.

**temp\_set\_point\_init:** Union[int, float] = 20.0

Default temperature set point

**upper\_temp:** Union[int, float] = 80.0

Upper temperature limit (safe for Galden HT135 cooling liquid)

**class LaudaProRp245eTcpCommunication**(*configuration*)

Bases: *hvl\_ccb.comm.tcp.Tcp*

Implements the Communication Protocol for Lauda PRO RP245E TCP connection.

**close()** → *None*

Close the Lauda PRO RP245E TCP connection.

**static config\_cls()** → Type[*hvl\_ccb.dev.lauda.LaudaProRp245eTcpCommunicationConfig*]

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**open()** → *None*

Open the Lauda PRO RP245E TCP connection.

**Raises** *LaudaProRp245eCommandError* – if the connection fails.

**query\_command**(*command*: *hvl\_ccb.dev.lauda.LaudaProRp245eCommand*, *param*: str = '') → str

Send and receive function. E.g. to be used when setting/changing device setting. :param command: first part of command string, defined in *LaudaProRp245eCommand* :param param: second part of command string, parameter (by default '') :return: None

**read()** → str

Receive value function. :return: reply from device as a string, the terminator, as well as the 'OK' stripped from the reply to make it directly useful as a value (e.g. in case the internal bath temperature is requested)

**write\_command**(*command*: *hvl\_ccb.dev.lauda.LaudaProRp245eCommand*, *param*: str = '') → *None*

Send command function. :param command: first part of command string, defined in *LaudaProRp245eCommand* :param param: second part of command string, parameter (by default '') :return: None

**class LaudaProRp245eTcpCommunicationConfig**(*host*: Union[str, *ipaddress.IPv4Address*,  
*ipaddress.IPv6Address*], *port*: int = 54321, *bufsize*: int =  
1024, *wait\_sec\_pre\_read\_or\_write*: Union[int, float] =  
0.005, *terminator*: str = '\r\n')

Bases: *hvl\_ccb.comm.tcp.TcpCommunicationConfig*

Configuration dataclass for *LaudaProRp245eTcpCommunication*.

**clean\_values()** → *None*

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.


**terminator: str = '\r\n'**

The terminator character

**wait\_sec\_pre\_read\_or\_write: Union[int, float] = 0.005**

Delay time between commands in seconds

## hvl\_ccb.dev.mbw973



Device class for controlling a MBW 973 SF6 Analyzer over a serial connection.

The MBW 973 is a gas analyzer designed for gas insulated switchgear and measures humidity, SF6 purity and SO2 contamination in one go. Manufacturer homepage: <https://www.mbw.ch/products/sf6-gas-analysis/973-sf6-analyzer/>

**class MBW973(com, dev\_config=None)**  
Bases: *hvl\_ccb.dev.base.SingleCommDevice*

MBW 973 dew point mirror device class.

**static config\_cls()**  
Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**  
Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**is\_done()** → bool

Poll status of the dew point mirror and return True, if all measurements are done.

**Returns** True, if all measurements are done; False otherwise.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**read**(*cast\_type*: typing.Type = <class 'str'>)

Read value from *self.com* and cast to *cast\_type*. Raises `ValueError` if read text (*str*) is not convertible to *cast\_type*, e.g. to `float` or to `int`.

**Returns** Read value of *cast\_type* type.

**read\_float()** → float

Convenience wrapper for *self.read()*, with typing hint for return value.

**Returns** Read `float` value.

**read\_int()** → int

Convenience wrapper for *self.read()*, with typing hint for return value.

**Returns** Read `int` value.

**read\_measurements()** → Dict[str, float]

Read out measurement values and return them as a dictionary.

**Returns** Dictionary with values.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**set\_measuring\_options**(*humidity*: bool = True, *sf6\_purity*: bool = False) → None

Send measuring options to the dew point mirror.

#### Parameters

- **humidity** – Perform humidity test or not?
- **sf6\_purity** – Perform SF6 purity test or not?

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**start()** → None

Start this device. Opens the communication protocol and retrieves the set measurement options from the device.

**Raises** `SerialCommunicationIOError` – when communication port cannot be opened.

**start\_control()** → None

Start dew point control to acquire a new value set.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**stop()** → None

Stop the device. Closes also the communication protocol.

**write**(*value*) → None

Send *value* to *self.com*.

**Parameters** **value** – Value to send, converted to *str*.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**class MBW973Config**(*polling\_interval*: Union[int, float] = 2)

Bases: object

Device configuration dataclass for MBW973.

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**polling\_interval: Union[int, float] = 2**

Polling period for *is\_done* status queries [in seconds].

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception MBW973ControlRunningException**

Bases: *hvl\_ccb.dev.mbw973.MBW973Error*

Error indicating there is still a measurement running, and a new one cannot be started.

**exception MBW973Error**

Bases: *hvl\_ccb.dev.base.DeviceException*

General error with the MBW973 dew point mirror device.

**exception MBW973PumpRunningException**

Bases: *hvl\_ccb.dev.mbw973.MBW973Error*

Error indicating the pump of the dew point mirror is still recovering gas, unable to start a new measurement.

**class MBW973SerialCommunication(configuration)**

Bases: *hvl\_ccb.comm.serial.SerialCommunication*

Specific communication protocol implementation for the MBW973 dew point mirror. Already predefines device-specific protocol parameters in config.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class MBW973SerialCommunicationConfig(terminator: bytes = b'\r', encoding: str = 'utf-8',
                                         encoding_error_handling: str = 'strict',
                                         wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                         default_n_attempts_read_text_nonempty: int = 10, port: Union[str,
                                         NoneType] = None, baudrate: int = 9600, parity: Union[str,
                                         hvl_ccb.comm.serial.SerialCommunicationParity] =
                                         <SerialCommunicationParity.NONE: 'N'>, stopbits: Union[int,
                                         hvl_ccb.comm.serial.SerialCommunicationStopbits] =
                                         <SerialCommunicationStopbits.ONE: 1>, bytesize: Union[int,
                                         hvl_ccb.comm.serial.SerialCommunicationBytesize] =
                                         <SerialCommunicationBytesize.EIGHTBITS: 8>, timeout:
                                         Union[int, float] = 3)
```

Bases: `hvl_ccb.comm.serial.SerialCommunicationConfig`

**baudrate: int = 9600**

Baudrate for MBW973 is 9600 baud

**bytesize: Union[int, hvl\_ccb.comm.serial.SerialCommunicationBytesize] = 8**

One byte is eight bits long

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity: Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = 'N'**

MBW973 does not use parity

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits: Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = 1**

MBW973 does use one stop bit


**terminator: bytes = b'\r'**

The terminator is only CR

**timeout: Union[int, float] = 3**

use 3 seconds timeout as default

## hvl\_ccb.dev.newport



Device class for Newport SMC100PP stepper motor controller with serial communication.

The SMC100PP is a single axis motion controller/driver for stepper motors up to 48 VDC at 1.5 A rms. Up to 31 controllers can be networked through the internal RS-485 communication link.

Manufacturer homepage: <https://www.newport.com/f/smci00-single-axis-dc-or-stepper-motion-controller>

```
class NewportConfigCommands(value=<no_arg>, names=None, module=None, type=None, start=1,
                             boundary=None)
```

Bases: `hvl_ccb.utils.enum.NameEnum`

Commands predefined by the communication protocol of the SMC100PP

```
AC = 'acceleration'
BA = 'backlash_compensation'
BH = 'hysteresis_compensation'
FRM = 'micro_step_per_full_step_factor'
FRS = 'motion_distance_per_full_step'
HT = 'home_search_type'
JR = 'jerk_time'
OH = 'home_search_velocity'
OT = 'home_search_timeout'
QIL = 'peak_output_current_limit'
SA = 'rs485_address'
SL = 'negative_software_limit'
SR = 'positive_software_limit'
VA = 'velocity'
VB = 'base_velocity'
ZX = 'stage_configuration'
```

```
exception NewportControllerError
```

Bases: `hvl_ccb.dev.newport.NewportError`

Error with the Newport controller.

```
exception NewportError
```

Bases: `hvl_ccb.dev.base.DeviceException`

General Exception for Newport Device

**exception NewportMotorError**

Bases: *hvl\_ccb.dev.newport.NewportError*

Error with the Newport motor.

**exception NewportMotorPowerSupplyWasCutError**

Bases: *hvl\_ccb.dev.newport.NewportError*

Error with the Newport motor after the power supply was cut and then restored, without interrupting the communication with the controller.

**class NewportSMC100PP(*com, dev\_config=None*)**

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Device class of the Newport motor controller SMC100PP

**class MotorErrors(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)**

Bases: *aenum.Enum*

Possible motor errors reported by the motor during `get_state()`.

**DC\_VOLTAGE\_TOO\_LOW = 3**

**FOLLOWING\_ERROR = 6**

**HOMING\_TIMEOUT = 5**

**NED\_END\_OF\_TURN = 11**

**OUTPUT\_POWER\_EXCEEDED = 2**

**PEAK\_CURRENT\_LIMIT = 9**

**POS\_END\_OF\_TURN = 10**

**RMS\_CURRENT\_LIMIT = 8**

**SHORT\_CIRCUIT = 7**

**WRONG\_ESP\_STAGE = 4**

**class StateMessages(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)**

Bases: *aenum.Enum*

Possible messages returned by the controller on `get_state()` query.

**CONFIG = '14'**

**DISABLE\_FROM\_JOGGING = '3E'**

**DISABLE\_FROM\_MOVING = '3D'**

**DISABLE\_FROM\_READY = '3C'**

**HOMING\_FROM\_RS232 = '1E'**

**HOMING\_FROM\_SMC = '1F'**

**JOGGING\_FROM\_DISABLE = '47'**

**JOGGING\_FROM\_READY = '46'**

**MOVING = '28'**

**NO\_REF\_ESP\_STAGE\_ERROR = '10'**

**NO\_REF\_FROM\_CONFIG = '0C'**

```

NO_REF_FROM_DISABLED = '0D'
NO_REF_FROM_HOMING = '0B'
NO_REF_FROM_JOGGING = '11'
NO_REF_FROM_MOVING = '0F'
NO_REF_FROM_READY = '0E'
NO_REF_FROM_RESET = '0A'
READY_FROM_DISABLE = '34'
READY_FROM_HOMING = '32'
READY_FROM_JOGGING = '35'
READY_FROM_MOVING = '33'

```

**States**

alias of `hvl_ccb.dev.newport.NewportStates`

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**exit\_configuration(*add*: Optional[int] = None) → None**

Exit the CONFIGURATION state and go back to the NOT REFERENCED state. All configuration parameters are saved to the device's memory.

**Parameters** *add* – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

**get\_acceleration(*add*: Optional[int] = None) → Union[int, float]**

Leave the configuration state. The configuration parameters are saved to the device's memory.

**Parameters** *add* – controller address (1 to 31)

**Returns** acceleration (preset units/s^2), value between 1e-6 and 1e12

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

**get\_controller\_information(*add*: Optional[int] = None) → str**

Get information on the controller name and driver version

**Parameters** *add* – controller address (1 to 31)

**Returns** controller information

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_motor\_configuration**(*add*: *Optional[int]* = *None*) → Dict[str, float]

Query the motor configuration and returns it in a dictionary.

**Parameters** *add* – controller address (1 to 31)

**Returns** dictionary containing the motor's configuration

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_move\_duration**(*dist*: *Union[int, float]*, *add*: *Optional[int]* = *None*) → float

Estimate the time necessary to move the motor of the specified distance.

**Parameters**

- **dist** – distance to travel
- **add** – controller address (1 to 31), defaults to self.address

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_negative\_software\_limit**(*add*: *Optional[int]* = *None*) → Union[int, float]

Get the negative software limit (the maximum position that the motor is allowed to travel to towards the left).

**Parameters** *add* – controller address (1 to 31)

**Returns** negative software limit (preset units), value between -1e12 and 0

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_position**(*add*: *Optional[int]* = *None*) → float

Returns the value of the current position.

**Parameters** *add* – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error
- *NewportUncertainPositionError* – if the position is ambiguous

**get\_positive\_software\_limit**(*add: Optional[int] = None*) → Union[int, float]

Get the positive software limit (the maximum position that the motor is allowed to travel to towards the right).

**Parameters** **add** – controller address (1 to 31)

**Returns** positive software limit (preset units), value between 0 and 1e12

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**get\_state**(*add: int = None*) → *StateMessages*

Check on the motor errors and the controller state

**Parameters** **add** – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error
- **NewportMotorError** – if the motor reports an error

**Returns** state message from the device (member of StateMessages)

**go\_home**(*add: Optional[int] = None*) → *None*

Move the motor to its home position.

**Parameters** **add** – controller address (1 to 31), defaults to self.address

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**go\_to\_configuration**(*add: Optional[int] = None*) → *None*

This method is executed during start(). It can also be executed after a reset(). The controller is put in CONFIG state, where configuration parameters can be changed.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**initialize**(*add: Optional[int] = None*) → *None*

Puts the controller from the NOT\_REF state to the READY state. Sends the motor to its “home” position.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed

- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**move\_to\_absolute\_position**(*pos: Union[int, float]*, *add: Optional[int] = None*) → *None*  
Move the motor to the specified position.

#### Parameters

- **pos** – target absolute position (affected by the configured offset)
- **add** – controller address (1 to 31), defaults to self.address

#### Raises

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**move\_to\_relative\_position**(*pos: Union[int, float]*, *add: Optional[int] = None*) → *None*  
Move the motor of the specified distance.

#### Parameters

- **pos** – distance to travel (the sign gives the direction)
- **add** – controller address (1 to 31), defaults to self.address

#### Raises

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**reset**(*add: Optional[int] = None*) → *None*

Resets the controller, equivalent to a power-up. This puts the controller back to NOT REFERENCED state, which is necessary for configuring the controller.

Parameters **add** – controller address (1 to 31)

#### Raises

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**set\_acceleration**(*acc: Union[int, float]*, *add: Optional[int] = None*) → *None*

Leave the configuration state. The configuration parameters are saved to the device's memory.

#### Parameters

- **acc** – acceleration (preset units/s^2), value between 1e-6 and 1e12
- **add** – controller address (1 to 31)

#### Raises

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**set\_motor\_configuration**(*add*: *Optional[int]* = *None*, *config*: *Optional[dict]* = *None*) → *None*

Set the motor configuration. The motor must be in CONFIG state.

**Parameters**

- **add** – controller address (1 to 31)
- **config** – dictionary containing the motor's configuration

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**set\_negative\_software\_limit**(*lim*: *Union[int, float]*, *add*: *Optional[int]* = *None*) → *None*

Set the negative software limit (the maximum position that the motor is allowed to travel to towards the left).

**Parameters**

- **lim** – negative software limit (preset units), value between -1e12 and 0
- **add** – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**set\_positive\_software\_limit**(*lim*: *Union[int, float]*, *add*: *Optional[int]* = *None*) → *None*

Set the positive software limit (the maximum position that the motor is allowed to travel to towards the right).

**Parameters**

- **lim** – positive software limit (preset units), value between 0 and 1e12
- **add** – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**start()**

Opens the communication protocol and applies the config.

**Raises** **SerialCommunicationIOError** – when communication port cannot be opened

**stop()** → *None*

Stop the device. Close the communication protocol.

**stop\_motion**(*add*: *Optional[int]* = *None*) → *None*

Stop a move in progress by decelerating the positioner immediately with the configured acceleration until it stops. If a controller address is provided, stops a move in progress on this controller, else stops the moves on all controllers.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**wait\_until\_motor\_initialized**(*add*: *Optional[int]* = *None*) → *None*

Wait until the motor leaves the HOMING state (at which point it should have arrived to the home position).

**Parameters** *add* – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

```
class NewportSMC100PPConfig(address: int = 1, user_position_offset: Union[int, float] = 23.987,
                             screw_scaling: Union[int, float] = 1, exit_configuration_wait_sec: Union[int,
                             float] = 5, move_wait_sec: Union[int, float] = 1, acceleration: Union[int, float]
                             = 10, backlash_compensation: Union[int, float] = 0, hysteresis_compensation:
                             Union[int, float] = 0.015, micro_step_per_full_step_factor: int = 100,
                             motion_distance_per_full_step: Union[int, float] = 0.01, home_search_type:
                             Union[int, hvl_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch] =
                             HomeSearch.HomeSwitch, jerk_time: Union[int, float] = 0.04,
                             home_search_velocity: Union[int, float] = 4, home_search_timeout: Union[int,
                             float] = 27.5, home_search_polling_interval: Union[int, float] = 1,
                             peak_output_current_limit: Union[int, float] = 0.4, rs485_address: int = 2,
                             negative_software_limit: Union[int, float] = - 23.5, positive_software_limit:
                             Union[int, float] = 25, velocity: Union[int, float] = 4, base_velocity: Union[int,
                             float] = 0, stage_configuration: Union[int,
                             hvl_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig] =
                             EspStageConfig.EnableEspStageCheck)
```

Bases: `object`

Configuration dataclass for the Newport motor controller SMC100PP.

```
class EspStageConfig(value=<no_arg>, names=None, module=None, type=None, start=1,
                     boundary=None)
```

Bases: `aenum.IntEnum`

Different configurations to check or not the motor configuration upon power-up.

```
DisableEspStageCheck = 1
EnableEspStageCheck = 3
UpdateEspStageInfo = 2
```

```
class HomeSearch(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
```

Bases: `aenum.IntEnum`

Different methods for the motor to search its home position during initialization.

```
CurrentPosition = 1
EndOfRunSwitch = 4
EndOfRunSwitch_and_Index = 3
HomeSwitch = 2
```

```
HomeSwitch_and_Index = 0
acceleration: Union[int, float] = 10
address: int = 1
backlash_compensation: Union[int, float] = 0
base_velocity: Union[int, float] = 0
clean_values()
exit_configuration_wait_sec: Union[int, float] = 5
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
home_search_polling_interval: Union[int, float] = 1
home_search_timeout: Union[int, float] = 27.5
home_search_type: Union[int, hvl_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch] = 2
home_search_velocity: Union[int, float] = 4
hysteresis_compensation: Union[int, float] = 0.015
is_configdataclass = True
jerk_time: Union[int, float] = 0.04
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
micro_step_per_full_step_factor: int = 100
motion_distance_per_full_step: Union[int, float] = 0.01
property motor_config: Dict[str, float]
```

Gather the configuration parameters of the motor into a dictionary.

**Returns** dict containing the configuration parameters of the motor

```
move_wait_sec: Union[int, float] = 1
negative_software_limit: Union[int, float] = -23.5
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
peak_output_current_limit: Union[int, float] = 0.4
positive_software_limit: Union[int, float] = 25
```

```
post_force_value(fieldname, value)
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

rs485_address: int = 2
screw_scaling: Union[int, float] = 1
stage_configuration: Union[int,
    hvl_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig] = 3
user_position_offset: Union[int, float] = 23.987
velocity: Union[int, float] = 4

class NewportSMC100PPSerialCommunication(configuration)
Bases: hvl_ccb.comm.serial.SerialCommunication

Specific communication protocol implementation for NewportSMC100 controller. Already predefines device-specific protocol parameters in config.

class ControllerErrors(value=<no_arg>, names=None, module=None, type=None, start=1,
    boundary=None)
Bases: aenum.Enum

Possible controller errors with values as returned by the device in response to sent commands.

ADDR_INCORRECT = 'B'
CMD_EXEC_ERROR = 'V'
CMD_NOT_ALLOWED = 'D'
CMD_NOT_ALLOWED_CC = 'X'
CMD_NOT_ALLOWED_CONFIGURATION = 'I'
CMD_NOT_ALLOWED_DISABLE = 'J'
CMD_NOT_ALLOWED_HOMING = 'L'
CMD_NOT_ALLOWED_MOVING = 'M'
CMD_NOT_ALLOWED_NOT_REFERENCED = 'H'
CMD_NOT_ALLOWED_PP = 'W'
CMD_NOT_ALLOWED_READY = 'K'
CODE_OR_ADDR_INVALID = 'A'
COM_TIMEOUT = 'S'
DISPLACEMENT_OUT_OF_LIMIT = 'G'
EEPROM_ACCESS_ERROR = 'U'
ESP_STAGE_NAME_INVALID = 'F'
HOME_STARTED = 'E'
NO_ERROR = '@'
PARAM_MISSING_OR_INVALID = 'C'
```

```
POSITION_OUT_OF_LIMIT = 'N'
```

**check\_for\_error**(*add: int*) → *None*

Ask the Newport controller for the last error it recorded.

This method is called after every command or query.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query**(*add: int, cmd: str, param: Optional[Union[int, float, str]] = None*) → *str*

Send a query to the controller, read the answer, and check for errors. The prefix add+cmd is removed from the answer.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **param** – optional parameter (int/float/str) appended to the command

**Returns** the answer from the device without the prefix

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**query\_multiple**(*add: int, cmd: str, prefixes: List[str]*) → *List[str]*

Send a query to the controller, read the answers, and check for errors. The prefixes are removed from the answers.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **prefixes** – prefixes of each line expected in the answer

**Returns** list of answers from the device without prefix

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**read\_text()** → str

Read one line of text from the serial port, and check for presence of a null char which indicates that the motor power supply was cut and then restored. The input buffer may hold additional data afterwards, since only one line is read.

This method uses *self.access\_lock* to ensure thread-safety.

**Returns** String read from the serial port; ‘’ if there was nothing to read.

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **NewportMotorPowerSupplyWasCutError** – if a null char is read

**send\_command**(*add*: int, *cmd*: str, *param*: Optional[Union[int, float, str]] = None) → None

Send a command to the controller, and check for errors.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **param** – optional parameter (int/float/str) appended to the command

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**send\_stop**(*add*: int) → None

Send the general stop ST command to the controller, and check for errors.

**Parameters** *add* – the controller address (1 to 31)

**Returns** ControllerErrors reported by Newport Controller

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained

**class NewportSMC100PPSerialCommunicationConfig**(*terminator*: bytes = b'\r\n', *encoding*: str = 'ascii', *encoding\_error\_handling*: str = 'replace', *wait\_sec\_read\_text\_nonempty*: Union[int, float] = 0.5, *default\_n\_attempts\_read\_text\_nonempty*: int = 10, *port*: Union[str, NoneType] = None, *baudrate*: int = 57600, *parity*: Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = <SerialCommunicationParity.NONE: 'N'>, *stopbits*: Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = <SerialCommunicationStopbits.ONE: 1>, *bytesize*: Union[int, hvl\_ccb.comm.serial.SerialCommunicationBytesize] = <SerialCommunicationBytesize.EIGHTBITS: 8>, *timeout*: Union[int, float] = 10)

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

**baudrate: int = 57600**  
Baudrate for NewportSMC100 controller is 57600 baud

**bytesize: Union[int, hvl\_ccb.comm.serial.SerialCommunicationBytesize] = 8**  
NewportSMC100 controller uses 8 bits for one data byte

**encoding: str = 'ascii'**  
use ASCII as de-/encoding, cf. the manual

**encoding\_error\_handling: str = 'replace'**  
replace bytes with instead of raising utf-8 exception when decoding fails

**force\_value(fieldname, value)**  
Forces a value to a dataclass field despite the class being frozen.  
NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys() → Sequence[str]**  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity: Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = 'N'**  
NewportSMC100 controller does not use parity

**classmethod required\_keys() → Sequence[str]**  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits: Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = 1**  
NewportSMC100 controller uses one stop bit

**terminator: bytes = b'\r\n'**  
The terminator is CR/LF

**timeout: Union[int, float] = 10**  
use 10 seconds timeout as default

**exception NewportSerialCommunicationError**  
Bases: *hvl\_ccb.dev.newport.NewportError*  
  
Communication error with the Newport controller.

**class NewportStates(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**  
Bases: *hvl\_ccb.utils.enum.AutoNumberNameEnum*  
  
States of the Newport controller. Certain commands are allowed only in certain states.

**CONFIG = 3**

**DISABLE = 6**


```

HOMING = 2
JOGGING = 7
MOVING = 5
NO_REF = 1
READY = 4

exception NewportUncertainPositionError
    Bases: hvl_ccb.dev.newport.NewportError
        Error with the position of the Newport motor.

```

## **[hvl\\_ccb.dev.pfeiffer\\_tpg](#)**



Device class for Pfeiffer TPG controllers.

The Pfeiffer TPG control units are used to control Pfeiffer Compact Gauges. Models: TPG 251 A, TPG 252 A, TPG 256A, TPG 261, TPG 262, TPG 361, TPG 362 and TPG 366.

Manufacturer homepage: <https://www.pfeiffer-vacuum.com/en/products/measurement-analysis/> measurement/activeline/controllers/

```

class PfeifferTPG(com, dev_config=None)
    Bases: hvl_ccb.dev.base.SingleCommDevice

```

Pfeiffer TPG control unit device class

```

class PressureUnits(value=<no_arg>, names=None, module=None, type=None, start=1,
                     boundary=None)
    Bases: hvl_ccb.utils.enum.NameEnum

```

Enum of available pressure units for the digital display. “0” corresponds either to bar or to mbar depending on the TPG model. In case of doubt, the unit is visible on the digital display.

**Micron** = 3

**Pascal** = 2

**Torr** = 1

```
Volt = 5
bar = 0
hPascal = 4
mbar = 0

class SensorStatus(value)
    Bases: enum.IntEnum
    An enumeration.

    Identification_error = 6
        No_sensor = 5
        Ok = 0
        Overrange = 2
        Sensor_error = 3
        Sensor_off = 4
        Underrange = 1

class SensorTypes(value)
    Bases: enum.Enum
    An enumeration.

    CMR = 4
    IKR = 2
    IKR11 = 2
    IKR9 = 2
    IMR = 5
    None = 7
    PBR = 6
    PKR = 3
    TPR = 1
    noSENSOR = 7
    noSen = 7

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

get_full_scale_mbar() → List[Union[int, float]]
    Get the full scale range of the attached sensors

    Returns full scale range values in mbar, like [0.01, 1, 0.1, 1000, 50000, 10]

    Raises
```

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**get\_full\_scale\_unitless()** → List[int]

Get the full scale range of the attached sensors. See lookup table between command and corresponding pressure in the device user manual.

**Returns** list of full scale range values, like [0, 1, 3, 3, 2, 0]

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**identify\_sensors()** → *None*

Send identification request TID to sensors on all channels.

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**measure(channel: int)** → Tuple[str, float]

Get the status and measurement of one sensor

**Parameters** **channel** – int channel on which the sensor is connected, with 1 <= channel <= number\_of\_sensors

**Returns** measured value as float if measurement successful, sensor status as string if not

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**measure\_all()** → List[Tuple[str, float]]

Get the status and measurement of all sensors (this command is not available on all models)

**Returns** list of measured values as float if measurements successful, and or sensor status as strings if not

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**property number\_of\_sensors**

**set\_display\_unit(unit: Union[str, hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits])** → *None*

Set the unit in which the measurements are shown on the display.

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**set\_full\_scale\_mbar(fsr: List[Union[int, float]])** → *None*

Set the full scale range of the attached sensors (in unit mbar)

**Parameters** **fsr** – full scale range values in mbar, for example [0.01, 1000]

**Raises**

- ***SerialCommunicationIOError*** – when communication port is not opened
- ***PfeifferTPGError*** – if command fails

**set\_full\_scale\_unitless**(*fsr: List[int]*) → *None*

Set the full scale range of the attached sensors. See lookup table between command and corresponding pressure in the device user manual.

**Parameters** **fsr** – list of full scale range values, like [0, 1, 3, 3, 2, 0]

**Raises**

- ***SerialCommunicationIOError*** – when communication port is not opened
- ***PfeifferTPGError*** – if command fails

**start()** → *None*

Start this device. Opens the communication protocol, and identify the sensors.

**Raises** ***SerialCommunicationIOError*** – when communication port cannot be opened

**stop()** → *None*

Stop the device. Closes also the communication protocol.

**property unit**

The pressure unit of readings is always mbar, regardless of the display unit.

**class PfeifferTPGConfig**(*model: Union[str, hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig.Model] = Model.TPG25xA*)

Bases: *object*

Device configuration dataclass for Pfeiffer TPG controllers.

**class Model**(*value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None*)

Bases: *hvl\_ccb.utils.enum.NameEnum*

An enumeration.

**TPG25xA** = {0.1: 8, 1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6, 50000: 7}

**TPGx6x** = {0.01: 0, 0.1: 1, 1: 2, 10: 3, 100: 4, 1000: 5, 2000: 6, 5000: 7, 10000: 8, 50000: 9}

**is\_valid\_scale\_range\_reversed\_str**(*v: str*) → *bool*

Check if given string represents a valid reversed scale range of a model.

**Parameters** **v** – Reversed scale range string.

**Returns** *True* if valid, *False* otherwise.

**clean\_values()**

**force\_value**(*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass** = *True*

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**model:** Union[str, [hvl\\_ccb.dev.pfeiffer\\_tpg.PfeifferTPGConfig.Model](#)] = {0.1: 8, 1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6, 50000: 7}

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception PfeifferTPGError**

Bases: [hvl\\_ccb.dev.base.DeviceException](#)

Error with the Pfeiffer TPG Controller.

**class PfeifferTPGSerialCommunication(configuration)**

Bases: [hvl\\_ccb.comm.serial.SerialCommunication](#)

Specific communication protocol implementation for Pfeiffer TPG controllers. Already predefines device-specific protocol parameters in config.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query(cmd: str)** → str

Send a query, then read and returns the first line from the com port.

**Parameters** cmd – query message to send to the device

**Returns** first line read on the com

**Raises**

- [SerialCommunicationIOError](#) – when communication port is not opened
- [PfeifferTPGError](#) – if the device does not acknowledge the command or if the answer from the device is empty

**send\_command(cmd: str)** → None

Send a command to the device and check for acknowledgement.

**Parameters** cmd – command to send to the device

**Raises**

- [SerialCommunicationIOError](#) – when communication port is not opened
- [PfeifferTPGError](#) – if the answer from the device differs from the expected acknowledgement character ‘chr(6)’.

```
class PfeifferTPGSerialCommunicationConfig(terminator: bytes = b'\r\n', encoding: str = 'utf-8',
                                             encoding_error_handling: str = 'strict',
                                             wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                             default_n_attempts_read_text_nonempty: int = 10, port:
                                             Union[str, NoneType] = None, baudrate: int = 9600, parity:
                                             Union[str,
                                               hvl_ccb.comm.serial.SerialCommunicationParity] =
                                               <SerialCommunicationParity.NONE: 'N'>, stopbits:
                                               Union[int,
                                                 hvl_ccb.comm.serial.SerialCommunicationStopbits] =
                                                 <SerialCommunicationStopbits.ONE: 1>, bytesize:
                                               Union[int,
                                                 hvl_ccb.comm.serial.SerialCommunicationBytesize] =
                                                 <SerialCommunicationBytesize.EIGHTBITS: 8>, timeout:
                                               Union[int, float] = 3)
```

Bases: `hvl_ccb.comm.serial.SerialCommunicationConfig`

**baudrate: int = 9600**

Baudrate for Pfeiffer TPG controllers is 9600 baud

**bytesize: Union[int, hvl\_ccb.comm.serial.SerialCommunicationBytesize] = 8**

One byte is eight bits long

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity: Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = 'N'**

Pfeiffer TPG controllers do not use parity

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits: Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = 1**

Pfeiffer TPG controllers use one stop bit

**terminator: bytes = b'\r\n'**

The terminator is <CR><LF>

```
timeout: Union[int, float] = 3
use 3 seconds timeout as default
```

## **hvl\_ccb.dev.picotech\_pt104**


**NOTE:** PicoSDK Python wrappers already on import attempt to load the PicoSDK library; thus, the API docs can only be generated in a system with the latter installed and are by default disabled.

To build the API docs for this submodule locally edit the `docs/hvl_ccb.dev.picotech_pt104.rst` file to remove the `.. code-block::` directive preceding the following directives:

```
.. inheritance-diagram:: hvl_ccb.dev.picotech_pt104
:parts: 1

.. automodule:: hvl_ccb.dev.picotech_pt104
:members:
:undoc-members:
:show-inheritance:
```

## **hvl\_ccb.dev.rs\_rto1024**



Python module for the Rhode & Schwarz RTO 1024 oscilloscope. The communication to the device is through VISA, type TCPIP / INSTR.

```
class RT01024(com: Union[hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunication,
hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig, dict], dev_config:
Union[hvl_ccb.dev.rs_rto1024.RTO1024Config, dict])
```

Bases: `hvl_ccb.dev.visa.VisaDevice`

Device class for the Rhode & Schwarz RTO 1024 oscilloscope.

```
class TriggerModes(value=<no_arg>, names=None, module=None, type=None, start=1,
boundary=None)
```

Bases: `hvl_ccb.utils.enum.AutoNumberNameEnum`

Enumeration for the three available trigger modes.

```
AUTO = 1  
FREE RUN = 3  
NORMAL = 2
```

#### **classmethod names()**

Returns a list of the available trigger modes. :return: list of strings

#### **activate\_measurements(meas\_n: int, source: str, measurements: List[str], category: str = 'AMPTime')**

Activate the list of ‘measurements’ of the waveform ‘source’ in the measurement box number ‘meas\_n’. The list ‘measurements’ starts with the main measurement and continues with additional measurements of the same ‘category’.

##### **Parameters**

- **meas\_n** – measurement number 1..8
- **source** – measurement source, for example C1W1
- **measurements** – list of measurements, the first one will be the main measurement.
- **category** – the category of measurements, by default AMPTime

#### **backup\_waveform(filename: str) → None**

Backup a waveform file from the standard directory specified in the device configuration to the standard backup destination specified in the device configuration. The filename has to be specified without .bin or path.

**Parameters filename** – The waveform filename without extension and path

#### **static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

#### **static default\_com\_cls()**

Return the default communication protocol for this device type, which is VisaCommunication.

**Returns** the VisaCommunication class

#### **file\_copy(source: str, destination: str) → None**

Copy a file from one destination to another on the oscilloscope drive. If the destination file already exists, it is overwritten without notice.

##### **Parameters**

- **source** – absolute path to the source file on the DSO filesystem
- **destination** – absolute path to the destination file on the DSO filesystem

**Raises RT01024Error** – if the operation did not complete

#### **get\_acquire\_length() → float**

Gets the time of one acquisition, that is the time across the 10 divisions of the diagram.

- Range: 250E-12 ... 500 [s]
- Increment: 1E-12 [s]

**Returns** the time for one acquisition. Range: 250e-12 ... 500 [s]

#### **get\_channel\_offset(channel: int) → float**

Gets the voltage offset of the indicated channel.

**Parameters** `channel` – is the channel number (1..4)

**Returns** channel offset voltage in V (value between -1 and 1)

**get\_channel\_position**(`channel: int`) → float

Gets the vertical position of the indicated channel.

**Parameters** `channel` – is the channel number (1..4)

**Returns** channel position in div (value between -5 and 5)

**get\_channel\_range**(`channel: int`) → float

Queries the channel range in V.

**Parameters** `channel` – is the input channel (1..4)

**Returns** channel range in V

**get\_channel\_scale**(`channel: int`) → float

Queries the channel scale in V/div.

**Parameters** `channel` – is the input channel (1..4)

**Returns** channel scale in V/div

**get\_channel\_state**(`channel: int`) → bool

Queries if the channel is active or not.

**Parameters** `channel` – is the input channel (1..4)

**Returns** True if active, else False

**get\_reference\_point**() → int

Gets the reference point of the time scale in % of the display. If the “Trigger offset” is zero, the trigger point matches the reference point. ReferencePoint = zero pint of the time scale

- Range: 0 … 100 [%]
- Increment: 1 [%]

**Returns** the reference in %

**get\_repetitions**() → int

Get the number of acquired waveforms with RUN Nx SINGLE. Also defines the number of waveforms used to calculate the average waveform.

- Range: 1 … 16777215
- Increment: 10
- \*RST = 1

**Returns** the number of waveforms to acquire

**get\_timestamps**() → List[float]

Gets the timestamps of all recorded frames in the history and returns them as a list of floats.

**Returns** list of timestamps in [s]

**Raises** `RT01024Error` – if the timestamps are invalid

**list\_directory**(`path: str`) → List[Tuple[str, str, int]]

List the contents of a given directory on the oscilloscope filesystem.

**Parameters** `path` – is the path to a folder

**Returns** a list of filenames in the given folder

**load\_configuration**(filename: str) → *None*

Load current settings from a configuration file. The filename has to be specified without base directory and ‘.dfl’ extension.

**Information from the manual** *ReCaLL* calls up the instrument settings from an intermediate memory identified by the specified number. The instrument settings can be stored to this memory using the command \*SAV with the associated number. It also activates the instrument settings which are stored in a file and loaded using *MMEMory:LOAD:STATE* .

**Parameters** **filename** – is the name of the settings file without path and extension

**local\_display**(state: bool) → *None*

Enable or disable local display of the scope.

**Parameters** **state** – is the desired local display state

**prepare\_ultra\_segmentation**() → *None*

Make ready for a new acquisition in ultra segmentation mode. This function does one acquisition without ultra segmentation to clear the history and prepare for a new measurement.

**read\_measurement**(meas\_n: int, name: str) → float

**Parameters**

- **meas\_n** – measurement number 1..8
- **name** – measurement name, for example “MAX”

**Returns** measured value

**run\_continuous\_acquisition**() → *None*

Start acquiring continuously.

**run\_single\_acquisition**() → *None*

Start a single or Nx acquisition.

**save\_configuration**(filename: str) → *None*

Save the current oscilloscope settings to a file. The filename has to be specified without path and ‘.dfl’ extension, the file will be saved to the configured settings directory.

**Information from the manual** *SAVe* stores the current instrument settings under the specified number in an intermediate memory. The settings can be recalled using the command \*RCL with the associated number. To transfer the stored instrument settings to a file, use *MMEMory:STORe:STATE* .

**Parameters** **filename** – is the name of the settings file without path and extension

**save\_waveform\_history**(filename: str, channel: int, waveform: int = 1) → *None*

Save the history of one channel and one waveform to a .bin file. This function is used after an acquisition using sequence trigger mode (with or without ultra segmentation) was performed.

**Parameters**

- **filename** – is the name (without extension) of the file
- **channel** – is the channel number
- **waveform** – is the waveform number (typically 1)

**Raises** *RT01024Error* – if storing waveform times out

**set\_acquire\_length**(timerange: float) → *None*

Defines the time of one acquisition, that is the time across the 10 divisions of the diagram.

- Range: 250E-12 ... 500 [s]
- Increment: 1E-12 [s]
- \*RST = 0.5 [s]

**Parameters** **timerange** – is the time for one acquisition. Range: 250e-12 ... 500 [s]

**set\_channel\_offset**(channel: int, offset: float) → *None*

Sets the voltage offset of the indicated channel.

- Range: Dependent on the channel scale and coupling [V]
- Increment: Minimum 0.001 [V], may be higher depending on the channel scale and coupling
- \*RST = 0

**Parameters**

- **channel** – is the channel number (1..4)
- **offset** – Offset voltage. Positive values move the waveform down, negative values move it up.

**set\_channel\_position**(channel: int, position: float) → *None*

Sets the vertical position of the indicated channel as a graphical value.

- Range: -5.0 ... 5.0 [div]
- Increment: 0.02
- \*RST = 0

**Parameters**

- **channel** – is the channel number (1..4)
- **position** – is the position. Positive values move the waveform up, negative values move it down.

**set\_channel\_range**(channel: int, v\_range: float) → *None*

Sets the voltage range across the 10 vertical divisions of the diagram. Use the command alternatively instead of set\_channel\_scale.

- Range for range: Depends on attenuation factors and coupling. With 1:1 probe and external attenuations and 50 input coupling, the range is 10 mV to 10 V. For 1 M input coupling, it is 10 mV to 100 V. If the probe and/or external attenuation is changed, multiply the range values by the attenuation factors.
- Increment: 0.01
- \*RST = 0.5

**Parameters**

- **channel** – is the channel number (1..4)
- **v\_range** – is the vertical range [V]

**set\_channel\_scale**(channel: int, scale: float) → *None*

Sets the vertical scale for the indicated channel. The scale value is given in volts per division.

- Range for scale: depends on attenuation factor and coupling. With 1:1 probe and external attenuations and 50 input coupling, the vertical scale (input sensitivity) is 1 mV/div to 1 V/div. For 1 M input coupling, it is 1 mV/div to 10 V/div. If the probe and/or external attenuation is changed, multiply the values by the attenuation factors to get the actual scale range.
- Increment: 1e-3
- \*RST = 0.05

See also: `set_channel_range`

#### Parameters

- **channel** – is the channel number (1..4)
- **scale** – is the vertical scaling [V/div]

**set\_channel\_state**(*channel: int, state: bool*) → *None*

Switches the channel signal on or off.

#### Parameters

- **channel** – is the input channel (1..4)
- **state** – is True for on, False for off

**set\_reference\_point**(*percentage: int*) → *None*

Sets the reference point of the time scale in % of the display. If the “Trigger offset” is zero, the trigger point matches the reference point. ReferencePoint = zero pint of the time scale

- Range: 0 … 100 [%]
- Increment: 1 [%]
- \*RST = 50 [%]

Parameters **percentage** – is the reference in %

**set\_repetitions**(*number: int*) → *None*

Set the number of acquired waveforms with RUN Nx SINGLE. Also defines the number of waveforms used to calculate the average waveform.

- Range: 1 … 16777215
- Increment: 10
- \*RST = 1

Parameters **number** – is the number of waveforms to acquire

**set\_trigger\_level**(*channel: int, level: float, event\_type: int = 1*) → *None*

Sets the trigger level for the specified event and source.

- Range: -10 to 10 V
- Increment: 1e-3 V
- \*RST = 0 V

#### Parameters

- **channel** – indicates the trigger source.
  - 1..4 = channel 1 to 4, available for all event types 1..3

- 5 = external trigger input on the rear panel for analog signals, available for A-event type  
= 1
- 6..9 = not available
- **level** – is the voltage for the trigger level in [V].
- **event\_type** – is the event type. 1: A-Event, 2: B-Event, 3: R-Event

**set\_trigger\_mode**(*mode*: Union[str, hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes]) → *None*  
Sets the trigger mode which determines the behavior of the instrument if no trigger occurs.

**Parameters** **mode** – is either auto, normal, or freerun.

**Raises** **RT01024Error** – if an invalid triggermode is selected

**set\_trigger\_source**(*channel*: int, *event\_type*: int = 1) → *None*  
Set the trigger (Event A) source channel.

**Parameters**

- **channel** – is the channel number (1..4)
- **event\_type** – is the event type. 1: A-Event, 2: B-Event, 3: R-Event

**start()** → *None*

Start the RTO1024 oscilloscope and bring it into a defined state and remote mode.

**stop()** → *None*

Stop the RTO1024 oscilloscope, reset events and close communication. Brings back the device to a state where local operation is possible.

**stop\_acquisition()** → *None*  
Stop any acquisition.

**class RT01024Config**(*waveforms\_path*: str, *settings\_path*: str, *backup\_path*: str, *spoll\_interval*: Union[int, float] = 0.5, *spoll\_start\_delay*: Union[int, float] = 2, *command\_timeout\_seconds*: Union[int, float] = 60, *wait\_sec\_short\_pause*: Union[int, float] = 0.1, *wait\_sec\_enable\_history*: Union[int, float] = 1, *wait\_sec\_post\_acquisition\_start*: Union[int, float] = 2)  
Bases: *hvl\_ccb.dev.visa.VisaDeviceConfig*, *hvl\_ccb.dev.rs\_rto1024.\_RT01024ConfigDefaultsBase*, *hvl\_ccb.dev.rs\_rto1024.\_RT01024ConfigBase*

Configdataclass for the RTO1024 device.

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception RT01024Error**

Bases: *hvl\_ccb.dev.base.DeviceException*

**class RT01024VisaCommunication(configuration)**

Bases: *hvl\_ccb.comm.visa.VisaCommunication*

Specialization of VisaCommunication for the RTO1024 oscilloscope

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class RT01024VisaCommunicationConfig(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address], interface\_type: Union[str, *hvl\_ccb.comm.visa.VisaCommunicationConfig.InterfaceType*] = InterfaceType.TCPIP\_INSTR, board: int = 0, port: int = 5025, timeout: int = 5000, chunk\_size: int = 204800, open\_timeout: int = 1000, write\_termination: str = '\n', read\_termination: str = '\n', visa\_backend: str = '')**

Bases: *hvl\_ccb.comm.visa.VisaCommunicationConfig*

Configuration dataclass for VisaCommunication with specifications for the RTO1024 device class.

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**interface\_type: Union[str, *hvl\_ccb.comm.visa.VisaCommunicationConfig.InterfaceType*] = 2**

Interface type of the VISA connection, being one of *InterfaceType*.

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.


**Returns** a list of strings containing all optional keys.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## [hvl\\_ccb.dev.se\\_il2t](#)



Device class for controlling a Schneider Electric ILS2T stepper drive over modbus TCP.

**class ILS2T(*com, dev\_config=None*)**  
Bases: [hvl\\_ccb.dev.base.SingleCommDevice](#)

Schneider Electric ILS2T stepper drive class.

**ACTION\_JOG\_VALUE = 0**

The single action value for *ILS2T.Mode.JOG*

**class ActionsPtp(*value*)**  
Bases: [enum.IntEnum](#)

Allowed actions in the point to point mode (*ILS2T.Mode.PTP*).

**ABSOLUTE\_POSITION = 0**

**RELATIVE\_POSITION\_MOTOR = 2**

**RELATIVE\_POSITION\_TARGET = 1**

**DEFAULT\_IO\_SCANNING\_CONTROL\_VALUES = {*'action': 2, 'continue\_after\_stop\_cu': 0, 'disable\_driver\_di': 0, 'enable\_driver\_en': 0, 'execute\_stop\_sh': 0, 'fault\_reset\_fr': 0, 'mode': 3, 'quick\_stop\_qs': 0, 'ref\_16': 1500, 'ref\_32': 0, 'reset\_stop\_ch': 0*}**

Default IO Scanning control mode values

**class Mode(*value*)**  
Bases: [enum.IntEnum](#)

ILS2T device modes

**JOG = 1**

**PTP = 3**

**class Ref16Jog(*value*)**  
Bases: [enum.Flag](#)

Allowed values for ILS2T ref\_16 register (the shown values are the integer representation of the bits), all in Jog mode = 1

**FAST = 4**

**NEG = 2**

**NEG\_FAST = 6**

**NONE = 0**

```

POS = 1
POS_FAST = 5

RegAddr
    Modbus Register Adresses
    alias of hvl\_ccb.dev.se\_ilst2t.ILS2TRegAddr

RegDatatype
    Modbus Register Datatypes
    alias of hvl\_ccb.dev.se\_ilst2t.ILS2TRegDatatype

class State(value)
    Bases: enum.IntEnum
    State machine status values
    ON = 6
    QUICKSTOP = 7
    READY = 4

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

disable(log_warn: bool = True, wait_sec_max: Optional[int] = None) → bool
    Disable the driver of the stepper motor and enable the brake.

    Note: the driver cannot be disabled if the motor is still running.

Parameters

- log_warn – if log a warning in case the motor cannot be disabled.
- wait_sec_max – maximal wait time for the motor to stop running and to disable it; by default, with None, use a config value

Returns True if disable request could and was sent, False otherwise.

do_iostressing_write(**kwargs: int) → None
    Perform a write operation using IO Scanning mode.

Parameters kwargs – Keyword-argument list with options to send, remaining are taken from the defaults.

enable() → None
    Enable the driver of the stepper motor and disable the brake.

execute_absolute_position(position: int) → bool
    Execute a absolute position change, i.e. enable motor, perform absolute position change, wait until done and disable motor afterwards.

    Check position at the end if wrong do not raise error; instead just log and return check result.

    Parameters position – absolute position of motor in user defined steps.

    Returns True if actual position is as expected, False otherwise.

```

**execute\_relative\_step**(*steps: int*) → bool

Execute a relative step, i.e. enable motor, perform relative steps, wait until done and disable motor afterwards.

Check position at the end if wrong do not raise error; instead just log and return check result.

**Parameters** **steps** – Number of steps.

**Returns** *True* if actual position is as expected, *False* otherwise.

**get\_dc\_volt**() → float

Read the DC supply voltage of the motor.

**Returns** DC input voltage.

**get\_error\_code**() → Dict[int, Dict[str, Any]]

Read all messages in fault memory. Will read the full error message and return the decoded values. At the end the fault memory of the motor will be deleted. In addition, *reset\_error* is called to re-enable the motor for operation.

**Returns** Dictionary with all information

**get\_position**() → int

Read the position of the drive and store into status.

**Returns** Position step value

**get\_status**() → Dict[str, int]

Perform an IO Scanning read and return the status of the motor.

**Returns** dict with status information.

**get\_temperature**() → int

Read the temperature of the motor.

**Returns** Temperature in degrees Celsius.

**jog\_run**(*direction: bool = True, fast: bool = False*) → *None*

Slowly turn the motor in positive direction.

**jog\_stop**() → *None*

Stop turning the motor in Jog mode.

**quickstop**() → *None*

Stops the motor with high deceleration rate and falls into error state. Reset with *reset\_error* to recover into normal state.

**reset\_error**() → *None*

Resets the motor into normal state after quick stop or another error occurred.

**set\_jog\_speed**(*slow: int = 60, fast: int = 180*) → *None*

Set the speed for jog mode. Default values correspond to startup values of the motor.

**Parameters**

- **slow** – RPM for slow jog mode.
- **fast** – RPM for fast jog mode.

**set\_max\_acceleration**(*rpm\_minute: int*) → *None*

Set the maximum acceleration of the motor.

**Parameters** **rpm\_minute** – revolution per minute per minute

**set\_max\_deceleration**(*rpm\_minute: int*) → *None*

Set the maximum deceleration of the motor.

**Parameters** `rpm_minute` – revolution per minute per minute

`set_max_rpm(rpm: int) → None`

Set the maximum RPM.

**Parameters** `rpm` – revolution per minute (  $0 < rpm \leq RPM\_MAX$  )

**Raises** `ILS2TException` – if RPM is out of range

`set_ramp_type(ramp_type: int = -1) → None`

**Set the ramp type. There are two options available:** 0: linear ramp -1: motor optimized ramp

**Parameters** `ramp_type` – 0: linear ramp | -1: motor optimized ramp

`start() → None`

Start this device.

`stop() → None`

Stop this device. Disables the motor (applies brake), disables access and closes the communication protocol.

`user_steps(steps: int = 16384, revolutions: int = 1) → None`

Define steps per revolution. Default is 16384 steps per revolution. Maximum precision is 32768 steps per revolution.

**Parameters**

- **steps** – number of steps in *revolutions*.
- **revolutions** – number of revolutions corresponding to *steps*.

`write_absolute_position(position: int) → None`

Write instruction to turn the motor until it reaches the absolute position. This function does not enable or disable the motor automatically.

**Parameters** `position` – absolute position of motor in user defined steps.

`write_relative_step(steps: int) → None`

Write instruction to turn the motor the relative amount of steps. This function does not enable or disable the motor automatically.

**Parameters** `steps` – Number of steps to turn the motor.

```
class ILS2TConfig(rpm_max_init: numbers.Integral = 1500, wait_sec_post_enable: Union[int, float] = 1,
                    wait_sec_max_disable: Union[int, float] = 10, wait_sec_post_cannot_disable: Union[int, float] = 1,
                    wait_sec_post_relative_step: Union[int, float] = 2,
                    wait_sec_post_absolute_position: Union[int, float] = 2)
```

Bases: `object`

Configuration for the ILS2T stepper motor device.

`clean_values()`

`force_value(fieldname, value)`

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field

- **value** – value to assign

**is\_configdataclass** = True

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**rpm\_max\_init: numbers.Integral = 1500**

initial maximum RPM for the motor, can be set up to 3000 RPM. The user is allowed to set a new max RPM at runtime using [ILS2T.set\\_max\\_rpm\(\)](#), but the value must never exceed this configuration setting.

**wait\_sec\_max\_disable: Union[int, float] = 10**

**wait\_sec\_post\_absolute\_position: Union[int, float] = 2**

**wait\_sec\_post\_CANNOT\_DISABLE: Union[int, float] = 1**

**wait\_sec\_post\_enable: Union[int, float] = 1**

**wait\_sec\_post\_relative\_step: Union[int, float] = 2**

**exception ILS2TException**

Bases: [hvl\\_ccb.dev.base.DeviceException](#)

Exception to indicate problems with the SE ILS2T stepper motor.

**class ILS2TModbusTcpCommunication(configuration)**

Bases: [hvl\\_ccb.comm.modbus\\_tcp.ModbusTcpCommunication](#)

Specific implementation of Modbus/TCP for the Schneider Electric ILS2T stepper motor.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class ILS2TModbusTcpCommunicationConfig(host: Union[str, ipaddress.IPv4Address,**

**ipaddress.IPv6Address], unit: int = 255, port: int = 502)**

Bases: [hvl\\_ccb.comm.modbus\\_tcp.ModbusTcpCommunicationConfig](#)

Configuration dataclass for Modbus/TCP communication specific for the Schneider Electric ILS2T stepper motor.

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field

- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**unit: int = 255**

The unit has to be 255 such that IO scanning mode works.

**class ILS2TRegAddr(value)**

Bases: enum.IntEnum

Modbus Register Adresses for Schneider Electric ILS2T stepper drive.

**ACCESS\_ENABLE = 282**

**FLT\_INFO = 15362**

**FLT\_MEM\_DEL = 15112**

**FLT\_MEM\_RESET = 15114**

**IO\_SCANNING = 6922**

**JOGN\_FAST = 10506**

**JOGN\_SLOW = 10504**

**POSITION = 7706**

**RAMP\_ACC = 1556**

**RAMP\_DECEL = 1558**

**RAMP\_N\_MAX = 1554**

**RAMP\_TYPE = 1574**

**SCALE = 1550**

**TEMP = 7200**

**VOLT = 7198**

**class ILS2TRegDatatype(value=<no\_arg>, names=None, module=None, type=None, start=1, boundary=None)**

Bases: aenum.Enum

Modbus Register Datatypes for Schneider Electric ILS2T stepper drive.

From the manual of the drive:

datatype	byte	min	max
INT8	1 Byte	-128	127
UINT8	1 Byte	0	255
INT16	2 Byte	-32_768	32_767
UINT16	2 Byte	0	65_535
INT32	4 Byte	-2_147_483_648	2_147_483_647
UINT32	4 Byte	0	4_294_967_295
BITS	just 32bits	N/A	N/A

**INT32 = (-2147483648, 2147483647)**

**is\_in\_range(value: int) → bool**

**exception IoScanningModeValueError**

Bases: *hvl\_ccb.dev.se\_ilst2t.ILS2TException*


Exception to indicate that the selected IO scanning mode is invalid.

**exception ScalingFactorValueError**

Bases: *hvl\_ccb.dev.se\_ilst2t.ILS2TException*

Exception to indicate that a scaling factor value is invalid.

## **hvl\_ccb.dev.sst\_luminox**



Device class for a SST Luminox Oxygen sensor. This device can measure the oxygen concentration between 0 % and 25 %.

Furthermore, it measures the barometric pressure and internal temperature. The device supports two operating modes: in streaming mode the device measures all parameters every second, in polling mode the device measures only after a query.

Technical specification and documentation for the device can be found at the manufacturer's page: <https://www.sstsensing.com/product/luminox-optical-oxygen-sensors-2/>

**class Luminox(com, dev\_config=None)**

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Luminox oxygen sensor device class.

**activate\_output(mode: hvl\_ccb.dev.sst\_luminox.LuminoxOutputMode) → None**

activate the selected output mode of the Luminox Sensor. :param mode: polling or streaming

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**query\_polling(measurement: Union[str, hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType]) → Union[Dict[Union[str, hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType], Union[float, int, str]], float, int, str]**

Query a value or values of Luminox measurements in the polling mode, according to a given measurement type.

**Parameters** **measurement** – type of measurement

**Returns** value of requested measurement

**Raises**

- **ValueError** – when a wrong key for LuminoxMeasurementType is provided
- **LuminoxOutputModeError** – when polling mode is not activated
- **LuminoxMeasurementTypeError** – when expected measurement value is not read

**read\_streaming() → Dict[Union[str, hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType], Union[float, int, str]]**

Read values of Luminox in the streaming mode. Convert the single string into separate values.

**Returns** dictionary with *LuminoxMeasurementType.all\_measurements\_types()* keys and accordingly type-parsed values.

**Raises**

- **LuminoxOutputModeError** – when streaming mode is not activated
- **LuminoxMeasurementTypeError** – when any of expected measurement values is not read

**start() → None**

Start this device. Opens the communication protocol.

**stop() → None**

Stop the device. Closes also the communication protocol.

**class LuminoxConfig(wait\_sec\_post\_activate: Union[int, float] = 0.5, wait\_sec\_trials\_activate: Union[int, float] = 0.1, nr\_trials\_activate: int = 5)**

Bases: object

Configuration for the SST Luminox oxygen sensor.

**clean\_values()****force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field

- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys()** → Sequence[str]  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**nr\_trials\_activate: int = 5**

**classmethod optional\_defaults()** → Dict[str, object]  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**wait\_sec\_post\_activate: Union[int, float] = 0.5**

**wait\_sec\_trials\_activate: Union[int, float] = 0.1**

**exception LuminoxException**  
Bases: *hvl\_ccb.dev.base.DeviceException*  
General Exception for Luminox Device.

**class LuminoxMeasurementType**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*,  
*boundary=None*)  
Bases: *hvl\_ccb.utils.enum.ValueEnum*  
Measurement types for *LuminoxOutputMode.polling*.

The *all\_measurements* type will read values for the actual measurement types as given in *LuminoxOutputMode.all\_measurements\_types()*; it parses multiple single values using regexp's for other measurement types, therefore, no regexp is defined for this measurement type.

**all\_measurements = 'A'**

**classmethod all\_measurements\_types()** → Tuple[*hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType*, ...]  
A tuple of *LuminoxMeasurementType* enum instances which are actual measurements, i.e. not date of manufacture or software revision.

**barometric\_pressure = 'P'**

**property command: str**

**date\_of\_manufacture = '# 0'**

**parse\_read\_measurement\_value**(*read\_txt: str*) → Union[Dict[Union[str, *hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType*], Union[float, int, str]], float, int, str]

**partial\_pressure\_o2 = '0'**

**percent\_o2 = '%'**

**sensor\_status = 'e'**

**serial\_number = '# 1'**

```

software_revision = '# 2'
temperature_sensor = 'T'

LuminoxMeasurementTypeDict
A typing hint for a dictionary holding LuminoxMeasurementType values. Keys are allowed as strings because
LuminoxMeasurementType is of a StrEnumBase type.

alias of Dict[Union[str, LuminoxMeasurementType], Union[float, int, str]]

exception LuminoxMeasurementTypeError
Bases: hvl_ccb.dev.sst_luminox.LuminoxException

Wrong measurement type for requested data

LuminoxMeasurementTypeValue
A typing hint for all possible LuminoxMeasurementType values as read in either streaming mode or in a polling
mode with LuminoxMeasurementType.all_measurements.

Beware: has to be manually kept in sync with LuminoxMeasurementType instances cast_type attribute values.

alias of Union[float, int, str]

class LuminoxOutputMode(value)
Bases: enum.Enum

output mode.

polling = 1
streaming = 0

exception LuminoxOutputModeError
Bases: hvl_ccb.dev.sst_luminox.LuminoxException

Wrong output mode for requested data

class LuminoxSerialCommunication(configuration)
Bases: hvl_ccb.comm.serial.SerialCommunication

Specific communication protocol implementation for the SST Luminox oxygen sensor. Already predefines
device-specific protocol parameters in config.

static config_cls()
Return the default configdataclass class.

    Returns a reference to the default configdataclass class

class LuminoxSerialCommunicationConfig(terminator: bytes = b'\r\n', encoding: str = 'utf-8',
                                             encoding_error_handling: str = 'strict',
                                             wait_sec_read_text_nonempty: Union[int, float] = 0.5,
                                             default_n_attempts_read_text_nonempty: int = 10, port:
                                             Union[str, NoneType] = None, baudrate: int = 9600, parity:
                                             Union[str, hvl_ccb.comm.serial.SerialCommunicationParity] =
                                             <SerialCommunicationParity.NONE: 'N'>, stopbits: Union[int,
                                             hvl_ccb.comm.serial.SerialCommunicationStopbits] =
                                             <SerialCommunicationStopbits.ONE: 1>, bytesize: Union[int,
                                             hvl_ccb.comm.serial.SerialCommunicationBytesize] =
                                             <SerialCommunicationBytesize.EIGHTBITS: 8>, timeout:
                                             Union[int, float] = 3)

Bases: hvl_ccb.comm.serial.SerialCommunicationConfig

baudrate: int = 9600
Baudrate for SST Luminox is 9600 baud

```

**bytesize:** Union[int, hvl\_ccb.comm.serial.SerialCommunicationBytesize] = 8  
One byte is eight bits long

**force\_value**(fieldname, value)  
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity:** Union[str, hvl\_ccb.comm.serial.SerialCommunicationParity] = 'N'  
SST Luminox does not use parity

**classmethod required\_keys()** → Sequence[str]  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits:** Union[int, hvl\_ccb.comm.serial.SerialCommunicationStopbits] = 1  
SST Luminox does use one stop bit

**terminator:** bytes = b'\r\n'  
The terminator is CR LF

**timeout:** Union[int, float] = 3  
use 3 seconds timeout as default

## hvl\_ccb.dev.utils

### Poller

**class Poller**(spoll\_handler: Callable, polling\_delay\_sec: Union[int, float] = 0, polling\_interval\_sec: Union[int, float] = 1, polling\_timeout\_sec: Optional[Union[int, float]] = None)  
Bases: object

Poller class wrapping *concurrent.futures.ThreadPoolExecutor* which enables passing of results and errors out of the polling thread.

**is\_polling()** → bool

Check if device status is being polled.

**Returns** *True* when polling thread is set and alive

**start\_polling()** → bool

Start polling.

**Returns** *True* if was not polling before, *False* otherwise

**stop\_polling()** → bool

Stop polling.

Wait for until polling function returns a result as well as any exception that might have been raised within a thread.

**Returns** *True* if was polling before, *False* otherwise, and last result of the polling function call.

**Raises** polling function exceptions


**wait\_for\_polling\_result()**

Wait for until polling function returns a result as well as any exception that might have been raised within a thread.

**Returns** polling function result

**Raises** polling function errors

## hvl\_ccb.dev.visa



```
class VisaDevice(com: Union[hvl_ccb.comm.visa.VisaCommunication,
hvl_ccb.comm.visa.VisaCommunicationConfig, dict], dev_config:
Optional[Union[hvl_ccb.dev.visa.VisaDeviceConfig, dict]] = None)
Bases: hvl_ccb.dev.base.SingleCommDevice

Device communicating over the VISA protocol using VisaCommunication.

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls() → Type[hvl_ccb.comm.visa.VisaCommunication]
    Return the default communication protocol for this device type, which is VisaCommunication.

    Returns the VisaCommunication class

get_error_queue() → str
    Read out error queue and logs the error.

    Returns Error string

get_identification() → str
    Queries “*IDN?” and returns the identification string of the connected device.

    Returns the identification string of the connected device

reset() → None
    Send “*RST” and “*CLS” to the device. Typically sets a defined state.

spoll_handler()
    Reads the status byte and decodes it. The status byte STB is defined in IEEE 488.2. It provides a rough
    overview of the instrument status.

    Returns

start() → None
    Start the VisaDevice. Sets up the status poller and starts it.

    Returns

stop() → None
    Stop the VisaDevice. Stops the polling thread and closes the communication protocol.

    Returns

wait_operation_complete(timeout: Optional[float] = None) → bool
    Waits for a operation complete event. Returns after timeout [s] has expired or the operation complete event
    has been caught.

    Parameters timeout – Time in seconds to wait for the event; None for no timeout.

    Returns True, if OPC event is caught, False if timeout expired

class VisaDeviceConfig(spoll_interval: Union[int, float] = 0.5, spoll_start_delay: Union[int, float] = 2)
Bases: hvl_ccb.dev.visa._VisaDeviceConfigDefaultsBase, hvl_ccb.dev.visa.
_VisaDeviceConfigBase

Configdataclass for a VISA device.

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing
    after value has been forced on fieldname.
```

## Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.


## Module contents

Devices subpackage.

### hvl\_ccb.utils

#### Submodules

##### hvl\_ccb.utils.conversion\_sensor



Sensors that are used by the devices implemented in the CCB

**class LEM4000S(shunt: float = 1.2, calibration\_factor: float = 1)**

Bases: [hvl\\_ccb.utils.conversion\\_sensor.Sensor](#)

Converts the output voltage (V) to the measured current (A) when using a LEM Current transducer LT 4000-S

```
CONVERSION: ClassVar[int] = 5000
calibration_factor: float = 1
convert(value, **kwargs)
shunt: float = 1.2

class LMT70A(temperature_unit: hvl_ccb.utils.conversion_unit.Temperature = Temperature.C)
Bases: hvl_ccb.utils.conversion_sensor.Sensor

Converts the output voltage (V) to the measured temperature (default °C) when using a TI Precision Analog
Temperature Sensor LMT70(A)

LUT: ClassVar[numpy.ndarray[Any,
numpy.dtype[numpy.typing._generic_alias.ScalarType]]] = array([[-55. , 1.375219],
[-50. , 1.350441], [-40. , 1.300593], [-30. , 1.250398], [-20. , 1.199884], [-10. ,
1.14907], [ 0. , 1.097987], [ 10. , 1.046647], [ 20. , 0.99505], [ 30. ,
0.943227], [ 40. , 0.891178], [ 50. , 0.838882], [ 60. , 0.78636], [ 70. ,
0.733608], [ 80. , 0.680654], [ 90. , 0.62749], [100. , 0.574117], [110. ,
0.520551], [120. , 0.46676], [130. , 0.412739], [140. , 0.358164], [150. ,
0.302785]])]


convert(value, **kwargs)
temperature_unit: hvl_ccb.utils.conversion_unit.Temperature = 'C'

class Sensor
Bases: abc.ABC

The BaseClass ‘Sensor’ is designed as a parent for all Sensors. Each attribute must be added to ‘__setattr__’, so
that the value is verified each time the value is changed. It is important to mark attributes that should be constant
with ‘typing.ClassVar[...]’. Together with ‘super().__setattr__(name, value)’, this guarantees that the values are
protected and cannot be altered by the user.

abstract convert(value, **kwargs)
```

## hvl\_ccb.utils.conversion\_unit



Unit conversion, within in the same group of units, for example Kelvin <-> Celsius

```
class Pressure(value)
Bases: hvl_ccb.utils.conversion_unit.Unit

An enumeration.

ATM = 'atm'
ATMOSPHERE = 'atm'
BAR = 'bar'
MILLIMETER_MERCURY = 'mmHg'
MMHG = 'mmHg'
PA = 'Pa'
PASCAL = 'Pa'
POUNDS_PER_SQUARE_INCH = 'psi'
PSI = 'psi'
TORR = 'torr'

class Temperature(value)
Bases: hvl_ccb.utils.conversion_unit.Unit

An enumeration.

C = 'C'
CELSIUS = 'C'
F = 'F'
FAHRENHEIT = 'F'
K = 'K'
KELVIN = 'K'


class Unit(value)
Bases: enum.Enum

An enumeration.

abstract classmethod convert(value, **kwargs)
```

**preserve\_type(func)**

This wrapper preserves the first order type of the input. Upto now the type of the data stored in a list, tuple, array or dict is not preserved. Integer will be converted to float!

**hvl\_ccb.utils.enum**

**class AutoNumberNameEnum**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)

Bases: *hvl\_ccb.utils.enum.NameEnum*, *aenum.AutoNumberEnum*

Auto-numbered enum with names used as string representation, and with lookup and equality based on this representation.

**class BoolEnum**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)

Bases: *hvl\_ccb.utils.enum.NameEnum*

**BoolEnum inherits from NameEnum and the type of the first value is** enforced to be ‘boolean’. For `bool()`-operation the `__bool__` is redefined here.

**class NameEnum**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)

Bases: *hvl\_ccb.utils.enum.StrEnumBase*

Enum with names used as string representation, and with lookup and equality based on this representation.

**class RangeEnum**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)

Bases: *float*, *hvl\_ccb.utils.enum.ValueEnum*

Range enumeration inherit from ValueEnum, find suitable voltage/current/resistance input range for devices such as multimeter and oscilloscope

**abstract classmethod unit()** → str

Returns the Unit of the values in the enumeration. :return: the unit of the values in the enumeration in string format

**class StrEnumBase**(*value=<no\_arg>*, *names=None*, *module=None*, *type=None*, *start=1*, *boundary=None*)

Bases: *aenum.Enum*

String representation-based equality and lookup.

```
class ValueEnum(value=<no_arg>, names=None, module=None, type=None, start=1, boundary=None)
Bases: hvl_ccb.utils.enum.StrEnumBase
```

Enum with string representation of values used as string representation, and with lookup and equality based on this representation.

Attention: to avoid errors, best use together with *unique* enum decorator.

## **hvl\_ccb.utils.typing**

Additional Python typing module utilities

### **ConvertibleTypes**

Typing hint for data type that can be used in conversion

alias of Union[int, float, List[Union[int, float]], Tuple[Union[int, float]], Dict[str, Union[int, float]], numpy.ndarray]

### **Number**

Typing hint auxiliary for a Python base number types: *int* or *float*.

alias of Union[int, float]

### **check\_generic\_type(value, type\_, name='instance')**

Check if *value* is of a generic type *type\_*. Raises *TypeError* if it's not.

#### **Parameters**

- **name** – name to report in case of an error
- **value** – value to check
- **type** – generic type to check against

### **is\_generic\_type\_hint(type\_)**

Check if class is a generic type, for example *Union[int, float]*, *List[int]*, or *List*.

#### **Parameters** **type** – type to check

## **hvl\_ccb.utils.validation**

```
validate_and_resolve_host(host: Union[str, ipaddress.IPv4Address, ipaddress.IPv6Address], logger:
    Optional[logging.Logger] = None) → str
```

```
validate_bool(x_name: str, x: object, logger: Optional[logging.Logger] = None) → None
```

Validate if given input *x* is a *bool*.

#### **Parameters**

- **x\_name** – string name of the validate input, use for the error message
- **x** – an input object to validate as boolean
- **logger** – logger of the calling submodule

**Raises** **TypeError** – when the validated input does not have boolean type

```
validate_number(x_name: str, x: object, limits: typing.Optional[typing.Tuple] = (None, None), number_type: typing.Union[typing.Type[typing.Union[int, float]], typing.Tuple[typing.Type[typing.Union[int, float]], ...]] = (<class 'int'>, <class 'float'>), logger: typing.Optional[logging.Logger] = None) → None
```

Validate if given input *x* is a number of given *number\_type* type, with value between given *limits[0]* and *limits[1]* (inclusive), if not *None*. For array-like objects (npt.NDArray, List, Tuple, Dict) it is checked if all elements are within the limits and have the correct type.

#### Parameters

- **x\_name** – string name of the validate input, use for the error message
- **x** – an input object to validate as number of given type within given range
- **logger** – logger of the calling submodule
- **limits** – [lower, upper] limit, with *None* denoting no limit: [-inf, +inf]
- **number\_type** – expected type or tuple of types of a number, by default (*int, float*)

#### Raises


- **TypeError** – when the validated input does not have expected type
- **ValueError** – when the validated input has correct number type but is not within given range

```
validate_tcp_port(port: int, logger: Optional[logging.Logger] = None)
```

### Module contents

#### 4.1.2 Submodules

##### hvl\_ccb.configuration



providing classes for handling configuration for communication protocols and devices.

Facilities

```
class ConfigurationMixin(configuration)
```

Bases: abc.ABC

Mixin providing configuration to a class.

```
property config
```

ConfigDataclass property.

```
Returns the configuration
```

```
abstract static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

configuration_save_json(path: str) → None
    Save current configuration as JSON file.

    Parameters path – path to the JSON file.

classmethod from_json(filename: str)
    Instantiate communication protocol using configuration from a JSON file.

    Parameters filename – Path and filename to the JSON configuration

exception ConfigurationValueWarning
    Bases: UserWarning

    User warnings category for values of @configdataclass fields.

class EmptyConfig
    Bases: object

    Empty configuration dataclass.

    clean_values()
        Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

    force_value(fieldname, value)
        Forces a value to a dataclass field despite the class being frozen.

        NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.

        Parameters
            • fieldname – name of the field
            • value – value to assign

    is_configdataclass = True

    classmethod keys() → Sequence[str]
        Returns a list of all configdataclass fields key-names.

        Returns a list of strings containing all keys.

    classmethod optional_defaults() → Dict[str, object]
        Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

        Returns a list of strings containing all optional keys.

    classmethod required_keys() → Sequence[str]
        Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.


        Returns a list of strings containing all required keys.

configdataclass(direct_decoration=None, frozen=True)
    Decorator to make a class a configdataclass. Types in these dataclasses are enforced. Implement a function clean_values(self) to do additional checking on value ranges etc.
```

It is possible to inherit from a configdataclass and re-decorate it with @configdataclass. In a subclass, default values can be added to existing fields. Note: adding additional non-default fields is prone to errors, since the order has to be respected through the whole chain (first non-default fields, only then default-fields).

**Parameters** `frozen` – defaults to True. False allows to later change configuration values. Attention:  
if configdataclass is not frozen and a value is changed, typing is not enforced anymore!


### `hvl_ccb.exception`



Introduce a common base exception for the CCB

**exception** `CCBException`  
Bases: `Exception`

### `hvl_ccb.experiment_manager`



Main module containing the top level `ExperimentManager` class. Inherit from this class to implement your own experiment functionality in another project and it will help you start, stop and manage your devices.

**exception** `ExperimentError`  
Bases: `Exception`

Exception to indicate that the current status of the experiment manager is on ERROR and thus no operations can be made until reset.

**class** `ExperimentManager(*args, **kwargs)`  
Bases: `hvl_ccb.dev.base.DeviceSequenceMixin`

Experiment Manager can start and stop communication protocols and devices. It provides methods to queue commands to devices and collect results.

**add\_device**(*name*: str, *device*: hvl\_ccb.dev.base.Device) → *None*

Add a new device to the manager. If the experiment is running, automatically start the device. If a device with this name already exists, raise an exception.

#### Parameters

- **name** – is the name of the device.
- **device** – is the instantiated Device object.

**Raises** *DeviceExistingException* –

**devices\_failed\_start**: Dict[str, hvl\_ccb.dev.base.Device]

Dictionary of named device instances from the sequence for which the most recent *start()* attempt failed.

Empty if *stop()* was called last; cf. *devices\_failed\_stop*.

**devices\_failed\_stop**: Dict[str, hvl\_ccb.dev.base.Device]

Dictionary of named device instances from the sequence for which the most recent *stop()* attempt failed.

Empty if *start()* was called last; cf. *devices\_failed\_start*.

**finish()** → *None*

Stop experimental setup, stop all devices.

**is\_error()** → bool

Returns true, if the status of the experiment manager is *error*.

**Returns** True if on error, false otherwise

**is\_finished()** → bool

Returns true, if the status of the experiment manager is *finished*.

**Returns** True if finished, false otherwise

**is\_running()** → bool

Returns true, if the status of the experiment manager is *running*.

**Returns** True if running, false otherwise

**run()** → *None*

Start experimental setup, start all devices.

**start()** → *None*

Alias for ExperimentManager.run()

**property status**: hvl\_ccb.experiment\_manager.ExperimentStatus

Get experiment status.

**Returns** experiment status enum code.

**stop()** → *None*

Alias for ExperimentManager.finish()

**class ExperimentStatus**(*value*)

Bases: enum.Enum

Enumeration for the experiment status

**ERROR** = 5

**FINISHED** = 4

**FINISHING** = 3

**INITIALIZED** = 0

```
INITIALIZING = -1
RUNNING = 2
STARTING = 1
```

#### 4.1.3 Module contents

Top-level package for HVL Common Code Base.

## **CONTRIBUTING**

Contributions are welcome, and they are greatly appreciated! Every little bit helps, and credit will always be given.

You can contribute in many ways:

### **5.1 Types of Contributions**

#### **5.1.1 Report Bugs**

Report bugs at [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/issues](https://gitlab.com/ethz_hvl/hvl_ccb/issues).

If you are reporting a bug, please include:

- Your operating system name and version.
- Any details about your local setup that might be helpful in troubleshooting.
- Detailed steps to reproduce the bug.

#### **5.1.2 Fix Bugs**

Look through the GitLab issues for bugs. Anything tagged with “bug” and “help wanted” is open to whoever wants to implement it.

#### **5.1.3 Implement Features**

Look through the GitLab issues for features. Anything tagged with “enhancement” and “help wanted” is open to whoever wants to implement it.

#### **5.1.4 Write Documentation**

HVL Common Code Base could always use more documentation, whether as part of the official HVL Common Code Base docs, in docstrings, or even on the web in blog posts, articles, and such.

### 5.1.5 Submit Feedback

The best way to send feedback is to file an issue at [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/issues](https://gitlab.com/ethz_hvl/hvl_ccb/issues).

If you are proposing a feature:

- Explain in detail how it would work.
- Keep the scope as narrow as possible, to make it easier to implement.
- Remember that this is a volunteer-driven project, and that contributions are welcome :)

## 5.2 Get Started!

Ready to contribute? Here's how to set up *hvl\_ccb* for local development.

1. Clone *hvl\_ccb* repo from GitLab.

```
$ git clone git@gitlab.com:ethz_hvl/hvl_ccb.git
```

2. Install your local copy into a virtualenv. Assuming you have `virtualenvwrapper` installed, this is how you set up your fork for local development:

```
$ mkvirtualenv hvl_ccb
$ cd hvl_ccb/
$ pip install -e .[all]
$ pip install -r requirements_dev.txt
```

3. Create a branch for local development:

```
$ git checkout -b name-of-your-bugfix-or-feature
```

Now you can make your changes locally.

4. When you're done making changes, check that your changes pass flake8 and the tests, including testing other Python versions with tox:

```
$ flake8 hvl_ccb tests
$ python setup.py test or py.test
$ tox
```

To get flake8 and tox, just pip install them into your virtualenv. You can also use the provided make-like shell script to run flake8 and tests:

```
$ ./make.sh lint
$ ./make.sh test
```

5. Commit your changes and push your branch to GitLab:

```
$ git add .
$ git commit -m "Your detailed description of your changes."
$ git push origin name-of-your-bugfix-or-feature
```

6. Submit a merge request through the GitLab website.

## 5.3 Merge Request Guidelines

Before you submit a merge request, check that it meets these guidelines:

1. The merge request should include tests.
2. If the merge request adds functionality, the docs should be updated. Put your new functionality into a function with a docstring, and add the feature to the list in README.rst.
3. The merge request should work for Python 3.7. Check [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/merge\\_requests](https://gitlab.com/ethz_hvl/hvl_ccb/merge_requests) and make sure that the tests pass for all supported Python versions.

## 5.4 Tips

- To run tests from a single file:

```
$ py.test tests/test_hvl_ccb.py
```

or a single test function:

```
$ py.test tests/test_hvl_ccb.py::test_command_line_interface
```

- If your tests are slow, profile them using the pytest-profiling plugin:

```
$ py.test tests/test_hvl_ccb.py --profile
```

or for a graphical overview (you need a SVG image viewer):

```
$ py.test tests/test_hvl_ccb.py --profile-svg  
$ open prof/combined.svg
```

- To add dependency, edit appropriate \*requirements variable in the setup.py file and re-run:

```
$ python setup.py develop
```

- To generate a PDF version of the Sphinx documentation instead of HTML use:

```
$ rm -rf docs/hvl_ccb.rst docs/modules.rst docs/_build && sphinx-apidoc -o docs/hvl_ccb && python -msphinx -M latexpdf docs/ docs/_build
```

This command can also be run through the make-like shell script:

```
$ ./make.sh docs-pdf
```

This requires a local installation of a LaTeX distribution, e.g. MikTeX.

## 5.5 Deploying

A reminder for the maintainers on how to deploy. Create `release-N.M.K` branch. Make sure all your changes are committed. Update or create entry in `HISTORY.rst` file, and, if applicable, update `AUTHORS.rst` file, update features tables in `README.rst` file, and update API docs:

```
$ make docs
```

Commit all of the above, except for the `docs/hvl_ccb.dev.picotech_pt104.rst`, and then run:

```
$ bumpversion patch # possible: major / minor / patch  
$ git push  
$ git push --tags
```

Go to <https://readthedocs.org/projects/hvlccb/builds/> and check if RTD docs build for the pushed tag passed.

Wait for the CI pipeline to finish successfully; afterwards, run a release check:

```
$ make release-check
```

Finally, prepare and push a release:

```
$ make release
```

Merge the release branch into master and devel branches with `--no-ff` flag and delete the release branch:

```
$ git switch master  
$ git merge --no-ff release-N.M.K  
$ git push  
$ git switch devel  
$ git merge --no-ff release-N.M.K  
$ git push  
$ git push --delete origin release-N.M.K  
$ git branch --delete release-N.M.K
```

Finally, prepare GitLab release and cleanup the corresponding milestone:

1. go to [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/-/tags/](https://gitlab.com/ethz_hvl/hvl_ccb/-/tags/), select the latest release tag, press “Edit release notes” and add the release notes (copy a corresponding entry from `HISTORY.rst` file with formatting adjusted from ReStructuredText to Markdown); press “Save changes”;
2. go to [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/-/releases](https://gitlab.com/ethz_hvl/hvl_ccb/-/releases), select the latest release, press “Edit this release” and under “Milestones” select the corresponding milestone; press “Save changes”;
3. go to [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/-/milestones](https://gitlab.com/ethz_hvl/hvl_ccb/-/milestones), make sure that it is 100% complete (otherwise, create a next patch-level milestone and assign it to the ongoing Issues and Merge Requests therein); press “Close Milestone”.

**CREDITS**

## 6.1 Maintainers

- Mikołaj Rybiński <mikolaj.rybinski@id.ethz.ch>
- Henrik Menne <henrik.menne@eeh.ee.ethz.ch>
- Henning Janssen <janssen@eeh.ee.ethz.ch>
- Maria Del <maria.del@id.ethz.ch>

## 6.2 Authors

- Mikołaj Rybiński <mikolaj.rybinski@id.ethz.ch>
- David Gruber <dev@davidgruber.ch>
- Henrik Menne <henrik.menne@eeh.ee.ethz.ch>
- Alise Chachereau <chachereau@eeh.ee.ethz.ch>
- Henning Janssen <janssen@eeh.ee.ethz.ch>
- David Taylor <dtaylor@ethz.ch>
- Joseph Engelbrecht <engelbrecht@eeh.ee.ethz.ch>

## 6.3 Contributors

- Luca Nembrini <lucane@student.ethz.ch>
- Maria Del <maria.del@id.ethz.ch>
- Raphael Faerber <raphael.faerber@eeh.ee.ethz.ch>
- Ruben Stadler <rstadler@student.ethz.ch>
- Hanut Vemulapalli <vemulapalli@eeh.ee.ethz.ch>
- Chi-Ching Hsu <hsu@ibi.baug.ethz.ch>



## HISTORY

### 7.1 0.10.1 (2022-01-24)

- several improvements and fixes for device cube:
  - privatize Alarms and AlarmsOverview
  - fix list of cube alarms
  - improve docs
  - fix bugs with earthing sticks
  - fix bug in config dataclass of cube
- introduction of BoolEnum
- introduction of RangeEnum
- bumpversion -> bump2version

### 7.2 0.10.0 (2022-01-17)

- Reimplementation of the Cube (before known as Supercube)
- new names:
  - Supercube Typ B -> BaseCube
  - Supercube Typ A -> PICube (power inverter Cube)
- new import:
  - `from hvl_ccb.dev.supercube import SupercubeB` -> `from hvl_ccb.dev.cube import BaseCube`
- new programming style:
  - getter / setter methods -> properties
  - e.g. get: `cube.get_support_output(port=1, contact=1)` -> `cube.support_1.output_1`
  - e.g. set: `cube.get_support_output(port=1, contact=1, state=True)` -> `cube.support_1.output_1 = True`
- unify Exceptions of Cube
- implement Fast Switch-Off of Cube

- remove method `support_output_impulse`
- all active alarms can now be queried `cube.active_alarms()`
- alarms will now result in different logging levels depending on the seriousness of the alarm.
- introduction of limits for slope and safety limit for RedReady
- during the startup the CCB will update the time of the cube.
- verification of inputs
- polarity of DC voltage
- Switch from `python-opcua` to `opcua-asyncio` (former package is no longer maintained)

## 7.3 0.9.0 (2022-01-07)

- New device: Highland T560 digital delay and pulse generator over Telnet.
- **Rework of the Technix Capacitor Charger.**
  - Moved into a separate sub-package
  - NEW import over `import hvl_ccb.dev.technix as XXX`
  - Slightly adapted behaviour
- Add `validate_tcp_port` to validate port number.
- **Add `validate_and_resolve_host` to validate and resolve host names and IPs.**
  - Remove requirement IPy
- Add a unified CCB Exception schema for all devices and communication protocols.
- Add data conversion functions to README.
- Update CI and devel images from Debian 10 buster to Debian 11 bullseye.
- Fix typing due to numpy update.
- Fix incorrect overloading of `clean_values()` in classes of type `XCommunicationConfig`.

## 7.4 0.8.5 (2021-11-05)

- Added arbitrary waveform for TiePie signal generation, configurable via `dev.tiepie.generator.TiePieGeneratorConfig.waveform` property.
- In `utils.conversion_sensor`: improvements for class constants; removed SciPy dependency.
- Added Python 3.10 support.

## 7.5 0.8.4 (2021-10-22)

- `utils.validation.validate_number` extension to handle NumPy arrays and array-like objects.
- `utils.conversion_unit` utility classes handle correctly `NamedTuple` instances.
- `utils.conversion_sensor` and `utils.conversion_unit` code simplification (no `transfer_function_order` attribute) and cleanups.
- Fixed incorrect error logging in `configuration.configdataclass`.
- `comm.telnet.TelnetCommunication` tests fixes for local run errors.

## 7.6 0.8.3 (2021-09-27)

- New data conversion functions in `utils.conversion_sensor` and `utils.conversion_unit` modules. Note: to use these functions you must install `hvl_ccb` with extra requirement, either `hvl_ccb[conversion]` or `hvl_ccb[all]`.
- Improved documentation with respect to installation of external libraries.

## 7.7 0.8.2 (2021-08-27)

- New functionality in `dev.labjack.LabJack`:
  - configure clock and send timed pulse sequences
  - set DAC/analog output voltage
- Bugfix: ignore random bits sent by to `dev.newport.NewportSMC100PP` controller during start-up/powering-up.

## 7.8 0.8.1 (2021-08-13)

- Add Python version check (min version error; max version warning).
- Daily checks for upstream dependencies compatibility and devel environment improvements.

## 7.9 0.8.0 (2021-07-02)

- TCP communication protocol.
- Lauda PRO RP 245 E circulation thermostat device over TCP.
- Pico Technology PT-104 Platinum Resistance Data Logger device as a wrapper of the Python bindings for the PicoSDK.
- In `com.visa.VisaCommunication`: periodic status polling when VISA/TCP keep alive connection is not supported by a host.

## 7.10 0.7.1 (2021-06-04)

- New `utils.validation` submodule with `validate_bool` and `validate_number` utilities extracted from internal use within a `dev.tiepie` subpackage.
- In `comm.serial.SerialCommunication`:
  - strict encoding errors handling strategy for subclasses,
  - user warning for a low communication timeout value.

## 7.11 0.7.0 (2021-05-25)

- The `dev.tiepie` module was splitted into a subpackage with, in particular, submodules for each of the device types – oscilloscope, generator, and i2c – and with backward-incompatible direct imports from the submodules.
- In `dev.technix`:
  - fixed communication crash on nested status byte query;
  - added enums for GET and SET register commands.
- Further minor logging improvements: added missing module level logger and removed some error logs in `except` blocks used for a flow control.
- In `examples/` folder renamed consistently all the examples.
- In API documentation: fix incorrect links mapping on inheritance diagrams.

## 7.12 0.6.1 (2021-05-08)

- In `dev.tiepie`:
  - dynamically set oscilloscope's channel limits in `OscilloscopeChannelParameterLimits`: `input_range` and `trigger_level_abs`, incl. update of latter on each change of `input_range` value of a `TiePieOscilloscopeChannelConfig` instances;
  - quick fix for opening of combined instruments by disabling `OscilloscopeParameterLimits`. `trigger_delay` (an advanced feature);
  - enable automatic devices detection to be able to find network devices with `TiePieOscilloscope.list_devices()`.
- Fix `examples/example_labjack.py`.
- Improved logging: consistently use module level loggers, and always log exception tracebacks.
- Improve API documentation: separate pages per modules, each with an inheritance diagram as an overview.

## 7.13 0.6.0 (2021-04-23)

- Technix capacitor charger using either serial connection or Telnet protocol.
- **Extensions, improvements and fixes in existing devices:**
  - In `dev.tiepie.TiePieOscilloscope`:
    - \* redesigned measurement start and data collection API, incl. time out argument, with no/infinite time out option;
    - \* trigger allows now a no/infinite time out;
    - \* record length and trigger level were fixed to accept, respectively, floating point and integer numbers;
    - \* fixed resolution validation bug;
  - `dev.heinzinger.HeinzingerDI` and `dev.rs_rto1024.RT01024` instances are now resilient to multiple `stop()` calls.
  - In `dev.crylas.CryLasLaser`: default configuration timeout and polling period were adjusted;
  - Fixed PSI9080 example script.
- **Package and source code improvements:**
  - Update to backward-incompatible `pyvisa-py>=0.5.2`. Developers, do update your local development environments!
  - External libraries, like LibTiePie SDK or LJM Library, are now not installed by default; they are now extra installation options.
  - Added Python 3.9 support.
  - Improved number formatting in logs.
  - Typing improvements and fixes for `mypy>=0.800`.

## 7.14 0.5.0 (2020-11-11)

- TiePie USB oscilloscope, generator and I2C host devices, as a wrapper of the Python bindings for the LibTiePie SDK.
- a FuG Elektronik Power Supply (e.g. Capacitor Charger HCK) using the built-in ADDAT controller with the Probus V protocol over a serial connection
- All devices polling status or measurements use now a `dev.utils.Poller` utility class.
- **Extensions and improvements in existing devices:**
  - In `dev.rs_rto1024.RT01024`: added Channel state, scale, range, position and offset accessors, and measurements activation and read methods.
  - In `dev.sst_luminox.Luminox`: added querying for all measurements in polling mode, and made output mode activation more robust.
  - In `dev.newport.NewportSMC100PP`: an error-prone `wait_until_move_finished` method of replaced by a fixed waiting time, device operations are now robust to a power supply cut, and device restart is not required to apply a start configuration.
- **Other minor improvements:**

- Single failure-safe starting and stopping of devices sequenced via `dev.base.DeviceSequenceMixin`.
- Moved `read_text_nonempty` up to `comm.serial.SerialCommunication`.
- Added development Dockerfile.
- Updated package and development dependencies: `pymodbus`, `pytest-mock`.

## 7.15 0.4.0 (2020-07-16)

- **Significantly improved new Supercube device controller:**
  - more robust error-handling,
  - status polling with generic Poller helper,
  - messages and status boards.
  - tested with a physical device,
- Improved OPC UA client wrapper, with better error handling, incl. re-tries on `concurrent.futures.TimeoutError`.
- SST Luminox Oxygen sensor device controller.
- **Backward-incompatible changes:**
  - `CommunicationProtocol.access_lock` has changed type from `threading.Lock` to `threading.RLock`.
  - `ILS2T.relative_step` and `ILS2T.absolute_position` are now called, respectively, `ILS2T.write_relative_step` and `ILS2T.write_absolute_position`.
- **Minor bugfixes and improvements:**
  - fix use of max resolution in `Labjack.set_ain_resolution()`,
  - resolve ILS2T devices relative and absolute position setters race condition,
  - added acoustic horn function in the 2015 Supercube.
- **Toolchain changes:**
  - add Python 3.8 support,
  - drop pytest-runner support,
  - ensure compatibility with `labjack_ljm` 2019 version library.

## 7.16 0.3.5 (2020-02-18)

- Fix issue with reading integers from LabJack LJM Library (device's product ID, serial number etc.)
- Fix development requirements specification (tox version).

## 7.17 0.3.4 (2019-12-20)

- New devices using serial connection:
  - Heinzinger Digital Interface I/II and a Heinzinger PNC power supply
  - Q-switched Pulsed Laser and a laser attenuator from CryLas
  - Newport SMC100PP single axis motion controller for 2-phase stepper motors
  - Pfeiffer TPG controller (TPG 25x, TPG 26x and TPG 36x) for Compact pressure Gauges
- PEP 561 compatibility and related corrections for static type checking (now in CI)
- Refactorings:
  - Protected non-thread safe read and write in communication protocols
  - Device sequence mixin: start/stop, add/rm and lookup
  - `.format()` to f-strings
  - more enumerations and a quite some improvements of existing code
- Improved error docstrings (`:raises:` annotations) and extended tests for errors.

## 7.18 0.3.3 (2019-05-08)

- Use PyPI labjack-ljm (no external dependencies)

## 7.19 0.3.2 (2019-05-08)

- INSTALLATION.rst with LJMPython prerequisite info

## 7.20 0.3.1 (2019-05-02)

- readthedocs.org support

## 7.21 0.3 (2019-05-02)

- Prevent an automatic close of VISA connection when not used.
- Rhode & Schwarz RTO 1024 oscilloscope using VISA interface over `TCP::INSTR`.
- Extended tests incl. messages sent to devices.
- Added Supercube device using an OPC UA client
- Added Supercube 2015 device using an OPC UA client (for interfacing with old system version)

## 7.22 0.2.1 (2019-04-01)

- Fix issue with LJMPython not being installed automatically with setuptools.

## 7.23 0.2.0 (2019-03-31)

- LabJack LJM Library communication wrapper and LabJack device.
- Modbus TCP communication protocol.
- Schneider Electric ILS2T stepper motor drive device.
- Elektro-Automatik PSI9000 current source device and VISA communication wrapper.
- Separate configuration classes for communication protocols and devices.
- Simple experiment manager class.

## 7.24 0.1.0 (2019-02-06)

- Communication protocol base and serial communication implementation.
- Device base and MBW973 implementation.

---

**CHAPTER  
EIGHT**

---

**INDICES AND TABLES**

- genindex
- modindex
- search



## PYTHON MODULE INDEX

### h

hvl\_ccb, 182

hvl\_ccb.comm, 33

hvl\_ccb.comm.base, 11

hvl\_ccb.comm.labjack\_ljm, 15

hvl\_ccb.comm.modbus\_tcp, 18

hvl\_ccb.comm.opc, 20

hvl\_ccb.comm.serial, 23

hvl\_ccb.comm.tcp, 26

hvl\_ccb.comm.telnet, 28

hvl\_ccb.comm.visa, 30

hvl\_ccb.configuration, 178

hvl\_ccb.dev, 173

hvl\_ccb.dev.base, 83

hvl\_ccb.dev.crylas, 85

hvl\_ccb.dev.cube, 44

hvl\_ccb.dev.cube.base, 33

hvl\_ccb.dev.cube.constants, 38

hvl\_ccb.dev.cube.picube, 41

hvl\_ccb.dev.ea\_psi9000, 95

hvl\_ccb.dev.fug, 99

hvl\_ccb.dev.heinzinger, 110

hvl\_ccb.dev.highland\_t560, 49

hvl\_ccb.dev.highland\_t560.base, 44

hvl\_ccb.dev.highland\_t560.channel, 46

hvl\_ccb.dev.highland\_t560.device, 47

hvl\_ccb.dev.labjack, 116

hvl\_ccb.dev.lauda, 123

hvl\_ccb.dev.mbw973, 129

hvl\_ccb.dev.newport, 133

hvl\_ccb.dev.pfeiffer\_tpg, 146

hvl\_ccb.dev.rs\_rto1024, 152

hvl\_ccb.dev.se\_ilst2t, 160

hvl\_ccb.dev.sst\_luminox, 166

hvl\_ccb.dev.supercube2015, 64

hvl\_ccb.dev.supercube2015.base, 49

hvl\_ccb.dev.supercube2015.constants, 55

hvl\_ccb.dev.supercube2015.typ\_a, 62

hvl\_ccb.dev.technix, 70

hvl\_ccb.dev.technix.base, 65

hvl\_ccb.dev.technix.device, 68

hvl\_ccb.dev.tiepie, 82

hvl\_ccb.dev.tiepie.base, 71

hvl\_ccb.dev.tiepie.channel, 73

hvl\_ccb.dev.tiepie.device, 75

hvl\_ccb.dev.tiepie.generator, 76

hvl\_ccb.dev.tiepie.i2c, 78

hvl\_ccb.dev.tiepie.oscilloscope, 79

hvl\_ccb.dev.tiepie.utils, 81

hvl\_ccb.dev.utils, 170

hvl\_ccb.dev.visa, 171

hvl\_ccb.exception, 180

hvl\_ccb.experiment\_manager, 180

hvl\_ccb.utils, 178

hvl\_ccb.utils.conversion\_sensor, 173

hvl\_ccb.utils.conversion\_unit, 174

hvl\_ccb.utils.enum, 176

hvl\_ccb.utils.typing, 177

hvl\_ccb.utils.validation, 177



# INDEX

## A

A (*HeinzingerPNC.UnitCurrent attribute*), 113  
A (*SupercubeOpcEndpoint attribute*), 62  
ABSOLUTE (*TiePieOscilloscopeTriggerLevelMode attribute*), 75  
ABSOLUTE\_POSITION (*ILS2T.ActionsPtp attribute*), 160  
AC (*NewportConfigCommands attribute*), 133  
AC\_100KV (*PowerSetup attribute*), 40  
AC\_150KV (*PowerSetup attribute*), 40  
AC\_200KV (*PowerSetup attribute*), 40  
AC\_50KV (*PowerSetup attribute*), 40  
AC\_DoubleStage\_150kV (*PowerSetup attribute*), 60  
AC\_DoubleStage\_200kV (*PowerSetup attribute*), 60  
AC\_SingleStage\_100kV (*PowerSetup attribute*), 60  
AC\_SingleStage\_50kV (*PowerSetup attribute*), 60  
ACA (*TiePieOscilloscopeChannelCoupling attribute*), 74  
acceleration (*NewportSMC100PPConfig attribute*), 141  
ACCESS\_ENABLE (*ILS2TRegAddr attribute*), 165  
access\_lock (*CommunicationProtocol attribute*), 14  
ACTION\_JOG\_VALUE (*ILS2T attribute*), 160  
activate\_clock\_output() (*T560 method*), 47  
activate\_measurements() (*RTO1024 method*), 153  
activate\_output() (*Luminox method*), 166  
activated (*BreakdownDetection attribute*), 56  
ACTIVE (*CryLasLaser.AnswersStatus attribute*), 88  
active\_alarms() (*BaseCube method*), 33  
ACTIVE\_HIGH (*Polarity attribute*), 44  
ACTIVE\_LOW (*Polarity attribute*), 45  
actualsetvalue (*FuGProbusVSetRegisters property*), 107  
ACV (*TiePieOscilloscopeChannelCoupling attribute*), 74  
adc\_mode (*FuGProbusVMonitorRegisters property*), 107  
add\_device() (*DeviceSequenceMixin method*), 83  
add\_device() (*ExperimentManager method*), 180  
ADDR\_INCORRECT (*NewportSMC100PPSerialCommunication.ControllerErrors attribute*), 142  
address (*NewportSMC100PPConfig attribute*), 141  
address (*VisaCommunicationConfig property*), 31  
address() (*VisaCommunicationConfig.InterfaceType method*), 31

ADMODE (*FuGProbusIVCommands attribute*), 104  
Alarm0 (*AlarmText attribute*), 55  
Alarm1 (*AlarmText attribute*), 55  
Alarm10 (*AlarmText attribute*), 55  
Alarm11 (*AlarmText attribute*), 55  
Alarm12 (*AlarmText attribute*), 55  
Alarm13 (*AlarmText attribute*), 55  
Alarm14 (*AlarmText attribute*), 55  
Alarm17 (*AlarmText attribute*), 55  
Alarm19 (*AlarmText attribute*), 55  
Alarm2 (*AlarmText attribute*), 55  
Alarm20 (*AlarmText attribute*), 55  
Alarm21 (*AlarmText attribute*), 55  
Alarm22 (*AlarmText attribute*), 55  
Alarm3 (*AlarmText attribute*), 55  
Alarm4 (*AlarmText attribute*), 55  
Alarm5 (*AlarmText attribute*), 55  
Alarm6 (*AlarmText attribute*), 55  
Alarm7 (*AlarmText attribute*), 55  
Alarm8 (*AlarmText attribute*), 55  
Alarm9 (*AlarmText attribute*), 55  
AlarmText (*class in hvl\_ccb.dev.supercube2015.constants*), 55  
ALL (*TiePieOscilloscopeAutoResolutionModes attribute*), 80  
all\_measurements (*LuminoxMeasurementType attribute*), 168  
all\_measurements\_types() (*LuminoxMeasurementType class method*), 168  
amplitude (*TiePieGeneratorConfig property*), 76  
ANALOG (*LaudaProRp245eConfig.ExtControlModeEnum attribute*), 126  
analog\_control (*FuGProbusVDIRegisters property*), 106  
ANY (*LabJack.DeviceType attribute*), 118  
ANY (*LJMCommunicationConfig.ConnectionType attribute*), 16  
ANY (*LJMCommunicationConfig.DeviceType attribute*), 16  
ANY (*TiePieOscilloscopeTriggerKind attribute*), 75  
ARBITRARY (*TiePieGeneratorSignalType attribute*), 77  
AsyncCommunicationProtocol (*class in*

*hvl\_ccb.comm.base), 11*

**AsyncCommunicationProtocolConfig** (class in *hvl\_ccb.comm.base*), 13

**ATM** (*Pressure attribute*), 175

**ATMOSPHERE** (*Pressure attribute*), 175

**attenuation** (*CryLasAttenuator property*), 86

**AUTO** (*EarthingStickOperatingStatus attribute*), 39

**AUTO** (*LaudaProRp245eConfig.OperationModeEnum attribute*), 127

**AUTO** (*RTO1024.TriggerModes attribute*), 153

**auto\_install\_mode** (*T560 property*), 47

**auto\_install\_mode** (*T560Config attribute*), 48

**auto\_laser\_on** (*CryLasLaserConfig attribute*), 91

**auto\_resolution\_mode** (*TiePieOscilloscopeConfig property*), 80

**AutoInstallMode** (class in *hvl\_ccb.dev.hIGHLAND\_T560.base*), 44

**AutoNumberNameEnum** (class in *hvl\_ccb.utils.enum*), 176

**B**

**B** (*SupercubeOpcEndpoint attribute*), 62

**BA** (*NewportConfigCommands attribute*), 133

**backlash\_compensation** (*NewportSMC100PPConfig attribute*), 141

**backup\_waveform()** (*RTO1024 method*), 153

**bar** (*PfeifferTPG.PressureUnits attribute*), 147

**BAR** (*Pressure attribute*), 175

**barometric\_pressure** (*LuminoxMeasurementType attribute*), 168

**base\_velocity** (*NewportSMC100PPConfig attribute*), 141

**BaseCube** (class in *hvl\_ccb.dev.cube.base*), 33

**BaseCubeConfiguration** (class in *hvl\_ccb.dev.cube.base*), 36

**BaseCubeOpcUaCommunication** (class in *hvl\_ccb.dev.cube.base*), 37

**BaseCubeOpcUaCommunicationConfig** (class in *hvl\_ccb.dev.cube.base*), 37

**BATH\_TEMP** (*LaudaProRp245eCommand attribute*), 125

**baudrate** (*CryLasAttenuatorSerialCommunicationConfig attribute*), 87

**baudrate** (*CryLasLaserSerialCommunicationConfig attribute*), 93

**baudrate** (*FuGSerialCommunicationConfig attribute*), 109

**baudrate** (*HeinzingerSerialCommunicationConfig attribute*), 115

**baudrate** (*LuminoxSerialCommunicationConfig attribute*), 169

**baudrate** (*MBW973SerialCommunicationConfig attribute*), 132

**baudrate** (*NewportSMC100PPSerialCommunicationConfig attribute*), 144

**baudrate** (*PfeifferTPGSerialCommunicationConfig attribute*), 151

**baudrate** (*SerialCommunicationConfig attribute*), 24

**BH** (*NewportConfigCommands attribute*), 133

**board** (*VisaCommunicationConfig attribute*), 31

**BoolEnum** (class in *hvl\_ccb.utils.enum*), 176

**breakdown\_detection\_active** (*BaseCube property*), 33

**breakdown\_detection\_reset()** (*BaseCube method*), 34

**breakdown\_detection\_triggered** (*BaseCube property*), 34

**BreakdownDetection** (class in *hvl\_ccb.dev.supercube2015.constants*), 56

**bufsize** (*TcpCommunicationConfig attribute*), 27

**build\_str()** (*LaudaProRp245eCommand method*), 126

**bytesize** (*CryLasAttenuatorSerialCommunicationConfig attribute*), 87

**bytesize** (*CryLasLaserSerialCommunicationConfig attribute*), 94

**bytesize** (*FuGSerialCommunicationConfig attribute*), 109

**bytesize** (*HeinzingerSerialCommunicationConfig attribute*), 115

**bytesize** (*LuminoxSerialCommunicationConfig attribute*), 169

**bytesize** (*MBW973SerialCommunicationConfig attribute*), 132

**bytesize** (*NewportSMC100PPSerialCommunicationConfig attribute*), 145

**bytesize** (*PfeifferTPGSerialCommunicationConfig attribute*), 151

**Bytesize** (*SerialCommunicationConfig attribute*), 24

**bytesize** (*SerialCommunicationConfig attribute*), 24

**C**

**C** (*LabJack.TemperatureUnit attribute*), 118

**C** (*LabJack.ThermocoupleType attribute*), 119

**C** (*Temperature attribute*), 175

**calibration\_factor** (*CryLasLaserConfig attribute*), 91

**calibration\_factor** (*LEM4000S attribute*), 174

**calibration\_mode** (*FuGProbusVDIRegisters property*), 106

**cc\_mode** (*FuGProbusVDIRegisters property*), 106

**CCBException**, 180

**cee16** (*GeneralSockets attribute*), 58

**cee16\_socket** (*BaseCube property*), 34

**CELSIUS** (*Temperature attribute*), 175

**ch\_a** (*T560 property*), 47

**ch\_b** (*T560 property*), 47

**ch\_c** (*T560 property*), 47

**ch\_d** (*T560 property*), 47

**channels\_enabled** (*TiePieOscilloscope property*), 79

**check\_for\_error()** (*New-portSMC100PPSerialCommunication method*),  
 143  
**check\_generic\_type()** (*in module hvl\_ccb.utils.typing*), 177  
**check\_master\_slave\_config()** (*PSI9000 method*),  
 95  
**chunk\_size** (*VisaCommunicationConfig attribute*), 31  
**clean\_amplitude()** (*TiePieGeneratorConfig method*),  
 76  
**clean\_auto\_resolution\_mode()** (*TiePieOscilloscopeConfig static method*), 80  
**clean\_coupling()** (*TiePieOscilloscopeChannelConfig static method*), 73  
**clean\_enabled()** (*TiePieGeneratorConfig static method*), 76  
**clean\_enabled()** (*TiePieOscilloscopeChannelConfig static method*), 73  
**clean\_frequency()** (*TiePieGeneratorConfig method*),  
 76  
**clean\_input\_range()** (*TiePieOscilloscopeChannelConfig method*), 73  
**clean\_offset()** (*TiePieGeneratorConfig method*), 76  
**clean\_pre\_sample\_ratio()** (*TiePieOscilloscopeConfig method*), 81  
**clean\_probe\_offset()** (*TiePieOscilloscopeChannelConfig method*), 73  
**clean\_record\_length()** (*TiePieOscilloscopeConfig method*), 81  
**clean\_resolution()** (*TiePieOscilloscopeConfig static method*), 81  
**clean\_sample\_frequency()** (*TiePieOscilloscopeConfig method*), 81  
**clean\_signal\_type()** (*TiePieGeneratorConfig static method*), 76  
**clean\_trigger\_enabled()** (*TiePieOscilloscopeChannelConfig static method*), 73  
**clean\_trigger\_hysteresis()** (*TiePieOscilloscopeChannelConfig method*), 73  
**clean\_trigger\_kind()** (*TiePieOscilloscopeChannelConfig static method*), 73  
**clean\_trigger\_level()** (*TiePieOscilloscopeChannelConfig method*), 73  
**clean\_trigger\_level\_mode()** (*TiePieOscilloscopeChannelConfig static method*), 73  
**clean\_trigger\_timeout()** (*TiePieOscilloscopeConfig method*), 81  
**clean\_values()** (*AsyncCommunicationProtocolConfig method*), 13  
**clean\_values()** (*BaseCubeConfiguration method*), 36  
**clean\_values()** (*CryLasAttenuatorConfig method*), 86  
**clean\_values()** (*CryLasLaserConfig method*), 91  
**clean\_values()** (*EmptyConfig method*), 84, 179  
**clean\_values()** (*FuGConfig method*), 101  
**clean\_values()** (*HeinzingerConfig method*), 111  
**clean\_values()** (*ILS2TConfig method*), 163  
**clean\_values()** (*LaudaProRp245eConfig method*),  
 127  
**clean\_values()** (*LaudaProRp245eTcpCommunicationConfig method*), 128  
**clean\_values()** (*LJMCommunicationConfig method*),  
 17  
**clean\_values()** (*LuminoxConfig method*), 167  
**clean\_values()** (*MBW973Config method*), 130  
**clean\_values()** (*ModbusTcpCommunicationConfig method*), 19  
**clean\_values()** (*NewportSMC100PPConfig method*),  
 141  
**clean\_values()** (*OpcUaCommunicationConfig method*), 21  
**clean\_values()** (*PfeifferTPGConfig method*), 149  
**clean\_values()** (*PICubeConfiguration method*), 42  
**clean\_values()** (*PSI9000Config method*), 97  
**clean\_values()** (*SerialCommunicationConfig method*), 24  
**clean\_values()** (*SupercubeConfiguration method*), 53  
**clean\_values()** (*T560Config method*), 48  
**clean\_values()** (*TcpCommunicationConfig method*),  
 27  
**clean\_values()** (*TechnixConfig method*), 69  
**clean\_values()** (*TelnetCommunicationConfig method*), 29  
**clean\_values()** (*TiePieDeviceConfig method*), 71  
**clean\_values()** (*VisaCommunicationConfig method*),  
 31  
**clean\_waveform()** (*TiePieGeneratorConfig method*),  
 76  
**Client** (*class in hvl\_ccb.comm.opc*), 20  
**CLOSE** (*EarthingStickOperation attribute*), 39  
**close()** (*CommunicationProtocol method*), 14  
**close()** (*LaudaProRp245eTcpCommunication method*),  
 128  
**close()** (*LJMCommunication method*), 15  
**close()** (*ModbusTcpCommunication method*), 18  
**close()** (*NullCommunicationProtocol method*), 14  
**close()** (*OpcUaCommunication method*), 20  
**close()** (*SerialCommunication method*), 23  
**close()** (*Tcp method*), 26  
**close()** (*TelnetCommunication method*), 28  
**close()** (*VisaCommunication method*), 30  
**close\_shutter()** (*CryLasLaser method*), 89  
**CLOSED** (*CryLasLaser.AnswersShutter attribute*), 88  
**CLOSED** (*CryLasLaserShutterStatus attribute*), 94  
**CLOSED** (*DoorStatus attribute*), 38  
**closed** (*DoorStatus attribute*), 56  
**CLOSED** (*EarthingStickStatus attribute*), 39  
**closed** (*EarthingStickStatus attribute*), 57

```

CMD_EXEC_ERROR                                (New- config_cls() (AsyncCommunicationProtocol static
                                                portSMC100PPSerialCommunication.ControllerErrors method), 11
                                                attribute), 142
                                                config_cls() (BaseCube static method), 34
CMD_NOT_ALLOWED                               (New- config_cls() (BaseCubeOpcUaCommunication static
                                                portSMC100PPSerialCommunication.ControllerErrors method), 37
                                                attribute), 142
                                                config_cls() (ConfigurationMixin static method), 178
CMD_NOT_ALLOWED_CC                          (New- config_cls() (CryLasAttenuator static method), 86
                                                portSMC100PPSerialCommunication.ControllerErrors method), 87
                                                attribute), 142
                                                config_cls() (CryLasAttenuatorSerialCommunication
                                                               static method), 87
CMD_NOT_ALLOWED_CONFIGURATION      (New- config_cls() (CryLasLaser static method), 89
                                                portSMC100PPSerialCommunication.ControllerErrors method), 93
                                                attribute), 142
                                                config_cls() (CryLasLaserSerialCommunication
                                                               static method), 93
CMD_NOT_ALLOWED_DISABLE        (New- config_cls() (Device static method), 83
                                                portSMC100PPSerialCommunication.ControllerErrors method), 104
                                                attribute), 142
                                                config_cls() (FuGProbusIV static method), 104
CMD_NOT_ALLOWED_HOMING          (New- config_cls() (HeinzingerDI static method), 112
                                                portSMC100PPSerialCommunication.ControllerErrors method), 109
                                                attribute), 142
                                                config_cls() (HeinzingerSerialCommunication static
                                                               method), 109
CMD_NOT_ALLOWED_MOVING          (New- config_cls() (ILS2T static method), 161
                                                portSMC100PPSerialCommunication.ControllerErrors method), 115
                                                attribute), 142
                                                config_cls() (ILS2TModbusTcpCommunication static
                                                               method), 115
CMD_NOT_ALLOWED_NOT_REFERENCED (New- config_cls() (LaudaProRp245e static method), 123
                                                portSMC100PPSerialCommunication.ControllerErrors method), 164
                                                attribute), 142
                                                config_cls() (LaudaProRp245eTcpCommunication
                                                               static method), 123
CMD_NOT_ALLOWED_PP              (New- config_cls() (LJMCommunication static method), 15
                                                portSMC100PPSerialCommunication.ControllerErrors method), 128
                                                attribute), 142
                                                config_cls() (Luminox static method), 166
CMD_NOT_ALLOWED_READY          (New- config_cls() (LuminoxSerialCommunication static
                                                portSMC100PPSerialCommunication.ControllerErrors method), 169
                                                attribute), 142
                                                config_cls() (MBW973 static method), 129
CMR (PfeifferTPG.SensorTypes attribute), 147
CODE_OR_ADDR_INVALID            (New- config_cls() (ModbusTcpCommunication static
                                                portSMC100PPSerialCommunication.ControllerErrors method), 131
                                                attribute), 142
                                                config_cls() (ModbusTcpCommunication static
                                                               method), 18
collect_measurement_data()    (TiePieOscilloscope config_cls() (NewportSMC100PP static method), 135
                                method), 79
                                config_cls() (NewportSMC100PPSerialCommunication
                                                               static method), 143
com (SingleCommDevice property), 85
COM_TIME_OUT (LaudaProRp245eCommand attribute), 125
com_time_out (LaudaProRp245eConfig attribute), 127
COM_TIMEOUT (NewportSMC100PPSerialCommunication.ControllerErrors
             attribute), 142
command (LuminoxMeasurementType property), 168
COMMAND (TriggerMode attribute), 46
command() (FuGProbusIV method), 104
communication_channel (TechnixConfig attribute), 69
CommunicationException, 14
CommunicationProtocol          (class      in config_cls() (PSI9000 static method), 95
                                hvl_ccb.comm.base), 14
config (ConfigurationMixin property), 178
CONFIG (FuGProbusVRegisterGroups attribute), 107
CONFIG (NewportSMC100PP.StateMessages attribute), 134
CONFIG (NewportStates attribute), 145
                                                config_cls() (PfeifferTPG static method), 147
                                                config_cls() (PfeifferTPGSerialCommunication static
                                                               method), 150
                                                config_cls() (PICube static method), 41
                                                config_cls() (PICubeOpcUaCommunication static
                                                               method), 43
                                                config_cls() (PSI9000VisaCommunication static
                                                               method), 98
                                                config_cls() (RTO1024 static method), 153
                                                config_cls() (RTO1024VisaCommunication static
                                                               method), 159
                                                config_cls() (SerialCommunication static method), 23

```

`config_cls()` (*Supercube2015Base static method*), 49  
`config_cls()` (*SupercubeAOpCuaCommunication static method*), 63  
`config_cls()` (*SupercubeOpcUaCommunication static method*), 53  
`config_cls()` (*SyncCommunicationProtocol static method*), 14  
`config_cls()` (*T560 static method*), 47  
`config_cls()` (*T560Communication static method*), 45  
`config_cls()` (*Tcp static method*), 26  
`config_cls()` (*Technix static method*), 68  
`config_cls()` (*TechnixSerialCommunication static method*), 65  
`config_cls()` (*TechnixTelnetCommunication static method*), 66  
`config_cls()` (*TelnetCommunication static method*), 28  
`config_cls()` (*TiePieOscilloscope static method*), 79  
`config_cls()` (*VisaCommunication static method*), 30  
`config_cls()` (*VisaDevice static method*), 172  
`config_gen()` (*TiePieGeneratorMixin attribute*), 77  
`config_high_pulse()` (*LabJack method*), 119  
`config_i2c()` (*TiePieI2CHostMixin attribute*), 78  
`config_osc()` (*TiePieOscilloscope attribute*), 79  
`config_osc_channel_dict` (*TiePieOscilloscope attribute*), 79  
`config_status` (*FuG property*), 99  
`configdataclass()` (*in module hvl\_ccb.configuration*), 179  
`configuration_save_json()` (*ConfigurationMixin method*), 179  
`ConfigurationMixin` (*class in hvl\_ccb.configuration*), 178  
`ConfigurationValueWarning`, 179  
`connection_type` (*LJMCommunicationConfig attribute*), 17  
`CONT_MODE` (*LaudaProRp245eCommand attribute*), 125  
`continue_ramp()` (*LaudaProRp245e method*), 123  
`control_mode` (*LaudaProRp245eConfig attribute*), 127  
`CONVERSION` (*LEM4000S attribute*), 173  
`convert()` (*LEM4000S method*), 174  
`convert()` (*LMT70A method*), 174  
`convert()` (*Sensor method*), 174  
`convert()` (*Unit class method*), 175  
`ConvertibleTypes` (*in module hvl\_ccb.utils.typing*), 177  
`COOLOFF` (*LaudaProRp245eConfig.OperationModeEnum attribute*), 127  
`COOLON` (*LaudaProRp245eConfig.OperationModeEnum attribute*), 127  
`coupling` (*TiePieOscilloscopeChannelConfig property*), 73  
`CR` (*FuGTerminators attribute*), 110  
`create_serial_port()` (*SerialCommunicationConfig method*), 24  
`create_telnet()` (*TelnetCommunicationConfig method*), 29  
`CRLF` (*FuGTerminators attribute*), 110  
`CryLasAttenuator` (*class in hvl\_ccb.dev.crylas*), 86  
`CryLasAttenuatorConfig` (*class in hvl\_ccb.dev.crylas*), 86  
`CryLasAttenuatorError`, 87  
`CryLasAttenuatorSerialCommunication` (*class in hvl\_ccb.dev.crylas*), 87  
`CryLasAttenuatorSerialCommunicationConfig` (*class in hvl\_ccb.dev.crylas*), 87  
`CryLasLaser` (*class in hvl\_ccb.dev.crylas*), 88  
`CryLasLaser.AnswersShutter` (*class in hvl\_ccb.dev.crylas*), 88  
`CryLasLaser.AnswersStatus` (*class in hvl\_ccb.dev.crylas*), 88  
`CryLasLaser.LaserStatus` (*class in hvl\_ccb.dev.crylas*), 89  
`CryLasLaser.RepetitionRates` (*class in hvl\_ccb.dev.crylas*), 89  
`CryLasLaserConfig` (*class in hvl\_ccb.dev.crylas*), 91  
`CryLasLaserError`, 92  
`CryLasLaserNotReadyError`, 92  
`CryLasLaserPoller` (*class in hvl\_ccb.dev.crylas*), 92  
`CryLasLaserSerialCommunication` (*class in hvl\_ccb.dev.crylas*), 92  
`CryLasLaserShutterStatus` (*class in hvl\_ccb.dev.crylas*), 94  
`CubeEarthingStickOperationError`, 38  
`CubeError`, 38  
`CubeRemoteControlError`, 38  
`CubeStatusChangeError`, 38  
`CubeStopError`, 38  
`current` (*FuG property*), 99  
`CURRENT` (*FuGProbusIVCommands attribute*), 104  
`CURRENT` (*FuGReadbackChannels attribute*), 108  
`current` (*Technix property*), 68  
`current_lower_limit` (*PSI9000Config attribute*), 97  
`current_monitor` (*FuG property*), 100  
`current_primary` (*PICube property*), 41  
`current_primary` (*Power attribute*), 60  
`current_upper_limit` (*PSI9000Config attribute*), 97  
`CurrentPosition` (*New-portSMC100PPConfig.HomeSearch attribute*), 140  
`cv_mode` (*FuGProbusVDIRegisters property*), 106

## D

`datachange_notification()` (*OpcUaSubHandler method*), 22  
`datachange_notification()` (*SupercubeSubscriptionHandler method*), 54

date\_of\_manufacture (*LuminoxMeasurementType attribute*), 168  
 datetime\_to\_opc() (*BaseCube class method*), 34  
 DC (*TiePieGeneratorSignalType attribute*), 77  
 DC\_140KV (*PowerSetup attribute*), 40  
 DC\_280KV (*PowerSetup attribute*), 40  
 DC\_DoubleStage\_280kV (*PowerSetup attribute*), 61  
 DC\_SingleStage\_140kV (*PowerSetup attribute*), 61  
 DC\_VOLTAGE\_TOO\_LOW (*NewportSMC100PP.MotorErrors attribute*), 134  
 DCA (*TiePieOscilloscopeChannelCoupling attribute*), 74  
 DCV (*TiePieOscilloscopeChannelCoupling attribute*), 74  
 default\_com\_cls() (*BaseCube static method*), 34  
 default\_com\_cls() (*CryLasAttenuator static method*), 86  
 default\_com\_cls() (*CryLasLaser static method*), 89  
 default\_com\_cls() (*FuGProbusIV static method*), 104  
 default\_com\_cls() (*HeinzingerDI static method*), 112  
 default\_com\_cls() (*ILS2T static method*), 161  
 default\_com\_cls() (*LabJack static method*), 119  
 default\_com\_cls() (*LaudaProRp245e static method*), 123  
 default\_com\_cls() (*Luminox static method*), 167  
 default\_com\_cls() (*MBW973 static method*), 129  
 default\_com\_cls() (*NewportSMC100PP static method*), 135  
 default\_com\_cls() (*PfeifferTPG static method*), 147  
 default\_com\_cls() (*PICube static method*), 41  
 default\_com\_cls() (*PSI9000 static method*), 95  
 default\_com\_cls() (*RTO1024 static method*), 153  
 default\_com\_cls() (*SingleCommDevice static method*), 85  
 default\_com\_cls() (*Supercube2015Base static method*), 50  
 default\_com\_cls() (*Supercube2015WithFU static method*), 62  
 default\_com\_cls() (*T560 static method*), 47  
 default\_com\_cls() (*Technix method*), 68  
 default\_com\_cls() (*TiePieOscilloscope static method*), 79  
 default\_com\_cls() (*VisaDevice static method*), 172  
 DEFAULT\_IO\_SCANNING\_CONTROL\_VALUES (*ILS2T attribute*), 160  
 default\_n\_attempts\_read\_text\_nonempty (*AsyncCommunicationProtocolConfig attribute*), 13  
 default\_n\_attempts\_read\_text\_nonempty (*FuGSerialCommunicationConfig attribute*), 109  
 default\_n\_attempts\_read\_text\_nonempty (*HeinzingerSerialCommunicationConfig attribute*), 115  
 default\_number\_of\_recordings (*HeinzingerConfig attribute*), 111  
 Device (*class in hvl\_ccb.dev.base*), 83  
 DEVICE\_TYPE (*LaudaProRp245eCommand attribute*), 125  
 device\_type (*LJMCommunicationConfig attribute*), 17  
 DeviceException, 83  
 DeviceExistingException, 83  
 DeviceFailuresException, 83  
 devices\_failed\_start (*DeviceSequenceMixin attribute*), 84  
 devices\_failed\_start (*ExperimentManager attribute*), 181  
 devices\_failed\_stop (*DeviceSequenceMixin attribute*), 84  
 devices\_failed\_stop (*ExperimentManager attribute*), 181  
 DeviceSequenceMixin (*class in hvl\_ccb.dev.base*), 83  
 di (*FuG property*), 100  
 digital\_control (*FuGProbusVDIRegisters property*), 106  
 DIOChannel (*LabJack attribute*), 118  
 DISABLE (*NewportStates attribute*), 145  
 disable() (*ILS2T method*), 161  
 DISABLE\_FROM\_JOGGING (*NewportSMC100PP.StateMessages attribute*), 134  
 DISABLE\_FROM\_MOVING (*NewportSMC100PP.StateMessages attribute*), 134  
 DISABLE\_FROM\_READY (*NewportSMC100PP.StateMessages attribute*), 134  
 disable\_pulses() (*LabJack method*), 119  
 DISABLED (*TiePieOscilloscopeAutoResolutionModes attribute*), 80  
 DisableEspStageCheck (*NewportSMC100PPConfig.EspStageConfig attribute*), 140  
 disarm\_trigger() (*T560 method*), 48  
 disconnect() (*Client method*), 20  
 DISPLACEMENT\_OUT\_OF\_LIMIT (*NewportSMC100PPSerialCommunication.ControllerErrors attribute*), 142  
 display\_message\_board() (*BaseCube method*), 34  
 display\_status\_board() (*BaseCube method*), 34  
 do\_iostressing\_write() (*ILS2T method*), 161  
 door\_1\_status (*BaseCube attribute*), 34  
 door\_2\_status (*BaseCube attribute*), 34  
 door\_3\_status (*BaseCube attribute*), 34  
 DoorStatus (*class in hvl\_ccb.dev.cube.constants*), 38  
 DoorStatus (*class in hvl\_ccb.dev.supercube2015.constants*), 56

## E

E (*LabJack.ThermocoupleType attribute*), 119  
 E0 (*FuGErrorcodes attribute*), 102

E1 (*FuGErrorcodes* attribute), 102  
E10 (*FuGErrorcodes* attribute), 102  
E100 (*FuGErrorcodes* attribute), 102  
E106 (*FuGErrorcodes* attribute), 102  
E11 (*FuGErrorcodes* attribute), 102  
E115 (*FuGErrorcodes* attribute), 102  
E12 (*FuGErrorcodes* attribute), 102  
E125 (*FuGErrorcodes* attribute), 102  
E13 (*FuGErrorcodes* attribute), 102  
E135 (*FuGErrorcodes* attribute), 102  
E14 (*FuGErrorcodes* attribute), 102  
E145 (*FuGErrorcodes* attribute), 102  
E15 (*FuGErrorcodes* attribute), 102  
E155 (*FuGErrorcodes* attribute), 102  
E16 (*FuGErrorcodes* attribute), 103  
E165 (*FuGErrorcodes* attribute), 103  
E2 (*FuGErrorcodes* attribute), 103  
E206 (*FuGErrorcodes* attribute), 103  
E306 (*FuGErrorcodes* attribute), 103  
E4 (*FuGErrorcodes* attribute), 103  
E5 (*FuGErrorcodes* attribute), 103  
E504 (*FuGErrorcodes* attribute), 103  
E505 (*FuGErrorcodes* attribute), 103  
E6 (*FuGErrorcodes* attribute), 103  
E666 (*FuGErrorcodes* attribute), 103  
E7 (*FuGErrorcodes* attribute), 103  
E8 (*FuGErrorcodes* attribute), 103  
E9 (*FuGErrorcodes* attribute), 103  
earthing\_rod\_1\_status (*BaseCube* attribute), 34  
earthing\_rod\_2\_status (*BaseCube* attribute), 34  
earthing\_rod\_3\_status (*BaseCube* attribute), 34  
EarthingRodStatus (class *hvl\_ccb.dev(cube.constants)*), 39  
EarthingStick (class *hvl\_ccb.dev.supercube2015.constants*), 56  
EarthingStickOperatingStatus (class *hvl\_ccb.dev(cube.constants)*), 39  
EarthingStickOperation (class *hvl\_ccb.dev(cube.constants)*), 39  
EarthingStickStatus (class *hvl\_ccb.dev(cube.constants)*), 39  
EarthingStickStatus (class *hvl\_ccb.dev.supercube2015.constants*), 57  
EEPROM\_ACCESS\_ERROR (New-*portSMC100PPSerialCommunication.ControllerErrors* attribute), 142  
EIGHT (*HeinzingerConfig.RecordingsEnum* attribute), 111  
EIGHT\_BIT (*TiePieOscilloscopeResolution* attribute), 81  
EIGHT\_HUNDRED\_MILLI\_VOLT (*TiePieOscilloscopeRange* attribute), 74  
EIGHT\_VOLT (*TiePieOscilloscopeRange* attribute), 74  
EIGHTBITS (*SerialCommunicationBytesize* attribute), 24  
EIGHTY\_VOLT (*TiePieOscilloscopeRange* attribute), 74  
EmptyConfig (class in *hvl\_ccb.configuration*), 179  
EmptyConfig (class in *hvl\_ccb.dev.base*), 84  
enable() (*ILS2T* method), 161  
enable\_clock() (*LabJack* method), 119  
enabled (*TiePieGeneratorConfig* property), 76  
enabled (*TiePieOscilloscopeChannelConfig* property), 73  
EnableEspStageCheck (New-*portSMC100PPConfig.EspStageConfig* attribute), 140  
encoding (*AsyncCommunicationProtocolConfig* attribute), 13  
encoding (*NewportSMC100PPSerialCommunicationConfig* attribute), 145  
encoding\_error\_handling (*AsyncCommunicationProtocolConfig* attribute), 13  
encoding\_error\_handling (New-*portSMC100PPSerialCommunicationConfig* attribute), 145  
EndOfRunSwitch (New-*portSMC100PPConfig.HomeSearch* attribute), 140  
EndOfRunSwitch\_and\_Index (New-*portSMC100PPConfig.HomeSearch* attribute), 140  
endpoint\_name (*BaseCubeOpcUaCommunicationConfig* attribute), 37  
endpoint\_name (*OpcUaCommunicationConfig* attribute), 21  
endpoint\_name (*PICubeOpcUaCommunicationConfig* attribute), 43  
in endpoint\_name (*SupercubeAOpcUaConfiguration* attribute), 64  
in ERROR (*DoorStatus* attribute), 38  
error (*DoorStatus* attribute), 56  
in ERROR (*EarthingStickStatus* attribute), 39  
error (*EarthingStickStatus* attribute), 58  
in ERROR (*ExperimentStatus* attribute), 181  
ERROR (*SafetyStatus* attribute), 40  
in Error (*SafetyStatus* attribute), 61  
errorcode (*FuGError* attribute), 102  
in Errors (class in *hvl\_ccb.dev.supercube2015.constants*), 58  
ESP\_STAGE\_NAME\_INVALID (New-*portSMC100PPSerialCommunication.ControllerErrors* attribute), 142  
ETH (*LaudaProRp245eConfig.ExtControlModeEnum* attribute), 126  
ETHERNET (*LJMCommunicationConfig.ConnectionType* attribute), 16  
EVEN (*SerialCommunicationParity* attribute), 25  
event\_notification() (*OpcUaSubHandler* method), 22  
EXECUTE (*FuGProbusIVCommands* attribute), 104

```

execute_absolute_position() (ILS2T method), 161
execute_on_x (FuGProbusVConfigRegisters property), 105
execute_relative_step() (ILS2T method), 161
EXECUTEONX (FuGProbusIVCommands attribute), 104
exit_configuration() (NewportSMC100PP method), 135
exit_configuration_wait_sec (NewportSMC100PPConfig attribute), 141
EXPERIMENT_BLOCKED (EarthingRodStatus attribute), 39
EXPERIMENT_READY (EarthingRodStatus attribute), 39
ExperimentError, 180
ExperimentManager (class) in hvl_ccb.experiment_manager, 180
ExperimentStatus (class) in hvl_ccb.experiment_manager, 181
EXPT100 (LaudaProRp245eConfig.ExtControlModeEnum attribute), 126
EXT_FALLING_EDGE (TriggerMode attribute), 46
EXT_RISING_EDGE (TriggerMode attribute), 46
External (PowerSetup attribute), 61
EXTERNAL_SOURCE (PowerSetup attribute), 40
EXTERNAL_TEMP (LaudaProRp245eCommand attribute), 125

F
F (LabJack.TemperatureUnit attribute), 118
F (Temperature attribute), 175
FAHRENHEIT (Temperature attribute), 175
failures (DeviceFailuresException attribute), 83
FALLING (TiePieOscilloscopeTriggerKind attribute), 75
FAST (ILS2T.Ref16Jog attribute), 160
file_copy() (RTO1024 method), 153
finish() (ExperimentManager method), 181
FINISHED (ExperimentStatus attribute), 181
FINISHING (ExperimentStatus attribute), 181
fire_trigger() (T560 method), 48
FIRMWARE (FuGReadbackChannels attribute), 108
FIVE_MHZ (LabJack.ClockFrequency attribute), 117
FIVEBITS (SerialCommunicationBytesize attribute), 24
FLT_INFO (ILS2TRegAddr attribute), 165
FLT_MEM_DEL (ILS2TRegAddr attribute), 165
FLT_MEM_RESET (ILS2TRegAddr attribute), 165
FOLLOWING_ERROR (NewportSMC100PP.MotorErrors attribute), 134
FOLLOWRAMP (FuGRampModes attribute), 108
force_trigger() (TiePieOscilloscope method), 80
force_value() (AsyncCommunicationProtocolConfig method), 13
force_value() (BaseCubeConfiguration method), 36
force_value() (BaseCubeOpcUaCommunicationConfig method), 37
force_value() (CryLasAttenuatorConfig method), 86
force_value() (CryLasAttenuatorSerialCommunicationConfig method), 87
force_value() (CryLasLaserConfig method), 91
force_value() (CryLasLaserSerialCommunicationConfig method), 94
force_value() (EmptyConfig method), 84, 179
force_value() (FuGConfig method), 101
force_value() (FuGSerialCommunicationConfig method), 109
force_value() (HeinzingerConfig method), 111
force_value() (HeinzingerSerialCommunicationConfig method), 115
force_value() (ILS2TConfig method), 163
force_value() (ILS2TModbusTcpCommunicationConfig method), 164
force_value() (LaudaProRp245eConfig method), 127
force_value() (LaudaProRp245eTcpCommunicationConfig method), 128
force_value() (LJMCommunicationConfig method), 17
force_value() (LuminoxConfig method), 167
force_value() (LuminoxSerialCommunicationConfig method), 170
force_value() (MBW973Config method), 131
force_value() (MBW973SerialCommunicationConfig method), 132
force_value() (ModbusTcpCommunicationConfig method), 19
force_value() (NewportSMC100PPConfig method), 141
force_value() (NewportSMC100PPSerialCommunicationConfig method), 145
force_value() (OpcUaCommunicationConfig method), 21
force_value() (PfeifferTPGConfig method), 149
force_value() (PfeifferTPGSerialCommunicationConfig method), 151
force_value() (PICubeConfiguration method), 42
force_value() (PICubeOpcUaCommunicationConfig method), 43
force_value() (PSI9000Config method), 97
force_value() (PSI9000VisaCommunicationConfig method), 98
force_value() (RTO1024Config method), 158
force_value() (RTO1024VisaCommunicationConfig method), 159
force_value() (SerialCommunicationConfig method), 24
force_value() (SupercubeAOpcUaConfiguration method), 64
force_value() (SupercubeConfiguration method), 53
force_value() (SupercubeOpcUaCommunicationConfig method), 54

```

**force\_value()** (*T560CommunicationConfig method*), 45  
**force\_value()** (*T560Config method*), 48  
**force\_value()** (*TcpCommunicationConfig method*), 27  
**force\_value()** (*TechnixConfig method*), 69  
**force\_value()** (*TechnixSerialCommunicationConfig method*), 66  
**force\_value()** (*TechnixTelnetCommunicationConfig method*), 67  
**force\_value()** (*TelnetCommunicationConfig method*), 29  
**force\_value()** (*TiePieDeviceConfig method*), 71  
**force\_value()** (*VisaCommunicationConfig method*), 31  
**force\_value()** (*VisaDeviceConfig method*), 172  
**FORTY\_MHZ** (*LabJack.ClockFrequency attribute*), 118  
**FORTY\_VOLT** (*TiePieOscilloscopeRange attribute*), 74  
**FOUR** (*HeinzingerConfig.RecordingsEnum attribute*), 111  
**FOUR\_HUNDRED\_MILLI\_VOLT** (*TiePieOscilloscopeRange attribute*), 74  
**FOUR\_VOLT** (*TiePieOscilloscopeRange attribute*), 74  
**FOURTEEN\_BIT** (*TiePieOscilloscopeResolution attribute*), 81  
**FREERUN** (*RTO1024.TriggerModes attribute*), 153  
**frequency** (*PICube property*), 41  
**frequency** (*Power attribute*), 60  
**frequency** (*T560 property*), 48  
**frequency** (*TiePieGeneratorConfig property*), 76  
**FRM** (*NewportConfigCommands attribute*), 133  
**from\_json()** (*ConfigurationMixin class method*), 179  
**FRS** (*NewportConfigCommands attribute*), 133  
**fso\_reset()** (*Supercube2015WithFU method*), 63  
**FuG** (*class in hvl\_ccb.dev.fug*), 99  
**FuGConfig** (*class in hvl\_ccb.dev.fug*), 101  
**FuGDigitalVal** (*class in hvl\_ccb.dev.fug*), 101  
**FuGError**, 102  
**FuGErrorcodes** (*class in hvl\_ccb.dev.fug*), 102  
**FuGMonitorModes** (*class in hvl\_ccb.dev.fug*), 103  
**FuGPolarities** (*class in hvl\_ccb.dev.fug*), 104  
**FuGProbusIV** (*class in hvl\_ccb.dev.fug*), 104  
**FuGProbusIVCommands** (*class in hvl\_ccb.dev.fug*), 104  
**FuGProbusV** (*class in hvl\_ccb.dev.fug*), 105  
**FuGProbusVConfigRegisters** (*class in hvl\_ccb.dev.fug*), 105  
**FuGProbusVDIRegisters** (*class in hvl\_ccb.dev.fug*), 106  
**FuGProbusVDORegisters** (*class in hvl\_ccb.dev.fug*), 106  
**FuGProbusVMonitorRegisters** (*class in hvl\_ccb.dev.fug*), 107  
**FuGProbusVRegisterGroups** (*class in hvl\_ccb.dev.fug*), 107  
**FuGProbusVSetRegisters** (*class in hvl\_ccb.dev.fug*), 107  
**FuGRampModes** (*class in hvl\_ccb.dev.fug*), 108  
**FuGReadbackChannels** (*class in hvl\_ccb.dev.fug*), 108  
**FuGSerialCommunication** (*class in hvl\_ccb.dev.fug*), 108  
**FuGSerialCommunicationConfig** (*class in hvl\_ccb.dev.fug*), 109  
**FuGTerminators** (*class in hvl\_ccb.dev.fug*), 110

## G

**gate\_mode** (*T560 property*), 48  
**gate\_polarity** (*T560 property*), 48  
**GateMode** (*class in hvl\_ccb.dev.hIGHLAND\_T560.base*), 44  
**GeneralSockets** (*class in hvl\_ccb.dev.supercube2015.constants*), 58  
**GeneralSupport** (*class in hvl\_ccb.dev.supercube2015.constants*), 58  
**GENERATOR** (*TiePieDeviceType attribute*), 72  
**generator\_start()** (*TiePieGeneratorMixin method*), 77  
**generator\_stop()** (*TiePieGeneratorMixin method*), 77  
**get()** (*AlarmText class method*), 55  
**get()** (*MeasurementsDividerRatio class method*), 59  
**get()** (*MeasurementsScaledInput class method*), 60  
**get\_acceleration()** (*NewportSMC100PP method*), 135  
**get\_acquire\_length()** (*RTO1024 method*), 153  
**get\_ain()** (*LabJack method*), 119  
**get\_bath\_temp()** (*LaudaProRp245e method*), 123  
**get\_by\_p\_id()** (*LabJack.DeviceType class method*), 118  
**get\_by\_p\_id()** (*LJMCommunicationConfig.DeviceType class method*), 16  
**get\_cal\_current\_source()** (*LabJack method*), 119  
**get\_cee16\_socket()** (*Supercube2015Base method*), 50  
**get\_channel\_offset()** (*RTO1024 method*), 153  
**get\_channel\_position()** (*RTO1024 method*), 154  
**get\_channel\_range()** (*RTO1024 method*), 154  
**get\_channel\_scale()** (*RTO1024 method*), 154  
**get\_channel\_state()** (*RTO1024 method*), 154  
**get\_clock()** (*LabJack method*), 120  
**get\_controller\_information()** (*NewportSMC100PP method*), 135  
**get\_current()** (*HeinzingerDI method*), 112  
**get\_dc\_volt()** (*ILS2T method*), 162  
**get\_device()** (*DeviceSequenceMixin method*), 84  
**get\_device\_by\_serial\_number()** (*in module hvl\_ccb.dev.tiepie.base*), 72  
**get\_device\_type()** (*LaudaProRp245e method*), 123  
**get\_devices()** (*DeviceSequenceMixin method*), 84  
**get\_digital\_input()** (*LabJack method*), 120  
**get\_door\_status()** (*Supercube2015Base method*), 50  
**get\_earthing\_manual()** (*Supercube2015Base method*), 50

<code>get_earthing_status()</code>	<i>(Supercube2015Base method)</i> , 50	<code>get_status()</code>	<i>(Supercube2015Base method)</i> , 50
<code>get_error_code()</code>	<i>(ILS2T method)</i> , 162	<code>get_support_input()</code>	<i>(Supercube2015Base method)</i> , 50
<code>get_error_queue()</code>	<i>(VisaDevice method)</i> , 172	<code>get_support_output()</code>	<i>(Supercube2015Base method)</i> , 51
<code>get_frequency()</code>	<i>(Supercube2015WithFU method)</i> , 63	<code>get_system_lock()</code>	<i>(PSI9000 method)</i> , 95
<code>get_fso_active()</code>	<i>(Supercube2015WithFU method)</i> , 63	<code>get_t13_socket()</code>	<i>(Supercube2015Base method)</i> , 51
<code>get_full_scale_mbar()</code>	<i>(PfeifferTPG method)</i> , 147	<code>get_target_voltage()</code>	<i>(Supercube2015WithFU method)</i> , 63
<code>get_full_scale_unitless()</code>	<i>(PfeifferTPG method)</i> , 148	<code>get_temperature()</code>	<i>(ILS2T method)</i> , 162
<code>get_identification()</code>	<i>(VisaDevice method)</i> , 172	<code>get_timestamps()</code>	<i>(RTO1024 method)</i> , 154
<code>get_interface_version()</code>	<i>(HeinzingerDI method)</i> , 112	<code>get_ui_lower_limits()</code>	<i>(PSI9000 method)</i> , 96
<code>get_max_voltage()</code>	<i>(Supercube2015WithFU method)</i> , 63	<code>get_uip_upper_limits()</code>	<i>(PSI9000 method)</i> , 96
<code>get_measurement_ratio()</code>	<i>(Supercube2015Base method)</i> , 50	<code>get_voltage()</code>	<i>(HeinzingerDI method)</i> , 112
<code>get_measurement_voltage()</code>	<i>(Supercube2015Base method)</i> , 50	<code>get_voltage_current_setpoint()</code>	<i>(PSI9000 method)</i> , 96
<code>get_motor_configuration()</code>	<i>(NewportSMC100PP method)</i> , 136	<code>go_home()</code>	<i>(NewportSMC100PP method)</i> , 137
<code>get_move_duration()</code>	<i>(NewportSMC100PP method)</i> , 136	<code>go_to_configuration()</code>	<i>(NewportSMC100PP method)</i> , 137
<code>get_negative_software_limit()</code>	<i>(NewportSMC100PP method)</i> , 136	<code>GREEN_NOT_READY</code>	<i>(SafetyStatus attribute)</i> , 40
<code>get_number_of_recordings()</code>	<i>(HeinzingerDI method)</i> , 112	<code>GREEN_READY</code>	<i>(SafetyStatus attribute)</i> , 41
<code>get_objects_node()</code>	<i>(Client method)</i> , 20	<code>GreenNotReady</code>	<i>(SafetyStatus attribute)</i> , 61
<code>get_objects_node()</code>	<i>(Server method)</i> , 22	<code>GreenReady</code>	<i>(SafetyStatus attribute)</i> , 61
<code>get_output()</code>	<i>(PSI9000 method)</i> , 95	<b>H</b>	
<code>get_position()</code>	<i>(ILS2T method)</i> , 162	<code>HARDWARE</code>	<i>(CryLasLaser.RepetitionRates attribute)</i> , 89
<code>get_position()</code>	<i>(NewportSMC100PP method)</i> , 136	<code>has_safe_ground</code>	<i>(TiePieOscilloscopeChannelConfig property)</i> , 73
<code>get_positive_software_limit()</code>	<i>(NewportSMC100PP method)</i> , 136	<code>HEAD</code>	<i>(CryLasLaser.AnswersStatus attribute)</i> , 88
<code>get_power_setup()</code>	<i>(Supercube2015WithFU method)</i> , 63	<code>HeinzingerConfig</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 110
<code>get_primary_current()</code>	<i>(Supercube2015WithFU method)</i> , 63	<code>HeinzingerConfig.RecordingsEnum</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 111
<code>get_primary_voltage()</code>	<i>(Supercube2015WithFU method)</i> , 63	<code>HeinzingerDI</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 111
<code>get_product_id()</code>	<i>(LabJack method)</i> , 120	<code>HeinzingerDI.OutputStatus</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 112
<code>get_product_name()</code>	<i>(LabJack method)</i> , 120	<code>HeinzingerPNC</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 113
<code>get_product_type()</code>	<i>(LabJack method)</i> , 120	<code>HeinzingerPNC.UnitCurrent</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 113
<code>get_pulse_energy_and_rate()</code>	<i>(CryLasLaser method)</i> , 89	<code>HeinzingerPNC.UnitVoltage</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 113
<code>get_reference_point()</code>	<i>(RTO1024 method)</i> , 154	<code>HeinzingerPNCDeviceNotRecognizedException</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 114
<code>get_register()</code>	<i>(FuGProbusV method)</i> , 105	<code>HeinzingerPNCError</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 114
<code>get_repetitions()</code>	<i>(RTO1024 method)</i> , 154	<code>HeinzingerPNCMaxCurrentExceededException</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 114
<code>get_sbus_rh()</code>	<i>(LabJack method)</i> , 120	<code>HeinzingerPNCMaxVoltageExceededException</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 114
<code>get_sbus_temp()</code>	<i>(LabJack method)</i> , 120	<code>HeinzingerSerialCommunication</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 114
<code>get_serial_number()</code>	<i>(HeinzingerDI method)</i> , 112	<code>HeinzingerSerialCommunicationConfig</code>	<i>(class in hvl_ccb.dev.heinzinger)</i> , 115
<code>get_serial_number()</code>	<i>(LabJack method)</i> , 120	<code>HIGH</code>	<i>(LabJack.DIOStatus attribute)</i> , 118
<code>get_state()</code>	<i>(NewportSMC100PP method)</i> , 137	<code>high_resolution</code>	<i>(FuGProbusVSetRegisters property)</i> , 107
<code>get_status()</code>	<i>(ILS2T method)</i> , 162		

home\_search\_polling\_interval (*NewportSMC100PPConfig attribute*), 141  
 home\_search\_timeout (*NewportSMC100PPConfig attribute*), 141  
 home\_search\_type (*NewportSMC100PPConfig attribute*), 141  
 home\_search\_velocity (*NewportSMC100PPConfig attribute*), 141  
 HOME\_STARTED (*NewportSMC100PPSerialCommunication attribute*), 142  
 HomeSwitch (*NewportSMC100PPConfig.HomeSearch attribute*), 140  
 HomeSwitch\_and\_Index (*NewportSMC100PPConfig.HomeSearch attribute*), 140  
 HOMING (*NewportStates attribute*), 145  
 HOMING\_FROM\_RS232 (*NewportSMC100PP.StateMessages attribute*), 134  
 HOMING\_FROM\_SMC (*NewportSMC100PP.StateMessages attribute*), 134  
 HOMING\_TIMEOUT (*NewportSMC100PP.MotorErrors attribute*), 134  
 horn (*Safety attribute*), 61  
 horn() (*Supercube2015Base method*), 51  
 host (*ModbusTcpCommunicationConfig attribute*), 19  
 host (*OpcUaCommunicationConfig attribute*), 21  
 host (*TcpCommunicationConfig attribute*), 27  
 host (*TelnetCommunicationConfig attribute*), 29  
 host (*VisaCommunicationConfig attribute*), 32  
 hPascal (*PfeifferTPG.PressureUnits attribute*), 147  
 HT (*NewportConfigCommands attribute*), 133  
 hvl\_ccb  
     module, 182  
 hvl\_ccb.comm  
     module, 33  
 hvl\_ccb.comm.base  
     module, 11  
 hvl\_ccb.comm.labjack\_ljm  
     module, 15  
 hvl\_ccb.comm.modbus\_tcp  
     module, 18  
 hvl\_ccb.comm.opc  
     module, 20  
 hvl\_ccb.comm.serial  
     module, 23  
 hvl\_ccb.comm.tcp  
     module, 26  
 hvl\_ccb.comm.telnet  
     module, 28  
 hvl\_ccb.comm.visa  
     module, 30  
 hvl\_ccb.configuration  
     module, 178  
 (New-  
     module, 173  
 hvl\_ccb.dev  
     module, 83  
 hvl\_ccb.dev.base  
     module, 44  
 hvl\_ccb.dev.crylas  
     module, 85  
 hvl\_ccb.dev.cube  
     module, 33  
 hvl\_ccb.dev.cube.constants  
     module, 38  
 hvl\_ccb.dev.cube.picube  
     module, 41  
 hvl\_ccb.dev.ea\_psi9000  
     module, 95  
 hvl\_ccb.dev.fug  
     module, 99  
 hvl\_ccb.dev.heinzinger  
     module, 110  
 hvl\_ccb.dev.highland\_t560  
     module, 49  
 hvl\_ccb.dev.highland\_t560.base  
     module, 44  
 hvl\_ccb.dev.highland\_t560.channel  
     module, 46  
 hvl\_ccb.dev.highland\_t560.device  
     module, 47  
 hvl\_ccb.dev.labjack  
     module, 116  
 hvl\_ccb.dev.lauda  
     module, 123  
 hvl\_ccb.dev.mbw973  
     module, 129  
 hvl\_ccb.dev.newport  
     module, 133  
 hvl\_ccb.dev.pfeiffer\_tpg  
     module, 146  
 hvl\_ccb.dev.rs\_rto1024  
     module, 152  
 hvl\_ccb.dev.se\_ilis2t  
     module, 160  
 hvl\_ccb.dev.sst\_luminox  
     module, 166  
 hvl\_ccb.dev.supercube2015  
     module, 64  
 hvl\_ccb.dev.supercube2015.base  
     module, 49  
 hvl\_ccb.dev.supercube2015.constants  
     module, 55  
 hvl\_ccb.dev.supercube2015.typ\_a  
     module, 62  
 hvl\_ccb.dev.technix  
     module, 70

```

hvl_ccb.dev.technix.base
    module, 65
hvl_ccb.dev.technix.device
    module, 68
hvl_ccb.dev.tiepie
    module, 82
hvl_ccb.dev.tiepie.base
    module, 71
hvl_ccb.dev.tiepie.channel
    module, 73
hvl_ccb.dev.tiepie.device
    module, 75
hvl_ccb.dev.tiepie.generator
    module, 76
hvl_ccb.dev.tiepie.i2c
    module, 78
hvl_ccb.dev.tiepie.oscilloscope
    module, 79
hvl_ccb.dev.tiepie.utils
    module, 81
hvl_ccb.dev.utils
    module, 170
hvl_ccb.dev.visa
    module, 171
hvl_ccb.exception
    module, 180
hvl_ccb.experiment_manager
    module, 180
hvl_ccb.utils
    module, 178
hvl_ccb.utils.conversion_sensor
    module, 173
hvl_ccb.utils.conversion_unit
    module, 174
hvl_ccb.utils.enum
    module, 176
hvl_ccb.utils.typing
    module, 177
hvl_ccb.utils.validation
    module, 177
hysteresis_compensation (New-
    portSMC100PPConfig attribute), 141
| I2C (TiePieDeviceType attribute), 72
ID (FuGProbusIVCommands attribute), 105
Identification_error (PfeifferTPG.SensorStatus attribute), 147
identifier (LJMCommunicationConfig attribute), 17
identify_device() (FuG method), 100
identify_device() (HeinzingerPNC method), 114
identify_sensors() (PfeifferTPG method), 148
IKR (PfeifferTPG.SensorTypes attribute), 147
IKR11 (PfeifferTPG.SensorTypes attribute), 147
IKR9 (PfeifferTPG.SensorTypes attribute), 147
ILS2T (class in hvl_ccb.dev.se_ilst), 160
ILS2T.ActionsPtp (class in hvl_ccb.dev.se_ilst), 160
ILS2T.Mode (class in hvl_ccb.dev.se_ilst), 160
ILS2T.Ref16Jog (class in hvl_ccb.dev.se_ilst), 160
ILS2T.State (class in hvl_ccb.dev.se_ilst), 161
ILS2TConfig (class in hvl_ccb.dev.se_ilst), 163
ILS2TException, 164
ILS2TModbusTcpCommunication (class in hvl_ccb.dev.se_ilst), 164
ILS2TModbusTcpCommunicationConfig (class in hvl_ccb.dev.se_ilst), 164
ILS2TRegAddr (class in hvl_ccb.dev.se_ilst), 165
ILS2TRegDatatype (class in hvl_ccb.dev.se_ilst), 165
IMMEDIATELY (FuGRampModes attribute), 108
IMPULSE_140KV (PowerSetup attribute), 40
IMR (PfeifferTPG.SensorTypes attribute), 147
in_1_1 (GeneralSupport attribute), 58
in_1_2 (GeneralSupport attribute), 58
in_2_1 (GeneralSupport attribute), 58
in_2_2 (GeneralSupport attribute), 58
in_3_1 (GeneralSupport attribute), 58
in_3_2 (GeneralSupport attribute), 58
in_4_1 (GeneralSupport attribute), 58
in_4_2 (GeneralSupport attribute), 58
in_5_1 (GeneralSupport attribute), 58
in_5_2 (GeneralSupport attribute), 58
in_6_1 (GeneralSupport attribute), 58
in_6_2 (GeneralSupport attribute), 59
INACTIVE (CryLasLaser.AnswersStatus attribute), 88
INACTIVE (DoorStatus attribute), 38
inactive (DoorStatus attribute), 56
INACTIVE (EarthingStickStatus attribute), 39
inactive (EarthingStickStatus attribute), 58
inhibit (Technix property), 68
init_attenuation (CryLasAttenuatorConfig attribute), 86
init_monitored_nodes() (OpcUaCommunication method), 20
init_shutter_status (CryLasLaserConfig attribute), 92
initialize() (NewportSMC100PP method), 137
INITIALIZED (ExperimentStatus attribute), 181
INITIALIZING (ExperimentStatus attribute), 181
INITIALIZING (SafetyStatus attribute), 41
Initializing (SafetyStatus attribute), 61
INPUT (FuGProbusVRegisterGroups attribute), 107
INPUT (GateMode attribute), 44
input() (GeneralSupport class method), 59
input_1 (MeasurementsDividerRatio attribute), 59
input_1 (MeasurementsScaledInput attribute), 60
input_2 (MeasurementsScaledInput attribute), 60
input_3 (MeasurementsScaledInput attribute), 60
input_4 (MeasurementsScaledInput attribute), 60

```

input\_range (*TiePieOscilloscopeChannelConfig* property), 74  
 INSTALL (*AutoInstallMode* attribute), 44  
 INT32 (*ILS2TRegDatatype* attribute), 166  
 INT\_SYNTHESIZER (*TriggerMode* attribute), 46  
 interface\_type (*PSI9000VisaCommunicationConfig* attribute), 98  
 interface\_type (*RTO1024VisaCommunicationConfig* attribute), 159  
 interface\_type (*VisaCommunicationConfig* attribute), 32  
 internal (*LabJack.CjcType* attribute), 117  
 INTERNAL (*LaudaProRp245eConfig.ExtControlModeEnum* attribute), 126  
 Internal (*PowerSetup* attribute), 61  
 InvalidSupercubeStatusError, 49  
 IO\_SCANNING (*ILS2TRegAddr* attribute), 165  
 IoScanningModeValueError, 166  
 is\_configdataclass (*AsyncCommunicationProtocolConfig* attribute), 13  
 is\_configdataclass (*BaseCubeConfiguration* attribute), 36  
 is\_configdataclass (*CryLasAttenuatorConfig* attribute), 86  
 is\_configdataclass (*CryLasLaserConfig* attribute), 92  
 is\_configdataclass (*EmptyConfig* attribute), 85, 179  
 is\_configdataclass (*FuGConfig* attribute), 101  
 is\_configdataclass (*HeinzingerConfig* attribute), 111  
 is\_configdataclass (*ILS2TConfig* attribute), 164  
 is\_configdataclass (*LaudaProRp245eConfig* attribute), 127  
 is\_configdataclass (*LJMCommunicationConfig* attribute), 17  
 is\_configdataclass (*LuminoxConfig* attribute), 168  
 is\_configdataclass (*MBW973Config* attribute), 131  
 is\_configdataclass (*ModbusTcpCommunicationConfig* attribute), 19  
 is\_configdataclass (*NewportSMC100PPConfig* attribute), 141  
 is\_configdataclass (*OpcUaCommunicationConfig* attribute), 21  
 is\_configdataclass (*PfeifferTPGConfig* attribute), 149  
 is\_configdataclass (*SupercubeConfiguration* attribute), 53  
 is\_configdataclass (*T560Config* attribute), 48  
 is\_configdataclass (*TcpCommunicationConfig* attribute), 27  
 is\_configdataclass (*TechnixConfig* attribute), 69  
 is\_configdataclass (*TiePieDeviceConfig* attribute), 71  
 is\_configdataclass (*VisaCommunicationConfig* attribute), 32  
 is\_data\_ready\_polling\_interval\_sec (*TiePieDeviceConfig* attribute), 71  
 is\_done() (*MBW973* method), 129  
 is\_error() (*ExperimentManager* method), 181  
 is\_finished() (*ExperimentManager* method), 181  
 is\_generic\_type\_hint() (in *module* *hvl\_ccb.utils.typing*), 177  
 is\_in\_range() (*ILS2TRegDatatype* method), 166  
 is\_inactive (*CryLasLaser.LaserStatus* property), 89  
 is\_measurement\_data\_ready() (*TiePieOscilloscope* method), 80  
 is\_open (*Client* property), 20  
 is\_open (*LJMCommunication* property), 15  
 is\_open (*OpcUaCommunication* property), 20  
 is\_open (*SerialCommunication* property), 23  
 is\_open (*TelnetCommunication* property), 28  
 is\_polling() (*Poller* method), 170  
 is\_ready (*CryLasLaser.LaserStatus* property), 89  
 is\_running() (*ExperimentManager* method), 181  
 is\_started (*Technix* property), 68  
 is\_triggered() (*TiePieOscilloscope* method), 80  
 is\_valid\_scale\_range\_reversed\_str() (*PfeifferTPGConfig.Model* method), 149

## J

J (*LabJack.ThermocoupleType* attribute), 119  
 jerk\_time (*NewportSMC100PPConfig* attribute), 141  
 JOG (*ILS2T.Mode* attribute), 160  
 jog\_run() (*ILS2T* method), 162  
 jog\_stop() (*ILS2T* method), 162  
 JOGGING (*NewportStates* attribute), 146  
 JOGGING\_FROM\_DISABLE (*NewportSMC100PP.StateMessages* attribute), 134  
 JOGGING\_FROM\_READY (*NewportSMC100PP.StateMessages* attribute), 134  
 JOGN\_FAST (*ILS2TRegAddr* attribute), 165  
 JOGN\_SLOW (*ILS2TRegAddr* attribute), 165  
 JR (*NewportConfigCommands* attribute), 133

## K

K (*LabJack.TemperatureUnit* attribute), 118  
 K (*LabJack.ThermocoupleType* attribute), 119  
 K (*Temperature* attribute), 175  
 KELVIN (*Temperature* attribute), 175  
 keys() (*AsyncCommunicationProtocolConfig* class method), 13  
 keys() (*BaseCubeConfiguration* class method), 36  
 keys() (*BaseCubeOpcUaCommunicationConfig* class method), 37  
 keys() (*CryLasAttenuatorConfig* class method), 86  
 keys() (*CryLasAttenuatorSerialCommunicationConfig* class method), 88

keys() (*CryLasLaserConfig class method*), 92  
 keys() (*CryLasLaserSerialCommunicationConfig class method*), 94  
 keys() (*EmptyConfig class method*), 85, 179  
 keys() (*FuGConfig class method*), 101  
 keys() (*FuGSerialCommunicationConfig class method*), 109  
 keys() (*HeinzingerConfig class method*), 111  
 keys() (*HeinzingerSerialCommunicationConfig class method*), 115  
 keys() (*ILS2TConfig class method*), 164  
 keys() (*ILS2TModbusTcpCommunicationConfig class method*), 165  
 keys() (*LaudaProRp245eConfig class method*), 127  
 keys() (*LaudaProRp245eTcpCommunicationConfig class method*), 128  
 keys() (*LJMCommunicationConfig class method*), 17  
 keys() (*LuminoxConfig class method*), 168  
 keys() (*LuminoxSerialCommunicationConfig class method*), 170  
 keys() (*MBW973Config class method*), 131  
 keys() (*MBW973SerialCommunicationConfig class method*), 132  
 keys() (*ModbusTcpCommunicationConfig class method*), 19  
 keys() (*NewportSMC100PPConfig class method*), 141  
 keys() (*NewportSMC100PPSerialCommunicationConfig class method*), 145  
 keys() (*OpcUaCommunicationConfig class method*), 21  
 keys() (*PfeifferTPGConfig class method*), 149  
 keys() (*PfeifferTPGSerialCommunicationConfig class method*), 151  
 keys() (*PICubeConfiguration class method*), 42  
 keys() (*PICubeOpcUaCommunicationConfig class method*), 43  
 keys() (*PSI9000Config class method*), 97  
 keys() (*PSI9000VisaCommunicationConfig class method*), 99  
 keys() (*RTO1024Config class method*), 158  
 keys() (*RTO1024VisaCommunicationConfig class method*), 159  
 keys() (*SerialCommunicationConfig class method*), 25  
 keys() (*SupercubeAOpcUaConfiguration class method*), 64  
 keys() (*SupercubeConfiguration class method*), 53  
 keys() (*SupercubeOpcUaCommunicationConfig class method*), 54  
 keys() (*T560CommunicationConfig class method*), 45  
 keys() (*T560Config class method*), 48  
 keys() (*TcpCommunicationConfig class method*), 27  
 keys() (*TechnixConfig class method*), 69  
 keys() (*TechnixSerialCommunicationConfig class method*), 66  
 keys() (*TechnixTelnetCommunicationConfig class method*), 67  
 keys() (*TelnetCommunicationConfig class method*), 29  
 keys() (*TiePieDeviceConfig class method*), 71  
 keys() (*VisaCommunicationConfig class method*), 32  
 keys() (*VisaDeviceConfig class method*), 173  
 kV (*HeinzingerPNC.UnitVoltage attribute*), 114

**L**

LabJack (*class in hvl\_ccb.dev.labjack*), 117  
 LabJack.AInRange (*class in hvl\_ccb.dev.labjack*), 117  
 LabJack.BitLimit (*class in hvl\_ccb.dev.labjack*), 117  
 LabJack.CalMicroAmpere (*class in hvl\_ccb.dev.labjack*), 117  
 LabJack.CjcType (*class in hvl\_ccb.dev.labjack*), 117  
 LabJack.ClockFrequency (*class in hvl\_ccb.dev.labjack*), 117  
 LabJack.DeviceType (*class in hvl\_ccb.dev.labjack*), 118  
 LabJack.DIOStatus (*class in hvl\_ccb.dev.labjack*), 118  
 LabJack.TemperatureUnit (*class in hvl\_ccb.dev.labjack*), 118  
 LabJack.ThermocoupleType (*class in hvl\_ccb.dev.labjack*), 118  
 LabJackError, 122  
 LabJackIdentifierDIOError, 122  
 laser\_off() (*CryLasLaser method*), 89  
 laser\_on() (*CryLasLaser method*), 90  
 LaudaProRp245e (*class in hvl\_ccb.dev.lauda*), 123  
 LaudaProRp245eCommand (*class in hvl\_ccb.dev.lauda*), 125  
 LaudaProRp245eCommandError, 126  
 LaudaProRp245eConfig (*class in hvl\_ccb.dev.lauda*), 126  
 LaudaProRp245eConfig.ExtControlModeEnum (*class in hvl\_ccb.dev.lauda*), 126  
 LaudaProRp245eConfig.OperationModeEnum (*class in hvl\_ccb.dev.lauda*), 127  
 LaudaProRp245eTcpCommunication (*class in hvl\_ccb.dev.lauda*), 128  
 LaudaProRp245eTcpCommunicationConfig (*class in hvl\_ccb.dev.lauda*), 128  
 LEM4000S (*class in hvl\_ccb.utils.conversion\_sensor*), 173  
 LF (*FuGTerminators attribute*), 110  
 LFCR (*FuGTerminators attribute*), 110  
 LINE\_1 (*MessageBoard attribute*), 39  
 LINE\_10 (*MessageBoard attribute*), 39  
 LINE\_11 (*MessageBoard attribute*), 39  
 LINE\_12 (*MessageBoard attribute*), 39  
 LINE\_13 (*MessageBoard attribute*), 39  
 LINE\_14 (*MessageBoard attribute*), 39  
 LINE\_15 (*MessageBoard attribute*), 39  
 LINE\_2 (*MessageBoard attribute*), 40  
 LINE\_3 (*MessageBoard attribute*), 40  
 LINE\_4 (*MessageBoard attribute*), 40

LINE\_5 (*MessageBoard* attribute), 40  
 LINE\_6 (*MessageBoard* attribute), 40  
 LINE\_7 (*MessageBoard* attribute), 40  
 LINE\_8 (*MessageBoard* attribute), 40  
 LINE\_9 (*MessageBoard* attribute), 40  
**list\_devices()** (*TiePieOscilloscope* static method), 80  
**list\_directory()** (*RTO1024* method), 154  
**LJMCommunication** (class in *hvl\_ccb.comm.labjack\_ljm*), 15  
**LJMCommunicationConfig** (class in *hvl\_ccb.comm.labjack\_ljm*), 16  
**LJMCommunicationConfig.ConnectionType** (class in *hvl\_ccb.comm.labjack\_ljm*), 16  
**LJMCommunicationConfig.DeviceType** (class in *hvl\_ccb.comm.labjack\_ljm*), 16  
**LJMCommunicationError**, 17  
**lm34** (*LabJack.CjcType* attribute), 117  
**LMT70A** (class in *hvl\_ccb.utils.conversion\_sensor*), 174  
**load\_configuration()** (*RTO1024* method), 155  
**load\_device\_configuration()** (*T560* method), 48  
**local\_display()** (*RTO1024* method), 155  
**LOCKED** (*DoorStatus* attribute), 38  
**locked** (*DoorStatus* attribute), 56  
**LOW** (*LabJack.DIStatus* attribute), 118  
**LOWER\_TEMP** (*LaudaProRp245eCommand* attribute), 125  
**lower\_temp** (*LaudaProRp245eConfig* attribute), 127  
**Luminox** (class in *hvl\_ccb.dev.sst\_luminox*), 166  
**LuminoxConfig** (class in *hvl\_ccb.dev.sst\_luminox*), 167  
**LuminoxException**, 168  
**LuminoxMeasurementType** (class in *hvl\_ccb.dev.sst\_luminox*), 168  
**LuminoxMeasurementTypeDict** (in module *hvl\_ccb.dev.sst\_luminox*), 169  
**LuminoxMeasurementTypeError**, 169  
**LuminoxMeasurementTypeValue** (in module *hvl\_ccb.dev.sst\_luminox*), 169  
**LuminoxOutputMode** (class in *hvl\_ccb.dev.sst\_luminox*), 169  
**LuminoxOutputModeError**, 169  
**LuminoxSerialCommunication** (class in *hvl\_ccb.dev.sst\_luminox*), 169  
**LuminoxSerialCommunicationConfig** (class in *hvl\_ccb.dev.sst\_luminox*), 169  
**LUT** (*LMT70A* attribute), 174

## M

**mA** (*HeinzingerPNC.UnitCurrent* attribute), 113  
**MANUAL** (*EarthingStickOperatingStatus* attribute), 39  
**manual()** (*EarthingStick* class method), 56  
**manual\_1** (*EarthingStick* attribute), 56  
**manual\_2** (*EarthingStick* attribute), 56  
**manual\_3** (*EarthingStick* attribute), 56  
**manual\_4** (*EarthingStick* attribute), 56  
**manual\_5** (*EarthingStick* attribute), 56  
**manual\_6** (*EarthingStick* attribute), 57  
**MARK** (*SerialCommunicationParity* attribute), 25  
**max\_current** (*FuG* property), 100  
**max\_current** (*HeinzingerPNC* property), 114  
**max\_current** (*Technix* property), 68  
**max\_current** (*TechnixConfig* attribute), 70  
**max\_current\_hardware** (*FuG* property), 100  
**max\_current\_hardware** (*HeinzingerPNC* property), 114  
**max\_pr\_number** (*LaudaProRp245eConfig* attribute), 127  
**max\_pump\_level** (*LaudaProRp245eConfig* attribute), 127  
**max\_timeout\_retry\_nr** (*OpcUaCommunicationConfig* attribute), 22  
**max\_voltage** (*FuG* property), 100  
**max\_voltage** (*HeinzingerPNC* property), 114  
**max\_voltage** (*Technix* property), 68  
**max\_voltage** (*TechnixConfig* attribute), 70  
**max\_voltage\_hardware** (*FuG* property), 100  
**max\_voltage\_hardware** (*HeinzingerPNC* property), 114  
**MAXIMUM** (*LabJack.ClockFrequency* attribute), 118  
**mbar** (*PfeifferTPG.PressureUnits* attribute), 147  
**MBW973** (class in *hvl\_ccb.dev.mbw973*), 129  
**MBW973Config** (class in *hvl\_ccb.dev.mbw973*), 130  
**MBW973ControlRunningException**, 131  
**MBW973Error**, 131  
**MBW973PumpRunningException**, 131  
**MBW973SerialCommunication** (class in *hvl\_ccb.dev.mbw973*), 131  
**MBW973SerialCommunicationConfig** (class in *hvl\_ccb.dev.mbw973*), 131  
**measure()** (*PfeifferTPG* method), 148  
**measure\_all()** (*PfeifferTPG* method), 148  
**measure\_current()** (*HeinzingerDI* method), 112  
**measure\_voltage()** (*HeinzingerDI* method), 112  
**measure\_voltage\_current()** (*PSI9000* method), 96  
**MeasurementsDividerRatio** (class in *hvl\_ccb.dev.supercube2015.constants*), 59  
**MeasurementsScaledInput** (class in *hvl\_ccb.dev.supercube2015.constants*), 59  
**MessageBoard** (class in *hvl\_ccb.dev.cube.constants*), 39  
**micro\_step\_per\_full\_step\_factor** (*NewportsSMC100PPConfig* attribute), 141  
**Micron** (*PfeifferTPG.PressureUnits* attribute), 146  
**MILLIMETER\_MERCURY** (*Pressure* attribute), 175  
**MINIMUM** (*LabJack.ClockFrequency* attribute), 118  
**MMHG** (*Pressure* attribute), 175  
**ModbusTcpCommunication** (class in *hvl\_ccb.comm.modbus\_tcp*), 18  
**ModbusTcpCommunicationConfig** (class in *hvl\_ccb.comm.modbus\_tcp*), 19  
**ModbusTcpConnectionFailedException**, 19  
**model** (*PfeifferTPGConfig* attribute), 150

```

module
    hvl_ccb, 182
    hvl_ccb.comm, 33
    hvl_ccb.comm.base, 11
    hvl_ccb.comm.labjack_ljm, 15
    hvl_ccb.comm.modbus_tcp, 18
    hvl_ccb.comm.opc, 20
    hvl_ccb.comm.serial, 23
    hvl_ccb.comm.tcp, 26
    hvl_ccb.comm.telnet, 28
    hvl_ccb.comm.visa, 30
    hvl_ccb.configuration, 178
    hvl_ccb.dev, 173
    hvl_ccb.dev.base, 83
    hvl_ccb.dev.crylas, 85
    hvl_ccb.dev.cube, 44
    hvl_ccb.dev.cube.base, 33
    hvl_ccb.dev.cube.constants, 38
    hvl_ccb.dev.cube.picube, 41
    hvl_ccb.dev.ea_psi9000, 95
    hvl_ccb.dev.fug, 99
    hvl_ccb.dev.heinzinger, 110
    hvl_ccb.dev.highland_t560, 49
    hvl_ccb.dev.highland_t560.base, 44
    hvl_ccb.dev.highland_t560.channel, 46
    hvl_ccb.dev.highland_t560.device, 47
    hvl_ccb.dev.labjack, 116
    hvl_ccb.dev.lauda, 123
    hvl_ccb.dev.mbw973, 129
    hvl_ccb.dev.newport, 133
    hvl_ccb.dev.pfeiffer_tpg, 146
    hvl_ccb.dev.rs_rto1024, 152
    hvl_ccb.dev.se_ils2t, 160
    hvl_ccb.dev.sst_luminox, 166
    hvl_ccb.dev.supercube2015, 64
    hvl_ccb.dev.supercube2015.base, 49
    hvl_ccb.dev.supercube2015.constants, 55
    hvl_ccb.dev.supercube2015.typ_a, 62
    hvl_ccb.dev.technix, 70
    hvl_ccb.dev.technix.base, 65
    hvl_ccb.dev.technix.device, 68
    hvl_ccb.dev.tiepie, 82
    hvl_ccb.dev.tiepie.base, 71
    hvl_ccb.dev.tiepie.channel, 73
    hvl_ccb.dev.tiepie.device, 75
    hvl_ccb.dev.tiepie.generator, 76
    hvl_ccb.dev.tiepie.i2c, 78
    hvl_ccb.dev.tiepie.oscilloscope, 79
    hvl_ccb.dev.tiepie.utils, 81
    hvl_ccb.dev.utils, 170
    hvl_ccb.dev.visa, 171
    hvl_ccb.exception, 180
    hvl_ccb.experiment_manager, 180
    hvl_ccb.utils, 178

    hvl_ccb.utils.conversion_sensor, 173
    hvl_ccb.utils.conversion_unit, 174
    hvl_ccb.utils.enum, 176
    hvl_ccb.utils.typing, 177
    hvl_ccb.utils.validation, 177
    MONITOR_I (FuGProbusVRegisterGroups attribute), 107
    MONITOR_V (FuGProbusVRegisterGroups attribute), 107
    most_recent_error (FuGProbusVConfigRegisters property), 105
    motion_distance_per_full_step (NewportSMC100PPConfig attribute), 141
    motor_config (NewportSMC100PPConfig property), 141
    move_to_absolute_position() (NewportSMC100PP method), 138
    move_to_relative_position() (NewportSMC100PP method), 138
    move_wait_sec (NewportSMC100PPConfig attribute), 141
    MOVING (NewportSMC100PP.StateMessages attribute), 134
    MOVING (NewportStates attribute), 146
    MS_NOMINAL_CURRENT (PSI9000 attribute), 95
    MS_NOMINAL_VOLTAGE (PSI9000 attribute), 95
    MULTI_COMMANDS_MAX (VisaCommunication attribute), 30
    MULTI_COMMANDS_SEPARATOR (VisaCommunication attribute), 30

N
n_channels (TiePieOscilloscope property), 80
n_max_try_get_device (TiePieDeviceConfig attribute), 71
NameEnum (class in hvl_ccb.utils.enum), 176
NAMES (SerialCommunicationParity attribute), 25
names() (RTO1024.TriggerModes class method), 153
namespace_index (BaseCubeConfiguration attribute), 36
namespace_index (SupercubeConfiguration attribute), 53
NATIVEONLY (TiePieOscilloscopeAutoResolutionModes attribute), 80
NED_END_OF_TURN (NewportSMC100PP.MotorErrors attribute), 134
NEG (ILS2T.Ref16Jog attribute), 160
NEG_FAST (ILS2T.Ref16Jog attribute), 160
NEGATIVE (FuGPolarities attribute), 104
NEGATIVE (Polarity attribute), 40
negative_software_limit (NewportSMC100PPConfig attribute), 141
NewportConfigCommands (class in hvl_ccb.dev.newport), 133
NewportControllerError, 133
NewportError, 133

```

NewportMotorError, 133  
NewportMotorPowerSupplyWasCutError, 134  
NewportSerialCommunicationError, 145  
NewportSMC100PP (*class in hvl\_ccb.dev.newport*), 134  
NewportSMC100PP.MotorErrors (*class in hvl\_ccb.dev.newport*), 134  
NewportSMC100PP.StateMessages (*class in hvl\_ccb.dev.newport*), 134  
NewportSMC100PPConfig (*class in hvl\_ccb.dev.newport*), 140  
NewportSMC100PPConfig.EspStageConfig (*class in hvl\_ccb.dev.newport*), 140  
NewportSMC100PPConfig.HomeSearch (*class in hvl\_ccb.dev.newport*), 140  
NewportSMC100PPSerialCommunication (*class in hvl\_ccb.dev.newport*), 142  
NewportSMC100PPSerialCommunication.ControllerError (*class in hvl\_ccb.dev.newport*), 142  
NewportSMC100PPSerialCommunicationConfig (*class in hvl\_ccb.dev.newport*), 144  
NewportStates (*class in hvl\_ccb.dev.newport*), 145  
NewportUncertainPositionError, 146  
NO (*FuGDigitalVal attribute*), 102  
NO\_ERROR (*NewportSMC100PPSerialCommunication.ControllerError attribute*), 142  
NO\_REF (*NewportStates attribute*), 146  
NO\_REF\_ESP\_STAGE\_ERROR  
  *portSMC100PP.StateMessages*  
  134  
NO\_REF\_FROM\_CONFIG  
  *portSMC100PP.StateMessages*  
  134  
NO\_REF\_FROM\_DISABLED  
  *portSMC100PP.StateMessages*  
  134  
NO\_REF\_FROM\_HOMING  
  *portSMC100PP.StateMessages*  
  135  
NO\_REF\_FROM\_JOGGING  
  *portSMC100PP.StateMessages*  
  135  
NO\_REF\_FROM\_MOVING  
  *portSMC100PP.StateMessages*  
  135  
NO\_REF\_FROM\_READY  
  *portSMC100PP.StateMessages*  
  135  
NO\_REF\_FROM\_RESET  
  *portSMC100PP.StateMessages*  
  135  
No\_sensor (*PfeifferTPG.SensorStatus attribute*), 147  
NO\_SOURCE (*PowerSetup attribute*), 40  
NOISE (*TiePieGeneratorSignalType attribute*), 77  
noise\_level\_measurement\_channel\_1 (*BaseCube-Configuration attribute*), 36  
noise\_level\_measurement\_channel\_2 (*BaseCube-Configuration attribute*), 36  
noise\_level\_measurement\_channel\_3 (*BaseCube-Configuration attribute*), 36  
noise\_level\_measurement\_channel\_4 (*BaseCube-Configuration attribute*), 36  
NONE (*ILS2T.Ref16Jog attribute*), 160  
NONE (*LabJack.ThermocoupleType attribute*), 119  
None (*PfeifferTPG.SensorTypes attribute*), 147  
NONE (*SerialCommunicationParity attribute*), 25  
NORMAL (*RTO1024.TriggerModes attribute*), 153  
noSen (*PfeifferTPG.SensorTypes attribute*), 147  
noSENSOR (*PfeifferTPG.SensorTypes attribute*), 147  
not\_defined (*AlarmText attribute*), 56  
nr\_trials\_activate (*LuminoxConfig attribute*), 168  
HvLCommunicationProtocol (*class in hvl\_ccb.comm.base*), 14  
Number (*in module hvl\_ccb.utils.typing*), 177  
number\_of\_decimals (*HeinzingerConfig attribute*), 111  
number\_of\_sensors (*PfeifferTPG property*), 148

## O

OFF (*SerialCommunicationParity attribute*), 25  
OFF (*AutoInstallMode attribute*), 44  
OFF (*FuGDigitalVal attribute*), 102  
OFF (*GateMode attribute*), 44  
OFF (*HeinzingerDI.OutputStatus attribute*), 112  
OFF (*TriggerMode attribute*), 46  
offset (*TiePieGeneratorConfig property*), 77  
OH (*NewportConfigCommands attribute*), 133  
Ok (*PfeifferTPG.SensorStatus attribute*), 147  
on (*FuG property*), 100  
ON (*FuGDigitalVal attribute*), 102  
on (*FuGProbusVDIRegisters property*), 106  
ON (*HeinzingerDI.OutputStatus attribute*), 112  
ON (*ILS2T.State attribute*), 161  
ONE (*HeinzingerConfig.RecordingsEnum attribute*), 111  
ONE (*LabJack.AInRange attribute*), 117  
ONE (*SerialCommunicationStopbits attribute*), 26  
ONE\_HUNDREDTH (*LabJack.AInRange attribute*), 117  
ONE\_POINT\_FIVE (*SerialCommunicationStopbits attribute*), 26  
ONE\_TENTH (*LabJack.AInRange attribute*), 117  
ONLYUPWARDSOFTZERO (*FuGRampModes attribute*), 108  
OPC\_MAX\_YEAR (*BaseCube attribute*), 33  
OPC\_MIN\_YEAR (*BaseCube attribute*), 33  
OpcUaCommunication (*class in hvl\_ccb.comm.opc*), 20  
OpcUaCommunicationConfig (*class in hvl\_ccb.comm.opc*), 21  
OpcUaCommunicationIOError, 22  
OpcUaCommunicationTimeoutError, 22  
OpcUaSubHandler (*class in hvl\_ccb.comm.opc*), 22

OPEN (*DoorStatus attribute*), 38  
open (*DoorStatus attribute*), 56  
OPEN (*EarthingStickOperation attribute*), 39  
OPEN (*EarthingStickStatus attribute*), 39  
open (*EarthingStickStatus attribute*), 58  
open() (*CommunicationProtocol method*), 14  
open() (*LaudaProRp245eTcpCommunication method*), 128  
open() (*LJMCommunication method*), 15  
open() (*ModbusTcpCommunication method*), 18  
open() (*NullCommunicationProtocol method*), 14  
open() (*OpcUaCommunication method*), 20  
open() (*SerialCommunication method*), 23  
open() (*Tcp method*), 26  
open() (*TelnetCommunication method*), 28  
open() (*VisaCommunication method*), 30  
open\_interlock (*Technix property*), 69  
open\_shutter() (*CryLasLaser method*), 90  
open\_timeout (*VisaCommunicationConfig attribute*), 32  
OPENED (*CryLasLaser.AnswersShutter attribute*), 88  
OPENED (*CryLasLaserShutterStatus attribute*), 94  
operate (*BaseCube property*), 35  
operate (*PICube property*), 42  
operate() (*Supercube2015Base method*), 51  
OPERATION\_MODE (*LaudaProRp245eCommand attribute*), 125  
operation\_mode (*LaudaProRp245eConfig attribute*), 127  
optional\_defaults() (*AsyncCommunicationProtocolConfig class method*), 13  
optional\_defaults() (*BaseCubeConfiguration class method*), 36  
optional\_defaults() (*BaseCubeOpcUaCommunicationConfig class method*), 37  
optional\_defaults() (*CryLasAttenuatorConfig class method*), 87  
optional\_defaults() (*CryLasAttenuatorSerialCommunicationConfig class method*), 88  
optional\_defaults() (*CryLasLaserConfig class method*), 92  
optional\_defaults() (*CryLasLaserSerialCommunicationConfig class method*), 94  
optional\_defaults() (*EmptyConfig class method*), 85, 179  
optional\_defaults() (*FuGConfig class method*), 101  
optional\_defaults() (*FuGSerialCommunicationConfig class method*), 109  
optional\_defaults() (*HeinzingerConfig class method*), 111  
optional\_defaults() (*HeinzingerSerialCommunicationConfig class method*), 115  
optional\_defaults() (*ILS2TConfig class method*), 164  
optional\_defaults()  
                  (*ILS2TModbusTcpCommunicationConfig class method*), 165  
optional\_defaults() (*LaudaProRp245eConfig class method*), 127  
optional\_defaults() (*LaudaProRp245eTcpCommunicationConfig class method*), 129  
optional\_defaults() (*LJMCommunicationConfig class method*), 17  
optional\_defaults() (*LuminoxConfig class method*), 168  
optional\_defaults() (*LuminoxSerialCommunicationConfig class method*), 170  
optional\_defaults() (*MBW973Config class method*), 131  
optional\_defaults()  
                  (*MBW973SerialCommunicationConfig class method*), 132  
optional\_defaults() (*ModbusTcpCommunicationConfig class method*), 19  
optional\_defaults() (*NewportSMC100PPConfig class method*), 141  
optional\_defaults()  
                  (*NewportSMC100PPSerialCommunicationConfig class method*), 145  
optional\_defaults() (*OpcUaCommunicationConfig class method*), 22  
optional\_defaults() (*PfeifferTPGConfig class method*), 150  
optional\_defaults() (*PfeifferTPGSerialCommunicationConfig class method*), 151  
optional\_defaults() (*PICubeConfiguration class method*), 42  
optional\_defaults() (*PICubeOpcUaCommunicationConfig class method*), 43  
optional\_defaults() (*PSI9000Config class method*), 97  
optional\_defaults()  
                  (*PSI9000VisaCommunicationConfig class method*), 99  
optional\_defaults() (*RTO1024Config class method*), 158  
optional\_defaults()  
                  (*RTO1024VisaCommunicationConfig class method*), 159  
optional\_defaults() (*SerialCommunicationConfig class method*), 25  
optional\_defaults() (*SupercubeAOpcUaConfiguration class method*), 64  
optional\_defaults() (*SupercubeConfiguration class method*), 53  
optional\_defaults() (*SupercubeOpcUaCommunicationConfig class method*), 54  
optional\_defaults() (*T560CommunicationConfig*

*class method*), 45  
**optional\_defaults()** (*T560Config class method*), 48  
**optional\_defaults()** (*TcpCommunicationConfig class method*), 27  
**optional\_defaults()** (*TechnixConfig class method*), 70  
**optional\_defaults()** (*TechnixSerialCommunicationConfig class method*), 66  
**optional\_defaults()** (*TechnixTelnetCommunicationConfig class method*), 67  
**optional\_defaults()** (*TelnetCommunicationConfig class method*), 29  
**optional\_defaults()** (*TiePieDeviceConfig class method*), 72  
**optional\_defaults()** (*VisaCommunicationConfig class method*), 32  
**optional\_defaults()** (*VisaDeviceConfig class method*), 173  
**OSCILLOSCOPE** (*TiePieDeviceType attribute*), 72  
**OT** (*NewportConfigCommands attribute*), 133  
**out** (*FuGProbusVDORegisters property*), 106  
**out\_1\_1** (*GeneralSupport attribute*), 59  
**out\_1\_2** (*GeneralSupport attribute*), 59  
**out\_2\_1** (*GeneralSupport attribute*), 59  
**out\_2\_2** (*GeneralSupport attribute*), 59  
**out\_3\_1** (*GeneralSupport attribute*), 59  
**out\_3\_2** (*GeneralSupport attribute*), 59  
**out\_4\_1** (*GeneralSupport attribute*), 59  
**out\_4\_2** (*GeneralSupport attribute*), 59  
**out\_5\_1** (*GeneralSupport attribute*), 59  
**out\_5\_2** (*GeneralSupport attribute*), 59  
**out\_6\_1** (*GeneralSupport attribute*), 59  
**out\_6\_2** (*GeneralSupport attribute*), 59  
**OUTPUT** (*FuGProbusIVCommands attribute*), 105  
**OUTPUT** (*GateMode attribute*), 44  
**output** (*Technix property*), 69  
**output()** (*GeneralSupport class method*), 59  
**output\_off()** (*FuGProbusIV method*), 104  
**output\_off()** (*HeinzingerDI method*), 112  
**output\_on()** (*HeinzingerDI method*), 113  
**OUTPUT\_POWER\_EXCEEDED** (*New-portSMC100PP.MotorErrors attribute*), 134  
**output\_status** (*HeinzingerDI property*), 113  
**OUTPUTONCMD** (*FuGProbusVRegisterGroups attribute*), 107  
**OUTPUTX0** (*FuGProbusVRegisterGroups attribute*), 107  
**OUTPUTX1** (*FuGProbusVRegisterGroups attribute*), 107  
**OUTPUTX2** (*FuGProbusVRegisterGroups attribute*), 107  
**OUTPUTXCMD** (*FuGProbusVRegisterGroups attribute*), 107  
**outX0** (*FuG property*), 100  
**outX1** (*FuG property*), 100  
**outX2** (*FuG property*), 100  
**outXCMD** (*FuG property*), 100  
**Overrange** (*PfeifferTPG.SensorStatus attribute*), 147  

## P

**PA** (*Pressure attribute*), 175  
**PARAM\_MISSING\_OR\_INVALID** (*New-portSMC100PPSerialCommunication.ControllerErrors attribute*), 142  
**parity** (*CryLasAttenuatorSerialCommunicationConfig attribute*), 88  
**parity** (*CryLasLaserSerialCommunicationConfig attribute*), 94  
**parity** (*FuGSerialCommunicationConfig attribute*), 109  
**parity** (*HeinzingerSerialCommunicationConfig attribute*), 115  
**parity** (*LuminoxSerialCommunicationConfig attribute*), 170  
**parity** (*MBW973SerialCommunicationConfig attribute*), 132  
**parity** (*NewportSMC100PPSerialCommunicationConfig attribute*), 145  
**parity** (*PfeifferTPGSerialCommunicationConfig attribute*), 151  
**Parity** (*SerialCommunicationConfig attribute*), 24  
**parity** (*SerialCommunicationConfig attribute*), 25  
**parse\_read\_measurement\_value()** (*LuminoxMeasurementType method*), 168  
**partial\_pressure\_o2** (*LuminoxMeasurementType attribute*), 168  
**Pascal** (*PfeifferTPG.PressureUnits attribute*), 146  
**PASCAL** (*Pressure attribute*), 175  
**pause()** (*LaudaProRp245e method*), 123  
**pause\_ramp()** (*LaudaProRp245e method*), 123  
**PBR** (*PfeifferTPG.SensorTypes attribute*), 147  
**PEAK\_CURRENT\_LIMIT** (*New-portSMC100PP.MotorErrors attribute*), 134  
**peak\_output\_current\_limit** (*New-portSMC100PPConfig attribute*), 141  
**percent\_o2** (*LuminoxMeasurementType attribute*), 168  
**period** (*T560 property*), 48  
**PfeifferTPG** (*class in hvl\_ccb.dev.pfeiffer\_tpg*), 146  
**PfeifferTPG.PressureUnits** (*class in hvl\_ccb.dev.pfeiffer\_tpg*), 146  
**PfeifferTPG.SensorStatus** (*class in hvl\_ccb.dev.pfeiffer\_tpg*), 147  
**PfeifferTPG.SensorTypes** (*class in hvl\_ccb.dev.pfeiffer\_tpg*), 147  
**PfeifferTPGConfig** (*class in hvl\_ccb.dev.pfeiffer\_tpg*), 149  
**PfeifferTPGConfig.Model** (*class in hvl\_ccb.dev.pfeiffer\_tpg*), 149  
**PfeifferTPGError**, 150

PfeifferTPGSerialCommunication (class in *hvl\_ccb.dev.pfeiffer\_tpg*), 150  
 PfeifferTPGSerialCommunicationConfig (class in *hvl\_ccb.dev.pfeiffer\_tpg*), 150  
 PICube (class in *hvl\_ccb.dev(cube.picube)*), 41  
 PICubeConfiguration (class in *hvl\_ccb.dev(cube.picube)*), 42  
 PICubeOpcUaCommunication (class in *hvl\_ccb.dev(cube.picube)*), 43  
 PICubeOpcUaCommunicationConfig (class in *hvl\_ccb.dev(cube.picube)*), 43  
 PICubeTestParameterError, 40  
 PKR (*PfeifferTPG.SensorTypes* attribute), 147  
 Polarity (class in *hvl\_ccb.dev(cube.constants)*), 40  
 Polarity (class in *hvl\_ccb.dev(highland\_t560.base)*), 44  
 POLARITY (*FuGProbusIVCommands* attribute), 105  
 polarity (*PICube* property), 42  
 Poller (class in *hvl\_ccb.dev.utils*), 170  
 polling (*LuminoxOutputMode* attribute), 169  
 polling\_delay\_sec (BaseCubeConfiguration attribute), 37  
 polling\_interval (*MBW973Config* attribute), 131  
 polling\_interval\_sec (BaseCubeConfiguration attribute), 37  
 polling\_interval\_sec (*TechnixConfig* attribute), 70  
 polling\_period (*CryLasLaserConfig* attribute), 92  
 polling\_timeout (*CryLasLaserConfig* attribute), 92  
 port (*ModbusTcpCommunicationConfig* attribute), 19  
 port (*OpcUaCommunicationConfig* attribute), 22  
 port (*SerialCommunicationConfig* attribute), 25  
 port (*SupercubeOpcUaCommunicationConfig* attribute), 54  
 port (*T560CommunicationConfig* attribute), 45  
 port (*TcpCommunicationConfig* attribute), 27  
 port (*TechnixTelnetCommunicationConfig* attribute), 67  
 port (*TelnetCommunicationConfig* attribute), 29  
 port (*VisaCommunicationConfig* attribute), 32  
 POS (*ILS2T.Ref16Jog* attribute), 161  
 POS\_END\_OF\_TURN (*NewportSMC100PP.MotorErrors* attribute), 134  
 POS\_FAST (*ILS2T.Ref16Jog* attribute), 161  
 POSITION (*ILS2TRegAddr* attribute), 165  
 POSITION\_OUT\_OF\_LIMIT (*NewportSMC100PPSerialCommunication.ControllerErrors* attribute), 142  
 POSITIVE (*FuGPolarities* attribute), 104  
 POSITIVE (*Polarity* attribute), 40  
 positive\_software\_limit (*NewportSMC100PPConfig* attribute), 141  
 post\_force\_value() (*NewportSMC100PPConfig* method), 141  
 post\_stop\_pause\_sec (*TechnixConfig* attribute), 70  
 POUNDS\_PER\_SQUARE\_INCH (*Pressure* attribute), 175  
 Power (class in *hvl\_ccb.dev.supercube2015.constants*), 60  
 POWER\_INVERTER\_220V (*PowerSetup* attribute), 40  
 power\_limit (*PSI9000Config* attribute), 98  
 power\_setup (*PICube* property), 42  
 PowerSetup (class in *hvl\_ccb.dev(cube.constants)*), 40  
 PowerSetup (class in *hvl\_ccb.dev.supercube2015.constants*), 60  
 pre\_sample\_ratio (*TiePieOscilloscopeConfig* property), 81  
 prepare\_ultra\_segmentation() (*RTO1024* method), 155  
 preserve\_type() (in module *hvl\_ccb.utils.conversion\_unit*), 175  
 Pressure (class in *hvl\_ccb.utils.conversion\_unit*), 174  
 probe\_offset (*TiePieOscilloscopeChannelConfig* property), 74  
 PSI (*Pressure* attribute), 175  
 PSI9000 (class in *hvl\_ccb.dev.ea\_psi9000*), 95  
 PSI9000Config (class in *hvl\_ccb.dev.ea\_psi9000*), 97  
 PSI9000Error, 98  
 PSI9000VisaCommunication (class in *hvl\_ccb.dev.ea\_psi9000*), 98  
 PSI9000VisaCommunicationConfig (class in *hvl\_ccb.dev.ea\_psi9000*), 98  
 PT100 (*LabJack.ThermocoupleType* attribute), 119  
 PT1000 (*LabJack.ThermocoupleType* attribute), 119  
 PT500 (*LabJack.ThermocoupleType* attribute), 119  
 PTP (*ILS2T.Mode* attribute), 160  
 PublicPropertiesReprMixin (class in *hvl\_ccb.dev.tiepie.utils*), 81  
 PULSE (*TiePieGeneratorSignalType* attribute), 77  
 pump\_init (*LaudaProRp245eConfig* attribute), 127  
 PUMP\_LEVEL (*LaudaProRp245eCommand* attribute), 125

## Q

QIL (*NewportConfigCommands* attribute), 133  
 QUERY (*FuGProbusIVCommands* attribute), 105  
 query() (*CryLasLaserSerialCommunication* method), 93  
 query() (*FuGSerialCommunication* method), 109  
 query() (*NewportSMC100PPSerialCommunication* method), 143  
 query() (*PfeifferTPGSerialCommunication* method), 150  
 query() (*SyncCommunicationProtocol* method), 14  
 query() (*T560Communication* method), 45  
 query() (*VisaCommunication* method), 30  
 query\_all() (*CryLasLaserSerialCommunication* method), 93  
 query\_command() (*LaudaProRp245eTcpCommunication* method), 128

**query\_multiple()** (*NewportSMC100PPSerialCommunication method*), 143  
**query\_polling()** (*Luminox method*), 167  
**query\_status()** (*Technix method*), 69  
**QUEUE** (*AutoInstallMode attribute*), 44  
**QUICK\_STOP** (*SafetyStatus attribute*), 41  
**QUICKSTOP** (*ILS2T.State attribute*), 161  
**QuickStop** (*SafetyStatus attribute*), 61  
**quickstop()** (*ILS2T method*), 162  
**quit** (*Errors attribute*), 58  
**quit\_error()** (*BaseCube method*), 35  
**quit\_error()** (*Supercube2015Base method*), 51

**R**

**R** (*LabJack.ThermocoupleType attribute*), 119  
**raise\_()** (*FuGErrorcodes method*), 103  
**RAMP\_ACC** (*ILS2TRegAddr attribute*), 165  
**RAMP\_CONTINUE** (*LaudaProRp245eCommand attribute*), 125  
**RAMP\_DECEL** (*ILS2TRegAddr attribute*), 165  
**RAMP\_DELETE** (*LaudaProRp245eCommand attribute*), 125  
**RAMP\_ITERATIONS** (*LaudaProRp245eCommand attribute*), 125  
**RAMP\_N\_MAX** (*ILS2TRegAddr attribute*), 165  
**RAMP\_PAUSE** (*LaudaProRp245eCommand attribute*), 125  
**RAMP\_SELECT** (*LaudaProRp245eCommand attribute*), 125  
**RAMP\_SET** (*LaudaProRp245eCommand attribute*), 125  
**RAMP\_START** (*LaudaProRp245eCommand attribute*), 126  
**RAMP\_STOP** (*LaudaProRp245eCommand attribute*), 126  
**RAMP\_TYPE** (*ILS2TRegAddr attribute*), 165  
**rampmode** (*FuGProbusVSetRegisters property*), 108  
**ramprate** (*FuGProbusVSetRegisters property*), 108  
**rampstate** (*FuGProbusVSetRegisters property*), 108  
**RAMPUPWARDS** (*FuGRampModes attribute*), 108  
**RangeEnum** (*class in hvl\_ccb.utils.enum*), 176  
**RATEDCURRENT** (*FuGReadbackChannels attribute*), 108  
**RATEDVOLTAGE** (*FuGReadbackChannels attribute*), 108  
**read()** (*AsyncCommunicationProtocol method*), 11  
**read()** (*BaseCube method*), 35  
**read()** (*CryLasLaserSerialCommunication method*), 93  
**read()** (*LaudaProRp245eTcpCommunication method*), 128  
**read()** (*MBW973 method*), 130  
**read()** (*OpcUaCommunication method*), 21  
**read()** (*Supercube2015Base method*), 51  
**read()** (*Tcp method*), 26  
**read\_all()** (*AsyncCommunicationProtocol method*), 11  
**read\_bytes()** (*AsyncCommunicationProtocol method*), 12  
**read\_bytes()** (*SerialCommunication method*), 23  
**read\_bytes()** (*TelnetCommunication method*), 28  
**read\_float()** (*MBW973 method*), 130  
**read\_holding\_registers()** (*ModbusTcpCommunication method*), 18  
**read\_input\_registers()** (*ModbusTcpCommunication method*), 18  
**read\_int()** (*MBW973 method*), 130  
**read\_measurement()** (*RTO1024 method*), 155  
**read\_measurements()** (*MBW973 method*), 130  
**read\_name()** (*LJMCommunication method*), 15  
**read\_nonempty()** (*AsyncCommunicationProtocol method*), 12  
**read\_output\_while\_polling** (*TechnixConfig attribute*), 70  
**read\_resistance()** (*LabJack method*), 120  
**read\_single\_bytes()** (*SerialCommunication method*), 23  
**read\_streaming()** (*Luminox method*), 167  
**read\_termination** (*VisaCommunicationConfig attribute*), 32  
**read\_text()** (*AsyncCommunicationProtocol method*), 12  
**read\_text()** (*NewportSMC100PPSerialCommunication method*), 143  
**read\_text\_nonempty()** (*AsyncCommunicationProtocol method*), 12  
**READ\_TEXT\_SKIP\_PREFIXES** (*CryLasLaserSerialCommunication attribute*), 92  
**read\_thermocouple()** (*LabJack method*), 121  
**readback\_data** (*FuGProbusVConfigRegisters property*), 105  
**READBACKCHANNEL** (*FuGProbusIVCommands attribute*), 105  
**ready** (*BaseCube property*), 35  
**READY** (*CryLasLaser.AnswersStatus attribute*), 88  
**READY** (*ILS2T.State attribute*), 161  
**READY** (*NewportStates attribute*), 146  
**ready()** (*Supercube2015Base method*), 51  
**READY\_ACTIVE** (*CryLasLaser.LaserStatus attribute*), 89  
**READY\_FROM\_DISABLE** (*NewportSMC100PP.StateMessages attribute*), 135  
**READY\_FROM\_HOMING** (*NewportSMC100PP.StateMessages attribute*), 135  
**READY\_FROM\_JOGGING** (*NewportSMC100PP.StateMessages attribute*), 135  
**READY\_FROM\_MOVING** (*NewportSMC100PP.StateMessages attribute*), 135  
**READY\_INACTIVE** (*CryLasLaser.LaserStatus attribute*), 89  
**record\_length** (*TiePieOscilloscopeConfig property*), 81

RED\_OPERATE (*SafetyStatus attribute*), 41  
RED\_READY (*SafetyStatus attribute*), 41  
RedOperate (*SafetyStatus attribute*), 62  
RedReady (*SafetyStatus attribute*), 62  
reg\_3 (*FuGProbusVDIRegisters property*), 106  
RegAddr (*ILS2T attribute*), 161  
RegDatatype (*ILS2T attribute*), 161  
register\_pulse\_time (*TechnixConfig attribute*), 70  
RELATIVE (*TiePieOscilloscopeTriggerLevelMode attribute*), 75  
RELATIVE\_POSITION\_MOTOR (*ILS2T.ActionsPtp attribute*), 160  
RELATIVE\_POSITION\_TARGET (*ILS2T.ActionsPtp attribute*), 160  
remote (*Technix property*), 69  
remove\_device() (*DeviceSequenceMixin method*), 84  
require\_block\_measurement\_support (*TiePieDeviceConfig attribute*), 72  
required\_keys() (*AsyncCommunicationProtocolConfig class method*), 13  
required\_keys() (*BaseCubeConfiguration class method*), 37  
required\_keys() (*BaseCubeOpcUaCommunicationConfig class method*), 38  
required\_keys() (*CryLasAttenuatorConfig class method*), 87  
required\_keys() (*CryLasAttenuatorSerialCommunicationConfig class method*), 88  
required\_keys() (*CryLasLaserConfig class method*), 92  
required\_keys() (*CryLasLaserSerialCommunicationConfig class method*), 94  
required\_keys() (*EmptyConfig class method*), 85, 179  
required\_keys() (*FuGConfig class method*), 101  
required\_keys() (*FuGSerialCommunicationConfig class method*), 110  
required\_keys() (*HeinzingerConfig class method*), 111  
required\_keys() (*HeinzingerSerialCommunicationConfig class method*), 115  
required\_keys() (*ILS2TConfig class method*), 164  
required\_keys() (*ILS2TModbusTcpCommunicationConfig class method*), 165  
required\_keys() (*LaudaProRp245eConfig class method*), 127  
required\_keys() (*LaudaProRp245eTcpCommunicationConfig class method*), 129  
required\_keys() (*LJMCommunicationConfig class method*), 17  
required\_keys() (*LuminoxConfig class method*), 168  
required\_keys() (*LuminoxSerialCommunicationConfig class method*), 170  
required\_keys() (*MBW973Config class method*), 131  
required\_keys() (*MBW973SerialCommunicationConfig class method*), 132  
required\_keys() (*ModbusTcpCommunicationConfig class method*), 19  
required\_keys() (*NewportSMC100PPConfig class method*), 142  
required\_keys() (*NewportSMC100PPSerialCommunicationConfig class method*), 145  
required\_keys() (*OpcUaCommunicationConfig class method*), 22  
required\_keys() (*PfeifferTPGConfig class method*), 150  
required\_keys() (*PfeifferTPGSerialCommunicationConfig class method*), 151  
required\_keys() (*PICubeConfiguration class method*), 43  
required\_keys() (*PICubeOpcUaCommunicationConfig class method*), 43  
required\_keys() (*PSI9000Config class method*), 98  
required\_keys() (*PSI9000VisaCommunicationConfig class method*), 99  
required\_keys() (*RTO1024Config class method*), 159  
required\_keys() (*RTO1024VisaCommunicationConfig class method*), 159  
required\_keys() (*SerialCommunicationConfig class method*), 25  
required\_keys() (*SupercubeAOpcUaConfiguration class method*), 64  
required\_keys() (*SupercubeConfiguration class method*), 53  
required\_keys() (*SupercubeOpcUaCommunicationConfig class method*), 54  
required\_keys() (*T560CommunicationConfig class method*), 45  
required\_keys() (*T560Config class method*), 49  
required\_keys() (*TcpCommunicationConfig class method*), 27  
required\_keys() (*TechnixConfig class method*), 70  
required\_keys() (*TechnixSerialCommunicationConfig class method*), 66  
required\_keys() (*TechnixTelnetCommunicationConfig class method*), 67  
required\_keys() (*TelnetCommunicationConfig class method*), 29  
required\_keys() (*TiePieDeviceConfig class method*), 72  
required\_keys() (*VisaCommunicationConfig class method*), 32  
required\_keys() (*VisaDeviceConfig class method*), 173  
reset (*BreakdownDetection attribute*), 56  
RESET (*FuGProbusIVCommands attribute*), 105  
reset() (*FuGProbusIV method*), 104

**reset()** (*NewportSMC100PP method*), 138  
**reset()** (*VisaDevice method*), 172  
**reset\_error()** (*ILS2T method*), 162  
**reset\_interface()** (*HeinzingerDI method*), 113  
**reset\_ramp()** (*LaudaProRp245e method*), 123  
**resolution** (*TiePieOscilloscopeConfig property*), 81  
**response\_sleep\_time** (*CryLasAttenuatorConfig attribute*), 87  
**RISING** (*TiePieOscilloscopeTriggerKind attribute*), 75  
**RISING\_OR\_FALLING** (*TiePieOscilloscopeTriggerKind attribute*), 75  
**RMS\_CURRENT\_LIMIT** (*NewportSMC100PP.MotorErrors attribute*), 134  
**rpm\_max\_init** (*ILS2TConfig attribute*), 164  
**rs485\_address** (*NewportSMC100PPConfig attribute*), 142  
**RT01024** (*class in hvl\_ccb.dev.rs\_rto1024*), 152  
**RT01024.TriggerModes** (*class in hvl\_ccb.dev.rs\_rto1024*), 152  
**RT01024Config** (*class in hvl\_ccb.dev.rs\_rto1024*), 158  
**RT01024Error**, 159  
**RT01024VisaCommunication** (*class in hvl\_ccb.dev.rs\_rto1024*), 159  
**RT01024VisaCommunicationConfig** (*class in hvl\_ccb.dev.rs\_rto1024*), 159  
**run()** (*ExperimentManager method*), 181  
**run()** (*LaudaProRp245e method*), 123  
**run\_continuous\_acquisition()** (*RT01024 method*), 155  
**run\_single\_acquisition()** (*RT01024 method*), 155  
**RUNNING** (*ExperimentStatus attribute*), 182

**S**

**S** (*LabJack.ThermocoupleType attribute*), 119  
**SA** (*NewportConfigCommands attribute*), 133  
**SafeGround** (*class in hvl\_ccb.dev.tiepie.channel*), 73  
**Safety** (*class in hvl\_ccb.dev.supercube2015.constants*), 61  
**SafetyStatus** (*class in hvl\_ccb.dev.cube.constants*), 40  
**SafetyStatus** (*class in hvl\_ccb.dev.supercube2015.constants*), 61  
**sample\_frequency** (*TiePieOscilloscopeConfig property*), 81  
**save\_configuration()** (*RT01024 method*), 155  
**save\_device\_configuration()** (*T560 method*), 48  
**save\_waveform\_history()** (*RT01024 method*), 155  
**SCALE** (*ILS2TRegAddr attribute*), 165  
**ScalingFactorValueError**, 166  
**screw\_scaling** (*NewportSMC100PPConfig attribute*), 142  
**send\_command()** (*NewportSMC100PPSerialCommunication method*), 144  
**send\_command()** (*PfeifferTPGSerialCommunication method*), 150  
**send\_hello()** (*Client method*), 20  
**send\_pulses()** (*LabJack method*), 121  
**send\_stop()** (*NewportSMC100PPSerialCommunication method*), 144  
**Sensor** (*class in hvl\_ccb.utils.conversion\_sensor*), 174  
**Sensor\_error** (*PfeifferTPG.SensorStatus attribute*), 147  
**Sensor\_off** (*PfeifferTPG.SensorStatus attribute*), 147  
**sensor\_status** (*LuminoxMeasurementType attribute*), 168  
**SERIAL** (*LaudaProRp245eConfig.ExtControlModeEnum attribute*), 126  
**serial\_number** (*LuminoxMeasurementType attribute*), 168  
**serial\_number** (*TiePieDeviceConfig attribute*), 72  
**SerialCommunication** (*class in hvl\_ccb.comm.serial*), 23  
**SerialCommunicationBytesize** (*class in hvl\_ccb.comm.serial*), 24  
**SerialCommunicationConfig** (*class in hvl\_ccb.comm.serial*), 24  
**SerialCommunicationIOError**, 25  
**SerialCommunicationParity** (*class in hvl\_ccb.comm.serial*), 25  
**SerialCommunicationStopbits** (*class in hvl\_ccb.comm.serial*), 25  
**Server** (*class in hvl\_ccb.comm.opc*), 22  
**set\_acceleration()** (*NewportSMC100PP method*), 138  
**set\_acquire\_length()** (*RT01024 method*), 155  
**set\_ain\_differential()** (*LabJack method*), 121  
**set\_ain\_range()** (*LabJack method*), 121  
**set\_ain\_resistance()** (*LabJack method*), 121  
**set\_ain\_resolution()** (*LabJack method*), 121  
**set\_ain\_thermocouple()** (*LabJack method*), 121  
**set\_analog\_output()** (*LabJack method*), 122  
**set\_attenuation()** (*CryLasAttenuator method*), 86  
**set\_cee16\_socket()** (*Supercube2015Base method*), 51  
**set\_channel\_offset()** (*RT01024 method*), 156  
**set\_channel\_position()** (*RT01024 method*), 156  
**set\_channel\_range()** (*RT01024 method*), 156  
**set\_channel\_scale()** (*RT01024 method*), 156  
**set\_channel\_state()** (*RT01024 method*), 157  
**set\_clock()** (*LabJack method*), 122  
**set\_control\_mode()** (*LaudaProRp245e method*), 123  
**set\_current()** (*HeinzingerDI method*), 113  
**set\_current()** (*HeinzingerPNC method*), 114  
**set\_digital\_output()** (*LabJack method*), 122  
**set\_display\_unit()** (*PfeifferTPG method*), 148  
**set\_earthing\_manual()** (*Supercube2015Base method*), 51  
**set\_external\_temp()** (*LaudaProRp245e method*), 124

set\_full\_scale\_mbar() (*PfeifferTPG method*), 148  
set\_full\_scale\_unitless() (*PfeifferTPG method*), 149  
set\_init\_attenuation() (*CryLasAttenuator method*), 86  
set\_init\_shutter\_status() (*CryLasLaser method*), 90  
set\_jog\_speed() (*ILS2T method*), 162  
set\_lower\_limits() (*PSI9000 method*), 96  
set\_max\_acceleration() (*ILS2T method*), 162  
set\_max\_deceleration() (*ILS2T method*), 162  
set\_max\_rpm() (*ILS2T method*), 163  
set\_measuring\_options() (*MBW973 method*), 130  
set\_message\_board() (*BaseCube method*), 35  
set\_motor\_configuration() (*NewportSMC100PP method*), 138  
set\_negative\_software\_limit() (*NewportSMC100PP method*), 139  
set\_number\_of\_recordings() (*HeinzingerDI method*), 113  
set\_output() (*PSI9000 method*), 96  
set\_positive\_software\_limit() (*NewportSMC100PP method*), 139  
set\_pulse\_energy() (*CryLasLaser method*), 90  
set\_pump\_level() (*LaudaProRp245e method*), 124  
set\_ramp\_iterations() (*LaudaProRp245e method*), 124  
set\_ramp\_program() (*LaudaProRp245e method*), 124  
set\_ramp\_segment() (*LaudaProRp245e method*), 124  
set\_ramp\_type() (*ILS2T method*), 163  
set\_reference\_point() (*RTO1024 method*), 157  
set\_register() (*FuGProbusV method*), 105  
set\_remote\_control() (*Supercube2015Base method*), 51  
set\_repetition\_rate() (*CryLasLaser method*), 90  
set\_repetitions() (*RTO1024 method*), 157  
set\_slope() (*Supercube2015WithFU method*), 63  
set\_status\_board() (*BaseCube method*), 35  
set\_support\_output() (*Supercube2015Base method*), 52  
set\_support\_output\_impulse() (*Supercube2015Base method*), 52  
set\_system\_lock() (*PSI9000 method*), 96  
set\_t13\_socket() (*Supercube2015Base method*), 52  
set\_target\_voltage() (*Supercube2015WithFU method*), 63  
set\_temp\_set\_point() (*LaudaProRp245e method*), 124  
set\_transmission() (*CryLasAttenuator method*), 86  
set\_trigger\_level() (*RTO1024 method*), 157  
set\_trigger\_mode() (*RTO1024 method*), 158  
set\_trigger\_source() (*RTO1024 method*), 158  
set\_upper\_limits() (*PSI9000 method*), 96  
set\_voltage() (*HeinzingerDI method*), 113  
set\_voltage() (*HeinzingerPNC method*), 114  
set\_voltage\_current() (*PSI9000 method*), 96  
SETCURRENT (*FuGProbusVRegisterGroups attribute*), 107  
setup (*Power attribute*), 60  
setvalue (*FuGProbusVSetRegisters property*), 108  
SETVOLTAGE (*FuGProbusVRegisterGroups attribute*), 107  
SEVENBITS (*SerialCommunicationBytesize attribute*), 24  
SHORT\_CIRCUIT (*NewportSMC100PP.MotorErrors attribute*), 134  
shunt (*LEM4000S attribute*), 174  
SHUTDOWN\_CURRENT\_LIMIT (*PSI9000 attribute*), 95  
SHUTDOWN\_VOLTAGE\_LIMIT (*PSI9000 attribute*), 95  
ShutterStatus (*CryLasLaser attribute*), 89  
ShutterStatus (*CryLasLaserConfig attribute*), 91  
signal\_type (*TiePieGeneratorConfig property*), 77  
SINE (*TiePieGeneratorSignalType attribute*), 77  
SingleCommDevice (class in *hvl\_ccb.dev.base*), 85  
SIXBITS (*SerialCommunicationBytesize attribute*), 24  
SIXTEEN (*HeinzingerConfig.RecordingsEnum attribute*), 111  
SIXTEEN\_BIT (*TiePieOscilloscopeResolution attribute*), 81  
SL (*NewportConfigCommands attribute*), 133  
SN (*FuGReadbackChannels attribute*), 108  
SOFTWARE\_INTERNAL\_SIXTY (*CryLasLaser.RepetitionRates attribute*), 89  
SOFTWARE\_INTERNAL\_TEN (*CryLasLaser.RepetitionRates attribute*), 89  
SOFTWARE\_INTERNAL\_TWENTY (*CryLasLaser.RepetitionRates attribute*), 89  
software\_revision (*LuminoxMeasurementType attribute*), 168  
SPACE (*SerialCommunicationParity attribute*), 25  
SPECIALRAMPUPWARDS (*FuGRampModes attribute*), 108  
spoll() (*VisaCommunication method*), 31  
spoll\_handler() (*VisaDevice method*), 172  
SQUARE (*TiePieGeneratorSignalType attribute*), 77  
SR (*NewportConfigCommands attribute*), 133  
srq\_mask (*FuGProbusVConfigRegisters property*), 105  
srq\_status (*FuGProbusVConfigRegisters property*), 106  
stage\_configuration (*NewportSMC100PPConfig attribute*), 142  
START (*LaudaProRp245eCommand attribute*), 126  
start() (*BaseCube method*), 35  
start() (*CryLasAttenuator method*), 86  
start() (*CryLasLaser method*), 90  
start() (*Device method*), 83  
start() (*DeviceSequenceMixin method*), 84  
start() (*ExperimentManager method*), 181  
start() (*FuG method*), 100  
start() (*FuGProbusIV method*), 104

**start()** (*HeinzingerDI method*), 113  
**start()** (*HeinzingerPNC method*), 114  
**start()** (*ILS2T method*), 163  
**start()** (*LabJack method*), 122  
**start()** (*LaudaProRp245e method*), 124  
**start()** (*Luminox method*), 167  
**start()** (*MBW973 method*), 130  
**start()** (*NewportSMC100PP method*), 139  
**start()** (*PfeifferTPG method*), 149  
**start()** (*PSI9000 method*), 97  
**start()** (*RTO1024 method*), 158  
**start()** (*SingleCommDevice method*), 85  
**start()** (*Supercube2015Base method*), 52  
**start()** (*Technix method*), 69  
**start()** (*TiePieGeneratorMixin method*), 77  
**start()** (*TiePieI2CHostMixin method*), 78  
**start()** (*TiePieOscilloscope method*), 80  
**start()** (*VisaDevice method*), 172  
**start\_control()** (*MBW973 method*), 130  
**start\_measurement()** (*TiePieOscilloscope method*), 80  
**start\_polling()** (*Poller method*), 171  
**start\_ramp()** (*LaudaProRp245e method*), 124  
**STARTING** (*ExperimentStatus attribute*), 182  
**States** (*NewportSMC100PP attribute*), 135  
**status** (*BaseCube property*), 35  
**status** (*ExperimentManager property*), 181  
**status** (*FuGProbusVConfigRegisters property*), 106  
**status** (*FuGProbusVDORegisters property*), 107  
**status** (*Technix property*), 69  
**status\_1\_closed** (*EarthingStick attribute*), 57  
**status\_1\_connected** (*EarthingStick attribute*), 57  
**status\_1\_open** (*EarthingStick attribute*), 57  
**status\_2\_closed** (*EarthingStick attribute*), 57  
**status\_2\_connected** (*EarthingStick attribute*), 57  
**status\_2\_open** (*EarthingStick attribute*), 57  
**status\_3\_closed** (*EarthingStick attribute*), 57  
**status\_3\_connected** (*EarthingStick attribute*), 57  
**status\_3\_open** (*EarthingStick attribute*), 57  
**status\_4\_closed** (*EarthingStick attribute*), 57  
**status\_4\_connected** (*EarthingStick attribute*), 57  
**status\_4\_open** (*EarthingStick attribute*), 57  
**status\_5\_closed** (*EarthingStick attribute*), 57  
**status\_5\_connected** (*EarthingStick attribute*), 57  
**status\_5\_open** (*EarthingStick attribute*), 57  
**status\_6\_closed** (*EarthingStick attribute*), 57  
**status\_6\_connected** (*EarthingStick attribute*), 57  
**status\_6\_open** (*EarthingStick attribute*), 57  
**status\_closed()** (*EarthingStick class method*), 57  
**status\_connected()** (*EarthingStick class method*), 57  
**status\_error** (*Safety attribute*), 61  
**status\_green** (*Safety attribute*), 61  
**status\_open()** (*EarthingStick class method*), 57  
**status\_ready\_for\_red** (*Safety attribute*), 61  
**status\_red** (*Safety attribute*), 61  
**STATUSBYTE** (*FuGReadbackChannels attribute*), 108  
**stop** (*Errors attribute*), 58  
**STOP** (*LaudaProRp245eCommand attribute*), 126  
**stop()** (*BaseCube method*), 35  
**stop()** (*CryLasLaser method*), 90  
**stop()** (*Device method*), 83  
**stop()** (*DeviceSequenceMixin method*), 84  
**stop()** (*ExperimentManager method*), 181  
**stop()** (*FuGProbusIV method*), 104  
**stop()** (*HeinzingerDI method*), 113  
**stop()** (*ILS2T method*), 163  
**stop()** (*LabJack method*), 122  
**stop()** (*LaudaProRp245e method*), 125  
**stop()** (*Luminox method*), 167  
**stop()** (*MBW973 method*), 130  
**stop()** (*NewportSMC100PP method*), 139  
**stop()** (*PfeifferTPG method*), 149  
**stop()** (*PSI9000 method*), 97  
**stop()** (*RTO1024 method*), 158  
**stop()** (*SingleCommDevice method*), 85  
**stop()** (*Supercube2015Base method*), 52  
**stop()** (*Technix method*), 69  
**stop()** (*TiePieGeneratorMixin method*), 77  
**stop()** (*TiePieI2CHostMixin method*), 78  
**stop()** (*TiePieOscilloscope method*), 80  
**stop()** (*VisaDevice method*), 172  
**stop\_acquisition()** (*RTO1024 method*), 158  
**stop\_motion()** (*NewportSMC100PP method*), 139  
**stop\_number** (*Errors attribute*), 58  
**stop\_polling()** (*Poller method*), 171  
**stop\_ramp()** (*LaudaProRp245e method*), 125  
**STOP\_SAFETY\_STATUSES** (in module *hvl\_ccb.dev(cube.constants)*), 40  
**stopbits** (*CryLasAttenuatorSerialCommunicationConfig attribute*), 88  
**stopbits** (*CryLasLaserSerialCommunicationConfig attribute*), 94  
**stopbits** (*FuGSerialCommunicationConfig attribute*), 110  
**stopbits** (*HeinzingerSerialCommunicationConfig attribute*), 116  
**stopbits** (*LuminoxSerialCommunicationConfig attribute*), 170  
**stopbits** (*MBW973SerialCommunicationConfig attribute*), 132  
**stopbits** (*NewportSMC100PPSerialCommunicationConfig attribute*), 145  
**stopbits** (*PfeifferTPGSerialCommunicationConfig attribute*), 151  
**Stopbits** (*SerialCommunicationConfig attribute*), 24  
**stopbits** (*SerialCommunicationConfig attribute*), 25  
**streaming** (*LuminoxOutputMode attribute*), 169  
**StrEnumBase** (class in *hvl\_ccb.utils.enum*), 176

```

sub_handler (BaseCubeOpcUaCommunicationConfig attribute), 38
sub_handler (OpcUaCommunicationConfig attribute), 22
sub_handler (SupercubeOpcUaCommunicationConfig attribute), 54
suitable_range() (TiePieOscilloscopeRange static method), 74
Supercube2015Base (class hvl_ccb.dev.supercube2015.base), 49
Supercube2015WithFU (class hvl_ccb.dev.supercube2015.typ_a), 62
SupercubeAOpCuaCommunication (class hvl_ccb.dev.supercube2015.typ_a), 63
SupercubeAOpCuaConfiguration (class hvl_ccb.dev.supercube2015.typ_a), 63
SupercubeConfiguration (class hvl_ccb.dev.supercube2015.base), 52
SupercubeOpcEndpoint (class hvl_ccb.dev.supercube2015.constants), 62
SupercubeOpcUaCommunication (class hvl_ccb.dev.supercube2015.base), 53
SupercubeOpcUaCommunicationConfig (class hvl_ccb.dev.supercube2015.base), 53
SupercubeSubscriptionHandler (class hvl_ccb.dev.supercube2015.base), 54
switchto_green (Safety attribute), 61
switchto_operate (Safety attribute), 61
switchto_ready (Safety attribute), 61
SyncCommunicationProtocol (class hvl_ccb.comm.base), 14
SyncCommunicationProtocolConfig (class hvl_ccb.comm.base), 14

T
T (LabJack.ThermocoupleType attribute), 119
t13_1 (GeneralSockets attribute), 58
t13_2 (GeneralSockets attribute), 58
t13_3 (GeneralSockets attribute), 58
t13_socket_1 (BaseCube attribute), 36
t13_socket_2 (BaseCube attribute), 36
t13_socket_3 (BaseCube attribute), 36
T13_SOCKET_PORTS (in module hvl_ccb.dev.supercube2015.constants), 62
T1MS (FuGMonitorModes attribute), 103
T200MS (FuGMonitorModes attribute), 103
T20MS (FuGMonitorModes attribute), 103
T256US (FuGMonitorModes attribute), 103
T4 (LabJack.DeviceType attribute), 118
T4 (LJMCommunicationConfig.DeviceType attribute), 16
T40MS (FuGMonitorModes attribute), 103
T4MS (FuGMonitorModes attribute), 103
T560 (class in hvl_ccb.dev.higland_t560.device), 47
T560Communication (class hvl_ccb.dev.higland_t560.base), 45
T560CommunicationConfig (class hvl_ccb.dev.higland_t560.base), 45
T560Config (class in hvl_ccb.dev.higland_t560.device), 48
T560Error, 46
T7 (LabJack.DeviceType attribute), 118
in T7 (LJMCommunicationConfig.DeviceType attribute), 16
T7_PRO (LabJack.DeviceType attribute), 118
in T7_PRO (LJMCommunicationConfig.DeviceType attribute), 16
T800MS (FuGMonitorModes attribute), 103
T80MS (FuGMonitorModes attribute), 104
in target_pulse_energy (CryLasLaser property), 91
Tcp (class in hvl_ccb.comm.tcp), 26
in TCP (LJMCommunicationConfig.ConnectionType attribute), 16
in TcpCommunicationConfig (class in hvl_ccb.comm.tcp), 27
in TCPIP_INSTR (VisaCommunicationConfig.InterfaceType attribute), 31
in TCPIP_SOCKET (VisaCommunicationConfig.InterfaceType attribute), 31
in TEC1 (CryLasLaser.AnswersStatus attribute), 88
TEC2 (CryLasLaser.AnswersStatus attribute), 89
Technix (class in hvl_ccb.dev.technix.device), 68
TechnixConfig (class in hvl_ccb.dev.technix.device), 69
TechnixError, 65
TechnixFaultError, 65
TechnixSerialCommunication (class hvl_ccb.dev.technix.base), 65
TechnixSerialCommunicationConfig (class hvl_ccb.dev.technix.base), 65
TechnixTelnetCommunication (class hvl_ccb.dev.technix.base), 66
TechnixTelnetCommunicationConfig (class hvl_ccb.dev.technix.base), 66
TelnetCommunication (class in hvl_ccb.comm.telnet), 28
TelnetCommunicationConfig (class hvl_ccb.comm.telnet), 28
TelnetError, 29
TEMP (ILS2TRegAddr attribute), 165
TEMP_SET_POINT (LaudaProRp245eCommand attribute), 126
temp_set_point_init (LaudaProRp245eConfig attribute), 128
Temperature (class in hvl_ccb.utils.conversion_unit), 175
temperature_sensor (LuminoxMeasurementType attribute), 169
temperature_unit (LMT70A attribute), 174
TEN (LabJack.AInRange attribute), 117

```

TEN (*LabJack.CalMicroAmpere attribute*), 117  
 TEN\_MHZ (*LabJack.ClockFrequency attribute*), 118  
 terminator (*AsyncCommunicationProtocolConfig attribute*), 14  
 terminator (*CryLasAttenuatorSerialCommunicationConfig attribute*), 88  
 terminator (*CryLasLaserSerialCommunicationConfig attribute*), 94  
 TERMINATOR (*FuGProbusIVCommands attribute*), 105  
 terminator (*FuGProbusVConfigRegisters property*), 106  
 terminator (*FuGSerialCommunicationConfig attribute*), 110  
 terminator (*HeinzingerSerialCommunicationConfig attribute*), 116  
 terminator (*LaudaProRp245eTcpCommunicationConfig attribute*), 129  
 terminator (*LuminoxSerialCommunicationConfig attribute*), 170  
 terminator (*MBW973SerialCommunicationConfig attribute*), 132  
 terminator (*NewportSMC100PPSerialCommunicationConfig attribute*), 145  
 terminator (*PfeifferTPGSerialCommunicationConfig attribute*), 151  
 terminator (*T560CommunicationConfig attribute*), 46  
 terminator\_str() (*SerialCommunicationConfig method*), 25  
 THIRTY\_TWO\_BIT (*LabJack.BitLimit attribute*), 117  
 TiePieDeviceConfig (*class in hvl\_ccb.dev.tiepie.base*), 71  
 TiePieDeviceType (*class in hvl\_ccb.dev.tiepie.base*), 72  
 TiePieError, 72  
 TiePieGeneratorConfig (*class in hvl\_ccb.dev.tiepie.generator*), 76  
 TiePieGeneratorConfigLimits (*class in hvl\_ccb.dev.tiepie.generator*), 77  
 TiePieGeneratorMixin (*class in hvl\_ccb.dev.tiepie.generator*), 77  
 TiePieGeneratorSignalType (*class in hvl\_ccb.dev.tiepie.generator*), 77  
 TiePieHS5 (*class in hvl\_ccb.dev.tiepie.device*), 75  
 TiePieHS6 (*class in hvl\_ccb.dev.tiepie.device*), 76  
 TiePieI2CHostConfig (*class in hvl\_ccb.dev.tiepie.i2c*), 78  
 TiePieI2CHostConfigLimits (*class in hvl\_ccb.dev.tiepie.i2c*), 78  
 TiePieI2CHostMixin (*class in hvl\_ccb.dev.tiepie.i2c*), 78  
 TiePieOscilloscope (*class in hvl\_ccb.dev.tiepie.oscilloscope*), 79  
 TiePieOscilloscopeAutoResolutionModes (*class in hvl\_ccb.dev.tiepie.oscilloscope*), 80  
 TiePieOscilloscopeChannelConfig (*class in hvl\_ccb.dev.tiepie.channel*), 73  
 TiePieOscilloscopeChannelConfigLimits (*class in hvl\_ccb.dev.tiepie.channel*), 74  
 TiePieOscilloscopeChannelCoupling (*class in hvl\_ccb.dev.tiepie.channel*), 74  
 TiePieOscilloscopeConfig (*class in hvl\_ccb.dev.tiepie.oscilloscope*), 80  
 TiePieOscilloscopeConfigLimits (*class in hvl\_ccb.dev.tiepie.oscilloscope*), 81  
 TiePieOscilloscopeRange (*class in hvl\_ccb.dev.tiepie.channel*), 74  
 TiePieOscilloscopeResolution (*class in hvl\_ccb.dev.tiepie.oscilloscope*), 81  
 TiePieOscilloscopeTriggerKind (*class in hvl\_ccb.dev.tiepie.channel*), 74  
 TiePieOscilloscopeTriggerLevelMode (*class in hvl\_ccb.dev.tiepie.channel*), 75  
 TiePieWS5 (*class in hvl\_ccb.dev.tiepie.device*), 76  
 timeout (*CryLasAttenuatorSerialCommunicationConfig attribute*), 88  
 timeout (*CryLasLaserSerialCommunicationConfig attribute*), 94  
 timeout (*FuGSerialCommunicationConfig attribute*), 110  
 timeout (*HeinzingerSerialCommunicationConfig attribute*), 116  
 timeout (*LuminoxSerialCommunicationConfig attribute*), 170  
 timeout (*MBW973SerialCommunicationConfig attribute*), 132  
 timeout (*NewportSMC100PPSerialCommunicationConfig attribute*), 145  
 timeout (*PfeifferTPGSerialCommunicationConfig attribute*), 151  
 timeout (*SerialCommunicationConfig attribute*), 25  
 timeout (*TelnetCommunicationConfig attribute*), 29  
 timeout (*VisaCommunicationConfig attribute*), 32  
 timeout\_interval (*BaseCubeConfiguration attribute*), 37  
 timeout\_status\_change (*BaseCubeConfiguration attribute*), 37  
 timeout\_test\_parameters (*PICubeConfiguration attribute*), 43  
 Torr (*PfeifferTPG.PressureUnits attribute*), 146  
 TORR (*Pressure attribute*), 175  
 TPG25xA (*PfeifferTPGConfig.Model attribute*), 149  
 TPGx6x (*PfeifferTPGConfig.Model attribute*), 149  
 TPR (*PfeifferTPG.SensorTypes attribute*), 147  
 transmission (*CryLasAttenuator property*), 86  
 TRIANGLE (*TiePieGeneratorSignalType attribute*), 77  
 trigger\_enabled (*TiePieOscilloscopeChannelConfig property*), 74  
 trigger\_hysteresis (*TiePieOscilloscopeChannel-*

*Config property), 74*  
**trigger\_kind** (*TiePieOscilloscopeChannelConfig property*), 74  
**trigger\_level** (*T560 property*), 48  
**trigger\_level** (*TiePieOscilloscopeChannelConfig property*), 74  
**trigger\_level\_mode** (*TiePieOscilloscopeChannelConfig property*), 74  
**trigger\_mode** (*T560 property*), 48  
**trigger\_timeout** (*TiePieOscilloscopeConfig property*), 81  
**triggered** (*BreakdownDetection attribute*), 56  
**TriggerMode** (*class in hvl\_ccb.dev.hIGHLAND\_t560.base*), 46  
**TWELVE\_BIT** (*TiePieOscilloscopeResolution attribute*), 81  
**TWELVE\_HUNDRED\_FIFTY\_KHZ** (*LabJack.ClockFrequency attribute*), 118  
**TWENTY\_FIVE\_HUNDRED\_KHZ** (*LabJack.ClockFrequency attribute*), 118  
**TWENTY\_MHZ** (*LabJack.ClockFrequency attribute*), 118  
**TWENTY\_VOLT** (*TiePieOscilloscopeRange attribute*), 74  
**TWO** (*HeinzingerConfig.RecordingsEnum attribute*), 111  
**TWO** (*SerialCommunicationStopbits attribute*), 26  
**TWO\_HUNDRED** (*LabJack.CalMicroAmpere attribute*), 117  
**TWO\_HUNDRED\_MILLI\_VOLT** (*TiePieOscilloscopeRange attribute*), 74  
**TWO\_VOLT** (*TiePieOscilloscopeRange attribute*), 74

**U**

**Underrange** (*PfeifferTPG.SensorStatus attribute*), 147  
**Unit** (*class in hvl\_ccb.utils.conversion\_unit*), 175  
**unit** (*ILS2TModbusTcpCommunicationConfig attribute*), 165  
**unit** (*ModbusTcpCommunicationConfig attribute*), 19  
**unit** (*PfeifferTPG property*), 149  
**unit()** (*RangeEnum class method*), 176  
**unit\_current** (*HeinzingerPNC property*), 114  
**unit\_voltage** (*HeinzingerPNC property*), 114  
**UNKNOWN** (*HeinzingerDI.OutputStatus attribute*), 112  
**UNKNOWN** (*HeinzingerPNC.UnitCurrent attribute*), 113  
**UNKNOWN** (*HeinzingerPNC.UnitVoltage attribute*), 114  
**UNKNOWN** (*TiePieGeneratorSignalType attribute*), 77  
**UNKNOWN** (*TiePieOscilloscopeAutoResolutionModes attribute*), 80  
**UNKNOWN** (*TiePieOscilloscopeTriggerLevelMode attribute*), 75  
**UNREADY\_INACTIVE** (*CryLasLaser.LaserStatus attribute*), 89  
**update\_laser\_status()** (*CryLasLaser method*), 91  
**update\_parameter** (*OpcUaCommunicationConfig attribute*), 22  
**update\_repetition\_rate()** (*CryLasLaser method*), 91  
**update\_shutter\_status()** (*CryLasLaser method*), 91  
**update\_target\_pulse\_energy()** (*CryLasLaser method*), 91  
**UpdateEspStageInfo** (*NewportSMC100PPConfig.EspStageConfig attribute*), 140  
**UPPER\_TEMP** (*LaudaProRp245eCommand attribute*), 126  
**upper\_temp** (*LaudaProRp245eConfig attribute*), 128  
**USB** (*LaudaProRp245eConfig.ExtControlModeEnum attribute*), 126  
**USB** (*LJMCommunicationConfig.ConnectionType attribute*), 16  
**use\_external\_clock()** (*T560 method*), 48  
**user\_position\_offset** (*NewportSMC100PPConfig attribute*), 142  
**user\_steps()** (*ILS2T method*), 163

**V**

**V** (*HeinzingerPNC.UnitVoltage attribute*), 114  
**VA** (*NewportConfigCommands attribute*), 133  
**validate\_and\_resolve\_host()** (*in module hvl\_ccb.utils.validation*), 177  
**validate\_bool()** (*in module hvl\_ccb.utils.validation*), 177  
**validate\_number()** (*in module hvl\_ccb.utils.validation*), 177  
**validate\_pump\_level()** (*LaudaProRp245e method*), 125  
**validate\_tcp\_port()** (*in module hvl\_ccb.utils.validation*), 178  
**value** (*FuGProbusVMonitorRegisters property*), 107  
**value** (*LabJack.AInRange property*), 117  
**value\_raw** (*FuGProbusVMonitorRegisters property*), 107  
**ValueEnum** (*class in hvl\_ccb.utils.enum*), 176  
**VB** (*NewportConfigCommands attribute*), 133  
**velocity** (*NewportSMC100PPConfig attribute*), 142  
**visa\_backend** (*VisaCommunicationConfig attribute*), 32  
**VisaCommunication** (*class in hvl\_ccb.comm.visa*), 30  
**VisaCommunicationConfig** (*class in hvl\_ccb.comm.visa*), 31  
**VisaCommunicationConfig.InterfaceType** (*class in hvl\_ccb.comm.visa*), 31  
**VisaCommunicationError**, 32  
**VisaDevice** (*class in hvl\_ccb.dev.visa*), 171  
**VisaDeviceConfig** (*class in hvl\_ccb.dev.visa*), 172  
**VOLT** (*ILS2TRegAddr attribute*), 165  
**Volt** (*PfeifferTPG.PressureUnits attribute*), 146  
**voltage** (*FuG property*), 101  
**VOLTAGE** (*FuGProbusIVCommands attribute*), 105  
**VOLTAGE** (*FuGReadbackChannels attribute*), 108  
**voltage** (*Technix property*), 69  
**voltage\_actual** (*PICube property*), 42  
**voltage\_lower\_limit** (*PSI9000Config attribute*), 98

**voltage\_max** (*PICube property*), 42  
**voltage\_max** (*Power attribute*), 60  
**voltage\_monitor** (*FuG property*), 101  
**voltage\_primary** (*PICube property*), 42  
**voltage\_primary** (*Power attribute*), 60  
**voltage\_regulation** (*Technix property*), 69  
**voltage\_slope** (*Power attribute*), 60  
**voltage\_target** (*Power attribute*), 60  
**voltage\_upper\_limit** (*PSI9000Config attribute*), 98

**W**

**WAIT\_AFTER\_WRITE** (*VisaCommunication attribute*), 30  
**wait\_for\_polling\_result()** (*Poller method*), 171  
**wait\_operation\_complete()** (*VisaDevice method*), 172  
**wait\_sec\_initialisation** (*PSI9000Config attribute*), 98  
**wait\_sec\_max\_disable** (*ILS2TConfig attribute*), 164  
**wait\_sec\_post\_absolute\_position** (*ILS2TConfig attribute*), 164  
**wait\_sec\_post\_activate** (*LuminoxConfig attribute*), 168  
**wait\_sec\_post\_cannot\_disable** (*ILS2TConfig attribute*), 164  
**wait\_sec\_post\_enable** (*ILS2TConfig attribute*), 164  
**wait\_sec\_post\_relative\_step** (*ILS2TConfig attribute*), 164  
**wait\_sec\_pre\_read\_or\_write** (*LaudaProRp245eTcpCommunicationConfig attribute*), 129  
**wait\_sec\_read\_text\_nonempty** (*AsyncCommunicationProtocolConfig attribute*), 14  
**wait\_sec\_read\_text\_nonempty** (*FuGSerialCommunicationConfig attribute*), 110  
**wait\_sec\_read\_text\_nonempty** (*HeinzingerSerialCommunicationConfig attribute*), 116  
**wait\_sec\_retry\_get\_device** (*TiePieDeviceConfig attribute*), 72  
**wait\_sec\_settings\_effect** (*PSI9000Config attribute*), 98  
**wait\_sec\_stop\_commands** (*FuGConfig attribute*), 101  
**wait\_sec\_stop\_commands** (*HeinzingerConfig attribute*), 111  
**wait\_sec\_system\_lock** (*PSI9000Config attribute*), 98  
**wait\_sec\_trials\_activate** (*LuminoxConfig attribute*), 168  
**wait\_timeout\_retry\_sec** (*OpcUaCommunicationConfig attribute*), 22  
**wait\_until\_motor\_initialized()** (*NewportSMC100PP method*), 140  
**wait\_until\_ready()** (*CryLasLaser method*), 91  
**waveform** (*TiePieGeneratorConfig property*), 77  
**WIFI** (*LJMCommunicationConfig.ConnectionType attribute*), 16

**wrap\_libtiepie\_exception()** (*in module hvl\_ccb.dev.tiepie.base*), 72  
**write()** (*AsyncCommunicationProtocol method*), 12  
**write()** (*BaseCube method*), 36  
**write()** (*MBW973 method*), 130  
**write()** (*OpcUaCommunication method*), 21  
**write()** (*Supercube2015Base method*), 52  
**write()** (*Tcp method*), 26  
**write()** (*VisaCommunication method*), 31  
**write\_absolute\_position()** (*ILS2T method*), 163  
**write\_bytes()** (*AsyncCommunicationProtocol method*), 12  
**write\_bytes()** (*SerialCommunication method*), 23  
**write\_bytes()** (*TelnetCommunication method*), 28  
**write\_command()** (*LaudaProRp245eTcpCommunication method*), 128  
**write\_name()** (*LJMCommunication method*), 16  
**write\_names()** (*LJMCommunication method*), 16  
**write\_registers()** (*ModbusTcpCommunication method*), 18  
**write\_relative\_step()** (*ILS2T method*), 163  
**write\_termination** (*VisaCommunicationConfig attribute*), 32  
**write\_text()** (*AsyncCommunicationProtocol method*), 13  
**WRONG\_ESP\_STAGE** (*NewportSMC100PP.MotorErrors attribute*), 134

**X**

**x\_stat** (*FuGProbusVDIRegisters property*), 106  
**XOUTPUTS** (*FuGProbusIVCommands attribute*), 105

**Y**

**YES** (*FuGDigitalVal attribute*), 102

**Z**

**ZX** (*NewportConfigCommands attribute*), 133