

---

# HVL Common Code Base Documentation

*Release 0.3.5*

**Mikolaj Rybiński, David Gruber**

**Feb 18, 2020**



---

## Contents:

---

<b>1</b>	<b>HVL Common Code Base</b>	<b>1</b>
1.1	Features . . . . .	1
1.2	Credits . . . . .	2
<b>2</b>	<b>Installation</b>	<b>3</b>
2.1	Stable release . . . . .	3
2.2	From sources . . . . .	3
<b>3</b>	<b>Usage</b>	<b>5</b>
<b>4</b>	<b>API Documentation</b>	<b>7</b>
4.1	hvl_ccb package . . . . .	7
<b>5</b>	<b>Contributing</b>	<b>119</b>
5.1	Types of Contributions . . . . .	119
5.2	Get Started! . . . . .	120
5.3	Merge Request Guidelines . . . . .	121
5.4	Tips . . . . .	121
5.5	Deploying . . . . .	121
<b>6</b>	<b>Credits</b>	<b>123</b>
6.1	Development Lead . . . . .	123
6.2	Contributors . . . . .	123
<b>7</b>	<b>History</b>	<b>125</b>
7.1	0.3.5 (2020-02-18) . . . . .	125
7.2	0.3.4 (2019-12-20) . . . . .	125
7.3	0.3.3 (2019-05-08) . . . . .	126
7.4	0.3.2 (2019-05-08) . . . . .	126
7.5	0.3.1 (2019-05-02) . . . . .	126
7.6	0.3 (2019-05-02) . . . . .	126
7.7	0.2.1 (2019-04-01) . . . . .	126
7.8	0.2.0 (2019-03-31) . . . . .	126
7.9	0.1.0 (2019-02-06) . . . . .	126
<b>8</b>	<b>Indices and tables</b>	<b>127</b>
	<b>Python Module Index</b>	<b>129</b>

---



# CHAPTER 1

---

## HVL Common Code Base

---

Python common code base to control devices high voltage research devices, in particular, as used in Christian Franck's High Voltage Lab (HVL), D-ITET, ETH.

- Free software: GNU General Public License v3

- **Documentation:**

- if you're planning to develop start w/ reading “CONTRIBUTING.rst”, otherwise either
  - read [HVL CCB documentation at RTD](#), or
  - install *Sphinx* and *sphinx\_rtd\_theme* Python packages and locally build docs on Windows in git-bash by running:

```
$ ./make.sh docs
```

from a shell with Make installed by running:

```
$ make docs
```

The target index HTML (“docs/\_build/html/index.html”) will open automatically in your Web browser.

## 1.1 Features

Manage experiments with `ExperimentManager` instance controlling one or more of the following devices:

- a MBW973 SF6 Analyzer / dew point mirror over a serial connection (COM-ports)
- a LabJack (T7-PRO) device using a LabJack LJM Library for communication

- a Schneider Electric ILS2T stepper motor drive over Modbus TCP
- a Elektro-Automatik PSI9000 DC power supply using VISA over TCP for communication
- a Rhode & Schwarz RTO 1024 oscilloscope using VISA interface over TCP : : INSTR
- a state-of-the-art HVL in-house Supercube device variants using an OPC UA client
- a Heinzinger Digital Interface I/II and a Heinzinger PNC power supply over a serial connection
- a passively Q-switched Pulsed Laser and a laser attenuator from CryLas over a serial connection
- a Newport SMC100PP single axis motion controller for 2-phase stepper motors over a serial connection
- a Pfeiffer TPG controller (TPG 25x, TPG 26x and TPG 36x) for Compact pressure Gauges

## 1.2 Credits

This package was created with [Cookiecutter](#) and the [audreyr/cookiecutter-pypackage](#) project template.

# CHAPTER 2

---

## Installation

---

### 2.1 Stable release

To install HVL Common Code Base, run this command in your terminal:

```
$ pip install hvl_ccb
```

This is the preferred method to install HVL Common Code Base, as it will always install the most recent stable release. If you don't have `pip` installed, this [Python installation guide](#) can guide you through the process.

### 2.2 From sources

The sources for HVL Common Code Base can be downloaded from the [GitLab](#) repo.

You can either clone the repository:

```
$ git clone git@gitlab.com:ethz_hvl/hvl_ccb.git
```

Or download the [tarball](#):

```
$ curl -OL https://gitlab.com/ethz_hvl/hvl_ccb/-/archive/master/hvl_ccb.tar.gz
```

Once you have a copy of the source, you can install it with:

```
$ python setup.py install
```



# CHAPTER 3

---

## Usage

---

To use HVL Common Code Base in a project:

```
import hvl_ccb
```



# CHAPTER 4

---

## API Documentation

---

### 4.1 hvl\_ccb package

#### 4.1.1 Subpackages

##### [hvl\\_ccb.comm package](#)

###### **Submodules**

###### [hvl\\_ccb.comm.base module](#)

Module with base classes for communication protocols.

**class** `hvl_ccb.comm.base.CommunicationProtocol(config)`  
Bases: `hvl_ccb.configuration.ConfigurationMixin, abc.ABC`

Communication protocol abstract base class.

Specifies the methods to implement for communication protocol, as well as implements some default settings and checks.

**access\_lock = None**

Access lock to use with context manager when accessing the communication protocol (thread safety)

**close()**

Close the communication protocol

**open()**

Open communication protocol

###### [hvl\\_ccb.comm.labjack\\_ljm module](#)

Communication protocol for LabJack using the LJM Library. Originally developed and tested for LabJack T7-PRO.

Makes use of the LabJack LJM Library Python wrapper. This wrapper needs an installation of the LJM Library for Windows, Mac OS X or Linux. Go to: <https://labjack.com/support/software/installers/ljm> and <https://labjack.com/support/software/examples/ljm/python>

**class** `hvl_ccb.comm.labjack_ljm.LJMCommunication`(*configuration*)

Bases: `hvl_ccb.comm.base.CommunicationProtocol`

Communication protocol implementing the LabJack LJM Library Python wrapper.

**close()** → None

Close the communication port.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**is\_open**

Flag indicating if the communication port is open.

**Returns** *True* if the port is open, otherwise *False*

**open()** → None

Open the communication port.

**read\_name**(\**names*, *return\_num\_type*: `Type[numbers.Real]`) = <class 'float'> →

`Union[numbers.Real, Sequence[numbers.Real]]`

Read one or more input numeric values by name.

### Parameters

- **names** – one or more names to read out from the LabJack
- **return\_num\_type** – optional numeric type specification for return values; by default `float`.

**Returns** answer of the LabJack, either single number or multiple numbers in a sequence, respectively, when one or multiple names to read were given

**Raises** `TypeError` – if read value of type not compatible with *return\_num\_type*

**write\_name**(*name*: str, *value*: `numbers.Real`) → None

Write one value to a named output.

### Parameters

- **name** – String or with name of LabJack IO
- **value** – is the value to write to the named IO port

**write\_names**(*name\_value\_dict*: `Dict[str, numbers.Real]`) → None

Write more than one value at once to named outputs.

**Parameters** **name\_value\_dict** – is a dictionary with string names of LabJack IO as keys and corresponding numeric values

**class** `hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig`(*device\_type*: `Union[str, hvl_ccb.dev.labjack.DeviceType]`) = `'ANY'`, *connection\_type*: `Union[str, hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig]` = `'ANY'`, *identifier*: str = `'ANY'`)

Bases: `object`

Configuration dataclass for *LJMCommunication*.

```
class ConnectionType(*args, **kwds)
    Bases: hvl_ccb.utils.enum.AutoNumberNameEnum

    LabJack connection type.

    ANY = 1
    ETHERNET = 4
    TCP = 3
    USB = 2
    WIFI = 5

class DeviceType(*args, **kwds)
    Bases: hvl_ccb.utils.enum.AutoNumberNameEnum

    LabJack device types.

    Can be also looked up by ambiguous Product ID (p_id) or by instance name: `python
    LabJackDeviceType(4) is LabJackDeviceType('T4')` 

    ANY = 1
    T4 = 2
    T7 = 3
    T7_PRO = 4

    get_by_p_id = <bound method DeviceType.get_by_p_id of <aenum 'DeviceType'>>
```

**clean\_values()** → None

Performs value checks on *device\_type* and *connection\_type*.

**connection\_type** = 'ANY'

Can be either string or of enum *ConnectionType*.

**device\_type** = 'ANY'

Can be either string 'ANY', 'T7\_PRO', 'T7', 'T4', or of enum *DeviceType*.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**identifier** = 'ANY'

The identifier specifies information for the connection to be used. This can be an IP address, serial number, or device name. See the LabJack docs ( <https://labjack.com/support/software/api/ljm/function-reference/ljmopens/identifier-parameter>) for more information.

**is\_configdataclass** = True

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception hvl\_ccb.comm.labjack\_ljm.LJMCommunicationError**

Bases: Exception

Errors coming from LJMCommunication.

### **hvl\_ccb.comm.modbus\_tcp module**

Communication protocol for modbus TCP ports. Makes use of the pymodbus library.

**class hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication(configuration)**  
Bases: *hvl\_ccb.comm.base.CommunicationProtocol*

Implements the Communication Protocol for modbus TCP.

**close()**

Close the Modbus TCP connection.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**open()** → None

Open the Modbus TCP connection.

**Raises ModbusTcpConnectionFailedException** – if the connection fails.

**read\_holding\_registers(address: int, count: int)** → List[int]

Read specified number of register starting with given address and return the values from each register.

#### Parameters

- **address** – address of the first register
- **count** – count of registers to read

**Returns** list of *int* values

**read\_input\_registers(address: int, count: int)** → List[int]

Read specified number of register starting with given address and return the values from each register in a list.

#### Parameters

- **address** – address of the first register
- **count** – count of registers to read

**Returns** list of *int* values

**write\_registers(address: int, values: Union[List[int], int])**

Write values from the specified address forward.

**Parameters**

- **address** – address of the first register
- **values** – list with all values

```
class hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig (host: str, unit: int,  
                                         port: int = 502)
```

Bases: object

Configuration dataclass for *ModbusTcpCommunication*.

**clean\_values()**

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**host** = None

Host is the IP address of the connected device.

**is\_configdataclass** = True

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port** = 502

TCP port

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**unit** = None

Unit number to be used when connecting with Modbus/TCP. Typically this is used when connecting to a relay having Modbus/RTU-connected devices.

```
exception hvl_ccb.comm.modbus_tcp.ModbusTcpConnectionFailedException (string=")
```

Bases: pymodbus.exceptions.ConnectionException

Exception raised when the connection failed.

## **hvl\_ccb.comm.opc module**

Communication protocol implementing an OPC UA connection. This protocol is used to interface with the “Super-cube” PLC from Siemens.

```
class hvl_ccb.comm.opc.OpcUaCommunication(config)
Bases: hvl_ccb.comm.base.CommunicationProtocol

Communication protocol implementing an OPC UA connection. Makes use of the package python-opcua.

close() → None
    Close the connection to the OPC UA server.

static config_cls()
    Return the default configdataclass class

    Returns a reference to the default configdataclass class

init_monitored_nodes(node_id: Union[str, Iterable[str]], ns_index: int) → None
    Initialize monitored nodes.

    Parameters
        • node_id – one or more strings of node IDs.
        • ns_index – the namespace index the nodes belong to.

    Raises OpcUaCommunicationIOError – when protocol was not opened or can't communicate with a OPC UA server

is_open
    Flag indicating if the communication port is open.

    Returns True if the port is open, otherwise False

open() → None
    Open the communication to the OPC UA server.

    Raises OpcUaCommunicationIOError – when communication port cannot be opened.

read(node_id, ns_index)
    Read a value from a node with id and namespace index.

    Parameters
        • node_id – the ID of the node to read the value from
        • ns_index – the namespace index of the node

    Returns the value of the node object.

    Raises OpcUaCommunicationIOError – when protocol was not opened or can't communicate with a OPC UA server

write(node_id, ns_index, value) → None
    Write a value to a node with name name.

    Parameters
        • node_id – the id of the node to write the value to.
        • ns_index – the namespace index of the node.
        • value – the value to write.

    Raises OpcUaCommunicationIOError – when protocol was not opened or can't communicate with a OPC UA server
```

```
class hvl_ccb.comm.opc.OpcUaCommunicationConfig(host: str, endpoint_name: str, port: int = 4840, sub_handler: hvl_ccb.comm.opc.OpcUaSubHandler = <hvl_ccb.comm.opc.OpcUaSubHandler object>, update_period: int = 500)
```

Bases: object

Configuration dataclass for OPC UA Communciation.

**clean\_values()**

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**endpoint\_name = None**

Endpoint of the OPC server, this is a path like ‘OPCUA/SimulationServer’

**force\_value(*fieldname, value*)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**host = None**

Hostname or IP-Address of the OPC UA server.

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port = 4840**

Port of the OPC UA server to connect to.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**sub\_handler = <hvl\_ccb.comm.opc.OpcUaSubHandler object>**

object to use for handling subscriptions.

**update\_period = 500**

Update period for generating datachange events in OPC UA [milli seconds]

**exception hvl\_ccb.comm.opc.OpcUaCommunicationIOError**

Bases: OSError

OPC-UA communication I/O error.

```
class hvl_ccb.comm.opc.OpcUaSubHandler
Bases: object
```

Base class for subscription handling of OPC events and data change events. Override methods from this class to add own handling capabilities.

To receive events from server for a subscription data\_change and event methods are called directly from receiving thread. Do not do expensive, slow or network operation there. Create another thread if you need to do such a thing.

```
datachange_notification(node, val, data)
event_notification(event)
```

## hvl\_ccb.comm.serial module

Communication protocol for serial ports. Makes use of the `pySerial` library.

```
class hvl_ccb.comm.serial.SerialCommunication(configuration)
Bases: hvl_ccb.comm.base.CommunicationProtocol
```

Implements the Communication Protocol for serial ports.

```
ENCODING = 'utf-8'
```

```
UNICODE_HANDLING = 'replace'
```

```
close()
```

Close the serial connection.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
is_open
```

Flag indicating if the serial port is open.

**Returns** `True` if the serial port is open, otherwise `False`

```
open()
```

Open the serial connection.

**Raises** `SerialCommunicationIOError` – when communication port cannot be opened.

```
read_bytes(size: int = 1) → bytes
```

Read the specified number of bytes from the serial port. The input buffer may hold additional data afterwards.

This method uses `self.access_lock` to ensure thread-safety.

**Parameters** `size` – number of bytes to read

**Returns** Bytes read from the serial port; `b''` if there was nothing to read.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

```
read_text() → str
```

Read one line of text from the serial port. The input buffer may hold additional data afterwards, since only one line is read.

This method uses `self.access_lock` to ensure thread-safety.

**Returns** String read from the serial port; `''` if there was nothing to read.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**write\_bytes** (`data: bytes`) → int  
Write bytes to the serial port.

This method uses `self.access_lock` to ensure thread-safety.

**Parameters** `data` – data to write to the serial port

**Returns** number of bytes written

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**write\_text** (`text: str`)

Write text to the serial port. The text is encoded and terminated by the configured terminator.

This method uses `self.access_lock` to ensure thread-safety.

**Parameters** `text` – Text to send to the port.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**class** `hvl_ccb.comm.serial.SerialCommunicationBytesize(*args, **kwds)`

Bases: `hvl_ccb.utils.enum.ValueEnum`

Serial communication bytesize.

**EIGHTBITS** = 8

**FIVEBITS** = 5

**SEVENBITS** = 7

**SIXBITS** = 6

**class** `hvl_ccb.comm.serial.SerialCommunicationConfig(port: str, baudrate: int, parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity], stopbits: Union[int, float, hvl_ccb.comm.serial.SerialCommunicationStopbits], bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize], terminator: bytes = b'\r\n', timeout: Union[int, float] = 2)`

Bases: `object`

Configuration dataclass for `SerialCommunication`.

**Bytesize**

alias of `SerialCommunicationBytesize`

**Parity**

alias of `SerialCommunicationParity`

**Stopbits**

alias of `SerialCommunicationStopbits`

**baudrate = None**

Baudrate of the serial port

**bytesize = None**

Size of a byte, 5 to 8

**clean\_values()**

```
create_serial_port () → serial.serialposix.Serial
    Create a serial port instance according to specification in this configuration

    Returns Closed serial port instance

force_value (fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing
    after value has been forced on fieldname.

    Parameters
        • fieldname – name of the field
        • value – value to assign

is_configdataclass = True

classmethod keys () → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.

classmethod optional_defaults () → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified
    on instantiation.

    Returns a list of strings containing all optional keys.

parity = None
    Parity to be used for the connection.

port = None
    Port is a string referring to a COM-port (e.g. 'COM3') or a URL. The full list of capabilities is found on the pyserial documentation.

classmethod required_keys () → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on
    instantiation.

    Returns a list of strings containing all required keys.

stopbits = None
    Stopbits setting, can be 1, 1.5 or 2.

terminator = b'\r\n'
    The terminator character. Typically this is b'\r\n' or b'\n', but can also be b'\r' or other combinations.

timeout = 2
    Timeout in seconds for the serial port

exception hvl_ccb.comm.serial.SerialCommunicationIOError
    Bases: OSError

    Serial communication related I/O errors.

class hvl_ccb.comm.serial.SerialCommunicationParity (*args, **kwds)
    Bases: hvl\_ccb.utils.enum.ValueEnum

    Serial communication parity.

EVEN = 'E'
MARK = 'M'
```

```

NAMEs = {'E': 'Even', 'M': 'Mark', 'N': 'None', 'O': 'Odd', 'S': 'Space'}
NONE = 'N'
ODD = 'O'
SPACE = 'S'

class hvl_ccb.comm.serial.SerialCommunicationStopbits(*args, **kwds)
    Bases: hvl_ccb.utils.enum.ValueEnum
        Serial communication stopbits.

    ONE = 1
    ONE_POINT_FIVE = 1.5
    TWO = 2

```

## **hvl\_ccb.comm.visa module**

Communication protocol for VISA. Makes use of the pyvisa library. The backend can be NI-Visa or pyvisa-py.

Information on how to install a VISA backend can be found here: [https://pyvisa.readthedocs.io/en/master/getting\\_nivisa.html](https://pyvisa.readthedocs.io/en/master/getting_nivisa.html)

So far only TCPIP SOCKET and TCPIP INSTR interfaces are supported.

```

class hvl_ccb.comm.visa.VisaCommunication(configuration)
    Bases: hvl_ccb.comm.base.CommunicationProtocol
        Implements the Communication Protocol for VISA / SCPI.

    MULTI_COMMANDS_MAX = 5
        The maximum of commands that can be sent in one round is 5 according to the VISA standard.

    MULTI_COMMANDS_SEPARATOR = ';'
        The character to separate two commands is ; according to the VISA standard.

    WAIT_AFTER_WRITE = 0.08
        Small pause in seconds to wait after write operations, allowing devices to really do what we tell them before continuing with further tasks.

    close() → None
        Close the VISA connection and invalidates the handle.

    static config_cls() → Type[hvl_ccb.comm.visa.VisaCommunicationConfig]
        Return the default configdataclass class.

            Returns a reference to the default configdataclass class

    open() → None
        Open the VISA connection and create the resource.

    query(*commands) → Union[str, Tuple[str, ...]]
        A combination of write(message) and read.

            Parameters commands – list of commands

            Returns list of values

            Raises VisaCommunicationError – when connection was not started, or when trying to issue too many commands at once.

```

**spoll()** → int

Execute serial poll on the device. Reads the status byte register STB. This is a fast function that can be executed periodically in a polling fashion.

**Returns** integer representation of the status byte

**Raises** `VisaCommunicationError` – when connection was not started

**write(\*commands)** → None

Write commands. No answer is read or expected.

**Parameters** `commands` – one or more commands to send

**Raises** `VisaCommunicationError` – when connection was not started

```
class hvl_ccb.comm.visa.VisaCommunicationConfig(host: str, interface_type: Union[str, hvl_ccb.comm.visa.VisaCommunicationConfig.InterfaceType], board: int = 0, port: int = 5025, timeout: int = 5000, chunk_size: int = 204800, open_timeout: int = 1000, write_termination: str = '\n', read_termination: str = '\n', visa_backend: str = '')
```

Bases: `object`

*VisaCommunication* configuration dataclass.

**class InterfaceType(\*args, \*\*kwds)**

Bases: `hvl_ccb.utils.enum.AutoNumberNameEnum`

Supported VISA Interface types.

**TCPIP\_INSTR = 2**

VXI-11 protocol

**TCPIP\_SOCKET = 1**

VISA-Raw protocol

**address(host: str, port: int = None, board: int = None) → str**

Address string specific to the VISA interface type.

**Parameters**

- `host` – host IP address
- `port` – optional TCP port
- `board` – optional board number

**Returns** address string

**address**

Address string depending on the VISA protocol's configuration.

**Returns** address string corresponding to current configuration

**board = 0**

Board number is typically 0 and comes from old bus systems.

**chunk\_size = 204800**

Chunk size is the allocated memory for read operations. The standard is 20kB, and is increased per default here to 200kB. It is specified in bytes.

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**host = None**

IP address of the VISA device. DNS names are currently unsupported.

**interface\_type = None**

Interface type of the VISA connection, being one of `InterfaceType`.

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**open\_timeout = 1000**

Timeout for opening the connection, in milli seconds.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port = 5025**

TCP port, standard is 5025.

**read\_termination = '\n'**

Read termination character.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**timeout = 5000**

Timeout for commands in milli seconds.

**visa\_backend = ''**

Specifies the path to the library to be used with PyVISA as a backend. Defaults to None, which is NI-VISA (if installed), or pyvisa-py (if NI-VISA is not found). To force the use of pyvisa-py, specify '@py' here.

**write\_termination = '\n'**

Write termination character.

**exception hvl\_ccb.comm.visa.VisaCommunicationError**

Bases: Exception

Base class for VisaCommunication errors.

## Module contents

Communication protocols subpackage.

## **hvl\_ccb.dev package**

### **Subpackages**

#### **hvl\_ccb.dev.supercube package**

### **Submodules**

#### **hvl\_ccb.dev.supercube.base module**

Base classes for the Supercube device.

```
class hvl_ccb.dev.supercube.base.SupercubeBase (com, dev_config=None)  
    Bases: hvl_ccb.dev.base.SingleCommDevice
```

Base class for Supercube variants.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
static default_com_cls()
```

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

```
get_cee16_socket() → bool
```

Read the on-state of the IEC CEE16 three-phase power socket.

**Returns** the on-state of the CEE16 power socket

```
get_door_status(door: int) → hvl_ccb.dev.supercube.constants.DoorStatus
```

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**Parameters** `door` – the door number (1..3)

**Returns** the door status

```
get_earthing_manual(number: int) → bool
```

**TODO:** Test `get_earthing_manual` with device

Get the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters** `number` – number of the earthing stick (1..6)

**Returns** earthing stick manual status

```
get_earthing_status(number: int) → hvl_ccb.dev.supercube.constants.EarthingStickStatus
```

Get the status of an earthing stick, whether it is closed, open or undefined (moving).

**Parameters** `number` – number of the earthing stick (1..6)

**Returns** earthing stick status

```
get_measurement_ratio(channel: int) → float
```

**TODO:** test `get_measurement_ratio` with device

Get the set measurement ratio of an AC/DC analog input channel. Every input channel has a divider ratio assigned during setup of the Supercube system. This ratio can be read out.

**Parameters** `channel` – number of the input channel (1..4)

**Returns** the ratio

**get\_measurement\_voltage** (*channel: int*) → float  
**TODO:** test **get\_measurement\_voltage** with device

Get the measured voltage of an analog input channel. The voltage read out here is already scaled by the configured divider ratio.

**Parameters** **channel** – number of the input channel (1..4)

**Returns** measured voltage

**get\_status** () → hvl\_ccb.dev.supercube.constants.SafetyStatus  
Get the safety circuit status of the Supercube. :return: the safety status of the supercube's state machine.

**get\_support\_input** (*port: int, contact: int*) → bool  
Get the state of a support socket input.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital input read state

**get\_support\_output** (*port: int, contact: int*) → bool  
Get the state of a support socket output.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital output read state

**get\_t13\_socket** (*port: int*) → bool  
Read the state of a SEV T13 power socket.

**Parameters** **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*

**Returns** on-state of the power socket

**operate** (*state: bool*) → None  
Set operate state. If the state is RedReady, this will turn on the high voltage and close the safety switches.

**Parameters** **state** – set operate state

**quit\_error** () → None  
Quits errors that are active on the Supercube.

**read** (*node\_id: str*)  
Local wrapper for the OPC UA communication protocol read method.

**Parameters** **node\_id** – the id of the node to read.

**Returns** the value of the variable

**ready** (*state: bool*) → None  
Set ready state. Ready means locket safety circuit, red lamps, but high voltage still off.

**Parameters** **state** – set ready state

**set\_cee16\_socket** (*state: bool*) → None  
Switch the IEC CEE16 three-phase power socket on or off.

**Parameters** **state** – desired on-state of the power socket

**Raises ValueError** – if state is not of type bool

**set\_earthing\_manual** (*number: int, manual: bool*) → None

**TODO:** Test set\_earthing\_manual with device

Set the manual status of an earthing stick. If an earthing stick is set to manual, it stays closed even if the system is in states RedReady or RedOperate.

#### Parameters

- **number** – number of the earthing stick (1..6)
- **manual** – earthing stick manual status (True or False)

**set\_remote\_control** (*state: bool*) → None

**TODO:** test set\_remote\_control with device

Enable or disable remote control for the Supercube. This will effectively display a message on the touch-screen HMI.

**Parameters state** – desired remote control state

**set\_support\_output** (*port: int, contact: int, state: bool*) → None

Set the state of a support output socket.

#### Parameters

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **state** – is the desired state of the support output

**set\_support\_output\_impulse** (*port: int, contact: int, duration: float = 0.2, pos\_pulse: bool = True*) → None

Issue an impulse of a certain duration on a support output contact. The polarity of the pulse (On-wait-Off or Off-wait-On) is specified by the pos\_pulse argument.

This function is blocking.

#### Parameters

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **duration** – is the length of the impulse in seconds
- **pos\_pulse** – is True, if the pulse shall be HIGH, False if it shall be LOW

**set\_t13\_socket** (*port: int, state: bool*) → None

Set the state of a SEV T13 power socket.

#### Parameters

- **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*
- **state** – is the desired on-state of the socket

**start ()** → None

Starts the device. Sets the root node for all OPC read and write commands to the Siemens PLC object node which holds all our relevant objects and variables.

**stop ()** → None

Stop the Supercube device. Deactivates the remote control and closes the communication protocol.

**write** (*node\_id, value*) → None

Local wrapper for the OPC UA communication protocol write method.

**Parameters**

- **node\_id** – the id of the node to read
- **value** – the value to write to the variable

```
class hvl_ccb.dev.supercube.base.SupercubeConfiguration (namespace_index: int = 3)
```

Bases: object

Configuration dataclass for the Supercube devices.

**clean\_values()**

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass** = True

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**namespace\_index** = 3

Namespace of the OPC variables, typically this is 3 (coming from Siemens)

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
class hvl_ccb.dev.supercube.base.SupercubeOpcUaCommunication (config)
```

Bases: *hvl\_ccb.comm.opc.OpcUaCommunication*

Communication protocol specification for Supercube devices.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig (host:
    str, end-
    point_name:
    str, port:
    int =
    4840,
    sub_handler:
    hvl_ccb.comm.opc.OpcUaSubHandler
    =
    <hvl_ccb.dev.supercube.base.SupercubeObject>,
    up-
    date_period:
    int = 500)
```

Bases: *hvl\_ccb.comm.opc.OpcUaCommunicationConfig*

Communication protocol configuration for OPC UA, specifications for the Supercube devices.

**force\_value** (*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**sub\_handler** = <*hvl\_ccb.dev.supercube.base.SupercubeSubscriptionHandler object*>

Subscription handler for data change events

**class** hvl\_ccb.dev.supercube.base.**SupercubeSubscriptionHandler**

Bases: *hvl\_ccb.comm.opc.OpcUaSubHandler*

OPC Subscription handler for datachange events and normal events specifically implemented for the Supercube devices.

**datachange\_notification** (*node: opcua.common.node.Node, val, data*)

In addition to the standard operation (debug logging entry of the datachange), alarms are logged at INFO level using the alarm text.

**Parameters**

- **node** – the node object that triggered the datachange event

- **val** – the new value
- **data** –

## **hvl\_ccb.dev.supercube.constants module**

Constants, variable names for the Supercube OPC-connected devices.

```
class hvl_ccb.dev.supercube.constants.AlarmText (*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum
```

This enumeration contains textual representations for all error classes (stop, warning and message) of the Supercube system. Use the [AlarmText.get\(\)](#) method to retrieve the enum of an alarm number.

```
Alarm1 = 'STOP Emergency Stop 1'
Alarm10 = 'STOP Earthing stick 2 error while opening'
Alarm11 = 'STOP Earthing stick 3 error while opening'
Alarm12 = 'STOP Earthing stick 4 error while opening'
Alarm13 = 'STOP Earthing stick 5 error while opening'
Alarm14 = 'STOP Earthing stick 6 error while opening'
Alarm15 = 'STOP Earthing stick 1 error while closing'
Alarm16 = 'STOP Earthing stick 2 error while closing'
Alarm17 = 'STOP Earthing stick 3 error while closing'
Alarm18 = 'STOP Earthing stick 4 error while closing'
Alarm19 = 'STOP Earthing stick 5 error while closing'
Alarm2 = 'STOP Emergency Stop 2'
Alarm20 = 'STOP Earthing stick 6 error while closing'
Alarm21 = 'STOP Safety fence 1'
Alarm22 = 'STOP Safety fence 2'
Alarm23 = 'STOP OPC connection error'
Alarm24 = 'STOP Grid power failure'
Alarm25 = 'STOP UPS failure'
Alarm26 = 'STOP 24V PSU failure'
Alarm3 = 'STOP Emergency Stop 3'
Alarm4 = 'STOP Safety Switch 1 error'
Alarm41 = 'WARNING Door 1: Use earthing rod!'
Alarm42 = 'MESSAGE Door 1: Earthing rod is still in setup.'
Alarm43 = 'WARNING Door 2: Use earthing rod!'
Alarm44 = 'MESSAGE Door 2: Earthing rod is still in setup.'
Alarm45 = 'WARNING Door 3: Use earthing rod!'
Alarm46 = 'MESSAGE Door 3: Earthing rod is still in setup.'
```

```
Alarm47 = 'MESSAGE UPS charge < 85%'  
Alarm48 = 'MESSAGE UPS running on battery'  
Alarm5 = 'STOP Safety Switch 2 error'  
Alarm6 = 'STOP Door 1 lock supervision'  
Alarm7 = 'STOP Door 2 lock supervision'  
Alarm8 = 'STOP Door 3 lock supervision'  
Alarm9 = 'STOP Earthing stick 1 error while opening'  
get = <bound method AlarmText.get of <aenum 'AlarmText'>>  
not_defined = 'NO ALARM TEXT DEFINED'  
  
class hvl_ccb.dev.supercube.constants.Alarms(*args, **kwds)  
Bases: hvl_ccb.dev.supercube.constants._AlarmsBase  
  
Alarms enumeration containing all variable NodeID strings for the alarm array.  
  
Alarm1 = '"DB_Alarm_HMI"."Alarm1"'  
Alarm10 = '"DB_Alarm_HMI"."Alarm10"'  
Alarm100 = '"DB_Alarm_HMI"."Alarm100"'  
Alarm101 = '"DB_Alarm_HMI"."Alarm101"'  
Alarm102 = '"DB_Alarm_HMI"."Alarm102"'  
Alarm103 = '"DB_Alarm_HMI"."Alarm103"'  
Alarm104 = '"DB_Alarm_HMI"."Alarm104"'  
Alarm105 = '"DB_Alarm_HMI"."Alarm105"'  
Alarm106 = '"DB_Alarm_HMI"."Alarm106"'  
Alarm107 = '"DB_Alarm_HMI"."Alarm107"'  
Alarm108 = '"DB_Alarm_HMI"."Alarm108"'  
Alarm109 = '"DB_Alarm_HMI"."Alarm109"'  
Alarm11 = '"DB_Alarm_HMI"."Alarm11"'  
Alarm110 = '"DB_Alarm_HMI"."Alarm110"'  
Alarm111 = '"DB_Alarm_HMI"."Alarm111"'  
Alarm112 = '"DB_Alarm_HMI"."Alarm112"'  
Alarm113 = '"DB_Alarm_HMI"."Alarm113"'  
Alarm114 = '"DB_Alarm_HMI"."Alarm114"'  
Alarm115 = '"DB_Alarm_HMI"."Alarm115"'  
Alarm116 = '"DB_Alarm_HMI"."Alarm116"'  
Alarm117 = '"DB_Alarm_HMI"."Alarm117"'  
Alarm118 = '"DB_Alarm_HMI"."Alarm118"'  
Alarm119 = '"DB_Alarm_HMI"."Alarm119"'  
Alarm12 = '"DB_Alarm_HMI"."Alarm12"'
```

```
Alarm120 = '"DB_Alarm_HMI"."Alarm120"'
Alarm121 = '"DB_Alarm_HMI"."Alarm121"'
Alarm122 = '"DB_Alarm_HMI"."Alarm122"'
Alarm123 = '"DB_Alarm_HMI"."Alarm123"'
Alarm124 = '"DB_Alarm_HMI"."Alarm124"'
Alarm125 = '"DB_Alarm_HMI"."Alarm125"'
Alarm126 = '"DB_Alarm_HMI"."Alarm126"'
Alarm127 = '"DB_Alarm_HMI"."Alarm127"'
Alarm128 = '"DB_Alarm_HMI"."Alarm128"'
Alarm129 = '"DB_Alarm_HMI"."Alarm129"'
Alarm13 = '"DB_Alarm_HMI"."Alarm13"'
Alarm130 = '"DB_Alarm_HMI"."Alarm130"'
Alarm131 = '"DB_Alarm_HMI"."Alarm131"'
Alarm132 = '"DB_Alarm_HMI"."Alarm132"'
Alarm133 = '"DB_Alarm_HMI"."Alarm133"'
Alarm134 = '"DB_Alarm_HMI"."Alarm134"'
Alarm135 = '"DB_Alarm_HMI"."Alarm135"'
Alarm136 = '"DB_Alarm_HMI"."Alarm136"'
Alarm137 = '"DB_Alarm_HMI"."Alarm137"'
Alarm138 = '"DB_Alarm_HMI"."Alarm138"'
Alarm139 = '"DB_Alarm_HMI"."Alarm139"'
Alarm14 = '"DB_Alarm_HMI"."Alarm14"'
Alarm140 = '"DB_Alarm_HMI"."Alarm140"'
Alarm141 = '"DB_Alarm_HMI"."Alarm141"'
Alarm142 = '"DB_Alarm_HMI"."Alarm142"'
Alarm143 = '"DB_Alarm_HMI"."Alarm143"'
Alarm144 = '"DB_Alarm_HMI"."Alarm144"'
Alarm145 = '"DB_Alarm_HMI"."Alarm145"'
Alarm146 = '"DB_Alarm_HMI"."Alarm146"'
Alarm147 = '"DB_Alarm_HMI"."Alarm147"'
Alarm148 = '"DB_Alarm_HMI"."Alarm148"'
Alarm149 = '"DB_Alarm_HMI"."Alarm149"'
Alarm15 = '"DB_Alarm_HMI"."Alarm15"'
Alarm150 = '"DB_Alarm_HMI"."Alarm150"'
Alarm151 = '"DB_Alarm_HMI"."Alarm151"'
Alarm16 = '"DB_Alarm_HMI"."Alarm16"'
```

```
Alarm17 = '"DB_Alarm_HMI"."Alarm17"'
Alarm18 = '"DB_Alarm_HMI"."Alarm18"'
Alarm19 = '"DB_Alarm_HMI"."Alarm19"'
Alarm2 = '"DB_Alarm_HMI"."Alarm2"'
Alarm20 = '"DB_Alarm_HMI"."Alarm20"'
Alarm21 = '"DB_Alarm_HMI"."Alarm21"'
Alarm22 = '"DB_Alarm_HMI"."Alarm22"'
Alarm23 = '"DB_Alarm_HMI"."Alarm23"'
Alarm24 = '"DB_Alarm_HMI"."Alarm24"'
Alarm25 = '"DB_Alarm_HMI"."Alarm25"'
Alarm26 = '"DB_Alarm_HMI"."Alarm26"'
Alarm27 = '"DB_Alarm_HMI"."Alarm27"'
Alarm28 = '"DB_Alarm_HMI"."Alarm28"'
Alarm29 = '"DB_Alarm_HMI"."Alarm29"'
Alarm3 = '"DB_Alarm_HMI"."Alarm3"'
Alarm30 = '"DB_Alarm_HMI"."Alarm30"'
Alarm31 = '"DB_Alarm_HMI"."Alarm31"'
Alarm32 = '"DB_Alarm_HMI"."Alarm32"'
Alarm33 = '"DB_Alarm_HMI"."Alarm33"'
Alarm34 = '"DB_Alarm_HMI"."Alarm34"'
Alarm35 = '"DB_Alarm_HMI"."Alarm35"'
Alarm36 = '"DB_Alarm_HMI"."Alarm36"'
Alarm37 = '"DB_Alarm_HMI"."Alarm37"'
Alarm38 = '"DB_Alarm_HMI"."Alarm38"'
Alarm39 = '"DB_Alarm_HMI"."Alarm39"'
Alarm4 = '"DB_Alarm_HMI"."Alarm4"'
Alarm40 = '"DB_Alarm_HMI"."Alarm40"'
Alarm41 = '"DB_Alarm_HMI"."Alarm41"'
Alarm42 = '"DB_Alarm_HMI"."Alarm42"'
Alarm43 = '"DB_Alarm_HMI"."Alarm43"'
Alarm44 = '"DB_Alarm_HMI"."Alarm44"'
Alarm45 = '"DB_Alarm_HMI"."Alarm45"'
Alarm46 = '"DB_Alarm_HMI"."Alarm46"'
Alarm47 = '"DB_Alarm_HMI"."Alarm47"'
Alarm48 = '"DB_Alarm_HMI"."Alarm48"'
Alarm49 = '"DB_Alarm_HMI"."Alarm49"'
```

```
Alarm5 = '"DB_Alarm_HMI"."Alarm5"'
Alarm50 = '"DB_Alarm_HMI"."Alarm50"'
Alarm51 = '"DB_Alarm_HMI"."Alarm51"'
Alarm52 = '"DB_Alarm_HMI"."Alarm52"'
Alarm53 = '"DB_Alarm_HMI"."Alarm53"'
Alarm54 = '"DB_Alarm_HMI"."Alarm54"'
Alarm55 = '"DB_Alarm_HMI"."Alarm55"'
Alarm56 = '"DB_Alarm_HMI"."Alarm56"'
Alarm57 = '"DB_Alarm_HMI"."Alarm57"'
Alarm58 = '"DB_Alarm_HMI"."Alarm58"'
Alarm59 = '"DB_Alarm_HMI"."Alarm59"'
Alarm6 = '"DB_Alarm_HMI"."Alarm6"'
Alarm60 = '"DB_Alarm_HMI"."Alarm60"'
Alarm61 = '"DB_Alarm_HMI"."Alarm61"'
Alarm62 = '"DB_Alarm_HMI"."Alarm62"'
Alarm63 = '"DB_Alarm_HMI"."Alarm63"'
Alarm64 = '"DB_Alarm_HMI"."Alarm64"'
Alarm65 = '"DB_Alarm_HMI"."Alarm65"'
Alarm66 = '"DB_Alarm_HMI"."Alarm66"'
Alarm67 = '"DB_Alarm_HMI"."Alarm67"'
Alarm68 = '"DB_Alarm_HMI"."Alarm68"'
Alarm69 = '"DB_Alarm_HMI"."Alarm69"'
Alarm7 = '"DB_Alarm_HMI"."Alarm7"'
Alarm70 = '"DB_Alarm_HMI"."Alarm70"'
Alarm71 = '"DB_Alarm_HMI"."Alarm71"'
Alarm72 = '"DB_Alarm_HMI"."Alarm72"'
Alarm73 = '"DB_Alarm_HMI"."Alarm73"'
Alarm74 = '"DB_Alarm_HMI"."Alarm74"'
Alarm75 = '"DB_Alarm_HMI"."Alarm75"'
Alarm76 = '"DB_Alarm_HMI"."Alarm76"'
Alarm77 = '"DB_Alarm_HMI"."Alarm77"'
Alarm78 = '"DB_Alarm_HMI"."Alarm78"'
Alarm79 = '"DB_Alarm_HMI"."Alarm79"'
Alarm8 = '"DB_Alarm_HMI"."Alarm8"'
Alarm80 = '"DB_Alarm_HMI"."Alarm80"'
Alarm81 = '"DB_Alarm_HMI"."Alarm81"'
```

```
Alarm82 = '"DB_Alarm_HMI"."Alarm82"'
Alarm83 = '"DB_Alarm_HMI"."Alarm83"'
Alarm84 = '"DB_Alarm_HMI"."Alarm84"'
Alarm85 = '"DB_Alarm_HMI"."Alarm85"'
Alarm86 = '"DB_Alarm_HMI"."Alarm86"'
Alarm87 = '"DB_Alarm_HMI"."Alarm87"'
Alarm88 = '"DB_Alarm_HMI"."Alarm88"'
Alarm89 = '"DB_Alarm_HMI"."Alarm89"'
Alarm9 = '"DB_Alarm_HMI"."Alarm9"'
Alarm90 = '"DB_Alarm_HMI"."Alarm90"'
Alarm91 = '"DB_Alarm_HMI"."Alarm91"'
Alarm92 = '"DB_Alarm_HMI"."Alarm92"'
Alarm93 = '"DB_Alarm_HMI"."Alarm93"'
Alarm94 = '"DB_Alarm_HMI"."Alarm94"'
Alarm95 = '"DB_Alarm_HMI"."Alarm95"'
Alarm96 = '"DB_Alarm_HMI"."Alarm96"'
Alarm97 = '"DB_Alarm_HMI"."Alarm97"'
Alarm98 = '"DB_Alarm_HMI"."Alarm98"'
Alarm99 = '"DB_Alarm_HMI"."Alarm99"'

class hvl_ccb.dev.supercube.constants.BreakdownDetection(*args, **kwds)
    Bases: hvl\_ccb.utils.enum.ValueEnum

    Node ID strings for the breakdown detection.

    TODO: these variable NodeIDs are not tested and/or correct yet.

    activated = 'Ix_Allg_Breakdown_activated'  

        Boolean read-only variable indicating whether breakdown detection and fast switchoff is enabled in the system or not.

    reset = 'Qx_Allg_Breakdown_reset'  

        Boolean writable variable to reset the fast switch-off. Toggle to re-enable.

    triggered = 'Ix_Allg_Breakdown_triggered'  

        Boolean read-only variable telling whether the fast switch-off has triggered. This can also be seen using the safety circuit state, therefore no method is implemented to read this out directly.

class hvl_ccb.dev.supercube.constants.Door(*args, **kwds)
    Bases: hvl\_ccb.utils.enum.ValueEnum

    Variable NodeID strings for doors.

    get = <bound method Door.get of <aenum 'Door'>>
    status_1 = '"DB_Safety_Circuit"."Türe 1"."si_HMI_status_door"'
    status_2 = '"DB_Safety_Circuit"."Türe 2"."si_HMI_status_door"'
    status_3 = '"DB_Safety_Circuit"."Türe 3"."si_HMI_status_door"'
```

---

```
class hvl_ccb.dev.supercube.constants.DoorStatus (*args, **kwds)
Bases: aenum.IntEnum

Possible status values for doors.

closed = 2
    Door is closed, but not locked.

error = 4
    Door has an error or was opened in locked state (either with emergency stop or from the inside).

inactive = 0
    not enabled in Supercube HMI setup, this door is not supervised.

locked = 3
    Door is closed and locked (safe state).

open = 1
    Door is open.

class hvl_ccb.dev.supercube.constants.EarthingStick (*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum

Variable NodeID strings for all earthing stick statuses (read-only integer) and writable booleans for setting the earthing in manual mode.

manual = <bound method EarthingStick.manual of <aenum 'EarthingStick'>>
manual_1 = '"DB_Safety_Circuit"."Erdpeitsche 1"."sx_earthing_manually"'
manual_2 = '"DB_Safety_Circuit"."Erdpeitsche 2"."sx_earthing_manually"'
manual_3 = '"DB_Safety_Circuit"."Erdpeitsche 3"."sx_earthing_manually"'
manual_4 = '"DB_Safety_Circuit"."Erdpeitsche 4"."sx_earthing_manually"'
manual_5 = '"DB_Safety_Circuit"."Erdpeitsche 5"."sx_earthing_manually"'
manual_6 = '"DB_Safety_Circuit"."Erdpeitsche 6"."sx_earthing_manually"'
status = <bound method EarthingStick.status of <aenum 'EarthingStick'>>
status_1 = '"DB_Safety_Circuit"."Erdpeitsche 1"."si_HMI_Status"'
status_2 = '"DB_Safety_Circuit"."Erdpeitsche 2"."si_HMI_Status"'
status_3 = '"DB_Safety_Circuit"."Erdpeitsche 3"."si_HMI_Status"'
status_4 = '"DB_Safety_Circuit"."Erdpeitsche 4"."si_HMI_Status"'
status_5 = '"DB_Safety_Circuit"."Erdpeitsche 5"."si_HMI_Status"'
status_6 = '"DB_Safety_Circuit"."Erdpeitsche 6"."si_HMI_Status"'

class hvl_ccb.dev.supercube.constants.EarthingStickStatus (*args, **kwds)
Bases: aenum.IntEnum

Status of an earthing stick. These are the possible values in the status integer e.g. in EarthingStick.status\_1.

closed = 1
    Earthing is closed (safe).

error = 3
    Earthing is in error, e.g. when the stick did not close correctly or could not open.
```

**inactive = 0**

Earthing stick is deselected and not enabled in safety circuit. To get out of this state, the earthing has to be enabled in the Supercube HMI setup.

**open = 2**

Earthing is open (not safe).

**class hvl\_ccb.dev.supercube.constants.Errors (\*args, \*\*kwds)**

Bases: *hvl\_ccb.utils.enum.ValueEnum*

Variable NodeID strings for information regarding error, warning and message handling.

**message = '"DB\_Meldepuffer"."Hinweis\_aktiv"'**

Boolean read-only variable telling if a message is active.

**quit = '"DB\_Meldepuffer"."Quittierbutton"'**

Writable boolean for the error quit button.

**stop = '"DB\_Meldepuffer"."Stop\_aktiv"'**

Boolean read-only variable telling if a stop is active.

**warning = '"DB\_Meldepuffer"."Warnung\_aktiv"'**

Boolean read-only variable telling if a warning is active.

**class hvl\_ccb.dev.supercube.constants.GeneralSockets (\*args, \*\*kwds)**

Bases: *hvl\_ccb.utils.enum.ValueEnum*

NodeID strings for the power sockets (3x T13 and 1xCEE16).

**cee16 = '"Qx\_Allg\_Socket\_CEE16"'**

CEE16 socket (writeable boolean).

**t13\_1 = '"Qx\_Allg\_Socket\_T13\_1"'**

SEV T13 socket No. 1 (writable boolean).

**t13\_2 = '"Qx\_Allg\_Socket\_T13\_2"'**

SEV T13 socket No. 2 (writable boolean).

**t13\_3 = '"Qx\_Allg\_Socket\_T13\_3"'**

SEV T13 socket No. 3 (writable boolean).

**class hvl\_ccb.dev.supercube.constants.GeneralSupport (\*args, \*\*kwds)**

Bases: *hvl\_ccb.utils.enum.ValueEnum*

NodeID strings for the support inputs and outputs.

**in\_1\_1 = '"Ix\_Allg\_Support1\_1"'**

**in\_1\_2 = '"Ix\_Allg\_Support1\_2"'**

**in\_2\_1 = '"Ix\_Allg\_Support2\_1"'**

**in\_2\_2 = '"Ix\_Allg\_Support2\_2"'**

**in\_3\_1 = '"Ix\_Allg\_Support3\_1"'**

**in\_3\_2 = '"Ix\_Allg\_Support3\_2"'**

**in\_4\_1 = '"Ix\_Allg\_Support4\_1"'**

**in\_4\_2 = '"Ix\_Allg\_Support4\_2"'**

**in\_5\_1 = '"Ix\_Allg\_Support5\_1"'**

**in\_5\_2 = '"Ix\_Allg\_Support5\_2"'**

**in\_6\_1 = '"Ix\_Allg\_Support6\_1"'**

```

in_6_2 = '"Ix_Allg_Support6_2"'
input = <bound method GeneralSupport.input of <aenum 'GeneralSupport'>>
out_1_1 = '"Qx_Allg_Support1_1"'
out_1_2 = '"Qx_Allg_Support1_2"'
out_2_1 = '"Qx_Allg_Support2_1"'
out_2_2 = '"Qx_Allg_Support2_2"'
out_3_1 = '"Qx_Allg_Support3_1"'
out_3_2 = '"Qx_Allg_Support3_2"'
out_4_1 = '"Qx_Allg_Support4_1"'
out_4_2 = '"Qx_Allg_Support4_2"'
out_5_1 = '"Qx_Allg_Support5_1"'
out_5_2 = '"Qx_Allg_Support5_2"'
out_6_1 = '"Qx_Allg_Support6_1"'
out_6_2 = '"Qx_Allg_Support6_2"'

output = <bound method GeneralSupport.output of <aenum 'GeneralSupport'>>

```

```

class hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio(*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum
```

Variable NodeID strings for the measurement input scaling ratios. These ratios are defined in the Supercube HMI setup and are provided in the python module here to be able to read them out, allowing further calculations.

TODO: these variable nodeIDs are not tested and/or correct yet.

```

get = <bound method MeasurementsDividerRatio.get of <aenum 'MeasurementsDividerRatio'>>
input_1 = 'Ir_Measure_DividerRatio_1'
input_2 = 'Ir_Measure_DividerRatio_2'
input_3 = 'Ir_Measure_DividerRatio_3'
input_4 = 'Ir_Measure_DividerRatio_4'

class hvl_ccb.dev.supercube.constants.MeasurementsScaledInput(*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum
```

Variable NodeID strings for the four analog BNC inputs for measuring voltage. The voltage returned in these variables is already scaled with the set ratio, which can be read using the variables in *MeasurementsDividerRatio*.

TODO: these variable NodeIDs are not tested and/or correct yet.

```

get = <bound method MeasurementsScaledInput.get of <aenum 'MeasurementsScaledInput'>>
input_1 = 'Qr_Measure_Input1_scaledVoltage'
input_2 = 'Qr_Measure_Input2_scaledVoltage'
input_3 = 'Qr_Measure_Input3_scaledVoltage'
input_4 = 'Qr_Measure_Input4_scaledVoltage'
```

```
class hvl_ccb.dev.supercube.constants.OpcControl (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings for supervision of the OPC connection from the controlling workstation to the Supercube.

TODO: this variable nodeID string is not tested and/or correct yet.

active = 'Ix_OPc_active'
writable boolean to enable OPC remote control and display a message window on the Supercube HMI.

class hvl_ccb.dev.supercube.constants.Power (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings concerning power data.

TODO: these variable NodeIDs are not tested and/or correct yet, they don't exist yet on Supercube side.

current_primary = 'Qr_Power_FU_actualCurrent'
Primary current in ampere, measured by the frequency converter. (read-only)

frequency = 'Ir_Power_FU_Frequency'
Frequency converter output frequency. (read-only)

setup = 'Qi_Power_Setup'
Power setup that is configured using the Supercube HMI. The value corresponds to the ones in PowerSetup. (read-only)

voltage_max = 'Iw_Power_maxVoltage'
Maximum voltage allowed by the current experimental setup. (read-only)

voltage_primary = 'Qr_Power_FU_actualVoltage'
Primary voltage in volts, measured by the frequency converter at its output. (read-only)

voltage_slope = 'Ir_Power_dUDt'
Voltage slope in V/s.

voltage_target = 'Ir_Power_TargetVoltage'
Target voltage setpoint in V.

class hvl_ccb.dev.supercube.constants.PowerSetup (*args, **kwds)
Bases: aenum.IntEnum

Possible power setups corresponding to the value of variable Power.setup.

AC_DoubleStage_150kV = 4
AC voltage with two MWB transformers, one at 100kV and the other at 50kV, resulting in a total maximum voltage of 150kV.

AC_DoubleStage_200kV = 5
AC voltage with two MWB transformers both at 100kV, resulting in a total maximum voltage of 200kV

AC_SingleStage_100kV = 3
AC voltage with MWB transformer set to 100kV maximum voltage.

AC_SingleStage_50kV = 2
AC voltage with MWB transformer set to 50kV maximum voltage.

DC_DoubleStage_280kV = 8
DC voltage with two AC transformers set to 100kV AC each, resulting in 280kV DC in total (or a single stage transformer with Greinacher voltage doubling rectifier)

DC_SingleStage_140kV = 7
DC voltage with one AC transformer set to 100kV AC, resulting in 140kV DC
```

**External = 1**

External power supply fed through blue CEE32 input using isolation transformer and safety switches of the Supercube, or using an external safety switch attached to the Supercube Type B.

**Internal = 6**

Internal usage of the frequency converter, controlling to the primary voltage output of the supercube itself (no measurement transformer used)

**NoPower = 0**

No safety switches, use only safety components (doors, fence, earthing...) without any power.

```
class hvl_ccb.dev.supercube.constants.Safety(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum
```

NodeID strings for the basic safety circuit status and green/red switches “ready” and “operate”.

**status = 'DB\_Safety\_Circuit"."si\_safe\_status'**

Status is a read-only integer containing the state number of the supercube-internal state machine. The values correspond to numbers in *SafetyStatus*.

**switchto\_operate = '"DB\_Safety\_Circuit"."sx\_safe\_switchto\_operate"**

Writable boolean for switching to Red Operate (locket, HV on) state.

**switchto\_ready = '"DB\_Safety\_Circuit"."sx\_safe\_switchto\_ready"**

Writable boolean for switching to Red Ready (locked, HV off) state.

```
class hvl_ccb.dev.supercube.constants.SafetyStatus(*args, **kwds)
Bases: aenum.IntEnum
```

Safety status values that are possible states returned from `hvl_ccb.dev.supercube.base.Supercube.get_status()`. These values correspond to the states of the Supercube’s safety circuit statemachine.

**Error = 6**

System is in error mode.

**GreenNotReady = 1**

System is safe, lamps are green and some safety elements are not in place such that it cannot be switched to red currently.

**GreenReady = 2**

System is safe and all safety elements are in place to be able to switch to *ready*.

**Initializing = 0**

System is initializing or booting.

**QuickStop = 5**

Fast turn off triggered and switched off the system. Reset FSO to go back to a normal state.

**RedOperate = 4**

System is locked in red state and in *operate* mode, i.e. high voltage on.

**RedReady = 3**

System is locked in red state and *ready* to go to *operate* mode.

```
class hvl_ccb.dev.supercube.constants.SupercubeOpcEndpoint(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum
```

OPC Server Endpoint strings for the supercube variants.

**A = 'Supercube Typ A'**

**B = 'Supercube Typ B'**

```
hvl_ccb.dev.supercube.constants.T13_SOCKET_PORTS = (1, 2, 3)
Port numbers of SEV T13 power socket
```

## **hvl\_ccb.dev.supercube.typ\_a module**

Supercube Typ A module.

```
class hvl_ccb.dev.supercube.typ_a.SuperCubeAOpcUaCommunication(config)
```

Bases: *hvl\_ccb.dev.supercube.base.SuperCubeOpcUaCommunication*

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.supercube.typ_a.SuperCubeAOpcUaConfiguration(host: str, end-
```

point\_name: str

= 'Supercube

Typ A', port:

int = 4840,

sub\_handler:

*hvl\_ccb.comm.opc.OpcUaSubHandler*

=

<*hvl\_ccb.dev.supercube.base.SuperCube*

object at

0x7f2f86c81c50>,

update\_period:

int = 500)

Bases: *hvl\_ccb.dev.supercube.base.SuperCubeOpcUaCommunicationConfig*

```
endpoint_name = 'Supercube Typ A'
```

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

### **Parameters**

- **fieldname** – name of the field
- **value** – value to assign

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

---

**class** hvl\_ccb.dev.supercube.typ\_a.**SupercubeWithFU**(*com, dev\_config=None*)  
 Bases: *hvl\_ccb.dev.supercube.base.SupercubeBase*

Variant A of the Supercube with frequency converter.

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**fso\_reset()** → None

**TODO: test fso\_reset with device**

Reset the fast switch off circuitry to go back into normal state and allow to re-enable operate mode.

**get\_frequency()** → float

**TODO: test get\_frequency with device**

Read the electrical frequency of the current Supercube setup.

**Returns** the frequency in Hz

**get\_fso\_active()** → bool

**TODO: test get\_fso\_active with device**

Get the state of the fast switch off functionality. Returns True if it is enabled, False otherwise.

**Returns** state of the FSO functionality

**get\_max\_voltage()** → float

**TODO: test get\_max\_voltage with device**

Reads the maximum voltage of the setup and returns in V.

**Returns** the maximum voltage of the setup in V.

**get\_power\_setup()** → hvl\_ccb.dev.supercube.constants.PowerSetup

**TODO: test get\_power\_setup with device**

Return the power setup selected in the Supercube's settings.

**Returns** the power setup

**get\_primary\_current()** → float

**TODO: get\_primary\_current with device**

Read the current primary current at the output of the frequency converter ( before transformer).

**Returns** primary current in A

**get\_primary\_voltage()** → float

**TODO: test get\_primary\_voltage with device**

Read the current primary voltage at the output of the frequency converter ( before transformer).

**Returns** primary voltage in V

**get\_target\_voltage()** → float

**TODO: test get\_target\_voltage with device**

Gets the current setpoint of the output voltage value in V. This is not a measured value but is the corresponding function to *set\_target\_voltage()*.

**Returns** the setpoint voltage in V.

**set\_slope** (*slope: float*) → None

**TODO:** test **set\_slope** with device

Sets the dV/dt slope of the Supercube frequency converter to a new value in V/s.

**Parameters** **slope** – voltage slope in V/s (0..15)

**set\_target\_voltage** (*volt\_v: float*) → None

**TODO:** test **set\_target\_voltage** with device

Set the output voltage to a defined value in V.

**Parameters** **volt\_v** – the desired voltage in V

## [hvl\\_ccb.dev.supercube.typ\\_b module](#)

Supercube Typ B module.

**class** hvl\_ccb.dev.supercube.typ\_b.**SupercubeB** (*com, dev\_config=None*)

Bases: [hvl\\_ccb.dev.supercube.base.SupercubeBase](#)

Variant B of the Supercube without frequency converter but external safety switches.

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**class** hvl\_ccb.dev.supercube.typ\_b.**SupercubeBOpcUaCommunication** (*config*)

Bases: [hvl\\_ccb.dev.supercube.base.SupercubeOpcUaCommunication](#)

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class** hvl\_ccb.dev.supercube.typ\_b.**SupercubeBOpcUaConfiguration** (*host: str, endpoint\_name: str = 'Supercube Typ B', port: int = 4840, sub\_handler: hvl\_ccb.comm.opc.OpcUaSubHandler = <hvl\_ccb.dev.supercube.base.SupercubeBOpcUaCommunication object at 0x7f2f86c81c50>, update\_period: int = 500*)

Bases: [hvl\\_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig](#)

**endpoint\_name = 'Supercube Typ B'**

**force\_value** (*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field

- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## Module contents

Supercube package with implementation for system versions from 2019 on (new concept with hard-PLC Siemens S7-1500 as CPU).

### hvl\_ccb.dev.supercube2015 package

#### Submodules

##### hvl\_ccb.dev.supercube2015.base module

Base classes for the Supercube device.

**exception hvl\_ccb.dev.supercube2015.base.InvalidSupercubeStatusError**

Bases: Exception

Exception raised when supercube has invalid status.

**class hvl\_ccb.dev.supercube2015.base.Supercube2015Base (com, dev\_config=None)**

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Base class for Supercube variants.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**get\_cee16\_socket()** → bool

Read the on-state of the IEC CEE16 three-phase power socket.

**Returns** the on-state of the CEE16 power socket

**get\_door\_status(door: int)** → hvl\_ccb.dev.supercube2015.constants.DoorStatus

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**Parameters** **door** – the door number (1..3)

**Returns** the door status

**get\_earthing\_manual** (*number: int*) → bool

Get the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters** **number** – number of the earthing stick (1..6)

**Returns** earthing stick manual status

**get\_earthing\_status** (*number: int*) → int

Get the status of an earthing stick, whether it is closed, open or undefined (moving).

**Parameters** **number** – number of the earthing stick (1..6)

**Returns** earthing stick status; see constants.EarthingStickStatus

**get\_measurement\_ratio** (*channel: int*) → float

Get the set measurement ratio of an AC/DC analog input channel. Every input channel has a divider ratio assigned during setup of the Supercube system. This ratio can be read out.

**Attention:** Supercube 2015 does not have a separate ratio for every analog input. Therefore there is only one ratio for *channel* = 1.

**Parameters** **channel** – number of the input channel (1..4)

**Returns** the ratio

**get\_measurement\_voltage** (*channel: int*) → float

Get the measured voltage of an analog input channel. The voltage read out here is already scaled by the configured divider ratio.

**Attention:** In contrast to the *new* Supercube, the old one returns here the input voltage read at the ADC. It is not scaled by a factor.

**Parameters** **channel** – number of the input channel (1..4)

**Returns** measured voltage

**get\_status** () → int

Get the safety circuit status of the Supercube. :return: the safety status of the supercube's state machine;

see constants.SafetyStatus.

**get\_support\_input** (*port: int, contact: int*) → bool

Get the state of a support socket input.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital input read state

**get\_support\_output** (*port: int, contact: int*) → bool

Get the state of a support socket output.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital output read state

**get\_t13\_socket** (*port: int*) → bool

Read the state of a SEV T13 power socket.

**Parameters** **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*

**Returns** on-state of the power socket

**operate** (*state: bool*) → None

Set operate state. If the state is RedReady, this will turn on the high voltage and close the safety switches.

**Parameters** **state** – set operate state

**quit\_error** () → None

Quits errors that are active on the Supercube.

**read** (*node\_id: str*)

Local wrapper for the OPC UA communication protocol read method.

**Parameters** **node\_id** – the id of the node to read.

**Returns** the value of the variable

**ready** (*state: bool*) → None

Set ready state. Ready means locket safety circuit, red lamps, but high voltage still off.

**Parameters** **state** – set ready state

**set\_cee16\_socket** (*state: bool*) → None

Switch the IEC CEE16 three-phase power socket on or off.

**Parameters** **state** – desired on-state of the power socket

**Raises** **ValueError** – if state is not of type bool

**set\_earthing\_manual** (*number: int, manual: bool*) → None

Set the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters**

- **number** – number of the earthing stick (1..6)
- **manual** – earthing stick manual status (True or False)

**set\_remote\_control** (*state: bool*) → None

Enable or disable remote control for the Supercube. This will effectively display a message on the touch-screen HMI.

**Parameters** **state** – desired remote control state

**set\_support\_output** (*port: int, contact: int, state: bool*) → None

Set the state of a support output socket.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **state** – is the desired state of the support output

**set\_support\_output\_impulse** (*port: int, contact: int, duration: float = 0.2, pos\_pulse: bool = True*) → None

Issue an impulse of a certain duration on a support output contact. The polarity of the pulse (On-wait-Off or Off-wait-On) is specified by the pos\_pulse argument.

This function is blocking.

### Parameters

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **duration** – is the length of the impulse in seconds
- **pos\_pulse** – is True, if the pulse shall be HIGH, False if it shall be LOW

**set\_t13\_socket** (*port: int, state: bool*) → None

Set the state of a SEV T13 power socket.

### Parameters

- **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*
- **state** – is the desired on-state of the socket

**start** () → None

Starts the device. Sets the root node for all OPC read and write commands to the Siemens PLC object node which holds all our relevant objects and variables.

**stop** () → None

Stop the Supercube device. Deactivates the remote control and closes the communication protocol.

**write** (*node\_id, value*) → None

Local wrapper for the OPC UA communication protocol write method.

### Parameters

- **node\_id** – the id of the node to read
- **value** – the value to write to the variable

**class** hvl\_ccb.dev.supercube2015.base.**SupercubeConfiguration** (*namespace\_index: int = 7*)

Bases: object

Configuration dataclass for the Supercube devices.

**clean\_values** ()

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**force\_value** (*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**namespace\_index = 7**

Namespace of the OPC variables, typically this is 3 (coming from Siemens)

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**class** hvl\_ccb.dev.supercube2015.base.**SupercubeOpcUaCommunication**(*config*)

Bases: *hvl\_ccb.comm.opc.OpcUaCommunication*

Communication protocol specification for Supercube devices.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class** hvl\_ccb.dev.supercube2015.base.**SupercubeOpcUaCommunicationConfig**(*host*:

*str*,

*end-*

*point\_name*:

*str*,

*port*:

*int*

=

4845,

*sub\_handler*:

*hvl\_ccb.comm.opc.OpcUaSu*

=

<*hvl\_ccb.dev.supercube2015*

*ob-*

*ject*>,

*up-*

*date\_period*:

*int*

=

500)

Bases: *hvl\_ccb.comm.opc.OpcUaCommunicationConfig*

Communication protocol configuration for OPC UA, specifications for the Supercube devices.

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
port = 4845
```

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
sub_handler = <hvl_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler object>
```

Subscription handler for data change events

```
class hvl_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler
```

Bases: *hvl\_ccb.comm.opc.OpcUaSubHandler*

OPC Subscription handler for datachange events and normal events specifically implemented for the Supercube devices.

```
datachange_notification(node: opcua.common.node.Node, val, data)
```

In addition to the standard operation (debug logging entry of the datachange), alarms are logged at INFO level using the alarm text.

### Parameters

- **node** – the node object that triggered the datachange event
- **val** – the new value
- **data** –

## **hvl\_ccb.dev.supercube2015.constants module**

Constants, variable names for the Supercube OPC-connected devices.

```
class hvl_ccb.dev.supercube2015.constants.AlarmText(*args, **kwds)
```

Bases: *hvl\_ccb.utils.enum.ValueEnum*

This enumeration contains textual representations for all error classes (stop, warning and message) of the Supercube system. Use the *AlarmText.get()* method to retrieve the enum of an alarm number.

```
Alarm0 = 'No Alarm.'  
Alarm1 = 'STOP Safety switch 1 error'  
Alarm10 = 'STOP Earthing stick 2 error'  
Alarm11 = 'STOP Earthing stick 3 error'  
Alarm12 = 'STOP Earthing stick 4 error'  
Alarm13 = 'STOP Earthing stick 5 error'  
Alarm14 = 'STOP Earthing stick 6 error'  
Alarm17 = 'STOP Source switch error'  
Alarm19 = 'STOP Fence 1 error'  
Alarm2 = 'STOP Safety switch 2 error'
```

```

Alarm20 = 'STOP Fence 2 error'
Alarm21 = 'STOP Control error'
Alarm22 = 'STOP Power outage'
Alarm3 = 'STOP Emergency Stop 1'
Alarm4 = 'STOP Emergency Stop 2'
Alarm5 = 'STOP Emergency Stop 3'
Alarm6 = 'STOP Door 1 lock supervision'
Alarm7 = 'STOP Door 2 lock supervision'
Alarm8 = 'STOP Door 3 lock supervision'
Alarm9 = 'STOP Earthing stick 1 error'
get = <bound method AlarmText.get of <aenum 'AlarmText'>>
not_defined = 'NO ALARM TEXT DEFINED'

class hvl_ccb.dev.supercube2015.constants.BreakdownDetection (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Node ID strings for the breakdown detection.

activated = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.Breakdowndetection.connect'
Boolean read-only variable indicating whether breakdown detection and fast switchoff is enabled in the system or not.

reset = 'hvl-ipc.WINAC.Support6OutA'
Boolean writable variable to reset the fast switch-off. Toggle to re-enable.

triggered = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.Breakdowndetection.triggered'
Boolean read-only variable telling whether the fast switch-off has triggered. This can also be seen using the safety circuit state, therefore no method is implemented to read this out directly.

class hvl_ccb.dev.supercube2015.constants.DoorStatus (*args, **kwds)
Bases: aenum.IntEnum

Possible status values for doors.

closed = 2
Door is closed, but not locked.

error = 4
Door has an error or was opened in locked state (either with emergency stop or from the inside).

inactive = 0
not enabled in Supercube HMI setup, this door is not supervised.

locked = 3
Door is closed and locked (safe state).

open = 1
Door is open.

class hvl_ccb.dev.supercube2015.constants.EarthingStick (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings for all earthing stick statuses (read-only integer) and writable booleans for setting the earthing in manual mode.

manual = <bound method EarthingStick.manual of <aenum 'EarthingStick'>>

```

```
manual_1 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.MANUAL'
manual_2 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.MANUAL'
manual_3 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.MANUAL'
manual_4 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.MANUAL'
manual_5 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.MANUAL'
manual_6 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.MANUAL'
status_1_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.CLOSE'
status_1_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.CONNECT'
status_1_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.OPEN'
status_2_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.CLOSE'
status_2_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.CONNECT'
status_2_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.OPEN'
status_3_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.CLOSE'
status_3_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.CONNECT'
status_3_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.OPEN'
status_4_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.CLOSE'
status_4_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.CONNECT'
status_4_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.OPEN'
status_5_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.CLOSE'
status_5_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.CONNECT'
status_5_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.OPEN'
status_6_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.CLOSE'
status_6_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.CONNECT'
status_6_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.OPEN'
status_closed = <bound method EarthingStick.status_closed of <aenum 'EarthingStick'>>
status_connected = <bound method EarthingStick.status_connected of <aenum 'EarthingStick'>>
status_open = <bound method EarthingStick.status_open of <aenum 'EarthingStick'>>

class hvl_ccb.dev.supercube2015.constants.EarthingStickStatus(*args, **kwds)
Bases: aenum.IntEnum

Status of an earthing stick. These are the possible values in the status integer e.g. in EarthingStick.
status_1.

closed = 1
    Earthing is closed (safe).

error = 3
    Earthing is in error, e.g. when the stick did not close correctly or could not open.

inactive = 0
    Earthing stick is deselected and not enabled in safety circuit. To get out of this state, the earthing has to be
    enabled in the Supercube HMI setup.
```

```

open = 2
    Earthing is open (not safe).

class hvl_ccb.dev.supercube2015.constants.Errors (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings for information regarding error, warning and message handling.

quit = 'hvl-ipc.WINAC.SYSTEMSTATE.Faultconfirmation'
    Writable boolean for the error quit button.

stop = 'hvl-ipc.WINAC.SYSTEMSTATE.ERROR'
    Boolean read-only variable telling if a stop is active.

stop_number = 'hvl-ipc.WINAC.SYSTEMSTATE.Errornumber'

class hvl_ccb.dev.supercube2015.constants.GeneralSockets (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

NodeID strings for the power sockets (3x T13 and 1xCEE16).

cee16 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.CEE16'
    CEE16 socket (writable boolean).

t13_1 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.T13_1'
    SEV T13 socket No. 1 (writable boolean).

t13_2 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.T13_2'
    SEV T13 socket No. 2 (writable boolean).

t13_3 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.T13_3'
    SEV T13 socket No. 3 (writable boolean).

class hvl_ccb.dev.supercube2015.constants.GeneralSupport (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

NodeID strings for the support inputs and outputs.

in_1_1 = 'hvl-ipc.WINAC.Support1InA'
in_1_2 = 'hvl-ipc.WINAC.Support1InB'
in_2_1 = 'hvl-ipc.WINAC.Support2InA'
in_2_2 = 'hvl-ipc.WINAC.Support2InB'
in_3_1 = 'hvl-ipc.WINAC.Support3InA'
in_3_2 = 'hvl-ipc.WINAC.Support3InB'
in_4_1 = 'hvl-ipc.WINAC.Support4InA'
in_4_2 = 'hvl-ipc.WINAC.Support4InB'
in_5_1 = 'hvl-ipc.WINAC.Support5InA'
in_5_2 = 'hvl-ipc.WINAC.Support5InB'
in_6_1 = 'hvl-ipc.WINAC.Support6InA'
in_6_2 = 'hvl-ipc.WINAC.Support6InB'

input = <bound method GeneralSupport.input of <aenum 'GeneralSupport'>>
out_1_1 = 'hvl-ipc.WINAC.Support1OutA'
out_1_2 = 'hvl-ipc.WINAC.Support1OutB'

```

```
out_2_1 = 'hvl-ipc.WINAC.Support2OutA'
out_2_2 = 'hvl-ipc.WINAC.Support2OutB'
out_3_1 = 'hvl-ipc.WINAC.Support3OutA'
out_3_2 = 'hvl-ipc.WINAC.Support3OutB'
out_4_1 = 'hvl-ipc.WINAC.Support4OutA'
out_4_2 = 'hvl-ipc.WINAC.Support4OutB'
out_5_1 = 'hvl-ipc.WINAC.Support5OutA'
out_5_2 = 'hvl-ipc.WINAC.Support5OutB'
out_6_1 = 'hvl-ipc.WINAC.Support6OutA'
out_6_2 = 'hvl-ipc.WINAC.Support6OutB'
output = <bound method GeneralSupport.output of <aenum 'GeneralSupport'>>

class hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio(*args,
**kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings for the measurement input scaling ratios. These ratios are defined in the Supercube HMI setup and are provided in the python module here to be able to read them out, allowing further calculations.

get = <bound method MeasurementsDividerRatio.get of <aenum 'MeasurementsDividerRatio'>>
input_1 = 'hvl-ipc.WINAC.SYSTEM_INTERN.DivididerRatio'

class hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInput(*args,
**kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings for the four analog BNC inputs for measuring voltage. The voltage returned in these variables is already scaled with the set ratio, which can be read using the variables in MeasurementsDividerRatio.

get = <bound method MeasurementsScaledInput.get of <aenum 'MeasurementsScaledInput'>>
input_1 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI1Volt'
input_2 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI2Volt'
input_3 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI3Volt'
input_4 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI4Volt'

class hvl_ccb.dev.supercube2015.constants.Power(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings concerning power data.

current_primary = 'hvl-ipc.WINAC.SYSTEM_INTERN.FUCurrentprim'
    Primary current in ampere, measured by the frequency converter. (read-only)

frequency = 'hvl-ipc.WINAC.FU.Frequency'
    Frequency converter output frequency. (read-only)

setup = 'hvl-ipc.WINAC.FU.TrafoSetup'
    Power setup that is configured using the Supercube HMI. The value corresponds to the ones in PowerSetup. (read-only)

voltage_max = 'hvl-ipc.WINAC.FU.maxVoltagekV'
    Maximum voltage allowed by the current experimental setup. (read-only)
```

```

voltage_primary = 'hvl-ipc.WINAC.SYSTEM_INTERN.FUVoltageprim'
    Primary voltage in volts, measured by the frequency converter at its output. (read-only)

voltage_slope = 'hvl-ipc.WINAC.FU.dUDt_-1'
    Voltage slope in V/s.

voltage_target = 'hvl-ipc.WINAC.FU.SOLL'
    Target voltage setpoint in V.

class hvl_ccb.dev.supercube2015.constants.PowerSetup(*args, **kwds)
Bases: aenum.IntEnum

    Possible power setups corresponding to the value of variable Power.setup.

AC_DoubleStage_150kV = 3
    AC voltage with two MWB transformers, one at 100kV and the other at 50kV, resulting in a total maximum voltage of 150kV.

AC_DoubleStage_200kV = 4
    AC voltage with two MWB transformers both at 100kV, resulting in a total maximum voltage of 200kV

AC_SingleStage_100kV = 2
    AC voltage with MWB transformer set to 100kV maximum voltage.

AC_SingleStage_50kV = 1
    AC voltage with MWB transformer set to 50kV maximum voltage.

DC_DoubleStage_280kV = 7
    DC voltage with two AC transformers set to 100kV AC each, resulting in 280kV DC in total (or a single stage transformer with Greinacher voltage doubling rectifier)

DC_SingleStage_140kV = 6
    DC voltage with one AC transformer set to 100kV AC, resulting in 140kV DC

External = 0
    External power supply fed through blue CEE32 input using isolation transformer and safety switches of the Supercube, or using an external safety switch attached to the Supercube Type B.

Internal = 5
    Internal usage of the frequency converter, controlling to the primary voltage output of the supercube itself (no measurement transformer used)

class hvl_ccb.dev.supercube2015.constants.Safety(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

    NodeID strings for the basic safety circuit status and green/red switches “ready” and “operate”.

status_error = 'hvl-ipc.WINAC.SYSTEMSTATE.ERROR'

status_green = 'hvl-ipc.WINAC.SYSTEMSTATE.GREEN'

status_ready_for_red = 'hvl-ipc.WINAC.SYSTEMSTATE.ReadyForRed'
    Status is a read-only integer containing the state number of the supercube-internal state machine. The values correspond to numbers in SafetyStatus.

status_red = 'hvl-ipc.WINAC.SYSTEMSTATE.RED'

switchto_green = 'hvl-ipc.WINAC.SYSTEMSTATE.GREEN_REQUEST'

switchto_operate = 'hvl-ipc.WINAC.SYSTEMSTATE.switchon'
    Writable boolean for switching to Red Operate (locket, HV on) state.

switchto_ready = 'hvl-ipc.WINAC.SYSTEMSTATE.RED_REQUEST'
    Writable boolean for switching to Red Ready (locked, HV off) state.

```

```
class hvl_ccb.dev.supercube2015.constants.SafetyStatus (*args, **kwds)
Bases: aenum.IntEnum

Safety status values that are possible states returned from hvl_ccb.dev.supercube.base.
Supercube.get_status(). These values correspond to the states of the Supercube's safety circuit
statemachine.

Error = 6
    System is in error mode.

GreenNotReady = 1
    System is safe, lamps are green and some safety elements are not in place such that it cannot be switched
    to red currently.

GreenReady = 2
    System is safe and all safety elements are in place to be able to switch to ready.

Initializing = 0
    System is initializing or booting.

QuickStop = 5
    Fast turn off triggered and switched off the system. Reset FSO to go back to a normal state.

RedOperate = 4
    System is locked in red state and in operate mode, i.e. high voltage on.

RedReady = 3
    System is locked in red state and ready to go to operate mode.

class hvl_ccb.dev.supercube2015.constants.SupercubeOpcEndpoint (*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

OPC Server Endpoint strings for the supercube variants.

A = 'OPC.SimaticNET.S7'
B = 'OPC.SimaticNET.S7'

hvl_ccb.dev.supercube2015.constants.T13_SOCKET_PORTS = (1, 2, 3)
Port numbers of SEV T13 power socket
```

## **hvl\_ccb.dev.supercube2015.typ\_a module**

Supercube Typ A module.

```
class hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU (com,
                                                               dev_config=None)
Bases: hvl_ccb.dev.supercube2015.base.Supercube2015Base

Variant A of the Supercube with frequency converter.

static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

fso_reset() → None
    Reset the fast switch off circuitry to go back into normal state and allow to re-enable operate mode.

get_frequency() → float
    Read the electrical frequency of the current Supercube setup.

    Returns the frequency in Hz
```

---

**get\_fso\_active()** → bool  
Get the state of the fast switch off functionality. Returns True if it is enabled, False otherwise.

**Returns** state of the FSO functionality

**get\_max\_voltage()** → float  
Reads the maximum voltage of the setup and returns in V.

**Returns** the maximum voltage of the setup in V.

**get\_power\_setup()** → hvl\_ccb.dev.supercube2015.constants.PowerSetup  
Return the power setup selected in the Supercube's settings.

**Returns** the power setup

**get\_primary\_current()** → float  
Read the current primary current at the output of the frequency converter ( before transformer).

**Returns** primary current in A

**get\_primary\_voltage()** → float  
Read the current primary voltage at the output of the frequency converter ( before transformer).

**Returns** primary voltage in V

**get\_target\_voltage()** → float  
Gets the current setpoint of the output voltage value in V. This is not a measured value but is the corresponding function to [`set\_target\_voltage\(\)`](#).

**Returns** the setpoint voltage in V.

**set\_slope(*slope: float*)** → None  
Sets the dV/dt slope of the Supercube frequency converter to a new value in V/s.

**Parameters** **slope** – voltage slope in V/s (0..15'000)

**set\_target\_voltage(*volt\_v: float*)** → None  
Set the output voltage to a defined value in V.

**Parameters** **volt\_v** – the desired voltage in V

**class** hvl\_ccb.dev.supercube2015.typ\_a.**SupercubeAOpcUaCommunication**(*config*)  
Bases: [`hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaCommunication`](#)

**static config\_cls()**  
Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.supercube2015.typ_a.SupercubeAOpcUaConfiguration(host:  
    str, end-  
    point_name:  
    str      =  
    'OPC.SimaticNET.S7',  
    port: int  
    = 4845,  
    sub_handler:  
    hvl_ccb.comm.opc.OpcUaSubHand  
    =  
    <hvl_ccb.dev.supercube2015.base.  
    object at  
    0x7f2f8698fbe0>,  
    up-  
    date_period:  
    int = 500)  
Bases: hvl_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig
```

**endpoint\_name** = 'OPC.SimaticNET.S7'

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

## Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## Module contents

Supercube package with implementation for the old system version from 2015 based on Siemens WinAC soft-PLC on an industrial 32bit Windows computer.

## Submodules

## **hvl\_ccb.dev.base module**

Module with base classes for devices.

**class** `hvl_ccb.dev.base.Device` (`dev_config=None`)

Bases: `hvl_ccb.configuration.ConfigurationMixin`, `abc.ABC`

Base class for devices. Implement this class for a concrete device, such as measurement equipment or voltage sources.

Specifies the methods to implement for a device.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**start()** → None

Start or restart this Device. To be implemented in the subclass.

**stop()** → None

Stop this Device. To be implemented in the subclass.

**exception** `hvl_ccb.dev.base.DeviceExistingException`

Bases: `Exception`

Exception to indicate that a device with that name already exists.

**class** `hvl_ccb.dev.base.DeviceSequenceMixin` (`devices: Dict[str, hvl_ccb.dev.base.Device]`)

Bases: `abc.ABC`

Mixin that can be used on a device or other classes to provide facilities for handling multiple devices in a sequence.

**add\_device** (`name: str, device: hvl_ccb.dev.base.Device`) → None

Add a new device to the device sequence.

### Parameters

- **name** – is the name of the device.
- **device** – is the instantiated Device object.

**Raises** `DeviceExistingException` –

**get\_device** (`name: str`) → `hvl_ccb.dev.base.Device`

Get a device by name.

**Parameters** `name` – is the name of the device.

**Returns** the device object from this sequence.

**get\_devices** () → `List[Tuple[str, hvl_ccb.dev.base.Device]]`

Get list of name, device pairs according to current sequence.

**Returns** A list of tuples with name and device each.

**remove\_device** (`name: str`) → `hvl_ccb.dev.base.Device`

Remove a device from this sequence and return the device object.

**Parameters** `name` – is the name of the device.

**Returns** device object or `None` if such device was not in the sequence.

**Raises** `ValueError` – when device with given name was not found

```
start() → None
    Start all devices in this sequence in their added order.

stop() → None
    Stop all devices in this sequence in their reverse order.

class hvl_ccb.dev.base.EmptyConfig
    Bases: object

    Empty configuration dataclass that is the default configuration for a Device.

clean_values()
    Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.
```

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
is_configdataclass = True

classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.
```

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.
```

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.
```

**Returns** a list of strings containing all required keys.

```
class hvl_ccb.dev.base.SingleCommDevice(com, dev_config=None)
    Bases: hvl\_ccb.dev.base.Device, abc.ABC
```

Base class for devices with a single communication protocol.

**com**

Get the communication protocol of this device.

**Returns** an instance of CommunicationProtocol subtype

```
static default_com_cls() → Type[hvl_ccb.comm.base.CommunicationProtocol]
    Get the class for the default communication protocol used with this device.
```

**Returns** the type of the standard communication protocol for this device

**start()** → None

Open the associated communication protocol.

**stop()** → None

Close the associated communication protocol.

## hvl\_ccb.dev.crylas module

Device classes for a CryLas pulsed laser controller and a CryLas laser attenuator, using serial communication.

There are three modes of operation for the laser 1. Laser-internal hardware trigger (default): fixed to 20 Hz and max energy per pulse. 2. Laser-internal software trigger (for diagnosis only). 3. External trigger: required for arbitrary pulse energy or repetition rate. Switch to “external” on the front panel of laser controller for using option 3.

After switching on the laser with `laser_on()`, the system must stabilize for some minutes. Do not apply abrupt changes of pulse energy or repetition rate.

Manufacturer homepage: [https://www.crylas.de/products/pulsed\\_laser.html](https://www.crylas.de/products/pulsed_laser.html)

**class** `hvl_ccb.dev.crylas.CryLasAttenuator`(*com, dev\_config=None*)

Bases: `hvl_ccb.dev.base.SingleCommDevice`

Device class for the CryLas laser attenuator.

**attenuation**

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**set\_attenuation**(*percent: Union[int, float]*) → None

Set the percentage of attenuated light (inverse of set\_transmission). :param percent: percentage of attenuation, number between 0 and 100 :raises ValueError: if param percent not between 0 and 100 :raises SerialCommunicationIOError: when communication port is not opened :raises CryLasAttenuatorError: if the device does not confirm success

**set\_init\_attenuation()**

Sets the attenuation to its configured initial/default value

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**set\_transmission**(*percent: Union[int, float]*) → None

Set the percentage of transmitted light (inverse of set\_attenuation). :param percent: percentage of transmitted light :raises ValueError: if param percent not between 0 and 100 :raises SerialCommunicationIOError: when communication port is not opened :raises CryLasAttenuatorError: if the device does not confirm success

**start()** → None

Open the com, apply the config value ‘init\_attenuation’

**Raises** `SerialCommunicationIOError` – when communication port cannot be opened

**transmission**

**class** `hvl_ccb.dev.crylas.CryLasAttenuatorConfig`(*init\_attenuation: Union[int, float] = 0, response\_sleep\_time: Union[int, float] = 1*)

Bases: `object`

Device configuration dataclass for CryLas attenuator.

**clean\_values()**

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**init\_attenuation** = 0

**is\_configdataclass** = True

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**response\_sleep\_time** = 1

**exception hvl\_ccb.dev.crylas.CryLasAttenuatorError**

Bases: Exception

General error with the CryLas Attenuator.

**class hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunication** (*configuration*)

Bases: *hvl\_ccb.comm.serial.SerialCommunication*

Specific communication protocol implementation for the CryLas attenuator. Already predefines device-specific protocol parameters in config.

**static config\_cls**()

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig (port: str,  
baudrate:  
int      =  
9600,  
parity:  
Union[str,  
hvl_ccb.comm.serial.SerialCommu  
= <Serial-  
Communi-  
cationPar-  
ity.NONE:  
'N',  
stopbits:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionStop-  
bits.ONE:  
1>, bytesize:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionByte-  
size.EIGHTBITS:  
8>, terminator:  
bytes  
= b", timeout:  
Union[int,  
float] =  
3)
```

Bases: [hvl\\_ccb.comm.serial.SerialCommunicationConfig](#)

**baudrate = 9600**

Baudrate for CryLas attenuator is 9600 baud

**bytesize = 8**

One byte is eight bits long

**force\_value (fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys () → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**  
CryLas attenuator does not use parity

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**  
CryLas attenuator uses one stop bit

**terminator = b''**  
No terminator

**timeout = 3**  
use 3 seconds timeout as default

**class hvl\_ccb.dev.crylas.CryLasLaser(com, dev\_config=None)**  
Bases: *hvl\_ccb.dev.base.SingleCommDevice*

CryLas laser controller device class.

**class AnswersShutter(\*args, \*\*kwds)**  
Bases: *aenum.Enum*

Standard answers of the CryLas laser controller to ‘Shutter’ command passed via *com*.

**CLOSED = 'Shutter inaktiv'**

**OPENED = 'Shutter aktiv'**

**class AnswersStatus(\*args, \*\*kwds)**  
Bases: *aenum.Enum*

Standard answers of the CryLas laser controller to ‘STATUS’ command passed via *com*.

**ACTIVE = 'STATUS: Laser active'**

**HEAD = 'STATUS: Head ok'**

**INACTIVE = 'STATUS: Laser inactive'**

**READY = 'STATUS: System ready'**

**TEC1 = 'STATUS: TEC1 Regulation ok'**

**TEC2 = 'STATUS: TEC2 Regulation ok'**

**class LaserStatus(\*args, \*\*kwds)**  
Bases: *aenum.Enum*

Status of the CryLas laser

**READY\_INACTIVE = 1**

**READ\_ACTIVE = 2**

```

UNREADY_INACTIVE = 0
is_inactive
is_ready

class RepetitionRates(*args, **kwds)
Bases: aenum.IntEnum

Repetition rates for the internal software trigger in Hz

HARDWARE = 0
SOFTWARE_INTERNAL_SIXTY = 60
SOFTWARE_INTERNAL_TEN = 10
SOFTWARE_INTERNAL_TWENTY = 20

ShutterStatus
alias of CryLasLaserShutterStatus

close_shutter() → None
Close the laser shutter.

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserError – if success is not confirmed by the device

static config_cls()
Return the default configdataclass class.

Returns a reference to the default configdataclass class

static default_com_cls()
Get the class for the default communication protocol used with this device.

Returns the type of the standard communication protocol for this device

get_pulse_energy_and_rate() → Tuple[int, int]
Use the debug mode, return the measured pulse energy and rate.

Returns (energy in micro joule, rate in Hz)

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserError – if the device does not answer the query

laser_off() → None
Turn the laser off.

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserError – if success is not confirmed by the device

laser_on() → None
Turn the laser on.

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserNotReadyError – if the laser is not ready to be turned on

```

- *CryLasLaserError* – if success is not confirmed by the device

**open\_shutter ()** → None  
Open the laser shutter.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**set\_init\_shutter\_status ()** → None  
Open or close the shutter, to match the configured shutter\_status.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**set\_pulse\_energy (energy: int)** → None  
Sets the energy of pulses (works only with external hardware trigger). Proceed with small energy steps, or the regulation may fail.

**Parameters** **energy** – energy in micro joule

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if the device does not confirm success

**set\_repetition\_rate (rate: Union[int, hvl\_ccb.dev.crylas.CryLasLaser.RepetitionRates])** → None  
Sets the repetition rate of the internal software trigger.

**Parameters** **rate** – frequency (Hz) as an integer

**Raises**

- *ValueError* – if rate is not an accepted value in RepetitionRates Enum
- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**start ()** → None  
Opens the communication protocol and configures the device.

**Raises** *SerialCommunicationIOError* – when communication port cannot be opened

**stop ()** → None  
Stops the device and closes the communication protocol.

**Raises**

- *SerialCommunicationIOError* – if com port is closed unexpectedly
- *CryLasLaserError* – if laser\_off() or close\_shutter() fail

**target\_pulse\_energy**

**update\_laser\_status ()** → None

Update the laser status to *LaserStatus.NOT\_READY* or *LaserStatus.INACTIVE* or *LaserStatus.ACTIVE*.

Note: laser never explicitly says that it is not ready (*LaserStatus.NOT\_READY*) in response to ‘STATUS’ command. It only says that it is ready (heated-up and implicitly inactive/off) or active (on). If it’s not

either of these then the answer is *Answers.HEAD*. Moreover, the only time the laser explicitly says that its status is inactive (*Answers.INACTIVE*) is after issuing a ‘LASER OFF’ command.

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**update\_repetition\_rate()** → None

Query the laser repetition rate.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**update\_shutter\_status()** → None

Update the shutter status (OPENED or CLOSED)

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**update\_target\_pulse\_energy()** → None

Query the laser pulse energy.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**wait\_until\_ready()** → None

Block execution until the laser is ready

**Raises** *CryLasLaserError* – if the polling thread stops before the laser is ready

```
class hvl_ccb.dev.crylas.CryLasLaserConfig(calibration_factor: Union[int, float]
                                              = 4.35, polling_period: Union[int, float] = 5, polling_timeout: Union[int, float] = 300, on_start_wait_until_ready: bool = False, auto_laser_on: bool = True, init_shutter_status: Union[int, hvl_ccb.dev.crylas.CryLasLaserShutterStatus] = <CryLasLaserShutterStatus.CLOSED: 0>)
```

Bases: object

Device configuration dataclass for the CryLas laser controller.

**ShutterStatus**

alias of *CryLasLaserShutterStatus*

**auto\_laser\_on = True**

**calibration\_factor = 4.35**

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field

- **value** – value to assign

```
init_shutter_status = 0
```

```
is_configdataclass = True
```

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
on_start_wait_until_ready = False
```

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
polling_period = 5
```

```
polling_timeout = 300
```

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
exception hvl_ccb.dev.crylas.CryLasLaserError
```

Bases: Exception

General error with the CryLas Laser.

```
exception hvl_ccb.dev.crylas.CryLasLaserNotReadyError
```

Bases: [hvl\\_ccb.dev.crylas.CryLasLaserError](#)

Error when trying to turn on the CryLas Laser before it is ready.

```
class hvl_ccb.dev.crylas.CryLasLaserSerialCommunication(configuration)
```

Bases: [hvl\\_ccb.comm.serial.SerialCommunication](#)

Specific communication protocol implementation for the CryLas laser controller. Already predefines device-specific protocol parameters in config.

```
READ_TEXT_SKIP_PREFIXES = ('>', 'MODE:')
```

Prefixes of lines that are skipped when read from the serial port.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
query(cmd: str, prefix: str, post_cmd: str = None) → str
```

Send a command, then read the com until a line starting with prefix, or an empty line, is found. Returns the line in question.

**Parameters**

- **cmd** – query message to send to the device
- **prefix** – start of the line to look for in the device answer
- **post\_cmd** – optional additional command to send after the query

**Returns** line in question as a string

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**query\_all** (`cmd: str, prefix: str`)

Send a command, then read the com until a line starting with prefix, or an empty line, is found. Returns a list of successive lines starting with prefix.

#### Parameters

- **cmd** – query message to send to the device
- **prefix** – start of the line to look for in the device answer

**Returns** line in question as a string

**Raises** `SerialCommunicationIOError` – when communication port is not opened

```
class hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig(port: str, baudrate: int = 19200, parity: Union[str, hvl_ccb.comm.serial.SerialCommunication] = <SerialCommunicationParity.NONE: 'N'>, stopbits: Union[int, hvl_ccb.comm.serial.SerialCommunication] = <SerialCommunicationStopbits.ONE: 1>, bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunication] = <SerialCommunicationByteSize.EIGHTBITS: 8>, terminator: bytes = b'n', timeout: Union[int, float] = 3)
```

Bases: `hvl_ccb.comm.serial.SerialCommunicationConfig`

**baudrate = 19200**

Baudrate for CryLas laser is 19200 baud

**bytesize = 8**

One byte is eight bits long

**force\_value** (`fieldname, value`)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.

parity = 'N'
    CryLas laser does not use parity

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

stopbits = 1
    CryLas laser uses one stop bit

terminator = b'\n'
    The terminator is LF

timeout = 3
    use 3 seconds timeout as default

class hvl_ccb.dev.crylas.CryLasLaserShutterStatus(*args, **kwds)
Bases: aenum.Enum

Status of the CryLas laser shutter

CLOSED = 0

OPENED = 1
```

### [hvl\\_ccb.dev.ea\\_psi9000 module](#)

Device class for controlling a Elektro Automatik PSI 9000 power supply over VISA.

It is necessary that a backend for pyvisa is installed. This can be NI-Visa oder pyvisa-py (up to know, all the testing was done with NI-Visa)

```
class hvl_ccb.dev.ea_psi9000.PSI9000(com: Union[hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunication,
                                             hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig,
                                             dict], dev_config: Union[hvl_ccb.dev.ea_psi9000.PSI9000Config,
                                             dict, None] = None)
Bases: hvl_ccb.dev.visa.VisaDevice

Elektro Automatik PSI 9000 power supply.

MS_NOMINAL_CURRENT = 2040

MS_NOMINAL_VOLTAGE = 80

SHUTDOWN_CURRENT_LIMIT = 0.1

SHUTDOWN_VOLTAGE_LIMIT = 0.1

check_master_slave_config() → None
    Checks if the master / slave configuration and initializes if successful
```

**Raises `PSI9000Error`** – if master-slave configuration failed

**static `config_cls()`**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static `default_com_cls()`**

Return the default communication protocol for this device type, which is VisaCommunication.

**Returns** the VisaCommunication class

**`get_output() → bool`**

Reads the current state of the DC output of the source. Returns True, if it is enabled, false otherwise.

**Returns** the state of the DC output

**`get_system_lock() → bool`**

Get the current lock state of the system. The lock state is true, if the remote control is active and false, if not.

**Returns** the current lock state of the device

**`get_ui_lower_limits() → Tuple[float, float]`**

Get the lower voltage and current limits. A lower power limit does not exist.

**Returns** Umin in V, Imin in A

**`get_uip_upper_limits() → Tuple[float, float, float]`**

Get the upper voltage, current and power limits.

**Returns** Umax in V, Imax in A, Pmax in W

**`get_voltage_current_setpoint() → Tuple[float, float]`**

Get the voltage and current setpoint of the current source.

**Returns** Uset in V, Iset in A

**`measure_voltage_current() → Tuple[float, float]`**

Measure the DC output voltage and current

**Returns** Umeas in V, Imeas in A

**`set_lower_limits(voltage_limit: float = None, current_limit: float = None) → None`**

Set the lower limits for voltage and current. After writing the values a check is performed if the values are set correctly.

**Parameters**

- **voltage\_limit** – is the lower voltage limit in V
- **current\_limit** – is the lower current limit in A

**Raises `PSI9000Error`** – if the limits are out of range

**`set_output(target_onstate: bool) → None`**

Enables / disables the DC output.

**Parameters** **target\_onstate** – enable or disable the output power

**Raises `PSI9000Error`** – if operation was not successful

**`set_system_lock(lock: bool) → None`**

Lock / unlock the device, after locking the control is limited to this class unlocking only possible when voltage and current are below the defined limits

**Parameters** **lock** – True: locking, False: unlocking

```
set_upper_limits (voltage_limit: float = None, current_limit: float = None, power_limit: float = None) → None
```

Set the upper limits for voltage, current and power. After writing the values a check is performed if the values are set. If a parameter is left blank, the maximum configurable limit is set.

#### Parameters

- **voltage\_limit** – is the voltage limit in V
- **current\_limit** – is the current limit in A
- **power\_limit** – is the power limit in W

**Raises** `PSI9000Error` – if limits are out of range

```
set_voltage_current (volt: float, current: float) → None
```

Set voltage and current setpoints.

After setting voltage and current, a check is performed if writing was successful.

#### Parameters

- **volt** – is the setpoint voltage: 0..81.6 V (1.02 \* 0-80 V) (absolute max, can be smaller if limits are set)
- **current** – is the setpoint current: 0..2080.8 A (1.02 \* 0 - 2040 A) (absolute max, can be smaller if limits are set)

**Raises** `PSI9000Error` – if the desired setpoint is out of limits

```
start () → None
```

Start this device.

```
stop () → None
```

Stop this device. Turns off output and lock, if enabled.

```
class hvl_ccb.dev.ea_psi9000.PSI9000Config (spoll_interval: Union[int, float] = 0.5, spoll_start_delay: Union[int, float] = 2, power_limit: Union[int, float] = 43500, voltage_lower_limit: Union[int, float] = 0.0, voltage_upper_limit: Union[int, float] = 10.0, current_lower_limit: Union[int, float] = 0.0, current_upper_limit: Union[int, float] = 2040.0, wait_sec_system_lock: Union[int, float] = 0.5, wait_sec_settings_effect: Union[int, float] = 1, wait_sec_initialisation: Union[int, float] = 2)
```

Bases: `hvl_ccb.dev.visa.VisaDeviceConfig`

Elektro Automatik PSI 9000 power supply device class. The device is communicating over a VISA TCP socket.

Using this power supply, DC voltage and current can be supplied to a load with up to 2040 A and 80 V (using all four available units in parallel). The maximum power is limited by the grid, being at 43.5 kW available through the CEE63 power socket.

```
clean_values () → None
```

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

```
current_lower_limit = 0.0
```

Lower current limit in A, depending on the experimental setup.

```
current_upper_limit = 2040.0
```

Upper current limit in A, depending on the experimental setup.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**power\_limit = 43500**

Power limit in W depending on the experimental setup. With 3x63A, this is 43.5kW. Do not change this value, if you do not know what you are doing. There is no lower power limit.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**voltage\_lower\_limit = 0.0**

Lower voltage limit in V, depending on the experimental setup.

**voltage\_upper\_limit = 10.0**

Upper voltage limit in V, depending on the experimental setup.

**wait\_sec\_initialisation = 2****wait\_sec\_settings\_effect = 1****wait\_sec\_system\_lock = 0.5****exception hvl\_ccb.dev.ea\_psi9000.PSI9000Error**

Bases: Exception

Base error class regarding problems with the PSI 9000 supply.

**class hvl\_ccb.dev.ea\_psi9000.PSI9000VisaCommunication** (*configuration*)

Bases: [hvl\\_ccb.comm.visa.VisaCommunication](#)

Communication protocol used with the PSI 9000 power supply.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig(host: str, interface_type: Union[str, hvl_ccb.comm.visa.VisaCommunicationConfig] = <Interface-Type.TCPIP_SOCKET: I>, board: int = 0, port: int = 5025, timeout: int = 5000, chunk_size: int = 204800, open_timeout: int = 1000, write_termination: str = '\n', read_termination: str = '\n', visa_backend: str = '')
```

Bases: `hvl_ccb.comm.visa.VisaCommunicationConfig`

Visa communication protocol config dataclass with specification for the PSI 9000 power supply.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**interface\_type** = 1

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## hvl\_ccb.dev.heinzinger module

Device classes for Heinzinger Digital Interface I/II and Heinzinger PNC power supply.

The Heinzinger Digital Interface I/II is used for many Heinzinger power units. Manufacturer homepage: <https://www.heinzinger.com/products/accessories-and-more/digital-interfaces/>

The Heinzinger PNC series is a series of high voltage direct current power supplies. The class HeinzingerPNC is tested with two PNChp 60000-1neg and a PNChp 1500-1neg. Check the code carefully before using it with other PNC devices, especially PNC3p or PNCCap. Manufacturer homepage: <https://www.heinzinger.com/products/high-voltage/universal-high-voltage-power-supplies/>

```
class hvl_ccb.dev.heinzinger.HeinzingerConfig (default_number_of_recordings:  
    Union[int,  
          hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum]  
    = 1, number_of_decimals: int = 6,  
    wait_sec_stop_commands: Union[int,  
        float] = 0.5)
```

Bases: object

Device configuration dataclass for Heinzinger power supplies.

```
class RecordingsEnum
```

Bases: enum.IntEnum

An enumeration.

```
EIGHT = 8
```

```
FOUR = 4
```

```
ONE = 1
```

```
SIXTEEN = 16
```

```
TWO = 2
```

```
clean_values()
```

```
default_number_of_recordings = 1
```

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
is_configdataclass = True
```

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
number_of_decimals = 6
```

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
wait_sec_stop_commands = 0.5
    Time to wait after subsequent commands during stop (in seconds)

class hvl_ccb.dev.heinzinger.HeinzingerDI (com, dev_config=None)
    Bases: hvl_ccb.dev.base.SingleCommDevice, abc.ABC
    Heinzinger Digital Interface I/II device class

    Sends basic SCPI commands and reads the answer. Only the standard instruction set from the manual is implemented.

    static config_cls()
        Return the default configdataclass class.

        Returns a reference to the default configdataclass class

    static default_com_cls()
        Get the class for the default communication protocol used with this device.

        Returns the type of the standard communication protocol for this device

    get_current() → float
        Queries the set current of the Heinzinger PNC (not the measured current!).
        Raises SerialCommunicationIOError – when communication port is not opened

    get_interface_version() → str
        Queries the version number of the digital interface.
        Raises SerialCommunicationIOError – when communication port is not opened

    get_number_of_recordings() → int
        Queries the number of recordings the device is using for average value calculation.

        Returns int number of recordings
        Raises SerialCommunicationIOError – when communication port is not opened

    get_serial_number() → str
        Ask the device for its serial number and returns the answer as a string.

        Returns string containing the device serial number
        Raises SerialCommunicationIOError – when communication port is not opened

    get_voltage() → float
        Queries the set voltage of the Heinzinger PNC (not the measured voltage!).
        Raises SerialCommunicationIOError – when communication port is not opened

    measure_current() → float
        Ask the Device to measure its output current and return the measurement result.

        Returns measured current as float
        Raises SerialCommunicationIOError – when communication port is not opened

    measure_voltage() → float
        Ask the Device to measure its output voltage and return the measurement result.

        Returns measured voltage as float
        Raises SerialCommunicationIOError – when communication port is not opened

    output_off() → None
        Switch DC voltage output off.
```

**Raises `SerialCommunicationIOError`** – when communication port is not opened

**`output_on()`** → None  
Switch DC voltage output on.

**Raises `SerialCommunicationIOError`** – when communication port is not opened

**`reset_interface()`** → None  
Reset of the digital interface; only Digital Interface I: Power supply is switched to the Local-Mode (Manual operation)

**Raises `SerialCommunicationIOError`** – when communication port is not opened

**`set_current(value: Union[int, float])`** → None  
Sets the output current of the Heinzinger PNC to the given value.

**Parameters `value`** – current expressed in `self.unit_current`

**Raises `SerialCommunicationIOError`** – when communication port is not opened

**`set_number_of_recordings(value: Union[int, hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum])`**  
→ None  
Sets the number of recordings the device is using for average value calculation. The possible values are 1, 2, 4, 8 and 16.

**Raises `SerialCommunicationIOError`** – when communication port is not opened

**`set_voltage(value: Union[int, float])`** → None  
Sets the output voltage of the Heinzinger PNC to the given value.

**Parameters `value`** – voltage expressed in `self.unit_voltage`

**Raises `SerialCommunicationIOError`** – when communication port is not opened

**`start()`**  
Opens the communication protocol.

**Raises `SerialCommunicationIOError`** – when communication port cannot be opened.

**`stop()`** → None  
Stop the device. Closes also the communication protocol.

**class `hvl_ccb.dev.heinzinger.HeinzingerPNC(com, dev_config=None)`**  
Bases: `hvl_ccb.dev.heinzinger.HeinzingerDI`

Heinzinger PNC power supply device class.

The power supply is controlled over a Heinzinger Digital Interface I/II

**class `UnitCurrent(*args, **kwds)`**  
Bases: `hvl_ccb.utils.enum.AutoNumberNameEnum`

```
A = 3
UNKNOWN = 1
mA = 2
```

**class `UnitVoltage(*args, **kwds)`**  
Bases: `hvl_ccb.utils.enum.AutoNumberNameEnum`

```
UNKNOWN = 1
V = 2
kV = 3
```

```
identify_device() → None
    Identify the device nominal voltage and current based on its serial number.

    Raises SerialCommunicationIOError – when communication port is not opened

max_current
max_current_hardware
max_voltage
max_voltage_hardware

set_current(value: Union[int, float]) → None
    Sets the output current of the Heinzinger PNC to the given value.

    Parameters value – current expressed in self.unit_current

    Raises SerialCommunicationIOError – when communication port is not opened

set_voltage(value: Union[int, float]) → None
    Sets the output voltage of the Heinzinger PNC to the given value.

    Parameters value – voltage expressed in self.unit_voltage

    Raises SerialCommunicationIOError – when communication port is not opened

start() → None
    Opens the communication protocol and configures the device.

unit_current
unit_voltage

exception hvl_ccb.dev.heinzinger.HeinzingerPNCDeviceNotRecognizedException
    Bases: hvl\_ccb.dev.heinzinger.HeinzingerPNCError

    Error indicating that the serial number of the device is not recognized.

exception hvl_ccb.dev.heinzinger.HeinzingerPNCError
    Bases: Exception

    General error with the Heinzinger PNC voltage source.

exception hvl_ccb.dev.heinzinger.HeinzingerPNCMaxCurrentExceededException
    Bases: hvl\_ccb.dev.heinzinger.HeinzingerPNCError

    Error indicating that program attempted to set the current to a value exceeding ‘max_current’.

exception hvl_ccb.dev.heinzinger.HeinzingerPNCMaxVoltageExceededException
    Bases: hvl\_ccb.dev.heinzinger.HeinzingerPNCError

    Error indicating that program attempted to set the voltage to a value exceeding ‘max_voltage’.

class hvl_ccb.dev.heinzinger.HeinzingerSerialCommunication(configuration)
    Bases: hvl\_ccb.comm.serial.SerialCommunication

    Specific communication protocol implementation for Heinzinger power supplies. Already predefines device-specific protocol parameters in config.

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

read_text_nonempty(n_attempts_max: int = 40) → str
    Reads from the serial port, until a non-empty line is found, or the number of attempts is exceeded.
```

**Parameters**

- **n\_attempts\_max** – maximum number of read attempts
- **attempt\_interval\_sec** – interval between subsequent attempts in seconds

**Returns** String read from the serial port; “” if number of attempts is exceeded or serial port is not opened.

```
class hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig(port:      str,
                                                                     baudrate:
                                                                     int = 9600,
                                                                     parity:
                                                                     Union[str,
                                                                           hvl_ccb.comm.serial.SerialCommunicationParity.NONE:
                                                                           'N',
                                                                           'N'>,
                                                                           stopbits:
                                                                           Union[int,
                                                                             hvl_ccb.comm.serial.SerialCommunicationStopbits.ONE:
                                                                             1,
                                                                             1>, bytesize:
                                                                             Union[int,
                                                                               hvl_ccb.comm.serial.SerialCommunicationBytesize.EIGHTBITS:
                                                                               8,
                                                                               8>, termina-
                                                                               tor: bytes =
                                                                               b'n', timeout:
                                                                               Union[int,
                                                                                 float] = 3,
                                                                                 wait_sec_read_text_nonempty:
                                                                                 Union[int,
                                                                                   float] = 0.5)
```

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

**baudrate = 9600**

Baudrate for Heinzinger power supplies is 9600 baud

**bytesize = 8**

One byte is eight bits long

**clean\_values()**

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**

Heinzinger does not use parity

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**

Heinzinger uses one stop bit

**terminator = b'\n'**

The terminator is LF

**timeout = 3**

use 3 seconds timeout as default

**wait\_sec\_read\_text\_nonempty = 0.5**

time to wait between attempts of reading a non-empty text

## **hvl\_ccb.dev.labjack module**

Labjack Device for hvl\_ccb. Originally developed and tested for LabJack T7-PRO.

Makes use of the LabJack LJM Library Python wrapper. This wrapper needs an installation of the LJM Library for Windows, Mac OS X or Linux. Go to: <https://labjack.com/support/software/installers/ljm> and <https://labjack.com/support/software/examples/ljm/python>

**class hvl\_ccb.dev.labjack.LabJack(com, dev\_config=None)**

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

LabJack Device.

This class is tested with a LabJack T7-Pro and should also work with T4 and T7 devices communicating through the LJM Library. Other or older hardware versions and variants of LabJack devices are not supported.

**class AInRange(\*args, \*\*kwds)**

Bases: *hvl\_ccb.utils.enum.StrEnumBase*

**ONE = 1.0**

**ONE\_HUNDREDTH = 0.01**

**ONE\_TENTH = 0.1**

**TEN = 10.0**

**value**

```

class CalMicroAmpere(*args, **kwds)
    Bases: aenum.Enum

    Pre-defined microampere (uA) values for calibration current source query.

    TEN = '10uA'

    TWO_HUNDRED = '200uA'

class CjcType(*args, **kwds)
    Bases: hvl\_ccb.utils.enum.NameEnum

    CJC slope and offset

    internal = (1, 0)

    lm34 = (55.56, 255.37)

DIOChannel
    alias of hvl_ccb._dev.labjack.TSeriesDIOChannel

class DIOStatus(*args, **kwds)
    Bases: aenum.IntEnum

    State of a digital I/O channel.

    HIGH = 1

    LOW = 0

class DeviceType(*args, **kwds)
    Bases: hvl\_ccb.utils.enum.AutoNumberNameEnum

    LabJack device types.

    Can be also looked up by ambiguous Product ID (p_id) or by instance name: `python
    LabJackDeviceType(4)` is LabJackDeviceType('T4') ``

    ANY = 1

    T4 = 2

    T7 = 3

    T7_PRO = 4

    get_by_p_id = <bound method DeviceType.get_by_p_id of <aenum 'DeviceType'>>

class TemperatureUnit(*args, **kwds)
    Bases: hvl\_ccb.utils.enum.NameEnum

    Temperature unit (to be returned)

    C = 1

    F = 2

    K = 0

class ThermocoupleType(*args, **kwds)
    Bases: hvl\_ccb.utils.enum.NameEnum

    Thermocouple type; NONE means disable thermocouple mode.

    C = 30

    E = 20

    J = 21

```

**K = 22**

**NONE = 0**

**PT100 = 40**

**PT1000 = 42**

**PT500 = 41**

**R = 23**

**S = 25**

**T = 24**

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**get\_ain(\*channels) → Union[float, Sequence[float]]**

Read currently measured value (voltage, resistance, ...) from one or more of analog inputs.

**Parameters** **channels** – AIN number or numbers (0..254)

**Returns** the read value (voltage, resistance, ...) as *float* or ‘tuple of them in case multiple channels given

**get\_cal\_current\_source(name: Union[str, CalMicroAmpere]) → float**

This function will return the calibration of the chosen current source, this ist not a measurement!

The value was stored during fabrication.

**Parameters** **name** – ‘200uA’ or ‘10uA’ current source

**Returns** calibration of the chosen current source in ampere

**get\_digital\_input(address: Union[str, hvl\_ccb.\_dev.labjack.TSeriesDIOChannel]) → hvl\_ccb.dev.labjack.LabJack.DIOStatus**

Get the value of a digital input.

allowed names for T7 (Pro): FIO0 - FIO7, EIO0 - EIO 7, CIO0- CIO3, MIO0 - MIO2 :param address: name of the output -> ‘FIO0’ :return: HIGH when *address* DIO is high, and LOW when *address* DIO is low

**get\_product\_id() → int**

This function returns the product ID reported by the connected device.

Attention: returns 7 for both T7 and T7-Pro devices!

**Returns** integer product ID of the device

**get\_product\_name(force\_query\_id=False) → str**

This function will return the product name based on product ID reported by the device.

Attention: returns “T7” for both T7 and T7-Pro devices!

**Parameters** **force\_query\_id** – boolean flag to force *get\_product\_id* query to device instead of using cached device type from previous queries.

**Returns** device name string, compatible with *LabJack.DeviceType*

**get\_product\_type(force\_query\_id: bool = False) → hvl\_ccb.\_dev.labjack.DeviceType**

This function will return the device type based on reported device type and in case of unambiguity based on configuration of device’s communication protocol (e.g. for “T7” and “T7\_PRO” devices), or, if not available first matching.

**Parameters** `force_query_id` – boolean flag to force `get_product_id` query to device instead of using cached device type from previous queries.

**Returns** `DeviceType` instance

**Raises** `LabJackIdentifierDIOError` – when read Product ID is unknown

**get\_sbush\_rh** (`number: int`) → float

Read the relative humidity value from a serial SBUS sensor.

**Parameters** `number` – port number (0..22)

**Returns** relative humidity in %RH

**get\_sbush\_temp** (`number: int`) → float

Read the temperature value from a serial SBUS sensor.

**Parameters** `number` – port number (0..22)

**Returns** temperature in Kelvin

**get\_serial\_number** () → int

Returns the serial number of the connected LabJack.

**Returns** Serial number.

**read\_resistance** (`channel: int`) → float

Read resistance from specified channel.

**Parameters** `channel` – channel with resistor

**Returns** resistance value with 2 decimal places

**read\_thermocouple** (`pos_channel: int`) → float

Read the temperature of a connected thermocouple.

**Parameters** `pos_channel` – is the AIN number of the positive pin

**Returns** temperature in specified unit

**set\_ain\_differential** (`pos_channel: int, differential: bool`) → None

Sets an analog input to differential mode or not. T7-specific: For base differential channels, positive must be even channel from 0-12 and negative must be positive+1. For extended channels 16-127, see Mux80 datasheet.

**Parameters**

- `pos_channel` – is the AIN number (0..12)
- `differential` – True or False

**Raises** `LabJackError` – if parameters are unsupported

**set\_ain\_range** (`channel: int, vrange: Union[Real, AInRange]`) → None

Set the range of an analog input port.

**Parameters**

- `channel` – is the AIN number (0..254)
- `vrange` – is the voltage range to be set

**set\_ain\_resistance** (`channel: int, vrange: Union[Real, AInRange], resolution: int`) → None

Set the specified channel to resistance mode. It utilized the 200uA current source of the LabJack.

**Parameters**

- `channel` – channel that should measure the resistance

- **vrange** – voltage range of the channel
- **resolution** – resolution index of the channel T4: 0-5, T7: 0-8, T7-Pro 0-12

**set\_ain\_resolution** (*channel: int, resolution: int*) → None  
Set the resolution index of an analog input port.

#### Parameters

- **channel** – is the AIN number (0..254)
- **resolution** – is the resolution index within 0...‘get\_product\_type().ain\_max\_resolution’ range; 0 will set the resolution index to default value.

**set\_ain\_thermocouple** (*pos\_channel: int, thermocouple: Union[None, str, ThermocoupleType], cjc\_address: int = 60050, cjc\_type: Union[str, CjcType] = <CjcType.internal: (1, 0)>, vrange: Union[Real, AInRange] = <AInRange.ONE\_HUNDREDTH: '0.01'>, resolution: int = 10, unit: Union[str, TemperatureUnit] = <TemperatureUnit.K: 0>*) → None  
Set the analog input channel to thermocouple mode.

#### Parameters

- **pos\_channel** – is the analog input channel of the positive part of the differential pair
- **thermocouple** – None to disable thermocouple mode, or string specifying the thermocouple type
- **cjc\_address** – modbus register address to read the CJC temperature
- **cjc\_type** – determines cjc slope and offset, ‘internal’ or ‘lm34’
- **vrange** – measurement voltage range
- **resolution** – resolution index (T7-Pro: 0-12)
- **unit** – is the temperature unit to be returned (‘K’, ‘C’ or ‘F’)

**Raises** `LabJackError` – if parameters are unsupported

**set\_digital\_output** (*address: str, state: Union[int, DIOStatus]*) → None  
Set the value of a digital output.

#### Parameters

- **address** – name of the output -> ‘FIO0’
- **state** – state of the output -> *DIOStatus* instance or corresponding *int* value

**start()** → None  
Start the Device.

**stop()** → None  
Stop the Device.

**exception** `hvl_ccb.dev.labjack.LabJackError`  
Bases: Exception

Errors of the LabJack device.

**exception** `hvl_ccb.dev.labjack.LabJackIdentifierDIOError`  
Bases: Exception

Error indicating a wrong DIO identifier

## hvl\_ccb.dev.mbw973 module

Device class for controlling a MBW 973 SF6 Analyzer over a serial connection.

The MBW 973 is a gas analyzer designed for gas insulated switchgear and measures humidity, SF6 purity and SO2 contamination in one go. Manufacturer homepage: <https://www.mbw.ch/products/sf6-gas-analysis/973-sf6-analyzer/>

**class** hvl\_ccb.dev.mbw973.**MBW973** (*com, dev\_config=None*)

Bases: hvl\_ccb.dev.base.SingleCommDevice

MBW 973 dew point mirror device class.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**is\_done() → bool**

Poll status of the dew point mirror and return True, if all measurements are done.

**Returns** True, if all measurements are done; False otherwise.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**read(*cast\_type: Type[CT\_co] = <class 'str'>*)**

Read value from *self.com* and cast to *cast\_type*. Raises `ValueError` if read text (*str*) is not convertible to *cast\_type*, e.g. to `float` or to `int`.

**Returns** Read value of *cast\_type* type.

**read\_float() → float**

Convenience wrapper for *self.read()*, with typing hint for return value.

**Returns** Read `float` value.

**read\_int() → int**

Convenience wrapper for *self.read()*, with typing hint for return value.

**Returns** Read `int` value.

**read\_measurements() → Dict[str, float]**

Read out measurement values and return them as a dictionary.

**Returns** Dictionary with values.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**set\_measuring\_options(*humidity: bool = True, sf6\_purity: bool = False*) → None**

Send measuring options to the dew point mirror.

### Parameters

- **humidity** – Perform humidity test or not?
- **sf6\_purity** – Perform SF6 purity test or not?

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**start ()** → None  
Start this device. Opens the communication protocol and retrieves the set measurement options from the device.

**Raises** `SerialCommunicationIOError` – when communication port cannot be opened.

**start\_control ()** → None  
Start dew point control to acquire a new value set.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**stop ()** → None  
Stop the device. Closes also the communication protocol.

**write (value)** → None  
Send `value` to `self.com`.

**Parameters** `value` – Value to send, converted to `str`.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**class** `hvl_ccb.dev.mbw973.MBW973Config` (`polling_interval: Union[int, float] = 2`)  
Bases: `object`

Device configuration dataclass for MBW973.

**clean\_values ()**

**force\_value (fieldname, value)**  
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

**Parameters**

- `fieldname` – name of the field
- `value` – value to assign

**is\_configdataclass = True**

**classmethod keys ()** → Sequence[str]  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults ()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**polling\_interval = 2**

Polling period for `is_done` status queries [in seconds].

**classmethod required\_keys ()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception** `hvl_ccb.dev.mbw973.MBW973ControlRunningException`  
Bases: `hvl_ccb.dev.mbw973.MBW973Error`

Error indicating there is still a measurement running, and a new one cannot be started.

```
exception hvl_ccb.dev.mbw973.MBW973Error
```

Bases: Exception

General error with the MBW973 dew point mirror device.

```
exception hvl_ccb.dev.mbw973.MBW973PumpRunningException
```

Bases: [hvl\\_ccb.dev.mbw973.MBW973Error](#)

Error indicating the pump of the dew point mirror is still recovering gas, unable to start a new measurement.

```
class hvl_ccb.dev.mbw973.MBW973SerialCommunication(configuration)
```

Bases: [hvl\\_ccb.comm.serial.SerialCommunication](#)

Specific communication protocol implementation for the MBW973 dew point mirror. Already predefines device-specific protocol parameters in config.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig(port: str, baudrate: int = 9600, parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity] = <SerialCommunicationParity.NONE: 'N'>, stopbits: Union[int, hvl_ccb.comm.serial.SerialCommunicationStopbits = <SerialCommunicationStopbits.ONE: 1>, bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize = <SerialCommunicationByte-size.EIGHTBITS: 8>, terminator: bytes = b'r', timeout: Union[int, float] = 3))
```

Bases: [hvl\\_ccb.comm.serial.SerialCommunicationConfig](#)

```
baudrate = 9600
```

Baudrate for MBW973 is 9600 baud

```
bytesize = 8
```

One byte is eight bits long

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.

parity = 'N'
    MBW973 does not use parity

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

stopbits = 1
    MBW973 does use one stop bit

terminator = b'\r'
    The terminator is only CR

timeout = 3
    use 3 seconds timeout as default

class hvl_ccb.dev.mbw973.Poller(period: float, callback: Callable)
Bases: object

Wrapper for the threading.Timer class to periodically poll data.

start() → None
    Start the polling timer.

stop() → None
    Stop the polling timer.

timer_callback() → None
    Callback method that is called every time the timer elapses. It calls the specified user callback function and restarts the timer.
```

### **hvl\_ccb.dev.newport module**

Device class for Newport SMC100PP stepper motor controller with serial communication.

The SMC100PP is a single axis motion controller/driver for stepper motors up to 48 VDC at 1.5 A rms. Up to 31 controllers can be networked through the internal RS-485 communication link.

Manufacturer homepage: <https://www.newport.com/f/smc100-single-axis-dc-or-stepper-motion-controller>

```
class hvl_ccb.dev.newport.NewportConfigCommands(*args, **kwds)
Bases: hvl_ccb.utils.enum.NameEnum

Commands predefined by the communication protocol of the SMC100PP

AC = 'acceleration'
BA = 'backlash_compensation'
BH = 'hysteresis_compensation'
FRM = 'micro_step_per_full_step_factor'
FRS = 'motion_distance_per_full_step'
HT = 'home_search_type'
```

```

JR = 'jerk_time'
OH = 'home_search_velocity'
OT = 'home_search_timeout'
QIL = 'peak_output_current_limit'
SA = 'rs485_address'
SL = 'negative_software_limit'
SR = 'positive_software_limit'
VA = 'velocity'
VB = 'base_velocity'
ZX = 'stage_configuration'

exception hvl_ccb.dev.newport.NewportControllerError
    Bases: Exception
    Error with the Newport controller.

exception hvl_ccb.dev.newport.NewportMotorError
    Bases: Exception
    Error with the Newport motor.

class hvl_ccb.dev.newport.NewportSMC100PP (com, dev_config=None)
    Bases: hvl_ccb.dev.base.SingleCommDevice
    Device class of the Newport motor controller SMC100PP

class MotorErrors (*args, **kwds)
    Bases: aenum.Enum
    Possible motor errors reported by the motor during get_state().

    DC_VOLTAGE_TOO_LOW = 3
    FOLLOWING_ERROR = 6
    HOMING_TIMEOUT = 5
    NED_END_OF_TURN = 11
    OUTPUT_POWER_EXCEEDED = 2
    PEAK_CURRENT_LIMIT = 9
    POS_END_OF_TURN = 10
    RMS_CURRENT_LIMIT = 8
    SHORT_CIRCUIT = 7
    WRONG_ESP_STAGE = 4

class StateMessages (*args, **kwds)
    Bases: aenum.Enum
    Possible messages returned by the controller on get_state() query.

    CONFIG = '14'
    DISABLE_FROM_JOGGING = '3E'
    DISABLE_FROM_MOVING = '3D'

```

```
DISABLE_FROM_READY = '3C'
HOMING_FROM_RS232 = '1E'
HOMING_FROM_SMC = '1F'
JOGGING_FROM_DISABLE = '47'
JOGGING_FROM_READY = '46'
MOVING = '28'
NO_REF_ESP_STAGE_ERROR = '10'
NO_REF_FROM_CONFIG = '0C'
NO_REF_FROM_DISABLED = '0D'
NO_REF_FROM_HOMING = '0B'
NO_REF_FROM_JOGGING = '11'
NO_REF_FROM_MOVING = '0F'
NO_REF_FROM_READY = '0E'
NO_REF_FROM_RESET = '0A'
READY_FROM_DISABLE = '34'
READY_FROM_HOMING = '32'
READY_FROM_JOGGING = '35'
READY_FROM_MOVING = '33'
```

**States**

alias of [NewportStates](#)

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**exit\_configuration (add: int = None) → None**

Exit the CONFIGURATION state and go back to the NOT REFERENCED state. All configuration parameters are saved to the device's memory.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- [SerialCommunicationIOError](#) – if the com is closed
- [NewportSerialCommunicationError](#) – if an unexpected answer is obtained
- [NewportControllerError](#) – if the controller reports an error

**get\_acceleration (add: int = None) → Union[int, float]**

Leave the configuration state. The configuration parameters are saved to the device's memory.

**Parameters** **add** – controller address (1 to 31)

**Returns** acceleration (preset units/s^2), value between 1e-6 and 1e12

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_controller\_information**(*add: int = None*) → str

Get information on the controller name and driver version

**Parameters** *add* – controller address (1 to 31)

**Returns** controller information

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_motor\_configuration**(*add: int = None*) → Dict[str, float]

Query the motor configuration and returns it in a dictionary.

**Parameters** *add* – controller address (1 to 31)

**Returns** dictionary containing the motor's configuration

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_move\_duration**(*dist: Union[int, float], add: int = None*) → float

Estimate the time necessary to move the motor of the specified distance.

**Parameters**

- **dist** – distance to travel
- **add** – controller address (1 to 31), defaults to self.address

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_negative\_software\_limit**(*add: int = None*) → Union[int, float]

Get the negative software limit (the maximum position that the motor is allowed to travel to towards the left).

**Parameters** *add* – controller address (1 to 31)

**Returns** negative software limit (preset units), value between -1e12 and 0

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

`get_position(add: int = None) → float`

Returns the value of the current position.

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

`get_positive_software_limit(add: int = None) → Union[int, float]`

Get the positive software limit (the maximum position that the motor is allowed to travel to towards the right).

**Parameters** `add` – controller address (1 to 31)

**Returns** positive software limit (preset units), value between 0 and 1e12

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

`get_state(add: int = None) → StateMessages`

Check on the motor errors and the controller state

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error
- `NewportMotorError` – if the motor reports an error

**Returns** state message from the device (member of StateMessages)

`go_home(add: int = None) → None`

Move the motor to its home position.

**Parameters** `add` – controller address (1 to 31), defaults to self.address

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

`go_to_configuration(add: int = None) → None`

This method is executed during start(). It can also be executed after a reset(). The controller is put in CONFIG state, where configuration parameters can be changed.

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed

- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**initialize**(*add*: int = None) → None

Puts the controller from the NOT\_REF state to the READY state. Sends the motor to its “home” position.

**Parameters** *add* – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**move\_to\_absolute\_position**(*pos*: Union[int, float], *add*: int = None) → None

Move the motor to the specified position.

**Parameters**

- **pos** – target absolute position (affected by the configured offset)
- **add** – controller address (1 to 31), defaults to self.address

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**move\_to\_relative\_position**(*pos*: Union[int, float], *add*: int = None) → None

Move the motor of the specified distance.

**Parameters**

- **pos** – distance to travel (the sign gives the direction)
- **add** – controller address (1 to 31), defaults to self.address

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**reset**(*add*: int = None) → None

Resets the controller, equivalent to a power-up. This puts the controller back to NOT REFERENCED state, which is necessary for configuring the controller.

**Parameters** *add* – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_acceleration**(*acc*: Union[int, float], *add*: int = None) → None

Leave the configuration state. The configuration parameters are saved to the device’s memory.

**Parameters**

- **acc** – acceleration (preset units/s<sup>2</sup>), value between 1e-6 and 1e12
- **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_motor\_configuration** (*add*: int = None, *config*: dict = None) → None

Set the motor configuration. The motor must be in CONFIG state.

**Parameters**

- **add** – controller address (1 to 31)
- **config** – dictionary containing the motor's configuration

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_negative\_software\_limit** (*lim*: Union[int, float], *add*: int = None) → None

Set the negative software limit (the maximum position that the motor is allowed to travel to towards the left).

**Parameters**

- **lim** – negative software limit (preset units), value between -1e12 and 0
- **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_positive\_software\_limit** (*lim*: Union[int, float], *add*: int = None) → None

Set the positive software limit (the maximum position that the motor is allowed to travel to towards the right).

**Parameters**

- **lim** – positive software limit (preset units), value between 0 and 1e12
- **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**start()**

Opens the communication protocol and applies the config.

**Raises** *SerialCommunicationIOError* – when communication port cannot be opened

**stop()** → None

Stop the device. Close the communication protocol.

**stop\_motion(add: int = None)** → None

Stop a move in progress by decelerating the positioner immediately with the configured acceleration until it stops. If a controller address is provided, stops a move in progress on this controller, else stops the moves on all controllers.

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

**wait\_until\_motor\_initialized(add: int = None)** → None

Wait until the motor leaves the HOMING state (at which point it should have arrived to the home position).

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

**wait\_until\_move\_finished(add: int = None)** → None

Wait until the motor leaves the MOVING state (at which point it should have arrived to the target position).

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

```
class hvl_ccb.dev.newport.NewportSMC100PPConfig(address: int = 1, user_position_offset: Union[int, float] = 23.987, screw_scaling: Union[int, float] = 1, exit_configuration_wait_sec: Union[int, float] = 5, move_finished_extra_wait_sec: Union[int, float] = 1, acceleration: Union[int, float] = 10, backlash_compensation: Union[int, float] = 0, hysteresis_compensation: Union[int, float] = 0.015, micro_step_per_full_step_factor: int = 100, motion_distance_per_full_step: Union[int, float] = 0.01, home_search_type: Union[int, hvl_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch] = <HomeSearch.HomeSwitch: 2>, jerk_time: Union[int, float] = 0.04, home_search_velocity: Union[int, float] = 4, home_search_timeout: Union[int, float] = 27.5, home_search_polling_interval: Union[int, float] = 1, peak_output_current_limit: Union[int, float] = 0.4, rs485_address: int = 2, negative_software_limit: Union[int, float] = -23.5, positive_software_limit: Union[int, float] = 25, velocity: Union[int, float] = 4, base_velocity: Union[int, float] = 0, stage_configuration: Union[int, hvl_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig] = <EspStageConfig.EnableEspStageCheck: 3>)
```

Bases: object

Configuration dataclass for the Newport motor controller SMC100PP.

```
class EspStageConfig(*args, **kwds)
```

Bases: aenum.IntEnum

Different configurations to check or not the motor configuration upon power-up.

```
DisableEspStageCheck = 1
```

```
EnableEspStageCheck = 3
```

```
UpdateEspStageInfo = 2
```

```
class HomeSearch(*args, **kwds)
```

Bases: aenum.IntEnum

Different methods for the motor to search its home position during initialization.

```
CurrentPosition = 1
```

```
EndOfRunSwitch = 4
```

```
EndOfRunSwitch_and_Index = 3
```

```
HomeSwitch = 2
HomeSwitch_and_Index = 0
acceleration = 10
address = 1
backlash_compensation = 0
base_velocity = 0
clean_values()
exit_configuration_wait_sec = 5
force_value(fieldname, value)
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define post_force_value method with same signature as this method to do extra processing
after value has been forced on fieldname.
```

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

```
home_search_polling_interval = 1
home_search_timeout = 27.5
home_search_type = 2
home_search_velocity = 4
hysteresis_compensation = 0.015
is_configdataclass = True
jerk_time = 0.04
classmethod keys() → Sequence[str]
Returns a list of all configdataclass fields key-names.

>Returns a list of strings containing all keys.

micro_step_per_full_step_factor = 100
motion_distance_per_full_step = 0.01
motor_config
Gather the configuration parameters of the motor into a dictionary.

>Returns dict containing the configuration parameters of the motor

move_finished_extra_wait_sec = 1
negative_software_limit = -23.5
classmethod optional_defaults() → Dict[str, object]
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified
on instantiation.

>Returns a list of strings containing all optional keys.

peak_output_current_limit = 0.4
positive_software_limit = 25
```

```
post_force_value (fieldname, value)
classmethod required_keys () → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

rs485_address = 2
screw_scaling = 1
stage_configuration = 3
user_position_offset = 23.987
velocity = 4

class hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication (configuration)
Bases: hvl_ccb.comm.serial.SerialCommunication

Specific communication protocol implementation Heinzinger power supplies. Already predefines device-specific protocol parameters in config.

class ControllerErrors (*args, **kwds)
Bases: aenum.Enum

Possible controller errors with values as returned by the device in response to sent commands.

ADDR_INCORRECT = 'B'
CMD_EXEC_ERROR = 'V'
CMD_NOT_ALLOWED = 'D'
CMD_NOT_ALLOWED_CC = 'X'
CMD_NOT_ALLOWED_CONFIGURATION = 'I'
CMD_NOT_ALLOWED_DISABLE = 'J'
CMD_NOT_ALLOWED_HOMING = 'L'
CMD_NOT_ALLOWED_MOVING = 'M'
CMD_NOT_ALLOWED_NOT_REFERENCED = 'H'
CMD_NOT_ALLOWED_PP = 'W'
CMD_NOT_ALLOWED_READY = 'K'
CODE_OR_ADDR_INVALID = 'A'
COM_TIMEOUT = 'S'
DISPLACEMENT_OUT_OF_LIMIT = 'G'
EEPROM_ACCESS_ERROR = 'U'
ESP_STAGE_NAME_INVALID = 'F'
HOME_STARTED = 'E'
NO_ERROR = '@'
PARAM_MISSING_OR_INVALID = 'C'
POSITION_OUT_OF_LIMIT = 'N'
```

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query (add: int, cmd: str, param: Union[int, float, str, None] = None) → str**

Send a query to the controller, read the answer, and check for errors. The prefix add+cmd is removed from the answer.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **param** – optional parameter (int/float/str) appended to the command

**Returns** the answer from the device without the prefix

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**query\_multiple (add: int, cmd: str, prefixes: List[str]) → List[str]**

Send a query to the controller, read the answers, and check for errors. The prefixes are removed from the answers.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **prefixes** – prefixes of each line expected in the answer

**Returns** list of answers from the device without prefix

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**send\_command (add: int, cmd: str, param: Union[int, float, str, None] = None) → None**

Send a command to the controller, and check for errors.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **param** – optional parameter (int/float/str) appended to the command

**Raises**

- **SerialCommunicationIOError** – if the com is closed
- **NewportSerialCommunicationError** – if an unexpected answer is obtained
- **NewportControllerError** – if the controller reports an error

**send\_stop** (*add: int*) → None  
Send the general stop ST command to the controller, and check for errors.

**Parameters** *add* – the controller address (1 to 31)

**Returns** ControllerErrors reported by Newport Controller

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained

```
class hvl_ccb.dev.newport.NewportSMC100PPSSerialCommunicationConfig (port: str,  
baudrate:  
int =  
57600,  
parity:  
Union[str,  
hvl_ccb.comm.serial.SerialCommu  
= <Serial-  
Commu-  
nica-  
tionPar-  
ity.NONE:  
'N'>,  
stopbits:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionStop-  
bits.ONE:  
1>, byte-  
size:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionByte-  
size.EIGHTBITS:  
8>, ter-  
minator:  
bytes =  
b'\r\n',  
timeout:  
Union[int,  
float] =  
10)  
  
Bases: hvl_ccb.comm.serial.SerialCommunicationConfig  
  
baudrate = 57600  
Baudrate for Heinzinger power supplies is 9600 baud  
  
bytesize = 8  
One byte is eight bits long
```

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys () → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults () → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**

Heinzinger does not use parity

**classmethod required\_keys () → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**

Heinzinger uses one stop bit

**terminator = b'\r\n'**

The terminator is CR/LF

**timeout = 10**

use 10 seconds timeout as default

**exception hvl\_ccb.dev.newport.NewportSerialCommunicationError**

Bases: Exception

Communication error with the Newport controller.

**class hvl\_ccb.dev.newport.NewportStates (\*args, \*\*kwds)**

Bases: [hvl\\_ccb.utils.enum.AutoNumberNameEnum](#)

States of the Newport controller. Certain commands are allowed only in certain states.

**CONFIG = 3****DISABLE = 6****HOMING = 2****JOGGING = 7****MOVING = 5****NO\_REF = 1****READY = 4**

## `hvl_ccb.dev.pfeiffer_tpg module`

Device class for Pfeiffer TPG controllers.

The Pfeiffer TPG control units are used to control Pfeiffer Compact Gauges. Models: TPG 251 A, TPG 252 A, TPG 256A, TPG 261, TPG 262, TPG 361, TPG 362 and TPG 366.

Manufacturer homepage: <https://www.pfeiffer-vacuum.com/en/products/measurement-analysis/> measurement/activeline/controllers/

```
class hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG(com, dev_config=None)
Bases: hvl_ccb.dev.base.SingleCommDevice

Pfeiffer TPG control unit device class

class PressureUnits(*args, **kwds)
Bases: hvl_ccb.utils.enum.NameEnum

Enum of available pressure units for the digital display. "0" corresponds either to bar or to mbar depending on the TPG model. In case of doubt, the unit is visible on the digital display.

Micron = 3
Pascal = 2
Torr = 1
Volt = 5
bar = 0
hPascal = 4
mbar = 0

class SensorStatus
Bases: enum.IntEnum

An enumeration.

Identification_error = 6
No_sensor = 5
Ok = 0
Overrange = 2
Sensor_error = 3
Sensor_off = 4
Underrange = 1

class SensorTypes
Bases: enum.Enum

An enumeration.

CMR = 4
IKR = 2
IKR11 = 2
IKR9 = 2
IMR = 5
```

```

None = 7
PBR = 6
PKR = 3
TPR = 1
noSENSOR = 7
noSen = 7

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

get_full_scale_mbar() → List[Union[int, float]]
    Get the full scale range of the attached sensors

    Returns full scale range values in mbar, like [0.01, 1, 0.1, 1000, 50000, 10]

    Raises
        • SerialCommunicationIOError – when communication port is not opened
        • PfeifferTPGError – if command fails

get_full_scale_unitless() → List[int]
    Get the full scale range of the attached sensors. See lookup table between command and corresponding pressure in the device user manual.

    Returns list of full scale range values, like [0, 1, 3, 3, 2, 0]

    Raises
        • SerialCommunicationIOError – when communication port is not opened
        • PfeifferTPGError – if command fails

identify_sensors() → None
    Send identification request TID to sensors on all channels.

    Raises
        • SerialCommunicationIOError – when communication port is not opened
        • PfeifferTPGError – if command fails

measure(channel: int) → Tuple[str, float]
    Get the status and measurement of one sensor

    Parameters channel – int channel on which the sensor is connected, with
    1 <= channel <= number_of_sensors :return: measured value as float if measurement successful, sensor status as string if not :raises SerialCommunicationIOError: when communication port is not opened :raises PfeifferTPGError: if command fails

measure_all() → List[Tuple[str, float]]
    Get the status and measurement of all sensors (this command is not available on all models)

    Returns list of measured values as float if measurements successful,

```

and or sensor status as strings if not :raises SerialCommunicationIOError: when communication port is not opened :raises PfeifferTPGError: if command fails

**number\_of\_sensors**

**set\_display\_unit** (*unit*: Union[str, hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits]) → None  
Set the unit in which the measurements are shown on the display.

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**set\_full\_scale\_mbar** (*fsr*: List[Union[int, float]]) → None

Set the full scale range of the attached sensors (in unit mbar)

**Parameters** **fsr** – full scale range values in mbar, for example [0.01, 1000]

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**set\_full\_scale\_unitless** (*fsr*: List[int]) → None

Set the full scale range of the attached sensors. See lookup table between command and corresponding pressure in the device user manual.

**Parameters** **fsr** – list of full scale range values, like [0, 1, 3, 3, 2, 0]

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if command fails

**start()** → None

Start this device. Opens the communication protocol, and identify the sensors.

**Raises** **SerialCommunicationIOError** – when communication port cannot be opened

**stop()** → None

Stop the device. Closes also the communication protocol.

**unit**

The pressure unit of readings is always mbar, regardless of the display unit.

```
class hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig(model: Union[str, hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig.Model] = <Model.TPG25xA: {1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6, 50000: 7, 0.1: 8}>)
```

Bases: object

Device configuration dataclass for Pfeiffer TPG controllers.

**class Model(\*args, \*\*kwargs)**

Bases: [hvl\\_ccb.utils.enum.NameEnum](#)

**TPG25xA** = {0.1: 8, 1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6,

**TPGx6x** = {0.01: 0, 0.1: 1, 1: 2, 10: 3, 100: 4, 1000: 5, 2000: 6, 5000: 7,

**is\_valid\_scale\_range\_reversed\_str(v: str)** → bool

Check if given string represents a valid reversed scale range of a model.

**Parameters** **v** – Reversed scale range string.

**Returns** *True* if valid, *False* otherwise.

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**model = {0.1: 8, 1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6, 50000: 7}**

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGError**

Bases: Exception

Error with the Pfeiffer TPG Controller.

**class hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunication(configuration)**

Bases: *hvl\_ccb.comm.serial.SerialCommunication*

Specific communication protocol implementation for Pfeiffer TPG controllers. Already predefines device-specific protocol parameters in config.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query(cmd: str) → str**

Send a query, then read and returns the first line from the com port.

**Parameters cmd** – query message to send to the device

**Returns** first line read on the com

#### Raises

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if the device does not acknowledge the command or if

the answer from the device is empty

**send\_command**(*cmd: str*) → None

Send a command to the device and check for acknowledgement.

**Parameters** **cmd** – command to send to the device

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if the answer from the device differs from the expected acknowledgement character ‘chr(6)’.

```
class hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig(port:  
    str, baudrate:  
    int = 9600,  
    parity:  
    Union[str,  
          hvl_ccb.comm.serial.SerialComm = <SerialCom-  
          munica-  
          tionPar-  
          ity.NONE:  
          'N'>,  
          stopbits:  
          Union[int,  
                hvl_ccb.comm.serial.SerialComm = <SerialCom-  
                munica-  
                tionStop-  
                bits.ONE:  
                1>,  
                bytesize:  
                Union[int,  
                      hvl_ccb.comm.serial.SerialComm = <SerialCom-  
                      munica-  
                      tionByte-  
                      size.EIGHTBITS:  
                      8>, terminator:  
                      bytes = b'\r\n',  
                      timeout:  
                      Union[int,  
                            float] = 3)
```

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

**baudrate = 9600**

Baudrate for Pfeiffer TPG controllers is 9600 baud

**bytesize = 8**

One byte is eight bits long

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**

Pfeiffer TPG controllers do not use parity

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**

Pfeiffer TPG controllers use one stop bit

**terminator = b'\r\n'**

The terminator is <CR><LF>

**timeout = 3**

use 3 seconds timeout as default

## **hvl\_ccb.dev.rs\_rto1024 module**

Python module for the Rhode & Schwarz RTO 1024 oscilloscope. The communication to the device is through VISA, type TCPIP / INSTR.

```
class hvl_ccb.dev.rs_rto1024.RTO1024 (com: Union[hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunication,  
                                              hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig,  
                                              dict], dev_config: Union[hvl_ccb.dev.rs_rto1024.RTO1024Config,  
                                              dict])
```

Bases: *hvl\_ccb.dev.visa.VisaDevice*

Device class for the Rhode & Schwarz RTO 1024 oscilloscope.

**class TriggerModes** (\**args*, \*\**kwds*)

Bases: *hvl\_ccb.utils.enum.AutoNumberNameEnum*

Enumeration for the three available trigger modes.

**AUTO = 1**

**FREERUN = 3**

```
NORMAL = 2

names = <bound method RTO1024.TriggerModes.names of <aenum 'TriggerModes'>>

backup_waveform(filename: str) → None
    Backup a waveform file from the standard directory specified in the device configuration to the standard
    backup destination specified in the device configuration. The filename has to be specified without .bin or
    path.

    Parameters filename – The waveform filename without extension and path

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls()
    Return the default communication protocol for this device type, which is VisaCommunication.

    Returns the VisaCommunication class

file_copy(source: str, destination: str) → None
    Copy a file from one destination to another on the oscilloscope drive. If the destination file already exists,
    it is overwritten without notice.

    Parameters
        • source – absolute path to the source file on the DSO filesystem
        • destination – absolute path to the destination file on the DSO filesystem

    Raises RTO1024Error – if the operation did not complete

get_timestamps() → List[float]
    Gets the timestamps of all recorded frames in the history and returns them as a list of floats.

    Returns list of timestamps in [s]

    Raises RTO1024Error – if the timestamps are invalid

list_directory(path: str) → List[Tuple[str, str, int]]
    List the contents of a given directory on the oscilloscope filesystem.

    Parameters path – is the path to a folder

    Returns a list of filenames in the given folder

load_configuration(filename: str) → None
    Load current settings from a configuration file. The filename has to be specified without base directory
    and '.dfl' extension.

    Information from the manual ReCall calls up the instrument settings from an intermediate memory
    identified by the specified number. The instrument settings can be stored to this memory using the com-
    mand *SAV with the associated number. It also activates the instrument settings which are stored in a file
    and loaded using MMEMory:LOAD:STATE .

    Parameters filename – is the name of the settings file without path and extension

local_display(state: bool) → None
    Enable or disable local display of the scope.

    Parameters state – is the desired local display state

prepare_ultra_segmentation() → None
    Make ready for a new acquisition in ultra segmentation mode. This function does one acquisition without
    ultra segmentation to clear the history and prepare for a new measurement.
```

**run\_continuous\_acquisition()** → None

Start acquiring continuously.

**run\_single\_acquisition()** → None

Start a single or Nx acquisition.

**save\_configuration(filename: str)** → None

Save the current oscilloscope settings to a file. The filename has to be specified without path and ‘.dfl’ extension, the file will be saved to the configured settings directory.

**Information from the manual** *SAVe* stores the current instrument settings under the specified number in an intermediate memory. The settings can be recalled using the command *\*RCL* with the associated number. To transfer the stored instrument settings to a file, use *MMEMemory:STORe:STATe*.

**Parameters** **filename** – is the name of the settings file without path and extension

**save\_waveform\_history(filename: str, channel: int, waveform: int = 1)** → None

Save the history of one channel and one waveform to a .bin file. This function is used after an acquisition using sequence trigger mode (with or without ultra segmentation) was performed.

**Parameters**

- **filename** – is the name (without extension) of the file
- **channel** – is the channel number
- **waveform** – is the waveform number (typically 1)

**Raises** *RTO1024Error* – if storing waveform times out

**set\_acquire\_length(timerange: float)** → None

Defines the time of one acquisition, that is the time across the 10 divisions of the diagram.

- Range: 250E-12 … 500 [s]
- Increment: 1E-12 [s]
- \*RST = 0.5 [s]

**Parameters** **timerange** – is the time for one acquisition. Range: 250e-12 … 500 [s]

**set\_channel\_position(channel: int, position: float)** → None

Sets the vertical position of the indicated channel as a graphical value.

- Range: -5.0 … 5.0 [div]
- Increment: 0.02
- \*RST = 0

**Parameters**

- **channel** – is the channel number (1..4)
- **position** – is the position. Positive values move the waveform up, negative values move it down.

**set\_channel\_range(channel: int, v\_range: float)** → None

Sets the voltage range across the 10 vertical divisions of the diagram. Use the command alternatively instead of set\_channel\_scale.

- Range for range: Depends on attenuation factors and coupling. With 1:1 probe and external attenuations and 50 Ω input coupling, the range is 10 mV to 10 V. For 1 MΩ input coupling, it is 10 mV to

100 V. If the probe and/or external attenuation is changed, multiply the range values by the attenuation factors.

- Increment: 0.01
- \*RST = 0.5

#### Parameters

- **channel1** – is the channel number (1..4)
- **v\_range** – is the vertical range [V]

**set\_channel\_scale** (*channel: int, scale: float*) → None

Sets the vertical scale for the indicated channel. The scale value is given in volts per division.

- Range for scale: depends on attenuation factor and coupling. With 1:1 probe and external attenuations and  $50\ \Omega$  input coupling, the vertical scale (input sensitivity) is 1 mV/div to 1 V/div. For  $1\ M\Omega$  input coupling, it is 1 mV/div to 10 V/div. If the probe and/or external attenuation is changed, multiply the values by the attenuation factors to get the actual scale range.

- Increment: 1e-3
- \*RST = 0.05

See also: [set\\_channel\\_range](#)

#### Parameters

- **channel1** – is the channel number (1..4)
- **scale** – is the vertical scaling [V/div]

**set\_channel\_state** (*channel: int, state: bool*) → None

Switches the channel signal on or off.

#### Parameters

- **channel** – is the input channel (1..4)
- **state** – is True for on, False for off

**set\_reference\_point** (*percentage: int*) → None

Sets the reference point of the time scale in % of the display. If the “Trigger offset” is zero, the trigger point matches the reference point. ReferencePoint = zero pint of the time scale

- Range: 0 … 100 [%]
- Increment: 1 [%]
- \*RST = 50 [%]

**Parameters percentage** – is the reference in %

**set\_repetitions** (*number: int*) → None

Set the number of acquired waveforms with RUN Nx SINGLE. Also defines the number of waveforms used to calculate the average waveform.

- Range: 1 … 16777215
- Increment: 10
- \*RST = 1

**Parameters number** – is the number of waveforms to acquire

---

**set\_trigger\_level** (*channel: int, level: float, event\_type: int = 1*) → None

Sets the trigger level for the specified event and source.

- Range: -10 to 10 V
- Increment: 1e-3 V
- \*RST = 0 V

#### Parameters

- **channel** – indicates the trigger source.
  - 1..4 = channel 1 to 4, available for all event types 1..3
  - 5 = external trigger input on the rear panel for analog signals, available for A-event type = 1
  - 6..9 = not available
- **level** – is the voltage for the trigger level in [V].
- **event\_type** – is the event type. 1: A-Event, 2: B-Event, 3: R-Event

**set\_trigger\_mode** (*mode: Union[str, hvl\_ccb.dev.rs\_rto1024.RTO1024TriggerModes]*) → None

Sets the trigger mode which determines the behavior of the instrument if no trigger occurs.

**Parameters** **mode** – is either auto, normal, or freerun.

**Raises** *RTO1024Error* – if an invalid triggermode is selected

**set\_trigger\_source** (*channel: int, event\_type: int = 1*) → None

Set the trigger (Event A) source channel.

#### Parameters

- **channel** – is the channel number (1..4)
- **event\_type** – is the event type. 1: A-Event, 2: B-Event, 3: R-Event

**start()** → None

Start the RTO1024 oscilloscope and bring it into a defined state and remote mode.

**stop()** → None

Stop the RTO1024 oscilloscope, reset events and close communication. Brings back the device to a state where local operation is possible.

**stop\_acquisition()** → None

Stop any acquisition.

```
class hvl_ccb.dev.rs_rto1024.RTO1024Config(waveforms_path: str, settings_path: str, backup_path: str, spoll_interval: Union[int, float] = 0.5, spoll_start_delay: Union[int, float] = 2, command_timeout_seconds: Union[int, float] = 60, wait_sec_short_pause: Union[int, float] = 0.1, wait_sec_enable_history: Union[int, float] = 1, wait_sec_post_acquisition_start: Union[int, float] = 2)
```

Bases: *hvl\_ccb.dev.visa.VisaDeviceConfig*, *hvl\_ccb.dev.rs\_rto1024.\_RTO1024ConfigDefaultsBase*, *hvl\_ccb.dev.rs\_rto1024.\_RTO1024ConfigBase*

Configdataclass for the RTO1024 device.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception hvl\_ccb.dev.rs\_rto1024.RTO1024Error**

Bases: Exception

**class hvl\_ccb.dev.rs\_rto1024.RTO1024VisaCommunication** (*configuration*)

Bases: [hvl\\_ccb.comm.visa.VisaCommunication](#)

Specialization of VisaCommunication for the RTO1024 oscilloscope

**static config\_cls** ()

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class hvl\_ccb.dev.rs\_rto1024.RTO1024VisaCommunicationConfig** (*host: str, interface\_type: Union[str, hvl\_ccb.comm.visa.VisaCommunicationConfig] = <InterfaceType.TCPPIP\_INSTR: 2>, board: int = 0, port: int = 5025, timeout: int = 5000, chunk\_size: int = 204800, open\_timeout: int = 1000, write\_termination: str = '\n', read\_termination: str = '\n', visa\_backend: str = ''*)

Bases: [hvl\\_ccb.comm.visa.VisaCommunicationConfig](#)

Configuration dataclass for VisaCommunication with specifications for the RTO1024 device class.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**interface\_type = 2**

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## hvl\_ccb.dev.se\_il2t module

Device class for controlling a Schneider Electric ILS2T stepper drive over modbus TCP.

**class hvl\_ccb.dev.se\_il2t.ILS2T** (*com*, *dev\_config=None*)

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Schneider Electric ILS2T stepper drive class.

**ACTION\_JOG\_VALUE = 0**

The single action value for *ILS2T.Mode.JOG*

**class ActionsPtp**

Bases: *enum.IntEnum*

Allowed actions in the point to point mode (*ILS2T.Mode.PTP*).

**ABSOLUTE\_POSITION = 0**

**RELATIVE\_POSITION\_MOTOR = 2**

**RELATIVE\_POSITION\_TARGET = 1**

**DEFAULT\_IO\_SCANNING\_CONTROL\_VALUES = {'action': 2, 'continue\_after\_stop\_cu': 0, 'dis**

Default IO Scanning control mode values

**class Mode**

Bases: *enum.IntEnum*

ILS2T device modes

**JOG = 1**

```
PTP = 3

class Ref16Jog
    Bases: enum.Flag
    Allowed values for ILS2T ref_16 register (the shown values are the integer representation of the bits), all
    in Jog mode = 1

        FAST = 4
        NEG = 2
        NEG_FAST = 6
        NONE = 0
        POS = 1
        POS_FAST = 5

class RegAddr
    Bases: enum.IntEnum
    ILS2T Modbus Register Adresses
        ACCESS_ENABLE = 282
        FLT_INFO = 15362
        FLT_MEM_DEL = 15112
        FLT_MEM_RESET = 15114
        IO_SCANNING = 6922
        JOGN_FAST = 10506
        JOGN_SLOW = 10504
        POSITION = 7706
        RAMP_ACC = 1556
        RAMP_DECEL = 1558
        RAMP_N_MAX = 1554
        RAMP_TYPE = 1574
        SCALE = 1550
        TEMP = 7200
        VOLT = 7198

class RegDatatype(*args, **kwds)
    Bases: aenum.Enum
    Modbus Register Datatypes
    From the manual of the drive:
```

datatype	byte	min	max
INT8	1 Byte	-128	127
UINT8	1 Byte	0	255
INT16	2 Byte	-32_768	32_767
UINT16	2 Byte	0	65_535
INT32	4 Byte	-2_147_483_648	2_147_483_647
UINT32	4 Byte	0	4_294_967_295
BITS	just 32bits	N/A	N/A

**INT32 = (-2147483648, 2147483647)**

**is\_in\_range**(*value: int*) → bool

**class State**

Bases: enum.IntEnum

State machine status values

**ON = 6**

**QUICKSTOP = 7**

**READY = 4**

**absolute\_position**(*position: int*) → None

Turn the motor until it reaches the absolute position. This function does not enable or disable the motor automatically.

**Parameters** **position** – absolute position of motor in user defined steps.

**absolute\_position\_and\_wait**(*position: int*) → None

Enable motor, perform absolute position and wait until done, disable.

**Parameters** **position** – absolute position of motor in user defined steps.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**disable()** → None

Disable the driver of the stepper motor and enable the brake.

**do\_ioscanning\_write**(\*\*kwargs) → None

Perform a write operation using IO Scanning mode.

**Parameters** **kwargs** – Keyword-argument list with options to send, remaining are taken from the defaults.

**enable()** → None

Enable the driver of the stepper motor and disable the brake.

**get\_dc\_volt()** → float

Read the DC supply voltage of the motor.

**Returns** DC input voltage.

**get\_error\_code()** → Dict[int, Dict[str, Any]]

Read all messages in fault memory. Will read the full error message and return the decoded values. At the end the fault memory of the motor will be deleted. In addition, `reset_error` is called to re-enable the motor for operation.

**Returns** Dictionary with all information

**get\_position()** → int

Read the position of the drive and store into status.

**Returns** Position step value

**get\_status()** → Dict[str, int]

Perform an IO Scanning read and return the status of the motor.

**Returns** dict with status information.

**get\_temperature()** → int

Read the temperature of the motor.

**Returns** Temperature in degrees Celsius.

**jog\_run** (*direction: bool = True, fast: bool = False*) → None

Slowly turn the motor in positive direction.

**jog\_stop()** → None

Stop turning the motor in Jog mode.

**quickstop()** → None

Stops the motor with high deceleration rate and falls into error state. Reset with `reset_error` to recover into normal state.

**relative\_step** (*steps: int*) → None

Turn the motor the relative amount of steps. This function does not enable or disable the motor automatically. positive numbers -> CW negative numbers -> CCW

**Parameters** **steps** – Number of steps to turn the motor.

**relative\_step\_and\_wait** (*steps: int*) → None

Enable motor, perform relative steps and wait until done, disable.

**Parameters** **steps** – Number of steps.

**reset\_error()** → None

Resets the motor into normal state after quick stop or another error occurred.

**set\_jog\_speed** (*slow: int = 60, fast: int = 180*) → None

Set the speed for jog mode. Default values correspond to startup values of the motor.

**Parameters**

- **slow** – RPM for slow jog mode.
- **fast** – RPM for fast jog mode.

**set\_max\_acceleration** (*rpm\_minute: int*) → None

Set the maximum acceleration of the motor.

**Parameters** **rpm\_minute** – revolution per minute per minute

**set\_max\_deceleration** (*rpm\_minute: int*) → None

Set the maximum deceleration of the motor.

**Parameters** **rpm\_minute** – revolution per minute per minute

**set\_max\_rpm**(*rpm: int*) → None

Set the maximum RPM.

**Parameters** **rpm** – revolution per minute ( 0 < rpm <= RPM\_MAX)

**Raises** *ILS2TException* – if RPM is out of range

**set\_ramp\_type**(*ramp\_type: int = -1*) → None

**Set the ramp type. There are two options available:** 0: linear ramp -1: motor optimized ramp

**Parameters** **ramp\_type** – 0: linear ramp | -1: motor optimized ramp

**start**() → None

Start this device.

**stop**() → None

Stop this device. Disables the motor (applies brake), disables access and closes the communication protocol.

**user\_steps**(*steps: int = 16384, revolutions: int = 1*) → None

Define steps per revolution. Default is 16384 steps per revolution. Maximum precision is 32768 steps per revolution.

#### Parameters

- **steps** – number of steps in *revolutions*.
- **revolutions** – number of revolutions corresponding to *steps*.

```
class hvl_ccb.dev.se_il2t.ILS2TConfig(rpm_max_init: numbers.Integral = 1500,
                                         wait_sec_post_enable: Union[float, int] = 1,
                                         wait_sec_post_relative_step: Union[float, int] = 2,
                                         wait_sec_post_absolute_position: Union[float, int]
                                         = 2)
```

Bases: object

Configuration for the ILS2T stepper motor device.

**clean\_values**()

**force\_value**(*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass** = True

**classmethod keys**() → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults**() → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
rpm_max_init = 1500
```

initial maximum RPM for the motor, can be set up to 3000 RPM. The user is allowed to set a new max RPM at runtime using `ILS2T.set_max_rpm()`, but the value must never exceed this configuration setting.

```
wait_sec_post_absolute_position = 2
```

```
wait_sec_post_enable = 1
```

```
wait_sec_post_relative_step = 2
```

```
exception hvl_ccb.dev.se_ilts2t.ILS2TException
```

Bases: Exception

Exception to indicate problems with the SE ILS2T stepper motor.

```
class hvl_ccb.dev.se_ilts2t.ILS2TModbusTcpCommunication(configuration)
```

Bases: `hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication`

Specific implementation of Modbus/TCP for the Schneider Electric ILS2T stepper motor.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.se_ilts2t.ILS2TModbusTcpCommunicationConfig(host: str, unit: int = 255, port: int = 502)
```

Bases: `hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig`

Configuration dataclass for Modbus/TCP communication specific for the Schneider Electric ILS2T stepper motor.

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**unit = 255**

The unit has to be 255 such that IO scanning mode works.

**exception hvl\_ccb.dev.se\_ils2t.IoScanningModeValueError**

Bases: *hvl\_ccb.dev.se\_ils2t.ILS2TException*

Exception to indicate that the selected IO scanning mode is invalid.

**exception hvl\_ccb.dev.se\_ils2t.ScalingFactorValueError**

Bases: *hvl\_ccb.dev.se\_ils2t.ILS2TException*

Exception to indicate that a scaling factor value is invalid.

## **hvl\_ccb.dev.visa module**

**class hvl\_ccb.dev.visa.VisaDevice**(*com:* Union[hvl\_ccb.comm.visa.VisaCommunication, hvl\_ccb.comm.visa.VisaCommunicationConfig, dict], *dev\_config:* Union[hvl\_ccb.dev.visa.VisaDeviceConfig, dict, None] = None)

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Device communicating over the VISA protocol using VisaCommunication.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls() → Type[hvl\_ccb.comm.visa.VisaCommunication]**

Return the default communication protocol for this device type, which is VisaCommunication.

**Returns** the VisaCommunication class

**get\_error\_queue() → str**

Read out error queue and logs the error.

**Returns** Error string

**get\_identification() → str**

Queries “\*IDN?” and returns the identification string of the connected device.

**Returns** the identification string of the connected device

**reset() → None**

Send “\*RST” and “\*CLS” to the device. Typically sets a defined state.

**spoll\_handler()**

Reads the status byte and decodes it. The status byte STB is defined in IEEE 488.2. It provides a rough overview of the instrument status.

**Returns**

**start() → None**

Start the VisaDevice. Sets up the status poller and starts it.

**Returns**

**stop() → None**

Stop the VisaDevice. Stops the polling thread and closes the communication protocol.

**Returns**

**wait\_operation\_complete**(*timeout: Optional[float] = None*) → bool

Waits for a operation complete event. Returns after timeout [s] has expired or the operation complete event has been caught.

**Parameters** **timeout** – Time in seconds to wait for the event; *None* for no timeout.

**Returns** True, if OPC event is caught, False if timeout expired

**class** hvl\_ccb.dev.visa.**VisaDeviceConfig**(*spoll\_interval: Union[int, float] = 0.5,*  
*spoll\_start\_delay: Union[int, float] = 2*)  
Bases: hvl\_ccb.dev.visa.\_VisaDeviceConfigDefaultsBase, hvl\_ccb.dev.visa.\_VisaDeviceConfigBase

Configdataclass for a VISA device.

**force\_value**(*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys**() → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults**() → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys**() → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**class** hvl\_ccb.dev.visa.**VisaStatusPoller**(*target: Callable, interval: float = 0.5, start\_delay:*  
*float = 5*)

Bases: threading.Thread

Thread to periodically poll the status byte of a VISA device.

**run**() → None

Threaded method.

**stop**() → None

Gracefully stop the poller.

## Module contents

Devices subpackage.

## hvl\_ccb.utils package

### Submodules

#### hvl\_ccb.utils.enum module

**class** hvl\_ccb.utils.enum.AutoNumberNameEnum (\*args, \*\*kwds)

Bases: hvl\_ccb.utils.enum.StrEnumBase, aenum.AutoNumberEnum

Auto-numbered enum with names used as string representation, and with lookup and equality based on this representation.

**class** hvl\_ccb.utils.enum.NameEnum (\*args, \*\*kwds)

Bases: hvl\_ccb.utils.enum.StrEnumBase

Enum with names used as string representation, and with lookup and equality based on this representation.

**class** hvl\_ccb.utils.enum.StrEnumBase (\*args, \*\*kwds)

Bases: aenum.Enum

String representation-based equality and lookup.

**class** hvl\_ccb.utils.enum.ValueEnum (\*args, \*\*kwds)

Bases: hvl\_ccb.utils.enum.StrEnumBase

Enum with string representation of values used as string representation, and with lookup and equality based on this representation.

Attention: to avoid errors, best use together with *unique* enum decorator.

#### hvl\_ccb.utils.typing module

Additional Python typing module utilities

hvl\_ccb.utils.typing.check\_generic\_type (value, type\_, name='instance')

Check if value is of a generic type type\_. Raises *TypeError* if it's not.

##### Parameters

- **name** – name to report in case of an error
- **value** – value to check
- **type** – generic type to check against

hvl\_ccb.utils.typing.is\_generic (type\_)

Check if class is a user-defined generic type, for example *Union[int, float]* but not *List*.

##### Parameters type – type to check

hvl\_ccb.utils.typing.is\_type\_hint (type\_)

Check if class is a generic type, for example *Union* or *List[int]*

##### Parameters type – type to check

### Module contents

#### 4.1.2 Submodules

### hvl\_ccb.configuration module

Facilities providing classes for handling configuration for communication protocols and devices.

**class** hvl\_ccb.configuration.ConfigurationMixin (*configuration*)

Bases: abc.ABC

Mixin providing configuration to a class.

**config**

ConfigDataclass property.

**Returns** the configuration

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**configuration\_save\_json** (*path: str*) → None

Save current configuration as JSON file.

**Parameters** **path** – path to the JSON file.

**classmethod from\_json** (*filename: str*)

Instantiate communication protocol using configuration from a JSON file.

**Parameters** **filename** – Path and filename to the JSON configuration

hvl\_ccb.configuration.configdataclass (*direct\_decoration=None, frozen=True*)

Decorator to make a class a configdataclass. Types in these dataclasses are enforced. Implement a function clean\_values(self) to do additional checking on value ranges etc.

It is possible to inherit from a configdataclass and re-decorate it with @configdataclass. In a subclass, default values can be added to existing fields. Note: adding additional non-default fields is prone to errors, since the order has to be respected through the whole chain (first non-default fields, only then default-fields).

**Parameters** **frozen** – defaults to True. False allows to later change configuration values. Attention: if configdataclass is not frozen and a value is changed, typing is not enforced anymore!

### hvl\_ccb.experiment\_manager module

Main module containing the top level ExperimentManager class. Inherit from this class to implement your own experiment functionality in another project and it will help you start, stop and manage your devices.

**exception** hvl\_ccb.experiment\_manager.ExperimentError

Bases: Exception

Exception to indicate that the current status of the experiment manager is on ERROR and thus no operations can be made until reset.

**class** hvl\_ccb.experiment\_manager.ExperimentManager (*devices: Dict[str, hvl\_ccb.dev.base.Device]*)

Bases: hvl\_ccb.dev.base.DeviceSequenceMixin

Experiment Manager can start and stop communication protocols and devices. It provides methods to queue commands to devices and collect results.

**add\_device** (*name: str, device: hvl\_ccb.dev.base.Device*) → None

Add a new device to the manager. If the experiment is running, automatically start the device. If a device with this name already exists, raise an exception.

**Parameters**

- **name** – is the name of the device.
- **device** – is the instantiated Device object.

**Raises** *DeviceExistingException* –

**finish()** → None

Stop experimental setup, stop all devices.

**is\_error()** → bool

Returns true, if the status of the experiment manager is *error*.

**Returns** True if on error, false otherwise

**is\_finished()** → bool

Returns true, if the status of the experiment manager is *finished*.

**Returns** True if finished, false otherwise

**is\_running()** → bool

Returns true, if the status of the experiment manager is *running*.

**Returns** True if running, false otherwise

**run()** → None

Start experimental setup, start all devices.

**start()** → None

Alias for ExperimentManager.run()

**status**

Get experiment status.

**Returns** experiment status enum code.

**stop()** → None

Alias for ExperimentManager.finish()

**class** hvl\_ccb.experiment\_manager.**ExperimentStatus**

Bases: enum.Enum

Enumeration for the experiment status

**ERROR** = 5

**FINISHED** = 4

**FINISHING** = 3

**INITIALIZED** = 0

**INITIALIZING** = -1

**RUNNING** = 2

**STARTING** = 1

### 4.1.3 Module contents

Top-level package for HVL Common Code Base.



# CHAPTER 5

---

## Contributing

---

Contributions are welcome, and they are greatly appreciated! Every little bit helps, and credit will always be given.

You can contribute in many ways:

### 5.1 Types of Contributions

#### 5.1.1 Report Bugs

Report bugs at [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/issues](https://gitlab.com/ethz_hvl/hvl_ccb/issues).

If you are reporting a bug, please include:

- Your operating system name and version.
- Any details about your local setup that might be helpful in troubleshooting.
- Detailed steps to reproduce the bug.

#### 5.1.2 Fix Bugs

Look through the GitLab issues for bugs. Anything tagged with “bug” and “help wanted” is open to whoever wants to implement it.

#### 5.1.3 Implement Features

Look through the GitLab issues for features. Anything tagged with “enhancement” and “help wanted” is open to whoever wants to implement it.

### 5.1.4 Write Documentation

HVL Common Code Base could always use more documentation, whether as part of the official HVL Common Code Base docs, in docstrings, or even on the web in blog posts, articles, and such.

### 5.1.5 Submit Feedback

The best way to send feedback is to file an issue at [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/issues](https://gitlab.com/ethz_hvl/hvl_ccb/issues).

If you are proposing a feature:

- Explain in detail how it would work.
- Keep the scope as narrow as possible, to make it easier to implement.
- Remember that this is a volunteer-driven project, and that contributions are welcome :)

## 5.2 Get Started!

Ready to contribute? Here's how to set up *hvl\_ccb* for local development.

1. Clone *hvl\_ccb* repo from GitLab.

```
$ git clone git@gitlab.com:ethz_hvl/hvl_ccb.git
```

2. Install your local copy into a virtualenv. Assuming you have virtualenvwrapper installed, this is how you set up your fork for local development:

```
$ mkvirtualenv hvl_ccb
$ cd hvl_ccb/
$ python setup.py develop
$ pip install -r requirements_dev.txt
```

3. Create a branch for local development:

```
$ git checkout -b name-of-your-bugfix-or-feature
```

Now you can make your changes locally.

4. When you're done making changes, check that your changes pass flake8 and the tests, including testing other Python versions with tox:

```
$ flake8 hvl_ccb tests
$ python setup.py test or py.test
$ tox
```

To get flake8 and tox, just pip install them into your virtualenv. You can also use the provided make-like shell script to run flake8 and tests:

```
$ ./make.sh lint
$ ./make.sh test
```

5. Commit your changes and push your branch to GitLab:

```
$ git add .
$ git commit -m "Your detailed description of your changes."
$ git push origin name-of-your-bugfix-or-feature
```

6. Submit a merge request through the GitLab website.

## 5.3 Merge Request Guidelines

Before you submit a merge request, check that it meets these guidelines:

1. The merge request should include tests.
2. If the merge request adds functionality, the docs should be updated. Put your new functionality into a function with a docstring, and add the feature to the list in README.rst.
3. The merge request should work for Python 3.7. Check [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/merge\\_requests](https://gitlab.com/ethz_hvl/hvl_ccb/merge_requests) and make sure that the tests pass for all supported Python versions.

## 5.4 Tips

- To run tests from a single file:

```
$ py.test tests/test_hvl_ccb.py
```

or a single test function:

```
$ py.test tests/test_hvl_ccb.py::test_command_line_interface
```

- To add dependency, edit appropriate `*requirements` variable in the `setup.py` file and re-run:

```
$ python setup.py develop
```

- To generate a PDF version of the Sphinx documentation instead of HTML use:

```
$ rm -rf docs/hvl_ccb.rst docs/modules.rst docs/_build && sphinx-apidoc -o docs/_  
↳hvl_ccb && python -msphinx -M latexpdf docs/ docs/_build
```

This command can also be run through the make-like shell script:

```
$ ./make.sh docs-pdf
```

This requires a local installation of a LaTeX distribution, e.g. MikTeX.

## 5.5 Deploying

A reminder for the maintainers on how to deploy. Create release-N.M.K branch. Make sure all your changes are committed (including an entry in HISTORY.rst). Then run:

```
$ bumpversion patch # possible: major / minor / patch  
$ git push  
$ git push --tags  
$ make release
```

Merge the release branch into master and devel branches with `--no-ff` flag.

Optionally, go to [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/tags/vM.N.P/release/edit](https://gitlab.com/ethz_hvl/hvl_ccb/tags/vM.N.P/release/edit) and add release notes (e.g. changes lists).



# CHAPTER 6

---

## Credits

---

### 6.1 Development Lead

- Mikołaj Rybiński <mikolaj.rybinski@id.ethz.ch>
- (previously) David Graber <graber@eeh.ee.ethz.ch>

### 6.2 Contributors

- Henrik Menne <henrik.menne@eeh.ee.ethz.ch>
- Alise Chachereau <chachereau@eeh.ee.ethz.ch>



## History

---

### 7.1 0.3.5 (2020-02-18)

- Fix issue with reading integers from LabJack LJM Library (device's product ID, serial number etc.)
- Fix development requirements specification (tox version).

### 7.2 0.3.4 (2019-12-20)

- **New devices using serial connection:**
  - Heinzinger Digital Interface I/II and a Heinzinger PNC power supply
  - Q-switched Pulsed Laser and a laser attenuator from CryLas
  - Newport SMC100PP single axis motion controller for 2-phase stepper motors
  - Pfeiffer TPG controller (TPG 25x, TPG 26x and TPG 36x) for Compact pressure Gauges
- PEP 561 compatibility and related corrections for static type checking (now in CI)
- **Refactorings:**
  - Protected non-thread safe read and write in communication protocols
  - Device sequence mixin: start/stop, add/rm and lookup
  - `.format()` to f-strings
  - more enumerations and a quite some improvements of existing code
- Improved error docstrings (`:raises:` annotations) and extended tests for errors.

## 7.3 0.3.3 (2019-05-08)

- Use PyPI labjack-ljm (no external dependencies)

## 7.4 0.3.2 (2019-05-08)

- INSTALLATION.rst with LJMPython prerequisite info

## 7.5 0.3.1 (2019-05-02)

- readthedocs.org support

## 7.6 0.3 (2019-05-02)

- Prevent an automatic close of VISA connection when not used.
- Rhode & Schwarz RTO 1024 oscilloscope using VISA interface over `TCP::INSTR`.
- Extended tests incl. messages sent to devices.
- Added Supercube device using an OPC UA client
- Added Supercube 2015 device using an OPC UA client (for interfacing with old system version)

## 7.7 0.2.1 (2019-04-01)

- Fix issue with LJMPython not being installed automatically with setuptools.

## 7.8 0.2.0 (2019-03-31)

- LabJack LJM Library communication wrapper and LabJack device.
- Modbus TCP communication protocol.
- Schneider Electric ILS2T stepper motor drive device.
- Elektro-Automatik PSI9000 current source device and VISA communication wrapper.
- Separate configuration classes for communication protocols and devices.
- Simple experiment manager class.

## 7.9 0.1.0 (2019-02-06)

- Communication protocol base and serial communication implementation.
- Device base and MBW973 implementation.

# CHAPTER 8

---

## Indices and tables

---

- genindex
- modindex
- search



---

## Python Module Index

---

### h

hvl\_ccb, 117  
hvl\_ccb.comm, 19  
hvl\_ccb.comm.base, 7  
hvl\_ccb.comm.labjack\_ljm, 7  
hvl\_ccb.comm.modbus\_tcp, 10  
hvl\_ccb.comm.opc, 11  
hvl\_ccb.comm.serial, 14  
hvl\_ccb.comm.visa, 17  
hvl\_ccb.configuration, 116  
hvl\_ccb.dev, 114  
hvl\_ccb.dev.base, 53  
hvl\_ccb.dev.crylas, 55  
hvl\_ccb.dev.ea\_psi9000, 64  
hvl\_ccb.dev.heinzinger, 68  
hvl\_ccb.dev.labjack, 74  
hvl\_ccb.dev.mbw973, 79  
hvl\_ccb.dev.newport, 82  
hvl\_ccb.dev.pfeiffer\_tpg, 96  
hvl\_ccb.dev.rs\_rto1024, 101  
hvl\_ccb.dev.se\_ilis2t, 107  
hvl\_ccb.dev.supercube, 39  
hvl\_ccb.dev.supercube.base, 20  
hvl\_ccb.dev.supercube.constants, 25  
hvl\_ccb.dev.supercube.typ\_a, 36  
hvl\_ccb.dev.supercube.typ\_b, 38  
hvl\_ccb.dev.supercube2015, 52  
hvl\_ccb.dev.supercube2015.base, 39  
hvl\_ccb.dev.supercube2015.constants, 44  
hvl\_ccb.dev.supercube2015.typ\_a, 50  
hvl\_ccb.dev.visa, 113  
hvl\_ccb.experiment\_manager, 116  
hvl\_ccb.utils, 115  
hvl\_ccb.utils.enum, 115  
hvl\_ccb.utils.typing, 115



### A

A (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitCurrent attribute*), 71  
A (*hvl\_ccb.dev.supercube.constants.SupercubeOpcEndpoint attribute*), 35  
A (*hvl\_ccb.dev.supercube2015.constants.SupercubeOpcEndpoint attribute*), 50  
ABSOLUTE\_POSITION (*hvl\_ccb.dev.se\_ils2t.ILS2T.ActionsPtp attribute*), 107  
absolute\_position () (*hvl\_ccb.dev.se\_ils2t.ILS2T method*), 109  
absolute\_position\_and\_wait () (*hvl\_ccb.dev.se\_ils2t.ILS2T method*), 109  
AC (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 82  
AC\_DoubleStage\_150kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 34  
AC\_DoubleStage\_150kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 49  
AC\_DoubleStage\_200kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 34  
AC\_DoubleStage\_200kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 49  
AC\_SingleStage\_100kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 34  
AC\_SingleStage\_100kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 49  
AC\_SingleStage\_50kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 34  
AC\_SingleStage\_50kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 49  
attribute), 49  
acceleration (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 91  
ACCESS\_ENABLE (*hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr attribute*), 108  
ACCESS\_LOCK (*hvl\_ccb.comm.base.CommunicationProtocol attribute*), 7  
ACTION\_JOGL\_VALUE (*hvl\_ccb.dev.se\_ils2t.ILS2T attribute*), 107  
activated (*hvl\_ccb.dev.supercube.constants.BreakdownDetection attribute*), 30  
activated (*hvl\_ccb.dev.supercube2015.constants.BreakdownDetection attribute*), 45  
ACTIVE (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute*), 58  
active (*hvl\_ccb.dev.supercube.constants.OpcControl attribute*), 34  
add\_device () (*hvl\_ccb.dev.base.DeviceSequenceMixin method*), 53  
add\_device () (*hvl\_ccb.experiment\_manager.ExperimentManager method*), 116  
ADDR\_INCORRECT (*hvl\_ccb.dev.newport.NewportSMC100PPSerialComm attribute*), 92  
address (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 18  
address (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 91  
address () (*hvl\_ccb.comm.visa.VisaCommunicationConfig.InterfaceType method*), 18  
Alarm0 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 44  
Alarm1 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 26  
Alarm1 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 25  
Alarm1 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 44  
Alarm10 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 26  
Alarm10 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 26







tribute), 26  
 Alarm7 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 45  
 Alarm70 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm71 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm72 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm73 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm74 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm75 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm76 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm77 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm78 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm79 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm8 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm8 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 26  
 Alarm8 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 45  
 Alarm80 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm81 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm82 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 29  
 Alarm83 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm84 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm85 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm86 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm87 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm88 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm89 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm9 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm9 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 26  
 Alarm9 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 45  
 Alarm90 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm91 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm92 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm93 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm94 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm95 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm96 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm97 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm98 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarm99 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 30  
 Alarms (*class in hvl\_ccb.dev.supercube.constants*), 26  
 AlarmText (*class in hvl\_ccb.dev.supercube.constants*), 25  
 AlarmText (*class in hvl\_ccb.dev.supercube2015.constants*), 44  
 ANY (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.ConnectionType attribute*), 9  
 ANY (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.DeviceType attribute*), 9  
 ANY (*hvl\_ccb.dev.labjack.LabJack.DeviceType attribute*), 75  
 attenuation (*hvl\_ccb.dev.crylas.CryLasAttenuator attribute*), 55  
 AUTO (*hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes attribute*), 101  
 auto\_laser\_on (*hvl\_ccb.dev.crylas.CryLasLaserConfig attribute*), 61  
 AutoNumberNameEnum (*class in hvl\_ccb.utils.enum*), 115

## B

B (*hvl\_ccb.dev.supercube.constants.SupercubeOpcEndpoint attribute*), 35  
 B (*hvl\_ccb.dev.supercube2015.constants.SupercubeOpcEndpoint attribute*), 50  
 BA (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 82  
 backlash\_compensation (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 91  
 backup\_waveform () (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 102

```

bar (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG.PressureUnits    cee16 (hvl_ccb.dev.supercube2015.constants.GeneralSockets
attribute), 96                                attribute), 47
base_velocity (hvl_ccb.dev.newport.NewportSMC100PPConfiggeneric_type ()      (in      module
attribute), 91                                hvl_ccb.utils.typing), 115
baudrate (hvl_ccb.comm.serial.SerialCommunicationConfigcheck_master_slave_config ()
attribute), 15                                (hvl_ccb.dev.ea_psi9000.PSI9000      method),
baudrate (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig64Config
attribute), 57                                chunk_size (hvl_ccb.comm.visa.VisaCommunicationConfig
baudrate (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfigAttribute), 18
attribute), 63                                clean_values () (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig
baudrate (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
attribute), 73                                clean_values () (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig
baudrate (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfigmethod), 11
attribute), 81                                clean_values () (hvl_ccb.comm.opc.OpcUaCommunicationConfig
baudrate (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig
attribute), 94                                clean_values () (hvl_ccb.comm.serial.SerialCommunicationConfig
baudrate (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfigChnflg), 15
attribute), 100                                clean_values () (hvl_ccb.comm.visa.VisaCommunicationConfig
BH (hvl_ccb.dev.newport.NewportConfigCommands
attribute), 82                                method), 18
board (hvl_ccb.comm.visa.VisaCommunicationConfig
attribute), 18                                clean_values () (hvl_ccb.dev.base.EmptyConfig
BreakdownDetection (class      in
hvl_ccb.dev.supercube.constants), 30
BreakdownDetection (class      in
hvl_ccb.dev.supercube2015.constants), 45                                method), 54
clean_values () (hvl_ccb.dev.crylas.CryLasAttenuatorConfig
method), 55
clean_values () (hvl_ccb.dev.crylas.CryLasLaserConfig
method), 61
clean_values () (hvl_ccb.dev.ea_psi9000.PSI9000Config
method), 66
clean_values () (hvl_ccb.dev.heinzinger.HeinzingerConfig
method), 69
clean_values () (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
method), 73
bytesize (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig
attribute), 57                                clean_values () (hvl_ccb.dev.mbw973.MBW973Config
bytesize (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig
attribute), 15                                clean_values () (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
bytesize (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig
attribute), 63                                method), 73
bytesize (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
attribute), 73                                clean_values () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig
bytesize (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfigmethod), 99
attribute), 81                                clean_values () (hvl_ccb.dev.se_ils2t.ILS2TConfig
bytesize (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig
attribute), 94                                clean_values () (hvl_ccb.dev.supercube.base.SupercubeConfiguration
bytesize (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfigChnflg), 23
attribute), 100                                clean_values () (hvl_ccb.dev.supercube2015.base.SupercubeConfig
method), 42
C
C (hvl_ccb.dev.labjack.LabJack.TemperatureUnit
attribute), 75                                close () (hvl_ccb.comm.base.CommunicationProtocol
method), 7
C (hvl_ccb.dev.labjack.LabJack.ThermocoupleType
attribute), 75                                close () (hvl_ccb.comm.labjack_ljm.LJMCommunication
method), 8
calibration_factor
(hvl_ccb.dev.crylas.CryLasLaserConfig
attribute), 61                                close () (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication
method), 10
cee16 (hvl_ccb.dev.supercube.constants.GeneralSockets
attribute), 32                                close () (hvl_ccb.comm.opc.OpcUaCommunication
method), 12
close () (hvl_ccb.comm.serial.SerialCommunication
method), 14
close () (hvl_ccb.comm.visa.VisaCommunication
method), 16

```

*method), 17*  
`close_shutter()` (*hvl\_ccb.dev.crylas.CryLasLaser* *method), 59*  
`CLOSED` (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersShutter* *attribute), 58*  
`CLOSED` (*hvl\_ccb.dev.crylas.CryLasLaserShutterStatus* *attribute), 64*  
`closed` (*hvl\_ccb.dev.supercube.constants.DoorStatus* *attribute), 31*  
`closed` (*hvl\_ccb.dev.supercube.constants.EarthingStickStatus* *attribute), 31*  
`closed` (*hvl\_ccb.dev.supercube2015.constants.DoorStatus* *attribute), 45*  
`closed` (*hvl\_ccb.dev.supercube2015.constants.EarthingStickStatus* *attribute), 46*  
`CMD_EXEC_ERROR` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerErrorConfigurationMixin* *attribute), 92*  
`CMD_NOT_ALLOWED` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerErrorService* *attribute), 92*  
`CMD_NOT_ALLOWED_CC` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMD_NOT_ALLOWED_CONFIGURATION` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMD_NOT_ALLOWED_DISABLE` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMD_NOT_ALLOWED_HOMING` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMD_NOT_ALLOWED_MOVING` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMD_NOT_ALLOWED_NOT_REFERENCED` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMD_NOT_ALLOWED_PP` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMD_NOT_ALLOWED_READY` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CMR` (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes* *attribute), 96*  
`CODE_OR_ADDR_INVALID` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`com` (*hvl\_ccb.dev.base.SingleCommDevice* *attribute), 54*  
`COM_TIMEOUT` (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationControllerError* *attribute), 92*  
`CommunicationProtocol` (*class* *in* *config\_cls()* (*hvl\_ccb.dev.se\_il2t.ILS2T* *static method), 109*  
`config` (*hvl\_ccb.configuration.ConfigurationMixin* *at-* *config\_cls()* (*hvl\_ccb.dev.se\_il2t.ILS2TModbusTcpCommunication*

```

    static method), 112
config_cls () (hvl_ccb.dev.supercube.base.SuperCubeBase current_lower_limit
    static method), 20 (hvl_ccb.dev.ea_psi9000.PSI9000Config
config_cls () (hvl_ccb.dev.supercube.base.SuperCubeOpcUaCom attribute), 66
    static method), 23 current_primary (hvl_ccb.dev.supercube.constants.Power
config_cls () (hvl_ccb.dev.supercube.typ_a.SuperCubeAOpcUaCom attribute), 34
    static method), 36 current_primary (hvl_ccb.dev.supercube2015.constants.Power
config_cls () (hvl_ccb.dev.supercube.typ_b.SuperCubeBOpcUaCom attribute), 38
    static method), 38 current_upper_limit
config_cls () (hvl_ccb.dev.supercube2015.base.SuperCube2015B (hvl_ccb.dev.ea_psi9000.PSI9000Config
    static method), 39 attribute), 66
config_cls () (hvl_ccb.dev.supercube2015.base.SuperCubeOpcUaCommunication (hvl_ccb.dev.newport.NewportSMC100PPConfig.Ho
    static method), 43 attribute), 90
config_cls () (hvl_ccb.dev.supercube2015.typ_a.SuperCubeAOpcUaCommunication D
    static method), 51
config_cls () (hvl_ccb.dev.visa.VisaDevice static datachange_notification()
    method), 113 (hvl_ccb.comm.opc.OpcUaSubHandler
configdataclass () (in module datachange_notification()
    hvl_ccb.configuration), 116 (hvl_ccb.dev.supercube.base.SuperCubeSubscriptionHandler
configuration_save_json () method), 24
    (hvl_ccb.configuration.ConfigurationMixin
method), 116 datachange_notification()
ConfigurationMixin (class in hvl_ccb.dev.supercube2015.base.SuperCubeSubscriptionHandler
    hvl_ccb.configuration), 116 method), 44
connection_type (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig DC_DoubleStage_280kV
    attribute), 9 (hvl_ccb.dev.supercube.constants.PowerSetup
create_serial_port () attribute), 34
    (hvl_ccb.comm.serial.SerialCommunicationConfig DC_DoubleStage_280kV
    method), 15 (hvl_ccb.dev.supercube2015.constants.PowerSetup
CryLasAttenuator (class in hvl_ccb.dev.crylas), 55 attribute), 49
CryLasAttenuatorConfig (class in DC_SingleStage_140kV
    hvl_ccb.dev.crylas), 55 (hvl_ccb.dev.supercube.constants.PowerSetup
CryLasAttenuatorError, 56 attribute), 34
CryLasAttenuatorSerialCommunication (class in hvl_ccb.dev.supercube2015.constants.PowerSetup
    (class in hvl_ccb.dev.crylas), 56 attribute), 49
CryLasAttenuatorSerialCommunicationConfig (class in DC_VOLTAGE_TOO_LOW
    (class in hvl_ccb.dev.crylas), 56 (hvl_ccb.dev.newport.NewportSMC100PP.MotorErrors
CryLasLaser (class in hvl_ccb.dev.crylas), 58 attribute), 83
CryLasLaser.AnswersShutter (class in default_com_cls()
    hvl_ccb.dev.crylas), 58 (hvl_ccb.dev.base.SingleCommDevice static
CryLasLaser.AnswersStatus (class in method), 54
    hvl_ccb.dev.crylas), 58 default_com_cls()
CryLasLaser.LaserStatus (class in (hvl_ccb.dev.crylas.CryLasAttenuator static
    hvl_ccb.dev.crylas), 58 method), 55
CryLasLaser.RepetitionRates (class in default_com_cls()
    hvl_ccb.dev.crylas), 59 (hvl_ccb.dev.crylas.CryLasLaser static
CryLasLaserConfig (class in hvl_ccb.dev.crylas), 61 method), 59
CryLasLaserError, 62 default_com_cls()
CryLasLaserNotReadyError, 62 (hvl_ccb.dev.ea_psi9000.PSI9000 static
CryLasLaserSerialCommunication (class in method), 65
    hvl_ccb.dev.crylas), 62 default_com_cls()
CryLasLaserSerialCommunicationConfig (class in (hvl_ccb.dev.heinzinger.HeinzingerDI static
    hvl_ccb.dev.crylas), 63 method), 70
CryLasLaserShutterStatus (class in
    hvl_ccb.dev.crylas), 63

```

```

default_com_cls() (hvl_ccb.dev.labjack.LabJack
    static method), 76
default_com_cls() (hvl_ccb.dev.mbw973.MBW973 static method),
    79
default_com_cls() (hvl_ccb.dev.newport.NewportSMC100PP
    static method), 84
default_com_cls() (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG static
    method), 97
default_com_cls() (hvl_ccb.dev.rs_rto1024.RTO1024 static
    method), 102
default_com_cls() (hvl_ccb.dev.se_il2t.ILS2T static
    method), 109
default_com_cls() (hvl_ccb.dev.supercube.base.SuperCubeBase
    static method), 20
default_com_cls() (hvl_ccb.dev.supercube.typ_a.SuperCubeWithFU
    static method), 37
default_com_cls() (hvl_ccb.dev.supercube.typ_b.SuperCubeB
    static method), 38
default_com_cls() (hvl_ccb.dev.supercube2015.base.SuperCube2015Base
    static method), 39
default_com_cls() (hvl_ccb.dev.supercube2015.typ_a.SuperCube2015WithFU
    static method), 50
default_com_cls() (hvl_ccb.dev.visa.VisaDevice
    static method), 113
DEFAULT_IO_SCANNING_CONTROL_VALUES
    (hvl_ccb.dev.se_il2t.ILS2T attribute), 107
default_number_of_recordings
    (hvl_ccb.dev.heinzinger.HeinzingerConfig
    attribute), 69
Device (class in hvl_ccb.dev.base), 53
device_type (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig
    attribute), 9
DeviceExistingException, 53
DeviceSequenceMixin (class in hvl_ccb.dev.base),
    53
DIOChannel (hvl_ccb.dev.labjack.LabJack attribute),
    75
DISABLE (hvl_ccb.dev.newport.NewportStates attribute),
    95
disable() (hvl_ccb.dev.se_il2t.ILS2T method), 109
DISABLE_FROM_JOGGING
    (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    attribute), 83
DISABLE_FROM_MOVING
    (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    attribute), 83
DISABLE_FROM_READY
    (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    attribute), 83
DisableEspStageCheck
    (hvl_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig
    attribute), 90
DISPLACEMENT_OUT_OF_LIMIT
    (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication.C
    attribute), 92
do_ioscanning_write()
    (hvl_ccb.dev.se_il2t.ILS2T method), 109
Door (class in hvl_ccb.dev.supercube.constants), 30
DoorStatus (class
    hvl_ccb.dev.supercube.constants), 30
DoorStatus (class
    hvl_ccb.dev.supercube2015.constants), 45

```

**E**

```

E (hvl_ccb.dev.labjack.LabJack.ThermocoupleType attribute),
    75
EarthingStick (class
    hvl_ccb.dev.supercube.constants), 31
EarthingStick (class
    hvl_ccb.dev.supercube2015.constants), 45
EarthingStickStatus (class
    hvl_ccb.dev.supercube.constants), 31
EarthingStickStatus (class
    hvl_ccb.dev.supercube2015.constants), 46
EEPROM_ACCESS_ERROR
    (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication.C
    attribute), 92
EIGHT (hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum
    attribute), 69
EIGHTBITS (hvl_ccb.comm.serial.SerialCommunicationBytesize
    attribute), 15
EmptyConfig (class in hvl_ccb.dev.base), 54
enable() (hvl_ccb.dev.se_il2t.ILS2T method), 109
EnableEspStageCheck
    (hvl_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig
    attribute), 90
ENCODING (hvl_ccb.comm.serial.SerialCommunication
    attribute), 14
EndOfRunSwitch (hvl_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch
    attribute), 90
EndOfRunSwitch_and_Index
    (hvl_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch
    attribute), 90
endpoint_name (hvl_ccb.comm.opc.OpcUaCommunicationConfig
    attribute), 13
endpoint_name (hvl_ccb.dev.supercube.typ_a.SuperCubeAOpcUaConfig
    attribute), 36
endpoint_name (hvl_ccb.dev.supercube.typ_b.SuperCubeBOpcUaConfig
    attribute), 38

```

```

endpoint_name (hvl_ccb.dev.supercube2015.typ_a.SuperCube2015DoorOpenerExperiment_manager.ExperimentStatus
attribute), 52
error (hvl_ccb.dev.supercube.constants.DoorStatus attribute), 31
error (hvl_ccb.dev.supercube.constants.EarthingStickStatus attribute), 31
error (hvl_ccb.dev.supercube2015.constants.DoorStatus attribute), 45
error (hvl_ccb.dev.supercube2015.constants.EarthingStickStatus attribute), 46
Error (hvl_ccb.dev.supercube2015.constants.SafetyStatus attribute), 35
Error (hvl_ccb.dev.supercube2015.constants.SafetyStatus attribute), 50
ERROR (hvl_ccb.experiment_manager.ExperimentStatus attribute), 117
Errors (class in hvl_ccb.dev.supercube.constants), 32
Errors (class in hvl_ccb.dev.supercube2015.constants), 47
ESP_STAGE_NAME_INVALID
    (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationController, 92)
ETHERNET (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig attribute), 9
EVEN (hvl_ccb.comm.serial.SerialCommunicationParity attribute), 16
event_notification ()
    (hvl_ccb.comm.opc.OpcUaSubHandler method), 14
exit_configuration ()
    (hvl_ccb.dev.newport.NewportSMC100PP method), 84
exit_configuration_wait_sec
    (hvl_ccb.dev.newport.NewportSMC100PPConfig attribute), 91
ExperimentError, 116
ExperimentManager (class
    hvl_ccb.experiment_manager), 116
ExperimentStatus (class
    hvl_ccb.experiment_manager), 117
External (hvl_ccb.dev.supercube.constants.PowerSetup attribute), 34
External (hvl_ccb.dev.supercube2015.constants.PowerSetup attribute), 49
F
F (hvl_ccb.dev.labjack.LabJack.TemperatureUnit attribute), 75
FAST (hvl_ccb.dev.se_ils2t.ILS2T.Ref16Jog attribute), 108
file_copy () (hvl_ccb.dev.rs_rto1024.RTO1024 method), 102
finish () (hvl_ccb.experiment_manager.ExperimentManager method), 117
FINISHING (hvl_ccb.experiment_manager.ExperimentStatus attribute), 117
FIVEBITS (hvl_ccb.comm.serial.SerialCommunicationBytesize attribute), 15
FLT_INFO (hvl_ccb.dev.se_ils2t.ILS2T.RegAddr attribute), 108
FLT_MEM_DEL (hvl_ccb.dev.se_ils2t.ILS2T.RegAddr attribute), 108
RESET (hvl_ccb.dev.se_ils2t.ILS2T.RegAddr attribute), 108
FOLLOWING_ERROR (hvl_ccb.dev.newport.NewportSMC100PP.MotorError attribute), 83
force_value () (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig method), 9
force_value () (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication method), 11
force_value () (hvl_ccb.comm.opc.OpcUaCommunicationConfig method), 13
force_value () (hvl_ccb.comm.serial.SerialCommunicationConfig method), 16
force_value () (hvl_ccb.dev.base.EmptyConfig method), 54
force_value () (hvl_ccb.dev.crylas.CryLasAttenuatorConfig method), 55
force_value () (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunication method), 57
force_value () (hvl_ccb.dev.crylas.CryLasLaserConfig method), 61
force_value () (hvl_ccb.dev.crylas.CryLasLaserSerialCommunication method), 63
force_value () (hvl_ccb.dev.ea_psi9000.PSI9000Config method), 66
force_value () (hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig method), 68
force_value () (hvl_ccb.dev.heinzinger.HeinzingerConfig method), 69
force_value () (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunication method), 73
force_value () (hvl_ccb.dev.mbw973.MBW973Config method), 80
force_value () (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig method), 81
force_value () (hvl_ccb.dev.newport.NewportSMC100PPConfig method), 91
force_value () (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication method), 94
force_value () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig method), 99
force_value () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunication method), 101

```

```

force_value() (hvl_ccb.dev.rs_rto1024.RTO1024Config get (hvl_ccb.dev.supercube.constants.Door attribute),
               method), 105                                         30
force_value() (hvl_ccb.dev.rs_rto1024.RTO1024VisaConfig (hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio
               method), 107                                         attribute), 33
force_value() (hvl_ccb.dev.se_ils2t.ILS2TConfig get (hvl_ccb.dev.supercube.constants.MeasurementsScaledInput
               method), 111                                         attribute), 33
force_value() (hvl_ccb.dev.se_ils2t.ILS2TModbusTcpConfig (hvl_ccb.dev.supercube2015.constants.AlarmText at-
               method), 112                                         tribute), 45
force_value() (hvl_ccb.dev.supercube.base.SupercubeConfig (hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio
               method), 23                                         attribute), 48
force_value() (hvl_ccb.dev.supercube.base.SupercubeOpcUaConfiguration (hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInput
               method), 24                                         attribute), 48
force_value() (hvl_ccb.dev.supercube.typ_a.SupercubeAOpcUaConfiguration () (hvl_ccb.dev.newport.NewportSMC100PP
               method), 36
force_value() (hvl_ccb.dev.supercube.typ_b.SupercubeBOpcUaConfiguration () (hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio
               method), 38                                         attribute), 48
force_value() (hvl_ccb.dev.supercube2015.base.SupercubeConfiguration (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig.De-
               method), 42                                         attribute), 9
force_value() (hvl_ccb.dev.supercube2015.base.SupercubeOpcUaConfiguration (hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio
               method), 43                                         attribute), 75
force_value() (hvl_ccb.dev.supercube2015.typ_a.SupercubeAOpcUaConfiguration () (hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio
               method), 52                                         attribute), 76
force_value() (hvl_ccb.dev.visa.VisaDeviceConfig get_cee16_socket () (hvl_ccb.dev.supercube.base.SupercubeBase
               method), 114
FOUR (hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum method), 20
attribute), 69                                         get_cee16_socket ()
FREERUN (hvl_ccb.dev.rs_rto1024.RTO1024.TriggerModes (hvl_ccb.dev.supercube2015.base.Supercube2015Base
attribute), 101                                         method), 39
frequency (hvl_ccb.dev.supercube.constants.Power at- get_controller_information ()
tribute), 34                                         (hvl_ccb.dev.newport.NewportSMC100PP
frequency (hvl_ccb.dev.supercube2015.constants.Power method), 85
attribute), 48                                         get_current () (hvl_ccb.dev.heinzinger.HeinzingerDI
FRM (hvl_ccb.dev.newport.NewportConfigCommands at- method), 70
tribute), 82                                         get_dc_volt () (hvl_ccb.dev.se_ils2t.ILS2T method),
from_json () (hvl_ccb.configuration.ConfigurationMixin 109
class method), 116                                         get_device () (hvl_ccb.dev.base.DeviceSequenceMixin
FRS (hvl_ccb.dev.newport.NewportConfigCommands at- method), 53
tribute), 82                                         get_devices () (hvl_ccb.dev.base.DeviceSequenceMixin
fso_reset () (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU method), 53
method), 37                                         get_digital_input ()
fso_reset () (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU (hvl_ccb.dev.labjack.LabJack method), 76
method), 50                                         get_door_status ()
                                         (hvl_ccb.dev.supercube.base.SupercubeBase
                                         method), 20
G
GeneralSockets (class in get_door_status ()
hvl_ccb.dev.supercube.constants), 32 (hvl_ccb.dev.supercube2015.base.Supercube2015Base
GeneralSockets (class in get_earthing_manual ()
hvl_ccb.dev.supercube2015.constants), 47 (hvl_ccb.dev.supercube.base.SupercubeBase
GeneralSupport (class in get_earthing_manual ()
hvl_ccb.dev.supercube.constants), 32 (hvl_ccb.dev.supercube2015.base.Supercube2015Base
GeneralSupport (class in get_earthing_manual ()
hvl_ccb.dev.supercube2015.constants), 47 (hvl_ccb.dev.supercube2015.base.Supercube2015Base
get (hvl_ccb.dev.supercube.constants.AlarmText at- at- get_earthing_status ()
tribute), 26                                         attribute), 40

```

```

(hvl_ccb.dev.supercube.base.SupercubeBase
    method), 20
get_earthing_status ()           (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 40
get_error_code ()               (hvl_ccb.dev.se_il2t.ILS2T
    method), 109
get_error_queue ()              (hvl_ccb.dev.visa.VisaDevice
    method), 113
get_frequency ()                (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 37
get_frequency ()                (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 50
get_fso_active ()               (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 37
get_fso_active ()               (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 50
get_full_scale_mbar ()          (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG
    method), 97
get_full_scale_unitless ()      (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG
    method), 97
get_identification ()           (hvl_ccb.dev.visa.VisaDevice method), 113
get_interface_version ()        (hvl_ccb.dev.heinzinger.HeinzingerDI
    method), 70
get_max_voltage ()              (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 37
get_max_voltage ()              (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 51
get_measurement_ratio ()         (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 20
get_measurement_ratio ()         (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 40
get_measurement_voltage ()       (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 21
get_measurement_voltage ()       (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 40
get_motor_configuration ()       (hvl_ccb.dev.newport.NewportSMC100PP
    method), 85
get_move_duration ()             (hvl_ccb.dev.newport.NewportSMC100PP
    method), 85
get_negative_software_limit ()   (hvl_ccb.dev.newport.NewportSMC100PP
    method), 85
get_number_of_recordings ()      (hvl_ccb.dev.heinzinger.HeinzingerDI
    method), 70
get_output ()                   (hvl_ccb.dev.ea_psi9000.PSI9000
    method), 65
get_position ()                 (hvl_ccb.dev.newport.NewportSMC100PP
    method), 85
get_position ()                 (hvl_ccb.dev.se_il2t.ILS2T
    method), 110
get_positive_software_limit ()  (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 86
get_power_setup ()              (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 37
get_primary_current ()           (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 37
get_primary_current ()           (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 51
get_primary_voltage ()           (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 37
get_primary_voltage ()           (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 51
get_product_id ()                (hvl_ccb.dev.labjack.LabJack
    method), 76
get_product_name ()              (hvl_ccb.dev.labjack.LabJack
    method), 76
get_product_type ()              (hvl_ccb.dev.labjack.LabJack
    method), 76
get_pulse_energy_and_rate ()     (hvl_ccb.dev.crylas.CryLasLaser
    method), 59
get_sbus_rh ()                  (hvl_ccb.dev.labjack.LabJack
    method), 77
get_sbus_temp ()                (hvl_ccb.dev.labjack.LabJack
    method), 77
get_serial_number ()             (hvl_ccb.dev.heinzinger.HeinzingerDI
    method), 70
get_serial_number ()             (hvl_ccb.dev.labjack.LabJack
    method), 77
get_state ()                     (hvl_ccb.dev.newport.NewportSMC100PP
    method), 86
get_status ()                    (hvl_ccb.dev.se_il2t.ILS2T
    method), 110
get_status ()                    (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 85

```

*method), 21*  
`get_status () (hvl_ccb.dev.supercube2015.base.Supercube2015Base attribute), 50`  
*method), 40*  
`get_support_input ()`  
*(hvl\_ccb.dev.supercube.base.SupercubeBase method), 21*  
`get_support_input ()`  
*(hvl\_ccb.dev.supercube2015.base.Supercube2015Base method), 40*  
`get_support_output ()`  
*(hvl\_ccb.dev.supercube.base.SupercubeBase method), 21*  
`get_support_output ()`  
*(hvl\_ccb.dev.supercube2015.base.Supercube2015Base method), 40*  
`get_system_lock ()`  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000 method), HeinzingerPNCDeviceNotRecognizedException, 65*  
`get_t13_socket () (hvl_ccb.dev.supercube.base.Supercube2015Base attribute), 21`  
`get_t13_socket () (hvl_ccb.dev.supercube2015.base.Supercube2015Base attribute), 40`  
`get_target_voltage ()`  
*(hvl\_ccb.dev.supercube.typ\_a.SupercubeWithFU method), HeinzingerPNCMaxCurrentExceededException, 37*  
`get_target_voltage ()`  
*(hvl\_ccb.dev.supercube2015.typ\_a.Supercube2015WithFU method), HeinzingerSerialCommunication (class in hvl\_ccb.dev.heinzinger), 51*  
`get_temperature () (hvl_ccb.dev.se_ilst.ILS2T attribute), 110`  
`get.timestamps () (hvl_ccb.dev.rs_rto1024.RTO1024 attribute), 102`  
`get_ui_lower_limits ()`  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000 method), HeinzingerSerialCommunicationConfig (class in hvl\_ccb.dev.heinzinger), 65*  
`get_uip_upper_limits ()`  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000 method), HeinzingerSerialCommunicationConfig (class in hvl\_ccb.dev.heinzinger), 65*  
`get_voltage () (hvl_ccb.dev.heinzinger.HeinzingerDI attribute), 70`  
`get_voltage_current_setpoint ()`  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000 method), HeinzingerSerialCommunicationConfig (class in hvl\_ccb.dev.heinzinger), 65*  
`go_home () (hvl_ccb.dev.newport.NewportSMC100PP attribute), 86`  
`go_to_configuration ()`  
*(hvl\_ccb.dev.newport.NewportSMC100PP method), HeinzingerSerialCommunicationConfig (class in hvl\_ccb.dev.heinzinger), 86*  
`GreenNotReady (hvl_ccb.dev.supercube.constants.SafetyStatus attribute), 35`  
`GreenNotReady (hvl_ccb.dev.supercube2015.constants.SafetyStatus attribute), 50`  
`GreenReady (hvl_ccb.dev.supercube.constants.SafetyStatus attribute), 35`

**H**

*GreenReady (hvl\_ccb.dev.supercube2015.constants.SafetyStatus attribute), 50*  
**HARDWARE (hvl\_ccb.dev.crylas.CryLasLaser.RepetitionRates attribute), 59**  
**HEAD (hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute), 58**  
**HeinzingerConfig (class in hvl\_ccb.dev.heinzinger), 69**  
**HeinzingerConfig.RecordingsEnum (class in hvl\_ccb.dev.heinzinger), 69**  
**HeinzingerDI (class in hvl\_ccb.dev.heinzinger), 70**  
**HeinzingerPNC (class in hvl\_ccb.dev.heinzinger), 71**  
**HeinzingerPNC.UnitCurrent (class in hvl\_ccb.dev.heinzinger), 71**  
**HeinzingerPNC.UnitVoltage (class in hvl\_ccb.dev.heinzinger), 71**  
**HeinzingerPNCDeviceNotRecognizedException, 72**  
**HeinzingerPNCError, 72**  
**HeinzingerPNCMaxCurrentExceededException, 72**  
**HeinzingerPNCMaxVoltageExceededException, 72**  
**HeinzingerSerialCommunication (class in hvl\_ccb.dev.heinzinger), 72**  
**HeinzingerSerialCommunicationConfig (class in hvl\_ccb.dev.heinzinger), 73**  
**HIGH (hvl\_ccb.dev.labjack.LabJack.DIOStatus attribute), 75**  
**home\_search\_polling\_interval (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 91**  
**home\_search\_timeout (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 91**  
**home\_search\_type (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 91**  
**home\_search\_velocity (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 91**  
**HOME\_STARTED (hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication attribute), 92**  
**HomeSwitch (hvl\_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch attribute), 90**  
**HomeSwitch\_and\_Index (hvl\_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch attribute), 91**  
**HOMING (hvl\_ccb.dev.newport.NewportStates attribute), 95**  
**HOMING\_FROM\_RS232 (hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84**

HOMING\_FROM\_SMC (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages*, 91  
     attribute), 84  
 HOMING\_TIMEOUT (*hvl\_ccb.dev.newport.NewportSMC100PP.MotorErrors*  
     attribute), 83  
 host (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus*  
     attribute), 11  
 host (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig* identifier (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig*  
     attribute), 13  
 host (*hvl\_ccb.comm.visa.VisaCommunicationConfig* at-  
     tribute), 19  
 hPascal (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits*  
     attribute), 96  
 HT (*hvl\_ccb.dev.newport.NewportConfigCommands* at-  
     tribute), 82  
*hvl\_ccb* (*module*), 117  
*hvl\_ccb.comm* (*module*), 19  
*hvl\_ccb.comm.base* (*module*), 7  
*hvl\_ccb.comm.labjack\_ljm* (*module*), 7  
*hvl\_ccb.comm.modbus\_tcp* (*module*), 10  
*hvl\_ccb.comm.opc* (*module*), 11  
*hvl\_ccb.comm.serial* (*module*), 14  
*hvl\_ccb.comm.visa* (*module*), 17  
*hvl\_ccb.configuration* (*module*), 116  
*hvl\_ccb.dev* (*module*), 114  
*hvl\_ccb.dev.base* (*module*), 53  
*hvl\_ccb.dev.crylas* (*module*), 55  
*hvl\_ccb.dev.ea\_psi9000* (*module*), 64  
*hvl\_ccb.dev.heinzinger* (*module*), 68  
*hvl\_ccb.dev.labjack* (*module*), 74  
*hvl\_ccb.dev.mbw973* (*module*), 79  
*hvl\_ccb.dev.newport* (*module*), 82  
*hvl\_ccb.dev.pfeiffer\_tpg* (*module*), 96  
*hvl\_ccb.dev.rs\_rto1024* (*module*), 101  
*hvl\_ccb.dev.se\_ilis2t* (*module*), 107  
*hvl\_ccb.dev.supercube* (*module*), 39  
*hvl\_ccb.dev.supercube.base* (*module*), 20  
*hvl\_ccb.dev.supercube.constants* (*module*),  
     25  
*hvl\_ccb.dev.supercube.typ\_a* (*module*), 36  
*hvl\_ccb.dev.supercube.typ\_b* (*module*), 38  
*hvl\_ccb.dev.supercube2015* (*module*), 52  
*hvl\_ccb.dev.supercube2015.base* (*module*),  
     39  
*hvl\_ccb.dev.supercube2015.constants*  
     (*module*), 44  
*hvl\_ccb.dev.supercube2015.typ\_a* (*module*),  
     50  
*hvl\_ccb.dev.visa* (*module*), 113  
*hvl\_ccb.experiment\_manager* (*module*), 116  
*hvl\_ccb.utils* (*module*), 115  
*hvl\_ccb.utils.enum* (*module*), 115  
*hvl\_ccb.utils.typing* (*module*), 115  
*hysteresis\_compensation*  
     (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
     attribute), 47  
     Identification\_error  
     attribute), 96  
     identifier (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig*  
     attribute), 9  
     identify\_device()  
         (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*  
         method), 71  
     identify\_sensors()  
         (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG* method),  
         97  
     IKR (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes*  
     attribute), 96  
     IKR11 (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes*  
     attribute), 96  
     IKR9 (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes*  
     attribute), 96  
     ILS2T (class in *hvl\_ccb.dev.se\_ilis2t*), 107  
     ILS2T.ActionsPtp (class in *hvl\_ccb.dev.se\_ilis2t*),  
     107  
     ILS2T.Mode (class in *hvl\_ccb.dev.se\_ilis2t*), 107  
     ILS2T.Ref16Jog (class in *hvl\_ccb.dev.se\_ilis2t*), 108  
     ILS2T.RegAddr (class in *hvl\_ccb.dev.se\_ilis2t*), 108  
     ILS2T.RegDatatype (class in *hvl\_ccb.dev.se\_ilis2t*),  
     108  
     ILS2T.State (class in *hvl\_ccb.dev.se\_ilis2t*), 109  
     ILS2TConfig (class in *hvl\_ccb.dev.se\_ilis2t*), 111  
     ILS2TException, 112  
     ILS2TModbusTcpCommunication (class in  
         *hvl\_ccb.dev.se\_ilis2t*), 112  
     ILS2TModbusTcpCommunicationConfig (class  
         in *hvl\_ccb.dev.se\_ilis2t*), 112  
     IMR (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes*  
     attribute), 96  
     in\_1\_1 (*hvl\_ccb.dev.supercube.constants.GeneralSupport*  
     attribute), 32  
     in\_1\_1 (*hvl\_ccb.dev.supercube2015.constants.GeneralSupport*  
     attribute), 47  
     in\_1\_2 (*hvl\_ccb.dev.supercube.constants.GeneralSupport*  
     attribute), 32  
     in\_1\_2 (*hvl\_ccb.dev.supercube2015.constants.GeneralSupport*  
     attribute), 47  
     in\_2\_1 (*hvl\_ccb.dev.supercube.constants.GeneralSupport*  
     attribute), 32  
     in\_2\_1 (*hvl\_ccb.dev.supercube2015.constants.GeneralSupport*  
     attribute), 47  
     in\_2\_2 (*hvl\_ccb.dev.supercube.constants.GeneralSupport*  
     attribute), 32  
     in\_2\_2 (*hvl\_ccb.dev.supercube2015.constants.GeneralSupport*  
     attribute), 47

```
in_3_1(hvl_ccb.dev.supercube.constants.GeneralSupport Initialize(hvl_ccb.dev.supercube.constants.SafetyStatus
attribute), 32 attribute), 35
in_3_1(hvl_ccb.dev.supercube2015.constants.GeneralSupport Initialize(hvl_ccb.dev.supercube2015.constants.SafetyStatus
attribute), 47 attribute), 50
in_3_2(hvl_ccb.dev.supercube.constants.GeneralSupport INITIALIZING(hvl_ccb.experiment_manager.ExperimentStatus
attribute), 32 attribute), 117
in_3_2(hvl_ccb.dev.supercube2015.constants.GeneralSupport put(hvl_ccb.dev.supercube.constants.GeneralSupport
attribute), 47 attribute), 33
in_4_1(hvl_ccb.dev.supercube.constants.GeneralSupport input(hvl_ccb.dev.supercube2015.constants.GeneralSupport
attribute), 32 attribute), 47
in_4_1(hvl_ccb.dev.supercube2015.constants.GeneralSupport put_1(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio
attribute), 47 attribute), 33
in_4_2(hvl_ccb.dev.supercube.constants.GeneralSupport input_1(hvl_ccb.dev.supercube.constants.MeasurementsScaledInput
attribute), 32 attribute), 33
in_4_2(hvl_ccb.dev.supercube2015.constants.GeneralSupport put_1(hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio
attribute), 47 attribute), 48
in_5_1(hvl_ccb.dev.supercube.constants.GeneralSupport input_1(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInput
attribute), 32 attribute), 48
in_5_1(hvl_ccb.dev.supercube2015.constants.GeneralSupport put_2(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio
attribute), 47 attribute), 33
in_5_2(hvl_ccb.dev.supercube.constants.GeneralSupport input_2(hvl_ccb.dev.supercube.constants.MeasurementsScaledInput
attribute), 32 attribute), 33
in_5_2(hvl_ccb.dev.supercube2015.constants.GeneralSupport put_2(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInput
attribute), 47 attribute), 48
in_6_1(hvl_ccb.dev.supercube.constants.GeneralSupport input_3(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio
attribute), 32 attribute), 33
in_6_1(hvl_ccb.dev.supercube2015.constants.GeneralSupport put_3(hvl_ccb.dev.supercube.constants.MeasurementsScaledInput
attribute), 47 attribute), 33
in_6_2(hvl_ccb.dev.supercube.constants.GeneralSupport input_3(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInput
attribute), 32 attribute), 48
in_6_2(hvl_ccb.dev.supercube2015.constants.GeneralSupport put_4(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio
attribute), 47 attribute), 33
INACTIVE(hvl_ccb.dev.crylas.CryLasLaser.AnswersStatus input_4(hvl_ccb.dev.supercube.constants.MeasurementsScaledInput
attribute), 58 attribute), 33
inactive(hvl_ccb.dev.supercube.constants.DoorStatus input_4(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInput
attribute), 31 attribute), 48
inactive(hvl_ccb.dev.supercube.constants.EarthingStickSM102(hvl_ccb.dev.se_ilst2.ILS2T.RegDatatype at-
attribute), 31 tribute), 109
inactive(hvl_ccb.dev.supercube2015.constants.DoorStatus interface_type(hvl_ccb.comm.visa.VisaCommunicationConfig
attribute), 45 attribute), 19
inactive(hvl_ccb.dev.supercube2015.constants.EarthingStickSM102 interface_type(hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunication
attribute), 46 attribute), 68
init_attenuation(hvl_ccb.dev.crylas.CryLasAttenuatorConfig interface_type(hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunication
attribute), 56 attribute), 107
init_monitored_nodes() internal(hvl_ccb.dev.labjack.LabJack.CjcType
(hvl_ccb.comm.opc.OpcUaCommunication attribute), 75
method), 12 Internal(hvl_ccb.dev.supercube.constants.PowerSetup
attribute), 35
init_shutter_status Internal(hvl_ccb.dev.supercube2015.constants.PowerSetup
(hvl_ccb.dev.crylas.CryLasLaserConfig attribute), 49
tribute), 62
initialize() (hvl_ccb.dev.newport.NewportSMC100PPInvalidSupercubeStatusError, 39
method), 87 IO_SCANNING(hvl_ccb.dev.se_ilst2.ILS2T.RegAddr at-
attribute), 108
INITIALIZED(hvl_ccb.experiment_manager.ExperimentStatus IoScanningModeValueError, 113
attribute), 117
```

```

is_configdataclass                                is_inactive(hvl_ccb.dev.crylas.CryLasLaser.LaserStatus
    (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig   attribute), 59
    attribute), 9
is_configdataclass                                is_open(hvl_ccb.comm.labjack_ljm.LJMCommunication
    (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig attribute), 10
    attribute), 11
is_configdataclass                                is_open(hvl_ccb.comm.serial.SerialCommunication
    (hvl_ccb.comm.opc.OpcUaCommunicationConfig         attribute), 14
    attribute), 13
is_configdataclass                                is_ready(hvl_ccb.dev.crylas.CryLasLaser.LaserStatus
    (hvl_ccb.comm.serial.SerialCommunicationConfig   attribute), 59
    attribute), 16
is_configdataclass                                is_running()(hvl_ccb.experiment_manager.ExperimentManager
    (hvl_ccb.comm.visa.VisaCommunicationConfig        method), 117
    attribute), 19
is_configdataclass                                is_type_hint() (in module hvl_ccb.utils.typing),
    (hvl_ccb.dev.base.EmptyConfig      attribute), 54
is_configdataclass                                is_valid_scale_range_reversed_str()
    (hvl_ccb.dev.crylas.CryLasAttenuatorConfig       attribute), 56
is_configdataclass                                jerk_time(hvl_ccb.dev.newport.NewportSMC100PPConfig
    (hvl_ccb.dev.crylas.CryLasLaserConfig      attribute), 62
    attribute), 91
is_configdataclass                                JOG(hvl_ccb.dev.se_ils2t.ILS2T.Mode attribute), 107
    (hvl_ccb.dev.heinzinger.HeinzingerConfig       attribute), 69
is_configdataclass                                jog_run()(hvl_ccb.dev.se_ils2t.ILS2T method), 110
    (hvl_ccb.dev.mbw973.MBW973Config      attribute), 80
is_configdataclass                                jog_stop()(hvl_ccb.dev.se_ils2t.ILS2T method), 110
    (hvl_ccb.dev.newport.NewportSMC100PPConfig     attribute), 91
is_configdataclass                                JOGGING(hvl_ccb.dev.newport.NewportStates      attribute), 95
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig   attribute), 99
is_configdataclass                                JOGGING_DISABLE
    (hvl_ccb.dev.se_ils2t.ILS2TConfig      attribute), 111
is_configdataclass                                (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    (hvl_ccb.dev.supercube.base.SupercubeConfiguration attribute), 23
    attribute), 84
is_configdataclass                                JOGGING_READY
    (hvl_ccb.dev.supercube2015.base.SupercubeConfiguration attribute), 42
is_done() (hvl_ccb.dev.mbw973.MBW973 method), 79
is_error() (hvl_ccb.experiment_manager.ExperimentManager
    method), 117
is_finished() (hvl_ccb.experiment_manager.ExperimentManager
    method), 117
is_generic() (in module hvl_ccb.utils.typing), 115
is_in_range() (hvl_ccb.dev.se_ils2t.ILS2T.RegDatatype
    method), 109
is_inactive(hvl_ccb.dev.crylas.CryLasLaser.LaserStatus
    attribute), 59
is_open(hvl_ccb.comm.labjack_ljm.LJMCommunication
    attribute), 10
is_ready(hvl_ccb.dev.crylas.CryLasLaser.LaserStatus
    attribute), 59
is_running()(hvl_ccb.experiment_manager.ExperimentManager
    method), 117
is_type_hint() (in module hvl_ccb.utils.typing),
    115
is_valid_scale_range_reversed_str()
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig.Model
    method), 98

```

## J

```

J (hvl_ccb.dev.labjack.LabJack.ThermocoupleType      attribute), 75
jerk_time(hvl_ccb.dev.newport.NewportSMC100PPConfig
    attribute), 91
JOG(hvl_ccb.dev.se_ils2t.ILS2T.Mode attribute), 107
jog_run()(hvl_ccb.dev.se_ils2t.ILS2T method), 110
jog_stop()(hvl_ccb.dev.se_ils2t.ILS2T method), 110
JOGGING(hvl_ccb.dev.newport.NewportStates      attribute), 95
JOGGING_DISABLE
    (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    attribute), 84
JOGGING_READY
    (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    attribute), 84
JOGN_FAST
    (hvl_ccb.dev.se_ils2t.ILS2T.RegAddr
    attribute), 108
JOGN_SLOW
    (hvl_ccb.dev.se_ils2t.ILS2T.RegAddr
    attribute), 108
JR(hvl_ccb.dev.newport.NewportConfigCommands      attribute), 82

```

## K

```

K (hvl_ccb.dev.labjack.LabJack.TemperatureUnit      attribute), 75
K (hvl_ccb.dev.labjack.LabJack.ThermocoupleType      attribute), 75
keys() (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig
    class method), 9
keys() (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig
    class method), 11
keys() (hvl_ccb.comm.opc.OpcUaCommunicationConfig
    class method), 13
keys() (hvl_ccb.comm.serial.SerialCommunicationConfig
    class method), 16

```

keys () (hvl\_ccb.comm.visa.VisaCommunicationConfig class method), 19  
 keys () (hvl\_ccb.dev.base.EmptyConfig class method), 54  
 keys () (hvl\_ccb.dev.crylas.CryLasAttenuatorConfig class method), 56  
 keys () (hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig class method), 57  
 keys () (hvl\_ccb.dev.crylas.CryLasLaserConfig class method), 62  
 keys () (hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig class method), 63  
 keys () (hvl\_ccb.dev.ea\_psi9000.PSI9000Config class method), 67  
 keys () (hvl\_ccb.dev.ea\_psi9000.PSI9000VisaCommunicationConfig class method), 68  
 keys () (hvl\_ccb.dev.heinzinger.HeinzingerConfig class method), 69  
 keys () (hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig class method), 74  
 keys () (hvl\_ccb.dev.mbw973.MBW973Config class method), 80  
 keys () (hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig class method), 81  
 keys () (hvl\_ccb.dev.newport.NewportSMC100PPConfig class method), 91  
 keys () (hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig class method), 95  
 keys () (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig class method), 99  
 keys () (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig class method), 101  
 keys () (hvl\_ccb.dev.rs\_rto1024.RTO1024Config class method), 106  
 keys () (hvl\_ccb.dev.rs\_rto1024.RTO1024VisaCommunicationConfig class method), 107  
 keys () (hvl\_ccb.dev.se\_ils2t.ILS2TConfig class method), 111  
 keys () (hvl\_ccb.dev.se\_ils2t.ILS2TModbusTcpCommunicationConfig class method), 112  
 keys () (hvl\_ccb.dev.supercube.base.SupercubeConfiguration class method), 23  
 keys () (hvl\_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig class method), 24  
 keys () (hvl\_ccb.dev.supercube.typ\_a.SupercubeAOpcUaConfiguration class method), 36  
 keys () (hvl\_ccb.dev.supercube.typ\_b.SupercubeBOpcUaConfiguration class method), 39  
 keys () (hvl\_ccb.dev.supercube2015.base.SupercubeConfiguration class method), 42  
 keys () (hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig class method), 43  
 keys () (hvl\_ccb.dev.supercube2015.typ\_a.SupercubeAOpcUaConfiguration class method), 52

keys () (hvl\_ccb.dev.visa.VisaDeviceConfig class method), 114  
 kV (hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitVoltage attribute), 71

**L**

LabJack.AInRange (class in hvl\_ccb.dev.labjack), 74  
 LabJack.CalMicroAmpere (class in hvl\_ccb.dev.labjack), 74  
 LabJack.Cjctype (class in hvl\_ccb.dev.labjack), 75  
 LabJack.DeviceType (class in hvl\_ccb.dev.labjack), 75  
 LabJack.DIOStatus (class in hvl\_ccb.dev.labjack), 75  
 LabJack.TemperatureUnit (class in hvl\_ccb.dev.labjack), 75  
 LabJack.ThermocoupleType (class in hvl\_ccb.dev.labjack), 75  
 LabJackError, 78  
 LabJackIdentifierDIOError, 78  
 laser\_off () (hvl\_ccb.dev.crylas.CryLasLaser method), 59  
 laser\_on () (hvl\_ccb.dev.crylas.CryLasLaser method), 59  
 list\_directory () (hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 102  
 LJMCommunication (class in hvl\_ccb.comm.labjack\_ljm), 8  
 LJMCommunicationConfig (class in hvl\_ccb.comm.labjack\_ljm), 8  
 LJMCommunicationConfig.ConnectionType (class in hvl\_ccb.comm.labjack\_ljm), 9  
 LJMCommunicationConfig.DeviceType (class in hvl\_ccb.comm.labjack\_ljm), 9  
 LJMCommunicationError, 10  
 lm34 (hvl\_ccb.dev.labjack.LabJack.Cjctype attribute), 75  
 RTO1024\_Configuration () (hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 102  
 local\_display () (hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 102  
 locked (hvl\_ccb.dev.supercube.constants.DoorStatus attribute), 31  
 locked (hvl\_ccb.dev.supercube2015.constants.DoorStatus attribute), 45  
 LOW (hvl\_ccb.dev.labjack.LabJack.DIOStatus attribute), 75  
 mA (hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitCurrent attribute), 71

```

manual (hvl_ccb.dev.supercube.constants.EarthingStick         method), 97
       attribute), 31                                         measure_current ()
manual (hvl_ccb.dev.supercube2015.constants.EarthingStick    (hvl_ccb.dev.heinzinger.HeinzingerDI method),
       attribute), 45                                         70
manual_1 (hvl_ccb.dev.supercube.constants.EarthingStickmeasure_voltage ()
       attribute), 31                                         (hvl_ccb.dev.heinzinger.HeinzingerDI method),
manual_1 (hvl_ccb.dev.supercube2015.constants.EarthingStick  70
       attribute), 45                                         measure_voltage_current ()
manual_2 (hvl_ccb.dev.supercube.constants.EarthingStick      (hvl_ccb.dev.ea_psi9000.PSI9000   method),
       attribute), 31                                         65
manual_2 (hvl_ccb.dev.supercube2015.constants.EarthingStickMeasurementsDividerRatio (class     in
       attribute), 46                                         hvl_ccb.dev.supercube.constants), 33
manual_3 (hvl_ccb.dev.supercube.constants.EarthingStickMeasurementsDividerRatio (class     in
       attribute), 31                                         hvl_ccb.dev.supercube2015.constants), 48
manual_3 (hvl_ccb.dev.supercube2015.constants.EarthingStickMeasurementsScaledInput (class     in
       attribute), 46                                         hvl_ccb.dev.supercube.constants), 33
manual_4 (hvl_ccb.dev.supercube.constants.EarthingStickMeasurementsScaledInput (class     in
       attribute), 31                                         hvl_ccb.dev.supercube2015.constants), 48
manual_4 (hvl_ccb.dev.supercube2015.constants.EarthingStickMessage (hvl_ccb.dev.supercube.constants.Errors at-
       attribute), 46                                         tribute), 32
manual_5 (hvl_ccb.dev.supercube.constants.EarthingStickmicro_step_per_full_step_factor
       attribute), 31                                         (hvl_ccb.dev.newport.NewportSMC100PPConfig
manual_5 (hvl_ccb.dev.supercube2015.constants.EarthingStickattribute), 91
       attribute), 46                                         Micron (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG.PressureUnits
manual_6 (hvl_ccb.dev.supercube.constants.EarthingStick        attribute), 96
       attribute), 31                                         ModbusTcpCommunication (class     in
manual_6 (hvl_ccb.dev.supercube2015.constants.EarthingStick    hvl_ccb.comm.modbus_tcp), 10
       attribute), 46                                         ModbusTcpCommunicationConfig (class     in
                                                               hvl_ccb.comm.modbus_tcp), 11
MARK (hvl_ccb.comm.serial.SerialCommunicationParity        ModbusTcpConnectionFailedException, 11
       attribute), 16
max_current (hvl_ccb.dev.heinzinger.HeinzingerPNC           model (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig at-
       attribute), 72                                         tribute), 99
max_current.hardware (hvl_ccb.dev.heinzinger.HeinzingerPNC motion_distance_per_full_step
       attribute), 72                                         (hvl_ccb.dev.newport.NewportSMC100PPConfig
                                                               attribute), 91
max_voltage (hvl_ccb.dev.heinzinger.HeinzingerPNC           motor_config (hvl_ccb.dev.newport.NewportSMC100PPConfig
       attribute), 72                                         attribute), 91
max_voltage.hardware (hvl_ccb.dev.heinzinger.HeinzingerPNC move_finished_extra_wait_sec
       attribute), 72                                         (hvl_ccb.dev.newport.NewportSMC100PPConfig
                                                               attribute), 91
mbar (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG.PressureUnits move_to_absolute_position ()
       attribute), 96                                         (hvl_ccb.dev.newport.NewportSMC100PP
                                                               method), 87
MBW973 (class in hvl_ccb.dev.mbw973), 79
MBW973Config (class in hvl_ccb.dev.mbw973), 80
MBW973ControlRunningException, 80
MBW973Error, 80
MBW973PumpRunningException, 81
MBW973SerialCommunication (class     in
       hvl_ccb.dev.mbw973), 81
MBW973SerialCommunicationConfig (class     in
       hvl_ccb.dev.mbw973), 81
measure () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG
       method), 97
measure_all () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG
               MS_NOMINAL_CURRENT
               (hvl_ccb.dev.ea_psi9000.PSI9000 attribute), 64
               MS_NOMINAL_VOLTAGE
               (hvl_ccb.dev.ea_psi9000.PSI9000 attribute), 64

```

MULTI\_COMMANDS\_MAX  
`(hvl_ccb.comm.visa.VisaCommunication attribute), 17`

MULTI\_COMMANDS\_SEPARATOR  
`(hvl_ccb.comm.visa.VisaCommunication attribute), 17`

**N**

NameEnum (*class in hvl\_ccb.utils.enum*), 115

NAMES (*hvl\_ccb.comm.serial.SerialCommunicationParity attribute*), 16

names (*hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes attribute*), 102

namespace\_index (*hvl\_ccb.dev.supercube.base.Supercube attribute*), 23

namespace\_index (*hvl\_ccb.dev.supercube2015.base.SupercubeConfiguration attribute*), 42

NED\_END\_OF\_TURN (*hvl\_ccb.dev.newport.NewportSMC100PP.MotorError attribute*), 83

NEG (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 108

NEG\_FAST (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 108

negative\_software\_limit  
`(hvl_ccb.dev.newport.NewportSMC100PPConfig attribute), 91`

NewportConfigCommands (*class in hvl\_ccb.dev.newport*), 82

NewportControllerError, 83

NewportMotorError, 83

NewportSerialCommunicationError, 95

NewportSMC100PP (*class in hvl\_ccb.dev.newport*), 83

NewportSMC100PP.MotorErrors (*class in hvl\_ccb.dev.newport*), 83

NewportSMC100PP.StateMessages (*class in hvl\_ccb.dev.newport*), 83

NewportSMC100PPConfig (*class in hvl\_ccb.dev.newport*), 89

NewportSMC100PPConfig.EspStageConfig  
`(class in hvl_ccb.dev.newport), 90`

NewportSMC100PPConfig.HomeSearch (*class in hvl\_ccb.dev.newport*), 90

NewportSMC100PPSerialCommunication (*class in hvl\_ccb.dev.newport*), 92

NewportSMC100PPSerialCommunication.ControllerError  
`(class in hvl_ccb.dev.newport), 92`

NewportSMC100PPSerialCommunicationConfig  
`(class in hvl_ccb.dev.newport), 94`

NewportStates (*class in hvl\_ccb.dev.newport*), 95

NO\_ERROR (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication attribute*), 92

NO\_REF (*hvl\_ccb.dev.newport.NewportStates attribute*), 95

NO\_REF\_ESP\_STAGE\_ERROR  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 17`

attribute), 84

NO\_REF\_FROM\_CONFIG  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84`

NO\_REF\_FROM\_DISABLED  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84`

NO\_REF\_FROM\_HOMING  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84`

NO\_REF\_FROM\_JOGGING  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84`

NO\_REF\_FROM\_MOVING  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84`

NO\_REF\_FROM\_READY  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84`

NO\_REF\_FROM\_RESET  
`(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 84`

No\_sensor (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute*), 96

NONE (*hvl\_ccb.comm.serial.SerialCommunicationParity attribute*), 17

NONE (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 76

None (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute*), 96

NONE (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 108

NoPower (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 35

NORMAL (*hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes attribute*), 101

noSen (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute*), 97

noSENSOR (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute*), 97

not\_defined (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 26

not\_defined (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 45

number\_of\_decimals

    (*hvl\_ccb.dev.heinzinger.HeinzingerConfig attribute*), 69

    number\_of\_sensors

        (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG attribute*), 98

O

ODD (*hvl\_ccb.comm.serial.SerialCommunicationParity attribute*), 17

```

OH (hvl_ccb.dev.newport.NewportConfigCommands attribute), 83
Ok (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG.SensorStatus attribute), 96
ON (hvl_ccb.dev.se_ilis2t.ILS2T.State attribute), 109
on_start_wait_until_ready (hvl_ccb.dev.crylas.CryLasLaserConfig attribute), 62
ONE (hvl_ccb.comm.serial.SerialCommunicationStopbits attribute), 17
ONE (hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum attribute), 69
ONE (hvl_ccb.dev.labjack.LabJack.AInRange attribute), 74
ONE_HUNDREDTH (hvl_ccb.dev.labjack.LabJack.AInRange attribute), 74
ONE_POINT_FIVE (hvl_ccb.comm.serial.SerialCommunicationStop attribute), 17
ONE_TENTH (hvl_ccb.dev.labjack.LabJack.AInRange attribute), 74
OpcControl (class in hvl_ccb.dev.supercube.constants), 33
OpcUaCommunication (class in hvl_ccb.comm.opc), 11
OpcUaCommunicationConfig (class in hvl_ccb.comm.opc), 12
OpcUaCommunicationIOError, 13
OpcUaSubHandler (class in hvl_ccb.comm.opc), 13
open (hvl_ccb.dev.supercube.constants.DoorStatus attribute), 31
open (hvl_ccb.dev.supercube.constants.EarthingStickStatus attribute), 32
open (hvl_ccb.dev.supercube2015.constants.DoorStatus attribute), 45
open (hvl_ccb.dev.supercube2015.constants.EarthingStickStatus attribute), 46
open () (hvl_ccb.comm.base.CommunicationProtocol method), 7
open () (hvl_ccb.comm.labjack_ljm.LJMCommunication method), 8
open () (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication method), 10
open () (hvl_ccb.comm.opc.OpcUaCommunication method), 12
open () (hvl_ccb.comm.serial.SerialCommunication method), 14
open () (hvl_ccb.comm.visa.VisaCommunication method), 17
open_shutter() (hvl_ccb.dev.crylas.CryLasLaser method), 60
open_timeout (hvl_ccb.comm.visa.VisaCommunicationConfig attribute), 19
OPENED (hvl_ccb.dev.crylas.CryLasLaser.AnswersShutter attribute), 58
OPENED (hvl_ccb.dev.crylas.CryLasLaserStatus attribute), 64
operate () (hvl_ccb.dev.supercube.base.SupercubeBase method), 21
operate () (hvl_ccb.dev.supercube2015.base.Supercube2015Base method), 41
optional_defaults () (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig class method), 9
optional_defaults () (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig class method), 11
optional_defaults () (hvl_ccb.comm.opc.OpcUaCommunicationConfig class method), 13
optional_defaults () (hvl_ccb.comm.visa.VisaCommunicationConfig class method), 16
optional_defaults () (hvl_ccb.comm.visa.VisaCommunicationConfig class method), 19
optional_defaults () (hvl_ccb.dev.base.EmptyConfig class method), 54
optional_defaults () (hvl_ccb.dev.crylas.CryLasAttenuatorConfig class method), 56
optional_defaults () (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig class method), 58
optional_defaults () (hvl_ccb.dev.crylas.CryLasLaserConfig class method), 62
optional_defaults () (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig class method), 64
optional_defaults () (hvl_ccb.dev.ea_psi9000.PSI9000Config class method), 67
optional_defaults () (hvl_ccb.dev.heinzinger.HeinzingerConfig class method), 69
optional_defaults () (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig class method), 74
optional_defaults () (hvl_ccb.dev.mbw973.MBW973Config class method), 80
optional_defaults () (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig class method), 81

```

```

optional_defaults()
    (hvl_ccb.dev.newport.NewportSMC100PPConfig
     class method), 91
optional_defaults()
    (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig
     class method), 95
optional_defaults()
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig
     class method), 99
optional_defaults()
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig
     class method), 101
optional_defaults()
    (hvl_ccb.dev.rs_rto1024.RTO1024Config
     class method), 106
optional_defaults()
    (hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig
     class method), 107
optional_defaults()
    (hvl_ccb.dev.se_ils2t.ILS2TConfig
     method), 111
optional_defaults()
    (hvl_ccb.dev.se_ils2t.ILS2TModbusTcpCommunicationConfig
     class method), 112
optional_defaults()
    (hvl_ccb.dev.supercube.base.SupercubeConfiguration
     class method), 23
optional_defaults()
    (hvl_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig
     class method), 24
optional_defaults()
    (hvl_ccb.dev.supercube.typ_a.SupercubeAOpcUaConfiguration
     class method), 36
optional_defaults()
    (hvl_ccb.dev.supercube.typ_b.SupercubeBOpcUaConfiguration
     class method), 39
optional_defaults()
    (hvl_ccb.dev.supercube2015.base.SupercubeConfiguration
     class method), 42
optional_defaults()
    (hvl_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig
     class method), 43
optional_defaults()
    (hvl_ccb.dev.supercube2015.typ_a.SupercubeAOpcUaConfiguration
     class method), 52
optional_defaults()
    (hvl_ccb.dev.visa.VisaDeviceConfig
     method), 114
OT (hvl_ccb.dev.newport.NewportConfigCommands
     attribute), 83
out_1_1 (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
out_1_1 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 47
out_1_2 (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
out_1_2 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 47
out_2_1 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 47
out_2_2 (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
out_3_1 (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
out_3_1 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 48
out_3_2 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 48
out_4_1 (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
out_4_2 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 33
out_4_2 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 48
out_5_1 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 48
out_5_2 (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
out_6_1 (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
out_6_1 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 48
out_6_2 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 33
out_6_2 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 48
output (hvl_ccb.dev.supercube.constants.GeneralSupport
     attribute), 33
output (hvl_ccb.dev.supercube2015.constants.GeneralSupport
     attribute), 48
output (hvl_ccb.dev.heinzinger.HeinzingerDI
     method), 70
output_on () (hvl_ccb.dev.heinzinger.HeinzingerDI
     method), 71
OUT_POWER_EXCEEDED
(hvl_ccb.dev.newport.NewportSMC100PP.MotorErrors
     attribute), 47

```

*attribute), 83*  
*Overrange (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute), 96*  
**P**  
*PARAM\_MISSING\_OR\_INVALID (hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig attribute), 92*  
*Parity (hvl\_ccb.comm.serial.SerialCommunicationConfig attribute), 15*  
*parity (hvl\_ccb.comm.serial.SerialCommunicationConfig attribute), 16*  
*parity (hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig attribute), 58*  
*parity (hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig attribute), 64*  
*parity (hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig attribute), 74*  
*parity (hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig attribute), 82*  
*parity (hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig attribute), 95*  
*parity (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig attribute), 101*  
*Pascal (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits attribute), 96*  
*PBR (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute), 97*  
*PEAK\_CURRENT\_LIMIT (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 83*  
*peak\_output\_current\_limit (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 91*  
*PfeifferTPG (class in hvl\_ccb.dev.pfeiffer\_tpg), 96*  
*PfeifferTPG.PressureUnits (class in hvl\_ccb.dev.pfeiffer\_tpg), 96*  
*PfeifferTPG.SensorStatus (class in hvl\_ccb.dev.pfeiffer\_tpg), 96*  
*PfeifferTPG.SensorTypes (class in hvl\_ccb.dev.pfeiffer\_tpg), 96*  
*PfeifferTPGConfig (class in hvl\_ccb.dev.pfeiffer\_tpg), 98*  
*PfeifferTPGConfig.Model (class in hvl\_ccb.dev.pfeiffer\_tpg), 98*  
*PfeifferTPGError, 99*  
*PfeifferTPGSerialCommunication (class in hvl\_ccb.dev.pfeiffer\_tpg), 99*  
*PfeifferTPGSerialCommunicationConfig (class in hvl\_ccb.dev.pfeiffer\_tpg), 100*  
*PKR (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute), 97*  
*Poller (class in hvl\_ccb.dev.mbw973), 82*  
*polling\_interval (hvl\_ccb.dev.mbw973.MBW973Config attribute), 80*  
*polling\_period (hvl\_ccb.dev.crylas.CryLasLaserConfig attribute), 62*  
*polling\_timeout (hvl\_ccb.dev.crylas.CryLasLaserConfig attribute), 62*  
*port (hvl\_ccb.comm.opc.OpcUaCommunicationConfig attribute), 13*  
*port (hvl\_ccb.comm.serial.SerialCommunicationConfig attribute), 16*  
*port (hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig attribute), 19*  
*port (hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig attribute), 44*  
*POS\_OF\_Coil (hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute), 108*  
*POS\_END\_OF\_TURN (hvl\_ccb.dev.newport.NewportSMC100PP.MotorError attribute), 108*  
*POSITION (hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr attribute), 108*  
*POSITION\_OUT\_OF\_LIMIT (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 91*  
*positive\_software\_limit (hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute), 91*  
*POSITION\_OF\_FORCE\_VALUE () (hvl\_ccb.dev.newport.NewportSMC100PPConfig method), 91*  
*Power (class in hvl\_ccb.dev.supercube.constants), 34*  
*Power (class in hvl\_ccb.dev.supercube2015.constants), 48*  
*power\_limit (hvl\_ccb.dev.ea\_psi9000.PSI9000Config attribute), 67*  
*PowerSetup (class in hvl\_ccb.dev.supercube.constants), 34*  
*PowerSetup (class in hvl\_ccb.dev.supercube2015.constants), 49*  
*prepare\_ultra\_segmentation () (hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 102*  
*PSI9000 (class in hvl\_ccb.dev.ea\_psi9000), 64*  
*PSI9000Config (class in hvl\_ccb.dev.ea\_psi9000), 66*  
*PSI9000Error, 67*  
*PSI9000VisaCommunication (class in hvl\_ccb.dev.ea\_psi9000), 67*  
*PSI9000VisaCommunicationConfig (class in hvl\_ccb.dev.ea\_psi9000), 67*  
*PT100 (hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute), 76*  
*PT1000 (hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute), 76*

*attribute), 76*  
**PT500** (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 76  
**PTP** (*hvl\_ccb.dev.se\_ils2t.ILS2T.Mode attribute*), 107

**Q**

**QIL** (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 83  
**query()** (*hvl\_ccb.comm.visa.VisaCommunication method*), 17  
**query()** (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunication method*), 62  
**query()** (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication method*), 93  
**query()** (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunication method*), 99  
**query\_all()** (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunication method*), 63  
**query\_multiple()** (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication method*), 93  
**QUICKSTOP** (*hvl\_ccb.dev.se\_ils2t.ILS2T.State attribute*), 109  
**QuickStop** (*hvl\_ccb.dev.supercube.constants.SafetyStatus attribute*), 35  
**QuickStop** (*hvl\_ccb.dev.supercube2015.constants.SafetyStatus attribute*), 50  
**quickstop()** (*hvl\_ccb.dev.se\_ils2t.ILS2T method*), 110  
**quit** (*hvl\_ccb.dev.supercube.constants.Errors attribute*), 32  
**quit** (*hvl\_ccb.dev.supercube2015.constants.Errors attribute*), 47  
**quit\_error()** (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 21  
**quit\_error()** (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 41

**R**

**R** (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 76  
**RAMP\_ACC** (*hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr attribute*), 108  
**RAMP\_DECEL** (*hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr attribute*), 108  
**RAMP\_N\_MAX** (*hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr attribute*), 108  
**RAMP\_TYPE** (*hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr attribute*), 108  
**read()** (*hvl\_ccb.comm.opc.OpcUaCommunication method*), 12  
**read()** (*hvl\_ccb.dev.mbw973.MBW973 method*), 79  
**read()** (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 21

**read()** (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 41  
**READ\_ACTIVE** (*hvl\_ccb.dev.crylas.CryLasLaser.LaserStatus attribute*), 58  
**read\_bytes()** (*hvl\_ccb.comm.serial.SerialCommunication method*), 14  
**read\_float()** (*hvl\_ccb.dev.mbw973.MBW973 method*), 79  
**read\_holding\_registers()** (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication method*), 10  
**read\_input\_registers()** (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication method*), 10  
**read\_int()** (*hvl\_ccb.dev.mbw973.MBW973 method*), 79  
**read\_measurements()** (*hvl\_ccb.dev.mbw973.MBW973 method*), 79  
**read\_name()** (*hvl\_ccb.comm.labjack\_ljm.LJMCommunication method*), 8  
**read\_resistance()** (*hvl\_ccb.dev.labjack.LabJack method*), 77  
**read\_termination** (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 19  
**read\_text()** (*hvl\_ccb.comm.serial.SerialCommunication method*), 14  
**read\_text\_nonempty()** (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunication method*), 72  
**READ\_TEXT\_SKIP\_PREFIXES** (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunication attribute*), 62  
**read\_thermocouple()** (*hvl\_ccb.dev.labjack.LabJack method*), 77  
**READY** (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute*), 58  
**READY** (*hvl\_ccb.dev.newport.NewportStates attribute*), 95  
**READY** (*hvl\_ccb.dev.se\_ils2t.ILS2T.State attribute*), 109  
**ready()** (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 21  
**ready()** (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 41  
**READY\_FROM\_DISABLE** (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 84  
**READY\_FROM\_HOMING** (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 84  
**READY\_FROM\_JOGGING** (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 84  
**READY\_FROM\_MOVING**

(*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessage* required\_keys () (*hvl\_ccb.dev.mbw973.MBW973SerialCommunication*  
*attribute*), 84  
*class method*), 82

READY\_INACTIVE (*hvl\_ccb.dev.crylas.CryLasLaser.LaserSetup* required\_keys () (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
*attribute*), 58  
*class method*), 92

RedOperate (*hvl\_ccb.dev.supercube.constants.SafetyStatus* required\_keys () (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCom-*  
*attribute*), 35  
*class method*), 95

RedOperate (*hvl\_ccb.dev.supercube2015.constants.SafetyStatus* required\_keys () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig*  
*attribute*), 50  
*class method*), 99

RedReady (*hvl\_ccb.dev.supercube.constants.SafetyStatus* required\_keys () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommuni-*  
*attribute*), 35  
*class method*), 101

RedReady (*hvl\_ccb.dev.supercube2015.constants.SafetyStatus* required\_keys () (*hvl\_ccb.dev.rs\_rto1024.RTO1024Config*  
*attribute*), 50  
*class method*), 106

RELATIVE\_POSITION\_MOTOR  
(*hvl\_ccb.dev.se\_ils2t.ILS2T.ActionsPtp* at-  
*tribute*), 107  
*required\_keys* () (*hvl\_ccb.dev.rs\_rto1024.RTO1024VisaCommunicati-*  
*class method*), 107

RELATIVE\_POSITION\_TARGET  
(*hvl\_ccb.dev.se\_ils2t.ILS2T.ActionsPtp* at-  
*tribute*), 107  
*required\_keys* () (*hvl\_ccb.dev.se\_ils2t.ILS2TConfig*  
*class method*), 111

relative\_step () (*hvl\_ccb.dev.se\_ils2t.ILS2T* required\_keys () (*hvl\_ccb.dev.supercube.base.SupercubeConfigura-*  
*method*), 110  
*class method*), 23

relative\_step\_and\_wait () required\_keys () (*hvl\_ccb.dev.supercube.base.SupercubeOpcUaCommuni-*  
*(hvl\_ccb.dev.se\_ils2t.ILS2T method)*, 110  
*class method*), 24

remove\_device () (*hvl\_ccb.dev.base.DeviceSequenceMix* required\_keys () (*hvl\_ccb.dev.supercube.typ\_a.SupercubeAOpcUaCom-*  
*method*), 53  
*class method*), 36

required\_keys () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfigs* () (*hvl\_ccb.dev.supercube.typ\_b.SupercubeBOpcUaCom-*  
*class method*), 10  
*class method*), 39

required\_keys () (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfigs* () (*hvl\_ccb.dev.supercube2015.base.SupercubeConfigura-*  
*class method*), 11  
*class method*), 43

required\_keys () (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig* keys () (*hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaCom-*  
*class method*), 13  
*class method*), 44

required\_keys () (*hvl\_ccb.comm.serial.SerialCommunicationConfig* keys () (*hvl\_ccb.dev.supercube2015.typ\_a.SupercubeAOpcUaCom-*  
*class method*), 16  
*class method*), 52

required\_keys () (*hvl\_ccb.comm.visa.VisaCommunicationConfig* keys () (*hvl\_ccb.dev.visa.VisaDeviceConfig*  
*class method*), 19  
*class method*), 114

required\_keys () (*hvl\_ccb.dev.base.EmptyConfig* reset (*hvl\_ccb.dev.supercube.constants.BreakdownDetection*  
*attribute*), 54  
*class method*), 30

required\_keys () (*hvl\_ccb.dev.crylas.CryLasAttenuatorConfig* (*hvl\_ccb.dev.supercube2015.constants.BreakdownDetection*  
*class method*), 56  
*attribute*), 45

required\_keys () (*hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig* (*hvl\_ccb.dev.heinzinger.NewportSMC100PP*  
*method*), 58  
*class method*), 87

required\_keys () (*hvl\_ccb.dev.crylas.CryLasLaserConfig* reset () (*hvl\_ccb.dev.visa.VisaDevice method*), 62  
*class method*), 113  
*reset\_error* () (*hvl\_ccb.dev.se\_ils2t.ILS2T* method), 71

required\_keys () (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig* reset\_interface ()  
*class method*), 64

required\_keys () (*hvl\_ccb.dev.ea\_psi9000.PSI9000Config* (*hvl\_ccb.dev.heinzinger.HeinzingerDI method*),  
*class method*), 67  
*71*

required\_keys () (*hvl\_ccb.dev.ea\_psi9000.PSI9000VisaCommunicationConfig* time  
*class method*), 68  
*(hvl\_ccb.dev.crylas.CryLasAttenuatorConfig* attribute), 56

required\_keys () (*hvl\_ccb.dev.heinzinger.HeinzingerConfig* RMS\_CURRENT\_LIMIT  
*class method*), 69

required\_keys () (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig* (*hvl\_ccb.dev.heinzinger.NewportSMC100PP.MotorErrors*  
*attribute*), 74  
*attribute*), 83

required\_keys () (*hvl\_ccb.dev.mbw973.MBW973Config* pm\_max\_init (*hvl\_ccb.dev.se\_ils2t.ILS2TConfig* at-  
*class method*), 80  
*tribute*), 112

rs485\_address (*hvl\_ccb.dev.newport.NewportSMC100PPConfig\_error* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute*), 92  
 RTO1024 (*class in hvl\_ccb.dev.rs\_rto1024*), 101  
 RTO1024.TriggerModes (*class in hvl\_ccb.dev.rs\_rto1024*), 101  
 RTO1024Config (*class in hvl\_ccb.dev.rs\_rto1024*), 105  
 RTO1024Error, 106  
 RTO1024VisaCommunication (*class in hvl\_ccb.dev.rs\_rto1024*), 106  
 RTO1024VisaCommunicationConfig (*class in hvl\_ccb.dev.rs\_rto1024*), 106  
 run () (*hvl\_ccb.dev.visa.VisaStatusPoller method*), 114  
 run () (*hvl\_ccb.experiment\_manager.ExperimentManager method*), 117  
 run\_continuous\_acquisition ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 102  
 run\_single\_acquisition ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 103  
 RUNNING (*hvl\_ccb.experiment\_manager.ExperimentStatus attribute*), 117

**S**

S (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 76  
 SA (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 83  
 Safety (*class in hvl\_ccb.dev.supercube.constants*), 35  
 Safety (*class in hvl\_ccb.dev.supercube2015.constants*), 49  
 SafetyStatus (*class in hvl\_ccb.dev.supercube.constants*), 35  
 SafetyStatus (*class in hvl\_ccb.dev.supercube2015.constants*), 49  
 save\_configuration ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 103  
 save\_waveform\_history ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 103  
 SCALE (*hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr attribute*), 108  
 ScalingFactorValueError, 113  
 screw\_scaling (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 92  
 send\_command () (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication* (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*)), 93  
 send\_command () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunication* (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*)), 99  
 send\_stop () (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication* (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*)), 93

PPGConfig\_error (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute*), 96  
 Sensor\_off (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute*), 96  
 SerialCommunication (*class in hvl\_ccb.comm.serial*), 14  
 SerialCommunicationBytesize (*class in hvl\_ccb.comm.serial*), 15  
 SerialCommunicationConfig (*class in hvl\_ccb.comm.serial*), 15  
 SerialCommunicationIOError, 16  
 SerialCommunicationParity (*class in hvl\_ccb.comm.serial*), 16  
 SerialCommunicationStopbits (*class in hvl\_ccb.comm.serial*), 17  
 set\_acceleration ()  
     (*hvl\_ccb.dev.newport.NewportSMC100PP method*), 87  
 set\_acquire\_length ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 103  
 set\_ain\_differential ()  
     (*hvl\_ccb.dev.labjack.LabJack method*), 77  
 set\_ain\_range ()  
     (*hvl\_ccb.dev.labjack.LabJack method*), 77  
 set\_ain\_resistance ()  
     (*hvl\_ccb.dev.labjack.LabJack method*), 77  
 set\_ain\_resolution ()  
     (*hvl\_ccb.dev.labjack.LabJack method*), 78  
 set\_ain\_thermocouple ()  
     (*hvl\_ccb.dev.labjack.LabJack method*), 78  
 set\_attenuation ()  
     (*hvl\_ccb.dev.crylas.CryLasAttenuator method*), 55  
 set\_cee16\_socket ()  
     (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 21  
 set\_cee16\_socket ()  
     (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 41  
 set\_channel\_position ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 103  
 set\_channel\_range ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 104  
 set\_channel\_scale ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 104  
 set\_current () (*hvl\_ccb.dev.heinzinger.HeinzingerDI method*), 71

```

set_current () (hvl_ccb.dev.heinzinger.HeinzingerPNC
    method), 72
set_digital_output ()          (hvl_ccb.dev.labjack.LabJack method), 78
set_display_unit ()           (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method),
    98
set_earthing_manual ()        (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 22
set_earthing_manual ()        (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 41
set_full_scale_mbar ()         (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method),
    98
set_full_scale_unitless ()     (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method),
    98
set_init_attenuation ()        (hvl_ccb.dev.crylas.CryLasAttenuator method),
    55
set_init_shutter_status ()     (hvl_ccb.dev.crylas.CryLasLaser      method),
    60
set_jog_speed ()              (hvl_ccb.dev.se_ilS2t.ILS2T
    method), 110
set_lower_limits ()            (hvl_ccb.dev.ea_psi9000.PSI9000      method),
    65
set_max_acceleration ()        (hvl_ccb.dev.se_ilS2t.ILS2T method), 110
set_max_deceleration ()        (hvl_ccb.dev.se_ilS2t.ILS2T method), 110
set_max_rpm ()                (hvl_ccb.dev.se_ilS2t.ILS2T method),
    110
set_measuring_options ()       (hvl_ccb.dev.mbw973.MBW973      method),
    79
set_motor_configuration ()     (hvl_ccb.dev.newport.NewportSMC100PP
    method), 88
set_negative_software_limit ()  (hvl_ccb.dev.newport.NewportSMC100PP
    method), 88
set_number_of_recordings ()    (hvl_ccb.dev.heinzinger.HeinzingerDI method),
    71
set_output ()                  (hvl_ccb.dev.ea_psi9000.PSI9000
    method), 65
set_positive_software_limit ()  (hvl_ccb.dev.newport.NewportSMC100PP
    method), 88
set_pulse_energy ()            (hvl_ccb.dev.crylas.CryLasLaser      method),
    60
set_ramp_type ()               (hvl_ccb.dev.se_ilS2t.ILS2T
    method), 111
set_reference_point ()         (hvl_ccb.dev.rs_rto1024.RTO1024      method),
    104
set_remote_control ()          (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 22
set_remote_control ()          (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 41
set_repetition_rate ()         (hvl_ccb.dev.crylas.CryLasLaser      method),
    60
set_repetitions ()             (hvl_ccb.dev.rs_rto1024.RTO1024      method),
    104
set_slope ()                   (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 37
set_slope ()                   (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 51
set_support_output ()           (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 22
set_support_output ()           (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 41
set_support_output_impulse ()   (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 22
set_support_output_impulse ()   (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 41
set_system_lock ()              (hvl_ccb.dev.ea_psi9000.PSI9000      method),
    65
set_t13_socket ()              (hvl_ccb.dev.supercube.base.SupercubeBase
    method), 22
set_t13_socket ()              (hvl_ccb.dev.supercube2015.base.Supercube2015Base
    method), 42
set_target_voltage ()            (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
    method), 38
set_target_voltage ()            (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
    method), 51
set_transmission ()              (hvl_ccb.dev.crylas.CryLasAttenuator method),
    55
set_trigger_level ()             (hvl_ccb.dev.rs_rto1024.RTO1024      method),
    105

```

set\_trigger\_mode ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024* method), 105

set\_trigger\_source ()  
     (*hvl\_ccb.dev.rs\_rto1024.RTO1024* method), 105

set\_upper\_limits ()  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000* method), 65

set\_voltage () (*hvl\_ccb.dev.heinzinger.HeinzingerDI* method), 71

set\_voltage () (*hvl\_ccb.dev.heinzinger.HeinzingerPNC* method), 72

set\_voltage\_current ()  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000* method), 66

setup (*hvl\_ccb.dev.supercube.constants.Power* attribute), 34

setup (*hvl\_ccb.dev.supercube2015.constants.Power* attribute), 48

SEVENBITS (*hvl\_ccb.comm.serial.SerialCommunicationBytes* attribute), 15

SHORT\_CIRCUIT (*hvl\_ccb.dev.newport.NewportSMC100PP.Motor* attribute), 88

SHUTDOWN\_CURRENT\_LIMIT  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000* attribute), 64

SHUTDOWN\_VOLTAGE\_LIMIT  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000* attribute), 64

ShutterStatus (*hvl\_ccb.dev.crylas.CryLasLaser* attribute), 59

ShutterStatus (*hvl\_ccb.dev.crylas.CryLasLaserConfig* attribute), 61

SingleCommDevice (class in *hvl\_ccb.dev.base*), 54

SIXBITS (*hvl\_ccb.comm.serial.SerialCommunicationBytes* attribute), 15

SIXTEEN (*hvl\_ccb.dev.heinzinger.HeinzingerConfig.Recording* attribute), 69

SL (*hvl\_ccb.dev.newport.NewportConfigCommands* attribute), 83

SOFTWARE\_INTERNAL\_SIXTY  
     (*hvl\_ccb.dev.crylas.CryLasLaser.RepetitionRates* attribute), 59

SOFTWARE\_INTERNAL\_TEN  
     (*hvl\_ccb.dev.crylas.CryLasLaser.RepetitionRates* attribute), 59

SOFTWARE\_INTERNAL\_TWENTY  
     (*hvl\_ccb.dev.crylas.CryLasLaser.RepetitionRates* attribute), 59

SPACE (*hvl\_ccb.comm.serial.SerialCommunicationParity* attribute), 17

spoll () (*hvl\_ccb.comm.visa.VisaCommunication* method), 17

spoll\_handler () (*hvl\_ccb.dev.visa.VisaDevice* method), 113

SR (*hvl\_ccb.dev.newport.NewportConfigCommands* attribute), 83

stage\_configuration  
     (*hvl\_ccb.dev.newport.NewportSMC100PPConfig* attribute), 92

start () (*hvl\_ccb.dev.base.Device* method), 53

start () (*hvl\_ccb.dev.base.DeviceSequenceMixin* method), 53

start () (*hvl\_ccb.dev.base.SingleCommDevice* method), 54

start () (*hvl\_ccb.dev.crylas.CryLasAttenuator* method), 55

start () (*hvl\_ccb.dev.crylas.CryLasLaser* method), 60

start () (*hvl\_ccb.dev.ea\_psi9000.PSI9000* method), 66

start () (*hvl\_ccb.dev.heinzinger.HeinzingerDI* method), 71

start () (*hvl\_ccb.dev.heinzinger.HeinzingerPNC* method), 72

start () (*hvl\_ccb.dev.labjack.LabJack* method), 78

start () (*hvl\_ccb.dev.mbw973.MBW973* method), 79

start () (*hvl\_ccb.dev.mbw973.Poller* method), 82

start () (*hvl\_ccb.dev.newport.NewportSMC100PP* method), 88

start () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG* method), 98

start () (*hvl\_ccb.dev.rs\_rto1024.RTO1024* method), 105

start () (*hvl\_ccb.dev.se\_ilst2.ILS2T* method), 111

start () (*hvl\_ccb.dev.supercube.base.SupercubeBase* method), 22

start () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base* method), 42

start () (*hvl\_ccb.dev.visa.VisaDevice* method), 113

start () (*hvl\_ccb.experiment\_manager.ExperimentManager* method), 117

start () (*hvl\_ccb.dev.mbw973.MBW973* method), 80

STARTING (*hvl\_ccb.experiment\_manager.ExperimentStatus* attribute), 117

States (*hvl\_ccb.dev.newport.NewportSMC100PP* attribute), 84

status (*hvl\_ccb.dev.supercube.constants.EarthingStick* attribute), 31

status (*hvl\_ccb.dev.supercube.constants.Safety* attribute), 35

status (*hvl\_ccb.experiment\_manager.ExperimentManager* attribute), 117

status\_1 (*hvl\_ccb.dev.supercube.constants.Door* attribute), 30

status\_1 (*hvl\_ccb.dev.supercube.constants.EarthingStick* attribute), 31

status\_1\_closed (*hvl\_ccb.dev.supercube2015.constants.EarthingStick* attribute), 46

status\_1\_connected

```

(hvl_ccb.dev.supercube2015.constants.EarthingStick      attribute), 46
attribute), 46
status_connected (hvl_ccb.dev.supercube2015.constants.EarthingStick
attribute), 46
status_error (hvl_ccb.dev.supercube2015.constants.Safety
attribute), 49
status_green (hvl_ccb.dev.supercube2015.constants.Safety
attribute), 49
status_open (hvl_ccb.dev.supercube2015.constants.EarthingStick
attribute), 49
status_ready_for_red
(hvl_ccb.dev.supercube2015.constants.Safety
attribute), 49
status_red (hvl_ccb.dev.supercube2015.constants.Safety
attribute), 49
stop (hvl_ccb.dev.supercube.constants.Errors attribute),
32
stop (hvl_ccb.dev.supercube2015.constants.Errors at-
tribute), 47
stop () (hvl_ccb.dev.base.Device method), 53
method), 54
stop () (hvl_ccb.dev.base.SingleCommDevice method),
54
stop () (hvl_ccb.dev.crylas.CryLasLaser method), 60
stop () (hvl_ccb.dev.heinzinger.HeinzingerDI method),
66
stop () (hvl_ccb.dev.labjack.LabJack method), 78
stop () (hvl_ccb.dev.mbw973.MBW973 method), 80
stop () (hvl_ccb.dev.mbw973.Poller method), 82
stop () (hvl_ccb.dev.newport.NewportSMC100PP
method), 88
stop () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method),
status_4_open (hvl_ccb.dev.supercube2015.constants.EarthingStick
attribute), 46
stop () (hvl_ccb.dev.rs_rto1024.RTO1024   method),
105
stop () (hvl_ccb.dev.se_il2t.ILS2T method), 111
stop () (hvl_ccb.dev.supercube.base.SupercubeBase
method), 22
stop () (hvl_ccb.dev.supercube2015.base.Supercube2015Base
method), 42
stop () (hvl_ccb.dev.visa.VisaDevice method), 113
stop () (hvl_ccb.dev.visa.VisaStatusPoller method), 114
stop () (hvl_ccb.experiment_manager.ExperimentManager
method), 117
stop_acquisition()
stop () (hvl_ccb.dev.rs_rto1024.RTO1024   method),
105
stop_motion () (hvl_ccb.dev.newport.NewportSMC100PP
method), 89
stop_number (hvl_ccb.dev.supercube2015.constants.Errors
Stopbits (hvl_ccb.comm.serial.SerialCommunicationConfig
status_closed (hvl_ccb.dev.supercube2015.constants.EarthingStick
attribute), 15

```

stopbits (*hvl\_ccb.comm.serial.SerialCommunicationConfig* attribute), 16

stopbits (*hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig* attribute), 58

stopbits (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig* attribute), 64

stopbits (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig* attribute), 74

stopbits (*hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig* attribute), 82

stopbits (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig* attribute), 95

stopbits (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig* attribute), 101

`StrEnumBase` (class in *hvl\_ccb.utils.enum*), 115

`sub_handler` (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig* attribute), 35

`sub_handler` (*hvl\_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig* attribute), 24

`sub_handler` (*hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig* attribute), 44

`Supercube2015Base` (class in *hvl\_ccb.dev.supercube2015.base*), 39

`Supercube2015WithFU` (class in *hvl\_ccb.dev.supercube2015.typ\_a*), 50

`SupercubeAOpcUaCommunication` (class in *hvl\_ccb.dev.supercube.typ\_a*), 36

`SupercubeAOpcUaCommunication` (class in *hvl\_ccb.dev.supercube2015.typ\_a*), 51

`SupercubeAOpcUaConfiguration` (class in *hvl\_ccb.dev.supercube.typ\_a*), 36

`SupercubeAOpcUaConfiguration` (class in *hvl\_ccb.dev.supercube2015.typ\_a*), 51

`SupercubeB` (class in *hvl\_ccb.dev.supercube.typ\_b*), 38

`SupercubeBase` (class in *hvl\_ccb.dev.supercube.base*), 20

`SupercubeBOpcUaCommunication` (class in *hvl\_ccb.dev.supercube.typ\_b*), 38

`SupercubeBOpcUaConfiguration` (class in *hvl\_ccb.dev.supercube.typ\_b*), 38

`SupercubeConfiguration` (class in *hvl\_ccb.dev.supercube.base*), 23

`SupercubeConfiguration` (class in *hvl\_ccb.dev.supercube2015.base*), 42

`SupercubeOpcEndpoint` (class in *hvl\_ccb.dev.supercube.constants*), 35

`SupercubeOpcEndpoint` (class in *hvl\_ccb.dev.supercube2015.constants*), 50

`SupercubeOpcUaCommunication` (class in *hvl\_ccb.dev.supercube.base*), 23

`SupercubeOpcUaCommunication` (class in *hvl\_ccb.dev.supercube2015.base*), 43

`SupercubeOpcUaCommunicationConfig` (class in *hvl\_ccb.dev.supercube.base*), 23

`supercubeOpcUaCommunicationConfig` (*hvl\_ccb.dev.supercube2015*.base), 43

`supercubeSubscriptionHandler` (*hvl\_ccb.dev.supercube.base*), 24

`supercubeSubscriptionHandler` (*hvl\_ccb.dev.supercube2015*.base), 44

`supercubeWithFU` (*hvl\_ccb.dev.supercube.typ\_a*), 36

`supercube2015constants.SafetyAttribute`, 49

`supercube2015constants.SafetyAttribute`, 35

`supercube2015constants.SafetyAttribute`, 49

`switchto_ready` (*hvl\_ccb.dev.supercube.constants.SafetyAttribute*), 49

`switchto_ready` (*hvl\_ccb.dev.supercube2015.constants.SafetyAttribute*), 49

`T` (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType* attribute), 76

`t13_1` (*hvl\_ccb.dev.supercube.constants.GeneralSockets* attribute), 32

`t13_1` (*hvl\_ccb.dev.supercube2015.constants.GeneralSockets* attribute), 47

`t13_2` (*hvl\_ccb.dev.supercube.constants.GeneralSockets* attribute), 32

`t13_2` (*hvl\_ccb.dev.supercube2015.constants.GeneralSockets* attribute), 47

`t13_3` (*hvl\_ccb.dev.supercube.constants.GeneralSockets* attribute), 32

`t13_3` (*hvl\_ccb.dev.supercube2015.constants.GeneralSockets* attribute), 47

`T13_SOCKET_PORTS` (in module *hvl\_ccb.dev.supercube.constants*), 35

`T13_SOCKET_PORTS` (in module *hvl\_ccb.dev.supercube2015.constants*), 50

`T4` (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.DeviceType* attribute), 9

`T4` (*hvl\_ccb.dev.labjack.LabJack.DeviceType* attribute), 75

`T7` (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.DeviceType* attribute), 9

`T7` (*hvl\_ccb.dev.labjack.LabJack.DeviceType* attribute), 75

`T7_PRO` (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.DeviceType* attribute), 9

`T7_PRO` (*hvl\_ccb.dev.labjack.LabJack.DeviceType* attribute), 75

`target_pulse_energy` (*hvl\_ccb.dev.crylas.CryLasLaser* attribute), 60

TCP (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig*)  
*attribute*, 9  
 TCP\_IP\_INSTR (*hvl\_ccb.comm.visa.VisaCommunicationConfig*)  
*InterfaceType* (*hvl\_ccb.dev.crylas.CryLasAttenuator*  
*attribute*), 55  
 TCP\_IP\_SOCKET (*hvl\_ccb.comm.visa.VisaCommunicationConfig*)  
*InterfaceType* (*hvl\_ccb.dev.supercube.constants.BreakdownDetection*  
*attribute*), 18  
 TEC1 (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus* at-  
*tribute*), 58  
 TEC2 (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus* at-  
*tribute*), 58  
 TEMP (*hvl\_ccb.dev.se\_ils2t.ILS2T.RegAddr* *attribute*), 108  
 TEN (*hvl\_ccb.dev.labjack.LabJack.AInRange* *attribute*), 74  
 TEN (*hvl\_ccb.dev.labjack.LabJack.CalMicroAmpere* at-  
*tribute*), 75  
 terminator (*hvl\_ccb.comm.serial.SerialCommunicationConfig*)  
*range* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus*  
*attribute*), 16  
 terminator (*hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig*)  
*range* (*hvl\_ccb.comm.serial.SerialCommunication*  
*attribute*), 58  
 terminator (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig*)  
*range* (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig*  
*attribute*), 64  
 terminator (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig*)  
*pfeiffer\_tpg.PfeifferTPG* *attribute*), 98  
 terminator (*hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig*)  
*attribute*, 82  
 terminator (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig*)  
*attribute*, 95  
 terminator (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig*)  
*attribute*, 101  
 timeout (*hvl\_ccb.comm.serial.SerialCommunicationConfig*)  
*attribute*, 16  
 timeout (*hvl\_ccb.comm.visa.VisaCommunicationConfig*)  
*attribute*, 19  
 timeout (*hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig*)  
*attribute*, 58  
 timeout (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig*)  
*attribute*, 64  
 timeout (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig*)  
*attribute*, 74  
 timeout (*hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig*)  
*attribute*, 82  
 timeout (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig*)  
*attribute*, 95  
 timeout (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig*)  
*attribute*, 101  
 timer\_callback () (*hvl\_ccb.dev.mbw973.Poller* *method*), 82  
 Torr (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits* *attribute*), 96  
 TPG25xA (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig.Model* *attribute*), 98  
 TPGx6x (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig.Model* *attribute*), 98

ConnectionType (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig*)  
*attribute*, 9

U

range (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus*  
*attribute*), 96  
 CommunicationConfig (*hvl\_ccb.comm.serial.SerialCommunication*  
*attribute*), 14  
 CommunicationConfig (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig*  
*attribute*), 11  
 CommunicationConfig (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig*)  
*pfeiffer\_tpg.PfeifferTPG* *attribute*), 98  
 CommunicationConfig (*hvl\_ccb.dev.se\_ils2t.ILS2TModbusTcpCommunicationConfig*)  
*attribute*, 113  
 CommunicationConfig (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*unit\_current* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*attribute*, 82  
 CommunicationConfig (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*unit\_voltage* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*attribute*, 95  
 CommunicationConfig (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*UnitCurrent* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*attribute*, 101  
 CommunicationConfig (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*UnitVoltage* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*attribute*, 71  
 CommunicationConfig (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*UNKNOWN* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*UnitCurrent* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*attribute*, 71  
 CommunicationConfig (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*UNKNOWN* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*UnitVoltage* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*)  
*attribute*, 71  
 CommunicationConfig (*hvl\_ccb.dev.crylas.CryLasLaser*)  
*update\_laser\_status* ()  
*method*, 60  
 CommunicationConfig (*hvl\_ccb.dev.crylas.CryLasLaser*)  
*update\_pulse\_iod* (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig*)  
*attribute*, 13  
 CommunicationConfig (*hvl\_ccb.dev.crylas.CryLasLaser*)  
*update\_repetition\_rate* ()  
*method*, 61  
 CommunicationConfig (*hvl\_ccb.dev.crylas.CryLasLaser*)  
*update\_shutter\_status* ()  
*method*, 61  
 CommunicationConfig (*hvl\_ccb.dev.crylas.CryLasLaser*)  
*update\_target\_pulse\_energy* ()  
*method*, 61  
 CommunicationConfig (*hvl\_ccb.dev.crylas.CryLasLaser*)  
*UpdateEspStageInfo* ()  
*method*, 61  
 CommunicationConfig (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*)  
*EspStageConfig* (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*)  
*attribute*, 90  
 CommunicationConfig (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig*)  
*USB* (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig*)  
*ConnectionType* (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig*)  
*attribute*, 9

user\_position\_offset  
     (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
         attribute), 92

user\_steps() (*hvl\_ccb.dev.se\_ilst.ILS2T* method),  
     111

**V**

V (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitVoltage*  
     attribute), 71

VA (*hvl\_ccb.dev.newport.NewportConfigCommands* at-  
     tribute), 83

value (*hvl\_ccb.dev.labjack.LabJack.AInRange* at-  
     tribute), 74

ValueEnum (class in *hvl\_ccb.utils.enum*), 115

VB (*hvl\_ccb.dev.newport.NewportConfigCommands* at-  
     tribute), 83

velocity (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
     attribute), 92

visa\_backend (*hvl\_ccb.comm.visa.VisaCommunicationConfig*  
     attribute), 19

VisaCommunication (class in *hvl\_ccb.comm.visa*),  
     17

VisaCommunicationConfig (class in  
     *hvl\_ccb.comm.visa*), 18

VisaCommunicationConfig.InterfaceType  
     (class in *hvl\_ccb.comm.visa*), 18

VisaCommunicationError, 19

VisaDevice (class in *hvl\_ccb.dev.visa*), 113

VisaDeviceConfig (class in *hvl\_ccb.dev.visa*), 114

VisaStatusPoller (class in *hvl\_ccb.dev.visa*), 114

Volt (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits*  
     attribute), 96

VOLT (*hvl\_ccb.dev.se\_ilst.ILS2T.RegAddr* attribute),  
     108

voltage\_lower\_limit  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
         attribute), 67

voltage\_max (*hvl\_ccb.dev.supercube.constants.Power*  
     attribute), 34

voltage\_max (*hvl\_ccb.dev.supercube2015.constants.Power*  
     attribute), 48

voltage\_primary (*hvl\_ccb.dev.supercube.constants.Power*  
     attribute), 34

voltage\_primary (*hvl\_ccb.dev.supercube2015.constants.Power*  
     attribute), 48

voltage\_slope (*hvl\_ccb.dev.supercube.constants.Power*  
     attribute), 34

voltage\_slope (*hvl\_ccb.dev.supercube2015.constants.Power*  
     attribute), 49

voltage\_target (*hvl\_ccb.dev.supercube.constants.Power*  
     attribute), 34

voltage\_target (*hvl\_ccb.dev.supercube2015.constants.Power*  
     attribute), 49

voltage\_upper\_limit

**W**

WAIT\_AFTER\_WRITE (*hvl\_ccb.comm.visa.VisaCommunication*  
     attribute), 17

wait\_operation\_complete()  
     (*hvl\_ccb.dev.visa.VisaDevice* method), 113

wait\_sec\_initialisation  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
         attribute), 67

wait\_sec\_post\_absolute\_position  
     (*hvl\_ccb.dev.se\_ilst.ILS2TConfig* attribute),  
     112

wait\_sec\_post\_enable  
     (*hvl\_ccb.dev.se\_ilst.ILS2TConfig* attribute),  
     112

wait\_sec\_post\_relative\_step  
     (*hvl\_ccb.dev.se\_ilst.ILS2TConfig* attribute),  
     112

wait\_sec\_read\_text\_nonempty  
     (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig*  
         attribute), 74

wait\_sec\_settings\_effect  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
         attribute), 67

wait\_sec\_stop\_commands  
     (*hvl\_ccb.dev.heinzinger.HeinzingerConfig*  
         attribute), 69

wait\_sec\_system\_lock  
     (*hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
         attribute), 67

wait\_until\_motor\_initialized()  
     (*hvl\_ccb.dev.newport.NewportSMC100PP*  
         method), 89

wait\_until\_move\_finished()  
     (*hvl\_ccb.dev.newport.NewportSMC100PP*  
         method), 89

wait\_until\_ready()  
     (*hvl\_ccb.dev.crylas.CryLasLaser* method),  
     61

warning (*hvl\_ccb.dev.supercube.constants.Errors* at-  
     tribute), 32

write () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.ConnectionType*  
     attribute), 9

write () (*hvl\_ccb.comm.opc.OpcUaCommunication*  
     method), 12

write () (*hvl\_ccb.comm.visa.VisaCommunication*  
     method), 18

write () (*hvl\_ccb.dev.mbw973.MBW973* method), 80

write () (*hvl\_ccb.dev.supercube.base.SupercubeBase*  
     method), 22

write () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base*  
     method), 42

write\_bytes () (*hvl\_ccb.comm.serial.SerialCommunication method*), 15  
write\_name () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunication method*), 8  
write\_names () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunication method*), 8  
write\_registers ()  
    (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication method*), 10  
write\_termination  
    (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 19  
write\_text () (*hvl\_ccb.comm.serial.SerialCommunication method*), 15  
WRONG\_ESP\_STAGE (*hvl\_ccb.dev.newport.NewportSMC100PP.MotorErrors attribute*), 83

## Z

zx (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 83