

---

# HVL Common Code Base Documentation

*Release 0.4.0*

**Mikolaj Rybiński, David Gruber**

**Jul 16, 2020**



---

## Contents:

---

<b>1</b>	<b>HVL Common Code Base</b>	<b>1</b>
1.1	Features . . . . .	1
1.2	Credits . . . . .	2
<b>2</b>	<b>Installation</b>	<b>3</b>
2.1	Stable release . . . . .	3
2.2	From sources . . . . .	3
2.3	Additional system libraries . . . . .	3
<b>3</b>	<b>Usage</b>	<b>5</b>
<b>4</b>	<b>API Documentation</b>	<b>7</b>
4.1	hvl_ccb package . . . . .	7
<b>5</b>	<b>Contributing</b>	<b>127</b>
5.1	Types of Contributions . . . . .	127
5.2	Get Started! . . . . .	128
5.3	Merge Request Guidelines . . . . .	129
5.4	Tips . . . . .	129
5.5	Deploying . . . . .	129
<b>6</b>	<b>Credits</b>	<b>131</b>
6.1	Development Lead . . . . .	131
6.2	Contributors . . . . .	131
<b>7</b>	<b>History</b>	<b>133</b>
7.1	0.4.0 (2020-07-16) . . . . .	133
7.2	0.3.5 (2020-02-18) . . . . .	134
7.3	0.3.4 (2019-12-20) . . . . .	134
7.4	0.3.3 (2019-05-08) . . . . .	134
7.5	0.3.2 (2019-05-08) . . . . .	134
7.6	0.3.1 (2019-05-02) . . . . .	134
7.7	0.3 (2019-05-02) . . . . .	134
7.8	0.2.1 (2019-04-01) . . . . .	135
7.9	0.2.0 (2019-03-31) . . . . .	135
7.10	0.1.0 (2019-02-06) . . . . .	135

<b>8 Indices and tables</b>	<b>137</b>
<b>Python Module Index</b>	<b>139</b>
<b>Index</b>	<b>141</b>

# CHAPTER 1

---

## HVL Common Code Base

---

Python common code base to control devices high voltage research devices, in particular, as used in Christian Franck's High Voltage Lab (HVL), D-ITET, ETH.

- Free software: GNU General Public License v3

- **Documentation:**

- if you're planning to develop start w/ reading “CONTRIBUTING.rst”, otherwise either
  - read [HVL CCB documentation at RTD](#), or
  - install *Sphinx* and *sphinx\_rtd\_theme* Python packages and locally build docs on Windows in git-bash by running:

```
$ ./make.sh docs
```

from a shell with Make installed by running:

```
$ make docs
```

The target index HTML (“docs/\_build/html/index.html”) will open automatically in your Web browser.

## 1.1 Features

Manage experiments with `ExperimentManager` instance controlling one or more of the following devices:

- a MBW973 SF6 Analyzer / dew point mirror over a serial connection (COM-ports)
- a LabJack (T7-PRO) device using a LabJack LJM Library for communication

- a Schneider Electric ILS2T stepper motor drive over Modbus TCP
- a Elektro-Automatik PSI9000 DC power supply using VISA over TCP for communication
- a Rhode & Schwarz RTO 1024 oscilloscope using VISA interface over TCP : : INSTR
- a state-of-the-art HVL in-house Supercube device variants using an OPC UA client
- a Heinzinger Digital Interface I/II and a Heinzinger PNC power supply over a serial connection
- a passively Q-switched Pulsed Laser and a laser attenuator from CryLas over a serial connection
- a Newport SMC100PP single axis motion controller for 2-phase stepper motors over a serial connection
- a Pfeiffer TPG controller (TPG 25x, TPG 26x and TPG 36x) for Compact pressure Gauges
- a SST Luminox Oxygen sensor device controller over a serial connection

## **1.2 Credits**

This package was created with [Cookiecutter](#) and the [audreyr/cookiecutter-pypackage](#) project template.

# CHAPTER 2

---

## Installation

---

### 2.1 Stable release

To install HVL Common Code Base, run this command in your terminal:

```
$ pip install hvl_ccb
```

This is the preferred method to install HVL Common Code Base, as it will always install the most recent stable release. If you don't have `pip` installed, this [Python installation guide](#) can guide you through the process.

### 2.2 From sources

The sources for HVL Common Code Base can be downloaded from the [GitLab](#) repo.

You can either clone the repository:

```
$ git clone git@gitlab.com:ethz_hvl/hvl_ccb.git
```

Or download the [tarball](#):

```
$ curl -OL https://gitlab.com/ethz_hvl/hvl_ccb/-/archive/master/hvl_ccb.tar.gz
```

Once you have a copy of the source, you can install it with:

```
$ python setup.py install
```

### 2.3 Additional system libraries

Please note that for some of the dependencies, like e.g. for `labjack-ljm`, you need to separately install additional system libraries.

For [LJM Library](#) please make sure that you have installed version 2019\_02\_14 or higher.

# CHAPTER 3

---

## Usage

---

To use HVL Common Code Base in a project:

```
import hvl_ccb
```



# CHAPTER 4

---

## API Documentation

---

### 4.1 hvl\_ccb package

#### 4.1.1 Subpackages

##### [hvl\\_ccb.comm package](#)

###### **Submodules**

###### [hvl\\_ccb.comm.base module](#)

Module with base classes for communication protocols.

**class** `hvl_ccb.comm.base.CommunicationProtocol(config)`  
Bases: `hvl_ccb.configuration.ConfigurationMixin, abc.ABC`

Communication protocol abstract base class.

Specifies the methods to implement for communication protocol, as well as implements some default settings and checks.

**access\_lock = None**

Access lock to use with context manager when accessing the communication protocol (thread safety)

**close()**

Close the communication protocol

**open()**

Open communication protocol

###### [hvl\\_ccb.comm.labjack\\_ljm module](#)

Communication protocol for LabJack using the LJM Library. Originally developed and tested for LabJack T7-PRO.

Makes use of the LabJack LJM Library Python wrapper. This wrapper needs an installation of the LJM Library for Windows, Mac OS X or Linux. Go to: <https://labjack.com/support/software/installers/ljm> and <https://labjack.com/support/software/examples/ljm/python>

**class** `hvl_ccb.comm.labjack_ljm.LJMCommunication`(*configuration*)

Bases: `hvl_ccb.comm.base.CommunicationProtocol`

Communication protocol implementing the LabJack LJM Library Python wrapper.

**close()** → None

Close the communication port.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**is\_open**

Flag indicating if the communication port is open.

**Returns** *True* if the port is open, otherwise *False*

**open()** → None

Open the communication port.

**read\_name**(\**names*, *return\_num\_type*: `Type[numbers.Real]`) = <class 'float'> →

`Union[numbers.Real, Sequence[numbers.Real]]`

Read one or more input numeric values by name.

## Parameters

- **names** – one or more names to read out from the LabJack
- **return\_num\_type** – optional numeric type specification for return values; by default `float`.

**Returns** answer of the LabJack, either single number or multiple numbers in a sequence, respectively, when one or multiple names to read were given

**Raises** `TypeError` – if read value of type not compatible with *return\_num\_type*

**write\_name**(*name*: str, *value*: `numbers.Real`) → None

Write one value to a named output.

## Parameters

- **name** – String or with name of LabJack IO
- **value** – is the value to write to the named IO port

**write\_names**(*name\_value\_dict*: `Dict[str, numbers.Real]`) → None

Write more than one value at once to named outputs.

**Parameters** **name\_value\_dict** – is a dictionary with string names of LabJack IO as keys and corresponding numeric values

**class** `hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig`(*device\_type*: `Union[str, hvl_ccb.dev.labjack.DeviceType]`) = `'ANY'`, *connection\_type*: `Union[str, hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig]` = `'ANY'`, *identifier*: str = `'ANY'`)

Bases: `object`

Configuration dataclass for *LJMCommunication*.

```
class ConnectionType(*args, **kwds)
    Bases: hvl_ccb.utils.enum.AutoNumberNameEnum

    LabJack connection type.

    ANY = 1
    ETHERNET = 4
    TCP = 3
    USB = 2
    WIFI = 5

class DeviceType(*args, **kwds)
    Bases: hvl_ccb.utils.enum.AutoNumberNameEnum

    LabJack device types.

    Can be also looked up by ambiguous Product ID (p_id) or by instance name: `python
    LabJackDeviceType(4) is LabJackDeviceType('T4')` 

    ANY = 1
    T4 = 2
    T7 = 3
    T7_PRO = 4

    get_by_p_id = <bound method DeviceType.get_by_p_id of <aenum 'DeviceType'>>
```

**clean\_values()** → None

Performs value checks on *device\_type* and *connection\_type*.

**connection\_type** = 'ANY'

Can be either string or of enum *ConnectionType*.

**device\_type** = 'ANY'

Can be either string 'ANY', 'T7\_PRO', 'T7', 'T4', or of enum *DeviceType*.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**identifier** = 'ANY'

The identifier specifies information for the connection to be used. This can be an IP address, serial number, or device name. See the LabJack docs ( <https://labjack.com/support/software/api/ljm/function-reference/ljmopens/identifier-parameter>) for more information.

**is\_configdataclass** = True

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception hvl\_ccb.comm.labjack\_ljm.LJMCommunicationError**

Bases: Exception

Errors coming from LJMCommunication.

### **hvl\_ccb.comm.modbus\_tcp module**

Communication protocol for modbus TCP ports. Makes use of the pymodbus library.

**class hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication(configuration)**  
Bases: *hvl\_ccb.comm.base.CommunicationProtocol*

Implements the Communication Protocol for modbus TCP.

**close()**

Close the Modbus TCP connection.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**open()** → None

Open the Modbus TCP connection.

**Raises ModbusTcpConnectionFailedException** – if the connection fails.

**read\_holding\_registers(address: int, count: int)** → List[int]

Read specified number of register starting with given address and return the values from each register.

#### Parameters

- **address** – address of the first register
- **count** – count of registers to read

**Returns** list of *int* values

**read\_input\_registers(address: int, count: int)** → List[int]

Read specified number of register starting with given address and return the values from each register in a list.

#### Parameters

- **address** – address of the first register
- **count** – count of registers to read

**Returns** list of *int* values

**write\_registers(address: int, values: Union[List[int], int])**

Write values from the specified address forward.



```
class hvl_ccb.comm.opc.OpcUaCommunication(config)
Bases: hvl_ccb.comm.base.CommunicationProtocol

Communication protocol implementing an OPC UA connection. Makes use of the package python-opcua.

close() → None
    Close the connection to the OPC UA server.

static config_cls()
    Return the default configdataclass class

    Returns a reference to the default configdataclass class

init_monitored_nodes(node_id: Union[object, Iterable[T_co]], ns_index: int) → None
    Initialize monitored nodes.

    Parameters
        • node_id – one or more strings of node IDs; node IDs are always casted via str() method here, hence do not have to be strictly string objects.
        • ns_index – the namespace index the nodes belong to.

    Raises OpcUaCommunicationIOError – when protocol was not opened or can't communicate with a OPC UA server

is_open
    Flag indicating if the communication port is open.

    Returns True if the port is open, otherwise False

open() → None
    Open the communication to the OPC UA server.

    Raises OpcUaCommunicationIOError – when communication port cannot be opened.

read(node_id, ns_index)
    Read a value from a node with id and namespace index.

    Parameters
        • node_id – the ID of the node to read the value from
        • ns_index – the namespace index of the node

    Returns the value of the node object.

    Raises OpcUaCommunicationIOError – when protocol was not opened or can't communicate with a OPC UA server

write(node_id, ns_index, value) → None
    Write a value to a node with name name.

    Parameters
        • node_id – the id of the node to write the value to.
        • ns_index – the namespace index of the node.
        • value – the value to write.

    Raises OpcUaCommunicationIOError – when protocol was not opened or can't communicate with a OPC UA server
```

```
class hvl_ccb.comm.opc.OpcUaCommunicationConfig(host: str, endpoint_name: str,  

                                                port: int = 4840, sub_handler:  

                                                hvl_ccb.comm.opc.OpcUaSubHandler  

                                                = <hvl_ccb.comm.opc.OpcUaSubHandler  

                                                object>, update_period: int = 500,  

                                                wait_timeout_retry_sec: Union[int,  

                                                float] = 1, max_timeout_retry_nr: int  

                                                = 5)
```

Bases: object

Configuration dataclass for OPC UA Communciation.

**clean\_values()**

**endpoint\_name = None**

Endpoint of the OPC server, this is a path like ‘OPCUA/SimulationServer’

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**host = None**

Hostname or IP-Address of the OPC UA server.

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**max\_timeout\_retry\_nr = 5**

Maximal number of call re-tries on underlying OPC UA client timeout error

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**port = 4840**

Port of the OPC UA server to connect to.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**sub\_handler = <hvl\_ccb.comm.opc.OpcUaSubHandler object>**

object to use for handling subscriptions.

**update\_period = 500**

Update period for generating datachange events in OPC UA [milli seconds]

**wait\_timeout\_retry\_sec = 1**

Wait time between re-trying calls on underlying OPC UA client timeout error

```
exception hvl_ccb.comm.opc.OpcUaCommunicationIOError
Bases: OSError
    OPC-UA communication I/O error.

exception hvl_ccb.comm.opc.OpcUaCommunicationTimeoutError
Bases: hvl_ccb.comm.opc.OpcUaCommunicationIOError
    OPC-UA communication timeout error.

class hvl_ccb.comm.opc.OpcUaSubHandler
Bases: object
    Base class for subscription handling of OPC events and data change events. Override methods from this class to add own handling capabilities.

    To receive events from server for a subscription data_change and event methods are called directly from receiving thread. Do not do expensive, slow or network operation there. Create another thread if you need to do such a thing.

        datachange_notification(node, val, data)
        event_notification(event)
```

## hvl\_ccb.comm.serial module

Communication protocol for serial ports. Makes use of the `pySerial` library.

```
class hvl_ccb.comm.serial.SerialCommunication(configuration)
Bases: hvl_ccb.comm.base.CommunicationProtocol
```

Implements the Communication Protocol for serial ports.

```
ENCODING = 'utf-8'
```

```
UNICODE_HANDLING = 'replace'
```

```
close()
```

Close the serial connection.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
is_open
```

Flag indicating if the serial port is open.

**Returns** `True` if the serial port is open, otherwise `False`

```
open()
```

Open the serial connection.

**Raises** `SerialCommunicationIOError` – when communication port cannot be opened.

```
read_bytes(size: int = 1) → bytes
```

Read the specified number of bytes from the serial port. The input buffer may hold additional data afterwards.

This method uses `self.access_lock` to ensure thread-safety.

**Parameters** `size` – number of bytes to read

**Returns** Bytes read from the serial port; `b''` if there was nothing to read.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

`read_text() → str`

Read one line of text from the serial port. The input buffer may hold additional data afterwards, since only one line is read.

This method uses `self.access_lock` to ensure thread-safety.

**Returns** String read from the serial port; ‘’ if there was nothing to read.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

`write_bytes(data: bytes) → int`

Write bytes to the serial port.

This method uses `self.access_lock` to ensure thread-safety.

**Parameters** `data` – data to write to the serial port

**Returns** number of bytes written

**Raises** `SerialCommunicationIOError` – when communication port is not opened

`write_text(text: str)`

Write text to the serial port. The text is encoded and terminated by the configured terminator.

This method uses `self.access_lock` to ensure thread-safety.

**Parameters** `text` – Text to send to the port.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

`class hvl_ccb.comm.serial.SerialCommunicationBytesize(*args, **kwds)`

Bases: `hvl_ccb.utils.enum.ValueEnum`

Serial communication bytesize.

`EIGHTBITS = 8`

`FIVEBITS = 5`

`SEVENBITS = 7`

`SIXBITS = 6`

`class hvl_ccb.comm.serial.SerialCommunicationConfig(port: str, baudrate: int, parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity], stopbits: Union[int, float, hvl_ccb.comm.serial.SerialCommunicationStopbits], bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize], terminator: bytes = b'\r\n', timeout: Union[int, float] = 2)`

Bases: `object`

Configuration dataclass for `SerialCommunication`.

**Bytesize**

alias of `SerialCommunicationBytesize`

**Parity**

alias of `SerialCommunicationParity`

**Stopbits**

alias of `SerialCommunicationStopbits`

**baudrate = None**  
Baudrate of the serial port

**bytesize = None**  
Size of a byte, 5 to 8

**clean\_values()**

**create\_serial\_port() → serial.serialposix.Serial**  
Create a serial port instance according to specification in this configuration

**Returns** Closed serial port instance

**force\_value(fieldname, value)**  
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = None**  
Parity to be used for the connection.

**port = None**  
Port is a string referring to a COM-port (e.g. 'COM3') or a URL. The full list of capabilities is found [on the pyserial documentation](#).

**classmethod required\_keys() → Sequence[str]**  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = None**  
Stopbits setting, can be 1, 1.5 or 2.

**terminator = b'\r\n'**  
The terminator character. Typically this is b'\r\n' or b'\n', but can also be b'\r' or other combinations.

**terminator\_str() → str**

**timeout = 2**  
Timeout in seconds for the serial port

**exception hvl\_ccb.comm.serial.SerialCommunicationIOError**  
Bases: OSError

Serial communication related I/O errors.

```
class hvl_ccb.comm.serial.SerialCommunicationParity (*args, **kwds)
    Bases: hvl_ccb.utils.enum.ValueEnum

    Serial communication parity.

    EVEN = 'E'
    MARK = 'M'
    NAMES = {'E': 'Even', 'M': 'Mark', 'N': 'None', 'O': 'Odd', 'S': 'Space'}
    NONE = 'N'
    ODD = 'O'
    SPACE = 'S'

class hvl_ccb.comm.serial.SerialCommunicationStopbits (*args, **kwds)
    Bases: hvl_ccb.utils.enum.ValueEnum

    Serial communication stopbits.

    ONE = 1
    ONE_POINT_FIVE = 1.5
    TWO = 2
```

## hvl\_ccb.comm.visa module

Communication protocol for VISA. Makes use of the pyvisa library. The backend can be NI-Visa or pyvisa-py.

Information on how to install a VISA backend can be found here: [https://pyvisa.readthedocs.io/en/master/getting\\_nivisa.html](https://pyvisa.readthedocs.io/en/master/getting_nivisa.html)

So far only TCPIP SOCKET and TCPIP INSTR interfaces are supported.

```
class hvl_ccb.comm.visa.VisaCommunication (configuration)
    Bases: hvl_ccb.comm.base.CommunicationProtocol

    Implements the Communication Protocol for VISA / SCPI.

    MULTI_COMMANDS_MAX = 5
        The maximum of commands that can be sent in one round is 5 according to the VISA standard.

    MULTI_COMMANDS_SEPARATOR = ';'
        The character to separate two commands is ; according to the VISA standard.

    WAIT_AFTER_WRITE = 0.08
        Small pause in seconds to wait after write operations, allowing devices to really do what we tell them before continuing with further tasks.

    close() → None
        Close the VISA connection and invalidates the handle.

    static config_cls() → Type[hvl_ccb.comm.visa.VisaCommunicationConfig]
        Return the default configdataclass class.

        Returns a reference to the default configdataclass class

    open() → None
        Open the VISA connection and create the resource.
```

**query** (\*commands) → Union[str, Tuple[str, ...]]

A combination of write(message) and read.

**Parameters** **commands** – list of commands

**Returns** list of values

**Raises** **VisaCommunicationError** – when connection was not started, or when trying to issue too many commands at once.

**spoll()** → int

Execute serial poll on the device. Reads the status byte register STB. This is a fast function that can be executed periodically in a polling fashion.

**Returns** integer representation of the status byte

**Raises** **VisaCommunicationError** – when connection was not started

**write** (\*commands) → None

Write commands. No answer is read or expected.

**Parameters** **commands** – one or more commands to send

**Raises** **VisaCommunicationError** – when connection was not started

```
class hvl_ccb.comm.visa.VisaCommunicationConfig(host: str, interface_type: Union[str, hvl_ccb.comm.visa.VisaCommunicationConfig.InterfaceType], board: int = 0, port: int = 5025, timeout: int = 5000, chunk_size: int = 204800, open_timeout: int = 1000, write_termination: str = '\n', read_termination: str = '\n', visa_backend: str = '')
```

Bases: object

VISACommunication configuration dataclass.

**class InterfaceType** (\*args, \*\*kwds)

Bases: [hvl\\_ccb.utils.enum.AutoNumberNameEnum](#)

Supported VISA Interface types.

**TCPIP\_INSTR** = 2

VXI-11 protocol

**TCPIP\_SOCKET** = 1

VISA-Raw protocol

**address** (host: str, port: int = None, board: int = None) → str

Address string specific to the VISA interface type.

**Parameters**

- **host** – host IP address
- **port** – optional TCP port
- **board** – optional board number

**Returns** address string

**address**

Address string depending on the VISA protocol's configuration.

**Returns** address string corresponding to current configuration

**board** = 0

Board number is typically 0 and comes from old bus systems.

---

```
chunk_size = 204800
    Chunk size is the allocated memory for read operations. The standard is 20kB, and is increased per default here to 200kB. It is specified in bytes.

clean_values()

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.
```

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

```
host = None
    IP address of the VISA device. DNS names are currently unsupported.

interface_type = None
    Interface type of the VISA connection, being one of InterfaceType.
```

```
is_configdataclass = True

classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.
```

```
open_timeout = 1000
    Timeout for opening the connection, in milli seconds.

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.
```

```
port = 5025
    TCP port, standard is 5025.

read_termination = '\n'
    Read termination character.

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.
```

```
timeout = 5000
    Timeout for commands in milli seconds.

visa_backend = ''
    Specifies the path to the library to be used with PyVISA as a backend. Defaults to None, which is NI-VISA (if installed), or pyvisa-py (if NI-VISA is not found). To force the use of pyvisa-py, specify '@py' here.

write_termination = '\n'
    Write termination character.
```

```
exception hvl_ccb.comm.visa.VisaCommunicationError
    Bases: Exception
```

Base class for VisaCommunication errors.

## Module contents

Communication protocols subpackage.

### **hvl\_ccb.dev package**

#### Subpackages

##### **hvl\_ccb.dev.supercube package**

#### Submodules

##### **hvl\_ccb.dev.supercube.base module**

Base classes for the Supercube device.

```
class hvl_ccb.dev.supercube.base.Poller(spoll_handler: Callable, polling_delay_sec: Union[int, float] = 0, polling_interval_sec: Union[int, float] = 1, polling_timeout_sec: Union[int, float, None] = None)
```

Bases: object

Poller class wrapping *concurrent.futures.ThreadPoolExecutor* which enables passing of results and errors out of the polling thread.

**is\_polling()** → bool

Check if device status is being polled.

**Returns** *True* when polling thread is set and alive

**start\_polling()** → bool

Start polling.

**Returns** *True* if was not polling before, *False* otherwise

**stop\_polling()** → Tuple[bool, Optional[object]]

Stop polling.

Wait for until polling function returns a result as well as any exception that might have been raised within a thread.

**Returns** *True* if was polling before, *False* otherwise, and last result of the polling function call.

**Raises** polling function exceptions

```
class hvl_ccb.dev.supercube.base.SupercubeBase(com, dev_config=None)
```

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Base class for Supercube variants.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**display\_message\_board()** → None

Display 15 newest messages

**display\_status\_board()** → None

Display status board.

**get\_cee16\_socket()** → bool

Read the on-state of the IEC CEE16 three-phase power socket.

**Returns** the on-state of the CEE16 power socket

**get\_door\_status(door: int)** → hvl\_ccb.dev.supercube.constants.DoorStatus

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**Parameters** `door` – the door number (1..3)

**Returns** the door status

**get\_earthing\_stick\_manual(number: int)** → hvl\_ccb.dev.supercube.constants.EarthingStickOperation

Get the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters** `number` – number of the earthing stick (1..6)

**Returns** operation of the earthing stick in a manual operating mode (open == 0, close == 1)

**Raises** `ValueError` – when earthing stick number is not valid

**get\_earthing\_stick\_operating\_status(number: int)** → hvl\_ccb.dev.supercube.constants.EarthingStickOperatingStatus

Get the operating status of an earthing stick.

**Parameters** `number` – number of the earthing stick (1..6)

**Returns** earthing stick operating status (auto == 0, manual == 1)

**Raises** `ValueError` – when earthing stick number is not valid

**get\_earthing\_stick\_status(number: int)** → hvl\_ccb.dev.supercube.constants.EarthingStickStatus

Get the status of an earthing stick, whether it is closed, open or undefined (moving).

**Parameters** `number` – number of the earthing stick (1..6)

**Returns** earthing stick status

**Raises** `ValueError` – when earthing stick number is not valid

**get\_measurement\_ratio(channel: int)** → float

Get the set measurement ratio of an AC/DC analog input channel. Every input channel has a divider ratio assigned during setup of the Supercube system. This ratio can be read out.

**Parameters** `channel` – number of the input channel (1..4)

**Returns** the ratio

**Raises** `ValueError` – when channel is not valid

**get\_measurement\_voltage(channel: int)** → float

Get the measured voltage of an analog input channel. The voltage read out here is already scaled by the configured divider ratio.

**Parameters** `channel` – number of the input channel (1..4)

**Returns** measured voltage

**Raises** `ValueError` – when channel is not valid

**get\_status ()** → hvl\_ccb.dev.supercube.constants.SafetyStatus

Get the safety circuit status of the Supercube. :return: the safety status of the supercube's state machine.

**get\_support\_input (port: int, contact: int)** → bool

Get the state of a support socket input.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital input read state

**Raises ValueError** – when port or contact number is not valid

**get\_support\_output (port: int, contact: int)** → bool

Get the state of a support socket output.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital output read state

**Raises ValueError** – when port or contact number is not valid

**get\_t13\_socket (port: int)** → bool

Read the state of a SEV T13 power socket.

**Parameters** **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*

**Returns** on-state of the power socket

**Raises ValueError** – when port is not valid

**operate (state: bool)** → None

Set operate state. If the state is RedReady, this will turn on the high voltage and close the safety switches.

**Parameters** **state** – set operate state

**operate\_earthing\_stick (number: int, operation: hvl\_ccb.dev.supercube.constants.EarthingStickOperation)**

→ None

Operation of an earthing stick, which is set to manual operation. If an earthing stick is set to manual, it stays closed even if the system is in states RedReady or RedOperate.

**Parameters**

- **number** – number of the earthing stick (1..6)
- **operation** – earthing stick manual status (close or open)

**Raises SupercubeEarthingStickOperationException** – when operating status of given number's earthing stick is not manual

**quit\_error ()** → None

Quits errors that are active on the Supercube.

**read (node\_id: str)**

Local wrapper for the OPC UA communication protocol read method.

**Parameters** **node\_id** – the id of the node to read.

**Returns** the value of the variable

**ready** (*state: bool*) → None

Set ready state. Ready means locket safety circuit, red lamps, but high voltage still off.

**Parameters** **state** – set ready state

**set\_cee16\_socket** (*state: bool*) → None

Switch the IEC CEE16 three-phase power socket on or off.

**Parameters** **state** – desired on-state of the power socket

**Raises** **ValueError** – if state is not of type bool

**set\_message\_board** (*msgs: List[str], display\_board: bool = True*) → None

Fills messages into message board that display that 15 newest messages with a timestamp.

**Parameters**

- **msgs** – list of strings

- **display\_board** – display 15 newest messages if *True* (default)

**Raises** **ValueError** – if there are too many messages or the positions indices are invalid.

**set\_remote\_control** (*state: bool*) → None

Enable or disable remote control for the Supercube. This will effectively display a message on the touch-screen HMI.

**Parameters** **state** – desired remote control state

**set\_status\_board** (*msgs: List[str], pos: List[int] = None, clear\_board: bool = True, display\_board: bool = True*) → None

Sets and displays a status board. The messages and the position of the message can be defined.

**Parameters**

- **msgs** – list of strings

- **pos** – list of integers [0...14]

- **clear\_board** – clear unspecified lines if *True* (default), keep otherwise

- **display\_board** – display new status board if *True* (default)

**Raises** **ValueError** – if there are too many messages or the positions indices are invalid.

**set\_support\_output** (*port: int, contact: int, state: bool*) → None

Set the state of a support output socket.

**Parameters**

- **port** – is the socket number (1..6)

- **contact** – is the contact on the socket (1..2)

- **state** – is the desired state of the support output

**Raises** **ValueError** – when port or contact number is not valid

**set\_support\_output\_impulse** (*port: int, contact: int, duration: float = 0.2, pos\_pulse: bool = True*) → None

Issue an impulse of a certain duration on a support output contact. The polarity of the pulse (On-wait-Off or Off-wait-On) is specified by the pos\_pulse argument.

This function is blocking.

**Parameters**

- **port** – is the socket number (1..6)

- **contact** – is the contact on the socket (1..2)
- **duration** – is the length of the impulse in seconds
- **pos\_pulse** – is True, if the pulse shall be HIGH, False if it shall be LOW

**Raises ValueError** – when port or contact number is not valid

**set\_t13\_socket** (*port: int, state: bool*) → None

Set the state of a SEV T13 power socket.

#### Parameters

- **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*
- **state** – is the desired on-state of the socket

**Raises ValueError** – when port is not valid or state is not of type bool

**start()** → None

Starts the device. Sets the root node for all OPC read and write commands to the Siemens PLC object node which holds all our relevant objects and variables.

**stop()** → None

Stop the Supercube device. Deactivates the remote control and closes the communication protocol.

**write** (*node\_id, value*) → None

Local wrapper for the OPC UA communication protocol write method.

#### Parameters

- **node\_id** – the id of the node to read
- **value** – the value to write to the variable

```
class hvl_ccb.dev.supercube.base.SupercubeConfiguration(namespace_index:      int
                                                       = 3, polling_delay_sec: Union[int, float] = 5.0,
                                                       polling_interval_sec: Union[int, float] = 1.0)
```

Bases: object

Configuration dataclass for the Supercube devices.

**clean\_values()**

**force\_value** (*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**namespace\_index = 3**

Namespace of the OPC variables, typically this is 3 (coming from Siemens)

---

**classmethod optional\_defaults()** → Dict[str, object]  
 Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
polling_delay_sec = 5.0
polling_interval_sec = 1.0
```

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception** hvl\_ccb.dev.supercube.base.**SupercubeEarthingStickOperationError**  
 Bases: Exception

**class** hvl\_ccb.dev.supercube.base.**SupercubeOpcUaCommunication** (config)  
 Bases: *hvl\_ccb.comm.opc.OpcUaCommunication*

Communication protocol specification for Supercube devices.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class** hvl\_ccb.dev.supercube.base.**SupercubeOpcUaCommunicationConfig** (host:  
*str, end-point\_name:*  
*str, port:*  
*int = 4840,*  
*sub\_handler:*  
*hvl\_ccb.comm.opc.OpcUaSubHandler = <hvl\_ccb.dev.supercube.base.SupercubeEarthingStickOperationError object>,*  
*update\_period:*  
*int = 500,*  
*wait\_timeout\_retry\_sec: Union[int, float] = 1,*  
*max\_timeout\_retry\_nr: int = 5*)

Bases: *hvl\_ccb.comm.opc.OpcUaCommunicationConfig*

Communication protocol configuration for OPC UA, specifications for the Supercube devices.

**force\_value** (*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
sub_handler = <hvl_ccb.dev.supercube.base.SupercubeSubscriptionHandler object>
```

Subscription handler for data change events

```
class hvl_ccb.dev.supercube.base.SupercubeSubscriptionHandler
```

Bases: [hvl\\_ccb.comm.opc.OpcUaSubHandler](#)

OPC Subscription handler for datachange events and normal events specifically implemented for the Supercube devices.

```
datachange_notification(node: opcua.common.node.Node, val, data)
```

In addition to the standard operation (debug logging entry of the datachange), alarms are logged at INFO level using the alarm text.

#### Parameters

- **node** – the node object that triggered the datachange event
- **val** – the new value
- **data** –

## [hvl\\_ccb.dev.supercube.constants module](#)

Constants, variable names for the Supercube OPC-connected devices.

```
class hvl_ccb.dev.supercube.constants.AlarmText(*args, **kwds)
```

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

This enumeration contains textual representations for all error classes (stop, warning and message) of the Supercube system. Use the [AlarmText.get\(\)](#) method to retrieve the enum of an alarm number.

```
Alarm1 = 'STOP Emergency Stop 1'
```

```
Alarm10 = 'STOP Earthing stick 2 error while opening'
```

```
Alarm11 = 'STOP Earthing stick 3 error while opening'
```

```
Alarm12 = 'STOP Earthing stick 4 error while opening'
```

```
Alarm13 = 'STOP Earthing stick 5 error while opening'
```

```
Alarm14 = 'STOP Earthing stick 6 error while opening'
```

```
Alarm15 = 'STOP Earthing stick 1 error while closing'
```

```
Alarm16 = 'STOP Earthing stick 2 error while closing'
```

```
Alarm17 = 'STOP Earthing stick 3 error while closing'
```

```

Alarm18 = 'STOP Earthing stick 4 error while closing'
Alarm19 = 'STOP Earthing stick 5 error while closing'
Alarm2 = 'STOP Emergency Stop 2'
Alarm20 = 'STOP Earthing stick 6 error while closing'
Alarm21 = 'STOP Safety fence 1'
Alarm22 = 'STOP Safety fence 2'
Alarm23 = 'STOP OPC connection error'
Alarm24 = 'STOP Grid power failure'
Alarm25 = 'STOP UPS failure'
Alarm26 = 'STOP 24V PSU failure'
Alarm3 = 'STOP Emergency Stop 3'
Alarm4 = 'STOP Safety Switch 1 error'
Alarm41 = 'WARNING Door 1: Use earthing rod!'
Alarm42 = 'MESSAGE Door 1: Earthing rod is still in setup.'
Alarm43 = 'WARNING Door 2: Use earthing rod!'
Alarm44 = 'MESSAGE Door 2: Earthing rod is still in setup.'
Alarm45 = 'WARNING Door 3: Use earthing rod!'
Alarm46 = 'MESSAGE Door 3: Earthing rod is still in setup.'
Alarm47 = 'MESSAGE UPS charge < 85%'
Alarm48 = 'MESSAGE UPS running on battery'
Alarm5 = 'STOP Safety Switch 2 error'
Alarm6 = 'STOP Door 1 lock supervision'
Alarm7 = 'STOP Door 2 lock supervision'
Alarm8 = 'STOP Door 3 lock supervision'
Alarm9 = 'STOP Earthing stick 1 error while opening'
get = <bound method AlarmText.get of <aenum 'AlarmText'>>
not_defined = 'NO ALARM TEXT DEFINED'

class hvl_ccb.dev.supercube.constants.Alarms(*args, **kwds)
    Bases: hvl_ccb.dev.supercube.constants._AlarmEnumBase
    Alarms enumeration containing all variable NodeID strings for the alarm array.

    Alarm1 = '"DB_Alarm_HMI"."Alarm1"'
    Alarm10 = '"DB_Alarm_HMI"."Alarm10"'
    Alarm100 = '"DB_Alarm_HMI"."Alarm100"'
    Alarm101 = '"DB_Alarm_HMI"."Alarm101"'
    Alarm102 = '"DB_Alarm_HMI"."Alarm102"'
    Alarm103 = '"DB_Alarm_HMI"."Alarm103"'

```

```
Alarm104 = '"DB_Alarm_HMI"."Alarm104"'
Alarm105 = '"DB_Alarm_HMI"."Alarm105"'
Alarm106 = '"DB_Alarm_HMI"."Alarm106"'
Alarm107 = '"DB_Alarm_HMI"."Alarm107"'
Alarm108 = '"DB_Alarm_HMI"."Alarm108"'
Alarm109 = '"DB_Alarm_HMI"."Alarm109"'
Alarm11 = '"DB_Alarm_HMI"."Alarm11"'
Alarm110 = '"DB_Alarm_HMI"."Alarm110"'
Alarm111 = '"DB_Alarm_HMI"."Alarm111"'
Alarm112 = '"DB_Alarm_HMI"."Alarm112"'
Alarm113 = '"DB_Alarm_HMI"."Alarm113"'
Alarm114 = '"DB_Alarm_HMI"."Alarm114"'
Alarm115 = '"DB_Alarm_HMI"."Alarm115"'
Alarm116 = '"DB_Alarm_HMI"."Alarm116"'
Alarm117 = '"DB_Alarm_HMI"."Alarm117"'
Alarm118 = '"DB_Alarm_HMI"."Alarm118"'
Alarm119 = '"DB_Alarm_HMI"."Alarm119"'
Alarm12 = '"DB_Alarm_HMI"."Alarm12"'
Alarm120 = '"DB_Alarm_HMI"."Alarm120"'
Alarm121 = '"DB_Alarm_HMI"."Alarm121"'
Alarm122 = '"DB_Alarm_HMI"."Alarm122"'
Alarm123 = '"DB_Alarm_HMI"."Alarm123"'
Alarm124 = '"DB_Alarm_HMI"."Alarm124"'
Alarm125 = '"DB_Alarm_HMI"."Alarm125"'
Alarm126 = '"DB_Alarm_HMI"."Alarm126"'
Alarm127 = '"DB_Alarm_HMI"."Alarm127"'
Alarm128 = '"DB_Alarm_HMI"."Alarm128"'
Alarm129 = '"DB_Alarm_HMI"."Alarm129"'
Alarm13 = '"DB_Alarm_HMI"."Alarm13"'
Alarm130 = '"DB_Alarm_HMI"."Alarm130"'
Alarm131 = '"DB_Alarm_HMI"."Alarm131"'
Alarm132 = '"DB_Alarm_HMI"."Alarm132"'
Alarm133 = '"DB_Alarm_HMI"."Alarm133"'
Alarm134 = '"DB_Alarm_HMI"."Alarm134"'
Alarm135 = '"DB_Alarm_HMI"."Alarm135"'
Alarm136 = '"DB_Alarm_HMI"."Alarm136"'
```

```
Alarm137 = '"DB_Alarm_HMI"."Alarm137"'
Alarm138 = '"DB_Alarm_HMI"."Alarm138"'
Alarm139 = '"DB_Alarm_HMI"."Alarm139"'
Alarm14 = '"DB_Alarm_HMI"."Alarm14"'
Alarm140 = '"DB_Alarm_HMI"."Alarm140"'
Alarm141 = '"DB_Alarm_HMI"."Alarm141"'
Alarm142 = '"DB_Alarm_HMI"."Alarm142"'
Alarm143 = '"DB_Alarm_HMI"."Alarm143"'
Alarm144 = '"DB_Alarm_HMI"."Alarm144"'
Alarm145 = '"DB_Alarm_HMI"."Alarm145"'
Alarm146 = '"DB_Alarm_HMI"."Alarm146"'
Alarm147 = '"DB_Alarm_HMI"."Alarm147"'
Alarm148 = '"DB_Alarm_HMI"."Alarm148"'
Alarm149 = '"DB_Alarm_HMI"."Alarm149"'
Alarm15 = '"DB_Alarm_HMI"."Alarm15"'
Alarm150 = '"DB_Alarm_HMI"."Alarm150"'
Alarm151 = '"DB_Alarm_HMI"."Alarm151"'
Alarm16 = '"DB_Alarm_HMI"."Alarm16"'
Alarm17 = '"DB_Alarm_HMI"."Alarm17"'
Alarm18 = '"DB_Alarm_HMI"."Alarm18"'
Alarm19 = '"DB_Alarm_HMI"."Alarm19"'
Alarm2 = '"DB_Alarm_HMI"."Alarm2"'
Alarm20 = '"DB_Alarm_HMI"."Alarm20"'
Alarm21 = '"DB_Alarm_HMI"."Alarm21"'
Alarm22 = '"DB_Alarm_HMI"."Alarm22"'
Alarm23 = '"DB_Alarm_HMI"."Alarm23"'
Alarm24 = '"DB_Alarm_HMI"."Alarm24"'
Alarm25 = '"DB_Alarm_HMI"."Alarm25"'
Alarm26 = '"DB_Alarm_HMI"."Alarm26"'
Alarm27 = '"DB_Alarm_HMI"."Alarm27"'
Alarm28 = '"DB_Alarm_HMI"."Alarm28"'
Alarm29 = '"DB_Alarm_HMI"."Alarm29"'
Alarm3 = '"DB_Alarm_HMI"."Alarm3"'
Alarm30 = '"DB_Alarm_HMI"."Alarm30"'
Alarm31 = '"DB_Alarm_HMI"."Alarm31"'
Alarm32 = '"DB_Alarm_HMI"."Alarm32"'
```

```
Alarm33 = '"DB_Alarm_HMI"."Alarm33"'
Alarm34 = '"DB_Alarm_HMI"."Alarm34"'
Alarm35 = '"DB_Alarm_HMI"."Alarm35"'
Alarm36 = '"DB_Alarm_HMI"."Alarm36"'
Alarm37 = '"DB_Alarm_HMI"."Alarm37"'
Alarm38 = '"DB_Alarm_HMI"."Alarm38"'
Alarm39 = '"DB_Alarm_HMI"."Alarm39"'
Alarm4 = '"DB_Alarm_HMI"."Alarm4"'
Alarm40 = '"DB_Alarm_HMI"."Alarm40"'
Alarm41 = '"DB_Alarm_HMI"."Alarm41"'
Alarm42 = '"DB_Alarm_HMI"."Alarm42"'
Alarm43 = '"DB_Alarm_HMI"."Alarm43"'
Alarm44 = '"DB_Alarm_HMI"."Alarm44"'
Alarm45 = '"DB_Alarm_HMI"."Alarm45"'
Alarm46 = '"DB_Alarm_HMI"."Alarm46"'
Alarm47 = '"DB_Alarm_HMI"."Alarm47"'
Alarm48 = '"DB_Alarm_HMI"."Alarm48"'
Alarm49 = '"DB_Alarm_HMI"."Alarm49"'
Alarm5 = '"DB_Alarm_HMI"."Alarm5"'
Alarm50 = '"DB_Alarm_HMI"."Alarm50"'
Alarm51 = '"DB_Alarm_HMI"."Alarm51"'
Alarm52 = '"DB_Alarm_HMI"."Alarm52"'
Alarm53 = '"DB_Alarm_HMI"."Alarm53"'
Alarm54 = '"DB_Alarm_HMI"."Alarm54"'
Alarm55 = '"DB_Alarm_HMI"."Alarm55"'
Alarm56 = '"DB_Alarm_HMI"."Alarm56"'
Alarm57 = '"DB_Alarm_HMI"."Alarm57"'
Alarm58 = '"DB_Alarm_HMI"."Alarm58"'
Alarm59 = '"DB_Alarm_HMI"."Alarm59"'
Alarm6 = '"DB_Alarm_HMI"."Alarm6"'
Alarm60 = '"DB_Alarm_HMI"."Alarm60"'
Alarm61 = '"DB_Alarm_HMI"."Alarm61"'
Alarm62 = '"DB_Alarm_HMI"."Alarm62"'
Alarm63 = '"DB_Alarm_HMI"."Alarm63"'
Alarm64 = '"DB_Alarm_HMI"."Alarm64"'
Alarm65 = '"DB_Alarm_HMI"."Alarm65"'
```

```
Alarm66 = '"DB_Alarm_HMI"."Alarm66"'
Alarm67 = '"DB_Alarm_HMI"."Alarm67"'
Alarm68 = '"DB_Alarm_HMI"."Alarm68"'
Alarm69 = '"DB_Alarm_HMI"."Alarm69"'
Alarm7 = '"DB_Alarm_HMI"."Alarm7"'
Alarm70 = '"DB_Alarm_HMI"."Alarm70"'
Alarm71 = '"DB_Alarm_HMI"."Alarm71"'
Alarm72 = '"DB_Alarm_HMI"."Alarm72"'
Alarm73 = '"DB_Alarm_HMI"."Alarm73"'
Alarm74 = '"DB_Alarm_HMI"."Alarm74"'
Alarm75 = '"DB_Alarm_HMI"."Alarm75"'
Alarm76 = '"DB_Alarm_HMI"."Alarm76"'
Alarm77 = '"DB_Alarm_HMI"."Alarm77"'
Alarm78 = '"DB_Alarm_HMI"."Alarm78"'
Alarm79 = '"DB_Alarm_HMI"."Alarm79"'
Alarm8 = '"DB_Alarm_HMI"."Alarm8"'
Alarm80 = '"DB_Alarm_HMI"."Alarm80"'
Alarm81 = '"DB_Alarm_HMI"."Alarm81"'
Alarm82 = '"DB_Alarm_HMI"."Alarm82"'
Alarm83 = '"DB_Alarm_HMI"."Alarm83"'
Alarm84 = '"DB_Alarm_HMI"."Alarm84"'
Alarm85 = '"DB_Alarm_HMI"."Alarm85"'
Alarm86 = '"DB_Alarm_HMI"."Alarm86"'
Alarm87 = '"DB_Alarm_HMI"."Alarm87"'
Alarm88 = '"DB_Alarm_HMI"."Alarm88"'
Alarm89 = '"DB_Alarm_HMI"."Alarm89"'
Alarm9 = '"DB_Alarm_HMI"."Alarm9"'
Alarm90 = '"DB_Alarm_HMI"."Alarm90"'
Alarm91 = '"DB_Alarm_HMI"."Alarm91"'
Alarm92 = '"DB_Alarm_HMI"."Alarm92"'
Alarm93 = '"DB_Alarm_HMI"."Alarm93"'
Alarm94 = '"DB_Alarm_HMI"."Alarm94"'
Alarm95 = '"DB_Alarm_HMI"."Alarm95"'
Alarm96 = '"DB_Alarm_HMI"."Alarm96"'
Alarm97 = '"DB_Alarm_HMI"."Alarm97"'
Alarm98 = '"DB_Alarm_HMI"."Alarm98"'
```

```
Alarm99 = '"DB_Alarm_HMI"."Alarm99"'  
  
class hvl_ccb.dev.supercube.constants.BreakdownDetection(*args, **kwds)  
Bases: hvl_ccb.utils.enum.ValueEnum  
  
Node ID strings for the breakdown detection.  
  
TODO: these variable NodeIDs are not tested and/or correct yet.  
  
activated = '"Ix_Allg_Breakdown_activated"'  
Boolean read-only variable indicating whether breakdown detection and fast switchoff is enabled in the system or not.  
  
reset = '"Qx_Allg_Breakdown_reset"'  
Boolean writable variable to reset the fast switch-off. Toggle to re-enable.  
  
triggered = '"Ix_Allg_Breakdown_triggered"'  
Boolean read-only variable telling whether the fast switch-off has triggered. This can also be seen using the safety circuit state, therefore no method is implemented to read this out directly.  
  
class hvl_ccb.dev.supercube.constants.Door(*args, **kwds)  
Bases: hvl_ccb.dev.supercube.constants._DoorEnumBase  
  
Variable NodeID strings for doors.  
  
status_1 = '"DB_Safety_Circuit"."Door_1"."si_HMI_status"'  
status_2 = '"DB_Safety_Circuit"."Door_2"."si_HMI_status"'  
status_3 = '"DB_Safety_Circuit"."Door_3"."si_HMI_status"'  
  
class hvl_ccb.dev.supercube.constants.DoorStatus(*args, **kwds)  
Bases: aenum.IntEnum  
  
Possible status values for doors.  
  
closed = 2  
Door is closed, but not locked.  
  
error = 4  
Door has an error or was opened in locked state (either with emergency stop or from the inside).  
  
inactive = 0  
not enabled in Supercube HMI setup, this door is not supervised.  
  
locked = 3  
Door is closed and locked (safe state).  
  
open = 1  
Door is open.  
  
class hvl_ccb.dev.supercube.constants.EarthingStick(*args, **kwds)  
Bases: hvl_ccb.utils.enum.ValueEnum  
  
Variable NodeID strings for all earthing stick statuses (read-only integer) and writable booleans for setting the earthing in manual mode.  
  
manual = <bound method EarthingStick.manual of <aenum 'EarthingStick'>>  
manual_1 = '"DB_Safety_Circuit"."Earthstick_1"."sx_earthing_manually"'  
manual_2 = '"DB_Safety_Circuit"."Earthstick_2"."sx_earthing_manually"'  
manual_3 = '"DB_Safety_Circuit"."Earthstick_3"."sx_earthing_manually"'  
manual_4 = '"DB_Safety_Circuit"."Earthstick_4"."sx_earthing_manually"'
```

```

manual_5 = '"DB_Safety_Circuit"."Earthstick_5"."sx_earthing_manually"'
manual_6 = '"DB_Safety_Circuit"."Earthstick_6"."sx_earthing_manually"'
manuals = <bound method EarthingStick.manuals of <aenum 'EarthingStick'>>
number
    Get corresponding earthing stick number.

    Returns earthing stick number (1..6)

operating_status = <bound method EarthingStick.operating_status of <aenum 'EarthingStick'>>
operating_status_1 = '"DB_Safety_Circuit"."Earthstick_1"."sx_manual_control_active"'
operating_status_2 = '"DB_Safety_Circuit"."Earthstick_2"."sx_manual_control_active"'
operating_status_3 = '"DB_Safety_Circuit"."Earthstick_3"."sx_manual_control_active"'
operating_status_4 = '"DB_Safety_Circuit"."Earthstick_4"."sx_manual_control_active"'
operating_status_5 = '"DB_Safety_Circuit"."Earthstick_5"."sx_manual_control_active"'
operating_status_6 = '"DB_Safety_Circuit"."Earthstick_6"."sx_manual_control_active"'
operating_statuses = <bound method EarthingStick.operating_statuses of <aenum 'EarthingStick'>>
range = <bound method EarthingStick.range of <aenum 'EarthingStick'>>
status = <bound method EarthingStick.status of <aenum 'EarthingStick'>>
status_1 = '"DB_Safety_Circuit"."Earthstick_1"."si_HMI_Status"'
status_2 = '"DB_Safety_Circuit"."Earthstick_2"."si_HMI_Status"'
status_3 = '"DB_Safety_Circuit"."Earthstick_3"."si_HMI_Status"'
status_4 = '"DB_Safety_Circuit"."Earthstick_4"."si_HMI_Status"'
status_5 = '"DB_Safety_Circuit"."Earthstick_5"."si_HMI_Status"'
status_6 = '"DB_Safety_Circuit"."Earthstick_6"."si_HMI_Status"'
statuses = <bound method EarthingStick.statuses of <aenum 'EarthingStick'>>

class hvl_ccb.dev.supercube.constants.EarthstickMeta(*args, **kwds)
    Bases: aenum.EnumMeta

class hvl_ccb.dev.supercube.constants.EarthstickOperatingStatus(*args,
                                                               **kwds)
    Bases: aenum.IntEnum

    Operating Status for an earthing stick. Stick can be used in auto or manual mode.

    auto = 0
    manual = 1

class hvl_ccb.dev.supercube.constants.EarthstickOperation(*args, **kwds)
    Bases: aenum.IntEnum

    Operation of the earthing stick in manual operating mode. Can be closed or opened.

    close = 1
    open = 0

```

```
class hvl_ccb.dev.supercube.constants.EarthingStickStatus(*args, **kwds)
```

Bases: aenum.IntEnum

Status of an earthing stick. These are the possible values in the status integer e.g. in `EarthingStick.status_1`.

```
closed = 1
```

Earthing is closed (safe).

```
error = 3
```

Earthing is in error, e.g. when the stick did not close correctly or could not open.

```
inactive = 0
```

Earthing stick is deselected and not enabled in safety circuit. To get out of this state, the earthing has to be enabled in the Supercube HMI setup.

```
open = 2
```

Earthing is open (not safe).

```
class hvl_ccb.dev.supercube.constants.Errors(*args, **kwds)
```

Bases: `hvl_ccb.utils.enum.ValueEnum`

Variable NodeID strings for information regarding error, warning and message handling.

```
message = '"DB_Message_Buffer"."Info_active"'
```

Boolean read-only variable telling if a message is active.

```
quit = '"DB_Message_Buffer"."Reset_button"'
```

Writable boolean for the error quit button.

```
stop = '"DB_Message_Buffer"."Stop_active"'
```

Boolean read-only variable telling if a stop is active.

```
warning = '"DB_Message_Buffer"."Warning_active"'
```

Boolean read-only variable telling if a warning is active.

```
class hvl_ccb.dev.supercube.constants.GeneralSockets(*args, **kwds)
```

Bases: `hvl_ccb.utils.enum.ValueEnum`

NodeID strings for the power sockets (3x T13 and 1xCEE16).

```
cee16 = '"Qx_Allg_Socket_CEE16"'
```

CEE16 socket (writeable boolean).

```
t13_1 = '"Qx_Allg_Socket_T13_1"'
```

SEV T13 socket No. 1 (writable boolean).

```
t13_2 = '"Qx_Allg_Socket_T13_2"'
```

SEV T13 socket No. 2 (writable boolean).

```
t13_3 = '"Qx_Allg_Socket_T13_3"'
```

SEV T13 socket No. 3 (writable boolean).

```
class hvl_ccb.dev.supercube.constants.GeneralSupport(*args, **kwds)
```

Bases: `hvl_ccb.utils.enum.ValueEnum`

NodeID strings for the support inputs and outputs.

```
contact_range = <bound method GeneralSupport.contact_range of <aenum 'GeneralSupport'>
```

```
in_1_1 = '"Ix_Allg_Support1_1"'
```

```
in_1_2 = '"Ix_Allg_Support1_2"'
```

```
in_2_1 = '"Ix_Allg_Support2_1"'
```

```

in_2_2 = '"Ix_Allg_Support2_2"'
in_3_1 = '"Ix_Allg_Support3_1"'
in_3_2 = '"Ix_Allg_Support3_2"'
in_4_1 = '"Ix_Allg_Support4_1"'
in_4_2 = '"Ix_Allg_Support4_2"'
in_5_1 = '"Ix_Allg_Support5_1"'
in_5_2 = '"Ix_Allg_Support5_2"'
in_6_1 = '"Ix_Allg_Support6_1"'
in_6_2 = '"Ix_Allg_Support6_2"'

input = <bound method GeneralSupport.input of <aenum 'GeneralSupport'>>
out_1_1 = '"Qx_Allg_Support1_1"'
out_1_2 = '"Qx_Allg_Support1_2"'
out_2_1 = '"Qx_Allg_Support2_1"'
out_2_2 = '"Qx_Allg_Support2_2"'
out_3_1 = '"Qx_Allg_Support3_1"'
out_3_2 = '"Qx_Allg_Support3_2"'
out_4_1 = '"Qx_Allg_Support4_1"'
out_4_2 = '"Qx_Allg_Support4_2"'
out_5_1 = '"Qx_Allg_Support5_1"'
out_5_2 = '"Qx_Allg_Support5_2"'
out_6_1 = '"Qx_Allg_Support6_1"'
out_6_2 = '"Qx_Allg_Support6_2"'

output = <bound method GeneralSupport.output of <aenum 'GeneralSupport'>>
port_range = <bound method GeneralSupport.port_range of <aenum 'GeneralSupport'>>

class hvl_ccb.dev.supercube.constants.GeneralSupportMeta(*args, **kwds)
Bases: aenum.EnumMeta

class hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio(*args, **kwds)
Bases: hvl_ccb.dev.supercube.constants._InputEnumBase

Variable NodeID strings for the measurement input scaling ratios. These ratios are defined in the Supercube HMI setup and are provided in the python module here to be able to read them out, allowing further calculations.

input_1 = "DB_Measurements"."si_Divider_Ratio_1"
input_2 = "DB_Measurements"."si_Divider_Ratio_2"
input_3 = "DB_Measurements"."si_Divider_Ratio_3"
input_4 = "DB_Measurements"."si_Divider_Ratio_4"

class hvl_ccb.dev.supercube.constants.MeasurementsScaledInput(*args, **kwds)
Bases: hvl_ccb.dev.supercube.constants._InputEnumBase

```

Variable NodeID strings for the four analog BNC inputs for measuring voltage. The voltage returned in these variables is already scaled with the set ratio, which can be read using the variables in [MeasurementsDividerRatio](#).

```
input_1 = '"DB_Measurements"."si_scaled_Voltage_Input_1"'
input_2 = '"DB_Measurements"."si_scaled_Voltage_Input_2"'
input_3 = '"DB_Measurements"."si_scaled_Voltage_Input_3"'
input_4 = '"DB_Measurements"."si_scaled_Voltage_Input_4"'

class hvl_ccb.dev.supercube.constants.MessageBoard(*args, **kwds)
Bases: hvl_ccb.dev.supercube.constants._LineEnumBase
```

Variable NodeID strings for message board lines.

```
line_1 = '"DB_OPC_Connection"."Is_status_Line_1"'
line_10 = '"DB_OPC_Connection"."Is_status_Line_10"'
line_11 = '"DB_OPC_Connection"."Is_status_Line_11"'
line_12 = '"DB_OPC_Connection"."Is_status_Line_12"'
line_13 = '"DB_OPC_Connection"."Is_status_Line_13"'
line_14 = '"DB_OPC_Connection"."Is_status_Line_14"'
line_15 = '"DB_OPC_Connection"."Is_status_Line_15"'
line_2 = '"DB_OPC_Connection"."Is_status_Line_2"'
line_3 = '"DB_OPC_Connection"."Is_status_Line_3"'
line_4 = '"DB_OPC_Connection"."Is_status_Line_4"'
line_5 = '"DB_OPC_Connection"."Is_status_Line_5"'
line_6 = '"DB_OPC_Connection"."Is_status_Line_6"'
line_7 = '"DB_OPC_Connection"."Is_status_Line_7"'
line_8 = '"DB_OPC_Connection"."Is_status_Line_8"'
line_9 = '"DB_OPC_Connection"."Is_status_Line_9"'

class hvl_ccb.dev.supercube.constants.OpcControl(*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum
```

Variable NodeID strings for supervision of the OPC connection from the controlling workstation to the Supercube.

```
active = '"DB_OPC_Connection"."sx_OPC_active"'
writable boolean to enable OPC remote control and display a message window on the Supercube HMI.
```

```
live = '"DB_OPC_Connection"."sx_OPC_lifebit"'
```

```
class hvl_ccb.dev.supercube.constants.Power(*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum
```

Variable NodeID strings concerning power data.

**TODO:** these variable NodeIDs are not tested and/or correct yet, they don't exist yet on Supercube side.

```
current_primary = 'Qr_Power_FU_actual_Current'
Primary current in ampere, measured by the frequency converter. (read-only)
```

```

frequency = 'Ir_Power_FU_Frequency'
    Frequency converter output frequency. (read-only)

setup = 'Qi_Power_Setup'
    Power setup that is configured using the Supercube HMI. The value corresponds to the ones in
    PowerSetup. (read-only)

voltage_max = 'Iw_Power_max_Voltage'
    Maximum voltage allowed by the current experimental setup. (read-only)

voltage_primary = 'Or_Power_FU_actual_Voltage'
    Primary voltage in volts, measured by the frequency converter at its output. (read-only)

voltage_slope = 'Ir_Power_dUDt'
    Voltage slope in V/s.

voltage_target = 'Ir_Power_Target_Voltage'
    Target voltage setpoint in V.

class hvl_ccb.dev.supercube.constants.PowerSetup(*args, **kwds)
Bases: aenum.IntEnum

Possible power setups corresponding to the value of variable Power.setup.

AC_DoubleStage_150kV = 4
    AC voltage with two MWB transformers, one at 100kV and the other at 50kV, resulting in a total maximum
    voltage of 150kV.

AC_DoubleStage_200kV = 5
    AC voltage with two MWB transformers both at 100kV, resulting in a total maximum voltage of 200kV

AC_SingleStage_100kV = 3
    AC voltage with MWB transformer set to 100kV maximum voltage.

AC_SingleStage_50kV = 2
    AC voltage with MWB transformer set to 50kV maximum voltage.

DC_DoubleStage_280kV = 8
    DC voltage with two AC transformers set to 100kV AC each, resulting in 280kV DC in total (or a single
    stage transformer with Greinacher voltage doubling rectifier)

DC_SingleStage_140kV = 7
    DC voltage with one AC transformer set to 100kV AC, resulting in 140kV DC

External = 1
    External power supply fed through blue CEE32 input using isolation transformer and safety switches of
    the Supercube, or using an external safety switch attached to the Supercube Type B.

Internal = 6
    Internal usage of the frequency converter, controlling to the primary voltage output of the supercube itself
    (no measurement transformer used)

NoPower = 0
    No safety switches, use only safety components (doors, fence, earthing...) without any power.

class hvl_ccb.dev.supercube.constants.Safety(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

NodeID strings for the basic safety circuit status and green/red switches “ready” and “operate”.

status = '"DB_Safety_Circuit"."si_safe_status"'
    Status is a read-only integer containing the state number of the supercube-internal state machine. The
    values correspond to numbers in SafetyStatus.

```

```
switch_to_operate = '"DB_Safety_Circuit"."sx_safe_switch_to_operate"'
    Writable boolean for switching to Red Operate (locket, HV on) state.

switch_to_ready = '"DB_Safety_Circuit"."sx_safe_switch_to_ready"'
    Writable boolean for switching to Red Ready (locked, HV off) state.

class hvl_ccb.dev.supercube.constants.SafetyStatus(*args, **kwargs)
Bases: aenum.IntEnum

Safety status values that are possible states returned from hvl_ccb.dev.supercube.base.Supercube.get_status(). These values correspond to the states of the Supercube's safety circuit statemachine.

Error = 6
    System is in error mode.

GreenNotReady = 1
    System is safe, lamps are green and some safety elements are not in place such that it cannot be switched to red currently.

GreenReady = 2
    System is safe and all safety elements are in place to be able to switch to ready.

Initializing = 0
    System is initializing or booting.

QuickStop = 5
    Fast turn off triggered and switched off the system. Reset FSO to go back to a normal state.

RedOperate = 4
    System is locked in red state and in operate mode, i.e. high voltage on.

RedReady = 3
    System is locked in red state and ready to go to operate mode.

class hvl_ccb.dev.supercube.constants.SupercubeOpcEndpoint(*args, **kwargs)
Bases: hvl_ccb.utils.enum.ValueEnum

OPC Server Endpoint strings for the supercube variants.

A = 'Supercube Typ A'
B = 'Supercube Typ B'

hvl_ccb.dev.supercube.constants.T13_SOCKET_PORTS = (1, 2, 3)
Port numbers of SEV T13 power socket
```

## [hvl\\_ccb.dev.supercube.typ\\_a module](#)

Supercube Typ A module.

```
class hvl_ccb.dev.supercube.typ_a.SupercubeAOpcUaCommunication(config)
Bases: hvl_ccb.dev.supercube.base.SupercubeOpcUaCommunication

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class
```

```
class hvl_ccb.dev.supercube.typ_a.SupercubeAopcUaConfiguration(host: str, endpoint_name: str = 'Supercube Typ A', port: int = 4840, sub_handler: hvl_ccb.comm.opc.OpcUaSubHandler = <hvl_ccb.dev.supercube.base.SupercubeS>, object at 0x7f155c83f8d0>, update_period: int = 500, wait_timeout_retry_sec: Union[int, float] = 1, max_timeout_retry_nr: int = 5)
```

Bases: *hvl\_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig*

**endpoint\_name** = 'Supercube Typ A'

**force\_value** (*fieldname*, *value*)  
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]  
Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
class hvl_ccb.dev.supercube.typ_a.SupercubeWithFU(com, dev_config=None)
```

Bases: *hvl\_ccb.dev.supercube.base.SupercubeBase*

Variant A of the Supercube with frequency converter.

**static default\_com\_cls()**  
Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**fso\_reset()** → None  
TODO: test **fso\_reset** with device

Reset the fast switch off circuitry to go back into normal state and allow to re-enable operate mode.

**get\_frequency()** → float

**TODO: test get\_frequency with device**

Read the electrical frequency of the current Supercube setup.

**Returns** the frequency in Hz

**get\_fso\_active()** → bool

**TODO: test get\_fso\_active with device**

Get the state of the fast switch off functionality. Returns True if it is enabled, False otherwise.

**Returns** state of the FSO functionality

**get\_max\_voltage()** → float

**TODO: test get\_max\_voltage with device**

Reads the maximum voltage of the setup and returns in V.

**Returns** the maximum voltage of the setup in V.

**get\_power\_setup()** → hvl\_ccb.dev.supercube.constants.PowerSetup

**TODO: test get\_power\_setup with device**

Return the power setup selected in the Supercube's settings.

**Returns** the power setup

**get\_primary\_current()** → float

**TODO: get\_primary\_current with device**

Read the current primary current at the output of the frequency converter ( before transformer).

**Returns** primary current in A

**get\_primary\_voltage()** → float

**TODO: test get\_primary\_voltage with device**

Read the current primary voltage at the output of the frequency converter ( before transformer).

**Returns** primary voltage in V

**get\_target\_voltage()** → float

**TODO: test get\_target\_voltage with device**

Gets the current setpoint of the output voltage value in V. This is not a measured value but is the corresponding function to [set\\_target\\_voltage\(\)](#).

**Returns** the setpoint voltage in V.

**set\_slope(slope: float)** → None

**TODO: test set\_slope with device**

Sets the dV/dt slope of the Supercube frequency converter to a new value in V/s.

**Parameters** **slope** – voltage slope in V/s (0..15)

**set\_target\_voltage(volt\_v: float)** → None

**TODO: test set\_target\_voltage with device**

Set the output voltage to a defined value in V.

**Parameters** **volt\_v** – the desired voltage in V

## `hvl_ccb.dev.supercube.typ_b module`

Supercube Typ B module.

**class** `hvl_ccb.dev.supercube.typ_b.SuperCubeB`(*com, dev\_config=None*)

Bases: `hvl_ccb.dev.supercube.base.SuperCubeBase`

Variant B of the Supercube without frequency converter but external safety switches.

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**class** `hvl_ccb.dev.supercube.typ_b.SuperCubeBOpcUaCommunication`(*config*)

Bases: `hvl_ccb.dev.supercube.base.SuperCubeOpcUaCommunication`

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class** `hvl_ccb.dev.supercube.typ_b.SuperCubeBOpcUaConfiguration`(*host: str, endpoint\_name: str*

*= 'Supercube*

*Typ B', port:*

*int = 4840,*

*sub\_handler:*

*hvl\_ccb.comm.opc.OpcUaSubHandler*

*=*

*<hvl\_ccb.dev.supercube.base.SuperCube*

*object at*

*0x7f155c83f8d0>,*

*update\_period:*

*int = 500,*

*wait\_timeout\_retry\_sec:*

*Union[int,*

*float] = 1,*

*max\_timeout\_retry\_nr:*

*int = 5)*

Bases: `hvl_ccb.dev.supercube.base.SuperCubeOpcUaCommunicationConfig`

**endpoint\_name = 'Supercube Typ B'**

**force\_value**(*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## Module contents

Supercube package with implementation for system versions from 2019 on (new concept with hard-PLC Siemens S7-1500 as CPU).

### [hvl\\_ccb.dev.supercube2015 package](#)

#### Submodules

##### [hvl\\_ccb.dev.supercube2015.base module](#)

Base classes for the Supercube device.

```
exception hvl_ccb.dev.supercube2015.base.InvalidSupercubeStatusError
```

Bases: Exception

Exception raised when supercube has invalid status.

```
class hvl_ccb.dev.supercube2015.base.Supercube2015Base(com, dev_config=None)
```

Bases: [hvl\\_ccb.dev.base.SingleCommDevice](#)

Base class for Supercube variants.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
static default_com_cls()
```

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

```
get_cee16_socket() → bool
```

Read the on-state of the IEC CEE16 three-phase power socket.

**Returns** the on-state of the CEE16 power socket

```
get_door_status(door: int) → hvl_ccb.dev.supercube2015.constants.DoorStatus
```

Get the status of a safety fence door. See `constants.DoorStatus` for possible returned door statuses.

**Parameters** `door` – the door number (1..3)

**Returns** the door status

```
get_earthing_manual(number: int) → bool
```

Get the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters** **number** – number of the earthing stick (1..6)

**Returns** earthing stick manual status

**get\_earthing\_status** (*number: int*) → int

Get the status of an earthing stick, whether it is closed, open or undefined (moving).

**Parameters** **number** – number of the earthing stick (1..6)

**Returns** earthing stick status; see constants.EarthlingStickStatus

**get\_measurement\_ratio** (*channel: int*) → float

Get the set measurement ratio of an AC/DC analog input channel. Every input channel has a divider ratio assigned during setup of the Supercube system. This ratio can be read out.

**Attention:** Supercube 2015 does not have a separate ratio for every analog input. Therefore there is only one ratio for channel = 1.

**Parameters** **channel** – number of the input channel (1..4)

**Returns** the ratio

**get\_measurement\_voltage** (*channel: int*) → float

Get the measured voltage of an analog input channel. The voltage read out here is already scaled by the configured divider ratio.

**Attention:** In contrast to the *new* Supercube, the old one returns here the input voltage read at the ADC. It is not scaled by a factor.

**Parameters** **channel** – number of the input channel (1..4)

**Returns** measured voltage

**get\_status** () → int

Get the safety circuit status of the Supercube.

**Returns** the safety status of the supercube's state machine; see *constants.SafetyStatus*.

**get\_support\_input** (*port: int, contact: int*) → bool

Get the state of a support socket input.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital input read state

**get\_support\_output** (*port: int, contact: int*) → bool

Get the state of a support socket output.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)

**Returns** digital output read state

**get\_t13\_socket** (*port: int*) → bool

Read the state of a SEV T13 power socket.

**Parameters** **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*

**Returns** on-state of the power socket

**horn** (*state: bool*) → None

Turns acoustic horn on or off.

**Parameters state** – Turns horn on (True) or off (False)

**operate** (*state: bool*) → None

Set operate state. If the state is RedReady, this will turn on the high voltage and close the safety switches.

**Parameters state** – set operate state

**quit\_error** () → None

Quits errors that are active on the Supercube.

**read** (*node\_id: str*)

Local wrapper for the OPC UA communication protocol read method.

**Parameters node\_id** – the id of the node to read.

**Returns** the value of the variable

**ready** (*state: bool*) → None

Set ready state. Ready means locket safety circuit, red lamps, but high voltage still off.

**Parameters state** – set ready state

**set\_cee16\_socket** (*state: bool*) → None

Switch the IEC CEE16 three-phase power socket on or off.

**Parameters state** – desired on-state of the power socket

**Raises ValueError** – if state is not of type bool

**set\_earthing\_manual** (*number: int, manual: bool*) → None

Set the manual status of an earthing stick. If an earthing stick is set to manual, it is closed even if the system is in states RedReady or RedOperate.

**Parameters**

- **number** – number of the earthing stick (1..6)
- **manual** – earthing stick manual status (True or False)

**set\_remote\_control** (*state: bool*) → None

Enable or disable remote control for the Supercube. This will effectively display a message on the touch-screen HMI.

**Parameters state** – desired remote control state

**set\_support\_output** (*port: int, contact: int, state: bool*) → None

Set the state of a support output socket.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **state** – is the desired state of the support output

**set\_support\_output\_impulse** (*port: int, contact: int, duration: float = 0.2, pos\_pulse: bool = True*) → None

Issue an impulse of a certain duration on a support output contact. The polarity of the pulse (On-wait-Off or Off-wait-On) is specified by the pos\_pulse argument.

This function is blocking.

**Parameters**

- **port** – is the socket number (1..6)
- **contact** – is the contact on the socket (1..2)
- **duration** – is the length of the impulse in seconds
- **pos\_pulse** – is True, if the pulse shall be HIGH, False if it shall be LOW

**set\_t13\_socket** (*port: int, state: bool*) → None

Set the state of a SEV T13 power socket.

#### Parameters

- **port** – is the socket number, one of *constants.T13\_SOCKET\_PORTS*
- **state** – is the desired on-state of the socket

**start()** → None

Starts the device. Sets the root node for all OPC read and write commands to the Siemens PLC object node which holds all our relevant objects and variables.

**stop()** → None

Stop the Supercube device. Deactivates the remote control and closes the communication protocol.

**write** (*node\_id, value*) → None

Local wrapper for the OPC UA communication protocol write method.

#### Parameters

- **node\_id** – the id of the node to read
- **value** – the value to write to the variable

**class** hvl\_ccb.dev.supercube2015.base.**SupercubeConfiguration** (*namespace\_index: int = 7*)

Bases: object

Configuration dataclass for the Supercube devices.

**clean\_values()**

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**force\_value** (*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**namespace\_index = 7**

Namespace of the OPC variables, typically this is 3 (coming from Siemens)

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**class** hvl\_ccb.dev.supercube2015.base.**SupercubeOpcUaCommunication**(config)

Bases: *hvl\_ccb.comm.opc.OpcUaCommunication*

Communication protocol specification for Supercube devices.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class** hvl\_ccb.dev.supercube2015.base.**SupercubeOpcUaCommunicationConfig**(host:

str,

end-

point\_name:

str,

port:

int

=

4845,

sub\_handler:

*hvl\_ccb.comm.opc.OpcUaSu*

=

<*hvl\_ccb.dev.supercube2015*

ob-

ject>,

up-

date\_period:

int

=

500,

wait\_timeout\_retry\_sec:

Union[int,

float]

= 1,

max\_timeout\_retry\_nr:

int

=

5)

Bases: *hvl\_ccb.comm.opc.OpcUaCommunicationConfig*

Communication protocol configuration for OPC UA, specifications for the Supercube devices.

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys () → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.**classmethod optional\_defaults () → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.**port = 4845****classmethod required\_keys () → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.**sub\_handler = <hvl\_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler object>**  
Subscription handler for data change events**class hvl\_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler**Bases: *hvl\_ccb.comm.opc.OpcUaSubHandler*

OPC Subscription handler for datachange events and normal events specifically implemented for the Supercube devices.

**datachange\_notification (node: opcua.common.node.Node, val, data)**

In addition to the standard operation (debug logging entry of the datachange), alarms are logged at INFO level using the alarm text.

**Parameters**

- **node** – the node object that triggered the datachange event
- **val** – the new value
- **data** –

**hvl\_ccb.dev.supercube2015.constants module**

Constants, variable names for the Supercube OPC-connected devices.

**class hvl\_ccb.dev.supercube2015.constants.AlarmText (\*args, \*\*kwds)**  
Bases: *hvl\_ccb.utils.enum.ValueEnum*This enumeration contains textual representations for all error classes (stop, warning and message) of the Supercube system. Use the *AlarmText.get()* method to retrieve the enum of an alarm number.

```
Alarm0 = 'No Alarm.'

Alarm1 = 'STOP Safety switch 1 error'
Alarm10 = 'STOP Earthing stick 2 error'
Alarm11 = 'STOP Earthing stick 3 error'
Alarm12 = 'STOP Earthing stick 4 error'
```

```
Alarm13 = 'STOP Earthing stick 5 error'
Alarm14 = 'STOP Earthing stick 6 error'
Alarm17 = 'STOP Source switch error'
Alarm19 = 'STOP Fence 1 error'
Alarm2 = 'STOP Safety switch 2 error'
Alarm20 = 'STOP Fence 2 error'
Alarm21 = 'STOP Control error'
Alarm22 = 'STOP Power outage'
Alarm3 = 'STOP Emergency Stop 1'
Alarm4 = 'STOP Emergency Stop 2'
Alarm5 = 'STOP Emergency Stop 3'
Alarm6 = 'STOP Door 1 lock supervision'
Alarm7 = 'STOP Door 2 lock supervision'
Alarm8 = 'STOP Door 3 lock supervision'
Alarm9 = 'STOP Earthing stick 1 error'
get = <bound method AlarmText.get of <aenum 'AlarmText'>>
not_defined = 'NO ALARM TEXT DEFINED'

class hvl_ccb.dev.supercube2015.constants.BreakdownDetection(*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum

Node ID strings for the breakdown detection.

activated = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.Breakdowndetection.connect'
    Boolean read-only variable indicating whether breakdown detection and fast switchoff is enabled in the system or not.

reset = 'hvl-ipc.WINAC.Support6OutA'
    Boolean writable variable to reset the fast switch-off. Toggle to re-enable.

triggered = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.Breakdowndetection.triggered'
    Boolean read-only variable telling whether the fast switch-off has triggered. This can also be seen using the safety circuit state, therefore no method is implemented to read this out directly.

class hvl_ccb.dev.supercube2015.constants.DoorStatus(*args, **kwds)
Bases: aenum.IntEnum

Possible status values for doors.

closed = 2
    Door is closed, but not locked.

error = 4
    Door has an error or was opened in locked state (either with emergency stop or from the inside).

inactive = 0
    not enabled in Supercube HMI setup, this door is not supervised.

locked = 3
    Door is closed and locked (safe state).
```

```

open = 1
Door is open.

class hvl_ccb.dev.supercube2015.constants.EarthingStick(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

Variable NodeID strings for all earthing stick statuses (read-only integer) and writable booleans for setting the earthing in manual mode.

manual = <bound method EarthingStick.manual of <aenum 'EarthingStick'>>
manual_1 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.MANUAL'
manual_2 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.MANUAL'
manual_3 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.MANUAL'
manual_4 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.MANUAL'
manual_5 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.MANUAL'
manual_6 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.MANUAL'
status_1_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.CLOSE'
status_1_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.CONNECT'
status_1_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_1.OPEN'
status_2_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.CLOSE'
status_2_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.CONNECT'
status_2_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_2.OPEN'
status_3_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.CLOSE'
status_3_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.CONNECT'
status_3_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_3.OPEN'
status_4_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.CLOSE'
status_4_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.CONNECT'
status_4_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_4.OPEN'
status_5_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.CLOSE'
status_5_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.CONNECT'
status_5_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_5.OPEN'
status_6_closed = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.CLOSE'
status_6_connected = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.CONNECT'
status_6_open = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.STICK_6.OPEN'
status_closed = <bound method EarthingStick.status_closed of <aenum 'EarthingStick'>>
status_connected = <bound method EarthingStick.status_connected of <aenum 'EarthingStick'>>
status_open = <bound method EarthingStick.status_open of <aenum 'EarthingStick'>>

class hvl_ccb.dev.supercube2015.constants.EarthingStickStatus(*args, **kwds)
Bases: aenum.IntEnum

```

Status of an earthing stick. These are the possible values in the status integer e.g. in EarthingStick. status\_1.

```
closed = 1
Earthing is closed (safe).

error = 3
Earthing is in error, e.g. when the stick did not close correctly or could not open.

inactive = 0
Earthing stick is deselected and not enabled in safety circuit. To get out of this state, the earthing has to be
enabled in the Supercube HMI setup.

open = 2
Earthing is open (not safe).

class hvl_ccb.dev.supercube2015.constants.Errors (*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum

Variable NodeID strings for information regarding error, warning and message handling.

quit = 'hvl-ipc.WINAC.SYSTEMSTATE.Faultconfirmation'
Writable boolean for the error quit button.

stop = 'hvl-ipc.WINAC.SYSTEMSTATE.ERROR'
Boolean read-only variable telling if a stop is active.

stop_number = 'hvl-ipc.WINAC.SYSTEMSTATE.Errornumber'

class hvl_ccb.dev.supercube2015.constants.GeneralSockets (*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum

NodeID strings for the power sockets (3x T13 and 1xCEE16).

cee16 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.CEE16'
CEE16 socket (writeable boolean).

t13_1 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.T13_1'
SEV T13 socket No. 1 (writable boolean).

t13_2 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.T13_2'
SEV T13 socket No. 2 (writable boolean).

t13_3 = 'hvl-ipc.WINAC.SYSTEM_COMPONENTS.T13_3'
SEV T13 socket No. 3 (writable boolean).

class hvl_ccb.dev.supercube2015.constants.GeneralSupport (*args, **kwds)
Bases: hvl\_ccb.utils.enum.ValueEnum

NodeID strings for the support inputs and outputs.

in_1_1 = 'hvl-ipc.WINAC.Support1InA'
in_1_2 = 'hvl-ipc.WINAC.Support1InB'
in_2_1 = 'hvl-ipc.WINAC.Support2InA'
in_2_2 = 'hvl-ipc.WINAC.Support2InB'
in_3_1 = 'hvl-ipc.WINAC.Support3InA'
in_3_2 = 'hvl-ipc.WINAC.Support3InB'
in_4_1 = 'hvl-ipc.WINAC.Support4InA'
in_4_2 = 'hvl-ipc.WINAC.Support4InB'
in_5_1 = 'hvl-ipc.WINAC.Support5InA'
in_5_2 = 'hvl-ipc.WINAC.Support5InB'
```

```

in_6_1 = 'hvl-ipc.WINAC.Support6InA'
in_6_2 = 'hvl-ipc.WINAC.Support6InB'
input = <bound method GeneralSupport.input of <aenum 'GeneralSupport'>>
out_1_1 = 'hvl-ipc.WINAC.Support1OutA'
out_1_2 = 'hvl-ipc.WINAC.Support1OutB'
out_2_1 = 'hvl-ipc.WINAC.Support2OutA'
out_2_2 = 'hvl-ipc.WINAC.Support2OutB'
out_3_1 = 'hvl-ipc.WINAC.Support3OutA'
out_3_2 = 'hvl-ipc.WINAC.Support3OutB'
out_4_1 = 'hvl-ipc.WINAC.Support4OutA'
out_4_2 = 'hvl-ipc.WINAC.Support4OutB'
out_5_1 = 'hvl-ipc.WINAC.Support5OutA'
out_5_2 = 'hvl-ipc.WINAC.Support5OutB'
out_6_1 = 'hvl-ipc.WINAC.Support6OutA'
out_6_2 = 'hvl-ipc.WINAC.Support6OutB'
output = <bound method GeneralSupport.output of <aenum 'GeneralSupport'>>

```

```
class hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio(*args,
                                                               **kwds)
```

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Variable NodeID strings for the measurement input scaling ratios. These ratios are defined in the Supercube HMI setup and are provided in the python module here to be able to read them out, allowing further calculations.

```
get = <bound method MeasurementsDividerRatio.get of <aenum 'MeasurementsDividerRatio'>>
input_1 = 'hvl-ipc.WINAC.SYSTEM_INTERN.DivididerRatio'
```

```
class hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInput(*args,
                                                               **kwds)
```

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Variable NodeID strings for the four analog BNC inputs for measuring voltage. The voltage returned in these variables is already scaled with the set ratio, which can be read using the variables in [MeasurementsDividerRatio](#).

```
get = <bound method MeasurementsScaledInput.get of <aenum 'MeasurementsScaledInput'>>
input_1 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI1Volt'
input_2 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI2Volt'
input_3 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI3Volt'
input_4 = 'hvl-ipc.WINAC.SYSTEM_INTERN.AI4Volt'
```

```
class hvl_ccb.dev.supercube2015.constants.Power(*args, **kwds)
```

Bases: [hvl\\_ccb.utils.enum.ValueEnum](#)

Variable NodeID strings concerning power data.

```
current_primary = 'hvl-ipc.WINAC.SYSTEM_INTERN.FUCurrentprim'
```

Primary current in ampere, measured by the frequency converter. (read-only)

```
frequency = 'hvl-ipc.WINAC.FU.Frequency'
Frequency converter output frequency. (read-only)

setup = 'hvl-ipc.WINAC.FU.TrafoSetup'
Power setup that is configured using the Supercube HMI. The value corresponds to the ones in
PowerSetup. (read-only)

voltage_max = 'hvl-ipc.WINAC.FU.maxVoltagekV'
Maximum voltage allowed by the current experimental setup. (read-only)

voltage_primary = 'hvl-ipc.WINAC.SYSTEM_INTERN.FUVoltageprim'
Primary voltage in volts, measured by the frequency converter at its output. (read-only)

voltage_slope = 'hvl-ipc.WINAC.FU.dUdt_-1'
Voltage slope in V/s.

voltage_target = 'hvl-ipc.WINAC.FU.SOLL'
Target voltage setpoint in V.

class hvl_ccb.dev.supercube2015.constants.PowerSetup(*args, **kwds)
Bases: aenum.IntEnum

Possible power setups corresponding to the value of variable Power.setup.

AC_DoubleStage_150kV = 3
AC voltage with two MWB transformers, one at 100kV and the other at 50kV, resulting in a total maximum
voltage of 150kV.

AC_DoubleStage_200kV = 4
AC voltage with two MWB transformers both at 100kV, resulting in a total maximum voltage of 200kV

AC_SingleStage_100kV = 2
AC voltage with MWB transformer set to 100kV maximum voltage.

AC_SingleStage_50kV = 1
AC voltage with MWB transformer set to 50kV maximum voltage.

DC_DoubleStage_280kV = 7
DC voltage with two AC transformers set to 100kV AC each, resulting in 280kV DC in total (or a single
stage transformer with Greinacher voltage doubling rectifier)

DC_SingleStage_140kV = 6
DC voltage with one AC transformer set to 100kV AC, resulting in 140kV DC

External = 0
External power supply fed through blue CEE32 input using isolation transformer and safety switches of
the Supercube, or using an external safety switch attached to the Supercube Type B.

Internal = 5
Internal usage of the frequency converter, controlling to the primary voltage output of the supercube itself
(no measurement transformer used)

class hvl_ccb.dev.supercube2015.constants.Safety(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

NodeID strings for the basic safety circuit status and green/red switches “ready” and “operate”.

horn = 'hvl-ipc.WINAC.SYSTEM_INTERN.hornen'
Writeable boolean to manually turn on or off the horn

status_error = 'hvl-ipc.WINAC.SYSTEMSTATE.ERROR'
status_green = 'hvl-ipc.WINAC.SYSTEMSTATE.GREEN'
```

```

status_ready_for_red = 'hvl-ipc.WINAC.SYSTEMSTATE.ReadyForRed'
    Status is a read-only integer containing the state number of the supercube-internal state machine. The values correspond to numbers in SafetyStatus.

status_red = 'hvl-ipc.WINAC.SYSTEMSTATE.RED'

switchto_green = 'hvl-ipc.WINAC.SYSTEMSTATE.GREEN_REQUEST'

switchto_operate = 'hvl-ipc.WINAC.SYSTEMSTATE.switchon'
    Writable boolean for switching to Red Operate (locket, HV on) state.

switchto_ready = 'hvl-ipc.WINAC.SYSTEMSTATE.RED_REQUEST'
    Writable boolean for switching to Red Ready (locked, HV off) state.

class hvl_ccb.dev.supercube2015.constants.SafetyStatus(*args, **kwds)
Bases: aenum.IntEnum

Safety status values that are possible states returned from hvl_ccb.dev.supercube.base.Supercube.get_status(). These values correspond to the states of the Supercube's safety circuit statemachine.

Error = 6
    System is in error mode.

GreenNotReady = 1
    System is safe, lamps are green and some safety elements are not in place such that it cannot be switched to red currently.

GreenReady = 2
    System is safe and all safety elements are in place to be able to switch to ready.

Initializing = 0
    System is initializing or booting.

QuickStop = 5
    Fast turn off triggered and switched off the system. Reset FSO to go back to a normal state.

RedOperate = 4
    System is locked in red state and in operate mode, i.e. high voltage on.

RedReady = 3
    System is locked in red state and ready to go to operate mode.

class hvl_ccb.dev.supercube2015.constants.SupercubeOpcEndpoint(*args, **kwds)
Bases: hvl_ccb.utils.enum.ValueEnum

OPC Server Endpoint strings for the supercube variants.

A = 'OPC.SimaticNET.S7'
B = 'OPC.SimaticNET.S7'

hvl_ccb.dev.supercube2015.constants.T13_SOCKET_PORTS = (1, 2, 3)
    Port numbers of SEV T13 power socket

```

## **hvl\_ccb.dev.supercube2015.typ\_a module**

Supercube Typ A module.

```

class hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU(com,
                                                               dev_config=None)
Bases: hvl_ccb.dev.supercube2015.base.Supercube2015Base

```

Variant A of the Supercube with frequency converter.

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**fso\_reset()** → None

Reset the fast switch off circuitry to go back into normal state and allow to re-enable operate mode.

**get\_frequency()** → float

Read the electrical frequency of the current Supercube setup.

**Returns** the frequency in Hz

**get\_fso\_active()** → bool

Get the state of the fast switch off functionality. Returns True if it is enabled, False otherwise.

**Returns** state of the FSO functionality

**get\_max\_voltage()** → float

Reads the maximum voltage of the setup and returns in V.

**Returns** the maximum voltage of the setup in V.

**get\_power\_setup()** → hvl\_ccb.dev.supercube2015.constants.PowerSetup

Return the power setup selected in the Supercube's settings.

**Returns** the power setup

**get\_primary\_current()** → float

Read the current primary current at the output of the frequency converter ( before transformer).

**Returns** primary current in A

**get\_primary\_voltage()** → float

Read the current primary voltage at the output of the frequency converter ( before transformer).

**Returns** primary voltage in V

**get\_target\_voltage()** → float

Gets the current setpoint of the output voltage value in V. This is not a measured value but is the corresponding function to [set\\_target\\_voltage\(\)](#).

**Returns** the setpoint voltage in V.

**set\_slope(slope: float)** → None

Sets the dV/dt slope of the Supercube frequency converter to a new value in V/s.

**Parameters** **slope** – voltage slope in V/s (0..15'000)

**set\_target\_voltage(volt\_v: float)** → None

Set the output voltage to a defined value in V.

**Parameters** **volt\_v** – the desired voltage in V

**class** hvl\_ccb.dev.supercube2015.typ\_a.**SupercubeAOpcUaCommunication**(config)

Bases: [hvl\\_ccb.dev.supercube2015.base.SupercubeOpcUaCommunication](#)

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.supercube2015.typ_a.SuperCubeAOpcUaConfiguration (host:
    str, end-
    point_name:
    str =
    'OPC.SimaticNET.S7',
    port: int
    = 4845,
    sub_handler:
    hvl_ccb.comm.opc.OpcUaSubHandler
    =
    <hvl_ccb.dev.supercube2015.base.
    object at
    0x7f155c314320>,
    up-
    date_period:
    int = 500,
    wait_timeout_retry_sec:
    Union[int,
    float] = 1,
    max_timeout_retry_nr:
    int = 5)
```

Bases: *hvl\_ccb.dev.supercube2015.base.SuperCubeOpcUaCommunicationConfig*

**endpoint\_name** = 'OPC.SimaticNET.S7'

**force\_value** (*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

## Module contents

Supercube package with implementation for the old system version from 2015 based on Siemens WinAC soft-PLC on an industrial 32bit Windows computer.

## Submodules

### hvl\_ccb.dev.base module

Module with base classes for devices.

**class** `hvl_ccb.dev.base.Device` (`dev_config=None`)  
Bases: `hvl_ccb.configuration.ConfigurationMixin`, `abc.ABC`

Base class for devices. Implement this class for a concrete device, such as measurement equipment or voltage sources.

Specifies the methods to implement for a device.

**static config\_cls()**  
Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**start()** → None  
Start or restart this Device. To be implemented in the subclass.

**stop()** → None  
Stop this Device. To be implemented in the subclass.

**exception** `hvl_ccb.dev.base.DeviceExistingException`  
Bases: `Exception`

Exception to indicate that a device with that name already exists.

**class** `hvl_ccb.dev.base.DeviceSequenceMixin` (`devices: Dict[str, hvl_ccb.dev.base.Device]`)  
Bases: `abc.ABC`

Mixin that can be used on a device or other classes to provide facilities for handling multiple devices in a sequence.

**add\_device** (`name: str, device: hvl_ccb.dev.base.Device`) → None  
Add a new device to the device sequence.

#### Parameters

- **name** – is the name of the device.
- **device** – is the instantiated Device object.

**Raises** `DeviceExistingException` –

**get\_device** (`name: str`) → `hvl_ccb.dev.base.Device`  
Get a device by name.

**Parameters** `name` – is the name of the device.

**Returns** the device object from this sequence.

**get\_devices** () → `List[Tuple[str, hvl_ccb.dev.base.Device]]`  
Get list of name, device pairs according to current sequence.

**Returns** A list of tuples with name and device each.

**remove\_device** (`name: str`) → `hvl_ccb.dev.base.Device`  
Remove a device from this sequence and return the device object.

**Parameters** `name` – is the name of the device.

**Returns** device object or `None` if such device was not in the sequence.

**Raises ValueError** – when device with given name was not found

**start()** → None

Start all devices in this sequence in their added order.

**stop()** → None

Stop all devices in this sequence in their reverse order.

**class hvl\_ccb.dev.base.EmptyConfig**

Bases: object

Empty configuration dataclass that is the default configuration for a Device.

**clean\_values()**

Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**class hvl\_ccb.dev.base.SingleCommDevice**(*com, dev\_config=None*)

Bases: *hvl\_ccb.dev.base.Device*, abc.ABC

Base class for devices with a single communication protocol.

**com**

Get the communication protocol of this device.

**Returns** an instance of CommunicationProtocol subtype

**static default\_com\_cls()** → Type[hvl\_ccb.comm.base.CommunicationProtocol]

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**start()** → None

Open the associated communication protocol.

**stop()** → None  
Close the associated communication protocol.

### **hvl\_ccb.dev.crylas module**

Device classes for a CryLas pulsed laser controller and a CryLas laser attenuator, using serial communication.

There are three modes of operation for the laser 1. Laser-internal hardware trigger (default): fixed to 20 Hz and max energy per pulse. 2. Laser-internal software trigger (for diagnosis only). 3. External trigger: required for arbitrary pulse energy or repetition rate. Switch to “external” on the front panel of laser controller for using option 3.

After switching on the laser with `laser_on()`, the system must stabilize for some minutes. Do not apply abrupt changes of pulse energy or repetition rate.

Manufacturer homepage: [https://www.crylas.de/products/pulsed\\_laser.html](https://www.crylas.de/products/pulsed_laser.html)

**class hvl\_ccb.dev.crylas.CryLasAttenuator**(*com, dev\_config=None*)  
Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Device class for the CryLas laser attenuator.

**attenuation**

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**set\_attenuation**(*percent: Union[int, float]*) → None

Set the percentage of attenuated light (inverse of set\_transmission). :param percent: percentage of attenuation, number between 0 and 100 :raises ValueError: if param percent not between 0 and 100 :raises SerialCommunicationIOError: when communication port is not opened :raises CryLasAttenuatorError: if the device does not confirm success

**set\_init\_attenuation()**

Sets the attenuation to its configured initial/default value

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**set\_transmission**(*percent: Union[int, float]*) → None

Set the percentage of transmitted light (inverse of set\_attenuation). :param percent: percentage of transmitted light :raises ValueError: if param percent not between 0 and 100 :raises SerialCommunicationIOError: when communication port is not opened :raises CryLasAttenuatorError: if the device does not confirm success

**start()** → None

Open the com, apply the config value ‘init\_attenuation’

**Raises** *SerialCommunicationIOError* – when communication port cannot be opened

**transmission**

**class hvl\_ccb.dev.crylas.CryLasAttenuatorConfig**(*init\_attenuation: Union[int, float] = 0, response\_sleep\_time: Union[int, float] = 1*)

Bases: object

Device configuration dataclass for CryLas attenuator.

```
clean_values()
force_value (fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing after value has been forced on fieldname.
```

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

```
init_attenuation = 0
is_configdataclass = True
classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.
```

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.
```

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.
```

**Returns** a list of strings containing all required keys.

```
response_sleep_time = 1
exception hvl_ccb.dev.crylas.CryLasAttenuatorError
    Bases: Exception

    General error with the CryLas Attenuator.
```

```
class hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunication (configuration)
    Bases: hvl\_ccb.comm.serial.SerialCommunication

    Specific communication protocol implementation for the CryLas attenuator. Already predefines device-specific protocol parameters in config.
```

```
static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class
```

```
class hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig (port: str,  
baudrate:  
int      =  
9600,  
parity:  
Union[str,  
hvl_ccb.comm.serial.SerialCommu  
= <Serial-  
Communi-  
cationPar-  
ity.NONE:  
'N',  
stopbits:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionStop-  
bits.ONE:  
1>, bytesize:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionByte-  
size.EIGHTBITS:  
8>, terminator:  
bytes  
= b", timeout:  
Union[int,  
float] =  
3)
```

Bases: [hvl\\_ccb.comm.serial.SerialCommunicationConfig](#)

**baudrate = 9600**

Baudrate for CryLas attenuator is 9600 baud

**bytesize = 8**

One byte is eight bits long

**force\_value (fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys () → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
parity = 'N'  
CryLas attenuator does not use parity
```

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
stopbits = 1  
CryLas attenuator uses one stop bit
```

```
terminator = b''  
No terminator
```

```
timeout = 3  
use 3 seconds timeout as default
```

```
class hvl_ccb.dev.crylas.CryLasLaser(com, dev_config=None)
```

Bases: `hvl_ccb.dev.base.SingleCommDevice`

CryLas laser controller device class.

```
class AnswersShutter(*args, **kwds)
```

Bases: `aenum.Enum`

Standard answers of the CryLas laser controller to ‘Shutter’ command passed via *com*.

```
CLOSED = 'Shutter inaktiv'  
OPENED = 'Shutter aktiv'
```

```
class AnswersStatus(*args, **kwds)
```

Bases: `aenum.Enum`

Standard answers of the CryLas laser controller to ‘STATUS’ command passed via *com*.

```
ACTIVE = 'STATUS: Laser active'  
HEAD = 'STATUS: Head ok'  
INACTIVE = 'STATUS: Laser inactive'  
READY = 'STATUS: System ready'  
TEC1 = 'STATUS: TEC1 Regulation ok'  
TEC2 = 'STATUS: TEC2 Regulation ok'
```

```
class LaserStatus(*args, **kwds)
```

Bases: `aenum.Enum`

Status of the CryLas laser

```
READY_INACTIVE = 1  
READ_ACTIVE = 2
```

```
UNREADY_INACTIVE = 0
is_inactive
is_ready

class RepetitionRates(*args, **kwds)
Bases: aenum.IntEnum

Repetition rates for the internal software trigger in Hz

HARDWARE = 0
SOFTWARE_INTERNAL_SIXTY = 60
SOFTWARE_INTERNAL_TEN = 10
SOFTWARE_INTERNAL_TWENTY = 20

ShutterStatus
alias of CryLasLaserShutterStatus

close_shutter() → None
Close the laser shutter.

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserError – if success is not confirmed by the device

static config_cls()
Return the default configdataclass class.

Returns a reference to the default configdataclass class

static default_com_cls()
Get the class for the default communication protocol used with this device.

Returns the type of the standard communication protocol for this device

get_pulse_energy_and_rate() → Tuple[int, int]
Use the debug mode, return the measured pulse energy and rate.

Returns (energy in micro joule, rate in Hz)

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserError – if the device does not answer the query

laser_off() → None
Turn the laser off.

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserError – if success is not confirmed by the device

laser_on() → None
Turn the laser on.

Raises
• SerialCommunicationIOError – when communication port is not opened
• CryLasLaserNotReadyError – if the laser is not ready to be turned on
```

- *CryLasLaserError* – if success is not confirmed by the device

**open\_shutter ()** → None  
Open the laser shutter.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**set\_init\_shutter\_status ()** → None  
Open or close the shutter, to match the configured shutter\_status.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**set\_pulse\_energy (energy: int)** → None  
Sets the energy of pulses (works only with external hardware trigger). Proceed with small energy steps, or the regulation may fail.

**Parameters** **energy** – energy in micro joule

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if the device does not confirm success

**set\_repetition\_rate (rate: Union[int, hvl\_ccb.dev.crylas.CryLasLaser.RepetitionRates])** → None  
Sets the repetition rate of the internal software trigger.

**Parameters** **rate** – frequency (Hz) as an integer

**Raises**

- *ValueError* – if rate is not an accepted value in RepetitionRates Enum
- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**start ()** → None  
Opens the communication protocol and configures the device.

**Raises** *SerialCommunicationIOError* – when communication port cannot be opened

**stop ()** → None  
Stops the device and closes the communication protocol.

**Raises**

- *SerialCommunicationIOError* – if com port is closed unexpectedly
- *CryLasLaserError* – if laser\_off() or close\_shutter() fail

**target\_pulse\_energy**

**update\_laser\_status ()** → None

Update the laser status to *LaserStatus.NOT\_READY* or *LaserStatus.INACTIVE* or *LaserStatus.ACTIVE*.

Note: laser never explicitly says that it is not ready (*LaserStatus.NOT\_READY*) in response to ‘STATUS’ command. It only says that it is ready (heated-up and implicitly inactive/off) or active (on). If it’s not

either of these then the answer is *Answers.HEAD*. Moreover, the only time the laser explicitly says that its status is inactive (*Answers.INACTIVE*) is after issuing a ‘LASER OFF’ command.

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**update\_repetition\_rate()** → None

Query the laser repetition rate.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**update\_shutter\_status()** → None

Update the shutter status (OPENED or CLOSED)

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**update\_target\_pulse\_energy()** → None

Query the laser pulse energy.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *CryLasLaserError* – if success is not confirmed by the device

**wait\_until\_ready()** → None

Block execution until the laser is ready

**Raises** *CryLasLaserError* – if the polling thread stops before the laser is ready

```
class hvl_ccb.dev.crylas.CryLasLaserConfig(calibration_factor: Union[int, float]
                                              = 4.35, polling_period: Union[int, float] = 5, polling_timeout: Union[int, float] = 300, on_start_wait_until_ready: bool = False, auto_laser_on: bool = True, init_shutter_status: Union[int, hvl_ccb.dev.crylas.CryLasLaserShutterStatus] = <CryLasLaserShutterStatus.CLOSED: 0>)
```

Bases: object

Device configuration dataclass for the CryLas laser controller.

**ShutterStatus**

alias of *CryLasLaserShutterStatus*

**auto\_laser\_on = True**

**calibration\_factor = 4.35**

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field

- **value** – value to assign

```
init_shutter_status = 0
```

```
is_configdataclass = True
```

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
on_start_wait_until_ready = False
```

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
polling_period = 5
```

```
polling_timeout = 300
```

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
exception hvl_ccb.dev.crylas.CryLasLaserError
```

Bases: Exception

General error with the CryLas Laser.

```
exception hvl_ccb.dev.crylas.CryLasLaserNotReadyError
```

Bases: [hvl\\_ccb.dev.crylas.CryLasLaserError](#)

Error when trying to turn on the CryLas Laser before it is ready.

```
class hvl_ccb.dev.crylas.CryLasLaserSerialCommunication(configuration)
```

Bases: [hvl\\_ccb.comm.serial.SerialCommunication](#)

Specific communication protocol implementation for the CryLas laser controller. Already predefines device-specific protocol parameters in config.

```
READ_TEXT_SKIP_PREFIXES = ('>', 'MODE:')
```

Prefixes of lines that are skipped when read from the serial port.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
query(cmd: str, prefix: str, post_cmd: str = None) → str
```

Send a command, then read the com until a line starting with prefix, or an empty line, is found. Returns the line in question.

**Parameters**

- **cmd** – query message to send to the device
- **prefix** – start of the line to look for in the device answer
- **post\_cmd** – optional additional command to send after the query

**Returns** line in question as a string

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**query\_all** (`cmd: str, prefix: str`)

Send a command, then read the com until a line starting with prefix, or an empty line, is found. Returns a list of successive lines starting with prefix.

#### Parameters

- `cmd` – query message to send to the device
- `prefix` – start of the line to look for in the device answer

**Returns** line in question as a string

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**read\_text** () → str

Read first line of text from the serial port that does not start with any of `self.READ_TEXT_SKIP_PREFIXES`.

**Returns** String read from the serial port; '' if there was nothing to read.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

```
class hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig(port: str, baudrate: int = 19200, parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity.NONE: 'N', stopbits: Union[int, hvl_ccb.comm.serial.SerialCommunicationStopbits.ONE: 1], bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationByteSize.EIGHTBITS: 8], terminator: bytes = b'n', timeout: Union[int, float] = 3)
```

Bases: `hvl_ccb.comm.serial.SerialCommunicationConfig`

**baudrate = 19200**

Baudrate for CryLas laser is 19200 baud

**bytesize = 8**

One byte is eight bits long

**force\_value** (`fieldname, value`)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**

CryLas laser does not use parity

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**

CryLas laser uses one stop bit

**terminator = b'\n'**

The terminator is LF

**timeout = 3**

use 3 seconds timeout as default

**class hvl\_ccb.dev.crylas.CryLasLaserShutterStatus(\*args, \*\*kwargs)**

Bases: aenum.Enum

Status of the CryLas laser shutter

**CLOSED = 0**

**OPENED = 1**

## hvl\_ccb.dev.ea\_psi9000 module

Device class for controlling a Elektro Automatik PSI 9000 power supply over VISA.

It is necessary that a backend for pyvisa is installed. This can be NI-Visa oder pyvisa-py (up to know, all the testing was done with NI-Visa)

**class hvl\_ccb.dev.ea\_psi9000.PSI9000(com: Union[hvl\_ccb.dev.ea\_psi9000.PSI9000VisaCommunication, hvl\_ccb.dev.ea\_psi9000.PSI9000VisaCommunicationConfig, dict], dev\_config: Union[hvl\_ccb.dev.ea\_psi9000.PSI9000Config, dict, None] = None)**

Bases: *hvl\_ccb.dev.visa.VisaDevice*

Elektro Automatik PSI 9000 power supply.

**MS\_NOMINAL\_CURRENT = 2040**

```
MS_NOMINAL_VOLTAGE = 80
SHUTDOWN_CURRENT_LIMIT = 0.1
SHUTDOWN_VOLTAGE_LIMIT = 0.1
check_master_slave_config() → None
    Checks if the master / slave configuration and initializes if successful
    Raises PSI9000Error – if master-slave configuration failed

static config_cls()
    Return the default configdataclass class.
    Returns a reference to the default configdataclass class

static default_com_cls()
    Return the default communication protocol for this device type, which is VisaCommunication.
    Returns the VisaCommunication class

get_output() → bool
    Reads the current state of the DC output of the source. Returns True, if it is enabled, false otherwise.
    Returns the state of the DC output

get_system_lock() → bool
    Get the current lock state of the system. The lock state is true, if the remote control is active and false, if not.
    Returns the current lock state of the device

get_ui_lower_limits() → Tuple[float, float]
    Get the lower voltage and current limits. A lower power limit does not exist.
    Returns Umin in V, Imin in A

get_uip_upper_limits() → Tuple[float, float, float]
    Get the upper voltage, current and power limits.
    Returns Umax in V, Imax in A, Pmax in W

get_voltage_current_setpoint() → Tuple[float, float]
    Get the voltage and current setpoint of the current source.
    Returns Uset in V, Iset in A

measure_voltage_current() → Tuple[float, float]
    Measure the DC output voltage and current
    Returns Umeas in V, Imeas in A

set_lower_limits(voltage_limit: float = None, current_limit: float = None) → None
    Set the lower limits for voltage and current. After writing the values a check is performed if the values are set correctly.

    Parameters
        • voltage_limit – is the lower voltage limit in V
        • current_limit – is the lower current limit in A
    Raises PSI9000Error – if the limits are out of range

set_output(target_onstate: bool) → None
    Enables / disables the DC output.
```

**Parameters** `target_onstate` – enable or disable the output power

**Raises** `PSI9000Error` – if operation was not successful

**set\_system\_lock** (`lock: bool`) → None

Lock / unlock the device, after locking the control is limited to this class unlocking only possible when voltage and current are below the defined limits

**Parameters** `lock` – True: locking, False: unlocking

**set\_upper\_limits** (`voltage_limit: float = None, current_limit: float = None, power_limit: float = None`) → None

Set the upper limits for voltage, current and power. After writing the values a check is performed if the values are set. If a parameter is left blank, the maximum configurable limit is set.

**Parameters**

- `voltage_limit` – is the voltage limit in V
- `current_limit` – is the current limit in A
- `power_limit` – is the power limit in W

**Raises** `PSI9000Error` – if limits are out of range

**set\_voltage\_current** (`volt: float, current: float`) → None

Set voltage and current setpoints.

After setting voltage and current, a check is performed if writing was successful.

**Parameters**

- `volt` – is the setpoint voltage: 0..81.6 V (1.02 \* 0-80 V) (absolute max, can be smaller if limits are set)
- `current` – is the setpoint current: 0..2080.8 A (1.02 \* 0 - 2040 A) (absolute max, can be smaller if limits are set)

**Raises** `PSI9000Error` – if the desired setpoint is out of limits

**start()** → None

Start this device.

**stop()** → None

Stop this device. Turns off output and lock, if enabled.

```
class hvl_ccb.dev.ea_psi9000.PSI9000Config(spoll_interval: Union[int, float] = 0.5,
                                             spoll_start_delay: Union[int, float] = 2,
                                             power_limit: Union[int, float] = 43500, voltage_lower_limit: Union[int, float] = 0.0, voltage_upper_limit: Union[int, float] = 10.0, current_lower_limit: Union[int, float] = 0.0, current_upper_limit: Union[int, float] = 2040.0, wait_sec_system_lock: Union[int, float] = 0.5, wait_sec_settings_effect: Union[int, float] = 1, wait_sec_initialisation: Union[int, float] = 2)
```

Bases: `hvl_ccb.dev.visa.VisaDeviceConfig`

Elektro Automatik PSI 9000 power supply device class. The device is communicating over a VISA TCP socket.

Using this power supply, DC voltage and current can be supplied to a load with up to 2040 A and 80 V (using all four available units in parallel). The maximum power is limited by the grid, being at 43.5 kW available through the CEE63 power socket.

**clean\_values ()** → None  
Cleans and enforces configuration values. Does nothing by default, but may be overridden to add custom configuration value checks.

**current\_lower\_limit = 0.0**  
Lower current limit in A, depending on the experimental setup.

**current\_upper\_limit = 2040.0**  
Upper current limit in A, depending on the experimental setup.

**force\_value (fieldname, value)**  
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys ()** → Sequence[str]  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults ()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**power\_limit = 43500**

Power limit in W depending on the experimental setup. With 3x63A, this is 43.5kW. Do not change this value, if you do not know what you are doing. There is no lower power limit.

**classmethod required\_keys ()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**voltage\_lower\_limit = 0.0**

Lower voltage limit in V, depending on the experimental setup.

**voltage\_upper\_limit = 10.0**

Upper voltage limit in V, depending on the experimental setup.

**wait\_sec\_initialisation = 2**

**wait\_sec\_settings\_effect = 1**

**wait\_sec\_system\_lock = 0.5**

**exception hvl\_ccb.dev.ea\_psi9000.PSI9000Error**

Bases: Exception

Base error class regarding problems with the PSI 9000 supply.

**class hvl\_ccb.dev.ea\_psi9000.PSI9000VisaCommunication (configuration)**

Bases: *hvl\_ccb.comm.visa.VisaCommunication*

Communication protocol used with the PSI 9000 power supply.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig(host: str, interface_type: Union[str, hvl_ccb.comm.visa.VisaCommunicationConfig] = <Interface-Type.TCPIP_SOCKET: I>, board: int = 0, port: int = 5025, timeout: int = 5000, chunk_size: int = 204800, open_timeout: int = 1000, write_termination: str = '\n', read_termination: str = '\n', visa_backend: str = '')
```

Bases: *hvl\_ccb.comm.visa.VisaCommunicationConfig*

Visa communication protocol config dataclass with specification for the PSI 9000 power supply.

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
interface_type = 1
```

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

```
classmethod required_keys() → Sequence[str]
```

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

### **hvl\_ccb.dev.heinzinger module**

Device classes for Heinzinger Digital Interface I/II and Heinzinger PNC power supply.

The Heinzinger Digital Interface I/II is used for many Heinzinger power units. Manufacturer homepage: <https://www.heinzinger.com/products/accessories-and-more/digital-interfaces/>

The Heinzinger PNC series is a series of high voltage direct current power supplies. The class HeinzingerPNC is tested with two PNChp 60000-1neg and a PNChp 1500-1neg. Check the code carefully before using it with other PNC devices, especially PNC3p or PNCcap. Manufacturer homepage: <https://www.heinzinger.com/products/high-voltage/universal-high-voltage-power-supplies/>

```
class hvl_ccb.dev.heinzinger.HeinzingerConfig(default_number_of_recordings:  
                                              Union[int,  
                                              hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum]  
                                              = 1, number_of_decimals: int = 6,  
                                              wait_sec_stop_commands: Union[int,  
                                              float] = 0.5)
```

Bases: object

Device configuration dataclass for Heinzinger power supplies.

```
class RecordingsEnum
```

Bases: enum.IntEnum

An enumeration.

```
EIGHT = 8
```

```
FOUR = 4
```

```
ONE = 1
```

```
SIXTEEN = 16
```

```
TWO = 2
```

```
clean_values()
```

```
default_number_of_recordings = 1
```

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
is_configdataclass = True
```

```
classmethod keys() → Sequence[str]
```

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

```
number_of_decimals = 6
```

```
classmethod optional_defaults() → Dict[str, object]
```

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys () → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**wait\_sec\_stop\_commands = 0.5**

Time to wait after subsequent commands during stop (in seconds)

**class hvl\_ccb.dev.heinzinger.HeinzingerDI (com, dev\_config=None)**

Bases: *hvl\_ccb.dev.base.SingleCommDevice*, abc.ABC

Heinzinger Digital Interface I/II device class

Sends basic SCPI commands and reads the answer. Only the standard instruction set from the manual is implemented.

**static config\_cls ()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls ()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**get\_current () → float**

Queries the set current of the Heinzinger PNC (not the measured current!).

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_interface\_version () → str**

Queries the version number of the digital interface.

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_number\_of\_recordings () → int**

Queries the number of recordings the device is using for average value calculation.

**Returns** int number of recordings

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_serial\_number () → str**

Ask the device for its serial number and returns the answer as a string.

**Returns** string containing the device serial number

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**get\_voltage () → float**

Queries the set voltage of the Heinzinger PNC (not the measured voltage!).

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**measure\_current () → float**

Ask the Device to measure its output current and return the measurement result.

**Returns** measured current as float

**Raises** *SerialCommunicationIOError* – when communication port is not opened

**measure\_voltage () → float**

Ask the Device to measure its output voltage and return the measurement result.

**Returns** measured voltage as float

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**output\_off()** → None  
Switch DC voltage output off.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**output\_on()** → None  
Switch DC voltage output on.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**reset\_interface()** → None  
Reset of the digital interface; only Digital Interface I: Power supply is switched to the Local-Mode (Manual operation)

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**set\_current** (`value: Union[int, float]`) → None  
Sets the output current of the Heinzinger PNC to the given value.

**Parameters** `value` – current expressed in `self.unit_current`

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**set\_number\_of\_recordings** (`value: Union[int, hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum]`) → None  
Sets the number of recordings the device is using for average value calculation. The possible values are 1, 2, 4, 8 and 16.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**set\_voltage** (`value: Union[int, float]`) → None  
Sets the output voltage of the Heinzinger PNC to the given value.

**Parameters** `value` – voltage expressed in `self.unit_voltage`

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**start()**  
Opens the communication protocol.

**Raises** `SerialCommunicationIOError` – when communication port cannot be opened.

**stop()** → None  
Stop the device. Closes also the communication protocol.

**class** `hvl_ccb.dev.heinzinger.HeinzingerPNC` (`com, dev_config=None`)  
Bases: `hvl_ccb.dev.heinzinger.HeinzingerDI`

Heinzinger PNC power supply device class.

The power supply is controlled over a Heinzinger Digital Interface I/II

**class** `UnitCurrent` (\*args, \*\*kwds)  
Bases: `hvl_ccb.utils.enum.AutoNumberNameEnum`

`A = 3`  
`UNKNOWN = 1`  
`mA = 2`

**class** `UnitVoltage` (\*args, \*\*kwds)  
Bases: `hvl_ccb.utils.enum.AutoNumberNameEnum`

```

UNKNOWN = 1
V = 2
kV = 3

identify_device() → None
Identify the device nominal voltage and current based on its serial number.

Raises SerialCommunicationIOError – when communication port is not opened

max_current
max_current_hardware
max_voltage
max_voltage_hardware

set_current (value: Union[int, float]) → None
Sets the output current of the Heinzinger PNC to the given value.

Parameters value – current expressed in self.unit_current
Raises SerialCommunicationIOError – when communication port is not opened

set_voltage (value: Union[int, float]) → None
Sets the output voltage of the Heinzinger PNC to the given value.

Parameters value – voltage expressed in self.unit_voltage
Raises SerialCommunicationIOError – when communication port is not opened

start() → None
Opens the communication protocol and configures the device.

unit_current
unit_voltage

exception hvl_ccb.dev.heinzinger.HeinzingerPNCDeviceNotRecognizedException
Bases: hvl_ccb.dev.heinzinger.HeinzingerPNCError
Error indicating that the serial number of the device is not recognized.

exception hvl_ccb.dev.heinzinger.HeinzingerPNCError
Bases: Exception
General error with the Heinzinger PNC voltage source.

exception hvl_ccb.dev.heinzinger.HeinzingerPNCMaxCurrentExceededException
Bases: hvl_ccb.dev.heinzinger.HeinzingerPNCError
Error indicating that program attempted to set the current to a value exceeding ‘max_current’.

exception hvl_ccb.dev.heinzinger.HeinzingerPNCMaxVoltageExceededException
Bases: hvl_ccb.dev.heinzinger.HeinzingerPNCError
Error indicating that program attempted to set the voltage to a value exceeding ‘max_voltage’.

class hvl_ccb.dev.heinzinger.HeinzingerSerialCommunication (configuration)
Bases: hvl_ccb.comm.serial.SerialCommunication
Specific communication protocol implementation for Heinzinger power supplies. Already predefines device-specific protocol parameters in config.

static config_cls()
Return the default configdataclass class.

```

**Returns** a reference to the default configdataclass class

**read\_text\_nonempty** (*n\_attempts\_max*: int = 40) → str

Reads from the serial port, until a non-empty line is found, or the number of attempts is exceeded.

**Parameters**

- **n\_attempts\_max** – maximum number of read attempts

- **attempt\_interval\_sec** – interval between subsequent attempts in seconds

**Returns** String read from the serial port; ‘’ if number of attempts is exceeded or serial port is not opened.

```
class hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig(port:      str,
                                                               baudrate:
                                                               int = 9600,
                                                               parity:
                                                               Union[str,
                                                               hvl_ccb.comm.serial.SerialCommunicationParity.NONE:
                                                               'N'],
                                                               stopbits:
                                                               Union[int,
                                                               hvl_ccb.comm.serial.SerialCommunicationStopbits.ONE:
                                                               1],
                                                               bytesize:
                                                               Union[int,
                                                               hvl_ccb.comm.serial.SerialCommunicationBytesize.EIGHTBITS:
                                                               8],
                                                               terminator:
                                                               bytes =
                                                               b'n',
                                                               timeout:
                                                               Union[int,
                                                               float] = 3,
                                                               wait_sec_read_text_nonempty:
                                                               Union[int,
                                                               float] = 0.5)
```

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

**baudrate = 9600**

Baudrate for Heinzinger power supplies is 9600 baud

**bytesize = 8**

One byte is eight bits long

**clean\_values()**

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**

Heinzinger does not use parity

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**

Heinzinger uses one stop bit

**terminator = b'\n'**

The terminator is LF

**timeout = 3**

use 3 seconds timeout as default

**wait\_sec\_read\_text\_nonempty = 0.5**

time to wait between attempts of reading a non-empty text

## hvl\_ccb.dev.labjack module

Labjack Device for hvl\_ccb. Originally developed and tested for LabJack T7-PRO.

Makes use of the LabJack LJM Library Python wrapper. This wrapper needs an installation of the LJM Library for Windows, Mac OS X or Linux. Go to: <https://labjack.com/support/software/installers/ljm> and <https://labjack.com/support/software/examples/ljm/python>

**class hvl\_ccb.dev.labjack.LabJack(com, dev\_config=None)**  
Bases: *hvl\_ccb.dev.base.SingleCommDevice*

LabJack Device.

This class is tested with a LabJack T7-Pro and should also work with T4 and T7 devices communicating through the LJM Library. Other or older hardware versions and variants of LabJack devices are not supported.

**class AInRange(\*args, \*\*kwds)**  
Bases: *hvl\_ccb.utils.enum.StrEnumBase*  
**ONE = 1.0**  
**ONE\_HUNDREDTH = 0.01**

```
ONE_TENTH = 0.1
TEN = 10.0
value

class CalMicroAmpere(*args, **kwds)
Bases: aenum.Enum

Pre-defined microampere (uA) values for calibration current source query.

TEN = '10uA'
TWO_HUNDRED = '200uA'

class CjcType(*args, **kwds)
Bases: hvl_ccb.utils.enum.NameEnum

CJC slope and offset

internal = (1, 0)
lm34 = (55.56, 255.37)

DIOChannel
alias of hvl_ccb._dev.labjack.TSeriesDIOChannel

class DIOStatus(*args, **kwds)
Bases: aenum.IntEnum

State of a digital I/O channel.

HIGH = 1
LOW = 0

class DeviceType(*args, **kwds)
Bases: hvl_ccb.utils.enum.AutoNumberNameEnum

LabJack device types.

Can be also looked up by ambiguous Product ID (p_id) or by instance name: ``python
LabJackDeviceType(4) is LabJackDeviceType('T4') ``

ANY = 1
T4 = 2
T7 = 3
T7_PRO = 4
get_by_p_id = <bound method DeviceType.get_by_p_id of <aenum 'DeviceType'>>

class TemperatureUnit(*args, **kwds)
Bases: hvl_ccb.utils.enum.NameEnum

Temperature unit (to be returned)

C = 1
F = 2
K = 0

class ThermocoupleType(*args, **kwds)
Bases: hvl_ccb.utils.enum.NameEnum

Thermocouple type; NONE means disable thermocouple mode.
```

```
C = 30
E = 20
J = 21
K = 22
NONE = 0
PT100 = 40
PT1000 = 42
PT500 = 41
R = 23
S = 25
T = 24
```

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**get\_ain(\*channels) → Union[float, Sequence[float]]**

Read currently measured value (voltage, resistance, ...) from one or more of analog inputs.

**Parameters** **channels** – AIN number or numbers (0..254)

**Returns** the read value (voltage, resistance, ...) as *float* or ‘tuple of them in case multiple channels given

**get\_cal\_current\_source(name: Union[str, CalMicroAmpere]) → float**

This function will return the calibration of the chosen current source, this ist not a measurement!

The value was stored during fabrication.

**Parameters** **name** – ‘200uA’ or ‘10uA’ current source

**Returns** calibration of the chosen current source in ampere

**get\_digital\_input(address: Union[str, hvl\_ccb.\_dev.labjack.TSeriesDIOChannel]) → hvl\_ccb.dev.labjack.LabJack.DIOStatus**

Get the value of a digital input.

allowed names for T7 (Pro): FIO0 - FIO7, EIO0 - EIO 7, CIO0- CIO3, MIO0 - MIO2 :param address: name of the output -> ‘FIO0’ :return: HIGH when *address* DIO is high, and LOW when *address* DIO is low

**get\_product\_id() → int**

This function returns the product ID reported by the connected device.

Attention: returns 7 for both T7 and T7-Pro devices!

**Returns** integer product ID of the device

**get\_product\_name(force\_query\_id=False) → str**

This function will return the product name based on product ID reported by the device.

Attention: returns “T7” for both T7 and T7-Pro devices!

**Parameters** **force\_query\_id** – boolean flag to force *get\_product\_id* query to device instead of using cached device type from previous queries.

**Returns** device name string, compatible with *LabJack.DeviceType*

**get\_product\_type** (*force\_query\_id*: bool = False) → hvl\_ccb.\_dev.labjack.DeviceType

This function will return the device type based on reported device type and in case of unambiguity based on configuration of device's communication protocol (e.g. for "T7" and "T7\_PRO" devices), or, if not available first matching.

**Parameters** **force\_query\_id** – boolean flag to force *get\_product\_id* query to device instead of using cached device type from previous queries.

**Returns** *DeviceType* instance

**Raises** *LabJackIdentifierDIOError* – when read Product ID is unknown

**get\_sbush\_rh** (*number*: int) → float

Read the relative humidity value from a serial SBUS sensor.

**Parameters** **number** – port number (0..22)

**Returns** relative humidity in %RH

**get\_sbush\_temp** (*number*: int) → float

Read the temperature value from a serial SBUS sensor.

**Parameters** **number** – port number (0..22)

**Returns** temperature in Kelvin

**get\_serial\_number** () → int

Returns the serial number of the connected LabJack.

**Returns** Serial number.

**read\_resistance** (*channel*: int) → float

Read resistance from specified channel.

**Parameters** **channel** – channel with resistor

**Returns** resistance value with 2 decimal places

**read\_thermocouple** (*pos\_channel*: int) → float

Read the temperature of a connected thermocouple.

**Parameters** **pos\_channel** – is the AIN number of the positive pin

**Returns** temperature in specified unit

**set\_ain\_differential** (*pos\_channel*: int, *differential*: bool) → None

Sets an analog input to differential mode or not. T7-specific: For base differential channels, positive must be even channel from 0-12 and negative must be positive+1. For extended channels 16-127, see Mux80 datasheet.

**Parameters**

- **pos\_channel** – is the AIN number (0..12)
- **differential** – True or False

**Raises** *LabJackError* – if parameters are unsupported

**set\_ain\_range** (*channel*: int, *vrange*: Union[Real, AInRange]) → None

Set the range of an analog input port.

**Parameters**

- **channel** – is the AIN number (0..254)
- **vrange** – is the voltage range to be set

**set\_ain\_resistance** (*channel: int, vrange: Union[Real, AInRange], resolution: int*) → None  
Set the specified channel to resistance mode. It utilized the 200uA current source of the LabJack.

#### Parameters

- **channel** – channel that should measure the resistance
- **vrange** – voltage range of the channel
- **resolution** – resolution index of the channel T4: 0-5, T7: 0-8, T7-Pro 0-12

**set\_ain\_resolution** (*channel: int, resolution: int*) → None  
Set the resolution index of an analog input port.

#### Parameters

- **channel** – is the AIN number (0..254)
- **resolution** – is the resolution index within 0...‘get\_product\_type().ain\_max\_resolution’ range; 0 will set the resolution index to default value.

**set\_ain\_thermocouple** (*pos\_channel: int, thermocouple: Union[None, str, ThermocoupleType], cjc\_address: int = 60050, cjc\_type: Union[str, CjcType] = <CjcType.internal: (1, 0)>, vrange: Union[Real, AInRange] = <AInRange.ONE\_HUNDREDTH: '0.01'>, resolution: int = 10, unit: Union[str, TemperatureUnit] = <TemperatureUnit.K: 0>*) → None  
Set the analog input channel to thermocouple mode.

#### Parameters

- **pos\_channel** – is the analog input channel of the positive part of the differential pair
- **thermocouple** – None to disable thermocouple mode, or string specifying the thermocouple type
- **cjc\_address** – modbus register address to read the CJC temperature
- **cjc\_type** – determines cjc slope and offset, ‘internal’ or ‘lm34’
- **vrange** – measurement voltage range
- **resolution** – resolution index (T7-Pro: 0-12)
- **unit** – is the temperature unit to be returned (‘K’, ‘C’ or ‘F’)

**Raises** `LabJackError` – if parameters are unsupported

**set\_digital\_output** (*address: str, state: Union[int, DIOStatus]*) → None  
Set the value of a digital output.

#### Parameters

- **address** – name of the output -> ‘FIO0’
- **state** – state of the output -> *DIOStatus* instance or corresponding *int* value

**start()** → None  
Start the Device.

**stop()** → None  
Stop the Device.

**exception** `hvl_ccb.dev.labjack.LabJackError`  
Bases: Exception

Errors of the LabJack device.

```
exception hvl_ccb.dev.labjack.LabJackIdentifierDIOError
Bases: Exception
Error indicating a wrong DIO identifier
```

## hvl\_ccb.dev.mbw973 module

Device class for controlling a MBW 973 SF6 Analyzer over a serial connection.

The MBW 973 is a gas analyzer designed for gas insulated switchgear and measures humidity, SF6 purity and SO2 contamination in one go. Manufacturer homepage: <https://www.mbw.ch/products/sf6-gas-analysis/973-sf6-analyzer/>

```
class hvl_ccb.dev.mbw973.MBW973 (com, dev_config=None)
Bases: hvl_ccb.dev.base.SingleCommDevice

MBW 973 dew point mirror device class.

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

is_done() → bool
    Poll status of the dew point mirror and return True, if all measurements are done.

    Returns True, if all measurements are done; False otherwise.

    Raises SerialCommunicationIOError – when communication port is not opened

read(cast_type: Type[CT_co] = <class 'str'>)
    Read value from self.com and cast to cast_type. Raises ValueError if read text (str) is not convertible to
    cast_type, e.g. to float or to int.

    Returns Read value of cast_type type.

read_float() → float
    Convenience wrapper for self.read(), with typing hint for return value.

    Returns Read float value.

read_int() → int
    Convenience wrapper for self.read(), with typing hint for return value.

    Returns Read int value.

read_measurements() → Dict[str, float]
    Read out measurement values and return them as a dictionary.

    Returns Dictionary with values.

    Raises SerialCommunicationIOError – when communication port is not opened

set_measuring_options(humidity: bool = True, sf6_purity: bool = False) → None
    Send measuring options to the dew point mirror.

Parameters
    • humidity – Perform humidity test or not?
```

- **sf6\_purity** – Perform SF6 purity test or not?

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**start()** → None

Start this device. Opens the communication protocol and retrieves the set measurement options from the device.

**Raises** `SerialCommunicationIOError` – when communication port cannot be opened.

**start\_control()** → None

Start dew point control to acquire a new value set.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**stop()** → None

Stop the device. Closes also the communication protocol.

**write(value)** → None

Send `value` to `self.com`.

**Parameters** `value` – Value to send, converted to `str`.

**Raises** `SerialCommunicationIOError` – when communication port is not opened

**class** `hvl_ccb.dev.mbw973.MBW973Config(polling_interval: Union[int, float] = 2)`

Bases: `object`

Device configuration dataclass for MBW973.

**clean\_values()**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

#### Parameters

- **fieldname** – name of the field

- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**polling\_interval = 2**

Polling period for `is_done` status queries [in seconds].

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

```
exception hvl_ccb.dev.mbw973.MBW973ControlRunningException
```

Bases: `hvl_ccb.dev.mbw973.MBW973Error`

Error indicating there is still a measurement running, and a new one cannot be started.

```
exception hvl_ccb.dev.mbw973.MBW973Error
```

Bases: `Exception`

General error with the MBW973 dew point mirror device.

```
exception hvl_ccb.dev.mbw973.MBW973PumpRunningException
```

Bases: `hvl_ccb.dev.mbw973.MBW973Error`

Error indicating the pump of the dew point mirror is still recovering gas, unable to start a new measurement.

```
class hvl_ccb.dev.mbw973.MBW973SerialCommunication(configuration)
```

Bases: `hvl_ccb.comm.serial.SerialCommunication`

Specific communication protocol implementation for the MBW973 dew point mirror. Already predefines device-specific protocol parameters in config.

```
static config_cls()
```

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig(port: str, baudrate: int =
```

`9600`, parity: `Union[str,`

`hvl_ccb.comm.serial.SerialCommunicationParity]`

`= <SerialCommuni-`

`cationParity.NONE: 'N',`

`stopbits: Union[int,`

`hvl_ccb.comm.serial.SerialCommunicationStopbit]`

`= <SerialCommuni-`

`cationStopbits.ONE:`

`1>, bytesize: Union[int,`

`hvl_ccb.comm.serial.SerialCommunicationBytesize]`

`= <SerialCom-`

`municationByte-`

`size.EIGHTBITS: 8>,`

`terminator: bytes = b'r',`

`timeout: Union[int, float]`

`= 3)`

Bases: `hvl_ccb.comm.serial.SerialCommunicationConfig`

```
baudrate = 9600
```

Baudrate for MBW973 is 9600 baud

```
bytesize = 8
```

One byte is eight bits long

```
force_value(fieldname, value)
```

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define `post_force_value` method with same signature as this method to do extra processing after `value` has been forced on `fieldname`.

## Parameters

- **fieldname** – name of the field
- **value** – value to assign

```
classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

Returns a list of strings containing all keys.

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

Returns a list of strings containing all optional keys.

parity = 'N'
    MBW973 does not use parity

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

Returns a list of strings containing all required keys.

stopbits = 1
    MBW973 does use one stop bit

terminator = b'\r'
    The terminator is only CR

timeout = 3
    use 3 seconds timeout as default

class hvl_ccb.dev.mbw973.Poller(period: float, callback: Callable)
Bases: object

Wrapper for the threading.Timer class to periodically poll data.

start() → None
    Start the polling timer.

stop() → None
    Stop the polling timer.

timer_callback() → None
    Callback method that is called every time the timer elapses. It calls the specified user callback function and restarts the timer.
```

## hvl\_ccb.dev.newport module

Device class for Newport SMC100PP stepper motor controller with serial communication.

The SMC100PP is a single axis motion controller/driver for stepper motors up to 48 VDC at 1.5 A rms. Up to 31 controllers can be networked through the internal RS-485 communication link.

Manufacturer homepage: <https://www.newport.com/f/smci00-single-axis-dc-or-stepper-motion-controller>

```
class hvl_ccb.dev.newport.NewportConfigCommands(*args, **kwds)
Bases: hvl_ccb.utils.enum.NameEnum

Commands predefined by the communication protocol of the SMC100PP

AC = 'acceleration'
BA = 'backlash_compensation'
BH = 'hysteresis_compensation'
```

```
FRM = 'micro_step_per_full_step_factor'
FRS = 'motion_distance_per_full_step'
HT = 'home_search_type'
JR = 'jerk_time'
OH = 'home_search_velocity'
OT = 'home_search_timeout'
QIL = 'peak_output_current_limit'
SA = 'rs485_address'
SL = 'negative_software_limit'
SR = 'positive_software_limit'
VA = 'velocity'
VB = 'base_velocity'
ZX = 'stage_configuration'

exception hvl_ccb.dev.newport.NewportControllerError
    Bases: Exception
    Error with the Newport controller.

exception hvl_ccb.dev.newport.NewportMotorError
    Bases: Exception
    Error with the Newport motor.

class hvl_ccb.dev.newport.NewportSMC100PP (com, dev_config=None)
    Bases: hvl_ccb.dev.base.SingleCommDevice
    Device class of the Newport motor controller SMC100PP

    class MotorErrors (*args, **kwds)
        Bases: aenum.Enum
        Possible motor errors reported by the motor during get_state().
        DC_VOLTAGE_TOO_LOW = 3
        FOLLOWING_ERROR = 6
        HOMING_TIMEOUT = 5
        NED_END_OF_TURN = 11
        OUTPUT_POWER_EXCEEDED = 2
        PEAK_CURRENT_LIMIT = 9
        POS_END_OF_TURN = 10
        RMS_CURRENT_LIMIT = 8
        SHORT_CIRCUIT = 7
        WRONG_ESP_STAGE = 4

    class StateMessages (*args, **kwds)
        Bases: aenum.Enum
        Possible messages returned by the controller on get_state() query.
```

```

CONFIG = '14'
DISABLE_FROM_JOGGING = '3E'
DISABLE_FROM_MOVING = '3D'
DISABLE_FROM_READY = '3C'
HOMING_FROM_RS232 = '1E'
HOMING_FROM_SMC = '1F'
JOGGING_FROM_DISABLE = '47'
JOGGING_FROM_READY = '46'
MOVING = '28'
NO_REF_ESP_STAGE_ERROR = '10'
NO_REF_FROM_CONFIG = '0C'
NO_REF_FROM_DISABLED = '0D'
NO_REF_FROM_HOMING = '0B'
NO_REF_FROM_JOGGING = '11'
NO_REF_FROM_MOVING = '0F'
NO_REF_FROM_READY = '0E'
NO_REF_FROM_RESET = '0A'
READY_FROM_DISABLE = '34'
READY_FROM_HOMING = '32'
READY_FROM_JOGGING = '35'
READY_FROM_MOVING = '33'

```

**States**

alias of *NewportStates*

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**exit\_configuration(add: int = None) → None**

Exit the CONFIGURATION state and go back to the NOT REFERENCED state. All configuration parameters are saved to the device's memory.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_acceleration** (*add: int = None*) → Union[int, float]

Leave the configuration state. The configuration parameters are saved to the device's memory.

**Parameters** **add** – controller address (1 to 31)

**Returns** acceleration (preset units/s^2), value between 1e-6 and 1e12

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_controller\_information** (*add: int = None*) → str

Get information on the controller name and driver version

**Parameters** **add** – controller address (1 to 31)

**Returns** controller information

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_motor\_configuration** (*add: int = None*) → Dict[str, float]

Query the motor configuration and returns it in a dictionary.

**Parameters** **add** – controller address (1 to 31)

**Returns** dictionary containing the motor's configuration

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_move\_duration** (*dist: Union[int, float], add: int = None*) → float

Estimate the time necessary to move the motor of the specified distance.

**Parameters**

- **dist** – distance to travel
- **add** – controller address (1 to 31), defaults to self.address

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_negative\_software\_limit** (*add: int = None*) → Union[int, float]

Get the negative software limit (the maximum position that the motor is allowed to travel to towards the left).

**Parameters** **add** – controller address (1 to 31)

**Returns** negative software limit (preset units), value between -1e12 and 0

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_position** (*add: int = None*) → float

Returns the value of the current position.

**Parameters** *add* – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_positive\_software\_limit** (*add: int = None*) → Union[int, float]

Get the positive software limit (the maximum position that the motor is allowed to travel to towards the right).

**Parameters** *add* – controller address (1 to 31)

**Returns** positive software limit (preset units), value between 0 and 1e12

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**get\_state** (*add: int = None*) → StateMessages

Check on the motor errors and the controller state

**Parameters** *add* – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error
- *NewportMotorError* – if the motor reports an error

**Returns** state message from the device (member of StateMessages)

**go\_home** (*add: int = None*) → None

Move the motor to its home position.

**Parameters** *add* – controller address (1 to 31), defaults to self.address

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

`go_to_configuration(add: int = None) → None`

This method is executed during start(). It can also be executed after a reset(). The controller is put in CONFIG state, where configuration parameters can be changed.

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

`initialize(add: int = None) → None`

Puts the controller from the NOT\_REF state to the READY state. Sends the motor to its “home” position.

**Parameters** `add` – controller address (1 to 31)

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

`move_to_absolute_position(pos: Union[int, float], add: int = None) → None`

Move the motor to the specified position.

**Parameters**

- `pos` – target absolute position (affected by the configured offset)
- `add` – controller address (1 to 31), defaults to self.address

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

`move_to_relative_position(pos: Union[int, float], add: int = None) → None`

Move the motor of the specified distance.

**Parameters**

- `pos` – distance to travel (the sign gives the direction)
- `add` – controller address (1 to 31), defaults to self.address

**Raises**

- `SerialCommunicationIOError` – if the com is closed
- `NewportSerialCommunicationError` – if an unexpected answer is obtained
- `NewportControllerError` – if the controller reports an error

`reset(add: int = None) → None`

Resets the controller, equivalent to a power-up. This puts the controller back to NOT REFERENCED state, which is necessary for configuring the controller.

**Parameters** `add` – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_acceleration** (*acc: Union[int, float], add: int = None*) → None

Leave the configuration state. The configuration parameters are saved to the device's memory.

#### Parameters

- **acc** – acceleration (preset units/s<sup>2</sup>), value between 1e-6 and 1e12
- **add** – controller address (1 to 31)

#### Raises

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_motor\_configuration** (*add: int = None, config: dict = None*) → None

Set the motor configuration. The motor must be in CONFIG state.

#### Parameters

- **add** – controller address (1 to 31)
- **config** – dictionary containing the motor's configuration

#### Raises

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_negative\_software\_limit** (*lim: Union[int, float], add: int = None*) → None

Set the negative software limit (the maximum position that the motor is allowed to travel to towards the left).

#### Parameters

- **lim** – negative software limit (preset units), value between -1e12 and 0
- **add** – controller address (1 to 31)

#### Raises

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**set\_positive\_software\_limit** (*lim: Union[int, float], add: int = None*) → None

Set the positive software limit (the maximum position that the motor is allowed to travel to towards the right).

#### Parameters

- **lim** – positive software limit (preset units), value between 0 and 1e12
- **add** – controller address (1 to 31)

#### Raises

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**start()**

Opens the communication protocol and applies the config.

**Raises** *SerialCommunicationIOError* – when communication port cannot be opened

**stop()** → None

Stop the device. Close the communication protocol.

**stop\_motion(add: int = None)** → None

Stop a move in progress by decelerating the positioner immediately with the configured acceleration until it stops. If a controller address is provided, stops a move in progress on this controller, else stops the moves on all controllers.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**wait\_until\_motor\_initialized(add: int = None)** → None

Wait until the motor leaves the HOMING state (at which point it should have arrived to the home position).

**Parameters** **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**wait\_until\_move\_finished(add: int = None)** → None

Wait until the motor leaves the MOVING state (at which point it should have arrived to the target position).

**Parameters** **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

```
class hvl_ccb.dev.newport.NewportSMC100PPConfig(address: int = 1, user_position_offset:
    Union[int, float] = 23.987,
    screw_scaling: Union[int, float] = 1, exit_configuration_wait_sec:
    Union[int, float] = 5,
    move_finished_extra_wait_sec:
    Union[int, float] = 1, acceleration: Union[int, float] = 10, backlash_compensation:
    Union[int, float] = 0, hysteresis_compensation:
    Union[int, float] = 0.015, micro_step_per_full_step_factor: int
    = 100, motion_distance_per_full_step:
    Union[int, float] = 0.01, home_search_type: Union[int,
    hvl_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch]
    = <HomeSearch.HomeSwitch: 2>, jerk_time: Union[int, float] = 0.04,
    home_search_velocity: Union[int, float] = 4, home_search_timeout:
    Union[int, float] = 27.5, home_search_polling_interval:
    Union[int, float] = 1, peak_output_current_limit: Union[int,
    float] = 0.4, rs485_address: int
    = 2, negative_software_limit:
    Union[int, float] = -23.5, positive_software_limit: Union[int, float]
    = 25, velocity: Union[int, float] = 4, base_velocity: Union[int, float]
    = 0, stage_configuration: Union[int,
    hvl_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig]
    = <EspStageConfig.EnableEspStageCheck: 3>)
```

Bases: object

Configuration dataclass for the Newport motor controller SMC100PP.

```
class EspStageConfig(*args, **kwds)
```

Bases: aenum.IntEnum

Different configurations to check or not the motor configuration upon power-up.

```
DisableEspStageCheck = 1
```

```
EnableEspStageCheck = 3
```

```
UpdateEspStageInfo = 2
```

```
class HomeSearch(*args, **kwds)
```

Bases: aenum.IntEnum

Different methods for the motor to search its home position during initialization.

```
CurrentPosition = 1
```

```
EndOfRunSwitch = 4
```

```
EndOfRunSwitch_and_Index = 3
```

```
HomeSwitch = 2
HomeSwitch_and_Index = 0
acceleration = 10
address = 1
backlash_compensation = 0
base_velocity = 0
clean_values()
exit_configuration_wait_sec = 5
force_value(fieldname, value)
Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define post_force_value method with same signature as this method to do extra processing
after value has been forced on fieldname.
```

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

```
home_search_polling_interval = 1
home_search_timeout = 27.5
home_search_type = 2
home_search_velocity = 4
hysteresis_compensation = 0.015
is_configdataclass = True
jerk_time = 0.04
classmethod keys() → Sequence[str]
Returns a list of all configdataclass fields key-names.

>Returns a list of strings containing all keys.
```

```
micro_step_per_full_step_factor = 100
motion_distance_per_full_step = 0.01
motor_config
Gather the configuration parameters of the motor into a dictionary.

>Returns dict containing the configuration parameters of the motor
```

```
move_finished_extra_wait_sec = 1
negative_software_limit = -23.5
classmethod optional_defaults() → Dict[str, object]
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified
on instantiation.

>Returns a list of strings containing all optional keys.
```

```
peak_output_current_limit = 0.4
positive_software_limit = 25
```

```

post_force_value(fieldname, value)
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

rs485_address = 2
screw_scaling = 1
stage_configuration = 3
user_position_offset = 23.987
velocity = 4

class hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication(configuration)
Bases: hvl_ccb.comm.serial.SerialCommunication

Specific communication protocol implementation Heinzinger power supplies. Already predefines device-specific protocol parameters in config.

class ControllerErrors(*args, **kwds)
Bases: aenum.Enum

Possible controller errors with values as returned by the device in response to sent commands.

ADDR_INCORRECT = 'B'
CMD_EXEC_ERROR = 'V'
CMD_NOT_ALLOWED = 'D'
CMD_NOT_ALLOWED_CC = 'X'
CMD_NOT_ALLOWED_CONFIGURATION = 'I'
CMD_NOT_ALLOWED_DISABLE = 'J'
CMD_NOT_ALLOWED_HOMING = 'L'
CMD_NOT_ALLOWED_MOVING = 'M'
CMD_NOT_ALLOWED_NOT_REFERENCED = 'H'
CMD_NOT_ALLOWED_PP = 'W'
CMD_NOT_ALLOWED_READY = 'K'
CODE_OR_ADDR_INVALID = 'A'
COM_TIMEOUT = 'S'
DISPLACEMENT_OUT_OF_LIMIT = 'G'
EEPROM_ACCESS_ERROR = 'U'
ESP_STAGE_NAME_INVALID = 'F'
HOME_STARTED = 'E'
NO_ERROR = '@'
PARAM_MISSING_OR_INVALID = 'C'
POSITION_OUT_OF_LIMIT = 'N'

```

**check\_for\_error** (*add: int*) → None  
Ask the Newport controller for the last error it recorded.

This method is called after every command or query.

**Parameters** **add** – controller address (1 to 31)

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query** (*add: int, cmd: str, param: Union[int, float, str, None] = None*) → str

Send a query to the controller, read the answer, and check for errors. The prefix add+cmd is removed from the answer.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **param** – optional parameter (int/float/str) appended to the command

**Returns** the answer from the device without the prefix

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**query\_multiple** (*add: int, cmd: str, prefixes: List[str]*) → List[str]

Send a query to the controller, read the answers, and check for errors. The prefixes are removed from the answers.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **prefixes** – prefixes of each line expected in the answer

**Returns** list of answers from the device without prefix

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**send\_command** (*add: int, cmd: str, param: Union[int, float, str, None] = None*) → None

Send a command to the controller, and check for errors.

**Parameters**

- **add** – the controller address (1 to 31)
- **cmd** – the command to be sent
- **param** – optional parameter (int/float/str) appended to the command

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained
- *NewportControllerError* – if the controller reports an error

**send\_stop** (*add: int*) → None

Send the general stop ST command to the controller, and check for errors.

**Parameters** **add** – the controller address (1 to 31)

**Returns** ControllerErrors reported by Newport Controller

**Raises**

- *SerialCommunicationIOError* – if the com is closed
- *NewportSerialCommunicationError* – if an unexpected answer is obtained

```
class hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig (port: str,  
baudrate:  
int      =  
57600,  
parity:  
Union[str,  
hvl_ccb.comm.serial.SerialCommu  
= <Serial-  
Communi-  
cationPar-  
ity.NONE:  
'N',  
stopbits:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionStop-  
bits.ONE:  
1>, bytes-  
size:  
Union[int,  
hvl_ccb.comm.serial.SerialCommu  
= <Seri-  
alCom-  
munica-  
tionByte-  
size.EIGHTBITS:  
8>, ter-  
minator:  
bytes     =  
b'rn',  
timeout:  
Union[int,  
float]    =  
10)
```

Bases: [hvl\\_ccb.comm.serial.SerialCommunicationConfig](#)

**baudrate = 57600**

Baudrate for Heinzinger power supplies is 9600 baud

**bytesize = 8**

One byte is eight bits long

**force\_value (fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys () → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults () → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**

Heinzinger does not use parity

**classmethod required\_keys () → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**

Heinzinger uses one stop bit

**terminator = b'\r\n'**

The terminator is CR/LF

**timeout = 10**

use 10 seconds timeout as default

**exception hvl\_ccb.dev.newport.NewportSerialCommunicationError**

Bases: Exception

Communication error with the Newport controller.

**class hvl\_ccb.dev.newport.NewportStates (\*args, \*\*kwds)**

Bases: *hvl\_ccb.utils.enum.AutoNumberNameEnum*

States of the Newport controller. Certain commands are allowed only in certain states.

**CONFIG = 3**

**DISABLE = 6**

**HOMING = 2**

**JOGGING = 7**

**MOVING = 5**

**NO\_REF = 1**

**READY = 4**

## **hvl\_ccb.dev.pfeiffer\_tpg module**

Device class for Pfeiffer TPG controllers.

The Pfeiffer TPG control units are used to control Pfeiffer Compact Gauges. Models: TPG 251 A, TPG 252 A, TPG 256A, TPG 261, TPG 262, TPG 361, TPG 362 and TPG 366.

Manufacturer homepage: <https://www.pfeiffer-vacuum.com/en/products/measurement-analysis/> measurement/activeline/controllers/

```
class hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG(com, dev_config=None)
Bases: hvl_ccb.dev.base.SingleCommDevice

Pfeiffer TPG control unit device class

class PressureUnits(*args, **kwds)
Bases: hvl_ccb.utils.enum.NameEnum

Enum of available pressure units for the digital display. “0” corresponds either to bar or to mbar depending on the TPG model. In case of doubt, the unit is visible on the digital display.

Micron = 3
Pascal = 2
Torr = 1
Volt = 5
bar = 0
hPascal = 4
mbar = 0

class SensorStatus
Bases: enum.IntEnum

An enumeration.

Identification_error = 6
No_sensor = 5
Ok = 0
Overrange = 2
Sensor_error = 3
Sensor_off = 4
Underrange = 1

class SensorTypes
Bases: enum.Enum

An enumeration.

CMR = 4
IKR = 2
IKR11 = 2
IKR9 = 2
IMR = 5
None = 7
PBR = 6
PKR = 3
TPR = 1
noSENSOR = 7
noSen = 7
```

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Get the class for the default communication protocol used with this device.

**Returns** the type of the standard communication protocol for this device

**get\_full\_scale\_mbar() → List[Union[int, float]]**

Get the full scale range of the attached sensors

**Returns** full scale range values in mbar, like [0.01, 1, 0.1, 1000, 50000, 10]

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**get\_full\_scale\_unitless() → List[int]**

Get the full scale range of the attached sensors. See lookup table between command and corresponding pressure in the device user manual.

**Returns** list of full scale range values, like [0, 1, 3, 3, 2, 0]

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**identify\_sensors() → None**

Send identification request TID to sensors on all channels.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**measure(channel: int) → Tuple[str, float]**

Get the status and measurement of one sensor

**Parameters** `channel` – int channel on which the sensor is connected, with  $1 \leq \text{channel} \leq \text{number\_of\_sensors}$

**Returns** measured value as float if measurement successful, sensor status as string if not

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**measure\_all() → List[Tuple[str, float]]**

Get the status and measurement of all sensors (this command is not available on all models)

**Returns** list of measured values as float if measurements successful, and or sensor status as strings if not

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**number\_of\_sensors**

**set\_display\_unit** (*unit: Union[str, hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits]*) → None  
Set the unit in which the measurements are shown on the display.

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**set\_full\_scale\_mbar** (*fsr: List[Union[int, float]]*) → None  
Set the full scale range of the attached sensors (in unit mbar)

**Parameters** *fsr* – full scale range values in mbar, for example [0.01, 1000]

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**set\_full\_scale\_unitless** (*fsr: List[int]*) → None

Set the full scale range of the attached sensors. See lookup table between command and corresponding pressure in the device user manual.

**Parameters** *fsr* – list of full scale range values, like [0, 1, 3, 3, 2, 0]

**Raises**

- *SerialCommunicationIOError* – when communication port is not opened
- *PfeifferTPGError* – if command fails

**start()** → None

Start this device. Opens the communication protocol, and identify the sensors.

**Raises** *SerialCommunicationIOError* – when communication port cannot be opened

**stop()** → None

Stop the device. Closes also the communication protocol.

**unit**

The pressure unit of readings is always mbar, regardless of the display unit.

**class** *hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig* (*model: Union[str, hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig.Model] = <Model.TPG25xA: {1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6, 50000: 7, 0.1: 8}>*)

Bases: *object*

Device configuration dataclass for Pfeiffer TPG controllers.

**class Model** (\**args*, \*\**kwargs*)

Bases: *hvl\_ccb.utils.enum.NameEnum*

**TPG25xA** = {0.1: 8, 1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6}

**TPGx6x** = {0.01: 0, 0.1: 1, 1: 2, 10: 3, 100: 4, 1000: 5, 2000: 6, 5000: 7, 10000: 8}

**is\_valid\_scale\_range\_reversed\_str** (*v: str*) → bool

Check if given string represents a valid reversed scale range of a model.

**Parameters** *v* – Reversed scale range string.

**Returns** *True* if valid, *False* otherwise.

**clean\_values()**

**force\_value**(*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True****classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**model = {0.1: 8, 1: 0, 10: 1, 100: 2, 1000: 3, 2000: 4, 5000: 5, 10000: 6, 50000: 7}****classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGError**

Bases: Exception

Error with the Pfeiffer TPG Controller.

**class hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunication(*configuration*)**

Bases: *hvl\_ccb.comm.serial.SerialCommunication*

Specific communication protocol implementation for Pfeiffer TPG controllers. Already predefines device-specific protocol parameters in config.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**query**(*cmd: str*) → str

Send a query, then read and returns the first line from the com port.

**Parameters** **cmd** – query message to send to the device

**Returns** first line read on the com

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if the device does not acknowledge the command or if the answer from the device is empty

**send\_command**(*cmd: str*) → None

Send a command to the device and check for acknowledgement.

**Parameters** **cmd** – command to send to the device

**Raises**

- **SerialCommunicationIOError** – when communication port is not opened
- **PfeifferTPGError** – if the answer from the device differs from the expected acknowledgement character ‘chr(6)’.

```
class hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig(port:  
                           str, baudrate:  
                           int      =  
                           9600,  
                           parity:  
                           Union[str,  
                           hvl_ccb.comm.serial.SerialComm  
                           = <SerialCom-  
                           munica-  
                           tionPar-  
                           ity.NONE:  
                           'N'>,  
                           stopbits:  
                           Union[int,  
                           hvl_ccb.comm.serial.SerialComm  
                           = <SerialCom-  
                           munica-  
                           tionStop-  
                           bits.ONE:  
                           1>,  
                           bytesize:  
                           Union[int,  
                           hvl_ccb.comm.serial.SerialComm  
                           = <SerialCom-  
                           munica-  
                           tionByte-  
                           size.EIGHTBITS:  
                           8>, terminator:  
                           bytes      =  
                           b'm', timeout:  
                           Union[int,  
                           float]     =  
                           3)
```

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

**baudrate = 9600**

Baudrate for Pfeiffer TPG controllers is 9600 baud

**bytesize = 8**

One byte is eight bits long

**force\_value (fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing

after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**

Pfeiffer TPG controllers do not use parity

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**stopbits = 1**

Pfeiffer TPG controllers use one stop bit

**terminator = b'\r\n'**

The terminator is <CR><LF>

**timeout = 3**

use 3 seconds timeout as default

## hvl\_ccb.dev.rs\_rto1024 module

Python module for the Rhode & Schwarz RTO 1024 oscilloscope. The communication to the device is through VISA, type TCPIP / INSTR.

```
class hvl_ccb.dev.rs_rto1024.RTO1024(com: Union[hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunication,
                                              hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig,
                                              dict], dev_config: Union[hvl_ccb.dev.rs_rto1024.RTO1024Config,
                                              dict])
```

Bases: [hvl\\_ccb.dev.visa.VisaDevice](#)

Device class for the Rhode & Schwarz RTO 1024 oscilloscope.

**class TriggerModes(\*args, \*\*kwds)**

Bases: [hvl\\_ccb.utils.enum.AutoNumberNameEnum](#)

Enumeration for the three available trigger modes.

**AUTO = 1**

**FREERUN = 3**

**NORMAL = 2**

**names = <bound method RTO1024.TriggerModes.names of <aenum 'TriggerModes'>>**

**backup\_waveform** (*filename: str*) → None

Backup a waveform file from the standard directory specified in the device configuration to the standard backup destination specified in the device configuration. The filename has to be specified without .bin or path.

**Parameters** **filename** – The waveform filename without extension and path

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**static default\_com\_cls()**

Return the default communication protocol for this device type, which is VisaCommunication.

**Returns** the VisaCommunication class

**file\_copy** (*source: str, destination: str*) → None

Copy a file from one destination to another on the oscilloscope drive. If the destination file already exists, it is overwritten without notice.

**Parameters**

- **source** – absolute path to the source file on the DSO filesystem
- **destination** – absolute path to the destination file on the DSO filesystem

**Raises** *RTO1024Error* – if the operation did not complete

**get\_timestamps()** → List[float]

Gets the timestamps of all recorded frames in the history and returns them as a list of floats.

**Returns** list of timestamps in [s]

**Raises** *RTO1024Error* – if the timestamps are invalid

**list\_directory** (*path: str*) → List[Tuple[str, str, int]]

List the contents of a given directory on the oscilloscope filesystem.

**Parameters** **path** – is the path to a folder

**Returns** a list of filenames in the given folder

**load\_configuration** (*filename: str*) → None

Load current settings from a configuration file. The filename has to be specified without base directory and ‘.dfl’ extension.

**Information from the manual** *ReCaLL* calls up the instrument settings from an intermediate memory identified by the specified number. The instrument settings can be stored to this memory using the command \*SAV with the associated number. It also activates the instrument settings which are stored in a file and loaded using *MMEMemory:LOAD:STATe*.

**Parameters** **filename** – is the name of the settings file without path and extension

**local\_display** (*state: bool*) → None

Enable or disable local display of the scope.

**Parameters** **state** – is the desired local display state

**prepare\_ultra\_segmentation()** → None

Make ready for a new acquisition in ultra segmentation mode. This function does one acquisition without ultra segmentation to clear the history and prepare for a new measurement.

**run\_continuous\_acquisition()** → None

Start acquiring continuously.

**run\_single\_acquisition()** → None

Start a single or Nx acquisition.

**save\_configuration** (*filename: str*) → None

Save the current oscilloscope settings to a file. The filename has to be specified without path and '.dfl' extension, the file will be saved to the configured settings directory.

**Information from the manual** *SAVe* stores the current instrument settings under the specified number in an intermediate memory. The settings can be recalled using the command *\*RCL* with the associated number. To transfer the stored instrument settings to a file, use *MMEMory:STORe:STATE*.

**Parameters** **filename** – is the name of the settings file without path and extension

**save\_waveform\_history** (*filename: str, channel: int, waveform: int = 1*) → None

Save the history of one channel and one waveform to a .bin file. This function is used after an acquisition using sequence trigger mode (with or without ultra segmentation) was performed.

**Parameters**

- **filename** – is the name (without extension) of the file
- **channel** – is the channel number
- **waveform** – is the waveform number (typically 1)

**Raises** *RTO1024Error* – if storing waveform times out

**set\_acquire\_length** (*timerange: float*) → None

Defines the time of one acquisition, that is the time across the 10 divisions of the diagram.

- Range: 250E-12 ... 500 [s]
- Increment: 1E-12 [s]
- \*RST = 0.5 [s]

**Parameters** **timerange** – is the time for one acquisition. Range: 250e-12 ... 500 [s]

**set\_channel\_position** (*channel: int, position: float*) → None

Sets the vertical position of the indicated channel as a graphical value.

- Range: -5.0 ... 5.0 [div]
- Increment: 0.02
- \*RST = 0

**Parameters**

- **channel** – is the channel number (1..4)
- **position** – is the position. Positive values move the waveform up, negative values move it down.

**set\_channel\_range** (*channel: int, v\_range: float*) → None

Sets the voltage range across the 10 vertical divisions of the diagram. Use the command alternatively instead of *set\_channel\_scale*.

- Range for range: Depends on attenuation factors and coupling. With 1:1 probe and external attenuations and 50 Ω input coupling, the range is 10 mV to 10 V. For 1 MΩ input coupling, it is 10 mV to 100 V. If the probe and/or external attenuation is changed, multiply the range values by the attenuation factors.
- Increment: 0.01

- \*RST = 0.5

#### Parameters

- **channel1** – is the channel number (1..4)
- **v\_range** – is the vertical range [V]

**set\_channel\_scale** (*channel: int, scale: float*) → None

Sets the vertical scale for the indicated channel. The scale value is given in volts per division.

- Range for scale: depends on attenuation factor and coupling. With 1:1 probe and external attenuations and  $50\ \Omega$  input coupling, the vertical scale (input sensitivity) is 1 mV/div to 1 V/div. For  $1\ M\Omega$  input coupling, it is 1 mV/div to 10 V/div. If the probe and/or external attenuation is changed, multiply the values by the attenuation factors to get the actual scale range.
- Increment: 1e-3
- \*RST = 0.05

See also: set\_channel\_range

#### Parameters

- **channel1** – is the channel number (1..4)
- **scale** – is the vertical scaling [V/div]

**set\_channel\_state** (*channel: int, state: bool*) → None

Switches the channel signal on or off.

#### Parameters

- **channel1** – is the input channel (1..4)
- **state** – is True for on, False for off

**set\_reference\_point** (*percentage: int*) → None

Sets the reference point of the time scale in % of the display. If the “Trigger offset” is zero, the trigger point matches the reference point. ReferencePoint = zero pint of the time scale

- Range: 0 … 100 [%]
- Increment: 1 [%]
- \*RST = 50 [%]

Parameters **percentage** – is the reference in %

**set\_repetitions** (*number: int*) → None

Set the number of acquired waveforms with RUN Nx SINGLE. Also defines the number of waveforms used to calculate the average waveform.

- Range: 1 … 16777215
- Increment: 10
- \*RST = 1

Parameters **number** – is the number of waveforms to acquire

**set\_trigger\_level** (*channel: int, level: float, event\_type: int = 1*) → None

Sets the trigger level for the specified event and source.

- Range: -10 to 10 V
- Increment: 1e-3 V
- \*RST = 0 V

#### Parameters

- **channel** – indicates the trigger source.
  - 1..4 = channel 1 to 4, available for all event types 1..3
  - 5 = external trigger input on the rear panel for analog signals, available for A-event type = 1
  - 6..9 = not available
- **level** – is the voltage for the trigger level in [V].
- **event\_type** – is the event type. 1: A-Event, 2: B-Event, 3: R-Event

**set\_trigger\_mode** (*mode*: Union[str, hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes]) → None  
 Sets the trigger mode which determines the behavior of the instrument if no trigger occurs.

**Parameters** **mode** – is either auto, normal, or freerun.

**Raises** **RTO1024Error** – if an invalid triggermode is selected

**set\_trigger\_source** (*channel*: int, *event\_type*: int = 1) → None  
 Set the trigger (Event A) source channel.

#### Parameters

- **channel** – is the channel number (1..4)
- **event\_type** – is the event type. 1: A-Event, 2: B-Event, 3: R-Event

**start()** → None  
 Start the RTO1024 oscilloscope and bring it into a defined state and remote mode.

**stop()** → None  
 Stop the RTO1024 oscilloscope, reset events and close communication. Brings back the device to a state where local operation is possible.

**stop\_acquisition()** → None  
 Stop any acquisition.

```
class hvl_ccb.dev.rs_rto1024.RTO1024Config(waveforms_path: str, settings_path: str, backup_path: str, spoll_interval: Union[int, float] = 0.5, spoll_start_delay: Union[int, float] = 2, command_timeout_seconds: Union[int, float] = 60, wait_sec_short_pause: Union[int, float] = 0.1, wait_sec_enable_history: Union[int, float] = 1, wait_sec_post_acquisition_start: Union[int, float] = 2)
```

Bases: [hvl\\_ccb.dev.visa.VisaDeviceConfig](#), [hvl\\_ccb.dev.rs\\_rto1024.\\_RTO1024ConfigDefaultsBase](#), [hvl\\_ccb.dev.rs\\_rto1024.\\_RTO1024ConfigBase](#)

Configdataclass for the RTO1024 device.

**force\_value** (*fieldname*, *value*)  
 Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**exception hvl\_ccb.dev.rs\_rto1024.RTO1024Error**

Bases: Exception

**class hvl\_ccb.dev.rs\_rto1024.RTO1024VisaCommunication(configuration)**

Bases: *hvl\_ccb.comm.visa.VisaCommunication*

Specialization of VisaCommunication for the RTO1024 oscilloscope

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**class hvl\_ccb.dev.rs\_rto1024.RTO1024VisaCommunicationConfig(host: str, interface\_type: Union[str, hvl\_ccb.comm.visa.VisaCommunicationConfig] = <InterfaceType.TCPIP\_INSTR: 2>, board: int = 0, port: int = 5025, timeout: int = 5000, chunk\_size: int = 204800, open\_timeout: int = 1000, write\_termination: str = '\n', read\_termination: str = '\n', visa\_backend: str = '')**

Bases: *hvl\_ccb.comm.visa.VisaCommunicationConfig*

Configuration dataclass for VisaCommunication with specifications for the RTO1024 device class.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**interface\_type = 2****classmethod keys()** → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults()** → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys()** → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**hvl\_ccb.dev.se\_il2t module**

Device class for controlling a Schneider Electric ILS2T stepper drive over modbus TCP.

**class hvl\_ccb.dev.se\_il2t.ILS2T** (*com*, *dev\_config=None*)

Bases: *hvl\_ccb.dev.base.SingleCommDevice*

Schneider Electric ILS2T stepper drive class.

**ACTION\_JOG\_VALUE = 0**

The single action value for *ILS2T.Mode.JOG*

**class ActionsPtp**

Bases: *enum.IntEnum*

Allowed actions in the point to point mode (*ILS2T.Mode.PTP*).

**ABSOLUTE\_POSITION = 0****RELATIVE\_POSITION\_MOTOR = 2****RELATIVE\_POSITION\_TARGET = 1****DEFAULT\_IO\_SCANNING\_CONTROL\_VALUES = {'action': 2, 'continue\_after\_stop\_cu': 0, 'dis**

Default IO Scanning control mode values

**class Mode**

Bases: *enum.IntEnum*

ILS2T device modes

**JOG = 1****PTP = 3**

```
class Ref16Jog
    Bases: enum.Flag

    Allowed values for ILS2T ref_16 register (the shown values are the integer representation of the bits), all
    in Jog mode = 1

    FAST = 4
    NEG = 2
    NEG_FAST = 6
    NONE = 0
    POS = 1
    POS_FAST = 5

RegAddr
    Modbus Register Adresses
    alias of ILS2TRegAddr

RegDatatype
    Modbus Register Datatypes
    alias of ILS2TRegDatatype

class State
    Bases: enum.IntEnum

    State machine status values

    ON = 6
    QUICKSTOP = 7
    READY = 4

static config_cls()
    Return the default configdataclass class.

    Returns a reference to the default configdataclass class

static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

disable(log_warn: bool = True, wait_sec_max: Optional[int] = None) → bool
    Disable the driver of the stepper motor and enable the brake.

    Note: the driver cannot be disabled if the motor is still running.

    Parameters
        • log_warn – if log a warning in case the motor cannot be disabled.
        • wait_sec_max – maximal wait time for the motor to stop running and to disable it; by
            default, with None, use a config value

    Returns True if disable request could and was sent, False otherwise.

do_ioscanning_write(**kwargs) → None
    Perform a write operation using IO Scanning mode.

    Parameters kwargs – Keyword-argument list with options to send, remaining are taken from
        the defaults.
```

**enable()** → None  
 Enable the driver of the stepper motor and disable the brake.

**execute\_absolute\_position(position: int)** → bool  
 Execute a absolute position change, i.e. enable motor, perform absolute position change, wait until done and disable motor afterwards.

Check position at the end if wrong do not raise error; instead just log and return check result.

**Parameters** **position** – absolute position of motor in user defined steps.

**Returns** *True* if actual position is as expected, *False* otherwise.

**execute\_relative\_step(steps: int)** → bool  
 Execute a relative step, i.e. enable motor, perform relative steps, wait until done and disable motor afterwards.

Check position at the end if wrong do not raise error; instead just log and return check result.

**Parameters** **steps** – Number of steps.

**Returns** *True* if actual position is as expected, *False* otherwise.

**get\_dc\_volt()** → float  
 Read the DC supply voltage of the motor.

**Returns** DC input voltage.

**get\_error\_code()** → Dict[int, Dict[str, Any]]  
 Read all messages in fault memory. Will read the full error message and return the decoded values. At the end the fault memory of the motor will be deleted. In addition, `reset_error` is called to re-enable the motor for operation.

**Returns** Dictionary with all information

**get\_position()** → int  
 Read the position of the drive and store into status.

**Returns** Position step value

**get\_status()** → Dict[str, int]  
 Perform an IO Scanning read and return the status of the motor.

**Returns** dict with status information.

**get\_temperature()** → int  
 Read the temperature of the motor.

**Returns** Temperature in degrees Celsius.

**jog\_run(direction: bool = True, fast: bool = False)** → None  
 Slowly turn the motor in positive direction.

**jog\_stop()** → None  
 Stop turning the motor in Jog mode.

**quickstop()** → None  
 Stops the motor with high deceleration rate and falls into error state. Reset with `reset_error` to recover into normal state.

**reset\_error()** → None  
 Resets the motor into normal state after quick stop or another error occurred.

**set\_jog\_speed(slow: int = 60, fast: int = 180)** → None  
 Set the speed for jog mode. Default values correspond to startup values of the motor.

**Parameters**

- **slow** – RPM for slow jog mode.
- **fast** – RPM for fast jog mode.

**set\_max\_acceleration** (*rpm\_minute: int*) → None

Set the maximum acceleration of the motor.

**Parameters** **rpm\_minute** – revolution per minute per minute

**set\_max\_deceleration** (*rpm\_minute: int*) → None

Set the maximum deceleration of the motor.

**Parameters** **rpm\_minute** – revolution per minute per minute

**set\_max\_rpm** (*rpm: int*) → None

Set the maximum RPM.

**Parameters** **rpm** – revolution per minute ( 0 < rpm <= RPM\_MAX)

**Raises** *ILS2TException* – if RPM is out of range

**set\_ramp\_type** (*ramp\_type: int = -1*) → None

Set the ramp type. There are two options available: 0: linear ramp -1: motor optimized ramp

**Parameters** **ramp\_type** – 0: linear ramp | -1: motor optimized ramp

**start** () → None

Start this device.

**stop** () → None

Stop this device. Disables the motor (applies brake), disables access and closes the communication protocol.

**user\_steps** (*steps: int = 16384, revolutions: int = 1*) → None

Define steps per revolution. Default is 16384 steps per revolution. Maximum precision is 32768 steps per revolution.

**Parameters**

- **steps** – number of steps in *revolutions*.
- **revolutions** – number of revolutions corresponding to *steps*.

**write\_absolute\_position** (*position: int*) → None

Write instruction to turn the motor until it reaches the absolute position. This function does not enable or disable the motor automatically.

**Parameters** **position** – absolute position of motor in user defined steps.

**write\_relative\_step** (*steps: int*) → None

Write instruction to turn the motor the relative amount of steps. This function does not enable or disable the motor automatically.

**Parameters** **steps** – Number of steps to turn the motor.

```
class hvl_ccb.dev.se_ilts2t.ILS2TConfig(rpm_max_init: numbers.Integral = 1500,
                                         wait_sec_post_enable: Union[float, int] = 1,
                                         wait_sec_max_disable: Union[float, int] = 10,
                                         wait_sec_post_cannot_disable: Union[float, int] =
                                         1, wait_sec_post_relative_step: Union[float, int] =
                                         2, wait_sec_post_absolute_position: Union[float,
                                         int] = 2)
```

Bases: object

Configuration for the ILS2T stepper motor device.

**clean\_values()**

**force\_value(fieldname, value)**

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**is\_configdataclass = True**

**classmethod keys() → Sequence[str]**

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults() → Dict[str, object]**

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys() → Sequence[str]**

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**rpm\_max\_init = 1500**

initial maximum RPM for the motor, can be set up to 3000 RPM. The user is allowed to set a new max RPM at runtime using *ILS2T.set\_max\_rpm()*, but the value must never exceed this configuration setting.

**wait\_sec\_max\_disable = 10**

**wait\_sec\_post\_absolute\_position = 2**

**wait\_sec\_post\_CANNOT\_DISABLE = 1**

**wait\_sec\_post\_enable = 1**

**wait\_sec\_post\_relative\_step = 2**

**exception hvl\_ccb.dev.se\_ilst2t.ILS2TException**

Bases: Exception

Exception to indicate problems with the SE ILS2T stepper motor.

**class hvl\_ccb.dev.se\_ilst2t.ILS2TModbusTcpCommunication(configuration)**

Bases: *hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication*

Specific implementation of Modbus/TCP for the Schneider Electric ILS2T stepper motor.

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

```
class hvl_ccb.dev.se_ilst2t.ILS2TModbusTcpCommunicationConfig(host: str, unit: int
= 255, port: int =
502)
```

Bases: *hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig*

Configuration dataclass for Modbus/TCP communication specific for the Schneider Electric ILS2T stepper motor.

**force\_value** (*fieldname*, *value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

#### Parameters

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys** () → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults** () → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys** () → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**unit = 255**

The unit has to be 255 such that IO scanning mode works.

```
class hvl_ccb.dev.se_ilst2t.ILS2TRegAddr
```

Bases: enum.IntEnum

Modbus Register Addresses for Schneider Electric ILS2T stepper drive.

**ACCESS\_ENABLE = 282**

**FLT\_INFO = 15362**

**FLT\_MEM\_DEL = 15112**

**FLT\_MEM\_RESET = 15114**

**IO\_SCANNING = 6922**

**JOGN\_FAST = 10506**

**JOGN\_SLOW = 10504**

**POSITION = 7706**

**RAMP\_ACC = 1556**

**RAMP\_DECEL = 1558**

**RAMP\_N\_MAX = 1554**

**RAMP\_TYPE = 1574**

```
SCALE = 1550
TEMP = 7200
VOLT = 7198

class hvl_ccb.dev.se_ilst2t.ILS2TRegDatatype(*args, **kwds)
Bases: aenum.Enum
```

Modbus Register Datatypes for Schneider Electric ILS2T stepper drive.

From the manual of the drive:

datatype	byte	min	max
INT8	1 Byte	-128	127
UINT8	1 Byte	0	255
INT16	2 Byte	-32_768	32_767
UINT16	2 Byte	0	65_535
INT32	4 Byte	-2_147_483_648	2_147_483_647
UINT32	4 Byte	0	4_294_967_295
BITS	just 32bits	N/A	N/A

```
INT32 = (-2147483648, 2147483647)
is_in_range(value: int) → bool

exception hvl_ccb.dev.se_ilst2t.IoScanningModeValueError
Bases: hvl_ccb.dev.se_ilst2t.ILS2TException

Exception to indicate that the selected IO scanning mode is invalid.

exception hvl_ccb.dev.se_ilst2t.ScalingFactorValueError
Bases: hvl_ccb.dev.se_ilst2t.ILS2TException

Exception to indicate that a scaling factor value is invalid.
```

## hvl\_ccb.dev.sst\_luminox module

Device class for a SST Luminox Oxygen sensor. This device can measure the oxygen concentration between 0 % and 25 %.

Furthermore, it measures the barometric pressure and internal temperature. The device supports two operating modes: in streaming mode the device measures all parameters every second, in polling mode the device measures only after a query.

Technical specification and documentation for the device can be found at the manufacturer's page: <https://www.sstsensing.com/product/luminox-optical-oxygen-sensors-2/>

```
class hvl_ccb.dev.sst_luminox.Luminox(com, dev_config=None)
Bases: hvl_ccb.dev.base.SingleCommDevice

Luminox oxygen sensor device class.

activate_output(mode: hvl_ccb.dev.sst_luminox.LuminoxOutputMode) → None
    activate the selected output mode of the Luminox Sensor. :param mode: polling or streaming

static config_cls()
    Return the default configdataclass class.

Returns a reference to the default configdataclass class
```

```
static default_com_cls()
    Get the class for the default communication protocol used with this device.

    Returns the type of the standard communication protocol for this device

output = None

query_single_measurement(measurement: Union[str, hvl_ccb.dev.sst_luminox.LuminoxMeasurementType])
    → Union[float, int, str]
    query one single measurement value in polling mode.

    Parameters measurement – type of measurement
    Returns value of requested measurement
    Raises
        • ValueError – when a wrong key for LuminoxMeasurementType is provided
        • LuminoxOutputModeError – when polling mode is not activated
        • LuminoxMeasurementTypeError – when expected measurement value is not read

read_stream() → dict
    read values of Luminox in stream mode. Convert the single string into sepearate numbers

    Returns dictionary with all values
    Raises
        • LuminoxOutputModeError – when streaming mode is not activated
        • LuminoxMeasurementTypeError – when any of expected measurement values is
            not read

start() → None
    Start this device. Opens the communication protocol.

stop() → None
    Stop the device. Closes also the communication protocol.

class hvl_ccb.dev.sst_luminox.LuminoxConfig(wait_sec_post_activate: Union[int, float] =
                                            0.5)
Bases: object

Configuration for the SST Luminox oxygen sensor.

clean_values()

force_value(fieldname, value)
    Forces a value to a dataclass field despite the class being frozen.

    NOTE: you can define post_force_value method with same signature as this method to do extra processing
    after value has been forced on fieldname.

    Parameters
        • fieldname – name of the field
        • value – value to assign

is_configdataclass = True

classmethod keys() → Sequence[str]
    Returns a list of all configdataclass fields key-names.

    Returns a list of strings containing all keys.
```

```

classmethod optional_defaults() → Dict[str, object]
    Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

    Returns a list of strings containing all optional keys.

classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

    Returns a list of strings containing all required keys.

wait_sec_post_activate = 0.5

class hvl_ccb.dev.sst_luminox.LuminoxMeasurementType (*args, **kwds)
    Bases: hvl\_ccb.utils.enum.ValueEnum

    measurement type.

    barometric_pressure = 'P'

    command

    date_of_manufacture = '# 0'

    parse_read_measurement_value(read_txt: str) → Union[float, int, str]

    partial_pressure_o2 = 'O'

    percent_o2 = '%'

    sensor_status = 'e'

    serial_number = '# 1'

    software_revision = '# 2'

    temperature_sensor = 'T'

exception hvl_ccb.dev.sst_luminox.LuminoxMeasurementTypeError
    Bases: Exception

    Wrong measurement type for requested data

class hvl_ccb.dev.sst_luminox.LuminoxOutputMode
    Bases: enum.Enum

    output mode.

    polling = 1

    streaming = 0

exception hvl_ccb.dev.sst_luminox.LuminoxOutputModeError
    Bases: Exception

    Wrong output mode for requested data

class hvl_ccb.dev.sst_luminox.LuminoxSerialCommunication (configuration)
    Bases: hvl\_ccb.comm.serial.SerialCommunication

    Specific communication protocol implementation for the SST Luminox oxygen sensor. Already predefines device-specific protocol parameters in config.

    static config_cls()
        Return the default configdataclass class.

        Returns a reference to the default configdataclass class

```

```
class hvl_ccb.dev.sst_luminox.LuminoxSerialCommunicationConfig(port: str, baudrate: int = 9600, parity: Union[str, hvl_ccb.comm.serial.SerialCommunicationParity.NONE: 'N', stopbits: Union[int, hvl_ccb.comm.serial.SerialCommunicationStopbits.ONE: 1], bytesize: Union[int, hvl_ccb.comm.serial.SerialCommunicationBytesize.EIGHTBITS: 8], terminator: bytes = b'\r\n', timeout: Union[int, float] = 3))
```

Bases: *hvl\_ccb.comm.serial.SerialCommunicationConfig*

**baudrate = 9600**  
Baudrate for SST Luminox is 9600 baud

**bytesize = 8**  
One byte is eight bits long

**force\_value (fieldname, value)**  
Forces a value to a dataclass field despite the class being frozen.  
NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys () → Sequence[str]**  
Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults () → Dict[str, object]**  
Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**parity = 'N'**  
SST Luminox does not use parity

---

```
classmethod required_keys() → Sequence[str]
    Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

Returns a list of strings containing all required keys.

stopbits = 1
    SST Luminox does use one stop bit

terminator = b'\r\n'
    The terminator is CR LF

timeout = 3
    use 3 seconds timeout as default
```

## hvl\_ccb.dev.visa module

```
class hvl_ccb.dev.visa.VisaDevice(com: Union[hvl_ccb.comm.visa.VisaCommunication,
                                              hvl_ccb.comm.visa.VisaCommunicationConfig, dict],
                                    dev_config: Union[hvl_ccb.dev.visa.VisaDeviceConfig, dict,
                                                      None] = None)
Bases: hvl_ccb.dev.base.SingleCommDevice

Device communicating over the VISA protocol using VisaCommunication.

static config_cls()
    Return the default configdataclass class.

Returns a reference to the default configdataclass class

static default_com_cls() → Type[hvl_ccb.comm.visa.VisaCommunication]
    Return the default communication protocol for this device type, which is VisaCommunication.

Returns the VisaCommunication class

get_error_queue() → str
    Read out error queue and logs the error.

Returns Error string

get_identification() → str
    Queries “*IDN?” and returns the identification string of the connected device.

Returns the identification string of the connected device

reset() → None
    Send “*RST” and “*CLS” to the device. Typically sets a defined state.

spoll_handler()
    Reads the status byte and decodes it. The status byte STB is defined in IEEE 488.2. It provides a rough overview of the instrument status.

Returns

start() → None
    Start the VisaDevice. Sets up the status poller and starts it.

Returns

stop() → None
    Stop the VisaDevice. Stops the polling thread and closes the communication protocol.

Returns
```

**wait\_operation\_complete**(*timeout: Optional[float] = None*) → bool

Waits for a operation complete event. Returns after timeout [s] has expired or the operation complete event has been caught.

**Parameters** **timeout** – Time in seconds to wait for the event; *None* for no timeout.

**Returns** True, if OPC event is caught, False if timeout expired

**class** hvl\_ccb.dev.visa.**VisaDeviceConfig**(*spoll\_interval: Union[int, float] = 0.5,*  
*spoll\_start\_delay: Union[int, float] = 2*)  
Bases: hvl\_ccb.dev.visa.\_VisaDeviceConfigDefaultsBase, hvl\_ccb.dev.visa.\_VisaDeviceConfigBase

Configdataclass for a VISA device.

**force\_value**(*fieldname, value*)

Forces a value to a dataclass field despite the class being frozen.

NOTE: you can define *post\_force\_value* method with same signature as this method to do extra processing after *value* has been forced on *fieldname*.

**Parameters**

- **fieldname** – name of the field
- **value** – value to assign

**classmethod keys**() → Sequence[str]

Returns a list of all configdataclass fields key-names.

**Returns** a list of strings containing all keys.

**classmethod optional\_defaults**() → Dict[str, object]

Returns a list of all configdataclass fields, that have a default value assigned and may be optionally specified on instantiation.

**Returns** a list of strings containing all optional keys.

**classmethod required\_keys**() → Sequence[str]

Returns a list of all configdataclass fields, that have no default value assigned and need to be specified on instantiation.

**Returns** a list of strings containing all required keys.

**class** hvl\_ccb.dev.visa.**VisaStatusPoller**(*target: Callable, interval: float = 0.5, start\_delay:*  
*float = 5*)

Bases: threading.Thread

Thread to periodically poll the status byte of a VISA device.

**run**() → None

Threaded method.

**stop**() → None

Gracefully stop the poller.

## Module contents

Devices subpackage.

## [hvl\\_ccb.utils package](#)

### Submodules

#### [hvl\\_ccb.utils.enum module](#)

**class** hvl\_ccb.utils.enum.**AutoNumberNameEnum**(\*args, \*\*kwds)

Bases: [hvl\\_ccb.utils.enum.StrEnumBase](#), aenum.AutoNumberEnum

Auto-numbered enum with names used as string representation, and with lookup and equality based on this representation.

**class** hvl\_ccb.utils.enum.**NameEnum**(\*args, \*\*kwds)

Bases: [hvl\\_ccb.utils.enum.StrEnumBase](#)

Enum with names used as string representation, and with lookup and equality based on this representation.

**class** hvl\_ccb.utils.enum.**StrEnumBase**(\*args, \*\*kwds)

Bases: aenum.Enum

String representation-based equality and lookup.

**class** hvl\_ccb.utils.enum.**ValueEnum**(\*args, \*\*kwds)

Bases: [hvl\\_ccb.utils.enum.StrEnumBase](#)

Enum with string representation of values used as string representation, and with lookup and equality based on this representation.

Attention: to avoid errors, best use together with *unique* enum decorator.

#### [hvl\\_ccb.utils.typing module](#)

Additional Python typing module utilities

`hvl_ccb.utils.typing.check_generic_type(value, type_, name='instance')`

Check if `value` is of a generic type `type_`. Raises `TypeError` if it's not.

##### Parameters

- **name** – name to report in case of an error
- **value** – value to check
- **type** – generic type to check against

`hvl_ccb.utils.typing.is_generic(type_)`

Check if class is a user-defined generic type, for example `Union[int, float]` but not `List`.

##### Parameters `type` – type to check

`hvl_ccb.utils.typing.is_type_hint(type_)`

Check if class is a generic type, for example `Union` or `List[int]`

##### Parameters `type` – type to check

### Module contents

#### 4.1.2 Submodules

### hvl\_ccb.configuration module

Facilities providing classes for handling configuration for communication protocols and devices.

**class** hvl\_ccb.configuration.ConfigurationMixin (*configuration*)

Bases: abc.ABC

Mixin providing configuration to a class.

**config**

ConfigDataclass property.

**Returns** the configuration

**static config\_cls()**

Return the default configdataclass class.

**Returns** a reference to the default configdataclass class

**configuration\_save\_json** (*path: str*) → None

Save current configuration as JSON file.

**Parameters** **path** – path to the JSON file.

**classmethod from\_json** (*filename: str*)

Instantiate communication protocol using configuration from a JSON file.

**Parameters** **filename** – Path and filename to the JSON configuration

hvl\_ccb.configuration.configdataclass (*direct\_decoration=None, frozen=True*)

Decorator to make a class a configdataclass. Types in these dataclasses are enforced. Implement a function clean\_values(self) to do additional checking on value ranges etc.

It is possible to inherit from a configdataclass and re-decorate it with @configdataclass. In a subclass, default values can be added to existing fields. Note: adding additional non-default fields is prone to errors, since the order has to be respected through the whole chain (first non-default fields, only then default-fields).

**Parameters** **frozen** – defaults to True. False allows to later change configuration values. Attention: if configdataclass is not frozen and a value is changed, typing is not enforced anymore!

### hvl\_ccb.experiment\_manager module

Main module containing the top level ExperimentManager class. Inherit from this class to implement your own experiment functionality in another project and it will help you start, stop and manage your devices.

**exception** hvl\_ccb.experiment\_manager.ExperimentError

Bases: Exception

Exception to indicate that the current status of the experiment manager is on ERROR and thus no operations can be made until reset.

**class** hvl\_ccb.experiment\_manager.ExperimentManager (*devices: Dict[str, hvl\_ccb.dev.base.Device]*)

Bases: hvl\_ccb.dev.base.DeviceSequenceMixin

Experiment Manager can start and stop communication protocols and devices. It provides methods to queue commands to devices and collect results.

**add\_device** (*name: str, device: hvl\_ccb.dev.base.Device*) → None

Add a new device to the manager. If the experiment is running, automatically start the device. If a device with this name already exists, raise an exception.

**Parameters**

- **name** – is the name of the device.
- **device** – is the instantiated Device object.

**Raises** *DeviceExistingException* –

**finish()** → None

Stop experimental setup, stop all devices.

**is\_error()** → bool

Returns true, if the status of the experiment manager is *error*.

**Returns** True if on error, false otherwise

**is\_finished()** → bool

Returns true, if the status of the experiment manager is *finished*.

**Returns** True if finished, false otherwise

**is\_running()** → bool

Returns true, if the status of the experiment manager is *running*.

**Returns** True if running, false otherwise

**run()** → None

Start experimental setup, start all devices.

**start()** → None

Alias for ExperimentManager.run()

**status**

Get experiment status.

**Returns** experiment status enum code.

**stop()** → None

Alias for ExperimentManager.finish()

**class** hvl\_ccb.experiment\_manager.ExperimentStatus

Bases: enum.Enum

Enumeration for the experiment status

**ERROR = 5**

**FINISHED = 4**

**FINISHING = 3**

**INITIALIZED = 0**

**INITIALIZING = -1**

**RUNNING = 2**

**STARTING = 1**

### 4.1.3 Module contents

Top-level package for HVL Common Code Base.



# CHAPTER 5

---

## Contributing

---

Contributions are welcome, and they are greatly appreciated! Every little bit helps, and credit will always be given.

You can contribute in many ways:

### 5.1 Types of Contributions

#### 5.1.1 Report Bugs

Report bugs at [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/issues](https://gitlab.com/ethz_hvl/hvl_ccb/issues).

If you are reporting a bug, please include:

- Your operating system name and version.
- Any details about your local setup that might be helpful in troubleshooting.
- Detailed steps to reproduce the bug.

#### 5.1.2 Fix Bugs

Look through the GitLab issues for bugs. Anything tagged with “bug” and “help wanted” is open to whoever wants to implement it.

#### 5.1.3 Implement Features

Look through the GitLab issues for features. Anything tagged with “enhancement” and “help wanted” is open to whoever wants to implement it.

### 5.1.4 Write Documentation

HVL Common Code Base could always use more documentation, whether as part of the official HVL Common Code Base docs, in docstrings, or even on the web in blog posts, articles, and such.

### 5.1.5 Submit Feedback

The best way to send feedback is to file an issue at [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/issues](https://gitlab.com/ethz_hvl/hvl_ccb/issues).

If you are proposing a feature:

- Explain in detail how it would work.
- Keep the scope as narrow as possible, to make it easier to implement.
- Remember that this is a volunteer-driven project, and that contributions are welcome :)

## 5.2 Get Started!

Ready to contribute? Here's how to set up *hvl\_ccb* for local development.

1. Clone *hvl\_ccb* repo from GitLab.

```
$ git clone git@gitlab.com:ethz_hvl/hvl_ccb.git
```

2. Install your local copy into a virtualenv. Assuming you have virtualenvwrapper installed, this is how you set up your fork for local development:

```
$ mkvirtualenv hvl_ccb
$ cd hvl_ccb/
$ python setup.py develop
$ pip install -r requirements_dev.txt
```

3. Create a branch for local development:

```
$ git checkout -b name-of-your-bugfix-or-feature
```

Now you can make your changes locally.

4. When you're done making changes, check that your changes pass flake8 and the tests, including testing other Python versions with tox:

```
$ flake8 hvl_ccb tests
$ python setup.py test or py.test
$ tox
```

To get flake8 and tox, just pip install them into your virtualenv. You can also use the provided make-like shell script to run flake8 and tests:

```
$ ./make.sh lint
$ ./make.sh test
```

5. Commit your changes and push your branch to GitLab:

```
$ git add .
$ git commit -m "Your detailed description of your changes."
$ git push origin name-of-your-bugfix-or-feature
```

6. Submit a merge request through the GitLab website.

## 5.3 Merge Request Guidelines

Before you submit a merge request, check that it meets these guidelines:

1. The merge request should include tests.
2. If the merge request adds functionality, the docs should be updated. Put your new functionality into a function with a docstring, and add the feature to the list in README.rst.
3. The merge request should work for Python 3.7. Check [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/merge\\_requests](https://gitlab.com/ethz_hvl/hvl_ccb/merge_requests) and make sure that the tests pass for all supported Python versions.

## 5.4 Tips

- To run tests from a single file:

```
$ py.test tests/test_hvl_ccb.py
```

or a single test function:

```
$ py.test tests/test_hvl_ccb.py::test_command_line_interface
```

- To add dependency, edit appropriate `*requirements` variable in the `setup.py` file and re-run:

```
$ python setup.py develop
```

- To generate a PDF version of the Sphinx documentation instead of HTML use:

```
$ rm -rf docs/hvl_ccb.rst docs/modules.rst docs/_build && sphinx-apidoc -o docs/ hvl_ccb && python -msphinx -M latexpdf docs/ docs/_build
```

This command can also be run through the make-like shell script:

```
$ ./make.sh docs-pdf
```

This requires a local installation of a LaTeX distribution, e.g. MikTeX.

## 5.5 Deploying

A reminder for the maintainers on how to deploy. Create `release-N.M.K` branch. Make sure all your changes are committed. Update or create entry in `HISTORY.rst` file, update features list in `README.rst` file and update API docs:

```
$ make docs
```

Commit all of the above and then run:

```
$ bumpversion patch # possible: major / minor / patch
$ git push
$ git push --tags
$ make release
```

Merge the release branch into master and devel branches with `--no-ff` flag and delete the release branch:

```
$ git checkout master
$ git merge --no-ff release-N.M.K
$ git checkout devel
$ git merge --no-ff release-N.M.K
$ git push --delete origin release-N.M.K
$ git branch --delete release-N.M.K
```

Finally, go to [https://gitlab.com/ethz\\_hvl/hvl\\_ccb/tags/](https://gitlab.com/ethz_hvl/hvl_ccb/tags/), select the latest release tag, press “Edit release notes” and add release notes (corresponding entry from `HISTORY.rst` file, but consider also additional brief header or synopsis if needed).

# CHAPTER 6

---

## Credits

---

### 6.1 Development Lead

- Mikołaj Rybiński <mikolaj.rybinski@id.ethz.ch>
- (previously) David Graber <graber@eeh.ee.ethz.ch>

### 6.2 Contributors

- Henrik Menne <henrik.menne@eeh.ee.ethz.ch>
- Alise Chachereau <chachereau@eeh.ee.ethz.ch>



## History

---

### 7.1 0.4.0 (2020-07-16)

- **Significantly improved new Supercube device controller:**
  - more robust error-handling,
  - status polling with generic Poller helper,
  - messages and status boards.
  - tested with a physical device,
- Improved OPC UA client wrapper, with better error handling, incl. re-tries on concurrent.futures.TimeoutError.
- SST Luminox Oxygen sensor device controller.
- **Backward-incompatible changes:**
  - CommunicationProtocol.access\_lock has changed type from threading.Lock to threading.RLock.
  - ILS2T.relative\_step and ILS2T.absolute\_position are now called, respectively, ILS2T.write\_relative\_step and ILS2T.write\_absolute\_position.
- **Minor bugfixes and improvements:**
  - fix use of max resolution in Labjack.set\_ain\_resolution(),
  - resolve ILS2T devices relative and absolute position setters race condition,
  - added acoustic horn function in the 2015 Supercube.
- **Toolchain changes:**
  - add Python 3.8 support,
  - drop pytest-runner support,
  - ensure compatibility with labjack\_ljm 2019 version library.

## 7.2 0.3.5 (2020-02-18)

- Fix issue with reading integers from LabJack LJM Library (device's product ID, serial number etc.)
- Fix development requirements specification (tox version).

## 7.3 0.3.4 (2019-12-20)

- **New devices using serial connection:**
  - Heinzinger Digital Interface I/II and a Heinzinger PNC power supply
  - Q-switched Pulsed Laser and a laser attenuator from CryLas
  - Newport SMC100PP single axis motion controller for 2-phase stepper motors
  - Pfeiffer TPG controller (TPG 25x, TPG 26x and TPG 36x) for Compact pressure Gauges
- PEP 561 compatibility and related corrections for static type checking (now in CI)
- **Refactorings:**
  - Protected non-thread safe read and write in communication protocols
  - Device sequence mixin: start/stop, add/rm and lookup
  - `.format()` to f-strings
  - more enumerations and a quite some improvements of existing code
- Improved error docstrings (`:raises:` annotations) and extended tests for errors.

## 7.4 0.3.3 (2019-05-08)

- Use PyPI labjack-ljm (no external dependencies)

## 7.5 0.3.2 (2019-05-08)

- INSTALLATION.rst with LJMPython prerequisite info

## 7.6 0.3.1 (2019-05-02)

- readthedocs.org support

## 7.7 0.3 (2019-05-02)

- Prevent an automatic close of VISA connection when not used.
- Rhode & Schwarz RTO 1024 oscilloscope using VISA interface over `TCP::INSTR`.
- Extended tests incl. messages sent to devices.

- Added Supercube device using an OPC UA client
- Added Supercube 2015 device using an OPC UA client (for interfacing with old system version)

## 7.8 0.2.1 (2019-04-01)

- Fix issue with LJMPython not being installed automatically with setuptools.

## 7.9 0.2.0 (2019-03-31)

- LabJack LJM Library communication wrapper and LabJack device.
- Modbus TCP communication protocol.
- Schneider Electric ILS2T stepper motor drive device.
- Elektro-Automatik PSI9000 current source device and VISA communication wrapper.
- Separate configuration classes for communication protocols and devices.
- Simple experiment manager class.

## 7.10 0.1.0 (2019-02-06)

- Communication protocol base and serial communication implementation.
- Device base and MBW973 implementation.



# CHAPTER 8

---

## Indices and tables

---

- genindex
- modindex
- search



---

## Python Module Index

---

### h

hvl\_ccb, 125  
hvl\_ccb.comm, 20  
hvl\_ccb.comm.base, 7  
hvl\_ccb.comm.labjack\_ljm, 7  
hvl\_ccb.comm.modbus\_tcp, 10  
hvl\_ccb.comm.opc, 11  
hvl\_ccb.comm.serial, 14  
hvl\_ccb.comm.visa, 17  
hvl\_ccb.configuration, 124  
hvl\_ccb.dev, 122  
hvl\_ccb.dev.base, 56  
hvl\_ccb.dev.crylas, 58  
hvl\_ccb.dev.ea\_psi9000, 67  
hvl\_ccb.dev.heinzinger, 72  
hvl\_ccb.dev.labjack, 77  
hvl\_ccb.dev.mbw973, 82  
hvl\_ccb.dev.newport, 85  
hvl\_ccb.dev.pfeiffer\_tpg, 99  
hvl\_ccb.dev.rs\_rto1024, 105  
hvl\_ccb.dev.se\_ilis2t, 111  
hvl\_ccb.dev.sst\_luminox, 117  
hvl\_ccb.dev.supercube, 42  
hvl\_ccb.dev.supercube.base, 20  
hvl\_ccb.dev.supercube.constants, 26  
hvl\_ccb.dev.supercube.typ\_a, 38  
hvl\_ccb.dev.supercube.typ\_b, 41  
hvl\_ccb.dev.supercube2015, 55  
hvl\_ccb.dev.supercube2015.base, 42  
hvl\_ccb.dev.supercube2015.constants, 47  
hvl\_ccb.dev.supercube2015.typ\_a, 53  
hvl\_ccb.dev.visa, 121  
hvl\_ccb.experiment\_manager, 124  
hvl\_ccb.utils, 123  
hvl\_ccb.utils.enum, 123  
hvl\_ccb.utils.typing, 123



---

## Index

---

### A

A (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitCurrent attribute*), 74  
A (*hvl\_ccb.dev.supercube.constants.SupercubeOpcEndpoint attribute*), 38  
A (*hvl\_ccb.dev.supercube2015.constants.SupercubeOpcEndpoint attribute*), 53  
ABSOLUTE\_POSITION (*hvl\_ccb.dev.se\_ilis2t.ILS2T.ActionsPtp attribute*), 111  
AC (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 85  
AC\_DoubleStage\_150kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 37  
AC\_DoubleStage\_150kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 52  
AC\_DoubleStage\_200kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 37  
AC\_DoubleStage\_200kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 52  
AC\_SingleStage\_100kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 37  
AC\_SingleStage\_100kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 52  
AC\_SingleStage\_50kV (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 37  
AC\_SingleStage\_50kV (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute*), 52  
acceleration (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
ACCESS\_ENABLE (*hvl\_ccb.dev.se\_ilis2t.ILS2TRegAddr attribute*), 116  
access\_lock (*hvl\_ccb.comm.base.CommunicationProtocol attribute*), 7  
ACTION\_JOG\_VALUE (*hvl\_ccb.dev.se\_ilis2t.ILS2T attribute*), 111  
activate\_output () (*hvl\_ccb.dev.sst\_luminox.Luminox method*), 117  
activated (*hvl\_ccb.dev.supercube.constants.BreakdownDetection attribute*), 32  
activated (*hvl\_ccb.dev.supercube2015.constants.BreakdownDetection attribute*), 48  
ACTIVE (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute*), 61  
active (*hvl\_ccb.dev.supercube.constants.OpcControl attribute*), 36  
add\_device () (*hvl\_ccb.dev.base.DeviceSequenceMixin method*), 56  
add\_device () (*hvl\_ccb.experiment\_manager.ExperimentManager method*), 124  
ADDR\_INCORRECT (*hvl\_ccb.dev.newport.NewportSMC100PPSerialComm attribute*), 95  
address (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 18  
address (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
address () (*hvl\_ccb.comm.visa.VisaCommunicationConfig.InterfaceType method*), 18  
Alarm0 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 47  
Alarm1 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 27  
Alarm1 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 26  
Alarm1 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 47  
Alarm10 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 27  
Alarm10 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 26







Alarm7 (*hvl\_ccb.dev.supercube2015.constants.AlarmText* Alarm90 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 48  
 Alarm70 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm71 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm72 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm73 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm74 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm75 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm76 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm77 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm78 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm79 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm8 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm8 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 27  
 Alarm8 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 48  
 Alarm80 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm81 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm82 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm83 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm84 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm85 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm86 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm87 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm88 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm89 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm9 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm9 (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 27  
 Alarm9 (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 48  
 Alarm90 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm91 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm92 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm93 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm94 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm95 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm96 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm97 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm98 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarm99 (*hvl\_ccb.dev.supercube.constants.Alarms attribute*), 31  
 Alarms (*class in hvl\_ccb.dev.supercube.constants*), 27  
 AlarmText (*class in hvl\_ccb.dev.supercube.constants*), 26  
 AlarmText (*class in hvl\_ccb.dev.supercube2015.constants*), 47  
 ANY (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.ConnectionType attribute*), 9  
 ANY (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.DeviceType attribute*), 9  
 ANY (*hvl\_ccb.dev.labjack.LabJack.DeviceType attribute*), 78  
 attenuation (*hvl\_ccb.dev.crylas.CryLasAttenuator attribute*), 58  
 AUTO (*hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes attribute*), 105  
 auto (*hvl\_ccb.dev.supercube.constants.EarthingStickOperatingStatus attribute*), 33  
 auto\_laser\_on (*hvl\_ccb.dev.crylas.CryLasLaserConfig attribute*), 64  
 AutoNumberNameEnum (*class in hvl\_ccb.utils.enum*), 123

**B**

B (*hvl\_ccb.dev.supercube.constants.SupercubeOpcEndpoint attribute*), 38  
 B (*hvl\_ccb.dev.supercube2015.constants.SupercubeOpcEndpoint attribute*), 53  
 BA (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 85  
 backlash\_compensation (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
 backup\_waveform () (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*),

```

    105
bar (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG.PressureUnits
     attribute), 100
barometric_pressure
    (hvl_ccb.dev.sst_luminox.LuminoxMeasurementType
     attribute), 119
base_velocity (hvl_ccb.dev.newport.NewportSMC100PPConfig
     attribute), 94
baudrate (hvl_ccb.comm.serial.SerialCommunicationConfig
     attribute), 15
baudrate (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig)
     (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication
      attribute), 60
baudrate (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig)
     generic_type() (in module
      hvl_ccb.utils.typing), 123
baudrate (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig)
     slave_config()
     (attribute), 76
baudrate (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig)
     attribute), 84
baudrate (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig)
     attribute), 98
baudrate (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig)
     attribute), 104
baudrate (hvl_ccb.dev.sst_luminox.LuminoxSerialCommunicationConfig)
     attribute), 120
BH (hvl_ccb.dev.newport.NewportConfigCommands
     attribute), 85
board (hvl_ccb.comm.visa.VisaCommunicationConfig
     attribute), 18
BreakdownDetection (class
     hvl_ccb.dev.supercube.constants), 32
BreakdownDetection (class
     hvl_ccb.dev.supercube2015.constants), 48
Bytesize (hvl_ccb.comm.serial.SerialCommunicationConfig
     attribute), 15
bytesize (hvl_ccb.comm.serial.SerialCommunicationConfig
     attribute), 16
bytesize (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig)
     attribute), 60
bytesize (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig)
     attribute), 66
bytesize (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig)
     attribute), 76
bytesize (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig)
     attribute), 76
bytesize (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig)
     attribute), 84
bytesize (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig)
     attribute), 84
bytesize (hvl_ccb.dev.newport.NewportSMC100PPConfig
     attribute), 98
bytesize (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig)
     attribute), 104
bytesize (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig)
     attribute), 104
C (hvl_ccb.dev.labjack.LabJack.ThermocoupleType
     attribute), 78
calibration_factor
    (hvl_ccb.dev.crylas.CryLasLaserConfig
     attribute), 64
cee16 (hvl_ccb.dev.supercube.constants.GeneralSockets
     attribute), 34
cee16 (hvl_ccb.dev.supercube2015.constants.GeneralSockets
     attribute), 94
check_for_error()
chunk_size (hvl_ccb.comm.visa.VisaCommunicationConfig
     attribute), 84
clean_values () (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig
     attribute), 98
clean_values () (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig
     attribute), 9
clean_values () (hvl_ccb.comm.opc.OpcUaCommunicationConfig
     attribute), 13
clean_values () (hvl_ccb.comm.serial.SerialCommunicationConfig
     attribute), 16
clean_values () (hvl_ccb.comm.visa.VisaCommunicationConfig
     attribute), 19
clean_values () (hvl_ccb.dev.base.EmptyConfig
     attribute), 57
clean_values () (hvl_ccb.dev.crylas.CryLasAttenuatorConfig
     attribute), 58
clean_values () (hvl_ccb.dev.crylas.CryLasLaserConfig
     attribute), 64
clean_values () (hvl_ccb.dev.ea_psi9000.PSI9000Config
     attribute), 64
clean_values () (hvl_ccb.dev.heinzinger.HeinzingerConfig
     attribute), 69
clean_values () (hvl_ccb.dev.heinzinger.HeinzingerConfig
     attribute), 72
clean_values () (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig)
     attribute), 66
clean_values () (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig)
     attribute), 76
clean_values () (hvl_ccb.dev.mbw973.MBW973Config
     attribute), 84
clean_values () (hvl_ccb.dev.newport.NewportSMC100PPConfig
     attribute), 84
clean_values () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig
     attribute), 98
clean_values () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig)
     attribute), 102
clean_values () (hvl_ccb.dev.se_ils2t.ILS2TConfig
     attribute), 104
clean_values () (hvl_ccb.dev.sst_luminox.LuminoxConfig
     attribute), 115
clean_values () (hvl_ccb.dev.supercube.base.SupercubeConfiguration
     attribute), 120
clean_values () (hvl_ccb.dev.supercube.base.SupercubeConfiguration
     attribute), 118
clean_values () (hvl_ccb.dev.supercube.base.SupercubeConfiguration
     attribute), 24
clean_values () (hvl_ccb.dev.supercube2015.base.SupercubeConfiguration
     attribute), 24

```

```

        method), 45
close(hvl_ccb.dev.supercube.constants.EarthingStickOperation CMD_NOT_ALLOWED_READY
      attribute), 33 (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication.C
close() (hvl_ccb.comm.base.CommunicationProtocol CMR (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG.SensorTypes
      method), 7 attribute), 100
close() (hvl_ccb.comm.labjack_ljm.LJMCommunication CODE_OR_ADDR_INVALID
      method), 8 (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication.C
close() (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication attribute), 95
      method), 10 com(hvl_ccb.dev.base.SingleCommDevice attribute), 57
close() (hvl_ccb.comm.opc.OpcUaCommunication COM_TIMEOUT (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication.C
      method), 12 attribute), 95
close() (hvl_ccb.comm.serial.SerialCommunication command(hvl_ccb.dev.sst_luminox.LuminoxMeasurementType
      method), 14 attribute), 119
close() (hvl_ccb.comm.visa.VisaCommunication CommunicationProtocol (class in
      method), 17 hvl_ccb.comm.base), 7
close_shutter() (hvl_ccb.dev.crylas.CryLasLaser config(hvl_ccb.configuration.ConfigurationMixin at-
      method), 62 tribute), 124
CLOSED(hvl_ccb.dev.crylas.CryLasLaser.AnswersShutter CONFIG(hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
      attribute), 61 attribute), 86
CLOSED(hvl_ccb.dev.crylas.CryLasLaserShutterStatus CONFIG(hvl_ccb.dev.newport.NewportStates attribute),
      attribute), 67 99
closed(hvl_ccb.dev.supercube.constants.DoorStatus config_cls() (hvl_ccb.comm.labjack_ljm.LJMCommunication
      attribute), 32 static method), 8
closed(hvl_ccb.dev.supercube.constants.EarthingStickStatus config_cls() (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication
      attribute), 34 static method), 10
closed(hvl_ccb.dev.supercube2015.constants.DoorStatus config_cls() (hvl_ccb.comm.opc.OpcUaCommunication
      attribute), 48 static method), 12
closed(hvl_ccb.dev.supercube2015.constants.EarthingStickStatus config_cls() (hvl_ccb.comm.serial.SerialCommunication
      attribute), 49 static method), 14
CMD_EXEC_ERROR(hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfigControllerErrorVisaCommunication
      attribute), 95 static method), 17
CMD_NOT_ALLOWED(hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfigControllerErrorConfigurationMixin
      attribute), 95 static method), 124
CMD_NOT_ALLOWED_CC config_cls() (hvl_ccb.dev.base.Device static
      attribute), 95 methodControllerErrors
CMD_NOT_ALLOWED_CONFIGURATION config_cls() (hvl_ccb.dev.crylas.CryLasAttenuator
      attribute), 95 static method), 58
CMD_NOT_ALLOWED_DISABLE config_cls() (hvl_ccb.dev.crylas.CryLasLaserSerialCommunication
      attribute), 95 static method), 59
CMD_NOT_ALLOWED_HOMING config_cls() (hvl_ccb.dev.crylas.CryLasLaser static
      attribute), 95 methodControllerErrors
CMD_NOT_ALLOWED_MOVING config_cls() (hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunication
      attribute), 95 static method), 65
CMD_NOT_ALLOWED_NOT_REFERENCED config_cls() (hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunication
      attribute), 95 static method), 68
CMD_NOT_ALLOWED_PP config_cls() (hvl_ccb.dev.heinzinger.HeinzingerDI
      attribute), 95 static method), 73
CMD_NOT_ALLOWED_NOT_REFERENCED config_cls() (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunication
      attribute), 95 static method), 75
CMD_NOT_ALLOWED_PP config_cls() (hvl_ccb.dev.mbw973.MBW973 static
      attribute), 95 methodControllerErrors
config_cls() (hvl_ccb.dev.mbw973.MBW973SerialCommunication
      attribute), 95

```

```

    static method), 84                               CryLasAttenuatorError, 59
config_cls () (hvl_ccb.dev.newport.NewportSMC100PPCryLasAttenuatorSerialCommunication
    static method), 87                               (class in hvl_ccb.dev.crylas), 59
config_cls () (hvl_ccb.dev.newport.NewportSMC100PRSerialCommunicationConfig
    static method), 96                               (class in hvl_ccb.dev.crylas), 59
config_cls () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG CryLasLaser (class in hvl_ccb.dev.crylas), 61
    static method), 100                             CryLasLaser.AnswersShutter (class in
config_cls () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig
    static method), 103                            CryLasLaser.AnswersStatus (class in
config_cls () (hvl_ccb.dev.rs_rto1024.RTO1024      hvl_ccb.dev.crylas), 61
    static method), 106                            CryLasLaser.LaserStatus (class in
config_cls () (hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunication
    static method), 110                            CryLasLaser.RepetitionRates (class in
config_cls () (hvl_ccb.dev.se_ils2t.ILS2T   static   hvl_ccb.dev.crylas), 62
    method), 112                                CryLasLaserConfig (class in hvl_ccb.dev.crylas), 64
config_cls () (hvl_ccb.dev.se_ils2t.ILS2TModbusTcpCommunication
    static method), 115                           CryLasLaserNotReadyError, 65
config_cls () (hvl_ccb.dev.sst_luminox.Luminox CryLasLaserSerialCommunication (class in
    static method), 117                           hvl_ccb.dev.crylas), 65
config_cls () (hvl_ccb.dev.sst_luminox.LuminoxSerialCommunicationConfig
    static method), 119                           (class in hvl_ccb.dev.crylas), 66
config_cls () (hvl_ccb.dev.supercube.base.SupercubeBaseCryLasLaserShutterStatus (class in
    static method), 20                            hvl_ccb.dev.crylas), 67
config_cls () (hvl_ccb.dev.supercube.base.SupercubeOpcaUaCommunicationLimit
    static method), 25                           current_primary (hvl_ccb.dev.supercube.constants.Power
config_cls () (hvl_ccb.dev.supercube.typ_a.SupercubeAOpcaUaAttributed)
    static method), 38                           current_primary (hvl_ccb.dev.supercube.constants.Power
config_cls () (hvl_ccb.dev.supercube.typ_b.SupercubeBOpcaUaAttributed)
    static method), 41                           current_primary (hvl_ccb.dev.supercube2015.constants.Power
config_cls () (hvl_ccb.dev.supercube2015.base.Supercube2015BaseAttribute), 51
    static method), 42                           current_upper_limit
config_cls () (hvl_ccb.dev.supercube2015.base.SupercubeOpcaUaAttributed
    static method), 46                           current_upper_limit
config_cls () (hvl_ccb.dev.supercube2015.typ_a.SupercubeAOpcaUaCommunication
    static method), 54                           attribute), 70
config_cls () (hvl_ccb.dev.visa.VisaDevice  static
    method), 121                                attribute), 93
configdataclass () (in module datachange_notification()
    hvl_ccb.configuration), 124                  (hvl_ccb.comm.opc.OpcUaSubHandler
configuration_save_json ()                         method), 14
    ConfigurationMixin (class in
    hvl_ccb.configuration), 124                  datachange_notification()
                                                (hvl_ccb.dev.supercube.base.SupercubeSubscriptionHandler
                                                method), 26
                                                datachange_notification()
connection_type (hvl_ccb.comm.labjack_ljm.LJMCommunication
    attribute), 9                                (hvl_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler
                                                method), 47
contact_range (hvl_ccb.dev.supercube.constants.GeneralSupport
    attribute), 34                                datachange_manufacture
                                                (hvl_ccb.dev.sst_luminox.LuminoxMeasurementType
create_serial_port ()                           attribute), 119
    (hvl_ccb.comm.serial.SerialCommunicationConfig
    method), 16                                DC_DoubleStage_280kV
                                                (hvl_ccb.dev.supercube.constants.PowerSetup
CryLasAttenuator (class in hvl_ccb.dev.crylas), 58
    attribute), 37
CryLasAttenuatorConfig (class in DC_DoubleStage_280kV
    hvl_ccb.dev.crylas), 58
                                                (hvl_ccb.dev.supercube2015.constants.PowerSetup

```

## D

```

datachange_notification()
(hvl_ccb.comm.opc.OpcUaSubHandler
method), 14
datachange_notification()
(hvl_ccb.dev.supercube.base.SupercubeSubscriptionHandler
method), 26
datachange_notification()
(hvl_ccb.dev.supercube2015.base.SupercubeSubscriptionHandler
method), 47
datachange_manufacture
(hvl_ccb.dev.sst_luminox.LuminoxMeasurementType
attribute), 119
DC_DoubleStage_280kV
(hvl_ccb.dev.supercube.constants.PowerSetup
attribute), 37
DC_DoubleStage_280kV
(hvl_ccb.dev.supercube2015.constants.PowerSetup

```

*attribute), 52*  
**DC\_SingleStage\_140kV** (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute), 37*  
**DC\_SingleStage\_140kV** (*hvl\_ccb.dev.supercube2015.constants.PowerSetup attribute), 52*  
**DC\_VOLTAGE\_TOO\_LOW** (*hvl\_ccb.dev.newport.NewportSMC100PP.MotorErrors attribute), 86*  
**default\_com\_cls ()** (*hvl\_ccb.dev.base.SingleCommDevice method), 57*  
**default\_com\_cls ()** (*hvl\_ccb.dev.crylas.CryLasAttenuator method), 58*  
**default\_com\_cls ()** (*hvl\_ccb.dev.crylas.CryLasLaser method), 62*  
**default\_com\_cls ()** (*hvl\_ccb.dev.ea\_psi9000.PSI9000 method), 68*  
**default\_com\_cls ()** (*hvl\_ccb.dev.heinzinger.HeinzingerDI method), 73*  
**default\_com\_cls ()** (*hvl\_ccb.dev.labjack.LabJack static method), 79*  
**default\_com\_cls ()** (*hvl\_ccb.dev.mbw973.MBW973 static method), 82*  
**default\_com\_cls ()** (*hvl\_ccb.dev.newport.NewportSMC100PP static method), 87*  
**default\_com\_cls ()** (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG method), 101*  
**default\_com\_cls ()** (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 106*  
**default\_com\_cls ()** (*hvl\_ccb.dev.se\_ilst2t.ILS2T static method), 112*  
**default\_com\_cls ()** (*hvl\_ccb.dev.sst\_luminox.Luminox method), 117*  
**default\_com\_cls ()** (*hvl\_ccb.dev.supercube.base.SupercubeBase static method), 20*  
**default\_com\_cls ()** (*hvl\_ccb.dev.supercube.typ\_a.SupercubeWithFU static method), 39*  
**default\_com\_cls ()** (*hvl\_ccb.dev.supercube.typ\_b.SupercubeB static method), 41*  
**default\_com\_cls ()** (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base static method), 42*  
**default\_com\_cls ()** (*hvl\_ccb.dev.supercube2015.typ\_a.Supercube2015WithFU static method), 54*  
**default\_com\_cls ()** (*hvl\_ccb.dev.visa.VisaDevice static method), 121*  
**DEFAULT\_IO\_SCANNING\_CONTROL\_VALUES**  
**Device** (*class in hvl\_ccb.dev.base), 56*  
**device\_type** (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig attribute), 9*  
**DeviceExistingException**, 56  
**DeviceSequenceMixin** (*class in hvl\_ccb.dev.base), 56*  
**DIOChannel** (*hvl\_ccb.dev.labjack.LabJack attribute), 78*  
**DISABLE** (*hvl\_ccb.dev.newport.NewportStates attribute), 99*  
**disable ()** (*hvl\_ccb.dev.se\_ilst2t.ILS2T method), 112*  
**DISABLE\_FROM\_JOGGING** (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 87*  
**DISABLE\_FROM\_MOVING** (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 87*  
**DISABLE\_FROM\_READY** (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute), 87*  
**DisableEspStageCheck** (*hvl\_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig attribute), 93*  
**DISPLACEMENT\_OUT\_OF\_LIMIT** (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication.C class in hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication), 95*  
**display\_message\_board ()** (*hvl\_ccb.dev.supercube.base.SupercubeBase method), 21*  
**display\_status\_board ()** (*hvl\_ccb.dev.supercube.base.SupercubeBase method), 21*  
**do\_ioscanning\_write ()** (*hvl\_ccb.dev.se\_ilst2t.ILS2T method), 112*  
**Door** (*class in hvl\_ccb.dev.supercube.constants), 32*  
**DoorStatus** (*class in hvl\_ccb.dev.supercube.constants), 32*  
**DoorStatus** (*class in hvl\_ccb.dev.supercube2015.constants), 48*  
**E**  
**E** (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute), 121*

<p>tribute), 79</p> <p>EarthingStick (class     <i>hvl_ccb.dev.supercube.constants</i>), 32</p> <p>EarthingStick (class     <i>hvl_ccb.dev.supercube2015.constants</i>), 49</p> <p>EarthingStickMeta (class     <i>hvl_ccb.dev.supercube.constants</i>), 33</p> <p>EarthingStickOperatingStatus (class     <i>hvl_ccb.dev.supercube.constants</i>), 33</p> <p>EarthingStickOperation (class     <i>hvl_ccb.dev.supercube.constants</i>), 33</p> <p>EarthingStickStatus (class     <i>hvl_ccb.dev.supercube.constants</i>), 33</p> <p>EarthingStickStatus (class     <i>hvl_ccb.dev.supercube2015.constants</i>), 49</p> <p>EPPROM_ACCESS_ERROR     (<i>hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication</i>.<i>hvl_CcbControllerError</i>.<i>ILS2T</i> method), 113</p> <p>EIGHT (<i>hvl_ccb.dev.heinzinger.HeinzingerConfig</i>.<i>RecordingsEnum</i> (<i>hvl_ccb.dev.se_ilst.ILS2T</i> method), 113     attribute), 72</p> <p>EIGHTBITS (<i>hvl_ccb.comm.serial.SerialCommunication</i>.<i>Bytesize</i> (<i>hvl_ccb.dev.newport.NewportSMC100PP</i>     attribute), 87</p> <p>EmptyConfig (class in <i>hvl_ccb.dev.base</i>), 57</p> <p>enable() (<i>hvl_ccb.dev.se_ilst.ILS2T</i> method), 112</p> <p>EnableEspStageCheck     (<i>hvl_ccb.dev.newport.NewportSMC100PPConfig</i>.<i>ExpStageConfig</i>.<i>Error</i>, 124     attribute), 93</p> <p>ENCODING (<i>hvl_ccb.comm.serial.SerialCommunication</i>     attribute), 14</p> <p>EndOfRunSwitch (<i>hvl_ccb.dev.newport.NewportSMC100PPConfig</i>.<i>HomeSearch</i> (<i>hvl_ccb.dev.supercube.constants.PowerSetup</i>     attribute), 93</p> <p>EndOfRunSwitch_and_Index     (<i>hvl_ccb.dev.newport.NewportSMC100PPConfig</i>.<i>HomeSearch</i> (<i>hvl_ccb.dev.supercube2015.constants.PowerSetup</i>     attribute), 93</p> <p>endpoint_name (<i>hvl_ccb.comm.opc.OpcUaCommunicationConfig</i>     attribute), 13</p> <p>endpoint_name (<i>hvl_ccb.dev.supercube.typ_a.SupercubeAOpquaConfiguration</i>.<i>LabJack.TemperatureUnit</i> at-     tribute), 39</p> <p>endpoint_name (<i>hvl_ccb.dev.supercube.typ_b.SupercubeBOpquaConfiguration</i>.<i>ils2t.ILS2T.Ref16Jog</i> attribute),     41</p> <p>endpoint_name (<i>hvl_ccb.dev.supercube2015.typ_a.SupercubeAOpquaConfiguration</i>.<i>hvl_ccb.dev.rs_rto1024.RTO1024</i>     attribute), 55</p> <p>error (<i>hvl_ccb.dev.supercube.constants.DoorStatus</i> at-     tribute), 32</p> <p>error (<i>hvl_ccb.dev.supercube.constants.EarthingStickStatus</i> attribute), 34</p> <p>Error (<i>hvl_ccb.dev.supercube.constants.SafetyStatus</i> at-     tribute), 38</p> <p>error (<i>hvl_ccb.dev.supercube2015.constants.DoorStatus</i> attribute), 48</p> <p>error (<i>hvl_ccb.dev.supercube2015.constants.EarthingStickStatus</i> attribute), 50</p> <p>Error (<i>hvl_ccb.dev.supercube2015.constants.SafetyStatus</i> attribute), 53</p>	<p>in     ERROR (<i>hvl_ccb.experiment_manager.ExperimentStatus</i>             attribute), 125</p> <p>in     Errors (class in <i>hvl_ccb.dev.supercube.constants</i>), 34</p> <p>in     Errors (class in <i>hvl_ccb.dev.supercube2015.constants</i>),             50</p> <p>in     ESP_STAGE_NAME_INVALID             (<i>hvl_ccb.dev.newport.NewportSMC100PPSerialCommunication</i>.<i>hvl_CcbControllerError</i>.<i>ILS2T</i> method), 113</p> <p>in     ETHERNET (<i>hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig</i>.<i>Connec-</i></p> <p>in     attribute), 9</p> <p>in     EVEN (<i>hvl_ccb.comm.serial.SerialCommunication</i>.<i>Parity</i>             attribute), 17</p> <p>in     event_notification ()             (<i>hvl_ccb.comm.opc.OpcUaSubHandler</i>             method), 14</p> <p>execute_absolute_position ()</p> <p>execute_relative_step ()</p> <p>exit_configuration ()</p> <p>exit_configuration_wait_sec     (<i>hvl_ccb.dev.newport.NewportSMC100PPConfig</i>     attribute), 94</p> <p>ExperimentManager (class     <i>hvl_ccb.experiment_manager</i>), 124</p> <p>ExperimentStatus (class     <i>hvl_ccb.experiment_manager</i>), 124</p> <p>External (<i>hvl_ccb.dev.supercube.constants.PowerSetup</i>     attribute), 37</p> <p>F</p> <p>FINISHED (<i>hvl_ccb.experiment_manager.ExperimentStatus</i>     attribute), 125</p> <p>FINISHING (<i>hvl_ccb.experiment_manager.ExperimentStatus</i>     attribute), 125</p> <p>FIVEBITS (<i>hvl_ccb.comm.serial.SerialCommunication</i>.<i>Bytesize</i>     attribute), 15</p> <p>FLT_STATUSINFO (<i>hvl_ccb.dev.se_ilst.ILS2TRegAddr</i> at-     tribute), 116</p> <p>FLT_MEM_DEL (<i>hvl_ccb.dev.se_ilst.ILS2TRegAddr</i> at-     tribute), 116</p>
--	--

FLT\_MEM\_RESET (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr* force\_value () (*hvl\_ccb.dev.sst\_luminox.LuminoxSerialCommunication attribute*), 116  
 FOLLOWING\_ERROR (*hvl\_ccb.dev.newport.NewportSMC100PPMotorErrors*) (*hvl\_ccb.dev.supercube.base.SupercubeConfiguration attribute*), 86  
*force\_value ()* (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig* ()) (*hvl\_ccb.dev.supercube.base.SupercubeOpcUaCommunication method*), 25  
*force\_value ()* (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig* (*hvl\_ccb.dev.supercube.typ\_a.SupercubeAOpcUaConfig attribute*)), 11  
*force\_value ()* (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig* value () (*hvl\_ccb.dev.supercube.typ\_b.SupercubeBOpcUaConfig method*)), 41  
*force\_value ()* (*hvl\_ccb.comm.serial.SerialCommunicationConfig* value () (*hvl\_ccb.dev.supercube2015.base.SupercubeConfigurations method*)), 16  
*force\_value ()* (*hvl\_ccb.comm.visa.VisaCommunicationConfig* value () (*hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaConfig method*)), 19  
*force\_value ()* (*hvl\_ccb.dev.base.EmptyConfig* force\_value () (*hvl\_ccb.dev.supercube2015.typ\_a.SupercubeAOpcUaConfig method*)), 55  
*force\_value ()* (*hvl\_ccb.dev.crylas.CryLasAttenuatorConfig* value () (*hvl\_ccb.dev.visa.VisaDeviceConfig method*)), 59  
*force\_value ()* (*hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig* value () (*hvl\_ccb.dev.visa.VisaDeviceConfig attribute*)), 72  
*force\_value ()* (*hvl\_ccb.dev.crylas.CryLasLaserConfig* FREERUN (*hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes method*)), 64  
*force\_value ()* (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig* (*hvl\_ccb.dev.supercube.constants.Power attribute*)), 36  
*force\_value ()* (*hvl\_ccb.dev.ea\_psi9000.PSI9000Config* frequency (*hvl\_ccb.dev.supercube2015.constants.Power method*)), 70  
*force\_value ()* (*hvl\_ccb.dev.ea\_psi9000.PSI9000VisaCommunicationConfig* (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*)), 71  
*force\_value ()* (*hvl\_ccb.dev.heinzinger.HeinzingerConfig* from\_json () (*hvl\_ccb.configuration.ConfigurationMixin method*)), 124  
*force\_value ()* (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig* (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*)), 76  
*force\_value ()* (*hvl\_ccb.dev.mbw973.MBW973Config* fso\_reset () (*hvl\_ccb.dev.supercube.typ\_a.SupercubeWithFU method*)), 83  
*force\_value ()* (*hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig* (*hvl\_ccb.dev.supercube2015.typ\_a.Supercube2015WithFU method*)), 84  
*force\_value ()* (*hvl\_ccb.dev.newport.NewportSMC100PPConfig* G)  
*force\_value ()* (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig* (class in *hvl\_ccb.dev.supercube.constants*), 34  
*force\_value ()* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig* GeneralSockets (class in *hvl\_ccb.dev.supercube2015.constants*), 50  
*force\_value ()* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig* GeneralSupport (class in *hvl\_ccb.dev.supercube.constants*), 34  
*force\_value ()* (*hvl\_ccb.dev.rs\_rto1024.RTO1024Config* GeneralSupport (class in *hvl\_ccb.dev.supercube2015.constants*), 50  
*force\_value ()* (*hvl\_ccb.dev.rs\_rto1024.RTO1024VisaCommunicationConfig* Meta (class in *hvl\_ccb.dev.supercube.constants*), 35  
*force\_value ()* (*hvl\_ccb.dev.se\_ils2t.ILS2TConfig* get (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 27  
*force\_value ()* (*hvl\_ccb.dev.se\_ils2t.ILS2TModbusTcpCommunicationConfig* get (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*)), 48  
*force\_value ()* (*hvl\_ccb.dev.sst\_luminox.LuminoxConfig* get (*hvl\_ccb.dev.supercube2015.constants.MeasurementsDividerRatio method*)), 51

```

get (hvl_ccb.dev.supercube2015.constants.MeasurementsS  
ettedInputOr_queue () (hvl_ccb.dev.visa.VisaDevice  
attribute), 51
get_acceleration ()  
    (hvl_ccb.dev.newport.NewportSMC100PP  
method), 87
get_ain () (hvl_ccb.dev.labjack.LabJack method), 79
get_by_p_id (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfigDeviceType(hvl_ccb.dev.supercube.typ_a.SupercubeWithFU  
attribute), 9
get_by_p_id (hvl_ccb.dev.labjack.LabJack.DeviceType  
attribute), 78
get_cal_current_source ()  
    (hvl_ccb.dev.labjack.LabJack method), 79
get_cee16_socket ()  
    (hvl_ccb.dev.supercube.base.SupercubeBase  
method), 21
get_cee16_socket ()  
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base  
method), 42
get_controller_information ()  
    (hvl_ccb.dev.newport.NewportSMC100PP  
method), 88
get_current () (hvl_ccb.dev.heinzinger.HeinzingerDI  
method), 73
get_dc_volt () (hvl_ccb.dev.se_ilst2t.ILS2T method),  
    113
get_device () (hvl_ccb.dev.base.DeviceSequenceMixin  
method), 56
get_devices () (hvl_ccb.dev.base.DeviceSequenceMixin  
method), 56
get_digital_input ()  
    (hvl_ccb.dev.labjack.LabJack method), 79
get_door_status ()  
    (hvl_ccb.dev.supercube.base.SupercubeBase  
method), 21
get_door_status ()  
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base  
method), 42
get_earthing_manual ()  
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base  
method), 42
get_earthing_status ()  
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base  
method), 43
get_earthing_stick_manual ()  
    (hvl_ccb.dev.supercube.base.SupercubeBase  
method), 21
get_earthing_stick_operating_status ()  
    (hvl_ccb.dev.supercube.base.SupercubeBase  
method), 21
get_earthing_stick_status ()  
    (hvl_ccb.dev.supercube.base.SupercubeBase  
method), 21
get_error_code ()      (hvl_ccb.dev.se_ilst2t.ILS2T  
method), 113
get_fso_active () (hvl_ccb.dev.supercube2015.typ_a.Supercube2015W  
method), 54
get_full_scale_mbar ()  
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method),  
    101
get_full_scale_unitless ()  
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method),  
    101
get_base_identification ()  
    (hvl_ccb.dev.visa.VisaDevice method), 121
get_interface_version ()  
    (hvl_ccb.dev.heinzinger.HeinzingerDI method),  
    73
get_max_voltage ()  
    (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU  
method), 40
get_max_voltage ()  
    (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU  
method), 54
get_measurement_ratio ()  
    (hvl_ccb.dev.supercube.base.SupercubeBase  
method), 21
get_measurement_ratio ()  
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base  
method), 43
get_measurement_voltage ()  
    (hvl_ccb.dev.supercube.base.SupercubeBase  
method), 21
get_measurement_voltage ()  
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base  
method), 43
get_motor_configuration ()  
    (hvl_ccb.dev.newport.NewportSMC100PP  
method), 88
get_move_duration ()  
    (hvl_ccb.dev.newport.NewportSMC100PP  
method), 88
get_negative_software_limit ()  
    (hvl_ccb.dev.newport.NewportSMC100PP  
method), 88
get_number_of_recordings ()  
    (hvl_ccb.dev.heinzinger.HeinzingerDI method),  
    73
get_output ()      (hvl_ccb.dev.ea_psi9000.PSI9000  
method), 68
get_position () (hvl_ccb.dev.newport.NewportSMC100PP

```

```

        method), 89
get_position()      (hvl_ccb.dev.se_ilst.ILS2T
        method), 113
get_positive_software_limit()
    (hvl_ccb.dev.newport.NewportSMC100PP
        method), 89
get_power_setup()
    (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
        method), 40
get_power_setup()
    (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
        method), 54
get_primary_current()
    (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
        method), 40
get_primary_current()
    (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
        method), 54
get_primary_voltage()
    (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
        method), 40
get_primary_voltage()
    (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
        method), 54
get_product_id()
    (hvl_ccb.dev.labjack.LabJack
        method), 79
get_product_name()
    (hvl_ccb.dev.labjack.LabJack
        method), 79
get_product_type()
    (hvl_ccb.dev.labjack.LabJack
        method), 79
get_pulse_energy_and_rate()
    (hvl_ccb.dev.crylas.CryLasLaser
        method), 62
get_sbus_rh()
    (hvl_ccb.dev.labjack.LabJack
        method), 80
get_sbus_temp()
    (hvl_ccb.dev.labjack.LabJack
        method), 80
get_serial_number()
    (hvl_ccb.dev.heinzinger.HeinzingerDI
        method), 73
get_serial_number()
    (hvl_ccb.dev.labjack.LabJack
        method), 80
get_state()
    (hvl_ccb.dev.newport.NewportSMC100PP
        method), 89
get_status()
    (hvl_ccb.dev.se_ilst.ILS2T
        method), 113
get_status()
    (hvl_ccb.dev.supercube.base.SupercubeBase
        method), 21
get_status()
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base
        method), 43
get_support_input()
    (hvl_ccb.dev.supercube.base.SupercubeBase
        method), 22
get_support_input()

(hvl_ccb.dev.supercube2015.base.Supercube2015Base
        method), 43
get_support_output()
    (hvl_ccb.dev.supercube.base.SupercubeBase
        method), 22
get_support_output()
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base
        method), 43
get_system_lock()
    (hvl_ccb.dev.ea_psi9000.PSI9000
        method),
get_t13_socket()
    (hvl_ccb.dev.supercube.base.SupercubeBase
        method), 22
get_t13_socket()
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base
        method), 43
get_target_voltage()
    (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
        method), 40
get_target_voltage()
    (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
        method), 54
get_temperature()
    (hvl_ccb.dev.se_ilst.ILS2T
        method), 113
get_timestamps()
    (hvl_ccb.dev.rs_rto1024.RTO1024
        method), 106
get_ui_lower_limits()
    (hvl_ccb.dev.ea_psi9000.PSI9000
        method),
    68
get_uip_upper_limits()
    (hvl_ccb.dev.ea_psi9000.PSI9000
        method),
    68
get_voltage()
    (hvl_ccb.dev.heinzinger.HeinzingerDI
        method), 73
get_voltage_current_setpoint()
    (hvl_ccb.dev.ea_psi9000.PSI9000
        method),
    68
go_home()
    (hvl_ccb.dev.newport.NewportSMC100PP
        method), 89
go_to_configuration()
    (hvl_ccb.dev.newport.NewportSMC100PP
        method), 89
GreenNotReady(hvl_ccb.dev.supercube.constants.SafetyStatus
        attribute), 38
GreenNotReady(hvl_ccb.dev.supercube2015.constants.SafetyStatus
        attribute), 53
GreenReady(hvl_ccb.dev.supercube.constants.SafetyStatus
        attribute), 38
GreenReady(hvl_ccb.dev.supercube2015.constants.SafetyStatus
        attribute), 53
H
HARDWARE(hvl_ccb.dev.crylas.CryLasLaser.RepetitionRates
        attribute), 62

```

HEAD (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute*), 61  
 HeinzingerConfig (*class in hvl\_ccb.dev.heinzinger*), 72  
 HeinzingerConfig.RecordingsEnum (*class in hvl\_ccb.dev.heinzinger*), 72  
 HeinzingerDI (*class in hvl\_ccb.dev.heinzinger*), 73  
 HeinzingerPNC (*class in hvl\_ccb.dev.heinzinger*), 74  
 HeinzingerPNC.UnitCurrent (*class in hvl\_ccb.dev.heinzinger*), 74  
 HeinzingerPNC.UnitVoltage (*class in hvl\_ccb.dev.heinzinger*), 74  
 HeinzingerPNCDeviceNotRecognizedException (*hvl\_ccb (module)*), 125  
 HeinzingerPNCError, 75  
 HeinzingerPNCMaxCurrentExceededException (*hvl\_ccb.comm.labjack\_ljm (module)*), 75  
 HeinzingerPNCMaxVoltageExceededException (*hvl\_ccb.comm.opc (module)*), 11  
 HeinzingerSerialCommunication (*class in hvl\_ccb.dev.heinzinger*), 75  
 HeinzingerSerialCommunicationConfig (*class in hvl\_ccb.dev.heinzinger*), 76  
 HIGH (*hvl\_ccb.dev.labjack.LabJack.DIOStatus attribute*), 78  
 home\_search\_polling\_interval (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
 home\_search\_timeout (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
 home\_search\_type (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
 home\_search\_velocity (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
 HOME\_STARTED (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication.ControllerErrors attribute*), 95  
 HomeSwitch (*hvl\_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch attribute*), 93  
 HomeSwitch\_and\_Index (*hvl\_ccb.dev.newport.NewportSMC100PPConfig.HomeSearch attribute*), 94  
 HOMING (*hvl\_ccb.dev.newport.NewportStates attribute*), 99  
 HOMING\_FROM\_RS232 (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 87  
 HOMING\_FROM\_SMC (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 87  
 HOMING\_TIMEOUT (*hvl\_ccb.dev.newport.NewportSMC100PPMotionEvents attribute*), 86  
 horn (*hvl\_ccb.dev.supercube2015.constants.Safety attribute*), 52  
 horn () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 43  
 host (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig attribute*), 11  
 host (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig attribute*), 13  
 host (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 19  
 hPascal (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits attribute*), 100  
 HT (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 86  
 hvl\_ccb.comm (i.e., 20  
 hvl\_ccb.comm.base (i.e., 7  
 hvl\_ccb.comm.labjack\_ljm (i.e., 7  
 hvl\_ccb.comm.modbus\_tcp (i.e., 10  
 hvl\_ccb.comm.opc (i.e., 11  
 hvl\_ccb.comm.serial (i.e., 14  
 hvl\_ccb.comm.visa (i.e., 17  
 hvl\_ccb.configuration (i.e., 124  
 hvl\_ccb.dev (i.e., 122  
 hvl\_ccb.dev.base (i.e., 56  
 hvl\_ccb.dev.crylas (i.e., 58  
 hvl\_ccb.dev.ea\_psi9000 (i.e., 67  
 hvl\_ccb.dev.heinzinger (i.e., 72  
 hvl\_ccb.dev.labjack (i.e., 77  
 hvl\_ccb.dev.mbw973 (i.e., 82  
 hvl\_ccb.dev.newport (i.e., 85  
 hvl\_ccb.dev.pfeiffer\_tpg (i.e., 99  
 hvl\_ccb.dev.rs\_rto1024 (i.e., 105  
 hvl\_ccb.dev.se\_ils2t (i.e., 111  
 hvl\_ccb.dev.sst\_luminox (i.e., 117  
 hvl\_ccb.dev.supercube (i.e., 42  
 hvl\_ccb.dev.supercube.base (i.e., 20  
 hvl\_ccb.dev.supercube.constants (i.e.,  
 HOME\_STARTED (i.e., 38  
 HomeSearch supercube.typ\_b (i.e., 41  
 hvl\_ccb.dev.supercube2015 (i.e., 55  
 hvl\_ccb.dev.supercube2015.base (i.e.,  
 hvl\_ccb.dev.supercube2015.constants (i.e., 47  
 hvl\_ccb.dev.supercube2015.typ\_a (i.e., 53  
 hvl\_ccb.dev.supercube2015.typ\_a (i.e., 121  
 hvl\_ccb.experiment\_manager (i.e., 124  
 hvl\_ccb.utils.enum (i.e., 123  
 hysteresis\_compensation (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94

| Identification\_error  
     (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus  
         attribute), 100  
 identifier (hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig  
         attribute), 9  
 identify\_device ()  
     (hvl\_ccb.dev.heinzinger.HeinzingerPNC  
         method), 75  
 identify\_sensors ()  
     (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG method),  
         101  
 IKR (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes  
         attribute), 100  
 IKR11 (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes  
         attribute), 100  
 IKR9 (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes  
         attribute), 100  
 ILS2T (class in hvl\_ccb.dev.se\_ilst), 111  
 ILS2T.ActionsPtp (class in hvl\_ccb.dev.se\_ilst),  
     111  
 ILS2T.Mode (class in hvl\_ccb.dev.se\_ilst), 111  
 ILS2T.Ref16Jog (class in hvl\_ccb.dev.se\_ilst), 111  
 ILS2T.State (class in hvl\_ccb.dev.se\_ilst), 112  
 ILS2TConfig (class in hvl\_ccb.dev.se\_ilst), 114  
 ILS2TException, 115  
 ILS2TModbusTcpCommunication (class in  
     hvl\_ccb.dev.se\_ilst), 115  
 ILS2TModbusTcpCommunicationConfig (class  
     in hvl\_ccb.dev.se\_ilst), 115  
 ILS2TRegAddr (class in hvl\_ccb.dev.se\_ilst), 116  
 ILS2TRegDatatype (class in hvl\_ccb.dev.se\_ilst),  
     117  
 IMR (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes  
         attribute), 100  
 in\_1\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 34  
 in\_1\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_1\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 34  
 in\_1\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_2\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 34  
 in\_2\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_2\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 34  
 in\_2\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_3\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_3\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_3\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_3\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_4\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_4\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_4\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_4\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_5\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_5\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_5\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_5\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_6\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_6\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 50  
 in\_6\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport  
         attribute), 35  
 in\_6\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport  
         attribute), 51  
 INACTIVE (hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus  
         attribute), 61  
 inactive (hvl\_ccb.dev.supercube.constants.DoorStatus  
         attribute), 32  
 inactive (hvl\_ccb.dev.supercube.constants.EarthingStickStatus  
         attribute), 34  
 inactive (hvl\_ccb.dev.supercube2015.constants.DoorStatus  
         attribute), 48  
 inactive (hvl\_ccb.dev.supercube2015.constants.EarthingStickStatus  
         attribute), 50  
 init\_attenuation (hvl\_ccb.dev.crylas.CryLasAttenuatorConfig  
         attribute), 59  
 init\_monitored\_nodes ()  
     (hvl\_ccb.comm.opc.OpcUaCommunication  
         method), 12  
 init\_shutter\_status  
     (hvl\_ccb.dev.crylas.CryLasLaserConfig  
         attribute), 65  
 initialize () (hvl\_ccb.dev.newport.NewportSMC100PP  
         method), 90  
 INITIALIZED (hvl\_ccb.experiment\_manager.ExperimentStatus  
         attribute), 125  
 Initializing (hvl\_ccb.dev.supercube.constants.SafetyStatus  
         attribute), 38  
 Initializing (hvl\_ccb.dev.supercube2015.constants.SafetyStatus  
         attribute), 38

```

        attribute), 53                                is_configdataclass
INITIALIZING(hvl_ccb.experiment_manager.ExperimentStatus (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunicationConfig
        attribute), 125                                attribute), 11
input(hvl_ccb.dev.supercube.constants.GeneralSupport is_configdataclass
        attribute), 35                                (hvl_ccb.comm.opc.OpcUaCommunicationConfig
input(hvl_ccb.dev.supercube2015.constants.GeneralSupport attribute), 13
        attribute), 51                                is_configdataclass
input_1(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio (hvl_ccb.comm.serial.SerialCommunicationConfig
        attribute), 35                                attribute), 16
input_1(hvl_ccb.dev.supercube.constants.MeasurementsScaledInput is_configdataclass
        attribute), 36                                (hvl_ccb.comm.visa.VisaCommunicationConfig
input_1(hvl_ccb.dev.supercube2015.constants.MeasurementsDividerRatio), 19
        attribute), 51                                is_configdataclass
input_1(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInputConfig attribute),
        attribute), 51                                57
input_2(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio is_configdataclass
        attribute), 35                                (hvl_ccb.dev.crylas.CryLasAttenuatorConfig
input_2(hvl_ccb.dev.supercube.constants.MeasurementsScaledInputAttribute), 59
        attribute), 36                                is_configdataclass
input_2(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInputConfig attribute),
        attribute), 51                                65
input_3(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio is_configdataclass
        attribute), 35                                (hvl_ccb.dev.heinzinger.HeinzingerConfig
input_3(hvl_ccb.dev.supercube.constants.MeasurementsScaledInputAttribute), 72
        attribute), 36                                is_configdataclass
input_3(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInputConfig attribute),
        attribute), 51                                83
input_4(hvl_ccb.dev.supercube.constants.MeasurementsDividerRatio is_configdataclass
        attribute), 35                                (hvl_ccb.dev.newport.NewportSMC100PPConfig
input_4(hvl_ccb.dev.supercube.constants.MeasurementsScaledInputAttribute), 94
        attribute), 36                                is_configdataclass
input_4(hvl_ccb.dev.supercube2015.constants.MeasurementsScaledInputConfig attribute),
        attribute), 51                                103
INT32 (hvl_ccb.dev.se_ils2t.ILS2TRegDatatype at- is_configdataclass
        tribute), 117                                (hvl_ccb.dev.se_ils2t.ILS2TConfig attribute),
interface_type(hvl_ccb.comm.visa.VisaCommunicationConfig 115
        attribute), 19                                is_configdataclass
interface_type(hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig (drivenConfig.sst_luminox.LuminoxConfig
        attribute), 71                                attribute), 118
interface_type(hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig
        attribute), 111                                (hvl_ccb.dev.supercube.base.SupercubeConfiguration
internal (hvl_ccb.dev.labjack.LabJack.CjcType attribute), 24
        attribute), 78                                is_configdataclass
Internal(hvl_ccb.dev.supercube.constants.PowerSetup (hvl_ccb.dev.supercube2015.base.SupercubeConfiguration
        attribute), 37                                attribute), 45
Internal(hvl_ccb.dev.supercube2015.constants.PowerSetup_done () (hvl_ccb.dev.mbw973.MBW973 method),
        attribute), 52                                82
InvalidSupercubeStatusError, 42 is_error () (hvl_ccb.experiment_manager.ExperimentManager
IO_SCANNING (hvl_ccb.dev.se_ils2t.ILS2TRegAddr attribute), 116 method), 125
        attribute), 116                                is_finished () (hvl_ccb.experiment_manager.ExperimentManager
IoScanningModeValueError, 117 method), 125
is_configdataclass                                is_generic () (in module hvl_ccb.utils.typing), 123
        (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig_n_range () (hvl_ccb.dev.se_ils2t.ILS2TRegDatatype
attribute), 9                                method), 117

```

```

is_inactive(hvl_ccb.dev.crylas.CryLasLaser.LaserStatus)
    attribute), 62
is_open(hvl_ccb.comm.labjack_ljm.LJMCommunication)
    attribute), 8
is_open(hvl_ccb.comm.opc.OpcUaCommunication at-
    tribute), 12
is_open(hvl_ccb.comm.serial.SerialCommunication
    attribute), 14
is_polling() (hvl_ccb.dev.supercube.base.Poller
    method), 20
is_ready(hvl_ccb.dev.crylas.CryLasLaser.LaserStatus
    attribute), 62
is_running() (hvl_ccb.experiment_manager.ExperimentManager
    method), 125
is_type_hint() (in module hvl_ccb.utils.typing),
    123
is_valid_scale_range_reversed_str()
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig.Model
    method), 102

```

**J**

```

J (hvl_ccb.dev.labjack.LabJack.ThermocoupleType
    attribute), 79
jerk_time(hvl_ccb.dev.newport.NewportSMC100PPConfig
    attribute), 94
JOG (hvl_ccb.dev.se_il2t.ILS2T.Mode attribute), 111
jog_run() (hvl_ccb.dev.se_il2t.ILS2T method), 113
jog_stop() (hvl_ccb.dev.se_il2t.ILS2T method), 113
JOGGING (hvl_ccb.dev.newport.NewportStates
    attribute), 99
JOGGING_FROM_DISABLE
    (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    attribute), 87
JOGGING_FROM_READY
    (hvl_ccb.dev.newport.NewportSMC100PP.StateMessages
    attribute), 87
JOGN_FAST (hvl_ccb.dev.se_il2t.ILS2TRegAddr
    attribute), 116
JOGN_SLOW (hvl_ccb.dev.se_il2t.ILS2TRegAddr
    attribute), 116
JR (hvl_ccb.dev.newport.NewportConfigCommands
    attribute), 86

```

**K**

```

K (hvl_ccb.dev.labjack.LabJack.TemperatureUnit
    attribute), 78
K (hvl_ccb.dev.labjack.LabJack.ThermocoupleType
    attribute), 79
keys() (hvl_ccb.comm.labjack_ljm.LJMCommunication
    class method), 9
keys() (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication
    class method), 11
keys() (hvl_ccb.comm.opc.OpcUaCommunicationConfig
    class method), 13

```

```

keys() (hvl_ccb.comm.serial.SerialCommunicationConfig
    class method), 16
keys() (hvl_ccb.comm.visa.VisaCommunicationConfig
    class method), 19
keys() (hvl_ccb.dev.base.EmptyConfig class method),
    57
keys() (hvl_ccb.dev.crylas.CryLasAttenuatorConfig
    class method), 59
keys() (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig
    class method), 60
keys() (hvl_ccb.dev.crylas.CryLasLaserConfig class
    method), 65
keys() (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig
    class method), 67
keys() (hvl_ccb.dev.ea_psi9000.PSI9000Config class
    method), 70
keys() (hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig
    class method), 71
keys() (hvl_ccb.dev.heinzinger.HeinzingerConfig class
    method), 72
keys() (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
    class method), 77
keys() (hvl_ccb.dev.mbw973.MBW973Config class
    method), 83
keys() (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig
    class method), 84
keys() (hvl_ccb.dev.newport.NewportSMC100PPConfig
    class method), 94
keys() (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig
    class method), 98
keys() (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig
    class method), 103
keys() (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig
    class method), 105
keys() (hvl_ccb.dev.rs_rto1024.RTO1024Config class
    method), 110
keys() (hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig
    class method), 111
keys() (hvl_ccb.dev.se_il2t.ILS2TConfig class
    method), 115
keys() (hvl_ccb.dev.se_il2t.ILS2TModbusTcpCommunicationConfig
    class method), 116
keys() (hvl_ccb.dev.sst_luminox.LuminoxConfig class
    method), 118
keys() (hvl_ccb.dev.sst_luminox.LuminoxSerialCommunicationConfig
    class method), 120
keys() (hvl_ccb.dev.supercube.base.SupercubeConfiguration
    class method), 24
keys() (hvl_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig
    class method), 26
keys() (hvl_ccb.dev.supercube.typ_a.SupercubeAOpcUaConfiguration
    class method), 39
keys() (hvl_ccb.dev.supercube.typ_b.SupercubeBOpcUaConfiguration
    class method), 41

```

keys () (hvl\_ccb.dev.supercube2015.base.SupercubeConfiguration (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
keys () (hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaCommunication (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
keys () (hvl\_ccb.dev.supercube2015.typ\_a.SupercubeAOpcUaConfig (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
keys () (hvl\_ccb.dev.visa.VisaDeviceConfig class line\_9 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
kV (hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitVoltage attribute), 75

**L**

LabJack (class in hvl\_ccb.dev.labjack), 77  
LabJack.AInRange (class in hvl\_ccb.dev.labjack), 77  
LabJack.CalMicroAmpere (class in hvl\_ccb.dev.labjack), 78  
LabJack.CjcType (class in hvl\_ccb.dev.labjack), 78  
LabJack.DeviceType (class in hvl\_ccb.dev.labjack), 78  
LabJack.DIOStatus (class in hvl\_ccb.dev.labjack), 78  
LabJack.TemperatureUnit (class in hvl\_ccb.dev.labjack), 78  
LabJack.ThermocoupleType (class in hvl\_ccb.dev.labjack), 78  
LabJackError, 81  
LabJackIdentifierDIOError, 81  
laser\_off () (hvl\_ccb.dev.crylas.CryLasLaser method), 62  
laser\_on () (hvl\_ccb.dev.crylas.CryLasLaser method), 62  
line\_1 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_10 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_11 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_12 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_13 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_14 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_15 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_2 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_3 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_4 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36  
line\_5 (hvl\_ccb.dev.supercube.constants.MessageBoard attribute), 36

list\_directory () (hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 106  
live (hvl\_ccb.dev.supercube.constants.OpcControl attribute), 36

**LJMCommunication** (class in hvl\_ccb.comm.labjack\_ljm), 8  
**LJMCommunicationConfig** (class in hvl\_ccb.comm.labjack\_ljm), 8  
**LJMCommunicationConfig.ConnectionType** (class in hvl\_ccb.comm.labjack\_ljm), 9  
**LJMCommunicationConfig.DeviceType** (class in hvl\_ccb.comm.labjack\_ljm), 9  
**LJMCommunicationError**, 10  
lm34 (hvl\_ccb.dev.labjack.LabJack.CjcType attribute), 78  
load\_configuration () (hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 106  
local\_display () (hvl\_ccb.dev.rs\_rto1024.RTO1024 method), 106  
locked (hvl\_ccb.dev.supercube.constants.DoorStatus attribute), 32  
locked (hvl\_ccb.dev.supercube2015.constants.DoorStatus attribute), 48  
LOW (hvl\_ccb.dev.labjack.LabJack.DIOStatus attribute), 78  
Luminox (class in hvl\_ccb.dev.sst\_luminox), 117  
LuminoxConfig (class in hvl\_ccb.dev.sst\_luminox), 118  
LuminoxMeasurementType (class in hvl\_ccb.dev.sst\_luminox), 119  
LuminoxMeasurementTypeError, 119  
LuminoxOutputMode (class in hvl\_ccb.dev.sst\_luminox), 119  
LuminoxOutputModeError, 119  
LuminoxSerialCommunication (class in hvl\_ccb.dev.sst\_luminox), 119  
LuminoxSerialCommunicationConfig (class in hvl\_ccb.dev.sst\_luminox), 119

**M**

mA (hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitCurrent attribute), 74  
manual (hvl\_ccb.dev.supercube.constants.EarthingStick attribute), 32

M

manual (*hvl\_ccb.dev.supercube.constants.EarthingStick*  
*attribute*), 33  
**MOVING** *SerialCommunicationConfig* (class in  
*hvl\_ccb.dev.mbw973*), 84  
manual (*hvl\_ccb.dev.supercube2015.constants.EarthingStick*  
*attribute*), 49  
*measure()* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG*  
*method*), 101  
*measure\_all()* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG*  
*attribute*), 32  
*measure\_current()* (*hvl\_ccb.dev.heinzinger.HeinzingerDI* method),  
*attribute*), 49  
*measure\_voltage()*  
*measure\_voltage()* (*hvl\_ccb.dev.heinzinger.HeinzingerDI* method),  
*attribute*, 49  
*measure\_voltage\_current()* (*hvl\_ccb.dev.ea\_psi9000.PSI9000* method),  
*attribute*, 32  
*MeasurementsDividerRatio* (class in  
*hvl\_ccb.dev.supercube.constants*), 35  
*MeasurementsDividerRatio* (class in  
*hvl\_ccb.dev.supercube2015.constants*), 51  
*MeasurementsScaledInput* (class in  
*hvl\_ccb.dev.supercube.constants*), 35  
*MeasurementsScaledInput* (class in  
*hvl\_ccb.dev.supercube2015.constants*), 51  
*message* (*hvl\_ccb.dev.supercube.constants.Errors* at-  
*tribute*), 34  
*MessageBoard* (class in  
*hvl\_ccb.dev.supercube.constants*), 36  
*micro\_step\_per\_full\_step\_factor*  
*manuals* (*hvl\_ccb.dev.supercube.constants.EarthingStick*  
*attribute*), 33  
*MARK* (*hvl\_ccb.comm.serial.SerialCommunicationParity*  
*attribute*), 17  
*max\_current* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*  
*attribute*), 75  
*max\_current\_hardware* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC* at-  
*tribute*), 75  
*max\_timeout\_retry\_nr* (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig*  
*attribute*), 13  
*max\_voltage* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC*  
*attribute*), 75  
*max\_voltage\_hardware* (*hvl\_ccb.dev.heinzinger.HeinzingerPNC* at-  
*tribute*), 75  
*mbar* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits*  
*attribute*), 100  
*MBW973* (class in *hvl\_ccb.dev.mbw973*), 82  
*MBW973Config* (class in *hvl\_ccb.dev.mbw973*), 83  
*MBW973ControlRunningException*, 83  
*MBW973Error*, 84  
*MBW973PumpRunningException*, 84  
*MBW973SerialCommunication* (class in  
*hvl\_ccb.dev.mbw973*), 84  
**Moving** *SerialCommunicationConfig* (class in  
*hvl\_ccb.dev.mbw973*), 84  
*ModbusTcpCommunication* (class in  
*hvl\_ccb.comm.modbus\_tcp*), 10  
*ModbusTcpCommunicationConfig* (class in  
*hvl\_ccb.comm.modbus\_tcp*), 11  
*ModbusTcpConnectionFailedException*, 11  
*model* (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig* at-  
*tribute*), 103  
*motion\_distance\_per\_full\_step* (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
*attribute*), 94  
*motor\_config* (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
*attribute*), 94  
*move\_finished\_extra\_wait\_sec* (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
*attribute*), 94  
*move\_to\_absolute\_position()* (*hvl\_ccb.dev.newport.NewportSMC100PP*  
*method*), 90  
*move\_to\_relative\_position()* (*hvl\_ccb.dev.newport.NewportSMC100PP*  
*method*), 90  
*MOVING* (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages*

*attribute), 87*  
 MOVING (*hvl\_ccb.dev.newport.NewportStates attribute*), 99  
 MS\_NOMINAL\_CURRENT  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000 attribute)*, 67  
 MS\_NOMINAL\_VOLTAGE  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000 attribute)*, 67  
 MULTI\_COMMANDS\_MAX  
*(hvl\_ccb.comm.visa.VisaCommunication attribute)*, 17  
 MULTI\_COMMANDS\_SEPARATOR  
*(hvl\_ccb.comm.visa.VisaCommunication attribute)*, 17

**N**

NameEnum (*class in hvl\_ccb.utils.enum*), 123  
 NAMES (*hvl\_ccb.comm.serial.SerialCommunicationParity attribute*), 17  
 names (*hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes attribute*), 105  
 namespace\_index (*hvl\_ccb.dev.supercube.base.Supercube attribute*), 24  
 namespace\_index (*hvl\_ccb.dev.supercube2015.base.Supercube attribute*), 45  
 NED\_END\_OF\_TURN (*hvl\_ccb.dev.newport.NewportSMC100PP.MotorErrors attribute*), 86  
 NEG (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 112  
 NEG\_FAST (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 112  
 negative\_software\_limit  
*(hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute)*, 94  
 NewportConfigCommands (*class in hvl\_ccb.dev.newport*), 85  
 NewportControllerError, 86  
 NewportMotorError, 86  
 NewportSerialCommunicationError, 99  
 NewportSMC100PP (*class in hvl\_ccb.dev.newport*), 86  
 NewportSMC100PP.MotorErrors (*class in hvl\_ccb.dev.newport*), 86  
 NewportSMC100PP.StateMessages (*class in hvl\_ccb.dev.newport*), 86  
 NewportSMC100PPConfig (*class in hvl\_ccb.dev.newport*), 92  
 NewportSMC100PPConfig.EspStageConfig  
*(class in hvl\_ccb.dev.newport)*, 93  
 NewportSMC100PPConfig.HomeSearch (*class in hvl\_ccb.dev.newport*), 93  
 NewportSMC100PPSerialCommunication (*class in hvl\_ccb.dev.newport*), 95  
 NewportSMC100PPSerialCommunication.ControllerError  
*(class in hvl\_ccb.dev.newport)*, 95  
 NewportSMC100PPSerialCommunicationConfig  
*(class in hvl\_ccb.dev.newport)*, 97

NewportStates (*class in hvl\_ccb.dev.newport*), 99  
 NO\_ERROR (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication attribute*), 95  
 NO\_REF (*hvl\_ccb.dev.newport.NewportStates attribute*), 99  
 NO\_REF\_ESP\_STAGE\_ERROR  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_FROM\_CONFIG  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_FROM\_DISABLED  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_FROM\_HOMING  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_FROM\_JOGGING  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_FROM\_MOVING  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_CONFIGURATION  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_READY  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 NO\_REF\_RESET  
*(hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute)*, 87  
 No\_sensor (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute*), 100  
 NONE (*hvl\_ccb.comm.serial.SerialCommunicationParity attribute*), 17  
 NONE (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 79  
 None (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute*), 100  
 NONE (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 112  
 NoPower (*hvl\_ccb.dev.supercube.constants.PowerSetup attribute*), 37  
 NORMAL (*hvl\_ccb.dev.rs\_rto1024.RTO1024.TriggerModes attribute*), 105  
 noSen (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute*), 100  
 noSENSOR (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute*), 100  
 not\_defined (*hvl\_ccb.dev.supercube.constants.AlarmText attribute*), 27  
 not\_defined (*hvl\_ccb.dev.supercube2015.constants.AlarmText attribute*), 48  
 number (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 number\_of\_decimals

(*hvl\_ccb.dev.heinzinger.HeinzingerConfig attribute*), 72  
**O**  
 number\_of\_sensors (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG attribute*), 101  
 ODD (*hvl\_ccb.comm.serial.SerialCommunicationParity attribute*), 17  
 OH (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 86  
 Ok (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute*), 100  
 ON (*hvl\_ccb.dev.se\_ilst2t.ILS2T.State attribute*), 112  
 on\_start\_wait\_until\_ready (*hvl\_ccb.dev.crylas.CryLasLaserConfig attribute*), 65  
 ONE (*hvl\_ccb.comm.serial.SerialCommunicationStopbits attribute*), 17  
 ONE (*hvl\_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum attribute*), 72  
 ONE (*hvl\_ccb.dev.labjack.LabJack.AInRange attribute*), 77  
 ONE\_HUNDREDTH (*hvl\_ccb.dev.labjack.LabJack.AInRange attribute*), 77  
 ONE\_POINT\_FIVE (*hvl\_ccb.comm.serial.SerialCommunicationStopbits attribute*), 17  
 ONE\_TENTH (*hvl\_ccb.dev.labjack.LabJack.AInRange attribute*), 77  
 OpcControl (*class in hvl\_ccb.dev.supercube.constants*), 36  
 OpcUaCommunication (*class in hvl\_ccb.comm.opc*), 11  
 OpcUaCommunicationConfig (*class in hvl\_ccb.comm.opc*), 12  
 OpcUaCommunicationIOError, 14  
 OpcUaCommunicationTimeoutError, 14  
 OpcUaSubHandler (*class in hvl\_ccb.comm.opc*), 14  
 open (*hvl\_ccb.dev.supercube.constants.DoorStatus attribute*), 32  
 open (*hvl\_ccb.dev.supercube.constants.EarthingStickOperation attribute*), 33  
 open (*hvl\_ccb.dev.supercube.constants.EarthingStickStatus attribute*), 34  
 open (*hvl\_ccb.dev.supercube2015.constants.DoorStatus attribute*), 48  
 open (*hvl\_ccb.dev.supercube2015.constants.EarthingStickStatus attribute*), 50  
 open () (*hvl\_ccb.comm.base.CommunicationProtocol method*), 7  
 open () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunication method*), 8  
 open () (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication method*), 10  
 open () (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig class method*), 13  
 open\_shutter () (*hvl\_ccb.dev.crylas.CryLasLaser method*), 63  
 open\_timeout (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 19  
 OPENED (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersShutter attribute*), 61  
 OPENED (*hvl\_ccb.dev.crylas.CryLasLaserShutterStatus attribute*), 67  
 operate () (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 22  
 operate () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 44  
 operate\_earthing\_stick ()  
 operating\_status (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 operating\_status\_1 (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 operating\_status\_2 (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 operating\_status\_3 (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 operating\_status\_4 (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 operating\_status\_5 (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 operating\_status\_6 (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 operating\_statuses (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 optional\_defaults () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig class method*), 9  
 optional\_defaults () (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig class method*), 11  
 optional\_defaults () (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig class method*), 13  
 optional\_defaults ()

```

(hvl_ccb.comm.serial.SerialCommunicationConfig
    class method), 16
optional_defaults()
    (hvl_ccb.comm.visa.VisaCommunicationConfig
        class method), 19
optional_defaults()
    (hvl_ccb.dev.base.EmptyConfig class method),
        57
optional_defaults()
    (hvl_ccb.dev.crylas.CryLasAttenuatorConfig
        class method), 59
optional_defaults()
    (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig
        class method), 61
optional_defaults()
    (hvl_ccb.dev.crylas.CryLasLaserConfig   class
        method), 65
optional_defaults()
    (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig
        class method), 67
optional_defaults()
    (hvl_ccb.dev.ea_psi9000.PSI9000Config
        class method), 70
optional_defaults()
    (hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig
        class method), 71
optional_defaults()
    (hvl_ccb.dev.heinzinger.HeinzingerConfig
        class method), 72
optional_defaults()
    (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
        class method), 77
optional_defaults()
    (hvl_ccb.dev.mbw973.MBW973Config     class
        method), 83
optional_defaults()
    (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig
        class method), 85
optional_defaults()
    (hvl_ccb.dev.newport.NewportSMC100PPConfig
        class method), 94
optional_defaults()
    (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig
        class method), 99
optional_defaults()
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig
        class method), 103
optional_defaults()
    (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig
        class method), 105
optional_defaults()
    (hvl_ccb.dev.rs_rto1024.RTO1024Config
        class method), 110
optional_defaults()
    (hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig
        class method), 111
optional_defaults()
    (hvl_ccb.dev.se_il2t.ILS2TConfig      class
        method), 115
optional_defaults()
    (hvl_ccb.dev.se_il2t.ILS2TModbusTcpCommunicationConfig
        class method), 116
optional_defaults()
    (hvl_ccb.dev.sst_luminox.LuminoxConfig
        class method), 118
optional_defaults()
    (hvl_ccb.dev.sst_luminox.LuminoxSerialCommunicationConfig
        class method), 120
optional_defaults()
    (hvl_ccb.dev.supercube.base.SupercubeConfiguration
        class method), 24
optional_defaults()
    (hvl_ccb.dev.supercube.typ_a.SupercubeAOpcUaConfiguration
        class method), 39
optional_defaults()
    (hvl_ccb.dev.supercube.typ_b.SupercubeBOpcUaConfiguration
        class method), 41
optional_defaults()
    (hvl_ccb.dev.supercube2015.base.SupercubeConfiguration
        class method), 45
optional_defaults()
    (hvl_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig
        class method), 47
optional_defaults()
    (hvl_ccb.dev.supercube2015.typ_a.SupercubeAOpcUaConfiguration
        class method), 55
optional_defaults()
    (hvl_ccb.dev.visa.VisaDeviceConfig      class
        method), 122
OT (hvl_ccb.dev.newport.NewportConfigCommands at-
tribute), 86
out_1_1 (hvl_ccb.dev.supercube.constants.GeneralSupport
attribute), 35
out_1_2 (hvl_ccb.dev.supercube.constants.GeneralSupport
attribute), 51
out_1_2 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
attribute), 51
out_2_1 (hvl_ccb.dev.supercube2015.constants.GeneralSupport
attribute), 51
out_2_2 (hvl_ccb.dev.supercube.constants.GeneralSupport
attribute), 35

```

out\_2\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_3\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_3\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_3\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_3\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_4\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_4\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_4\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_4\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_5\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_5\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_5\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_5\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_6\_1 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_6\_1 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
out\_6\_2 (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), 35  
out\_6\_2 (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
output (hvl\_ccb.dev.sst\_luminox.Luminox attribute), percent\_o2 (hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType attribute), 118  
output (hvl\_ccb.dev.supercube.constants.GeneralSupport attribute), PfeifferTPG (class in hvl\_ccb.dev.pfeiffer\_tpg), 99  
output (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 35  
PfeifferTPG.PressureUnits (class in hvl\_ccb.dev.pfeiffer\_tpg), 100  
output (hvl\_ccb.dev.supercube2015.constants.GeneralSupport attribute), 51  
PfeifferTPG.SensorStatus (class in hvl\_ccb.dev.pfeiffer\_tpg), 100  
output\_off () (hvl\_ccb.dev.heinzinger.HeinzingerDI method), 74  
PfeifferTPG.SensorTypes (class in hvl\_ccb.dev.pfeiffer\_tpg), 100  
output\_on () (hvl\_ccb.dev.heinzinger.HeinzingerDI method), 74  
PfeifferTPGConfig (class in hvl\_ccb.dev.pfeiffer\_tpg), 100  
OUTPUT\_POWER\_EXCEEDED  
PfeifferTPGConfig.Model (class in hvl\_ccb.dev.pfeiffer\_tpg), 102  
(hvl\_ccb.dev.newport.NewportSMC100PP.MotorErrors attribute), 86  
Overrange (hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus attribute), 100  
PfeifferTPGError, 103  
PfeifferTPGSerialCommunication (class in hvl\_ccb.dev.pfeiffer\_tpg), 103  
PfeifferTPGSerialCommunicationConfig (class in hvl\_ccb.dev.pfeiffer\_tpg), 104  
PARAM\_MISSING\_OR\_INVALID  
(hvl\_ccb.dev.newport.NewportSMC100PP.SerialCommunication attribute), 95  
PfeifferTPGCommunicationConfig (class in hvl\_ccb.dev.pfeiffer\_tpg), 104  
PfeifferTPGSensorTypes (class in hvl\_ccb.dev.pfeiffer\_tpg), 100

P

`PARAM_MISSING_OR_INVALID` (*class in hvl\_ccb.dev.pfeiffer\_tpg*), 104  
`(hvl_ccb.dev.newport.NewportSMC100PPSSerialCommunicationCommitter.PfeifferTPG.SensorTypes attribute)`, 95  
`(hvl_ccb.dev.newport.NewportSMC100PPSSerialCommunicationCommitter.PfeifferTPG.SensorTypes attribute)`, 100

Poller (*class in hvl\_ccb.dev.mbw973*), 85  
 Poller (*class in hvl\_ccb.dev.supercube.base*), 20  
 polling (*hvl\_ccb.dev.sst\_luminox.LuminoxOutputMode attribute*), 119  
 polling\_delay\_sec (*hvl\_ccb.dev.supercube.base.SupercubeConfiguration attribute*), 25  
 polling\_interval (*hvl\_ccb.dev.mbw973.MBW973Config attribute*), 83  
 polling\_interval\_sec (*hvl\_ccb.dev.supercube.base.SupercubeConfiguration attribute*), 25  
 polling\_period (*hvl\_ccb.dev.crylas.CryLasLaserConfig attribute*), 65  
 polling\_timeout (*hvl\_ccb.dev.crylas.CryLasLaserConfig attribute*), 65  
 port (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig attribute*), 11  
 port (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig attribute*), 13  
 port (*hvl\_ccb.comm.serial.SerialCommunicationConfig attribute*), 16  
 port (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 19  
 port (*hvl\_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig attribute*), 47  
 port\_range (*hvl\_ccb.dev.supercube.constants.GeneralSupport attribute*), 35  
 POS (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 112  
 POS\_END\_OF\_TURN (*hvl\_ccb.dev.newport.NewportSMC100PP.MotorMethods attribute*), 86  
 POS\_FAST (*hvl\_ccb.dev.se\_ils2t.ILS2T.Ref16Jog attribute*), 112  
 POSITION (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr attribute*), 116  
 POSITION\_OUT\_OF\_LIMIT (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 95  
 positive\_software\_limit (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 94  
 post\_force\_value () (*hvl\_ccb.dev.newport.NewportSMC100PPConfig method*), 94  
 Power (*class in hvl\_ccb.dev.supercube.constants*), 36  
 Power (*class in hvl\_ccb.dev.supercube2015.constants*), 51  
 power\_limit (*hvl\_ccb.dev.ea\_psi9000.PSI9000Config attribute*), 70  
 PowerSetup (*class hvl\_ccb.dev.supercube.constants*), 37  
 PowerSetup (*class hvl\_ccb.dev.supercube2015.constants*), 52  
 prepare\_ultra\_segmentation ()  
**R**  
 R (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType method*), 106  
 PSI9000 (*class in hvl\_ccb.dev.ea\_psi9000*), 67  
 PSI9000Config (*class in hvl\_ccb.dev.ea\_psi9000*), 69  
 PSI9000Error, 70  
 PT1000VisaCommunication (*class in hvl\_ccb.dev.ea\_psi9000*), 70  
 PT1000VisaCommunicationConfig (*class in hvl\_ccb.dev.ea\_psi9000*), 71  
 PT100 (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 79  
 PT1000 (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 79  
 PT500 (*hvl\_ccb.dev.labjack.LabJack.ThermocoupleType attribute*), 79  
 PTP (*hvl\_ccb.dev.se\_ils2t.ILS2T.Mode attribute*), 111  
**Q**  
 QIL (*hvl\_ccb.dev.newport.NewportConfigCommands attribute*), 86  
 query () (*hvl\_ccb.comm.visa.VisaCommunication method*), 17  
 query () (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunication method*), 65  
 query () (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication method*), 96  
 query\_endy () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunication method*), 103  
 query\_all () (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunication methods*), 66  
 query\_multiple () (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication method*), 96  
 query\_single\_measurement () (*hvl\_ccb.dev.sst\_luminox.Luminox method*), 118  
 QUICKSTOP (*hvl\_ccb.dev.se\_ils2t.ILS2T.State attribute*), 34  
 QuickStop (*hvl\_ccb.dev.supercube.constants.SafetyStatus attribute*), 38  
 QuickStop (*hvl\_ccb.dev.supercube2015.constants.SafetyStatus attribute*), 53  
 quickstop () (*hvl\_ccb.dev.se\_ils2t.ILS2T method*), 113  
 quit (*hvl\_ccb.dev.supercube.constants.Errors attribute*), 34  
 quit (*hvl\_ccb.dev.supercube2015.constants.Errors attribute*), 50  
 quit\_error () (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 22  
 quit\_error () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 44

*tribute), 79*  
 RAMP\_ACC (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr attribute*), 116  
 RAMP\_DECEL (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr attribute*), 116  
 RAMP\_N\_MAX (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr attribute*), 116  
 RAMP\_TYPE (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr attribute*), 116  
 range (*hvl\_ccb.dev.supercube.constants.EarthingStick attribute*), 33  
 read () (*hvl\_ccb.comm.opc.OpcUaCommunication method*), 12  
 read () (*hvl\_ccb.dev.mbw973.MBW973 method*), 82  
 read () (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 22  
 read () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 44  
 READ\_ACTIVE (*hvl\_ccb.dev.crylas.CryLasLaser.LaserStatus attribute*), 61  
 read\_bytes () (*hvl\_ccb.comm.serial.SerialCommunication method*), 14  
 read\_float () (*hvl\_ccb.dev.mbw973.MBW973 method*), 82  
 read\_holding\_registers () (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication method*), 10  
 read\_input\_registers () (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunication method*), 10  
 read\_int () (*hvl\_ccb.dev.mbw973.MBW973 method*), 82  
 read\_measurements () (*hvl\_ccb.dev.mbw973.MBW973 method*), 82  
 read\_name () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunication method*), 8  
 read\_resistance () (*hvl\_ccb.dev.labjack.LabJack method*), 80  
 read\_stream () (*hvl\_ccb.dev.sst\_luminox.Luminox method*), 118  
 read\_termination (*hvl\_ccb.comm.visa.VisaCommunication attribute*), 19  
 read\_text () (*hvl\_ccb.comm.serial.SerialCommunication method*), 15  
 read\_text () (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunication method*), 66  
 read\_text\_nonempty () (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunication method*), 76  
 READ\_TEXT\_SKIP\_PREFIXES (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunication attribute*), 65  
 read\_thermocouple ()

READY (*hvl\_ccb.dev.labjack.LabJack method*), 80  
 READY (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute*), 61  
 READY (*hvl\_ccb.dev.newport.NewportStates attribute*), 99  
 READY (*hvl\_ccb.dev.se\_ils2t.ILS2T.State attribute*), 112  
 ready () (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 22  
 ready () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 44  
 READY\_FROM\_DISABLE (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 87  
 READY\_FROM\_HOMING (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 87  
 READY\_FROM\_JOGGING (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 87  
 READY\_FROM\_MOVING (*hvl\_ccb.dev.newport.NewportSMC100PP.StateMessages attribute*), 87  
 READY\_INACTIVE (*hvl\_ccb.dev.crylas.CryLasLaser.LaserStatus attribute*), 61  
 RedOperate (*hvl\_ccb.dev.supercube.constants.SafetyStatus attribute*), 38  
 RedOperate (*hvl\_ccb.dev.supercube2015.constants.SafetyStatus attribute*), 53  
 Ready (*hvl\_ccb.dev.supercube.constants.SafetyStatus attribute*), 38  
 Ready (*hvl\_ccb.dev.supercube2015.constants.SafetyStatus attribute*), 53  
 RegAddr (*hvl\_ccb.dev.se\_ils2t.ILS2T attribute*), 112  
 RegDatatype (*hvl\_ccb.dev.se\_ils2t.ILS2T attribute*), 112  
 RELATIVE\_POSITION\_MOTOR (*hvl\_ccb.dev.se\_ils2t.ILS2T.ActionsPtp attribute*), 111  
 RELATIVE\_POSITION\_TARGET (*hvl\_ccb.dev.se\_ils2t.ILS2T.ActionsPtp attribute*), 111  
 ConfigDevice () (*hvl\_ccb.dev.base.DeviceSequenceMixin method*), 56  
 required\_keys () (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig class method*), 10  
 Communication\_keys () (*hvl\_ccb.comm.modbus\_tcp.ModbusTcpCommunicationConfig class method*), 11  
 required\_keys () (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig class method*), 13  
 required\_keys () (*hvl\_ccb.comm.serial.SerialCommunicationConfig class method*), 16  
 required\_keys () (*hvl\_ccb.comm.visa.VisaCommunicationConfig class method*), 19  
 required\_keys () (*hvl\_ccb.dev.base.EmptyConfig*)

```

    class method), 57
required_keys() (hvl_ccb.dev.crylas.CryLasAttenuatorConfigured_keys () (hvl_ccb.dev.visa.VisaDeviceConfig
    class method), 59
required_keys() (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfigurable.constants.BreakdownDetection
    class method), 61
required_keys() (hvl_ccb.dev.crylas.CryLasLaserConfigset (hvl_ccb.dev.supercube2015.constants.BreakdownDetection
    class method), 65
required_keys() (hvl_ccb.dev.crylas.CryLasLaserSerialCommunicationConfigurable.constants.BreakdownDetection
    class method), 67
required_keys() (hvl_ccb.dev.ea_psi9000.PSI9000Configset () (hvl_ccb.dev.visa.VisaDevice method), 121
    class method), 70
required_keys() (hvl_ccb.dev.ea_psi9000.PSI9000VisaCommunicationConfig
    class method), 71
required_keys() (hvl_ccb.dev.heinzinger.HeinzingerConfig      (hvl_ccb.dev.heinzinger.HeinzingerDI method),
    class method), 73
required_keys() (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
    class method), 77
required_keys() (hvl_ccb.dev.mbw973.MBW973Config           (hvl_ccb.dev.crylas.CryLasAttenuatorConfig
    class method), 83
required_keys() (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfignewport.NewportSMC100PP.MotorErrors
    class method), 85
required_keys() (hvl_ccb.dev.newport.NewportSMC100PPConfiginit (hvl_ccb.dev.se_ils2t.ILS2TConfig at-
    class method), 95
required_keys() (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfignewport.NewportSMC100PPConfig
    class method), 99
required_keys() (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGConfig24 (class in hvl_ccb.dev.rs_rto1024), 105
    class method), 103
required_keys() (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig024), 105
    class method), 105
required_keys() (hvl_ccb.dev.rs_rto1024.RTO1024Config      109
    class method), 110
required_keys() (hvl_ccb.dev.rs_rto1024.RTO1024VisaCommunicationConfig      (class      in
    class method), 111
required_keys() (hvl_ccb.dev.se_ils2t.ILS2TConfig          RTO1024VisaCommunicationConfig (class      in
    class method), 115
hvl_ccb.dev.rs_rto1024), 110
required_keys() (hvl_ccb.dev.se_ils2t.ILS2TModbusTcpGetIn(hvl_ccb.dev.rs_rto1024), 112
    class method), 116
run () (hvl_ccb.experiment_manager.ExperimentManager
required_keys() (hvl_ccb.dev.sst_luminox.LuminoxConfig      method), 125
    class method), 119
run_continuous_acquisition()
required_keys() (hvl_ccb.dev.sst_luminox.LuminoxSerialCommunicationConfigrto1024.RTO1024      method),
    class method), 120
106
required_keys() (hvl_ccb.dev.supercube.base.SupercubeConfigurationacquisition ()
    class method), 25
(hvl_ccb.dev.rs_rto1024.RTO1024      method),
required_keys() (hvl_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig
    class method), 26
RUNNING (hvl_ccb.experiment_manager.ExperimentStatus
required_keys() (hvl_ccb.dev.supercube.typ_a.SupercubeAOpcUaConfiguration
    class method), 39
required_keys() (hvl_ccb.dev.supercube.typ_b.SupercubeSeBOpcUaConfiguration
    class method), 42
S (hvl_ccb.dev.labjack.LabJack.ThermocoupleType at-
required_keys() (hvl_ccb.dev.supercube2015.base.SupercubeConfiguration
    class method), 46
SA (hvl_ccb.dev.newport.NewportConfigCommands at-
required_keys() (hvl_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig
    class method), 47
Safety (class in hvl_ccb.dev.supercube.constants), 37
required_keys() (hvl_ccb.dev.supercube2015.typ_a.SupercubeAOpcUaConfiguration

```

Safety (*class in hvl\_ccb.dev.supercube2015.constants*), 52  
 SafetyStatus (*class hvl\_ccb.dev.supercube.constants*), 38  
 SafetyStatus (*class hvl\_ccb.dev.supercube2015.constants*), 53  
 save\_configuration () (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 107  
 save\_waveform\_history () (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 107  
 SCALE (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr attribute*), 116  
 ScalingFactorValueError, 117  
 screw\_scaling (*hvl\_ccb.dev.newport.NewportSMC100PPConfig attribute*), 95  
 send\_command () (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication rto1024.RTO1024 method*), 96  
 send\_command () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunication method*), 103  
 send\_stop () (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunication method*), 97  
 Sensor\_error (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus (hvl\_ccb.dev.rs\_rto1024.RTO1024 method) attribute*), 100  
 Sensor\_off (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus current (hvl\_ccb.dev.heinzinger.HeinzingerDI method) attribute*), 100  
 sensor\_status (*hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType (hvl\_ccb.dev.heinzinger.HeinzingerPNC method) attribute*), 119  
 serial\_number (*hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType digital\_output (hvl\_ccb.dev.labjack.LabJack method) attribute*), 119  
 SerialCommunication (*class hvl\_ccb.comm.serial*), 14  
 SerialCommunicationBytesize (*class hvl\_ccb.comm.serial*), 15  
 SerialCommunicationConfig (*class hvl\_ccb.comm.serial*), 15  
 SerialCommunicationIOError, 16  
 SerialCommunicationParity (*class hvl\_ccb.comm.serial*), 17  
 SerialCommunicationStopbits (*class hvl\_ccb.comm.serial*), 17  
 set\_acceleration () (*hvl\_ccb.dev.newport.NewportSMC100PP method*), 91  
 set\_acquire\_length () (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 107  
 set\_ain\_differential () (*hvl\_ccb.dev.labjack.LabJack method*), 80  
 set\_ain\_range () (*hvl\_ccb.dev.labjack.LabJack method*), 80  
 set\_ain\_resistance () (*hvl\_ccb.dev.labjack.LabJack method*), 80  
 set\_ain\_resolution () (*hvl\_ccb.dev.labjack.LabJack method*), 81  
 set\_ain\_thermocouple () (*hvl\_ccb.dev.labjack.LabJack method*), 81  
 set\_attenuation () (*hvl\_ccb.dev.crylas.CryLasAttenuator method*), 58  
 set\_cee16\_socket () (*hvl\_ccb.dev.supercube.base.SupercubeBase method*), 23  
 set\_channel\_position () (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 107  
 set\_channel\_range () (*hvl\_ccb.dev.rs\_rto1024.RTO1024 method*), 107  
 set\_channel\_state ()  
 Sensor\_error (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorStatus (hvl\_ccb.dev.rs\_rto1024.RTO1024 method) attribute*), 108  
 Sensor\_Statuscurrent (*hvl\_ccb.dev.heinzinger.HeinzingerDI method*), 74  
 sensor\_Statuscurrent (*hvl\_ccb.dev.heinzinger.HeinzingerPNC method*), 75  
 serial\_number (*hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType digital\_output (hvl\_ccb.dev.labjack.LabJack method) attribute*), 81  
 set\_display\_unit () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG method*), 102  
 set\_earthing\_manual () (*hvl\_ccb.dev.supercube2015.base.Supercube2015Base method*), 44  
 set\_full\_scale\_mbar () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG method*), 102  
 set\_full\_scale\_unitless () (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG method*), 102  
 set\_init\_attenuation () (*hvl\_ccb.dev.crylas.CryLasAttenuator method*), 58  
 set\_init\_shutter\_status () (*hvl\_ccb.dev.crylas.CryLasLaser method*), 63  
 set\_jog\_speed () (*hvl\_ccb.dev.se\_ils2t.ILS2T method*), 113  
 set\_lower\_limits () (*hvl\_ccb.dev.ea\_psi9000.PSI9000 method*), 68

```

set_max_acceleration()
    (hvl_ccb.dev.se_ilst2t.ILS2T method), 114
set_max_deceleration()
    (hvl_ccb.dev.se_ilst2t.ILS2T method), 114
set_max_rpm() (hvl_ccb.dev.se_ilst2t.ILS2T method),
    114
set_measuring_options()
    (hvl_ccb.dev.mbw973.MBW973 method),
    82
set_message_board()
    (hvl_ccb.dev.supercube.base.SupercubeBase
method), 23
set_motor_configuration()
    (hvl_ccb.dev.newport.NewportSMC100PP
method), 91
set_negative_software_limit()
    (hvl_ccb.dev.newport.NewportSMC100PP
method), 91
set_number_of_recordings()
    (hvl_ccb.dev.heinzinger.HeinzingerDI method),
    74
set_output() (hvl_ccb.dev.ea_psi9000.PSI9000
method), 68
set_positive_software_limit()
    (hvl_ccb.dev.newport.NewportSMC100PP
method), 91
set_pulse_energy()
    (hvl_ccb.dev.crylas.CryLasLaser     method),
    63
set_ramp_type() (hvl_ccb.dev.se_ilst2t.ILS2T
method), 114
set_reference_point()
    (hvl_ccb.dev.rs_rto1024.RTO1024   method),
    108
set_remote_control()
    (hvl_ccb.dev.supercube.base.SupercubeBase
method), 23
set_remote_control()
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base
method), 44
set_repetition_rate()
    (hvl_ccb.dev.crylas.CryLasLaser     method),
    63
set_repetitions()
    (hvl_ccb.dev.rs_rto1024.RTO1024   method),
    108
set_slope() (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
method), 40
set_slope() (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
method), 54
set_status_board()
    (hvl_ccb.dev.supercube.base.SupercubeBase
method), 23
set_support_output()

```

```

    (hvl_ccb.dev.supercube.base.SupercubeBase
method), 23
set_support_output()
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base
method), 44
set_support_output_impulse()
    (hvl_ccb.dev.supercube.base.SupercubeBase
method), 23
set_support_output_impulse()
    (hvl_ccb.dev.supercube2015.base.Supercube2015Base
method), 44
set_system_lock()
    (hvl_ccb.dev.ea_psi9000.PSI9000     method),
    69
set_t13_socket() (hvl_ccb.dev.supercube.base.SupercubeBase
method), 24
set_t13_socket() (hvl_ccb.dev.supercube2015.base.Supercube2015Base
method), 45
set_target_voltage()
    (hvl_ccb.dev.supercube.typ_a.SupercubeWithFU
method), 40
set_target_voltage()
    (hvl_ccb.dev.supercube2015.typ_a.Supercube2015WithFU
method), 54
set_transmission()
    (hvl_ccb.dev.crylas.CryLasAttenuator method),
    58
set_trigger_level()
    (hvl_ccb.dev.rs_rto1024.RTO1024   method),
    108
set_trigger_mode()
    (hvl_ccb.dev.rs_rto1024.RTO1024   method),
    109
set_trigger_source()
    (hvl_ccb.dev.rs_rto1024.RTO1024   method),
    109
set_upper_limits()
    (hvl_ccb.dev.ea_psi9000.PSI9000   method),
    69
set_voltage() (hvl_ccb.dev.heinzinger.HeinzingerDI
method), 74
set_voltage() (hvl_ccb.dev.heinzinger.HeinzingerPNC
method), 75
set_voltage_current()
    (hvl_ccb.dev.ea_psi9000.PSI9000   method),
    69
set_voltage_up() (hvl_ccb.dev.supercube.constants.Power
attribute), 37
set_voltage2015WithFU() (hvl_ccb.dev.supercube2015.constants.Power
attribute), 52
SEVENBITS (hvl_ccb.comm.serial.SerialCommunicationBytesize
attribute), 15
SHORT_CIRCUIT (hvl_ccb.dev.newport.NewportSMC100PP.MotorErrors
attribute), 86

```

SHUTDOWN_CURRENT_LIMIT ( <i>hvl_ccb.dev.ea_psi9000.PSI9000 attribute</i> ), 68	start ()     ( <i>hvl_ccb.dev.newport.NewportSMC100PP method</i> ), 92
SHUTDOWN_VOLTAGE_LIMIT ( <i>hvl_ccb.dev.ea_psi9000.PSI9000 attribute</i> ), 68	start ()     ( <i>hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method</i> ), 102
ShutterStatus ( <i>hvl_ccb.dev.crylas.CryLasLaser attribute</i> ), 62	start ()     ( <i>hvl_ccb.dev.rs_rto1024.RTO1024 method</i> ), 109
ShutterStatus ( <i>hvl_ccb.dev.crylas.CryLasLaserConfig attribute</i> ), 64	start ()     ( <i>hvl_ccb.dev.se_ilst2.ILS2T method</i> ), 114
SingleCommDevice (class in <i>hvl_ccb.dev.base</i> ), 57	start ()     ( <i>hvl_ccb.dev.sst_luminox.Luminox method</i> ), 118
SIXBITS ( <i>hvl_ccb.comm.serial.SerialCommunicationBytes attribute</i> ), 15	start ()     ( <i>hvl_ccb.dev.supercube.base.SupercubeBase method</i> ), 24
SIXTEEN ( <i>hvl_ccb.dev.heinzinger.HeinzingerConfig.RecordingStateEnum attribute</i> ), 72	start ()     ( <i>hvl_ccb.dev.supercube2015.base.Supercube2015Base method</i> ), 45
SL ( <i>hvl_ccb.dev.newport.NewportConfigCommands attribute</i> ), 86	start ()     ( <i>hvl_ccb.dev.visa.VisaDevice method</i> ), 121
SOFTWARE_INTERNAL_SIXTY ( <i>hvl_ccb.dev.crylas.CryLasLaser.RepetitionRates attribute</i> ), 62	start ()     ( <i>hvl_ccb.experiment_manager.ExperimentManager method</i> ), 125
SOFTWARE_INTERNAL_TEN ( <i>hvl_ccb.dev.crylas.CryLasLaser.RepetitionRates attribute</i> ), 62	start_control ()     ( <i>hvl_ccb.dev.mbw973.MBW973 method</i> ), 83
SOFTWARE_INTERNAL_TWENTY ( <i>hvl_ccb.dev.crylas.CryLasLaser.RepetitionRates attribute</i> ), 62	start_polling ()     ( <i>hvl_ccb.dev.supercube.base.Poller method</i> ), 20
software_revision ( <i>hvl_ccb.dev.sst_luminox.LuminoxMeasurementType attribute</i> ), 119	STARTING ( <i>hvl_ccb.experiment_manager.ExperimentStatus attribute</i> ), 125
SPACE ( <i>hvl_ccb.comm.serial.SerialCommunicationParity attribute</i> ), 17	States     ( <i>hvl_ccb.dev.newport.NewportSMC100PP attribute</i> ), 87
spoll ()     ( <i>hvl_ccb.comm.visa.VisaCommunication method</i> ), 18	status     ( <i>hvl_ccb.dev.supercube.constants.EarthingStick attribute</i> ), 33
spoll_handler ()     ( <i>hvl_ccb.dev.visa.VisaDevice method</i> ), 121	status     ( <i>hvl_ccb.dev.supercube.constants.Safety attribute</i> ), 37
SR ( <i>hvl_ccb.dev.newport.NewportConfigCommands attribute</i> ), 86	status     ( <i>hvl_ccb.experiment_manager.ExperimentManager attribute</i> ), 125
stage_configuration ( <i>hvl_ccb.dev.newport.NewportSMC100PPConfig attribute</i> ), 95	status_1     ( <i>hvl_ccb.dev.supercube.constants.Door attribute</i> ), 32
start ()     ( <i>hvl_ccb.dev.base.Device method</i> ), 56	status_1     ( <i>hvl_ccb.dev.supercube.constants.EarthingStick attribute</i> ), 33
start ()     ( <i>hvl_ccb.dev.base.DeviceSequenceMixin method</i> ), 57	status_1_closed ( <i>hvl_ccb.dev.supercube2015.constants.EarthingStick attribute</i> ), 49
start ()     ( <i>hvl_ccb.dev.base.SingleCommDevice method</i> ), 57	status_1_connected ( <i>hvl_ccb.dev.supercube2015.constants.EarthingStick attribute</i> ), 49
start ()     ( <i>hvl_ccb.dev.crylas.CryLasAttenuator method</i> ), 58	status_1_open ( <i>hvl_ccb.dev.supercube2015.constants.EarthingStick attribute</i> ), 49
start ()     ( <i>hvl_ccb.dev.crylas.CryLasLaser method</i> ), 63	status_2     ( <i>hvl_ccb.dev.supercube.constants.Door attribute</i> ), 32
start ()     ( <i>hvl_ccb.dev.ea_psi9000.PSI9000 method</i> ), 69	status_2     ( <i>hvl_ccb.dev.supercube.constants.EarthingStick attribute</i> ), 33
start ()     ( <i>hvl_ccb.dev.heinzinger.HeinzingerDI method</i> ), 74	status_2_closed ( <i>hvl_ccb.dev.supercube2015.constants.EarthingStick attribute</i> ), 49
start ()     ( <i>hvl_ccb.dev.heinzinger.HeinzingerPNC method</i> ), 75	status_2_connected ( <i>hvl_ccb.dev.supercube2015.constants.EarthingStick attribute</i> ), 49
start ()     ( <i>hvl_ccb.dev.labjack.LabJack method</i> ), 81	status_2_open ( <i>hvl_ccb.dev.supercube2015.constants.EarthingStick attribute</i> ), 49
start ()     ( <i>hvl_ccb.dev.mbw973.MBW973 method</i> ), 83	status_3     ( <i>hvl_ccb.dev.supercube.constants.Door attribute</i> ), 32
start ()     ( <i>hvl_ccb.dev.mbw973.Poller method</i> ), 85	

```

status_3 (hvl_ccb.dev.supercube.constants.EarthingStick           34
         attribute), 33
stop (hvl_ccb.dev.supercube2015.constants.Errors at-
status_3_closed (hvl_ccb.dev.supercube2015.constants.EarthingStick 50
                attribute), 49
stop () (hvl_ccb.dev.base.Device method), 56
status_3_connected (hvl_ccb.dev.supercube2015.constants.EarthingStick 57
                   attribute), 49
stop () (hvl_ccb.dev.base.DeviceSequenceMixin
status_3_open (hvl_ccb.dev.supercube2015.constants.EarthingStick
               attribute), 49
stop () (hvl_ccb.dev.crylas.CryLasLaser method), 63
status_4 (hvl_ccb.dev.supercube.constants.EarthingStick stop () (hvl_ccb.dev.ea_psi9000.PSI9000 method), 69
         attribute), 33
stop () (hvl_ccb.dev.heinzinger.HeinzingerDI method),
status_4_closed (hvl_ccb.dev.supercube2015.constants.EarthingStick
                 attribute), 49
stop () (hvl_ccb.dev.labjack.LabJack method), 81
status_4_connected (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hvl_ccb.dev.mbw973.MBW973 method), 83
                   attribute), 49
stop () (hvl_ccb.dev.mbw973.Poller method), 85
status_4_open (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hvl_ccb.dev.newport.NewportSMC100PP
               attribute), 49
stop () (hvl_ccb.dev.pfeiffer_tpg.PfeifferTPG method),
status_5 (hvl_ccb.dev.supercube.constants.EarthingStick 102
         attribute), 33
stop () (hvl_ccb.dev.rs_rto1024.RTO1024 method),
status_5_closed (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hvl_ccb.dev.se_ilst.ILS2T method), 114
                  attribute), 49
stop () (hvl_ccb.dev.sst_luminox.Luminox method), 118
status_5_connected (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hvl_ccb.dev.supercube.base.SupercubeBase
                   attribute), 49
method), 24
status_5_open (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hkl_ccb.dev.supercube2015.base.Supercube2015Base
               attribute), 49
method), 45
status_6 (hvl_ccb.dev.supercube.constants.EarthingStick stop () (hvl_ccb.dev.visa.VisaDevice method), 121
         attribute), 33
stop () (hvl_ccb.dev.visa.VisaStatusPoller method), 122
status_6_closed (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hlcckcb.experiment_manager.ExperimentManager
                    attribute), 49
method), 125
status_6_connected (hvl_ccb.dev.supercube2015.constants.EarthingStick stop_acquisition()
                   attribute), 49
(hvl_ccb.dev.supercube2015.constants.EarthingStick) (hvl_ccb.dev.rs_rto1024.RTO1024 method), 109
status_6_open (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hvl_ccb.dev.newport.NewportSMC100PP
               attribute), 49
method), 92
status_closed (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hvl_ccb.dev.supercube2015.constants.Errors
                     attribute), 50
attribute), 50
status_connected (hvl_ccb.dev.supercube2015.constants.EarthingStick stop () (hvl_ccb.dev.supercube.base.Poller
                   attribute), 49
method), 20
status_error (hvl_ccb.dev.supercube2015.constants.Safety stopbits (hvl_ccb.comm.serial.SerialCommunicationConfig
                   attribute), 52
attribute), 15
status_green (hvl_ccb.dev.supercube2015.constants.Safety stopbits (hvl_ccb.comm.serial.SerialCommunicationConfig
                   attribute), 52
attribute), 16
status_open (hvl_ccb.dev.supercube2015.constants.EarthingStick stopbits (hvl_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig
                    attribute), 49
attribute), 61
status_ready_for_red (hvl_ccb.dev.supercube2015.constants.Safety stopbits (hvl_ccb.dev.crylas.CryLasSerialCommunicationConfig
                   attribute), 52
attribute), 67
status_red (hvl_ccb.dev.supercube2015.constants.Safety stopbits (hvl_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig
                   attribute), 53
attribute), 77
stopbits (hvl_ccb.dev.mbw973.MBW973SerialCommunicationConfig
statuses (hvl_ccb.dev.supercube.constants.EarthingStick stopbits (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig
                   attribute), 33
attribute), 85
stopbits (hvl_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig
stop (hvl_ccb.dev.supercube.constants.Errors attribute), 99
attribute), 99

```

stopbits ( <i>hvl_ccb.dev.pfeiffer_tpg.PfeifferTPGSerialCommunicationConfig</i> )	(class <i>hvl_ccb.dev.supercube.typ_a</i> )	39	<i>in</i>
stopbits ( <i>hvl_ccb.dev.sst_luminox.LuminoxSerialCommunicationConfig</i> )	(operate <i>hvl_ccb.dev.supercube.constants.Safety</i> )	at-	
streaming ( <i>hvl_ccb.dev.sst_luminox.LuminoxOutputMode</i> )	(attribute, 37)		
StrEnumBase (class in <i>hvl_ccb.utils.enum</i> )	(attribute, 38)		
sub_handler ( <i>hvl_ccb.comm.opc.OpcUaCommunicationConfig</i> )	(to_green <i>hvl_ccb.dev.supercube2015.constants.Safety</i> )	attribute), 53	
sub_handler ( <i>hvl_ccb.dev.supercube.base.SupercubeOpcUaCommunicationConfig</i> )	(attribute, 53)		
sub_handler ( <i>hvl_ccb.dev.supercube2015.base.SupercubeOpcUaCommunicationConfig</i> )	(attribute, 53)		
Supercube2015Base (class <i>hvl_ccb.dev.supercube2015.base</i> )	(attribute, 42)		
Supercube2015WithFU (class <i>hvl_ccb.dev.supercube2015.typ_a</i> )	(attribute, 53)		
SupercubeAOpcUaCommunication (class <i>hvl_ccb.dev.supercube.typ_a</i> )	(attribute, 38)		
SupercubeAOpcUaCommunication (class <i>hvl_ccb.dev.supercube2015.typ_a</i> )	(attribute, 54)		
SupercubeAOpcUaConfiguration (class <i>hvl_ccb.dev.supercube.typ_a</i> )	(attribute, 38)		
SupercubeAOpcUaConfiguration (class <i>hvl_ccb.dev.supercube2015.typ_a</i> )	(attribute, 54)		
SupercubeB (class in <i>hvl_ccb.dev.supercube.typ_b</i> )	(attribute, 41)		
SupercubeBase (class <i>hvl_ccb.dev.supercube.base</i> )	(attribute, 20)		
SupercubeBOpcUaCommunication (class <i>hvl_ccb.dev.supercube.typ_b</i> )	(attribute, 41)		
SupercubeBOpcUaConfiguration (class <i>hvl_ccb.dev.supercube.typ_b</i> )	(attribute, 41)		
SupercubeConfiguration (class <i>hvl_ccb.dev.supercube.base</i> )	(attribute, 24)		
SupercubeConfiguration (class <i>hvl_ccb.dev.supercube2015.base</i> )	(attribute, 45)		
SupercubeEarthlingStickOperationError,			
25			
SupercubeOpcEndpoint (class <i>hvl_ccb.dev.supercube.constants</i> )	(attribute, 38)		
SupercubeOpcEndpoint (class <i>hvl_ccb.dev.supercube2015.constants</i> )	(attribute, 53)		
SupercubeOpcUaCommunication (class <i>hvl_ccb.dev.supercube.base</i> )	(attribute, 25)		
SupercubeOpcUaCommunication (class <i>hvl_ccb.dev.supercube2015.base</i> )	(attribute, 46)		
SupercubeOpcUaCommunicationConfig (class in <i>hvl_ccb.dev.supercube.base</i> )	(attribute, 25)		
SupercubeOpcUaCommunicationConfig (class in <i>hvl_ccb.dev.supercube2015.base</i> )	(attribute, 46)		
SupercubeSubscriptionHandler (class <i>hvl_ccb.dev.supercube.base</i> )	(attribute, 26)		
SupercubeSubscriptionHandler (class <i>hvl_ccb.dev.supercube2015.base</i> )	(attribute, 47)		
		T	
		T (attribute, 79)	
		t13_1 (attribute, 34)	
		t13_1 (attribute, 50)	
		t13_2 (attribute, 34)	
		t13_2 (attribute, 50)	
		t13_3 (attribute, 34)	
		t13_3 (attribute, 50)	
		T13_SOCKET_PORTS (module <i>hvl_ccb.dev.supercube.constants</i> )	
		T13_SOCKET_PORTS (module <i>hvl_ccb.dev.supercube2015.constants</i> )	
		T4 (attribute, 9)	
		T4 (attribute, 78)	
		T7 (attribute, 9)	
		T7 (attribute, 78)	
		T7_PRO (attribute, 9)	
		T7_PRO (attribute, 78)	
		target_pulse_energy (attribute, 63)	
		TCP (attribute, 9)	
		TCPIP_INSTR (attribute, 18)	
		TCPIP_SOCKET (attribute, 18)	

TEC1 (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute*), 61

TEC2 (*hvl\_ccb.dev.crylas.CryLasLaser.AnswersStatus attribute*), 61

TEMP (*hvl\_ccb.dev.se\_ilS2t.ILS2TRegAddr attribute*), 117

temperature\_sensor (*hvl\_ccb.dev.sst\_luminox.LuminoxMeasurementType attribute*), 119

TEN (*hvl\_ccb.dev.labjack.LabJack.AInRange attribute*), 78

TEN (*hvl\_ccb.dev.labjack.LabJack.CalMicroAmpere attribute*), 78

terminator (*hvl\_ccb.comm.serial.SerialCommunicationConfig attribute*), 16

terminator (*hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig attribute*), 61

terminator (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig attribute*), 67

terminator (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig attribute*), 77

terminator (*hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig attribute*), 85

terminator (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig attribute*), 99

terminator (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig attribute*), 105

terminator (*hvl\_ccb.dev.sst\_luminox.LuminoxSerialCommunicationConfig attribute*), 121

terminator\_str () (*hvl\_ccb.comm.serial.SerialCommunicationConfig attribute*), 16

timeout (*hvl\_ccb.comm.serial.SerialCommunicationConfig attribute*), 16

timeout (*hvl\_ccb.comm.visa.VisaCommunicationConfig attribute*), 19

timeout (*hvl\_ccb.dev.crylas.CryLasAttenuatorSerialCommunicationConfig attribute*), 61

timeout (*hvl\_ccb.dev.crylas.CryLasLaserSerialCommunicationConfig attribute*), 67

timeout (*hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig attribute*), 77

timeout (*hvl\_ccb.dev.mbw973.MBW973SerialCommunicationConfig attribute*), 85

timeout (*hvl\_ccb.dev.newport.NewportSMC100PPSerialCommunicationConfig attribute*), 99

timeout (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGSerialCommunicationConfig attribute*), 105

timeout (*hvl\_ccb.dev.sst\_luminox.LuminoxSerialCommunicationConfig attribute*), 121

timer\_callback () (*hvl\_ccb.dev.mbw973.Poller method*), 85

Torr (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits attribute*), 100

TPG25xA (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig.Model attribute*), 102

TPGx6x (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPGConfig.Model attribute*), 102

TPR (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.SensorTypes attribute*), 100

transmission (*hvl\_ccb.dev.crylas.CryLasAttenuator attribute*), 58

triggered (*hvl\_ccb.dev.supercube.constants.BreakdownDetection attribute*), 32

triggered (*hvl\_ccb.dev.supercube2015.constants.BreakdownDetection attribute*), 48

TWO (*hvl\_ccb.comm.serial.SerialCommunicationStopbits attribute*), 17

UNHANDLING (*hvl\_ccb.dev.labjack.LabJack.CalMicroAmpere attribute*), 78

UNKNOWN (*hvl\_ccb.dev.heinzinger.HeinzingerConfig.RecordingsEnum attribute*), 72

UNKNOWN (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitCurrent attribute*), 75

UNKNOWN (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitVoltage attribute*), 74

UNKNOWN (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.Unreadiness attribute*), 74

UNREADINESS (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.Unreadiness attribute*), 61

update\_laser\_status () (*hvl\_ccb.dev.crylas.CryLasLaser method*), 64

update\_period (*hvl\_ccb.comm.opc.OpcUaCommunicationConfig attribute*), 3

update\_repetition\_rate ()

update\_shutter\_status () (*hvl\_ccb.dev.crylas.CryLasLaser method*), 64

update\_target\_pulse\_energy () (*hvl\_ccb.dev.crylas.CryLasLaser method*), 64

updateEspStageInfo (*hvl\_ccb.dev.newport.NewportSMC100PPConfig.EspStageConfig attribute*), 102

*attribute), 93*  
 USB (*hvl\_ccb.comm.labjack\_ljm.LJMCommunicationConfig.ConnectivityType*), 52  
*attribute), 9*  
 user\_position\_offset  
*(hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
*attribute), 95*  
 user\_steps () (*hvl\_ccb.dev.se\_ils2t.ILS2T method*),  
*114*

**V**  
 V (*hvl\_ccb.dev.heinzinger.HeinzingerPNC.UnitVoltage*  
*attribute), 75*  
 VA (*hvl\_ccb.dev.newport.NewportConfigCommands* at-  
*tribute), 86*  
 value (*hvl\_ccb.dev.labjack.LabJack.AInRange* at-  
*tribute), 78*  
 ValueEnum (*class in hvl\_ccb.utils.enum*), 123  
 VB (*hvl\_ccb.dev.newport.NewportConfigCommands* at-  
*tribute), 86*  
 velocity (*hvl\_ccb.dev.newport.NewportSMC100PPConfig*  
*attribute), 95*  
 visa\_backend (*hvl\_ccb.comm.visa.VisaCommunicationConfig*  
*attribute), 19*  
 VisaCommunication (*class in hvl\_ccb.comm.visa*),  
*17*  
 VisaCommunicationConfig (*class in*  
*hvl\_ccb.comm.visa*), 18  
 VisaCommunicationConfig.InterfaceType  
*(class in hvl\_ccb.comm.visa), 18*  
 VisaCommunicationError, 19  
 VisaDevice (*class in hvl\_ccb.dev.visa*), 121  
 VisaDeviceConfig (*class in hvl\_ccb.dev.visa*), 122  
 VisaStatusPoller (*class in hvl\_ccb.dev.visa*), 122  
 Volt (*hvl\_ccb.dev.pfeiffer\_tpg.PfeifferTPG.PressureUnits*  
*attribute), 100*  
 VOLT (*hvl\_ccb.dev.se\_ils2t.ILS2TRegAddr attribute*), 117  
 voltage\_lower\_limit  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
*attribute), 70*  
 voltage\_max (*hvl\_ccb.dev.supercube.constants.Power*  
*attribute), 37*  
 voltage\_max (*hvl\_ccb.dev.supercube2015.constants.Power*  
*attribute), 52*  
 voltage\_primary (*hvl\_ccb.dev.supercube.constants.Power*  
*attribute), 37*  
 voltage\_primary (*hvl\_ccb.dev.supercube2015.constants.Power*  
*attribute), 52*  
 voltage\_slope (*hvl\_ccb.dev.supercube.constants.Power*  
*attribute), 37*  
 voltage\_slope (*hvl\_ccb.dev.supercube2015.constants.Power*  
*attribute), 52*  
 voltage\_target (*hvl\_ccb.dev.supercube.constants.Power*  
*attribute), 37*

*voltage\_target (hvl\_ccb.dev.supercube2015.constants.Power*  
*attribute), 52*  
*voltage\_upper\_limit*  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
*attribute), 70*

**W**  
 WAIT\_AFTER\_WRITE (*hvl\_ccb.comm.visa.VisaCommunication*  
*attribute), 17*  
 wait\_operation\_complete ()  
*(hvl\_ccb.dev.visa.VisaDevice method), 121*  
 wait\_sec\_initialisation  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
*attribute), 70*  
 wait\_sec\_max\_disable  
*(hvl\_ccb.dev.se\_ils2t.ILS2TConfig attribute),*  
*115*  
 wait\_sec\_post\_absolute\_position  
*(hvl\_ccb.dev.se\_ils2t.ILS2TConfig attribute),*  
*115*  
 wait\_sec\_post\_activate  
*(hvl\_ccb.dev.sst\_luminox.LuminoxConfig*  
*attribute), 119*  
 wait\_sec\_post\_CANNOT\_DISABLE  
*(hvl\_ccb.dev.se\_ils2t.ILS2TConfig attribute),*  
*115*  
 wait\_sec\_post\_enable  
*(hvl\_ccb.dev.se\_ils2t.ILS2TConfig attribute),*  
*115*  
 wait\_sec\_post\_relative\_step  
*(hvl\_ccb.dev.se\_ils2t.ILS2TConfig attribute),*  
*115*  
 wait\_sec\_read\_text\_nonempty  
*(hvl\_ccb.dev.heinzinger.HeinzingerSerialCommunicationConfig*  
*attribute), 77*  
 wait\_sec\_settings\_effect  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
*attribute), 70*  
 wait\_sec\_stop\_commands  
*(hvl\_ccb.dev.heinzinger.HeinzingerConfig*  
*attribute), 73*  
 wait\_sec\_system\_lock  
*(hvl\_ccb.dev.ea\_psi9000.PSI9000Config*  
*attribute), 70*  
 wait\_timeout\_retry\_sec  
*(hvl\_ccb.comm.opc.OpcUaCommunicationConfig*  
*attribute), 13*  
 wait\_until\_motor\_initialized()  
*(hvl\_ccb.dev.newport.NewportSMC100PP*  
*method), 92*  
 wait\_until\_move\_finished()  
*(hvl\_ccb.dev.newport.NewportSMC100PP*  
*method), 92*  
 wait\_until\_ready()

```
(hvl_ccb.dev.crylas.CryLasLaser      method),  
64  
warning (hvl_ccb.dev.supercube.constants.Errors  attribute), 34  
WIFI (hvl_ccb.comm.labjack_ljm.LJMCommunicationConfig.ConnectionType  
      attribute), 9  
write () (hvl_ccb.comm.opc.OpcUaCommunication  
      method), 12  
write () (hvl_ccb.comm.visa.VisaCommunication  
      method), 18  
write () (hvl_ccb.dev.mbw973.MBW973 method), 83  
write () (hvl_ccb.dev.supercube.base.SupercubeBase  
      method), 24  
write () (hvl_ccb.dev.supercube2015.base.Supercube2015Base  
      method), 45  
write_absolute_position()  
    (hvl_ccb.dev.se_ilst2t.ILS2T method), 114  
write_bytes () (hvl_ccb.comm.serial.SerialCommunication  
      method), 15  
write_name () (hvl_ccb.comm.labjack_ljm.LJMCommunication  
      method), 8  
write_names () (hvl_ccb.comm.labjack_ljm.LJMCommunication  
      method), 8  
write_registers()  
  (hvl_ccb.comm.modbus_tcp.ModbusTcpCommunication  
    method), 10  
write_relative_step()  
  (hvl_ccb.dev.se_ilst2t.ILS2T method), 114  
write_termination  
  (hvl_ccb.comm.visa.VisaCommunicationConfig  
    attribute), 19  
write_text () (hvl_ccb.comm.serial.SerialCommunication  
      method), 15  
WRONG_ESP_STAGE (hvl_ccb.dev.newport.NewportSMC100PP.MotorErrors  
      attribute), 86
```

## Z

```
ZX (hvl_ccb.dev.newport.NewportConfigCommands  attribute), 86
```